

HALYARD* PUREZERO* NITRILE CLEANROOM GLOVES

HALYARD* offers an accelerator-free¹ portfolio of three non-sterile and two sterile cleanroom gloves to address the needs of pharmaceutical, medical device, microelectronics and semi-conductor manufacturing industries.

GLOVE SELECTION GUIDE

	Description	Designed for	Max Particle Count	ISO Class	Finish	Fingertips	Double Donning	Size Range	Fingertip Thickness	Case Count	
NON-STERILE	HALYARD* PUREZERO* HG3 White Nitrile Gloves, 12", 305 mm	Semiconductor, Pharmaceutical, Medical Device Manufacturing	<950	ISO Class 3 or higher and Grade B/C/D cleanrooms	Tacky Grip	Textured Fingertips	Recommended for Outer Glove	XS-XL	0.16mm (6.3 mil)	1000/cs Ambidextrous	
	HALYARD* PUREZERO* HG3 Light Blue Nitrile Gloves, 12", 305 mm	Semiconductor, Pharmaceutical, Medical Device Manufacturing	<1200	ISO Class 3 or higher and Grade B/C/D cleanrooms	Tacky Grip	Textured Fingertips	Recommended for Outer Glove	XS-XL	0.10mm (3.9 mil)	1500/cs Ambidextrous	
	HALYARD* PUREZERO* HG3 White SGX* Nitrile Gloves, 12", 305 mm	Semiconductor, Pharmaceutical, Medical Device Manufacturing	<950	ISO Class 3 or higher and Grade B/C/D cleanrooms	Smooth Grip	Textured Fingertips	Outer or Under	XS-XL	0.16mm (6.3 mil)	1000/cs Ambidextrous	
STERILE	HALYARD* PUREZERO* HG3 Light Blue Sterile Nitrile Gloves, 12", 305 mm	Pharmaceutical, Biotechnology, Sterile Compounding, Aseptic Processing	<1200	ISO Class 3 or higher and Grade A/B cleanrooms	Smooth Grip	Textured Fingertips and Palms	Outer or Under	6.0 6.5 7.0 7.5 8.0 8.5 9.0 10.0	0.10mm (3.9 mil)	300 pairs/cs Hand Specific	
	HALYARD* PUREZERO* HG3 White Sterile Nitrile Gloves, 12", 305 mm	Pharmaceutical, Biotechnology, Sterile Compounding, Aseptic Processing	<950	ISO Class 3 or higher and Grade A/B cleanrooms	Smooth Grip	Textured Fingertips and Palms	Outer or Under	6.0 6.5 7.0 7.5 8.0 8.5 9.0 10.0	0.16mm (6.3 mil)	200 pairs/cs Hand Specific	
Applies to all PUREZERO* Gloves:		Accelerator-Free ¹ Static Dissipative in Use ² Compliant with the REACH regulation		AQL 1.0	CE 2797 (PPE Cat. III)	ISO 374-5:2016 VIRUS					
 White Non-Sterile and White Sterile		 Light Blue Non-Sterile		 White SGX* and Light Blue Sterile							

For more information or samples, contact your distributor or visit: www.purezerogloves.com

¹ Not formulated with these commonly used vulcanizing chemicals: Sulfur, Thiurams, Thioxoles, Guanidines, Carbamates, Amines, Thiophosphates, and Sulfenamides.

² Tested per ASTM D6319, EN 455-2

*Registered Trademark or Trademark of O&M Halyard or its affiliates. © 2024. All rights reserved. 008735 COPY-04864 Jan 2024