




HALYARD* PUREZERO*
HG3 Cream Latex Gloves
 For the Cleanroom Environment
 For Industrial Use Only

TECHNICAL DATA SHEET

Description

HALYARD* **PUREZERO*** HG3 Cream Latex Gloves are designed for pharmaceutical and medical device manufacturing applications. They have undetectable residual accelerator and other contact sensitizers. These gloves are clean processed (washed repeatedly in deionized water) to ensure consistent control of low particles and extractables and are recommended for use in ISO Class 3 or higher and Grade B/C/D cleanrooms.



Manufactured at OUR OWN
 Safeskin Facility in Thailand.

Meets
 USP <800>
 Guidelines.

Tested
 Against
 Fentanyl

Cleanliness Properties

Max Particle Count (>0.5 µm)	<1200 Particles/cm ²	IEST RP-CC005
Ionic Content (Extractable ions)	Max Level (ug/g)	IEST RP-CC005
Calcium	15	
Chloride	100	
Magnesium	5	
Nitrate	15	
Potassium	5	
Sodium	7	
Sulfate	15	
Zinc	100	
Ammonium	7	

Physical Properties

AQL	1.0
Non-Sterile	✓
Ambidextrous	✓
Textured Fingertips & Palms	✓
Powder-Free	✓
Silicone-Free	✓
Tensile Strength ¹	32 MPa (Target)
Ultimate Elongation ¹	875% (Target)
Stress at 500% Elongation	<= 5.5 MPa
Protein Content	<= 50 ug/g
Shelf Life	3 Years

Glove Dimensions

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
Glove Length (inch/cm)	12"/30.5	12"/30.5	12"/30.5	12"/30.5	12"/30.5
Middle Finger Length (mm)	71	78.5	80.5	88.1	91.6
Width of Palm (mm)	70	85	95	110	120
Finger Tip Thickness	0.20 mm (7.9 mil)	0.20 mm (7.9 mil)	0.20 mm (7.9 mil)	0.20 mm (7.9 mil)	0.20 mm (7.9 mil)
Palm Thickness	0.18 mm (7.1 mil)	0.18 mm (7.1 mil)	0.18 mm (7.1 mil)	0.18 mm (7.1 mil)	0.18 mm (7.1 mil)
Cuff Thickness	0.15 mm (5.9 mil)	0.15 mm (5.9 mil)	0.15 mm (5.9 mil)	0.15 mm (5.9 mil)	0.15 mm (5.9 mil)



HALYARD* PUREZERO*

HG3 Cream Latex Gloves

For the Cleanroom Environment
For Industrial Use Only

TECHNICAL DATA SHEET

Packaging Data

Double bag, plus case liner
1000 gloves per case: 100 gloves per PE bag X 10 PE bags per lined carton
Packaged in ISO Class 5 Cleanroom

Quality & Regulatory Standards

Compliant to these regulatory standards:

- ISO 9001
- ISO 13485
- ISO 14001
- ISO 10993-10
- ISO 10993-11

Additional Glove Information

Recommended for use in ISO Class 3 or higher and Grade B/C/D cleanrooms.

Tested against 35 chemicals, 14 chemotherapy drugs and Fentanyl.

Made in Thailand

Declaration of Conformity (DoC) and Certificates of Analysis (COA) for every production lot available online at

halyardhealth.com/information

Storage Instructions

HALYARD* PUREZERO* HG3 Cream Latex Gloves shall be stored in conditions where the product is kept dry (away from moisture), away from direct sunlight, away from sources of heat, and away from radioactive sources.

Ordering Information

**HALYARD* PUREZERO* HG3 CREAM LATEX GLOVES,
NON-STERILE, AMBIDEXTROUS, TEXTURED
FINGERTIPS AND PALMS**

Size	Code
XS	CLN8231XS
SM	CLN8231SM
MD	CLN8231MD
LG	CLN8231LG
XL	CLN8231XL

For additional information or samples,
contact your local distributor

¹ Tested per ASTM D412

This fact sheet has been created using the most recent information. In the interest of continuous improvement, the characteristics of the products may change without prior notice.