

MEDI CHOICE®

Product Catalog

Connecting the world of medical products to the point of care™.

Owens & Minor is pleased to present the 2015 MediChoice product catalog, an indispensable resource for medical professionals seeking top-quality products at affordable prices. Since rolling out MediChoice in 2002, we have enhanced and expanded the portfolio each year, and this issue is our best yet. We've worked hard to make sure this catalog features the most products and the most diverse set of suppliers to date, all while keeping an eye on the evolving demands of today's healthcare professionals.

Innovation and Collaboration

Since putting out the last issue, we have been busy gathering feedback from administrators, technicians, nurses, and doctors in hospitals and clinics across the United States to see what more we could do to help them do their jobs effectively and efficiently. The voices have been diverse, but they all basically say the same thing: We need solutions that balance the needs of patient care *and* cost-saving budget initiatives.

We shared the input we received from you with our product development team, which in turn shared it with our manufacturers, all of whom were chosen because they meet our stringent requirements for quality, service, and value. We work to develop close relationships with those manufacturers, and that helps ensure we can keep delivering clinically accepted products that compare with leading brands.

From admissions to discharge and everything in between, know that each product was carefully selected by our trained product development team to ensure it's something you can purchase with confidence.

Want to know more about how MediChoice can work for you? Turn to the inside back cover of this catalog for a complete list of Owens & Minor locations, each with a team of representatives eager to put the power of MediChoice to work for you.

*Solutions that balance the needs of patient care
and cost-saving budget initiatives*

MEDI CHOICE®

*Your choice for quality,
value and unmatched service*

CONTENTS



Admissions

ID Cards	6
Specimen Pans	6
Triangular Graduates	6
Patient Supplies	6, 7
Medicine Cups	7
Personal Care	7, 9
Pillows	9, 10
Pillow Cases	10
Bedside Bags	12
Patient Belongings Bags	12



Anesthesia

Anesthesia Masks	14
Berman Airways	14
Color-Coded Guedel Airways	15
Nasopharyngeal Airways	15
Pressure Infuser Bags	15



Apparel—Patient

Bibs	16
Exam Capes	16
Exam Gowns	16
Disposable Exam Shorts	16
Slippers	17, 18
Telemetry Pouches	18



Apparel—Staff

Bouffant Caps	20
Disposable Aprons	20
Protective Eyewear	20
Surgeon Caps	20
Disposable Scrubs	21
Disposable Coveralls	21
Disposable Lab Coats	21
Disposable Staff Jackets	21
Face Masks	22, 23
Fluid-Resistant Gowns	23
Isolation Gowns	23

Shoe Covers	24
Personal Protective Gowns	24, 25



Central Supply

Cart Covers	26
Sterilization Pouches	26
Utility Scissors	26
Instruments	27
Glove Box Holders	31
Infection Control Cabinets	31
Enzyme Spray Pre-Cleaner	34
Enzyme Cleaner Pre-Soak	34
Instrument Lubricants	34
Specialty Products	34
Cart Cleaners	35
Neutral pH Instrument Cleaners	35
Rinse Aids	35



Equipment

Scales & Measuring Devices	37-44
IV Stands	45
Hampers	44, 45
IV Stand Accessories	45, 46
Basin Stands	47
Mayo Stands	47
Step Stools	47, 48
Bassinets	48
Blanket Warmers	49
Kick Buckets	50
Waste Receptacles	50
Instrument Tables	50
Lighting—Exam & Procedure	50
Lighting—Minor	50
Surgery & Diagnostic	52, 53
Exam Stools	53-56
Manual Exam Tables	56-64
Powered Exam Tables	57
Powered Podiatry Chairs	58
Medical Carts	59
Overbed Tables	59
Medical Storage Carts	64, 65
O.R. & Stretcher Pads	66
Mattress Pads	67, 68
Mobile Work Stations	68
Treatment Tables	71-77
Sleeper Chairs	75, 76
Side Chairs	76
Bariatrics	77
Blood Drawing Chairs	79
Powered Procedure Chairs	80, 81
Utility Carts	81, 82
Utility Tables	82
Computer Accessories	83
Recliners	92



Exam Room Products

Digital Thermometers	84
Cotton-Tipped Applicators	89
OB/GYN & Procto Applicators	89
Ultrasound Gel	89
Electrodes	90
Lubricating Jelly	90
Vaginal Specula	90
Exam Table Paper	91
Patient Drapes	91
Professional Towels	91
Blood Pressure Cuffs	84-87
Dopplers and Probes	87, 88
Hand Held Sphygmomanometers	89-95
Chart Paper	92
Stethoscopes	92-95



Gloves

Latex Exam Gloves	96
Nitrile Exam Gloves	96
Stretch Vinyl Exam Gloves	97
Vinyl Exam Gloves	97
High-Risk Gloves	98
Sterile Procedure Gloves	98



Housekeeping

Recloseable Bags	99
Biohazard Waste Bags	99
Can Liners	100-103
Linen Bags	104



Incontinence

Foaming Skin Cleansers	105
Odor Eliminators	105
Seamless Knit Pants	105
Underpads	105
Adult Briefs	106
Protective Underwear	106
Washcloths	106

NOT MADE WITH
natural rubber latex

Indicates Latex-Free

diverse
supplier

Indicates Diverse Supplier

sterile EO

Indicates Ethylene Oxide Sterilized

sterile R

Indicates Gamma Radiation Sterilized

sterile A

Indicates Aseptic
Processing Sterilized

made in the
USA

Indicates made in the
United State of America

Adult Wipes
Cleansing Creams
Skin Protectant
Skin Recovery
Skin Treatment



Kits & Trays

Soft Tissue Biopsy Trays
Lumbar Puncture Trays
Bone Marrow Trays
Incision & Drainage Trays
Skin Staple Remover Trays
Central Line Dressing Trays
IV Start Kits
Tracheostomy Trays
Laceration Trays
Suture Removal Trays
Suture Trays
Irrigation Trays



Laboratory

Lab Markers
Task Wipes
Tourniquets
Safety Lancets
Specimen Containers
Specimen Bags
Pregnancy Tests
Strep A Tests
Mono Tests
Fecal Occult Tests
H. pylori Tests
Laboratory Accessories
Differential Counters
Slides and Cover Glass
Microscopes
Incubators
Centrifuges
Mixers



Maternity & Baby Care

Footprinters
Umbilical Cord Clamps
Umbilical Cord Clamp Cutters
Fetal Monitoring Belts
Abdominal Fetal Monitoring Bands
Ear and Ulcer Syringes
Infant Heel Warmers
Baby Scale Liners
Baby Bonnets

107 Baby Booties
107 Baby Lotion
107 Talc-Free Baby Powder
107 Baby Shampoo
107 Baby Wipes
Mesh Briefs
Obstetric Pads
Perineal Cold Packs
Perineal Wash Spray



Mobility Aids

112 Bathroom Safety Aids
109 Commode Chairs
112 Crutches
109 Canes
113 Walkers
111 Wheelchairs
111 Wheelchair Accessories



O.R. Products

118 Arm Board Overlays
118 Abduction Pillows
118 Foam Head Positioners
122 Cradle Positioners
122 Foam Heel & Elbow Protectors
125 Foam Positioning Wedges
125 Anti-Fog Solution
124 Biopsy Punches
124 Cordless Cauteries
124 Decanting Devices
123 Disposable O.R. Towels
124 Disposable Safety Scalpels
118 Electrosurgical Pencils/Tips
119 Electrosurgical Tip Cleaners
121 Esmark Bandages
20, 121 Eye Protectors
119 Hypo-Hyperthermia Blankets
119 Insufflation Tubing
120 Kick Bucket Liners
Light Handle Covers
Scope Warmers
Skin Markers
Skin Staplers
Solidifiers
Sponge Counters
Bias-Cut Stockinettes
Impervious Stockinettes
Sterile 1-Ply Stockinettes
Sterile 2-Ply Stockinettes
Tubular Stockinettes
Suction Tubing
Surgical Bras
Tubing Connectors
Utility Markers
Vessel Loops
Yankauers




Ortho/Rehab

Abdominal Binders
Ankle Braces
Ankle Walkers
Cast Boots
Cervical Collars
Gait Belts
Knee Immobilizers
Night Splints
Post-Op Shoes
Arm Slings
Shoulder Immobilizers
Wrist Splints
Heel and Elbow Protectors



Patient Care

Alcohol and Peroxide
Alcohol Prep Pads
Towelettes
Nail Polish Remover Pads
Instant Hand Sanitizers
Instant Cold Packs
Instant Warm Packs
Reusable Hot or Cold Gel Packs
Refillable Ice Bags
Heat Therapy Pads
Heat Therapy Pumps
Lemon Glycerin Swabs
Oral Swabs
Multipurpose Tube Holders
Tongue Depressors
Penlights
Pill Crushers & Splitters
Enema Bag Kits
Post-Mortem Bags & Kits



Patient Monitoring

Cables and Leads 162-164



Respiratory

Aerosol Masks
CO₂ Cannulas
Face Tents
Nebulizers
Oxygen Cannulas
Oxygen Masks
Oxygen Tubing
Trach Masks
Respiratory Set-Up Bags
Sodium Chloride Solution
Specimen Traps
Suction Catheters
Suction Catheter Kits
Suction Catheter Trays
Trach Tube Holders



Skin Care

Mattress Overlays
Pressure-Relieving Heel Protectors
Wheelchair Cushions
Heel Pillows



Urology

Irrigation Syringes
Calculi Strainers
Catheter Leg Straps
Pediatric Urine Collectors
Leg Bags
Foley Insertion Trays
Urethral Catheterization Trays



Wound Care

Povidone Iodine Solutions
Cotton Balls
Adhesive Bandages
Wound Closure Strips
Transparent Film Dressings
Cohesive Bandages
Elastic Bandages
Tubular Net Bandages
Arm Sleeves
Tubular Elastic Support Bandages
Paste Bandages
Saline Wound Flush
Compression Bandage Systems



ID Cards



- Sturdy Polyvinyl Chloride material
- Available in numerous sizes and thicknesses to fit all standard embossers
- Matte or polished finish options
- Variety of colors to meet specific needs
- Special options available including hole-punches, notches and scoring

Item No.	Description	Quantity
----------	-------------	----------

Card Identification Patient Sturdy Polyvinyl Chloride Matte Polished

MC2451YL12	Yellow Super Stressed CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
M3081BL18	Blue CR80 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC2451BL	Blue CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC2451CR	Cranberry CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC2451WH	White CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC2451YL	Yellow CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC2481BL	Blue CR80 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC2481CR	Cranberry CR80 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC2481WH	White CR80 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3051BL	Blue CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3051CR	Cranberry CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3051DB	Dark Blue CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3051GR	Green CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3051OR	Orange CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3051RD	Red CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3051YL	Yellow CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3081BL	Blue CR80 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3081CR	Cranberry CR80 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3081GR	Green CR80 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA

Card Identification Patient Sturdy Polyvinyl Chloride Matte Polished with Hole

M2451BL1	Blue CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
M2451BL10	Blue CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
M2481BL1	Blue CR80 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
M3051BL1	Blue CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
M3081BL1	Blue CR80 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MCCR530B	Blue CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA

Card Identification Patient Sturdy Polyvinyl Chloride Both Sides Matte

MC2453WH	White CR50 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3053BL	Blue CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MC3083BL	Blue CR80 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
MCCR524M	Blue CR50 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA

Card Identification Patient Sturdy Polyvinyl Chloride Matte

MC2483BL	Blue CR80 × 0.024 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
----------	------------------------------	---------------------------

Card Identification Patient Sturdy Polyvinyl Chloride Both Sides Polished

MC3082WH	White CR80 × 0.030 Millimeter	1/EA, 500 EA/BX, 10 BX/CA
----------	-------------------------------	---------------------------

PATIENT SUPPLIES

Flexible Drinking Straws

NOT MADE WITH natural rubber latex

- Disposable, white plastic straw with flexible bend to conveniently assist patients while reclined
- Individually paper wrapped
- 500 straws per box
- 7.75 Inch length × 0.25 diameter
- Withstands temperatures up to 392 degrees Fahrenheit



Item No.	Description	Quantity
DS0500	Straw Drinking Flexible Plastic 7.75 × 0.25 Inch White Individually Wrapped Disposable	500/BX, 20 BX/CA

Souffle Cups

NOT MADE WITH natural rubber latex

diverse supplier

- Titanium White 0.75 ounce plastic cup
- Sturdier and more stable than paper portion cups
- Won't absorb liquids and weaken like paper cups
- Better for pill crushing; powder doesn't stick to inside of the cup
- Flexible, smooth edges are easy on sensitive lips
- Green product—100 percent recyclable
- Microwavable (unlike paper cups)
- Single patient use



Item No.	Description	Quantity
1060300	Cup Souffle Plastic White 0.75 Ounce	1/EA, 100 EA/SL, 50 SL/CA

Please see our full line of Pill Crushers on page 161 in the Patient Care Section.



PATIENT SUPPLIES

Disposable Utensils

NOT MADE WITH
natural rubber latex

- Beige Spoon constructed of durable polystyrene construction
- Dispenser Box of 100 spoons per box; 10 boxes per case
- White Knife and Fork constructed of durable polypropylene
- Knives and Forks are both conveniently packaged in resealable bag of 100; 10 bags per case
- All utensils are disposable and unwrapped



Item No.	Description	Quantity
Patient Plastic		
UT1000	Spoon Heavy Duty Beige 5.94 Inch Unwrapped Disposable	100/BX, 10 BX/CA
UT1020	Fork Medium Weight White 3.2 Gram Polypropylene Packaged In Zip-Lock Bag	100/BG, 10 BG/CA
UT1010	Knife Medium Weight White 3.0 Gram Polystyrene Packaged In Zip-Lock Bag	100/BG, 10 BG/CA

Triangular Graduates

NOT MADE WITH
natural rubber latex

- 32 fluid ounce, easy-to-hold, three-sided container
- Available in clear polystyrene or transparent polypropylene with or without black graduations
- Graduated in ounces and cubic centimeters (32 ounce/1000 cubic centimeter)
- Frosted area for recording patient information



Item No.	Description	Quantity
Container Laboratory Graduated Triangular Three-Sided 1000 Cubic Centimeter Disposable		
33302	Translucent Polypropylene	1/EA, 10 EA/SL, 20 SL/CA
33303	Clear Polystyrene	1/EA, 10 EA/SL, 20 SL/CA
33302BL	Translucent Polypropylene with Black Marking	1/EA, 10 EA/SL, 20 SL/CA

Medicine Cups

NOT MADE WITH
natural rubber latex

diverse
supplier

- Sturdy 1 ounce plastic clear cup
- Graduated measurements:
 - Drams: 1-8 in 1 dram increments
 - Ounces: 0.125, 0.25, 0.50, 0.75, 1
 - Tablespoons: 0.50 and 1
 - CCA: 2.5, 5, 7.5, 10, 15, 20, 25, 30
- Single patient use



Item No.	Description	Quantity
33304	Cup Medicine Graduated 1 Ounce Plastic Clear Single-Patient Use	1/EA, 100 EA/SL, 50 SL/CA

PERSONAL CARE

Specimen Pans

NOT MADE WITH
natural rubber latex

diverse
supplier

- Molded tray of sturdy, white plastic
- Easy-to-read measurements in ounces and cubic centimeters
- Capacity to measure 27 ounces and 800 cubic centimeters
- Easy handling with grooved pour spout for transferring specimen samples



Item No.	Description	Quantity
33305	Pan Specimen 27 Ounce Sturdy Plastic White	1/EA, 100 EA/CA

Comb, Hairbrush, Toothbrush

NOT MADE WITH
natural rubber latex

- Convenient inner packaging for efficient par stocking
- Black plastic comb; half fine-tooth and half coarse-tooth
- Hairbrush features seven rows of nylon bristles
- Toothbrush has an angle-head with clear polypropylene bristles
- Premium toothbrush has 40-Tuft angled head



Item No.	Description	Quantity
Toothbrush Angled Head Clear Polypropylene Bristles Individually Wrapped		
PTB100	Premium Assorted Colors 40-Tuft Bristles	1/EA, 144 EA/BX, 10 BX/CA
PC7000	Standard Ivory 30-Tuft Bristles	1/EA, 144 EA/BX, 10 BX/CA
Comb Hair Plastic Black Half Fine-Tooth Half Coarse-Tooth		
PC2000	7 Inch	1/EA, 12 EA/BG, 120 BG/CA
Brush Hair Nylon Plastic Ivory Handle		
PC3000	Black 86-Tuft Bristles	1/EA, 12 EA/BX, 24 BX/CA

Shampoo and Body Bath

made in the
USA

NOT MADE WITH
natural rubber latex

diverse
supplier

- Gentle non-irritating formula
- Can be used on hair and body
- Flip-top dispensing cap for easy application
- pH-balanced
- Pleasant, mild fragrance



Item No.	Description	Quantity
PC8020	Shampoo And Body Bath With Flip-Top Cap 2 Ounce	1/EA, 144 EA/CA

PERSONAL CARE

Shaving Cream

NOT MADE WITH
natural rubber latex

- Thick, rich foaming lather
- Gentle formula
- Pleasant fragrance
- Alcohol free



Item No.	Description	Quantity
PC2515	Cream Shaving Alcohol Free 1.5 Ounce	1/EA, 144 EA/CA

Hand and Body Lotion

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Available in cocoa butter fragrance or fragrance-free
- Contains aloe vera and Vitamin E to promote skin health
- Quick-absorbing, non-greasy formula
- Emollients and moisturizers to replenish natural skin oils
- Convenient flip-top cap for easy dispensing
- Extra mild, pH-balanced



Item No.	Description	Quantity
----------	-------------	----------

Lotion Hand and Body Fragrance-Free with Flip-Top Cap

PC102FF	Fragrance-Free 2 Ounce	1/EA, 144 EA/CA
PC104FF	Fragrance-Free 4 Ounce	1/EA, 96 EA/CA
PC4020	Scented 2 Ounce	1/EA, 144 EA/CA
PC4035	Scented 4 Ounce	1/EA, 96 EA/CA

Roll-on Deodorant

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Helps reduce wetness
- Controls odor
- Alcohol free
- Light, fresh scent



Item No.	Description	Quantity
----------	-------------	----------

PC5015	Deodorant Antiperspirant Roll-On Lightly Scented 1.5 Ounce	1/EA, 96 EA/CA
--------	--	----------------

No-Rinse Shampoo and Body Bath

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Gentle shampoo cleans without rinsing
- Can be used on hair or skin
- Flip-top cap for easy dispensing
- Gentle cleansing, less irritating than soap
- pH-balanced



Item No.	Description	Quantity
----------	-------------	----------

Shampoo Hair and Body No-Rinse with Flip-Top Cap

PC9020	2 Ounce	1/EA, 144 EA/CA
PC9080	8 Ounce	1/EA, 48 EA/CA

PERSONAL CARE

Daily Moisturizing Lotion

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Daily conditioning lotion
- Does not interfere with the skin antiseptic chlorhexidine gluconate (CHG)
- Fragrance- and dye-free
- Hypoallergenic



Item No.	Description	Quantity
----------	-------------	----------

PC6002	Lotion Moisturizing Chlorhexidine-Gluconate Fragrance-Free Sunflower Seed Oil 2 Ounce	1/EA, 96 EA/CA
--------	---	----------------

PC6004	Lotion Moisturizing Chlorhexidine-Gluconate Fragrance-Free Sunflower Seed Oil 4 Ounce	1/EA, 60 EA/CA
--------	---	----------------

Mouthwash

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Alcohol-Free, Gluten-Free
- Cinnamon flavor
- Refreshing, cleansing mouth rinse
- Contains Methyl Paraben as a preservative
- Safety sealed for your protection



Item No.	Description	Quantity
----------	-------------	----------

PC7500	Mouthwash Oral Care Alcohol Free Cinnamon Burst Flavor 4 Ounce	1/EA, 60 EA/CA
--------	--	----------------

Lip Balm

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Helps moisturize, nourish and hydrate
- Prevents dry, chapped lips
- SPF 15 Sunscreen



Item No.	Description	Quantity
----------	-------------	----------

Balm Lip SPF15

LBS100	Box Of 600	1/EA, 600 EA/CA
LBS100C	Box Of 100	1/EA, 100 EA/BX, 4 BX/CA

Shampoo and Body Wash

made in the USA **NOT MADE WITH** natural rubber latex

- Two great scents: Apricot and Cucumber Melon
- Contains gentle conditioners to keep hair soft and manageable
- pH balanced—less irritating to sensitive skin
- Two sizes: 4 Ounce and 8.5 Ounce
- Contains aloe vera and proteins to promote skin health



Item No.	Description	Quantity
----------	-------------	----------

Shampoo Hair and Body

PC28202	Apricot Scent 4 Ounce	1/EA, 48 EA/CA
PC28203	Apricot Scent 8.5 Ounce	1/EA, 48 EA/CA
PC29022	Cucumber Melon Scent 4 Ounce	1/EA, 48 EA/CA
PC29023	Cucumber Melon Scent 8.5 Ounce	1/EA, 48 EA/CA

PERSONAL CARE

Denture Products

NOT MADE WITH
natural rubber latex

- Denture Cream Adhesive (MCDA20 and MCDA24):
 - Aids in a great fit and a strong, long-lasting hold of dentures
 - Extra strength holding power
 - Enables patient to confidently eat, drink, and speak while keeping dentures in place
 - Comparable to store brands
 - Easy to open cap
- Denture Cleaning Tablets (MCDENT90):
 - Powerful effervescent cleansing action to keep dentures fresh while soaking and cleaning
 - Antibacterial
 - Fight odor-causing bacteria and plaque build-up while dissolving away stains
 - Individually foil wrapped



Item No.	Description	Quantity
Cream Denture Adhesive		
MCDA20	2.0 Ounce	1/EA, 6 EA/BX, 12 BX/CA
MCDA24	2.4 Ounce	1/EA, 12 EA/BX, 6 BX/CA
Denture Cleanser		
MCDENT90	Tablets	90 EA/BX, 24 BX/CA

Disposable Razors

NOT MADE WITH
natural rubber latex

- High-technology twin blades with microscopic serrated edges to produce a closer, safer shave
- Grooved handle for improved grip to help keep from slipping in hand
- Strong, fixed head reduces breakage
- Packaged as 10 each per bag for easy restocking and inventory



Item No.	Description	Quantity
RZR2000	Razor Shaving Twin Blade Disposable	1/EA, 10 EA/BG, 72 BG/CA

Toothpaste

NOT MADE WITH
natural rubber latex

- Contains Fluoride for protection against cavities
- Gentle formula for sensitive teeth and gums
- Fresh mint flavor and smooth texture
- Two personal care sizes; 0.6 ounce and .85 ounce



Item No.	Description	Quantity
Toothpaste Fluoride Mint Flavor Not Made with Natural Rubber Latex		
PC7306	0.6 Ounce	1/EA, 144 EA/BX, 5 BX/CA
PC7385	0.85 Ounce	1/EA, 144 EA/BX, 5 BX/CA

PERSONAL CARE

Nail Clippers

NOT MADE WITH
natural rubber latex

diverse
supplier

- Constructed of durable stainless steel for long-lasting performance
- Includes a nail file with a cleaning tip for shaping and smoothing
- Cutting edge provides a clean, even cut
- Toenail clippers include a thumb pad lever to help prevent slipping



Item No.	Description	Quantity
Clipper with File		
PC2600	Fingernail	1/EA, 6 EA/BX, 48 BX/CA
PC2610	Toenail	1/EA, 6 EA/BX, 48 BX/CA

Disposable Pillows

made in the
USA

NOT MADE WITH
natural rubber latex

- Soft and comfortable spunbond polypropylene cover
- 100 percent garneted polyester fiber fill
- Non-allergenic, non-irritating, non-linting
- Single patient use
- Compressed packaging saves storage space



Item No.	Description	Quantity
Pillow Disposable Spunbond Polypropylene Cover Polyester Fiber Fill White Single-Patient Use		
089DP10	8 Ounce 13 × 17 Inch	1/EA, 24 EA/CA
089DP20	10 Ounce 16 × 22 Inch	1/EA, 12 EA/CA
089DP30	12 Ounce 18 × 24 Inch	1/EA, 12 EA/CA
089DP40	15 Ounce 18 × 24 Inch	1/EA, 12 EA/CA
089DP50	17 Ounce 20 × 26 Inch	1/EA, 12 EA/CA
089DP55	20 Ounce 18 × 24 Inch	1/EA, 12 EA/CA
089DP60	20 Ounce 20 × 26 Inch	1/EA, 12 EA/CA
089DP65	24 Ounce 18 × 24 Inch	1/EA, 12 EA/CA
089EC100	15 Ounce 18 × 24 Inch	1/EA, 12 EA/CA
089EC10A	Vacuum-Packed 18 × 24 Inch	1/EA, 12 EA/CA

PILLOWS

Limited Reusable Pillows



- Cost effective limited reuse pillow
- Moisture and stain resistant
- Multiple patient uses



Item No.	Description	Quantity
Pillow Limited Reusable Moisture- and Stain-Resistant		
089LP20	18 Ounce Blue 18 × 24 Inch	1/EA, 12 EA/CA
089LP30	18 Ounce Blue 24 × 26 Inch	1/EA, 12 EA/CA

Reusable Pillows



- Soft and comfortable vinyl ticking wipes clean
- Model 89RP5 has easy-clean nylon ticking
- 1 percent garneted polyester fiber fill
- Fluid proof, flame retardant material
- Anti-microbial, anti-fungal
- Non-allergenic, non-irritating, non-linting
- Compressed packaging to save space



Item No.	Description	Quantity
Pillow Reusable Vinyl Ticking Cover Polyester Fiber Fill		
089RP10	7 Ounce Blue 13 × 17 Inch	1/EA, 24 EA/CA
089RP20	17 Ounce Blue 18 × 24 Inch	1/EA, 12 EA/CA
089RP30	19 Ounce Blue 20 × 26 Inch	1/EA, 12 EA/CA
089RP40	24 Ounce Blue 20 × 26 Inch	1/EA, 12 EA/CA
089RP40W	24 Ounce White 20 × 26 Inch	1/EA, 12 EA/CA
089RP50	22 Ounce Yellow 21 × 27 Inch	1/EA, 12 EA/CA
089RP60	18 Ounce Yellow 20 × 26 Inch	1/EA, 12 EA/CA
089RP65	16 Ounce Yellow 18 × 24 Inch	1/EA, 12 EA/CA

Pillow Cases



- Premium:
 - Made of Fabricel®: Superior comfort and strength, soft, crepe-like draping
 - Leak proof moisture barrier protects pillow
- Standard:
 - Durable embossed tissue over polyethylene construction
 - Leak proof moisture barrier protects pillow
- Both products are disposable and intended for single patient use



Item No.	Description	Quantity
Pillowcase Single-Patient Use Disposable		
MCDPC001	Standard Tissue Over Polyethylene White 21 × 30 Inch	1/EA, 100 EA/CA
MCDPC002	Premium Fabricel White 21 × 30 Inch	1/EA, 100 EA/CA

PATIENT BEDSIDE

Patient Mugs



- Double-walled for insulation and decreased condensation
- For use with hot or cold liquids
- Clear markings in both ounce and millimeter, 32 ounce capacity
- Can be marked with patient information
- Translucent to easily evaluate fluid level
- Sturdy handle for stable lifting, even when full
- Custom options available



Item No.	Description	Quantity
Mug Patient Insulated with Straw Clear		
PPB032C	With Graduation 32 Ounce	1/EA, 50 EA/CA



Bedpans



- Stackable Commode Bedpan
 - Stackable design to maximize storage space
 - Fits most standard toilets
- Pontoon Bedpan
 - Rolled reinforced edges for stability and easy placement
 - Saddle-style design
- Fracture Bedpan
 - Ideal for immobile patients
 - Used also as a female urinal and pediatric bedpan
 - Features tapered front and loop handle for easy positioning



Item No.	Description	Quantity
Bedpan		
PPB1010	Fracture Tapered Front And Loop Handle Mauve Size 12	1/EA, 50 EA/CA
PPB1110	Stackable Mauve Size 12	1/EA, 50 EA/CA
PPB1210	Pontoon Rolled Reinforced Edges Saddle Style Design Mauve Size 12	1/EA, 20 EA/CA

PATIENT BEDSIDE

Carafes and Tumblers



- Carafe with Lid
 - Large 32 ounce, 1000 milliliter capacity
 - Includes clear lid with ice retention guard
 - Sturdy reinforced ergonomic handle for easy grip and use
 - Foam inner liner (sold separately) to reduce condensation
- Clear Plastic Carafe Liner (sold separately)
 - Ideal for refills and reducing condensation
- Tumblers
 - Available in clear and opaque mauve
 - Dual scale graduated markings
 - 7 ounce size ideal for patient grasping
 - Conveniently packaged in sleeves of 25 each



Item No.	Description	Quantity
Tumbler Patient Dual Scale Graduated		
PPB2010	7 Ounce Mauve	1/EA, 500 EA/CA
PPB2020	7 Ounce Clear	1/EA, 500 EA/CA
Carafe Patient with Lid Plastic		
PPB2110	32 Ounce Mauve	1/EA, 500 EA/CA
Liner Carafe		
PPB2120	Foam	1/EA, 500 EA/CA
PPB2122	Clear Plastic Graduated 28 Fluid Ounces Or 800 Cubic Centimeters	1/EA, 25 EA/SL, 20 SL/CA

Wash Basins and Soap Dishes



- Wash Basins
 - Large Mauve 7-Quart, 6.5 liter size
 - Dual graduated measurement markings; U.S. and Metric
 - Reinforced, rolled, textured rim to prevent gloved hands from slipping
 - Prefabricated slot for MediChoice® soap dish
- Soap Dish
 - Can be clipped to MediChoice® wash basin or used separately
 - Bottom ridges keep soap dry
 - Mauve color to match basin



Item No.	Description	Quantity
PPB3010	Dish Soap Snap-On For Wash Basin Mauve	1/EA, 500 EA/CA
PPB3110	Basin Wash Rectangular Dual Graduated 6.5 Liter 7 Quart Mauve	1/EA, 50 EA/CA

PATIENT BEDSIDE

Emesis Basins



- Basin has two separate graduated scales: ounce and milliliter
- Contoured basins are excellent for oral hygiene and ear irrigation
- Sturdy construction for stability
- Smooth surface for easy cleaning



Item No.	Description	Quantity
Basin Emesis Contoured		
PPB4010	Large 700 Milliliter Mauve 9.75 Inch	1/EA, 250 EA/CA
PPB4110	Small 500 Milliliter Mauve 9 Inch	1/EA, 250 EA/CA

Male Urinals



- Sturdy translucent plastic resists accidental collapsing of sides
- Wide base design provides stability
- Secure, snap-on lid over rounded edge opening to prevent spilling and reduce odors
- Ergonomic handle designed to easily attach to most hospital beds
- Easy to read dual graduate milliliter and ounce measurement markings on two sides



Item No.	Description	Quantity
PPB6100	Urinal Male Premium 1000 Milliliter Sturdy Plastic	1/EA, 50 EA/CA

Emesis Bags and Dispensers



- Emesis Bags
 - Aids in collecting emesis and prevents spills and contact with liquid
 - Excellent for patient use during transport and while in emergency waiting areas
 - Hard plastic ring around the top keeps the bag open during use
 - Twisted top hooks into one of the notches in the funnel to seal the bag for disposal
 - Graduated markings to measure ounce, cubic centimeters, and milliliters; ounce measures in 4 ounce increments
 - PP4201 is biodegradable and clear for easy viewing of contents
 - Single patient use
- Emesis Bag Dispenser
 - Made of durable polyethylene rubber
 - Top-lid covers emesis bags, aiding in infection control
 - Securely holds up to 24 folded Emesis Bags
 - Easy to assemble, wall-mounted dispenser



Item No.	Description	Quantity
PP4200	Bag Emesis Blue 6.5 Inch x 15 Inch x 1.35 Millimeter	1/EA, 24 EA/SL, 6 SL/CA
PP4201	Bag Emesis Biodegradable Clear 6.5 Inch x 15 Inch x 1.35 Millimeter	1/EA, 24 EA/SL, 6 SL/CA
PP4200D	Dispenser For Emesis Bag Polyethylene	1/EA, 12 EA/CA

PATIENT BEDSIDE

Bedside Bag

NOT MADE WITH
natural rubber latex

- Available in plastic
- Size: 6.5 Inch x 12 Inch
- Pleated sides for more capacity
- Adhesive straps attach bag to bed or nightstand



Item No.	Description	Quantity
PC7181	Bag Patient Bedside Plastic with Straps 6.5 x 12 Inch	1/EA, 100 EA/PK, 20 PK/CA

Eye Glass Cases

NOT MADE WITH
natural rubber latex

- Stores and protects glasses from breakage and scratching
- Hard durable storage case is made of sturdy polypropylene for protecting glasses
- Hard case is translucent neon green for easy visibility and location
- Hard case has adjustable, expandable depth of 1.5–2.5 Inches to accommodate various eyeglass sizes
- Expandable closure locks in place to allow for larger sized glasses and secure closure
- Soft sleeve case is made of polyurethane 7 Inch x 3.5 Inch with 2 millimeter microfiber drawstring closure



Item No.	Description	Quantity
Case Eyeglasses		
EGC1000	Hard Plastic Neon Green 3 Inch x 6 Inch with Expandable Depth 1.5 To 2.5 Inches	1/EA, 36 EA/CA
EGC1010	Soft Black 7 Inch x 3.5 Inch with 2 Millimeter Drawstring Closure	1/EA, 10 EA/BG, 10 BG/BX, 10 BX/CA

Patient Belonging Bags

NOT MADE WITH
natural rubber latex

- "Patient Belongings" printed in blue on one side and a place for name and room number
- 1.6 millimeter thickness
- Choose from full or standard size options
- Drawstring or plastic handles
- Durable polypropylene construction
- Customized options available



Item No.	Description	Quantity
Bag Patient Belonging Polypropylene		
PC7190	With Handle Opaque White 20 Inch x 18.5 Inch	1/EA, 25 EA/PK, 10 PK/CA
PC7191	With Handle Clear 20 Inch x 18.5 Inch	1/EA, 25 EA/PK, 10 PK/CA
PC7200	Drawstring Opaque White 20 Inch x 18.5 Inch	1/EA, 50 EA/PK, 5 PK/CA
PC7201	Drawstring Clear 20 Inch x 18.5 Inch	1/EA, 50 EA/PK, 5 PK/CA
PC7202C1	Drawstring Single Patient Use 20 Inch x 18.5 Inch Green Non Sterile	1/EA, 50 EA/PK, 5 PK/CA

PATIENT BEDSIDE

Denture Cups

NOT MADE WITH
natural rubber latex

- Teal cup with opaque white lid or fluorescent green with translucent clear lid
- Lid surface friendly for writing patient's name
- Convenient 8 ounce size
- Durable polypropylene construction
- 25 cups per pack; 10 packs per case



Item No.	Description	Quantity
Cup Denture 8 Ounce		
DC5000	Teal with Opaque Lid	1/EA, 25 EA/SL, 10 SL/CA
DC5001	Fluorescent Green with Opaque Lid	1/EA, 25 EA/SL, 10 SL/CA

Personal Belonging Box

- Helps prevent patients' personal items from being misplaced or damaged
- Separate compartments for personal items such as, denture cup, wallet, smartphones, glasses, prescriptions or hearing aids
- Ideal for use in hospital ER, nursing floors, Long Term Care facilities, Outpatient Surgery Centers and Rehabilitation Units
- Three compartment, one piece design with hinged locking lid to prevent missing lids or leftover pieces
- Clear design ensures patient can see where personal aids are first thing in the morning or at night eliminating the confusion associated with missing items; improving patient satisfaction through reducing lost patient valuables
- Designed to be nested or stacked for storage in supply closets or nursing stations
- Patient information and room number can be written on the box with a permanent marker or skin marker



Item No.	Description	Quantity
Box Patient Belonging Polypropylene		
PBBX100	Box Personal Belonging Clear Plastic 3 Compartment	1/EA, 130 EA/CA

QUIET KITS



Creating a quiet and relaxing environment for your patients is an important part of the healing process. We understand your need to provide the best experience for your patients, their families, and their care team.

MediChoice Quiet Kits were designed with this in mind. Each kit provides patient's amenities that will help them feel comfortable, calm and well-rested. Among the amenities in the MediChoice Quiet Kits are a soft, black eye mask with elastic headband to reduce brightness levels, ear plugs to muffle unwanted noise, "quiet please" door hanger for patients to get rest, SPF-15 lip balm to prevent cracking, a notebook for care questions and instructions, and crosswords and Sudoku puzzles.

Customization Options

MediChoice offers customizable Quiet Kits to help improve the patient care experience and strengthen your healthcare facility's brand name. The kits can be tailored to meet your specific needs.

- Items in the kit come pre-packed by MediChoice.
- Contents can be customized to meet your patients' needs.
- Kits can be designed with your name, logo and colors.

Quiet Kits



- Quiet kit to help patient rest quietly and add to patient comfort
- ADKT101 Kit includes:
 - Black soft eye mask with elastic band to stay in place
 - Yellow, non-toxic ear plugs
 - Moisturizing Lip Balm SPF 15; made in USA
 - 4 x 6 Inch Notebook for patient care questions and instructions
 - 4 x 6 Inch puzzle booklet with Sudoku and Crossword activities
 - 3 Inch Pencil
- ADKT102 Kit includes:
 - All of the above
 - "Quiet Please" Door hang tag



ADKT101 Kit



ADKT102 Kit

Item No.	Description	Quantity
ADKT101	With Eye Mask Lip Balm Ear Plugs Sudoku Book Question Book Pencil	1/EA, 50 EA/CA
ADKT102	Containing Eye Mask Quiet Please Door Hanger Lip Balm Sudoku Book Question Book and Pencil	1/EA, 50 EA/CA

PATIENT BEDSIDE

Facial Tissues



- Standard bedside pack (FT0040)
 - 5.5 Inch x 7 Inch
 - 40 sheets per box
 - Soft, post-consumer recycled paper
- Premium bedside pack (FT0060)
 - 8 Inch x 6.5 Inch
 - 60 sheets per box
 - Soft, strong virgin fibers



Item No.	Description	Quantity
Tissue Facial 2 Ply		
FT0040	Standard Bedside Pack 5.7 x 7 Inch	40/BX, 200 BX/CA
FT0060	Premium Bedside Pack 8 x 5.6 Inch	60/BX, 48 BX/CA

Earbuds and Headsets

- Black
- Disposable
- 6 foot cord
- ROHS compliant; free of hazardous materials such as Lead, Cadmium, Mercury, Hexavalent Chromium, Polybrominated Biphenyls and Polybrominated Diphenyl Ethers



Item No.	Description	Quantity
EBD500	Ear Bud Black Disposable 3.5 Mil Adapter 6 Foot Cord	1/EA, 500 EA/CA
HDS100	Headset Black Disposable 3.5 Millimeter Adapter with 6 Foot Cord	1/EA, 100 EA/CA

Please see our full line of Nursery Kits on page 128 in the Maternity & Baby Care Section.



ANESTHESIA

Designed to Provide Maximum Patient Comfort



Anesthesia Masks

NOT MADE WITH
natural rubber latex

- Features clear, flexible faceplates and soft elastomeric cushions for maximum patient comfort
- Unique ergonomic shape fits all facial configurations for improved sealing capabilities
- Soft, pliable air-filled cushion with an inflation valve; the cushion can be adjusted (inflated/deflated) to suit the patient's anatomical contour
- ISO standard fitting which is compatible with all standard catheter mounts and breathing systems with 22mm (adult) and 15mm (pediatric) connectors
- Offers additional protection against the possibility of corneal abrasion due to the contour of the mask
- All masks are individually bagged with color-coded labels for size identification
- Non-sterile
- Single patient use

Item No.	Description	Quantity
Mask Anesthesia With Hook Ring Inflation Valve		
AM1000	Adult Large Non-Sterile	1/EA, 30 EA/CA
AM1001	Adult Medium Non-Sterile	1/EA, 30 EA/CA
AM1002	Adult Small Non-Sterile	1/EA, 20 EA/CA
AM1003	Pediatric Non-Sterile	1/EA, 20 EA/CA
Mask Anesthesia		
AM1004	With Inflation Valve Infant Non-Sterile	1/EA, 20 EA/CA
AM1005	Clear Flexible Faceplate With Inflation Valve Neonate Non-Sterile	1/EA, 20 EA/CA

Berman Airways

NOT MADE WITH
natural rubber latex

- Satin-smooth finish is safe and comfortable
- Size marked on bite block for easy identification
- Available in seven sizes, packaged in bulk or individually
- Manufactured in polypropylene
- Non-sterile
- Single patient use

Item No.	Description	Quantity
Airway Berman Satin-Smooth Finish Polypropylene Bulk		
BA1240	40 Millimeter Non-Sterile	1/EA, 12 EA/BX, 4 BX/CA
BA1250	50 Millimeter Non-Sterile	1/EA, 12 EA/BX, 4 BX/CA
BA1260	60 Millimeter Non-Sterile	1/EA, 12 EA/BX, 4 BX/CA
BA1270	70 Millimeter Non-Sterile	1/EA, 12 EA/BX, 4 BX/CA
BA1280	80 Millimeter Non-Sterile	1/EA, 12 EA/BX, 4 BX/CA
BA1290	90 Millimeter Non-Sterile	1/EA, 12 EA/BX, 4 BX/CA
BA12100	100 Millimeter Non-Sterile	1/EA, 12 EA/BX, 4 BX/CA
Airway Berman Satin-Smooth Finish Individually Bagged Polypropylene		
BA5040	40 Millimeter Non-Sterile	1/EA, 50 EA/CA
BA5050	50 Millimeter Non-Sterile	1/EA, 50 EA/CA
BA5060	60 Millimeter Non-Sterile	1/EA, 50 EA/CA
BA5070	70 Millimeter Non-Sterile	1/EA, 50 EA/CA
BA5080	80 Millimeter Non-Sterile	1/EA, 50 EA/CA
BA5090	90 Millimeter Non-Sterile	1/EA, 50 EA/CA
BA50100	100 Millimeter Non-Sterile	1/EA, 50 EA/CA



Color Coded Guedel Airways

NOT MADE WITH
natural rubber latex

- Smooth, plastic finish allows simple insertion, minimizing tissue trauma
- Rigid bite block avoids occlusion and tongue bites
- Color coded for size identification
- Size printed on airway
- Packaged in bulk or individually
- Automation
- Non-sterile
- Single patient use

Item No.	Description	Quantity
Airway Guedel Plastic Finish Bulk		
GA1050	Blue 50 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA1060	Black 60 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA1070	White 70 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA1080	Green 80 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA1090	Yellow 90 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA10100	Red 100 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA

Airway Guedel Plastic Finish Individually Bagged Single Patient Use

GA5050	Blue 50 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA5060	Black 60 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA5070	White 70 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA5080	Green 80 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA5090	Yellow 90 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA
GA50100	Red 100 Millimeter Non-Sterile	1/EA, 10 EA/BX, 10 BX/CA

Nasopharyngeal Airways

sterile R

- Soft bevel-tip airway is comfortable and easy to insert
- Available in six sizes, including pediatric
- Packaged in sterile, individual pouches
- Available in latex PVC or free of natural rubber latex PVC



Item No.	Description	Quantity
Airway Nasopharyngeal Soft Bevel-Tip Polyvinyl Chloride Contains Natural Rubber Latex		
123526	26 French Sterile	1/EA, 10 EA/BX
123528	28 French Sterile	1/EA, 10 EA/BX
123530	30 French Sterile	1/EA, 10 EA/BX
123532	32 French Sterile	1/EA, 10 EA/BX
123534	34 French Sterile	1/EA, 10 EA/BX
123536	36 French Sterile	1/EA, 10 EA/BX
Airway Nasopharyngeal Soft Bevel-Tip Polyvinyl Chloride Not Made With Natural Rubber Latex		
123626	26 French Sterile	1/EA, 10 EA/BX
123628	28 French Sterile	1/EA, 10 EA/BX
123630	30 French Sterile	1/EA, 10 EA/BX
123632	32 French Sterile	1/EA, 10 EA/BX
123634	34 French Sterile	1/EA, 10 EA/BX
123636	36 French Sterile	1/EA, 10 EA/BX

Pressure Infuser Bags

NOT MADE WITH
natural rubber latex

- Unique design maintains consistent infuser pressure
- Nylon netting makes loading easy and permits a clear view of fluid
- Large bore stopcock with quick release for rapid inflation and deflation
- Safety valve prevents over-inflation of bag
- Color-coded gauge and sizing is easy to read and maintains reliable pressure
- Pressure valve releases at 325 mm Hg
- Available in three sizes: 500cc, 1000cc and 3000cc
- Non-sterile
- Single Patient Use



Item No.	Description	Quantity
Bag Infuser Pressure		
230500	500 Cubic Centimeter Nylon Disposable	1/EA, 5 EA/BX, 5 BX/CA
231000	1000 Cubic Centimeter Nylon Disposable	1/EA, 5 EA/BX, 5 BX/CA
30500	500 Cubic Centimeter Nylon Disposable	1/EA, 2 EA/BX, 5 BX/CA
31000	1000 Cubic Centimeter Nylon Disposable	1/EA, 5 EA/BX, 5 BX/CA
33000	3000 Cubic Centimeter Nylon Disposable	1/EA, 5 EA/BX, 5 BX/CA

Exam Capes

made in the USA **NOT MADE WITH natural rubber latex**

- Capes are built to resist rips and tears and are intended for multi-phasic and extended exams while providing modesty
- Reversible front/back opening
- Available in absorbent 3-ply tissue or tissue/polyethylene/tissue (TPT) for added strength and soak-through protection
- T-P-T available in three colors: blue, white or mauve
- Unisize: 30 Inch × 21 Inch



Item No.	Description	Quantity
Cape Examination		
MC910516	Tissue/Poly/Tissue Mauve 30 × 21 Inch	1/EA, 100 EA/CA
MC11105	Tissue/Poly/Tissue Blue 30 × 21 Inch	1/EA, 100 EA/CA
MC11106	3 Ply Tissue White 30 × 21 Inch	1/EA, 100 EA/CA
MC11115	Tissue/Poly/Tissue White 30 × 21 Inch	1/EA, 100 EA/CA

Exam Gowns

made in the USA **NOT MADE WITH natural rubber latex**

- Available in absorbent 3-ply tissue, tissue/polyethylene/tissue (TPT), or SCRIM reinforced material
- Reversible front/back opening with stretchable waist tie
- Unisize: 30 Inch × 42 Inch provides patient comfort, hygiene and modesty with knee length coverage
- Also available: 2XL, Scrim Reinforced, 72 Inch × 45 Inch
- Premium Tissue-Poly-Tissue available in blue, white or mauve
- SCRIM reinforced gown best for extended wear during multi-phasic procedures



Item No.	Description	Quantity
Gown Examination, Knee Length Coverage		
MC910536	Tissue/Poly/Tissue Mauve 30 × 42 Inch	1/EA, 50 EA/CA
MC980830	Pediatric Tissue/Poly/Tissue Blue 21 × 36 Inch	1/EA, 50 EA/CA
MC11108	Tissue/Poly/Tissue Blue 30 × 42 Inch	1/EA, 50 EA/CA
MC18838	Universal Scrim Reinforced Blue 30 × 42 Inch	1/EA, 50 EA/CA
MC80831	Glued Scrim Reinforced Blue 30 × 42 Inch	1/EA, 50 EA/CA
MC18846	Premium Tissue/Poly/Tissue White 30 × 42 Inch	1/EA, 50 EA/CA
MC18156	2X Large Scrim Reinforced Blue 72 × 45 Inch	1/EA, 25 EA/CA
MC11107	3 Ply Tissue/Poly White 30 × 42 Inch	1/EA, 50 EA/CA
MC18842	Premium Tissue/Poly/Tissue Blue 30 × 42 Inch	1/EA, 50 EA/CA
MC18844	X-Ray 3 Ply Tissue/Poly Blue 30 × 42 Inch	1/EA, 50 EA/CA

Disposable Exam Shorts

NOT MADE WITH natural rubber latex

- Durable SMS nonwoven material
- Superior softness and opacity
- Elastic waist band ensures secure fit
- Ideal as a modesty cover or an activity short



Item No.	Description	Quantity
Short Examination Elastic Waist Band		
181101	Medium SMS Blue Disposable	1/EA, 25 EA/BG, 4 BG/CA
181102	Large SMS Blue Disposable	1/EA, 25 EA/BG, 4 BG/CA
181103	Extra-Large SMS Blue Disposable	1/EA, 25 EA/BG, 4 BG/CA
181104	2XL SMS Blue Disposable	1/EA, 25 EA/BG, 4 BG/CA



Bibs

made in the USA **NOT MADE WITH natural rubber latex**

- Slips over head to protect patient's apparel during feeding
- Available in two styles
 - 18 Inch × 32 Inch is constructed of tissue-poly for durability, absorption and a moisture barrier
 - 16 Inch × 24 Inch is constructed of thick poly providing durability and a moisture barrier with a debris-catching pocket along bottom edge
- Solid white or blue in color
- Bulk-packed

Item No.	Description	Quantity
Bib Patient Tissue Poly		
MBIB1003	Blue 20 × 29 Inch Extra-Large	1/EA, 500 EA/CA
MBIB1001	White 18 × 32 Inch Adult	1/EA, 300 EA/CA
MC18964	White 16 × 32 Inch Adult	1/EA, 300 EA/CA
MBIB1002	White 16 × 24 Inch Adult With Pocket	1/EA, 500 EA/CA

SLIPPERS

Textured Side Out Slippers



- Soft terry construction
- Skid-resistant, large-coverage treads
- Color coded by size
- Available in double tread (top and bottom) to ensure grip if sock rotates on the foot



Item No.	Description	Quantity
Slipper Patient Slip-Resistant Sole Single Tread		
66600	Small Child Soft Terry Yellow	1/PR, 48 PR/CA
66601	Medium Soft Terry Blue	1/PR, 48 PR/CA
66602	Large Soft Terry Beige	1/PR, 48 PR/CA
66603	Extra-Large Soft Terry Gray	1/PR, 48 PR/CA
66604	2X Large Soft Terry Blue	1/PR, 48 PR/CA
Slipper Patient Slip-Resistant Double Tread		
66699	Infant Toddler Soft Terry White	1/PR, 48 PR/CA
66608	Small Child Soft Terry Yellow	1/PR, 48 PR/CA
66609	Medium Soft Terry Blue	1/PR, 48 PR/CA
66605	Large Soft Terry Beige	1/PR, 48 PR/CA
66607	Extra-Large Soft Terry Gray	1/PR, 48 PR/CA
66610	2X Large Soft Terry Blue	1/PR, 48 PR/CA
66611	Bariatric Soft Terry Beige	1/PR, 48 PR/CA

Smooth Side Out Slippers



- Soft terry pile side of the slipper sock is on the inside for patient comfort and safety
- Large coverage, skid resistant, puff ink tread
- Available with double tread (top and bottom) to ensure grip regardless of how slipper is worn
- Color coded by size



Item No.	Description	Quantity
Slipper Patient Single Tread Smooth Side Out		
SLP1001	Toddler Soft Terry Pile White	1/PR, 48 PR/CA
SLP1002	Small Child Soft Terry Pile Yellow	1/PR, 48 PR/CA
SLP1003	Medium Soft Terry Pile Blue	1/PR, 48 PR/CA
SLP1004	Large Soft Terry Pile Beige	1/PR, 48 PR/CA
SLP1005	Extra-Large Soft Terry Pile Gray	1/PR, 48 PR/CA
SLP4005	Extra-Large Soft Terry Pile Red	1/PR, 48 PR/CA
SLP1006	2X Large Soft Terry Pile Blue	1/PR, 48 PR/CA
Slipper Patient Double Tread Smooth Side Out		
SLP2001	Toddler Soft Terry Pile White	1/PR, 48 PR/CA
SLP2002	Small Child Soft Terry Pile Yellow	1/PR, 48 PR/CA
SLP2003	Medium Soft Terry Pile Blue	1/PR, 48 PR/CA
SLP2004	Large Soft Terry Pile Beige	1/PR, 48 PR/CA
SLP4004	Large Soft Terry Pile Red	1/PR, 48 PR/CA
SLP2005	Extra-Large Soft Terry Pile Gray	1/PR, 48 PR/CA
SLP2006	2X Large Soft Terry Pile Blue	1/PR, 48 PR/CA
SLP2007	Bariatric Soft Terry Pile	1/PR, 48 PR/CA

SLIPPERS

Fall Management Slippers



- Bright colors visually alert staff to fall risk patients
- Sizes Small to Bariatric
- Double tread (top and bottom) to reduce fall risk if slipper rotates on the foot
- Large coverage, skid resistant, puff ink tread
- Soft terry cloth construction for optimal comfort



Item No.	Description	Quantity
Slipper Patient Fall Management Double-Treaded Smooth Side Out		
SLP6DY	Toddler Soft Terrycloth Yellow	1/PR, 48 PR/CA
SLP7DG	Small Soft Terrycloth Green	1/PR, 48 PR/CA
SLP7DO	Small Soft Terrycloth Orange	1/PR, 48 PR/CA
SLP7DR	Small Soft Terrycloth Red	1/PR, 48 PR/CA
SLP7DY	Small Soft Terrycloth Yellow	1/PR, 48 PR/CA
SLP8DG	Medium Soft Terrycloth Green	1/PR, 48 PR/CA
SLP8DO	Medium Soft Terrycloth Orange	1/PR, 48 PR/CA
SLP8DR	Medium Soft Terrycloth Red	1/PR, 48 PR/CA
SLP8DY	Medium Soft Terrycloth Yellow	1/PR, 48 PR/CA
SLP12DG	Large Soft Terrycloth Green	1/PR, 48 PR/CA
SLP12DO	Large Soft Terrycloth Orange	1/PR, 48 PR/CA
SLP12DR	Large Soft Terrycloth Red	1/PR, 48 PR/CA
SLP12DY	Large Soft Terrycloth Yellow	1/PR, 48 PR/CA
SLP14DG	Extra-Large Soft Terrycloth Green	1/PR, 48 PR/CA
SLP14DO	Extra-Large Soft Terrycloth Orange	1/PR, 48 PR/CA
SLP14DP	Extra-Large Soft Terrycloth Purple	1/PR, 48 PR/CA
SLP14DR	Extra-Large Soft Terrycloth Red	1/PR, 48 PR/CA
SLP14DY	Extra-Large Soft Terrycloth Yellow	1/PR, 48 PR/CA
SLP16DG	2X Large Soft Terrycloth Green	1/PR, 48 PR/CA
SLP16DO	2X Large Soft Terrycloth Orange	1/PR, 48 PR/CA
SLP16DR	2X Large Soft Terrycloth Red	1/PR, 48 PR/CA
SLP16DY	2X Large Soft Terrycloth Yellow	1/PR, 48 PR/CA
SLP18DG	Bariatric Soft Terrycloth Green	1/PR, 48 PR/CA
SLP18DO	Bariatric Soft Terrycloth Orange	1/PR, 48 PR/CA
SLP18DR	Bariatric Soft Terrycloth Red	1/PR, 48 PR/CA
SLP18DY	Bariatric Soft Terrycloth Yellow	1/PR, 48 PR/CA

Universal Slippers



- Unique Weave fits most adult feet
- Available in blue and yellow
- Double tread design ensures tread contact with floor even if the slipper rotates on the foot



Item No.	Description	Quantity
Slipper Patient Double-Treaded Universal		
SLP01UB	Blue	1/PR, 48 PR/CA
SLP01UY	Yellow	1/PR, 48 PR/CA

Please see additional Fall Management products, including our full line of Bathroom Safety Aids on pages 131-132 in the Mobility Aids Section.

SLIPPERS

Hard Sole Slippers

NOT MADE WITH
natural rubber latex

- Soft, terry cloth construction provides optimal comfort
- Firm, flexible sole helps protect against accidental falls
- Above-ankle profile with elastic band to help ensure a snug fit
- Child and adult sizes
- Each pair is individually packaged



Item No.	Description	Quantity
Slipper Patient Hard Sided Sole Soft Terrycloth Upper		
SLPH3001	Child Size 4-5 White	1/PR, 48 PR/CA
SLPH3002	Small Size 5-6 Blue	1/PR, 48 PR/CA
SLPH3003	Medium Size 7-8 Blue	1/PR, 48 PR/CA
SLPH3004	Large Size 9-10 Blue	1/PR, 48 PR/CA
SLPH3005	Extra-Large Size 11-12 Blue	1/PR, 48 PR/CA

TELEMETRY POUCHES

Clear Polyethylene Telemetry Pouches

NOT MADE WITH
natural rubber latex

diverse
supplier

- Clear polyethylene pouch offers convenience, patient comfort and protection for the transmitter.
- 5.25 Inch × 8.125 Inch pouch with clip helps protect transmitter from fluids
- Hook and loop closure easily allows transmitter wires to pass through top of pouch
- Pre-attached patient ID label



Item No.	Description	Quantity
MC500	Pouch Telemetry Clear Velcro Hook and Loop Closure Pre-Attached Patient ID Label Polyethylene 5.25 Inch × 8.125 Inch	1/EA, 50 EA/CA

Non-Woven Telemetry Pouches with Neck Straps

NOT MADE WITH
natural rubber latex

diverse
supplier

- Non-woven pouch offers convenience, patient comfort and protection for the transmitter.
- Straps that tie at the neck only
- 6 Inch × 9.25 Inch pouch with non-woven back and polyethylene front
- Hook and loop closure
- Pre-attached patient ID label



Item No.	Description	Quantity
MC901	Pouch Telemetry Non-Woven Hook and Loop Closure with Neck Strap 6 Inch × 9.25 Inch	1/EA, 50 EA/CA

TELEMETRY POUCHES

Foam Telemetry Pouches

NOT MADE WITH
natural rubber latex

diverse
supplier

- Foam pouch offers convenience, patient comfort and protection for the transmitter.
- Soft foam is comfortable for patient
- 4.75 Inch × 6 Inch fits most standard size monitors
- Ideal for small digital “beeper” size monitors
- Pre-attached patient ID label



Item No.	Description	Quantity
MC902	Pouch Telemetry Soft Foam Velcro Hook Loop Waist Belt Pre-Attached Patient ID Label 4.75 Inch × 6 Inch	1/EA, 50 EA/CA

Non-Woven Telemetry Pouches with Neck and Waist Straps

NOT MADE WITH
natural rubber latex

diverse
supplier

- Non-woven pouch offers convenience, patient comfort and protection for the transmitter.
- Straps that tie at the neck and waist
- 6 Inch × 9.25 Inch pouch with non-woven back and polyethylene front
- Hook and loop closure
- Pre-attached patient ID label



Item No.	Description	Quantity
MC903	Pouch Telemetry Non-Woven Hook Loop with Neck and Waist Strap Pre-Attached Patient ID Label 6 Inch × 9.25 Inch	1/EA, 50 EA/CA

FALL MANAGEMENT BLANKETS

Fall Management Blanket

NOT MADE WITH
natural rubber latex

diverse
supplier

- Bright Yellow color to help identify patients at high risk for falls.
- 100 percent Brushed Polyester
- Size: 30 Inch × 40 Inch
- Provides warmth and comfort to patient
- For use in bed or wheelchair
- Machine Washable



Item No.	Description	Quantity
FPB040	Blanket Fall Management 30 Inch × 40 Inch Polyester Red	1/EA, 60 EA/CA
FPB060	Blanket Fall Management 30 Inch × 40 Inch Polyester Yellow	1/EA, 60 EA/CA

Please see additional Fall Management products, including our Gait Belts on page 150 in the Ortho/Rehab Section and our Turn and Positioning Slider Sheets on page 169 in the Skin Care Section.



APPAREL—STAFF



At Owens & Minor we understand your need to have a full-line of Disposable Isolation Gowns and Face Masks that you can depend on to keep your staff and patients safe, without compromising quality, comfort, or durability. Our comprehensive line of MediChoice gowns come in a variety of back and wrist closures, with generous sizing, color options, and five levels of protection to match all levels of risk—from our least-expensive general cover gown, to the popular non-AAMI isolation gown, all the way up to the poly-coated, fluid-resistant gowns offering AAMI 3 and AAMI 4 protection. MediChoice face masks range from standard procedure masks to surgical masks with non-fogging visors, offering fluid resistance up to 160mm Hg.

Comfort Drives Compliance

All MediChoice gowns and face masks are made for comfort and ease of donning and removal, helping to reduce HAIs. With non-reimbursement, these savings go directly to your bottom line.

Savings Without Sacrifice

MediChoice offers a range of gowns and face mask options at a low cost—but you don't sacrifice what matters most: protection and comfort. By standardizing to MediChoice gowns and face masks, you simplify your supply chain with reduced SKUs and inventory to manage.

Disposable Aprons

NOT MADE WITH
natural rubber latex

- Bib-style aprons have die-cut head openings and extra-long ties
- Good utility apron; variety of uses
- Light but strong embossed polyethylene
- Convenient dispenser box (100 aprons/box)
- 30 Inch × 44 Inch Size



Item No.	Description	Quantity
A7779	Utility Apron Embossed Bib-Style Extra Long Tie Polyethylene White 30 × 44 Inch Disposable	1/EA, 100 EA/BX, 5 BX/CA

Bouffant Caps

NOT MADE WITH
natural rubber latex

- Fluid-resistant, lightweight spunbond polypropylene
- Soft, non-binding, latex-free elastic band
- 21 Inch and 24 Inch sizes
- Conveniently packaged in dispenser box (100 each/box)



Item No.	Description	Quantity
Cap Bouffant Lightweight Spunbond Polypropylene Non-Sterile Disposable		
77701	Blue 21 Inch	1/EA, 100 EA/BX, 5 BX/CA
77702	Blue 24 Inch	1/EA, 100 EA/BX, 5 BX/CA
77704	Red 24 Inch	1/EA, 100 EA/BX, 5 BX/CA
77703	White 24 Inch	1/EA, 100 EA/BX, 5 BX/CA

Surgeon Caps

NOT MADE WITH
natural rubber latex

- Made of absorbent material
- Soft, breathable spunbond polypropylene top
- Secures with full-length ties
- Conveniently packaged in dispenser box (100 each/box)



Item No.	Description	Quantity
77706	Cap Surgeon with Ties Scrim Material with Polyspun Bond Blue	1/EA, 100/BX, 500/CA

Surgeon Hoods

NOT MADE WITH
natural rubber latex

- Made of absorbent SMS material
- Soft, breathable spunbond polypropylene top
- Secures with full-length ties
- Conveniently packaged in dispenser box (100 each/box)



Item No.	Description	Quantity
77707	Hood Surgeon SMS Blue Ties About Neck	1/EA, 100 EA/BX, 3 BX/CA

Complementary Products

Item No.	Description
77702	Cap Bouffant Lightweight Spunbond Polypropylene Blue 24 Inch Non-Sterile Disposable
77724	Cover Shoe Fluid-Resistant Slip-Resistant Durable Elastic Band Universal Spunbond Meltblown Spunbond Polypropylene Blue Non-Sterile
77727	Cover Shoe Fluid-Resistant Non-Skid Tread Durable Elastic Band Upper Extra-Large Spun Bond Polypropylene Blue
77728	Cover Shoe Fluid-Resistant Durable Elastic Band Upper Universal Spun Bond Polypropylene Blue
77701	Cap Bouffant Lightweight Spunbond Polypropylene Blue 21 Inch Non-Sterile Disposable
77725	Cover Shoe Fluid-Resistant Slip-Resistant Durable Elastic Band Extra-Large Spunbond Meltblown Spunbond Polypropylene Blue Non-Sterile
77726	Cover Shoe Non-Skid Tread Fluid-Resistant Durable Elastic Band Upper Universal Spun Bond Polypropylene Blue
77729	Cover Shoe Fluid-Resistant Durable Elastic Band Upper Extra-Large Spun Bond Polypropylene Blue

Protective Eyewear

NOT MADE WITH
natural rubber latex

- Lightweight and comfortable for long term wear; weighs only .3 ounces
- Available in pre-assembled, ready-to-wear option or as two separate components (lens and frame)
- Wrap-around design provides both front and side eye protection
- Reusable frames available in 10 stylish colors
- Optical-grade clear replacement lenses
- Fits comfortably with a surgical mask and over prescription eyewear



Item No.	Description	Quantity
Eyewear Protective Non-Sterile		
PEW050	Assembled Multi Color	1/EA, 10 EA/BG, 5 BG/CA
PEW100	Frame Reusable Multi Color	1/EA, 10 EA/BG, 10 BG/CA
PEW250	Replacement Lenses Single Use	1/EA, 25 EA/PK, 10 PK/CA
FS1100F	Shield Face Full Length Visor with Anti-Fog and Anti-Glare Coating Elastic Headband Foam Header Plastic Clear	1/EA, 24 EA/BX, 10 BX/CA

Please see our full line of Personal Protection Equipment (PPE) Organizers on pages 28-29 in the Central Supply Section.



Please see our full line of Specialty Dispensers and Accessories on page 30 in the Central Supply Section.



Disposable Scrubs

NOT MADE WITH
natural rubber latex

- Made of heavy durable, dark blue SMS material
- Tops have hemmed, V-Neck opening and three pockets
- Pants have drawstring or elastic waist and open leg cuffs
- Tops and pants sold separately, sizes M–2XL
- Individually packaged



Item No.	Description	Quantity
----------	-------------	----------

Pant Scrub Drawstring Waist Open Leg Cuff Spunbond Meltblown Spunbond Dark Blue Disposable

DSP01M	Medium	1/EA, 48 EA/CA
DSP01L	Large	1/EA, 48 EA/CA
DSP1XL	Extra-Large	1/EA, 48 EA/CA
DSP2XL	2XL	1/EA, 48 EA/CA

Top Scrub Hemmed Three-Pocket V-Neck Open Spunbond Meltblown Spunbond Dark Blue Disposable

DST01M	Medium	1/EA, 48 EA/CA
DST01L	Large	1/EA, 48 EA/CA
DST1XL	Extra-Large	1/EA, 48 EA/CA
DST2XL	2X	1/EA, 48 EA/CA

Pant Scrub Elastic Waist Open Leg Cuff Spunbond Meltblown Spunbond Blue Disposable

EDSP01M	Medium	1/EA, 48 EA/CA
EDSP01L	Large	1/EA, 48 EA/CA
EDSP1XL	Extra-Large	1/EA, 48 EA/CA
EDSP2XL	2XL	1/EA, 48 EA/CA

Disposable Staff Jackets

NOT MADE WITH
natural rubber latex

- Made of Fluid Resistant SMS nonwoven material
- Two color options: blue or white
- Knit collar, knit cuffs, snap front
- Available in sizes Small to 2XL
- Packaged CA/25 (5 bags of 5 ea)



Item No.	Description	Quantity
----------	-------------	----------

Jacket Staff Open Collar Cuff Snap Front Fluid-Resistant Spunbond Meltblown Spunbond Disposable

SJ0602	Small Blue	1/EA, 25 EA/CA
SJ0603	Medium Blue	1/EA, 25 EA/CA
SJ0604	Large Blue	1/EA, 25 EA/CA
SJ0605	Extra-Large Blue	1/EA, 25 EA/CA
SJ0606	2XL Blue	1/EA, 25 EA/CA
SJ0502	Small White	1/EA, 25 EA/CA
SJ0503	Medium White	1/EA, 25 EA/CA
SJ0504	Large White	1/EA, 25 EA/CA
SJ0505	Extra-Large White	1/EA, 25 EA/CA
SJ0506	2X Large White	1/EA, 25 EA/CA

DISPOSABLE COVERALLS

Premium Disposable Coveralls

NOT MADE WITH
natural rubber latex

- Made of 3-layer blue or white Anti-Static SMS material
- Sizes Large to 4XL
- Knit cuffs, full zip front, open collar



Item No.	Description	Quantity
----------	-------------	----------

Coverall Anti-Static Premium Knit Cuff Full Front Zip Open Collar Open Ankles Spunbond Meltblown Spunbond Disposable

05503L	Large White	1/EA, 24 EA/CA
05503X	Extra-Large White	1/EA, 24 EA/CA
5503XX	2X Large White	1/EA, 24 EA/CA
5503XXX	3X Large White	1/EA, 24 EA/CA
5503XXXX	4X Large White	1/EA, 24 EA/CA
5503BL	Large Blue	1/EA, 24 EA/CA
5503BX	Extra-Large Blue	1/EA, 24 EA/CA
5503BXX	2X Large Blue	1/EA, 24 EA/CA
5503BXXX	3X Large Blue	1/EA, 24 EA/CA
5503BXXXX	4X Large Blue	1/EA, 24 EA/CA

Standard Disposable Coveralls

NOT MADE WITH
natural rubber latex

- Made of white PSB material
- Sizes Large (universal) to 4XL
- Knit cuffs, full zip front, open collar and legs



Item No.	Description	Quantity
----------	-------------	----------

Coverall Protective Standard Knit Cuff Full Front Zip Open Collar Open Ankles Polypropylene Disposable

055803L	Large White	1/EA, 24 EA/CA
055803X	Extra-Large White	1/EA, 24 EA/CA
55803XX	2X Large White	1/EA, 24 EA/CA
55803XXX	3X Large White	1/EA, 24 EA/CA

DISPOSABLE LAB COATS

Standard Disposable Lab Coats

NOT MADE WITH
natural rubber latex

- Made of white PSB material
- Open collar and knit cuffs
- 3 pockets, 5-snap front
- Sizes Small to 2XL
- White only



Item No.	Description	Quantity
----------	-------------	----------

Coat Laboratory Three-Pocket Knit Cuff 5 Snap-Front Standard Polypropylene White Disposable

03304S	Small White	1/EA, 25 EA/CA
03304M	Medium White	1/EA, 25 EA/CA
03304L	Large White	1/EA, 25 EA/CA
03304X	Extra-Large White	1/EA, 25 EA/CA
3304XX	2X Large White	1/EA, 25 EA/CA

DISPOSABLE LAB COATS

Premium Disposable Lab Coats

NOT MADE WITH
natural rubber latex

- Made of fluid resistant Antistatic SMS material
- Two color options: blue and white
- Open collar and knit cuffs
- 3 pockets, 5-snap front



Item No.	Description	Quantity
Coat Laboratory Anti-Static Three-Pocket Knit Cuff 5-Snap Premium Spunbond Meltblown Spunbond Disposable		
03504XS	Extra-Small White	1/EA, 25 EA/CA
03504S	Small White	1/EA, 25 EA/CA
03504M	Medium White	1/EA, 25 EA/CA
03504L	Large White	1/EA, 25 EA/CA
03504X	Extra-Large White	1/EA, 25 EA/CA
3504XX	2X Large White	1/EA, 25 EA/CA
3504XXX	3X Large White	1/EA, 25 EA/CA
03504BS	Small Blue	1/EA, 25 EA/CA
03504BM	Medium Blue	1/EA, 25 EA/CA
03504BL	Large Blue	1/EA, 25 EA/CA
03504BX	Extra-Large Blue	1/EA, 25 EA/CA
3504BXX	2X Large Blue	1/EA, 25 EA/CA
3504BXXX	3X Large Blue	1/EA, 25 EA/CA

FACE MASKS

Procedure Masks

NOT MADE WITH
natural rubber latex

- A variety of procedure masks offered at all levels of fluid resistance including 80 mm Hg, 120 mm Hg and 160 mm Hg as well as non-rated cone mask and child face cover.
- Soft, stretchy ear loops comfortably fit a wide variety of face sizes
- Full-width, aluminum nose piece will not kink, ensures a secure fit over the nose and cheekbones
- PRIMAGARD 160 procedure mask has anti-fog foam
- PRIMASOFT version available for those with skin sensitivities



Item No.	Description	Quantity
Mask Procedure Face		
PG4-1001	PRIMAGARD80 Isolation with Earloop Yellow	1/EA, 50 EA/BX, 10 BX/CA
PG4-1211	PRIMAGARD80 Non Sterile Blue	1/EA, 50 EA/BX, 10 BX/CA
PG4-1331	PRIMASOFT80 for Sensitive Skin White	1/EA, 50 EA/BX, 10 BX/CA
PG4-1092	PRIMAGARD120 Denim Blue	1/EA, 50 EA/BX, 10 BX/CA
PG4-1073	PRIMAGARD160 with Anti-Fog Foam Indigo	1/EA, 50 EA/BX, 10 BX/CA
PG4-0220	PRIMAGARD Cone Style with Elastic Strap Blue	1/EA, 50 EA/BX, 10 BX/CA
PG4-1680	Cover Face Procedure PRIMASOURUS with Dinosaur Print Pediatric	1/EA, 50 EA/BX, 10 BX/CA

Please see our full line of
Face Mask Dispensers & Respiratory
Hygiene Stations on pages 32-33
in the Central Supply Section.



FACE MASKS

Procedure Masks with Visor

NOT MADE WITH
natural rubber latex

- Highest protection from splashes at 160 mm Hg
- Soft, stretchy ear loops comfortably fit a variety of face sizes
- Full-width, heavy duty aluminum nose piece will not kink, and ensures a secure fit to keep visor steady and prevent glare
- Anti-fog foam helps reduce fogging, and anti-glare strip reduces reflections from below
- Visor material is optically clear with anti-fog, anti-static and anti-glare properties



Item No.	Description	Quantity
Mask Procedure with Face Shield EarLoops and Anti-Fog		
PG4-1563	PRIMAGARD160 Blue	1/EA, 25 EA/BX, 4 BX/CA
PG4-1573	PRIMAGARD160 Indigo	1/EA, 25 EA/BX, 4 BX/CA

Surgical Tie Masks

NOT MADE WITH
natural rubber latex

- A variety of surgical masks offered at two levels of fluid resistance, including 80 mm Hg and 160 mm Hg
- Full-width, heavy duty aluminum nose piece will not kink, and ensures a secure fit over the nose and cheekbones
- PRIMASOFT version available for those with skin sensitivities



Item No.	Description	Quantity
Mask Surgical Face Tie		
PG4-2001	PRIMAGARD80 Blue	1/EA, 50 EA/BX, 6 BX/CA
PG4-2321	PRIMAGARD80 Light And Cool Blue	1/EA, 50 EA/BX, 6 BX/CA
PG4-2331	PRIMASOFT80 for Sensitive Skin White	1/EA, 50 EA/BX, 6 BX/CA
PG4-2121	PRIMAGARD80 Soft Bill Blue	1/EA, 50 EA/BX, 6 BX/CA
PG4-2161	PRIMAGARD80 Soft Bill Green Large	1/EA, 50 EA/BX, 6 BX/CA
PG4-2561	PRIMAGARD80 with Horizontal Tie Blue	1/EA, 50 EA/BX, 6 BX/CA
PG4-3363	PRIMAGARD160 Soft Bill Large White	1/EA, 50 EA/BX, 6 BX/CA
PG4-2073	PRIMAGARD160 Indigo	1/EA, 50 EA/BX, 6 BX/CA
PG4-2333	PRIMASOFT160 for Sensitive Skin White	1/EA, 50 EA/BX, 6 BX/CA
PG4-6403	PRIMAGARD160 Plume Protector Anti-fog Foam Gray	1/EA, 50 EA/BX, 6 BX/CA
PG4-4051	PRIMAGARD80 with Adhesive Strip Green	1/EA, 50 EA/BX, 6 BX/CA

Surgical Masks with Visor

NOT MADE WITH
natural rubber latex

- Highest protection from splashes at 160 mm Hg
- Full-width, heavy duty aluminum nose piece will not kink, and ensures a secure fit to keep visor steady and prevent glare
- Anti-fog foam helps to reduce fogging, and anti-glare strip reduces reflections from below
- Visor material is optically clear with anti-fog, anti-static and anti-glare properties
- PG4-5033 available for those with skin sensitivities



Item No.	Description	Quantity
Mask Surgical Face Anti-Fog Foam Anti Glare Strip with Visor Tie Aluminum Nose Piece		
PG4-5033	PRIMASOFT160 for Sensitive Skin White	1/EA, 25 EA/BX, 4 BX/CA
PG4-5063	PRIMAGARD160 Blue	1/EA, 25 EA/BX, 4 BX/CA
PG4-5073	PRIMAGARD160 Blue	1/EA, 25 EA/BX, 4 BX/CA

FACE MASKS

Surgical Masks with Anti-Fog

NOT MADE WITH natural rubber latex

- A variety of surgical masks offered at two levels of fluid resistance, including 80 mm Hg and 160 mm Hg
- Full-width, heavy duty aluminum alloy nose piece will not kink, and ensures a secure fit over the nose and under the cheeks
- Anti-fog foam and film helps reduce condensation on eyewear
- PRIMASOFT version available for those with sensitive skin



Item No.	Description	Quantity
Mask Surgical Anti-Fog with Tie Aluminum Nose Piece		
PG4-4021	PRIMAGARD80 Denim Blue	1/EA, 50 EA/BX, 6 BX/CA
PG4-4073	PRIMAGARD160 Indigo	1/EA, 50 EA/BX, 6 BX/CA
PG4-4241	PRIMAGARD80 Green	1/EA, 50 EA/BX, 6 BX/CA
PG4-4331	PRIMASOFT80 For Sensitive Skin White	1/EA, 50 EA/BX, 6 BX/CA

SHOE COVERS

Standard Shoe Covers

NOT MADE WITH natural rubber latex

- Fluid-resistant, durable polypropylene
- Fully elasticized and nonconductive
- Universal and XL sizes
- Available with or without tread
- Conveniently packaged in 100 each dispenser box



Item No.	Description	Quantity
Cover Shoe Blue		
77726	Non-Skid Tread Universal	1/EA, 100 EA/BX, 3 BX/CA
77727	Non-Skid Tread Extra-Large	1/EA, 100 EA/BX, 3 BX/CA
77728	Universal	1/EA, 100 EA/BX, 3 BX/CA
77729	Extra-Large	1/EA, 100 EA/BX, 3 BX/CA

SMS Shoe Covers

NOT MADE WITH natural rubber latex

- Fluid-resistant, durable SMS Material
- Fully elasticized and nonconductive
- Universal and XL sizes
- Conveniently packaged in 100 each dispenser box



Item No.	Description	Quantity
Cover Shoe Slip-Resistant Blue Non-Sterile		
77724	Universal	1/EA, 100 EA/BX, 3 BX/CA
77725	Extra-Large	1/EA, 100 EA/BX, 3 BX/CA

Boot Covers

NOT MADE WITH natural rubber latex

- Fluid-resistant, durable polypropylene coated foot with breathable SMS top
- Knee High with elastic top and ankle to maximize fluid protection
- Universal and XL sizes
- Non-skid tread
- Conveniently packaged in 50 each dispenser box



Item No.	Description	Quantity
Boot Full Coverage Impervious Shoe Cover Blue Non-Sterile		
77721	Universal	1/EA, 50 EA/BX, 3 BX/CA
77722	Extra-Large	1/EA, 50 EA/BX, 3 BX/CA

FLUID RESISTANT GOWNS

Fluid-Resistant Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations requiring a moderate to high level of fluid protection
- Made from poly coated spunbond polypropylene material has been tested against ASTM F1670 Synthetic Blood Penetration Resistance
- Full coverage isolation gown
- Knit or elastic cuff options provide a secure fit under gloves
- Neck tie closure for individual size adjustment
- Available in yellow, blue or white
- Packaged 10 gowns/bag, 10 bags/case
- Universal and XL sizes
- Single use; nonsterile



Item No.	Description	Quantity
Gown Fluid Resistant PSBPE Full Back Tie Neck and Waist		
77785	Elastic Cuff Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77795	Knit Cuff Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77885	Elastic Cuff Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77895	Knit Cuff Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77685	Elastic Cuff White Universal	1/EA, 10 EA/BG, 10 BG/CA
77695	Knit Cuff White Universal	1/EA, 10 EA/BG, 10 BG/CA
77785XL	Elastic Cuff Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77795XL	Knit Cuff Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77895XL	Knit Cuff Blue XL	1/EA, 10 EA/BG, 10 BG/CA
77695XL	Knit Cuff White XL	1/EA, 10 EA/BG, 10 BG/CA

ISOLATION GOWNS

Cover Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations when little fluid exposure is expected
- Made from strong, comfortable and breathable spunbond polypropylene (PSB) material
- Full coverage cover gown
- Elastic cuffs provide a secure fit under gloves
- Neck tie closure for individual size adjustment
- Available in yellow or blue
- Packaged 10 gowns/bag, 10 bags/case
- Universal and XL sizes
- Single use, non-sterile



Item No.	Description	Quantity
Gown Cover PSB Full Back Elastic Cuff Tie Neck and Waist		
77760	Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77760L	Yellow Large	1/EA, 10 EA/BG, 10 BG/CA
77760XL	Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77860	Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77860XL	Blue XL	1/EA, 10 EA/BG, 10 BG/CA

ISOLATION GOWNS

AAMI Rated Isolation Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations requiring a low to moderate level of fluid protection
- Made from strong, fluid resistant and breathable spunbond-meltblown-spunbond (SMS) material
- Tested in accordance with the ANSI/AAMI PB70:2012 standard
- Ultrasonically welded seams and waist ties to prevent fluid penetration
- Full coverage isolation gowns
- Knit, elastic or thumbloop cuff options provide secure fit under gloves
- Neck tie, tape tab or velcro closure for individual size adjustment or apron style neck for rapid over the head donning ensures gown compliance
- Available in yellow or blue
- Packaged 10 gowns/bag, 10 bags/case
- Universal and XL sizes
- Single use, non-sterile



Item No.	Description	Quantity
----------	-------------	----------

Gown Isolation SMS Full Back Elastic Cuff

77751	AAMI Level 1 Tie Neck and Waist Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77751XL	AAMI Level 1 Tie Neck and Waist Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77761	AAMI Level 1 Tab Neck Tie Waist Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77761XL	AAMI Level 1 Tab Neck Tie Waist Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77752	AAMI Level 2 Tie Neck and Waist Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77752XL	AAMI Level 2 Tie Neck and Waist Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
55562	AAMI Level 2 Tab Neck Tie Waist Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
55562XL	AAMI Level 2 Tab Neck Tie Waist Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
55582	AAMI Level 2 Overhead Thumbloop with Side Tie Waist Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
55582XL	AAMI Level 2 Overhead Thumbloop with Side Tie Waist Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
55882	AAMI Level 2 Overhead Thumbloop with Side Tie Waist Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
55882XL	AAMI Level 2 Overhead Thumbloop with Side Tie Waist Blue XL	1/EA, 10 EA/BG, 10 BG/CA
77753	AAMI Level 3 Tie Neck and Waist Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77753XL	AAMI Level 3 Tie Neck and Waist Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
55583	AAMI Level 3 Overhead Thumbloop with Side Tie Waist Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
55583XL	AAMI Level 3 Overhead Thumbloop with Side Tie Waist Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
55883	AAMI Level 3 Overhead Thumbloop with Side Tie Waist Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
55883XL	AAMI Level 3 Overhead Thumbloop with Side Tie Waist Blue XL	1/EA, 10 EA/BG, 10 BG/CA

Gown SMS Isolation Procedure AAMI Level 3 Full Back Knit Cuff Hook and Loop Neck Closure Waist

77873	Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77873XL	Blue XL	1/EA, 10 EA/BG, 10 BG/CA

Non-AAMI Isolation Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations requiring a low level of fluid protection
- Made from strong, fluid resistant and breathable spunbond-meltblown-spunbond (SMS) material
- Appropriate for extended wear situations due to high level of comfort
- Full coverage isolation gown
- Elastic cuffs provide a secure fit under gloves
- Neck tie closure for individual size adjustment
- Available in yellow or blue
- Packaged 10 gowns/bag, 10 bags/case
- Universal and XL sizes
- Single use, non-sterile



Item No.	Description	Quantity
----------	-------------	----------

Gown Isolation SMS Full Back Elastic Cuff Tie Neck and Waist

77780	Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77780XL	Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77880	Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77880XL	Blue XL	1/EA, 10 EA/BG, 10 BG/CA

PERSONAL PROTECTIVE GOWNS

Polyethylene Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations requiring a high level of fluid protection
- Made from polyethylene film that has been tested against ASTM F1670 Synthetic Blood Penetration Resistance and ASTM F1671 Viral Penetration Resistance
- Open back design allows the user to remain cooler without sacrificing protection
- Thumbloop style cuff ensures complete wrist coverage under gloves
- Apron style neck for rapid over the head donning ensures gown compliance
- Perforated tear away neck design for clean, easy removal
- Available in blue
- Individually wrapped gowns packed in a convenient dispenser box
- Universal, Large and XL sizes
- Single use, non-sterile



Item No.	Description	Quantity
----------	-------------	----------

77755	Mid Weight Universal	1/EA, 15EA/BX, 5/CA
77755XL	Mid Weight Extra Large	1/EA, 15EA/BX, 5/CA
77765	Light Weight Universal	1/EA, 15EA/BX, 5/CA
77765L	Light Weight Large	1/EA, 20EA/BX, 6/CA
77765XL	Light Weight Extra Large	1/EA, 15EA/BX, 5/CA

Please see our full line of Personal Protection Equipment (PPE) Organizers & Apparel Dispensers on page 29-30 in the Central Supply Section.



PERSONAL PROTECTIVE GOWNS



AAMI Rated Material Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations requiring a low to high level of fluid protection
- Available in a variety of protective materials
- Material tested in accordance with the ANSI/AAMI PB70:2012 Standard
- Open back design allows the wearer to remain cooler without sacrificing protection
- Thumbloop style cuff provides complete wrist coverage under gloves
- Apron style neck for rapid over the head donning ensures gown compliance
- Perforated tear away neck design for clean, easy removal
- Available in yellow or blue
- Packaged 10 gowns/bag, 10 bags/case
- Universal and XL sizes
- Single use, non-sterile

Item No.	Description	Quantity
Gown Fluid Resistant AAMI Level 3 Material Open Back Overhead Thumbloop Cuff Tie Waist		
77757	Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77757XL	Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77857	Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77857XL	Blue XL	1/EA, 10 EA/BG, 10 BG/CA
Gown SMS Protective AAMI Level 2 Material Open Back Overhead Thumbloop Cuff Tie Waist		
77772	Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77772XL	Yellow XL	1/EA, 10 EA/BG, 10 BG/CA
77872	Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77872XL	Blue XL	1/EA, 10 EA/BG, 10 BG/CA

PERSONAL PROTECTIVE GOWNS

Light-Weight Cover Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations when little fluid exposure is expected
- Made from strong, comfortable and breathable spunbond polypropylene (PSB) material
- Open back design allows wearer to remain cooler without sacrificing protection
- Thumbloop style cuff provides complete wrist coverage under gloves
- Apron style neck for rapid over the head donning ensures gown compliance
- Perforated tear away neck design for clean, easy removal
- Available in yellow or blue
- Packaged 10 gowns/bag, 10 bags/case
- Universal and XL sizes
- Single use, non-sterile



Item No.	Description	Quantity
Gown Isolation Cover Mid Weight PSB Open Back Overhead Thumbloop Cuff Tie Waist		
77876	Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77876XL	Blue XL	1/EA, 10 EA/BG, 10 BG/CA
Gown Isolation Cover Light Weight PSB Open Back Overhead Thumbloop Cuff Tie Waist		
77776	Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77776XL	Yellow XL	1/EA, 10 EA/BG, 10 BG/CA

Mid-Weight Cover Gowns

NOT MADE WITH natural rubber latex

- Ideal for situations requiring a low level of fluid protection
- Made from strong, fluid resistant and breathable spunbond-meltblown-spunbond (SMS) material
- Appropriate for extended wear situations due to high level of comfort
- Open back design allows the wearer to remain cooler without sacrificing protection
- Thumbloop cuffs provide a secure fit under gloves
- Apron style neck for rapid over the head donning ensures gown compliance
- Perforated tear away neck design for clean, easy removal
- Available in yellow or blue
- Packaged 10 gowns/bag, 10 bags/case
- Universal and XL sizes
- Single use, non-sterile



Item No.	Description	Quantity
Gown Isolation Non AAMI SMS Open Back Overhead Thumbloop Cuff Tie Waist		
77879	Blue Universal	1/EA, 10 EA/BG, 10 BG/CA
77879XL	Blue XL	1/EA, 10 EA/BG, 10 BG/CA
Gown Isolation Non AAMI Light Weight SMS Open Back Overhead Thumbloop Cuff Tie Waist		
77779	Yellow Universal	1/EA, 10 EA/BG, 10 BG/CA
77779XL	Yellow XL	1/EA, 10 EA/BG, 10 BG/CA

Utility Scissor

NOT MADE WITH
natural rubber latex

- Razor-sharp, super-hardened surgical stainless steel blades with one serrated edge to cut through the toughest materials
- Stainless steel safety bandage tip
- Ergonomic polypropylene handles with one large ring for maximum control and comfort
- Extra large rivet assures smooth cutting action and rugged durability
- Fully autoclaveable to 132 degrees Celsius (270 degrees Fahrenheit)
- Traditional 7-1/4 Inch length for maximum leverage and cutting strength



Item No.	Description	Quantity
IST1001	Scissor Utility Polypropylene Handle 7.25 Inch	1/EA, 10 EA/BX, 12 BX/CA



Cart Covers

NOT MADE WITH
natural rubber latex

- Provide a dust and dirt barrier for cart contents that are stored or transported between procedures or patient use
- Available on rolls or individually folded for easy dispensing
- Puncture resistant
- Moisture resistant
- Constructed of high strength resins

Item No.	Description	Quantity
MC4807	Cover Cart Roll Clear 38 × 26 × 48 Inch 0.90 MM	1/EA, 50 EA/RL
MC4305	Cover Cart Folded Clear 28 × 22 × 60 Inch 1.35 MM	1/EA, 50 EA/CA
CC3848F	Cover Cart Folded Clear 38 × 26 × 48 Inch 1.85 MM	1/EA, 50 EA/CA
CC3848R	Cover Cart Roll Clear 38 × 26 × 48 Inch 1.85 MM	1/EA, 50 EA/RL
CC4548F	Cover Cart Folded Clear 45 × 27 × 48 Inch 1.85 MM	1/EA, 50 EA/CA
CC4548R	Cover Cart Roll Clear 45 × 27 × 48 Inch 1.8 MM	1/EA, 50 EA/RL

Sterilization Pouches

NOT MADE WITH
natural rubber latex

- Triple heat-sealed seams resist blowouts and provide uniform peel when opening
- High quality medical grade Kraft paper resists tearing—avoids fiber tears in opening process
- Blue-tinted, see-through front allows visual inspection of contents
- Heal seal and self seal designs
- Perforated fold allows for accurate folding of adhesive strip on self seal pouches
- Easy-to-read color changing process indicator arrows; compatible with steam, chemical vapor or ETO sterilization
- Easy-open thumb notch on heat seal and self seal pouches
- Easy-open paper lip on self seal pouches
- Tubing pouches are designed for long or odd-shaped instruments



Item No.	Description	Quantity
----------	-------------	----------

Pouch Sterilization Heat-Seal Ethylene Oxide Steam Dual Indicator Blue-Tinted Kraft Paper

HR3100	3 Inch × 100 Feet	1/RL, 8 RL/CA
HR4100	4 Inch × 100 Feet	1/RL, 16 RL/CA
HR6100	6 Inch × 100 Feet	1/RL, 8 RL/CA
HR1000	10 Inch × 100 Feet	1/RL, 4 RL/CA
HS0610	6 × 10 Inch	1/EA, 200 EA/BX, 10 BX/CA
HS0615	6 × 15 Inch	1/EA, 200 EA/BX, 10 BX/CA
HS0813	8 × 13 Inch	1/EA, 200 EA/BX, 5 BX/CA
HS1215	12 × 15 Inch	1/EA, 200 EA/BX, 5 BX/CA

Pouch Sterilization Self-Seal Ethylene Oxide Steam Dual Indicator Blue-Tinted Kraft Paper

SS3509	3.5 × 9 Inch	1/EA, 200 EA/BX, 10 BX/CA
SS5210	5.25 × 10 Inch	1/EA, 200 EA/BX, 10 BX/CA
SS5214	5.25 × 14 Inch	1/EA, 200 EA/BX, 10 BX/CA
SS7513	7.5 × 13 Inch	1/EA, 200 EA/BX, 5 BX/CA
SS1215	12 × 15 Inch	1/EA, 200 EA/BX, 5 BX/CA
SS1218	12 × 18 Inch	1/EA, 200 EA/BX, 5 BX/CA

INSTRUMENTS

Forceps

NOT MADE WITH
natural rubber latex

sterile EO

- Single-Use
- High quality, stainless steel, floor-grade instruments that meet FDA requirements
- Eliminates time and expense of cleaning, reprocessing and sterilizing instruments
- Helps control cross-contamination with single-use packaging
- Instruments are individually packaged in sterile pouch or non-sterile clear poly bag



Item No.	Description	Quantity
----------	-------------	----------

FG1001	Forceps Adson Dressing Straight Single-Patient Use Stainless Steel 4.75 Inch Sterile	1/EA, 25 EA/CA
FG1002	Forceps Adson Dressing Straight Single-Patient Use Stainless Steel 4.75 Inch Non-Sterile	1/EA, 12 EA/CA

INSTRUMENTS

Towel Clamps

NOT MADE WITH natural rubber latex **sterile EO**

- High quality, stainless steel, floor-grade instruments that meet FDA requirements
- Single-Use
- Eliminates time and expense of cleaning, reprocessing and sterilizing instruments
- Helps control cross-contamination with single-use packaging
- Instruments are individually packaged in sterile pouch



Item No.	Description	Quantity
FG14921	Clamp Backhaus Towel Floor Grade Stainless Steel 3.5 Inch Disposable	1/EA, 50 EA/CA
FG14922	Clamp Backhaus Towel Floor Grade Stainless Steel 5.25 Inch Disposable	1/EA, 50 EA/CA

Hemostats

NOT MADE WITH natural rubber latex **sterile EO**

- Single-Use
- High quality, stainless steel, floor-grade instruments that meet FDA requirements
- Eliminates time and expense of cleaning, reprocessing and sterilizing instruments
- Helps control cross-contamination with single-use packaging
- Instruments are individually packaged in sterile pouch or non-sterile clear poly bag



Item No.	Description	Quantity
Hemostat Kelly 5.5 Inch Stainless Steel Disposable		
FG2001	Curved Sterile	1/EA, 25 EA/CA
FG2002	Curved Non-Sterile	1/EA, 12 EA/CA
FG2003	Straight Sterile	1/EA, 25 EA/CA
FG2004	Straight Non-Sterile	1/EA, 12 EA/CA
Hemostat Surgical 5 Inch Stainless Steel Disposable		
FG2005	Mosquito-Curved Sterile	1/EA, 25 EA/CA
FG2007	Mosquito-Straight Sterile	1/EA, 25 EA/CA
Hemostat Halstead 5 Inch Stainless Steel Disposable		
FG2006	Mosquito-Curved Non-Sterile	1/EA, 12 EA/CA
FG2008	Mosquito-Straight Non-Sterile	1/EA, 12 EA/CA
Hemostat Rochester-Pean 6.25 Inch Stainless Steel Disposable		
FG2009	Curved Sterile	1/EA, 25 EA/CA
FG2010	Clean Curved Non-Sterile	1/EA, 12 EA/CA

INSTRUMENTS

Needle Holders

NOT MADE WITH natural rubber latex **sterile EO**

- High quality, stainless steel, floor-grade instruments that meet FDA requirements
- Single-Use
- Eliminates time and expense of cleaning, reprocessing and sterilizing instruments
- Helps control cross-contamination with single-use packaging
- Instruments are individually packaged in sterile pouch or non-sterile clear poly bag



Item No.	Description	Quantity
FG3001	Holder Needle Webster 5 Inch Disposable Sterile	1/EA, 25 EA/CA
FG3002	Holder Needle Webster 5 Inch Disposable Non-Sterile	1/EA, 12 EA/CA

Scissors

NOT MADE WITH natural rubber latex **sterile EO**

- High quality, stainless steel, floor-grade instruments that meet FDA requirements
- Single-Use
- Eliminates time and expense of cleaning, reprocessing and sterilizing instruments
- Helps control cross-contamination with single-use packaging
- Instruments are individually packaged in sterile pouch or non-sterile clear poly bag



Item No.	Description	Quantity
Scissor Iris Straight 4.5 Inch Stainless Steel		
FG4001	Sterile	1/EA, 25 EA/CA
FG4002	Non-Sterile	1/EA, 12 EA/CA
FG4101	Sterile	1/EA, 25 EA/CA
Scissor Lister Stainless Steel		
FG5001	Bandage 5.5 Inch Sterile	1/EA, 25 EA/CA
FG5002	Bandage 5.5 Inch Non-Sterile	1/EA, 12 EA/CA
FG5104	4.5 Inch Non-Sterile	1/EA, 12 EA/CA
Scissor Operating Room Straight 5.5 Inch Stainless Steel		
FG6001	Blunt/Sharp Sterile	1/EA, 25 EA/CA
FG6002	Blunt/Sharp Non-Sterile	1/EA, 12 EA/CA
FG6002BB	Blunt/Blunt Non-Sterile	1/EA, 12 EA/CA
FG6004	Sharp/Sharp Non-Sterile	1/EA, 12 EA/CA
FG6005	Sharp/Sharp Sterile	1/EA, 25 EA/CA

PERSONAL PROTECTION EQUIPMENT

MediChoice Personal Protection Equipment Organizers help protect staff and patients by providing a solution to keep PPE items visible, stocked, and easily accessible. They are available for holding MediChoice gowns, gloves, face masks, and additional supplies.

Versatility

Increase compliance by mounting the Protection System to the wall outside of each patient room and high traffic areas so protective equipment is visible and available to protect staff at all times.

- **Permanent**—Designed to be used for standardization of location of protective equipment for everyday protection by placing units permanently on walls outside of patient rooms.
- **Portable**—The Protection Systems can be hung intermittently on patient doors and are a cost effective adjunct to existing isolation carts.
- **Essential Infection Control Signage**—Clip-on sign holders can be added to units to ensure proper isolation is initiated and maintained.

Convenience & Accessibility

MediChoice Protection Systems provide a versatile option for standardizing the location of PPE (Personal Protective Equipment) supplies in patient care areas.

- **Visible Cue**—Supplies are highly visible to staff and visitors, helping to improve compliance with isolation policies related to PPE donning, thereby helping to decrease the potential for costly HAIs (Healthcare Associated Infections).
- **Easy Access**—The unit has no doors or drawers allowing for easy access and visualization of consumable products.
- **Save time and money**—Protection Systems provide organization of PPE supplies, reducing clutter in high traffic areas by mounting units on walls or doors located near patient rooms, minimizing the potential for delaying the restocking process and blocking entry and exit ways.



PS016-0233

Protection System

Keyholes enable permanent mounting to a wall
(Cherry Fauxwood ABS Plastic shown)



PS016-0212

Protection System

INFECTION CONTROL CABINETS

Personal Protection Equipment (PPE) Organizers



made in the USA NOT MADE WITH natural rubber latex

- Provides convenient point-of-care access to Personal Protection Equipment (PPE) products
- Allows storage of Personal Protection Equipment in areas with limited space
- Eliminates the need for bulky carts and other obstructive storage units normally found in the hallway
- Wall or door mounted (Hangers sold separately)
- PPE Visible; Easy Restocking; Increased Isolation Compliance
- Organizers can be customized to hospital protocol

Item No.	Description	Quantity
GD-011	Caddy Glove Single Protection Organizer Quartz ABS Plastic 10 x 5.25 x 4 Inch	1/EA
MB-450	Fire Door For Protection Organizer Quartz Metal 34 x 17.6875 x 3.0625 Inch	1/EA
RG-006	Quick Reference Cards (How to Wear PPE in the Hospital)	1/EA, 25 EA/PK
PS016-0512	Protection System Powder Coated Aluminum	1/EA
Cabinet For Personal Protection Equipment		
LD-007	3 Glove Box Gown 1 Box Face Mask Quartz ABS Plastic 16.625 x 25.25 x 4.5 Inch	1/EA
MCLD010	3 Glove Box Gown 2 Box Face Mask	1/EA
PS060512	2 Glove Box Holders 1 Box Face Mask 1 Box Face Shields Bulk Gowns Quartz Powder Coated Aluminum 22.5W x 20.5H x 4.5D	1/EA
Holder		
MW-001	Biohazard Bag 1 Gallon Metal 9.5 x 12.875 x 5.25 Inch	1/EA
MW-003	Biohazard Bag 3 Gallon Metal 11.75 x 15.875 x 7.5 Inch	1/EA
MW-005	Biohazard Bag 5 Gallon Metal 13.3125 x 19 x 7.5 Inch	1/EA
Dispenser		
BB-001	Cap Bouffant Metal White 11.25 x 11 x 5 Inch	1/EA
BB-034	Cap Bouffant Metal White 10 x 10 x 5.25 Inch	1/EA
BG001-0212	Bag 2 Roll Quartz ABS Plastic 10.125 x 7.125 x 3.9375 Inch	1/EA
BG002-0111	Bag 2 Roll Polyethylene Terephthalate Glycol Plastic Clear 10.10 x 7.1 x 3.95 Inch	1/EA
BG002-0212	Bag 2 Roll Quartz ABS Plastic 10.10 x 7.1 x 3.95 Inch	1/EA
BP-007	Cap Bouffant Plastic Clear 11 x 11.25 x 5 Inch	1/EA
EP-105	Ear Plug Bulk Single Bin With Hinged Lid Plastic Clear 9.25 x 11.5625 x 6.75 Inch	1/EA
HP-010	Hair Net Plastic Clear 4.4375 x 8 x 3.9375 Inch	1/EA
MP-015	Caddy For Wet Wipe Canister Clip On Plastic Clear 6 x 7 x 6 Inch	1/EA
MP-042	Caddy Clip On Plastic Clear 4.25 x 4 x 4.5 Inch	1/EA
MP-085	Caddy Clip On Plastic Clear 4.8125 x 8.125 x 5.0625 Inch	1/EA
NC001-0512	Bedpan Urinal Steel 12.33 x 12 x 4.63 Inch	1/EA
RP-011	Gown Bulk Hinged Lid With Bottom Opening Plastic Clear White 14.5 x 15 x 5 Inch	1/EA
RP-030	Gown Bulk Plastic Clear 7.75 x 9.75 x 4.5 Inch	1/EA
SP-573	Shoe Cover Plastic Clear 12.375 x 10.375 x 2.4375 Inch	1/EA
BG001-0111	Bag 1 Roll Polyethylene Terephthalate Glycol Plastic Clear 10.1 x 3.6 x 3.95 Inch	1/EA
FM002-0213	Unit Face Combo White ABS Plastic 10.5 x 23.5 x 5 Inch	1/EA
Eyewear Dispenser		
CP-071	For 10 Pair of Safety Glass with Open Top Plastic Clear 9.0625 x 4.75 x 6.6875 Inch	1/EA
CP-072	For 10 Pair of Safety Glass with Hinged Lid Plastic Clear 9.25 x 3.4375 x 6.6875 Inch	1/EA
CP-073	For 10+ Pair of Safety Glass with Hinged Lid Plastic Clear 9.25 x 5.6875 x 9.3125 Inch	1/EA

Item No.	Description	Quantity
Eyewear Cabinet		
CP-075	For 10 Pair of Safety Glass with Clear Locking Door Quartz ABS Plastic White 13.5 x 20 x 4.125 Inch	1/EA
Protective Apparel Bulk Dispenser		
BP-011	Tall Single Bin with Hinged Lid Plastic Clear 12.25 x 16.9375 x 5.625 Inch	1/EA
BP-012	Tall Double Bin with Hinged Lid Plastic Clear 12.25 x 17 x 5.5 Inch	1/EA
BP-021	Single Bin with Hinged Lid Plastic Clear 12.25 x 11 x 5.5 Inch	1/EA
BP-022	Double Bin with Hinged Lid Plastic Clear 12.25 x 11 x 5.5 Inch	1/EA
BP-023	Tall Triple Bin with Hinged Lid Plastic Clear 17.5 x 11.125 x 5.8125 Inch	1/EA
BP-090	Triple Bin with Hinged Lid Plastic Clear 18.25 x 17 x 5.5 Inch	1/EA
Protection System Dispenser		
PS014-0212	Quartz ABS Plastic Gloveless Holds Gowns 1 Box Earloop Face Mask 1 Box Shield Face Mask Wall Mount Door Hanger or Floor Stand Options 23.375 x 15 x 4.75 Inch	1/EA
PS016-0111	Holds Gowns 3 Boxes Gloves 1 Box Earloop Face Mask 1 Box Shield Face Mask Built In Storage Box Wall Mount Plastic 23.38 x 26.63 x 4.75 Inch Clear Non-Sterile Reusable	1/EA
PS016-0212	Quartz ABS Plastic Holds Gowns 3 Boxes Gloves 1 Box Earloop Face Mask 1 Box Shield Face Mask Built-In Storage Box Wall Mount or Door Hanger Option 23.38 x 26.63 x 4.75	1/EA
PS016-0223	Maple Fauxwood	1/EA
PS016-0233	Cherry Fauxwood	1/EA
PS017-0212	Surgical Quartz ABS and Clear PETG Plastic 4 Garment Bins Built-In Storage Wall Mount 25 x 30.19 x 6.75 Inch	1/EA
Protection System Dispenser All Quartz ABS Plastic		
PS020212	Right Side Load 15L x 22H x 4.5D Inch	1/EA
PS030212	Left Side Load 15L x 22H x 4.5D Inch	1/EA
PS015-0212	Double Gown Holds 2 Gown Types 3 Boxes Gloves 2 Boxes Earloop Face Mask Built-In Storage Box Wall Mount or Door Hanger Option 25.125 x 30.75 x 4.75 Inch	1/EA
Infection Prevention Station		
KS122-0529	Gray Aluminum	1/EA
KS123-0529	Accessory Gray Aluminum	1/EA
Personal Protection Equipment Organizer		
ED-097	Hand Sanitizer Tissue Wet Wipe with Sign Holder Quartz ABS Plastic 15 x 11.375 x 5.75 Inch	1/EA
LD-037	2 Glove Box Gown Face Mask Shield Quartz ABS Plastic 22.75 x 20.5 x 4.75 Inch	1/EA
LD-064	3 Glove Box Gown Face Mask Built Inch Storage Caddy Quartz ABS Plastic 16.75 x 25.75 x 4 Inch	1/EA
LD-067	2 Glove Box Laundered Gown 1 Box Face Mask Quartz ABS Plastic 19.6875 x 28.625 x 4.75 Inch	1/EA
LD-070	3 Glove Box Gown Face Mask Shield Quartz ABS Plastic 20.0625 x 25.25 x 4.75 Inch	1/EA
LP-004	3 Glove Box Bulk Apparel 1 Box Face Mask Quartz ABS Plastic 22 x 19.75 x 5 Inch	1/EA
PS012-0512	4 Glove Box Gown 2 Box Face Mask Quartz ABS Plastic	1/EA
PS007-0212	Gown 2 Built In Caddy 1 Box Face Shield Quartz ABS Plastic 20 x 20.9375 x 4.9375 Inch	1/EA
Wall Hanger for Protection Equipment Organizer		
MB-500	Glass J Hook Quartz Metal 22 x 4.625 x .5625 Inch	1/EA
MB-510	Wall Door J Hook Quartz Metal 5.5 x 12 x 2.5 Inch	1/EA
Door Hanger for Protection Equipment Organizer		
MB-600	Quartz Metal 5.5 x 12 x 2.5 Inch	1/EA
MB-601	Quartz Metal 5.5 x 19 x 2.5 Inch	1/EA

Please see our full line of Gowns on page 23-25 in the Apparel-Staff Section.

INFECTION CONTROL CABINETS

Specialty Dispensers and Accessories



- Provides convenient point-of-care access to personal protective products
- Allows storage of personal protective equipment in areas with limited space
- Eliminates the need for bulky carts and other obstructive storage units normally found in the hallway
- Wall or door mounted (Hangers sold separately)
- Cabinets provided in seven different configurations
- PPE Visible; Easy Restocking; Increased Isolation Compliance
- Organizers can be customized to hospital protocol



Item No.	Description	Quantity
CP-076	Cabinet Eyewear for 12 Pair of Safety Glass with Clear Locking Door Quartz ABS Plastic White 20 × 20 × 5 Inch	1/EA
SD-010	Organizer Surgical Apparel Gown Polyethylene Terephthalate Glycol Plastic Clear 12.375 × 30.6875 × 6.5 Inch	1/EA
MCMA205	Pull Up Shelf Powder Coated White Aluminum 11 × 7.5 × 15 Inch	1/EA
14-182	Key for Bowman Dispensers	1/EA
MC14182	KEY-SET of Two Each for Locking Dispensers	1/ST
RC010111	Spray Perineal Wash without Rinse 8.5 Ounce	1/EA
Bio Shield		
LB001-0111	Large Clear PETG 11.5 × 15 × 7.5 Inch	1/EA
LB002-0111	Small Plastic Clear 8 × 9 × 5 Inch	1/EA
Dental Organizer		
MB-085	1 Glove Box 1 Box Face Mask Bronze Metal 9.875 × 10.125 × 4.1875 Inch	1/EA,
MCMB095	22 Gauge Metal Gloves Mask and Bulk Towels White Powder Coated 10L × 17.5H × 4.5D	1/EA
MD-100	3 Glove Box 2 Box Face Mask Towel Bib Quartz ABS Plastic Bronze 17.625 × 23.5 × 5.0625 Inch	1/EA
Holder		
BW100-0212	Hand Sanitizer Bottle Beige 4.375 × 4 × 4.5 Inch	1/EA
MD-015	Bottle Hand Sanitizer ABS Plastic White 4 × 6.5 × 3.25 Inch	1/EA
Reference Card for Transmission Based Precautions		
MCRG006	Horizontal Attaches to Name Tag	1/EA, 25/CS
MCRG008	Vertical Attaches to Name Tag	1/EA, 25/CS
Dispenser		
BWX01-0412	Wipe Disposable Universal Powder Coated Steel 3.5 × 3.8 × 5 Inch	1/EA
CL002-0111	Task Wipe Holds 1 Box Clear PETG	1/EA
CL003-0111	Task Wipe Large 5.1875 × 15.875 × 3.9375 Inch	1/EA
FB-007	Face Mask Cone Style White Steel	1/EA
FM001-0213	Face Mask Surgical Holds up to 6 Boxes ABS Plastic	1/EA
MCQB001	Cup 22 Gauge Metal White Powder Coated	1/EA, 12/CS
Bag Dispenser		
BG008-0111	Flat Pack Biohazard Bag Clear Plastic Holds 2 Boxes	1/EA
RC001-0111	Valve Mask Ambu-bag Dispenser Clear PETG Plastic Quartz Powder Coated Aluminum 21.25W × 18.25H × 6.25D	1/EA
Injector Epinephrine Dispenser		
ED-750	For up to 5 Door Sleeve Quartz ABS Plastic 9 × 11.1875 × 2.5 Inch	1/EA
ED-760	For up to 10 Plastic Divider Quartz ABS Plastic Yellow 9.5 × 11 × 4.625 Inch	1/EA
Towel Dispenser		
TB-002	Towel Dispenser C or Z Fold White Steel	1/EA
TB-003	Towel Dispenser C or Z Fold White Steel	1/EA

INFECTION CONTROL CABINETS



MCRP012




MCRP013

Apparel Dispensers



- Provides convenient point-of-care access to apparel
- Allows storage of apparel in areas with limited space
- Eliminates the need for bulky carts and other obstructive storage units normally found in the hallway
- Wall or door mounted (Hangers sold separately)
- Manufactured with Quartz ABS plastic and Clear PETG plastic
- Apparel Visible; Easy Restocking; Increased Isolation Compliance
- Organizers can be customized to hospital protocol

Item No.	Description	Quantity
Gown Dispenser		
GN001-0212	Bulk with Lid Quartz ABS Plastic	1/EA,
GN100-0212	Universal Boxed Keyholes for Wall Mounting Quartz ABS Plastic	1/EA
BK001-0300	Bulk Single Tall with Lid Stainless Steel	1/EA
BK002-0300	Bulk Single Short with Lid Stainless Steel	1/EA
Shoe Cover Dispenser		
SP-003	Polyethylene Terephthalate Glycol Plastic Clear 10.5 × 10.25 × 6.75 Inch	1/EA
Shoe Cover/Cap Dispenser		
PD100-0212	Universal Boxed Holds a wide range of Box Sizes Keyholes for Wall Mounting Made of Quartz ABS Plastic	1/EA
PD200-0212	Universal Quartz Lightly-textured Finish ABS Plastic 12.50W × 10.13H × 9.88D Keyholes for Wall Mounting	1/EA
PD300-0212	Universal Quartz Lightly-textured Finish ABS Plastic 13.70W × 5.73H × 6.95D Three-way Keyholes for Horizontal or Vertical Wall Mounting	1/EA



INFECTION CONTROL CABINETS

Glove Box Holders

made in the USA NOT MADE WITH natural rubber latex

- Keeps gloves clean and available
- Multiple types and configurations available
- Industrial strength materials are easy to clean and provide years of trouble-free service
- Quick and easy loading
- Horizontal or vertical mounting available, depending on space constraints



Item No.	Description	Quantity
Glove Box Dispenser		
GB-001	Single Metal White 5.5 × 10 × 3.75 Inch	1/EA
GB-002	Double Metal White 10.5 × 10 × 3.75 Inch	1/EA
GB-003	Triple Metal White 15.75 × 10 × 3.75 Inch	1/EA
GB-027	For Dietary Powder-Coated 22 Gauge Metal White 7.75 × 6.25 × 1 Inch	1/EA
GB-059	Quadruple Wall Mountable Metal White 22.3125 × 10 × 3.8125 Inch	1/EA
GB-067	Double End-To-End Metal White 5.625 × 20.0625 × 3.9375 Inch	1/EA
GB-068	Triple End-To-End Metal White 5.5 × 30.0625 × 3.9375 Inch	1/EA
GB-074	Single with Adjustable Spring Metal White 6.5 × 9.9375 × 3.8125 Inch	1/EA
GD-111	Single with Additional Shelf Quartz ABS Plastic 22.5 × 5.5 × 4 Inch	1/EA
GB-144	Powder-Coated Single with Adjustable Spring 22 Gauge Metal White 5.5 × 10 × 3.75 Inch	1/EA
GD-012	Double Plastic White 22.5 × 5.25 × 4 Inch	1/EA
GD-042	Double with Divider Quartz ABS Plastic 11.375 × 9.9375 × 3.75 Inch	1/EA
GL100214	Single Frosted Polycarbonate Plastic	1/EA, 8/CS
GL003-0413	Triple with Dividers White Steel	1/EA
GL011-0613	Single White Coated Wire 5.67W × 8.16H × 3.75D Slot-style Screw Holes For Vertical or Horizontal Wall Mounting	1/EA
GL022-0613	Double White Coated Wire 10.97W × 8.16H × 3.75D Slot-style Screw Holes for Vertical or Horizontal Wall Mounting	1/EA
GL033-0613	Triple White Coated Wire 16.32W × 8.16H × 3.75D Slot-style Screw Holes for Vertical or Horizontal Wall Mounting	1/EA
GL040-0212	Quad Quartz ABS Plastic	1/EA
GL042-0212	Quadruple Horizontal Mountable with Shelf ABS Plastic 20 × 11.5 × 3.5 Inch	1/EA
GL044-0613	Quad White Coated Wire 21.67W × 8.16H × 3.75D Slot-style Screw Holes for Vertical or Horizontal Wall Mounting	1/EA
GL100-0212	Single Opaque Quartz Injection Molded ABS Plastic 5.8125 × 10.125 × 4.25 Inch Vertical or Horizontal Wall Mounting	1/EA
GL020-0212	Double Holds 2 Boxes Quartz ABS Plastic	1/EA
GL010-0223	Single Maple Faux Wood Signature Series	1/EA
GL020-0223	Double Maple Faux Wood Signature Series	1/EA
GL030-0223	Triple Maple Faux Wood Signature Series	1/EA
GL040-0223	Quad Maple Faux Wood Signature Series	1/EA
GL010-0233	Single Cherry Faux Wood Signature Series	1/EA
GL020-0233	Double Cherry Faux Wood Signature Series	1/EA
GL040-0233	Quad Cherry Faux Wood Signature Series	1/EA
GL030-0233	Triple Cherry Faux Wood Signature Series	1/EA
GL100-1214	Single Semi-transparent Injection Molded Polycarbonate Plastic 5.8125 × 10.125 × 4.25 Inch Vertical or Horizontal Wall Mounting	1/EA, 12 EA/CA

Item No.	Description	Quantity
GL300-0212	Triple Box with Dividers Opaque Quartz Injection Molded ABS Plastic 17.125 × 10.125 × 4.25 Inch Vertical or Horizontal Wall Mounting	1/EA
GL300-1214	Triple with Dividers Semi-transparent Injection Molded Polycarbonate Plastic 17.125 × 10.125 × 4.25 Inch Vertical or Horizontal Wall Mounting	1/EA, 6 EA/CA
GP-013	Single Clear Plastic 5.5 × 9.75 × 3.75 Inch	1/EA
GP-014	Double Plastic Clear 10.5 × 10 × 3.75 Inch	1/EA
GP-015	Triple Clear Polyethylene Terephthalate Plastic 15.75 × 10 × 3.75 Inch Vertical or Horizontal Wall Mounting	1/EA
GP-020	Single Plastic Clear 6.5 × 10 × 4 Inch	1/EA
GP-072	Double Narrow Clear PETG	1/EA
GP-073	Triple Narrow Clear PETG	1/EA
GP-025	Poly Food Service Wall Mountable Plastic Clear 8.3125 × 7 × 1.125 Inch	1/EA
GP-042	Double with Divider Plastic Clear 11.3125 × 9.875 × 3.6875 Inch	1/EA
GP-061	Quadruple Plastic Clear 22.25 × 10 × 3.75 Inch	1/EA
GP-106	Double Plastic Clear 5.5 × 20 × 3.75 Inch	1/EA
GP-107	Triple Plastic Clear 5.75 × 30 × 3.75 Inch	1/EA
GP-120	Chemotherapy Bulk Single Plastic Clear 5.75 × 30 × 3.75 Inch	1/EA
GP-310	Single Horizontal Mountable Plastic Clear 11.1875 × 4.875 × 3.6875 Inch	1/EA
GP-320	Double Horizontal Mountable Plastic Clear 10.25 × 10.5 × 4.25 Inch	1/EA
GP-330	Triple Horizontal Mountable Plastic Clear 11 × 14.25 × 3.5 Inch	1/EA
GP-340	Quadruple Horizontal Mountable Plastic Clear 11 × 19 × 3.5 Inch	1/EA
GS-004	Single with Adjustable Spring Stainless Steel 5.5 × 10 × 3.75 Inch	1/EA
GS-005	Double Stainless Steel 10.5 × 10 × 3.75 Inch	1/EA
GS-006	Triple Stainless Steel 15.75 × 10 × 3.75 Inch	1/EA
GS-108	Double End-To-End Horizontal Mountable Stainless Steel 5.625 × 20 × 3.9375 Inch	1/EA
GS-109	Triple 20 Gauge Stainless Steel 5.75 × 30 × 4 Inch	1/EA
GS-122	Double with Divider Stainless Steel 10.5625 × 10 × 3.8125 Inch	1/EA
GS-123	Triple with Divider Stainless Steel 15.875 × 10 × 3.8125 Inch	1/EA
GS-124	Quadruple with Divider Stainless Steel 21.1875 × 10 × 3.8125 Inch	1/EA
GS-060	Quadruple Stainless Steel 22.25 × 10 × 3.75 Inch	1/EA
GS-073	Single with Adjustable Spring Stainless Steel 6.5 × 9.9375 × 3.8125 Inch	1/EA
GW-055	Coated Wire Triple Vertical Mountable 16 × 7.25 × 3.5 Inch	1/EA, 6/CS
PS004-0111	Triple with Plastic Bag Dispenser Holds 3 Boxes Of Gloves 10.1 × 20.4 × 3.75 Inch	1/EA
PS004-0212	Triple with Biohazard Bag Dispenser Quartz ABS Plastic 10 × 20.5 × 4 Inch	1/EA

Please see our full line of Exam Gloves on page 96-98 in the Gloves Section.

RESPIRATORY HYGIENE STATIONS

MediChoice Respiratory Hygiene Stations help facilities protect patients, visitors, and staff by providing a solution to stock items used in preventing the transmission of infections. They are available for holding MediChoice face masks, hand sanitizing gel, and facial tissues.

Versatility

Most units may be mounted to a floor stand, hung on a wall, or set on a counter. Essential infection control signage can be added to units to ensure proper isolation is initiated and maintained.

Convenience & Accessibility

Hygiene Stations can be placed where infection prevention supplies are needed most—all ports of entry or near an elevator, hanging on a wall.

RS001-0212

Accommodates either one bottle of hand sanitizing gel or one container of hand sanitizing wipes.

Optional sign holder (MP-075 shown, or MP-070) sold separately.



MP-070

Sign Holder Vertical

Optional sign holder (MP-070 shown, or MP-075) sold separately.



RS001-0233

Respiratory Hygiene Station

Color Options

RS001-0233

Maple Fauxwood Plastic

RS001-0212

Quartz ABS Plastic

KS010-0420

Floor Stand


Color Options

KS010-0434

White Cold Rolled Steel

KS010-0412

Quartz Cold Rolled Steel



Please see our full line of Face Masks on pages 22-23 in the Apparel–Staff Section, and Facial Tissues on page 13 in the Admissions Section.

INFECTION CONTROL CABINETS

Respiratory Hygiene Stations

made in the USA NOT MADE WITH natural rubber latex



- Provides protection in public areas
- Available in multiple configurations
- Wall Mount, Countertop Display, Kiosk Stand
- Locking or non-locking

Item No.	Description	Quantity
RG-008	Transmission Based Precautions Quick Reference Guide	1/EA, 25 EA/PK
Kiosk Stand		
KS010-0412	Powder-Coated Holds Respiratory Hygiene Station 16.2 Pound Steel 18 x 36 x 18 Inch	1/EA
KS101-0029	Hand Sanitizer Plate with Stand Split Post and Base Metal	1/EA
Floor Stand		
KS010-0420	Black Steel	1/EA
KS010-0434	White Steel	1/EA
KS102-0029	Hand Sanitizer	1/EA
KS201-0029	Mounting Plate Gray Steel	1/EA
Face Mask Dispenser		
FB-037	Metal White 7.5625 x 4.25 x 4.125 Inch	1/EA
FB-040	Metal White 8 x 6 x 5.25 Inch	1/EA
FB-041	Bulk Metal White 8.5 x 5 x 3.25 Inch	1/EA
FB-090	Metal White 9.0625 x 7.5 x 4.0625 Inch	1/EA
FB-091	Metal White 7.875 x 7.6875 x 5 Inch	1/EA
FD-005	Bulk with Hinged Lid Plastic White Clear 8.6875 x 5 x 4.375 Inch	1/EA
FD-010	Bulk with Hinged Lid Plastic White Clear 8.6875 x 9.75 x 4.375 Inch	1/EA
FP-017	Cone Bulk with Hinged Lid Polyethylene Terephthalate Glycol Plastic Clear 6.6875 x 5 x 5.6875 Inch	1/EA
FP-038	Cone Polyethylene Terephthalate Glycol Plastic Clear 7.6875 x 4.3125 x 4.25 Inch	1/EA
FP-084	Cone Polyethylene Terephthalate Glycol Plastic Clear 5.1875 x 6.125 x 5.3125 Inch	1/EA
Sign Holder		
MP-047	Clip-On Plastic Clear 8.5 x 11.1875 x 0.625 Inch	1/EA
MP-070	Attaches to Respiratory Hygiene Station 8.5 x 11 Inch Vertically Plastic Clear 8.625 x 13.8125 x 0.3125 Inch	1/EA
MP-075	Attaches to Respiratory Hygiene Station 11 x 8.5 Inch Horizontally Plastic Clear 11.25 x 11.3125 x 0.3125 Inch	1/EA
SN003-0212	Mounts Vertical/Horizontal 5.5 x 8.5 Inch Quartz ABS Plastic	1/EA
Respiratory Hygiene Station		
FD-062	Locking Lid Mask Tissue Hand Sanitizer Quartz ABS Plastic 17.5 x 10 x 5 Inch	1/EA
FD-068	Organizer Infection Prevention 2 Glove Box Face Mask Hand Sanitizer Quartz ABS Plastic 18.5 x 11.625 x 4.5 Inch	1/EA
FD-111	Non-Locking Lid Mask Tissue Hand Sanitizer Quartz ABS Plastic 19.75 x 13.875 x 5 Inch	1/EA
FD-112	Non-Locking Lid Mask Tissue Hand Sanitizer Quartz ABS Plastic 19.75 x 13.875 x 5.4375 Inch	1/EA
RS001-0212	2 Box Face Mask Tissue Wipe Hand Sanitizer Quartz ABS Plastic 17.875 x 10.5 x 4.4375 Inch	1/EA
RS001-0223	Maple Faux Wood Signature Series	1/EA
RS001-0233	Cherry Faux Wood Signature Series	1/EA
RS001-0512	Holds 2 Boxes Face Masks 1 or 2 Boxes Facial Tissues and 1 Hand Sanitizer Bottle or Hand Sanitizing Wipe Canister Quartz Powder-Coated Aluminum 17.9 x 11.8 x 5.1 Inches	1/EA
RS002-0212	Station Respiratory Hygiene Hinged Locking Lid	1/EA
RS004-0512	Non-Locking Holds 1 Hand Sanitizer Dispenser 2 Boxes Face Masks 1 or 2 Boxes Facial Tissues or 1 Box Gloves Quartz Powder-Coated Aluminum 22.3 x 14.9 x 5.2 Inches	1/EA
RS005-0512	Locking Holds 1 Hand Sanitizer Dispenser 2 Boxes Face Masks 1 or 2 Boxes Facial Tissues or 1 Box Gloves Quartz Powder-Coated Aluminum 17.9 x 11.8 x 5.1 Inches	1/EA

INFECTION CONTROL CABINETS

Face Mask Dispensers

made in the USA NOT MADE WITH natural rubber latex



- Provides convenient point-of-care access to face masks
- Allows storage of face masks in areas with limited space
- Eliminates the need for bulky carts and other obstructive storage units normally found in the hallway
- Wall or door mounted (Hangers sold separately)
- Face Masks Visible; Easy Restocking; Increased Isolation Compliance
- Organizers can be customized to hospital protocol

Item No.	Description	Quantity
MP-035	Sign Holder for FM001-0214 Clear PETG Plastic	1/EA
Powder-Coated 22 Gauge Metal White		
FB-001	5.6875 x 6.4375 x 5.0625 Inch	1/EA
FB-002	6.4375 x 8.5 x 5.1875 Inch	1/EA
FB-003	8.0625 x 9 x 5.1875 Inch	1/EA
FB-004	5.3125 x 7.75 x 5.3125 Inch	1/EA
FB-005	9.5625 x 7.1875 x 3.5625 Inch	1/EA
FB-010	10.4375 x 4.125 x 4.8125 Inch	1/EA
FB-021	7.4275 x 4.75 x 4.8125 Inch	1/EA
FB-039	7.8125 x 7.5 x 5.25 Inch	1/EA
FB-063	7.8125 x 7.5 x 5.25 Inch	1/EA
FB-064	7.5625 x 4.25 x 4.125 Inch	1/EA
Polyethylene Terephthalate Glycol Plastic Clear		
FP-009	Keyhole for Wall-Mount 5.5 x 7.75 x 5.5 Inch	1/EA
FP-010	9.6875 x 7.1875 x 3.6875 Inch	1/EA
FP-015	10.5625 x 4.125 x 4.9375 Inch	1/EA
FP-021	Keyhole for Wall-Mount 7.9375 x 4.25 x 5.1875 Inch	1/EA
FP-022	Keyhole for Wall-Mount 8.125 x 7.75 x 4.9375 Inch	1/EA
FP-028	7.6875 x 5.625 x 5.1875 Inch	1/EA
FP-119	Keyhole for Wall-Mount 4.6875 x 6.5 x 4.0625 Inch	1/EA
FP-120	Keyhole for Wall-Mount 9.1875 x 7.5 x 4.1875 Inch	1/EA
FP-121	Keyhole for Wall-Mount 7.6875 x 4.375 x 4.125 Inch	1/EA
FP-122	Keyhole for Wall-Mount 8.0625 x 7.625 x 5.375 Inch	1/EA
Universal		
FM100-0212	Earloop Style 8.25 x 4.625 x 5.75 Inch Keyholes for Wall Mounting Quartz ABS Plastic	1/EA
FM200-0212	Shield Style 9.75 x 8.0625 x 4.75 Inch Keyholes for Wall Mounting Quartz ABS Plastic	1/EA
FM300-0212	Cone Style 9.75 x 8.0625 x 4.75 Inch Keyholes for Wall Mounting Quartz ABS Plastic	1/EA

Please see our full line of Hand Sanitizer on page 158 in the Patient Care Section.

INSTRUMENT CARE & DECONTAMINATION

Enzyme Cleaner Pre-Soak



- Neutral pH, non-corrosive
- Three enzymes work together to remove soils quickly:
 - Protease breaks down proteins
 - Lipase breaks down fat
 - Amylase breaks down starch and carbohydrates
- Disintegrates proteinaceous matter in hard to reach crevices
- Significantly reduces the need for manual cleaning
- Highly concentrated: Super Concentrate .1 ounce per gallon; Premium .5 ounce per gallon; Standard 1 ounce per gallon
- Controlled foaming
- Pleasant fragrance



Item No.	Description	Quantity
Pre Soak Non-Corrosive		
ENZ001	1 Gallon	1/EA, 4 EA/CA
ENZ005	5 Gallon	1/EA
ENZ015	15 Gallon	1/EA
ENZ501	Fragrance-Free 1 Gallon	1/EA, 4 EA/CA
Premium Pre Soak Non-Corrosive		
PEZ001	1 Gallon	1/EA, 4 EA/CA
PEZ005	5 Gallon	1/EA
PEZ015	15 Gallon	1/EA
Super Concentrate Pre Soak Non-Corrosive		
PEZ101	1 Gallon	1/EA, 4 EA/CA
PEZ110	2.5 Gallon	1/EA, 2EA/CA
PEZ301	High-Power 1 Gallon	1/EA, 4 EA/CA
MCPUN1G	Pump Enzymatic Cleaner Manual Non-Corrosive for Gallon Container	1/EA

Enzyme Spray Pre-Cleaner



- Initiates enzymatic pre-cleaning
- Reduces heavy instrument transport weight compared to liquid-filled trays
- Free-rinsing and compatible with all washers
- Minimizes potential spill exposures due to splashing



Item No.	Description	Quantity
Spray Pre-Cleaner Instrument Enzyme		
EFS201	Foam 2 Ounce	1/EA, 48 EA/CA
EFS200	Foam 32 Ounce	1/EA, 6 EA/CA
EGS100	Gel 32 Ounce	1/EA, 6 EA/CA

INSTRUMENT CARE & DECONTAMINATION

Specialty Products



- Rust & Stain Remover
 - Concentrate removes rust, stains, and biofilm from metal surfaces
 - Cost effective way to improve surgical instrument performance and appearance
- Liquid Delimer-Descaler Concentrate
 - Removes lime and other hardwater deposits from washers/decontaminators and sterilizers
 - Extends life of washers and sterilizers
- Adhesive Remover
 - Ready-to-use spray removes stickers/labels, tape residue, adhesives, oil, grease, candle wax, chewing gum, ink/markers, glue etc.
 - Safe for use on a variety of metals, glass, plastic and rubber



Item No.	Description	Quantity
MSP001	Delimer Descaler Concentrate 1 Gallon	1/EA, 4 EA/CA
ISP001	Rust & Stain Remover Concentrate 1 Gallon	1/EA, 4 EA/CA
ISP080	Adhesive Remover Ready to Use Spray 8 Ounce	1/EA, 12 EA/CA

Instrument Lubricants



- Neutral pH
- Contains anti-microbial ingredients
- Extends the life of the instruments and ensures proper operation of hinges and joints
- Easier cleaning: no rinsing or wiping
- Reduces rusting, discoloration, and corrosion
- Water soluble and steam permeable
- Super concentrated formula dilution ratio is 1/10th ounce per gallon of water



Item No.	Description	Quantity
Lubricant Instrument		
INL001	Concentrate 1 Gallon	1/EA, 4 EA/CA
INL005	Concentrate 5 Gallon	1/CS
INL015	Concentrate 15 Gallon	1/DR
INL101	Super Concentrate 1 Gallon	1/EA, 4 EA/CA
INL110	Super Concentrate 2.5 Gallon	1/EA, 2 EA/CA

INSTRUMENT CARE & DECONTAMINATION



In our MediChoice® line, Owens & Minor offers an extensive selection of items to serve Central Supply departments in hospitals and surgery centers. From cart covers to cleaners to sterile processed, patient-ready instruments, consider MediChoice for a quality and cost-effective alternative for the storage and protection you require. The MediChoice central supply product offering includes:

- Cart Covers
- Enzyme Cleaners
- Forceps
- Hemostats
- Instrument Lubricants
- Needle Holders
- Scissors
- Specimen Bags
- Sterilization Pouches

INSTRUMENT CARE & DECONTAMINATION

Neutral pH Instrument Cleaners



- Neutral pH
- Effective and safe way to maximize instrument cleanliness and service life
- Safe to use on brass, copper and aluminum
- Contains no phosphates
- Non-corrosive
- Rinse aid additive
- Super concentrated formula dilution ratio is 1/10th ounce per gallon of water



Item No.	Description	Quantity
Detergent Instrument Cleaning Neutral Ph		
NPH001	1 Gallon	1/EA, 4 EA/CA
NPH005	5 Gallon	1/EA
NPH015	15 Gallon	1/EA
NPH101	Super Concentrate 1 Gallon	1/EA, 4 EA/CA
NPH110	Super Concentrate 2.5 Gallon	1/EA, 2 EA/CA

INSTRUMENT CARE & DECONTAMINATION

Cart Cleaners

- Specially formulated for all makes of cart washing units
- Will not harm copper, stainless steel, tungsten carbide, brass or alloys
- Contains blood cutting agents and rust inhibitors
- Special surfactant blend prevents spotting
- Super concentrated to save money



Item No.	Description	Quantity
Detergent Cart Cleaning		
MCW005	Concentrate 5 Gallon Drum	1/EA
MCW015	15 Gallon Drum	1/EA

Rinse Aids



- Special surfactants provide supplementary cleaning
- Specially formulated for all makes of cart washing equipment
- Helps rinse away spots and films
- Creates a powerful "sheeting" action to reduce drying time
- Highly concentrated



Item No.	Description	Quantity
Aid Rinse Cart Washing Highly Concentrated Dry Shine Drum		
MRA005	5 Gallon	1/EA
MRA015	15 Gallon	1/EA

EQUIPMENT



Whether you are planning a new facility, an expansion of an existing one, or simply replacing equipment items over time, consider the more than 500 products on the following pages. When you include MediChoice equipment in your plans, you can trust that our development staff has selected products to ensure we deliver on the high quality, outstanding service, and superior value that you have come to expect from Owens & Minor.

SCALES AND MEASURING DEVICES

Neo-Natal and Intensive Care Infant Scales

NOT MADE WITH natural rubber latex

- Measuring rod included
- Integrated digital display
- Shielded from electronic interference
- Power supply: 6 AA rechargeable batteries and A/C adapter/recharger included
- Capacity: 44 pounds/20 kilograms in 0.2 ounce/5 gram increments
- Unique tray design
- Protected against variation in temperature and humidity
- Weight: 11 pounds
- Dimensions: 6.5 Inch height × 22 Inch width × 12.5 Inch depth
- 3-year manufacturer's warranty




Item No.	Description	Quantity
SCAL1603	Scale Infant with Fine Graduation and Measuring Rod 44 Pound Capacity 6.5 Inch × 22 Inch × 12.5 Inch	1/EA

Measuring Rods

NOT MADE WITH natural rubber latex

- Telescopes to allow shorter practitioners the ability to measure height of very tall patients
- Headpiece can be folded down for safety
- Floor mount eliminates need for calibration
- Drywall anchors and 16 spacers for easy installation
- Anodized aluminum to withstand heavy use in clinical setting
- Scale mounts on wall and directly against floor—no calibration necessary
- Foot plate for increased stability
- 4 Inch medical grade plastic head plate to ensure accuracy
- Broad measuring slide and heel positioner
- Mounts easily with 2 screws
- Measuring range: 2 Inch to 90 Inch/6 centimeters to 230 centimeters in .125 Inch/1 millimeter increments
- Dimensions: 52 Inch height × 4.5 Inch width × 1.5 Inch depth
- 3-year manufacturer's warranty




Item No.	Description	Quantity
MEAS1608	Rod Measuring Mechanical Telescopic with Large Measuring Range Wall-Mount with Head Foot Plate Aluminum 52 Inch × 4.5 Inch × 1.5 Inch	1/EA

SCALES AND MEASURING DEVICES

Baby Measuring Mats

NOT MADE WITH natural rubber latex

- Heavy-duty “soft touch” material
- Easy to transport—can be rolled up and stored in a drawer unlike rigid boards
- Increments in Inches and centimeters
- Measuring range: 4 Inches to 39 Inches or 10 centimeters to 99 centimeters
- Washable
- Can be hung for use in labor and delivery
- Dimensions: 53 Inch length × 12 Inch length × 5.5 Inch depth
- 3-year manufacturer's warranty



Item No.	Description	Quantity
MEAS1609	Mat Mobile Measuring Baby and Toddler Washable Heavy Duty Soft Touch Material 53 Inch x 12 inch x 5.5 Inch	1/EA

Scales with BMI, Digital Height Rod and Wireless Transmission

NOT MADE WITH natural rubber latex

- Seca 360 wireless technology directly to EMR patient files
- Capacity: 660 pounds/300 kilograms in 0.1 pound/50 gram Increments
- Measuring range: 11 Inch - 86 Inch height/30-200 centimeters in 0.05 Inch/1 millimeter increments
- Planar load cell technology
- Functions: TARE, Pre-TARE, HOLD, auto-HOLD, auto-BMI, and pound/kilogram/switch over
- Immediate printout of measurements compared to statistical norms
- Power supply: A/C UL adapter, batteries
- Low profile glass platform with casters: 17 Inch width × 2.5 in height × 14.5 Inch length
- Dimensions: 17 Inch width × 94 Inch height × 18 Inch depth
- 3-year manufacturer's warranty



Item No.	Description	Quantity
SCAL1284	Scale Digital Measuring Station with Wireless Transmission of Height and Weight 660 Pound Capacity	1/EA

Please see our full line of Baby Scale Liners on page 129 in the Maternity and Baby Care Section.



SCALES AND MEASURING DEVICES

Infant and Baby Scales

NOT MADE WITH
natural rubber latex




- Shell-shaped tray designed to cradle patient and provide safety
- Switches from pounds to kilograms with touch of a button next to the display, but can be configured for metric-only use
- Wireless EMR Connectivity
- Dimensions: 24.5 Inch length × 7.5 Inch height × 14 Inch width
- 3-year manufacturer's warranty
- AUTOHOLD feature allows nursing staff to care for patient first then return to the scale to note the weight
- Immediate printout of measurements compared to statistical norms
- TARE function allows tray liner to be ignored when weighing patient
- Power Supply: 6 AA batteries included, optional A/C adapter available
- Capacity: 44 pound /20 kilograms in 0.2 ounce/5 gram increments
- Damping function ensures precise weighing even with patient movement
- Easy-to-read display monitor
- Sanitary, medical grade plastic construction

Item No.	Description	Quantity
SCAL1374	Scale Baby Digital with Wireless Transmission Cradle Design Plastic 44 Pound Capacity 24.5 Inch × 7.5 Inch × 14 Inch	1/EA

Mechanical Physician Scales

NOT MADE WITH
natural rubber latex



- Large, slip-resistant, heavy-duty, cast iron base provides for super solid weighing mechanism
- Telescoping height rod and wheels included
- Easy assembly
- Programmable to pounds or kilograms
- Capacity: 500 pounds in 0.125 pound increments (225 kilogram in 100 gram increments)
- Heavy-duty lever system
- Weight: 52 pounds
- Dimensions: 61.375 Inch height × 20.5 Inch width × 20.5 Inch depth
- 6-year manufacturer's warranty

Item No.	Description	Quantity
SCAL1601	Scale Physician Mechanical Column with Eye-Level Beam and Height Rod 500 Pound Capacity 61.375 Inch × 20.5 Inch × 20.5 Inch Pound Model	1/EA

SCALES AND MEASURING DEVICES

Combination Bariatric Scales

NOT MADE WITH
natural rubber latex



- Four-function combination scale with BMI can be used as a handrail, wheelchair, chair, or platform scale ("The scale that does everything")
- Flip-seat mountable on right or left with integrated damping system
- Weight: 120 pounds
- RS232 compatible
- Swivel display
- Power Supply: A/C UL adapter
- Capacity: 800 pounds/360 Kilograms in 0.2 pound/100 gram increments
- TARE, pre-TARE, and HOLD functions
- Slip-resistant mat included
- Easily moved with integrated wheels
- Base Dimensions: 2 Inch height × 29.9 Inch width × 32.3 Inch depth
- Dimensions: 43.75 Inch height × 35.75 Inch width × 41.9 Inch depth
- 3-year manufacturer's warranty

Item No.	Description	Quantity
SCAL1602	Scale Digital Multifunctional Bariatric 800 Pound Capacity 43.75 Inch × 35.75 Inch × 41.9 Inch	1/EA

Please see our full line of Wheelchairs on page 137-138 in the Mobility Aids Section.

SCALES AND MEASURING DEVICES

Wheelchair Scales

NOT MADE WITH natural rubber latex

- Wide base will hold most bariatric wheelchairs
- Immediate printout of measurements compared to statistical norms
- Integrated pre-TARE function determines patient's net weight in seconds
- Folds up in just a few easy steps, with transport casters for easy mobility
- Power Supply: Internal rechargeable batteries and A/C UL adapter/recharger included
- Capacity: 800 pounds/360 kilograms in 0.2 pound/100 gram increments
- 3-year manufacturer's warranty
- SCAL1604:
 - Base Dimensions: 2 Inch height x 31.375 Inch width x 38 Inch depth
 - Side barriers ensure maximum safety
 - Rechargeable batteries
 - Overall Dimensions: 36 Inch height x 43.5 Inch width x 45.5 Inch depth
 - Weight: 71 pounds
- SCAL1676:
 - Wireless technology directly to EMR patient files
 - Graduation: 0.1 pounds | 50 g
 - Dimensions: 36.2 Inch W x 44.1 Inch H x 45.3 Inch D
 - Base Dimensions: 29.9 Inch W x 2.2 Inch H x 36.6 Inch D
- RS232 compatible



Item No.	Description	Quantity
Scale Wheelchair Digital 800 Pound Capacity Wireless Transmission		
SCAL1604	With Slip-Resistant Mat 36 Inch x 43.5 Inch x 45.5 Inch	1/EA
SCAL1676	With Handrail	1/EA

High Capacity Digital Column Scales with BMI and Digital Height Rod

NOT MADE WITH natural rubber latex

- Integrated Body Mass Index function calculates BMI at the touch of a button
- Digital height rod with LCD display
- Wheels included
- Power Supply: A/C UL adapter
- Capacity: 550 pounds/250 kilograms in 0.1 pound/50 gram increments
- Functions: TARE, HOLD, auto-BMI and pounds/kilograms switch over
- EMR compatible: RS 232 Adapter (seca 460)
- Dimensions: 14 Inch width x 52 Inch height x 27 Inch length
- Measuring range: 43 Inch to 79 Inch height/110 centimeters to 200 centimeters in 0.05 Inch/1 millimeter increments
- Large heavy duty base: 14 Inch width x 3½ Inch height x 14 Inch length
- 3-year manufacturer's warranty



Item No.	Description	Quantity
SCAL1763	Scale Digital High Capacity with Body Mass Index and Measuring Station 550 Pound Capacity 14 Inch x 52 Inch x 27 Inch	1/EA

SCALES AND MEASURING DEVICES

Digital Chair Scales

NOT MADE WITH natural rubber latex

- Seca 360 wireless technology directly to EMR patient files
- Capacity: 600 pounds/300 Kilograms in 0.2 pound/100 gram increments
- Easy weighing of non-ambulatory patients
- Locking wheels make chair maneuverable
- Tare/hold feature
- Swivel armrests allow patients to get in and out easily, accommodates larger patients
- 4-wheel design prevents nurses straining to move scale
- Sanitary, medical grade plastic seat
- Foot rests extend to weigh patients with limited knee flexibility
- Power Supply: Internal rechargeable batteries and A/C adapter/recharger included
- Made of tubular steel
- May be used to move patients in limited cases
- Dimensions: 36 Inch height x 22 Inch width x 38.5 Inch depth
- 3-year manufacturer's warranty



Item No.	Description	Quantity
SCAL1607	Scale Chair Digital with Wireless Transmission 600 Pound Capacity 36 Inch x 22 Inch x 38.5 Inch	1/EA

High Capacity Digital Column Scales with BMI

NOT MADE WITH natural rubber latex

- Robust and durable cast-iron construction
- Integrated Body Mass Index (BMI), Hold and Tare functions
- Altitude calibration
- Weight: 35 pounds
- Wheels and height rod included
- Power Supply: 6 AA batteries included that last for up to 16,000 uses, optional A/C adapter available
- Capacity: 660 pounds/300 kilograms in 0.2 pound/100 gram increments
- Immediate printout of measurements compared to statistical norms
- Programmable to pounds or kilograms
- Shielded from electronic interference and humidity
- Wireless EMR Compatibility
- Cast iron base platform dimensions: 3.5 Inch height x 14 Inch width x 14 Inch depth
- Dimensions: 36.5 Inch height x 14 Inch width x 20.5 Inch depth
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Scale Digital High Capacity with Body Mass Index and Wireless Transmission 650 Pound Capacity 36.5 Inch x 14 Inch x 20.5 Inch		
SCAL1610	Pound Model	1/EA

SCALES AND MEASURING DEVICES

Organ and Diaper Scales

NOT MADE WITH natural rubber latex



- Stainless steel surface with black plastic base is easy to clean to meet clinical hygiene demands.
- Electronics and battery compartment are protected to prevent damage from leaking fluids
- Graduations: 0.05 ounce < 600 ounce > 0.1 ounce
- Capacity: 11 pounds/5 kilograms
- Battery operated for use anywhere
- Automatic switch-off function guarantees longer battery life
- Dimensions: 10.5 Inch width × 2.2 Inch height × 10.5 Inch depth, 4.2 pounds overall weight
- 3 Year Warranty

Item No.	Description	Quantity
SCAL1856	Scale Organ Diaper Digital with Cover Stainless Steel 11 Pound Capacity 10.5 Inch × 2.2 Inch × 10.5 Inch	1/EA

Digital Floor Scales with Cable Remote Control

NOT MADE WITH natural rubber latex



- Remote display can be hand-held or conveniently mounted to wall at eye level
- Built-in Body-Mass-Index function, no separate calculator needed
- HOLD function to make weighing and charting easier
- TARE function, perfect for weighing a toddler while being held by mother
- Power Supply: 4 AA batteries included that will last thousands of uses, optional A/C adapter available
- Capacity: 550 pounds/250 kilograms in 0.2 pounds/100 gram increments
- Lightweight, easy to transport
- Auto switch-off conserves batteries
- Optional padded transport case with shoulder strap available
- 3-year manufacturer's warranty

Item No.	Description	Quantity
SCAL1882	Scale Floor Portable with Digital Cable Remote Display with Built In Body Mass Index 550 Pound Capacity	1/EA

SCALES AND MEASURING DEVICES

Digital BMI Scales with Handrail

made in the USA

NOT MADE WITH natural rubber latex

- MotionLock Technology automatically calculates accurate weight even with unsteady patients
- Ergonomically designed handrail for accommodating smallest to tallest patient reads correctly even when patient applies weight to the handrails
- Extremely durable, all metal construction for long lasting performance
- Tilt and roll design for easy transport
- Arrives fully assembled—just put in leveling feet and batteries
- SCAL81DH features integrated digital height rod that quickly and accurately measures height with no buttons to push
- SCAL31MC offers optional manual height rod
- 3-year manufacturer's warranty
- Large platform with heavy-duty built-in 3 Inch casters for easy movement
- Console displays weight, height, and BMI on a single screen
- Features include BMI, auto zero, weight recall, programmable auto-off from 12 seconds to 25 minutes, tare, and more
- Pound/Kilogram Lockout: quickly change setup to disable LB/KG button so scale will only weigh in selected unit of measure
- Body Mass Index feature with programmable clothing allowance
- EMR connectivity via serial output or optional wireless module
- Powered by 4 Dcell alkaline batteries (included) that last 100,000 weighings
- Capacity: 1000 pounds/450 kilograms in 0.1 pound/kilogram increments



SCAL81DH



SCAL31MC

Item No.	Description	Quantity
Scale Digital 1000 Pound Capacity 22 Inch × 25 Inch		
SCAL31MC	With Body Mass Index Handrail	1/EA
SCAL81DH	With Body Mass Index Handrail and Integrated Digital Height Rod	1/EA
Adapter		
ADPTRMC	Adapter Power Medical Grade Fits Scales SCAL17MC SCAL20MC SCAL31MC SCAL38CH SCAL38MC SCAL41MC SCAL42MC SCAL45MC SCAL46MC SCAL49MC SCAL81DH SCAL88DH	1/EA

SCALES AND MEASURING DEVICES

**Neo-Natal Scales with Handles
1 Gram Accuracy**

made in the USA NOT MADE WITH natural rubber latex



- Precise 1 gram accuracy with large 45 pound capacity for all infant care
- MotionLock Technology accurately calculates weight of active infants
- Displays in all popular unit of measure formats: grams, kilograms, and pounds/ounce
- Ability to lock-in the preferred display format (i.e. kilograms) insuring scale always displays in chosen format
- All metal construction with replaceable 25 Inch width x 14 Inch depth acrylic cradle with length measure
- Built-in heavy-duty handles for safe transport
- Power Supply: 6 C-cell alkaline batteries (included) that last over 80,000 weighings
- Capacity: 45 pounds/20 kilograms in 1 gram (0 kilogram - 10 kilogram), 5 gram (10 kilogram - 20 kilogram), 0.1 ounce (0 pound - 45 pound) increments over full range
- 3-year manufacturer's warranty

Item No.	Description	Quantity
SCAL20MC	Scale Neonatal 1 Gram Accuracy 45 Pound Capacity 25 Inch x 14 Inch	1/EA
CRT202	Cart Scale Fits MediChoice Neonatal Scale SCAL20MC	1/EA

SCALES AND MEASURING DEVICES

Portable Bariatric Scales with BMI

made in the USA NOT MADE WITH natural rubber latex



- Generous 20 Inch Width x 24 Inch depth platform with built-in handle for easy transport
- MotionLock Technology for accurate weight even with unsteady patients
- 3-year manufacturer's warranty
- Console displays weight, height, and BMI on a single screen
- Features include BMI calculator, auto zero, weight recall, programmable "auto-off" from 12 seconds to 25 minutes, tare, and more
- Pound/Kilogram Lockout: quickly change setup to disable LB/KG button so scale will only read in selected unit of measure
- EMR connectivity via serial output or optional wireless module
- Powered by 4 D-cell alkaline batteries (included) that last 100,000 weighings
- Capacity: 1000 pounds in 0.1 pound increments
- Easy transport with lightweight aluminum construction and built-in handle
- Console mounting brackets (attached) allow positioning on wall or angled on tabletop

Item No.	Description	Quantity
SCAL17MC	Scale Portable Bariatric with Body Mass Index Built-in Handle 1000 Pound Capacity Aluminum 20 Inch x 24 Inch	1/EA

Digital Scales with Handrail

NOT MADE WITH natural rubber latex



- Wireless EMR Compatibility
- Capacity: 800 pounds/360 kilograms
- Graduation: 0.1 pounds/50 grams
- Functions: TARE, pre-TARE, HOLD, BMI, auto-CLEAR, kg/lbs switch-over, damping, automatic switch-off, SEND/PRINT
- Provides an immediate printout of measurements
- Optional: seca 233 measuring rod
- RS232 compatible

Item No.	Description	Quantity
SCAL1644	Scale Digital Multifunctional Handrail with Wireless Transmission 800 Pound Capacity	1/EA

Digital Portable Scales

NOT MADE WITH natural rubber latex



- Capacity: 440 pounds /200 kilograms
- Graduation: 0.2 pounds /100 grams
- Functions: auto-HOLD, step-off, pounds/kg switch-over, tap-on automatic switch-on, automatic switch-off
- Recommended for moderate use

Item No.	Description	Quantity
SCAL1813	Scale Digital Flat Floor for Individual Patient Use	1/EA

SCALES AND MEASURING DEVICES

Convertible Chair Scales

made in the USA NOT MADE WITH natural rubber latex

- Two great scales in one—a portable handrail scale, as well as a wide-body chair scale which easily converts between modes in seconds, with no tools required
- Sturdy wrap-around handrail provides maximum security, as well as a full 25 Inch span between handgrips
- Powered by 4 D-cell alkaline batteries (included) that last 100,000 weighings
- Slide-out footrest for patient comfort
- Durable, all metal construction for years of dependable use
- Lightweight with large 3 Inch wheels that tilt and roll effortlessly on all flooring surfaces for maximum portability
- MotionLock Technology automatically calculates average weight reading for patients that are unable to remain still or are unsteady
- Console displays weight, height, and BMI on a single screen
- Pound/Kilogram Lockout: quickly change setup to disable Pound/Kilogram button so scale will only read in selected unit of measure
- Features include tare, BMI calculator, auto-zero, weight recall, autoweigh, an always-on option available with no buttons to press for normal weighing
- Quick two second weight measurement accurate to within 0.1 pounds (0.1 kilograms)
- Has 750 pound capacity in Handrail Scale mode, or 500 pound capacity in Chair mode
- EMR connectivity via serial output or optional wireless module
- Handrail scale platform dimensions: 16 Inch width × 18 Inch depth × 1.75 Inch height/Chair scale mode has 21 Inch wide seat bottom with 3 Inch padded cushion
- 3-year manufacturer's warranty



SCAL38CH in Chair Scale Mode



SCAL38CH in Scale Mode

Item No.	Description	Quantity
SCAL38CH	Scale Chair Convertible with Handrail 500 Pound Capacity Metal 16 Inch × 18 Inch × 1.75 Inch	1/EA

SCALES AND MEASURING DEVICES

Portable Handrail Scales

made in the USA NOT MADE WITH natural rubber latex

- Sturdy wraparound handrail provides maximum security, as well as a full 25 Inch span between grips
- Large, lightweight platform with heavy-duty 3 Inch wheels that tilt and roll effortlessly on all flooring surfaces for maximum portability
- SCAL88DH scale features integrated digital height rod that quickly and accurately measures height with no buttons to push
- 3-year manufacturer's warranty
- MotionLock Technology automatically calculates average weight reading for patients that are unable to remain still or are unsteady
- Console displays weight, height, and BMI on a single screen
- Features include BMI calculator, auto zero, weight recall, auto weigh, always-on option available with no buttons to press for normal weighing, tare, and more
- Programmable auto-off: 12 seconds to 25 minutes
- Kilogram/Pound Lockout: quickly change setup to disable LB/KG button so scale will only read in selected unit of measure
- Quick two second weight measurement accurate within 0.1 pounds (0.1 kilograms.) and 750 pounds (340 kilograms) capacity
- EMR connectivity via serial output or optional wireless module
- Powered by 4 D-cell alkaline batteries (included) that last 100,000 weighings



SCAL38MC



SCAL38MC

Item No.	Description	Quantity
Scale Portable 750 Pound Capacity		
SCAL38MC	With Wraparound Handrail 16 Inch × 18 Inch × 1.75 Inch	1/EA
SCAL88DH	With Wraparound Handrail and Integrated Digital Height Rod	1/EA

SCALES AND MEASURING DEVICES

Folding Wheelchair Scales with Handrail

made in the **USA** NOT MADE WITH natural rubber latex

- Ergonomic folding design (12 Inch x 32 Inch footprint) with four wheels is self-supporting for safe transport and storage
- Two platform sizes available to accommodate all wheelchair widths
- Display console pivots 320 degree for optimal viewing
- All metal construction for long lasting performance
- Scale arrives fully assembled—just put in leveling feet and batteries
- 3-year manufacturer’s warranty
- MotionLock Technology calculates accurate weight even with unsteady patients
- Console displays weight, height, and BMI on a single screen
- Features include stored tare, auto zero, weight recall, programmable auto-off from 12 seconds to 25 minutes, BMI calculator, and more
- Pound/Kilogram Lockout: quickly change setup to disable Pound/Kilogram button so scale will only read in selected unit of measure
- EMR connectivity via serial output or optional wireless module
- Powered by 4 D-cell alkaline batteries (included) that last 100,000 weighings
- Capacity: 1000 pounds in 0.1 pound increments
- Accurate weighing even when patient applies weight to heavy-duty handrail



Item No.	Description	Quantity
Scale Wheelchair Folding 1000 Pound Capacity		
SCAL41MC	With Handrail Single Ramp 32 Inch x 34 Inch	1/EA
SCAL42MC	With Handrail Dual Ramp 32 Inch x 34 Inch	1/EA
SCAL49MC	With Handrail Dual Ramp 32 Inch x 41 Inch	1/EA

SCALES AND MEASURING DEVICES

Tilt and Roll Wheelchair Scales with Handrail

made in the **USA** NOT MADE WITH natural rubber latex

- Versatile scale that can be used as a wheelchair, handrail, chair, or dialysis scale (Single and Dual Ramp models available)
- Extra-large platform: 28 Inch length x 32 Inch width x 2 Inch height with fold-down ramp(s) to accommodate wheelchairs, recliners, and patients
- MotionLock Technology automatically calculates average weight reading for patients that are unable to remain still or are unsteady
- Console displays weight, height, and BMI on a single screen
- Features include stored tare, auto zero, weight recall, programmable auto-off from 12 seconds to 25 minutes, BMI calculator, and more
- Pound/Kilogram Lockout: quickly change setup to disable Pound/Kilogram button so scale will only read in selected unit of measure
- EMR connectivity via serial output or optional wireless module
- Powered by 4 D-cell alkaline batteries (included) that last 100,000 weighings
- Heavy-duty aluminum and steel construction with sturdy handrail that allows for greater stability when used as a stand-on scale and easy access on and off wheelchair ramp
- Capacity: 1000 pounds/450 kilograms in 0.1 pound/kilogram increments
- Tilt-and-roll design and large casters allows for easy maneuverability
- Display console pivots 320 degrees for optimal viewing
- 3-year manufacturer’s warranty



Item No.	Description	Quantity
Scale Wheelchair Tilt and Roll 1000 Pound Capacity Aluminum Steel 28 Inch x 32 Inch x 2 Inch		
SCAL45MC	Single Ramp with Handrail	1/EA
SCAL46MC	Dual Ramp with Handrail	1/EA

SCALES AND MEASURING DEVICES

Take Home Scale

NOT MADE WITH natural rubber latex **diverse supplier**



- Triple measurement in pounds, kilograms, or stones
- Default setting is in pounds
- Switch to convert in kilograms or stones
- Strain gauge sensor system
- Auto Step-On technology, no tapping activation
- Automatic Zero Reset and automatic Turn-Off
- Lithium battery powered ONLY - one (1) of 3 volt CR2032
- Anti-Tipping mechanism
- Weight Capacity: 440 pound/200 kilogram/31.43 stone
- Resolution: 0.2 pound/0.1 kilogram
- Platform Size: 13 Inch width × 12.6 Inch length (33 centimeters × 32 centimeters)
- Scale Dimensions: 13 Inch width × 12.6 Inch length × 0.9 Inch height (33 centimeters × 32 centimeters × 2.3 centimeters)
- Display Size: 3.1 Inch width × 2 Inch length (8 centimeters × 5.15 centimeters)
- 1-year limited manufacturer's warranty

Item No.	Description	Quantity
SCAL101	Scale Take Home Digital for Patient Use White 440 Pound Capacity 13 Inch × 12.6 Inch	1/EA, 6EA/CA

HAMPERS

Chrome Linen Hampers with Optional Labels

NOT MADE WITH natural rubber latex



- Solid 1 Inch chrome-plated steel tubing with lower support platform for heavyweight applications
- Ivory ABS Plastic Lid
- Hands-free foot pedal actuation opens lid wide for easy access and closes automatically to help contain odors
- Upper locking ring holds bags of at least 30 Inch wide securely in place
- Open front design provides convenient storage for additional bags
- Four 3 Inch, free-turning casters—two with locking foot brakes—provide easy maneuverability
- Soiled Linen label and Isolation Linen label included
- Dimensions: 19 Inch width × 20 Inch depth × 38 Inch height

Item No.	Description	Quantity
HAMP9000	Hamper Linen Square Foot-Operated Ivory ABS Plastic Lid Chrome-Plated Steel Frame 3 Inch Casters 19 Inch × 20 Inch × 38 Inch Includes Lid Labels for Soiled Linen and Isolation Linen	1/EA

HAMPERS

Adjustable Folding Chrome Hampers

NOT MADE WITH natural rubber latex **diverse supplier**



- Chrome-plated steel frame
- Spring-loaded, foot-operated lid
- M-style base with bag support bars
- 4 ball-bearing casters
- Chrome-plated or biohazard red lid available
- Dimensions: 19.5 Inch depth × 36.5 Inch height
- 1-year manufacturer's warranty

Item No.	Description	Quantity
HAMP8001	Hamper Linen Adjustable Folding Chrome-Plated Steel Frame Fits 18 Inch, 22 Inch, or 25 Inch Bags 19.75 Inch × 42 Inch × 37 Inch	1/EA

Square Tilt-Lid Chrome Hampers

NOT MADE WITH natural rubber latex



- Chrome-plated steel frame
- Pedal-activated lid
- 4 Inch casters with 2 locks
- Heavy-duty bag support and holder
- White epoxy-coated steel lid
- Prop to keep lid open
- Dimensions: 18.625 Inch × 37.5 Inch × 19.5 Inch

Item No.	Description	Quantity
HAMP8002	Hamper Linen Square Foot-Operated White Epoxy Steel Tilt-Lid Chrome-Plated Steel Frame 18.625 Inch × 37.5 Inch × 19.5 Inch	1/EA

18" Round Tilt-Lid Chrome Hampers

NOT MADE WITH natural rubber latex **diverse supplier**



- Chrome-plated steel frame
- Spring-loaded, foot-operated lid
- M-style base with bag support bars
- 4 ball-bearing casters
- Chrome-plated
- Dimensions: 19.5 Inch depth × 36.5 Inch height
- 1-year manufacturer's warranty

Item No.	Description	Quantity
HAMP8003	Hamper Linen Round Foot-Operated Tilt-Lid Chrome-Plated Steel Frame 4 Ball-bearing Casters 19.5 Inch × 36.5 Inch	1/EA

Please see our full line of Linen Bags on page 104 in the Housekeeping Section.



HAMPERS

18" Round Chrome Hampers

NOT MADE WITH natural rubber latex **diverse supplier**

- Chrome plated steel frame
- Bag-support bars
- 3 ball-bearing casters
- Dimensions: 19 Inch depth x 35 Inch height
- 1-year manufacturer's warranty



Item No.	Description	Quantity
HAMP8004	Hamper Linen Round Bag-Support Bar Chrome-Plated Steel Frame 3 Ball-bearing Casters 19 Inch x 35 Inch	1/EA

HAMPER ACCESSORIES

Hamper Accessories

NOT MADE WITH natural rubber latex

- Replacement bag for HAMP8002 and HAMP9000
- 14 Inch x 14 Inch
- White



Item No.	Description	Quantity
HAMP1000	Bag Hamper Replacement for Hamp8002 and Hamp9000 14 Inch x 14 Inch White Non Sterile	1/EA

IV STANDS

Chrome IV Stands

NOT MADE WITH natural rubber latex

- Chrome-plated steel pole with single knob height adjustment
- Available with either two or four detachable ram's horn hooks with plastic tips
- Models 1010 and 1012 with 5-leg cast aluminum base and dual-wheel casters
- Models 1011 and 1013 with epoxy-coated steel 6-leg base and 3 Inch ball bearing casters
- Load capacity: 45 pounds
- Height adjustment:
 - Models IVCR1010 and ICVR1012: 52 Inches to 90 Inches
 - Models IVCR1011 and ICVR1013: 58 Inches to 99 Inches



Item No.	Description	Quantity
Stand IV Chrome Plated Pole 45 Pound Steel		
IVCR1010	2 Hooks 5 Legs	1/EA
IVCR1011	2-Hook 6-Leg Epoxy	1/EA
IVCR1012	4-Hook 5-Leg	1/EA
IVCR1013	4-Hook 5-Leg Epoxy	1/EA

IV STANDS

Heavy Base Chrome IV Stands

NOT MADE WITH natural rubber latex

- Chrome-plated steel pole/knob-lock adjustment
- 4-leg epoxy coated steel base with bumper guards
- Detachable ram's horn hooks with plastic tips
- 2 Inch rubber wheel ball bearing casters
- Capacity: 45 pounds
- Available in 2 or 4 hooks
- Base diameter: 22 Inch
- Height Range: 52 Inch - 91 Inch
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
Stand IV Heavy Base Chrome-Plated Steel 22 Inch Base Diameter Height Adjustable From 52 Inch to 91 Inch 45 Pound Capacity		
IVCR24H	2 Hooks 4 Legs	1/EA
IVCR44H	4 Hooks 4 Legs	1/EA

Stainless IV Stands

NOT MADE WITH natural rubber latex

- Sturdy, functional design mounted on 5-leg non-tip 24 Inch diameter reinforced plastic base for stability
- Fitted with 2 Inch swivel casters
- Telescopic, stainless steel rod with positive-locking hand wheel adjustments on 1 Inch diameter tubing
- 92.5 Inch extended height
- 2 or 4-ram's horn hooks
- Minimal assembly required



Item No.	Description	Quantity
Stand IV Stainless Steel		
IVSS1050	2 Hooks 5 Legs	1/EA
IVSS1051	4 Hooks 5 Legs	1/EA

Please see our full line of Sponge Counter Bags on page 145 in the OR Products Section.



IV STAND ACCESSORIES

IV Pole Outlet Strips

NOT MADE WITH
natural rubber latex

- Four (4) hospital grade outlets with covers
- 15 foot heavy-duty cord
- Stainless steel frame with hood
- Mounting requires 1 universal clamp - IVUC41 (sold separately)
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
IVOS50	Outlet Strip IV Pole Four Outlets 15 Foot Heavy-Duty Cord Stainless Steel Frame	1/EA

Organizer Hooks

NOT MADE WITH
natural rubber latex

- Drainage Bag Hooks (IVBH56)
 - Chrome plated with 2-hook positions
 - Clamps to pole with thumb screws, no universal clamp needed
- Extra Hooks and Organizer Hooks (IVOH55)
 - 4 ram's horn hooks and 4 open hooks
 - All stainless steel hooks
 - Not M.R.I. compatible
- Single Ram's Horn Accessory Hook
 - Requires 1 universal clamp - IVUC41 (sold separately)
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
Hooks Organizer		
IVBH56	Drainage Bag 2 Hooks Chrome	1/EA
IVOH55	4 Ram's Horn Hooks and 4 Open Hooks Stainless Steel	1/EA
IVSH57	Single Ram's Horn Accessory Hook	1/EA

Stability Weights

NOT MADE WITH
natural rubber latex

- Adds 2.5 or 5 pounds to base to improve stability
- Easily mounts to center of base
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
Weight Base Stability For IV Stand		
IVBW48	2.5 Pound	1/EA
IVBW485	5 Pound	1/EA

IV STAND ACCESSORIES

Oxygen Cylinder Holders

NOT MADE WITH
natural rubber latex

- Chrome plated steel frame
- Holds D or E size oxygen cylinders
- Requires 1 universal clamp - IVUC41 (sold separately)
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
IVCH44	Holder Cylinder Oxygen Chrome-Plated Steel for D or E Size	1/EA

Infusion Pump Supports

NOT MADE WITH
natural rubber latex

- Infusion pump frame holds up to four IV pumps
 - Frame constructed of 1 inch chrome plated steel
 - Requires 2 universal clamps - IVUC41 (sold separately)
 - Pump frame only suitable for 6-leg IV stand models
- Pump support tray provides easy pump attachment
 - Requires 1 universal clamp - IVUC41 (sold separately)
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
IVPF42	Frame Infusion Pump Chrome Plated	1/EA
IVST46	Tray Infusion Pump Support	1/EA

Steering Handles

NOT MADE WITH
natural rubber latex

- Easy-grip black vinyl coating
- Aids movement and control of poles
- Requires 1 universal clamp - IVUC41 (sold separately)
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
IVSH43	Handle Steering Easy Grip Black Vinyl Coating	1/EA

Universal Clamps

NOT MADE WITH
natural rubber latex

- Aluminum clamp used to mount most accessories without disassembling pole
- Easy position adjustments
- Set-screw locks accessories in place
- Not included with accessories
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
IVUC41	Clamp Universal Pole Mount Fits All MediChoice IV Stands Aluminum	1/EA

IV STAND ACCESSORIES

Heavy-Duty Wire Baskets

NOT MADE WITH natural rubber latex

- Available in either epoxy-coated steel or stainless steel
- Two sizes to choose from:
14 Inch x 8 Inch x 8.5 Inch or
12 Inch x 6 Inch x 6.5 Inch
- Requires 1 universal clamp - IVUC41 (sold separately)
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
Basket IV Stand Wire Heavy-Duty Fits All MediChoice IV Stands		
IWWB51B	Large Epoxy-Coated Steel 14 Inch x 8 Inch x 8.5 Inch	1/EA
IWWB51S	Large Stainless Steel 14 Inch x 8 Inch x 8.5 Inch	1/EA
IWWB52B	Small Epoxy-Coated Steel 12 Inch x 6 Inch x 6.5 Inch	1/EA
IWWB52S	Small Stainless Steel 12 Inch x 6 Inch x 6.5 Inch	1/EA

MAYO STANDS

Stainless Mayo Stands

made in the USA NOT MADE WITH natural rubber latex diverse supplier

- Non-tip weighted flat base slides easily under OR table
- Removable stainless steel tray for easy cleaning
- Mounted on 2 Inch ball-bearing swivel casters
- Foot-pedal adjustment
- Tray size: 16 Inch width x 21 Inch length
- Height: 38.25 Inch to 63.875 Inch
- 3-year manufacturer's warranty



Item No.	Description	Quantity
MAY01251	Stainless Steel Mayo Stand, 16 Inch x 21 Inch Tray	1/EA

Mobile Instrument Stands

NOT MADE WITH natural rubber latex

- Knob height adjustment for comfortable use while sitting or standing, 31 Inch-50 Inch height range
- Stable 5-leg cast aluminum base
- 1 Inch chrome-plated steel tubing
- Designed for limited space environments
- Removable stainless steel tray:
19 Inch length x 12.75 Inch width x .66 Inch height
- Dual wheel soft roll casters
- Some assembly required
- Shipping weight: 15 pounds



Item No.	Description	Quantity
Stand Instrument Mobile Five-Leg Cast Aluminum Base		
MAY02000	Stainless Steel Tray 19 Inch x 12.75 Inch x 0.66 Inch	1/EA

MAYO STANDS

Chrome Double Post Mayo Stands

NOT MADE WITH natural rubber latex

- Removable stainless steel tray
- 2 ball-bearing rubber casters
- 1 Inch chrome-plated poles and base
- Heavy, rectangular easy position base:
19.75 width x 19.24 Inch length
- Tray size: 12.75 Inch width x 19 Inch length
- Adjustable height: 34 Inch - 54 Inch



Item No.	Description	Quantity
Stand Mayo Double Post Stainless Steel Tray		
MAY01201	Double Post with Tray Chrome 2 Ball-Bearing Casters Height Adjustable from 34 Inch to 54 Inch	1/EA

BASIN STANDS

Basin Stands

made in the USA NOT MADE WITH natural rubber latex diverse supplier

- Available in single-basin and double-basin stainless designs
- Lower shelf for extra storage (edges bent and hemmed for stability and safety)
- Open frame design allows easy access from any side
- 16 gauge, type 304 stainless steel frame
- 2 Inch diameter ball-bearing swivel casters
- Easy to clean
- Dimensions:
 - Single:
15 Inch depth x 15 Inch width x 32 Inch height
 - Double:
15 Inch depth x 29 Inch width x 32 Inch height
 - Basin Diameter: 13.25 Inch deep
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Stand Basin Single with Shelf 7 Quart Stainless Steel		
BAST1001	15 Inch x 15 Inch x 32 Inch	1/EA
BAST1002	15 Inch x 29 Inch x 32 Inch	1/EA

STEP STOOLS

Chrome Step Stools

NOT MADE WITH
natural rubber latex

- Chrome powder-coated finish with all welded steel base and top, no screws or rivets
- Rubber tread top for safety
- Ideal for exam rooms and x-ray areas
- Available with .875 Inch chrome plated steel handrail 35.5 Inch height
- Capacity: 350 pounds
- Reinforced rubber feet
- 1 Inch steel tube with X-frame top support
- Some assembly required
- Dimensions: 9 Inch height × 11.25 Inch width × 14.25 Inch length



Item No.	Description	Quantity
Stool Step 350 Pound Chrome 9 Inch × 1.25 Inch × 14.25 Inch		
FT3500	Stool Step	1/EA
FT3500H	Stool Step with Handrail	1/EA

Stainless Stackable Step Stools

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Slip-resistant top
- Heavy gauge legs and handles
- Stacks safely for increased height
- Can interlock to create larger platforms
- Stools sold individually; pictured with optional mobile cart
- Capacity: 400 Pounds
- Extra-large platform dimensions: 12 Inch depth × 18 Inch width × 6 Inch height
- 3-year manufacturer's warranty



Item No.	Description	Quantity
STEP1750	Stool Step Stackable Stainless Steel 12 Inch × 18 Inch × 6 Inch 400 Pound Capacity	1/EA
CART1750	Cart Transport for MediChoice Stackable Step Stool STEP1750 18 Inch Wide × 12 Inch Deep × 10 Inch High	1/EA

STEP STOOLS

Stainless Step Stools

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Strong welded stainless steel
- Tip-resistant, rugged, and easy to maintain
- Slip-resistant top
- Available with 38.5 Inch handrail
- Capacity: 400 pounds
- Slip-resistant floor tips
- Platform dimensions: 12 Inch depth × 18 Inch width × 8 Inch height
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Stool Step 12 Inch × 18 Inch × 8 Inch 400 Pound Capacity		
STEP1751	Stainless Steel	1/EA
STEP1752	With Handrail Stainless Steel	1/EA

BASSINETS

Bassinets

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Basket may be rested horizontally or Trendelenburg can be achieved by elevating foot section of basket
- Full extension drawer mounted on ball-bearing slides
- Lower shelf for extra storage (edges bent and hemmed for stability and safety)
- 3 Inch ball-bearing swivel casters
- Drawer size: 24 Inch length × 12 Inch depth × 5 Inch height
- Overall dimensions: 32 Inch width × 17 Inch depth × 35.25 Inch length
- 3-year manufacturer's warranty



Item No.	Description	Quantity
BASN2050	Bassinet Medical with Drawer Basket Pad Stainless Steel 32 Inch × 17 Inch × 35.25 Inch Reusable	1/EA

BLANKET WARMERS



BLANKET WARMERS



Counter Top Blanket Warmers



- Single compartment stainless steel blanket warmer with digital temperature control to ensure fast and easy operation
- Cuverro copper alloy antimicrobial handles for infection prevention and control
- Glass or solid self-closing, energy efficient door
- 1 fully adjustable perforated solid stainless steel shelf
- UL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140° F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 5.8 amps
- Dimensions:
 - Overall Size: 30 Inch width × 24 Inch depth × 24 Inch height
 - Compartment Size: 28 Inch width × 17.5 Inch depth × 16 Inch height (4.5 cubic feet)
- 5-year manufacturer's warranty on parts, 1-year on labor

Item No.	Description	Quantity
Warmer Blanket Countertop Single Door 30 Inch Wide × 24 Inch Deep × 24 Inch High		
BW3024G	Glass Door	1/EA
BW3024S	Solid Door	1/EA
Blanket Warmer Stand Countertop		
BW3530ST	Non-Mobile 36 × 24 × 29 Inch	1/EA
BW3534ST	Mobile 36 × 24 × 34 Inch	1/EA

Under Counter Blanket Warmers



- Single compartment stainless steel blanket warmer with digital temperature control to ensure fast and easy operation
- Glass or solid self-closing, energy efficient door
- Cuverro copper alloy antimicrobial handles for infection prevention and control
- 3 fully adjustable perforated solid stainless steel shelves
- UL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140° F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 7.3 amps
- Dimensions:
 - Overall Size: 30 Inch width × 24 Inch depth × 72 Inch height
 - Compartment Size: 28 Inch width × 17.5 Inch depth × 65 Inch height (17.7 cubic feet)
- 5-year manufacturer's warranty on parts, 1-year on labor

Item No.	Description	Quantity
Warmer Blanket Undercounter Single Door 30 Inch Wide × 24 Inch Deep × 36 Inch High		
BW3036G	Glass Door	1/EA
BW3036S	Solid Door	1/EA

BLANKET WARMERS

Double Door Blanket Warmers



- Dual compartment stainless steel blanket warmer with independent digital temperature controls to ensure fast and easy operation
- Glass or solid self-closing, energy efficient door
- Cuverro copper alloy antimicrobial handles for infection prevention and control
- Upper and lower compartment each has 2 adjustable perforated solid stainless steel shelves
- UL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140° F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 7.3 amps
- Dimensions:
 - Overall Size: 30 Inch width × 24 Inch depth × 72 Inch height
 - Compartment Size: 28 Inch width × 17.5 Inch depth × 30 Inch height upper/35 Inch height lower (17 total cubic feet)
- 5-year manufacturer's warranty on parts, 1-year on labor



Item No.	Description	Quantity
Warmer Blanket Double Doors 30 Inch Wide × 24 Inch Deep × 72 Inch High		
BW3072DG	Glass Doors	1/EA
BW3072DS	Solid Doors	1/EA

Single Door Blanket Warmers



- Single compartment stainless steel blanket warmer with digital temperature control to ensure fast and easy operation
- Glass or solid self-closing, energy efficient door
- Cuverro copper alloy antimicrobial handles for infection prevention and control
- 3 fully adjustable perforated solid stainless steel shelves
- UL approved and meets AORN and ECRI guidelines for patient safety
- Temperature range 90°–140° F; AORN and ECRI recommended temperature 120°F
- Power Supply: 120V, 50/60 cycles, and 7.3 amps
- Dimensions:
 - Overall Size: 30 Inch width × 24 Inch depth × 72 Inch height
 - Compartment Size: 28 Inch width × 17.5 Inch depth × 65 Inch height (17.7 cubic feet)
- 5-year manufacturer's warranty on parts, 1-year on labor



Item No.	Description	Quantity
Warmer Blanket Single Door 30 Inch Wide × 24 Inch Deep × 72 Inch High		
BW3072G	Glass Door	1/EA
BW3072S	Solid Door	1/EA
Mobile Dolly		
BW3504	Mobile Dolly for Blanket Warmers Stainless Steel	1/EA

Instrument Tables



- 20 Inch width and 24 Inch width tables have stiffeners of 20-gauge #304 stainless steel laminated and welded to tops
- Tops and shelves are fabricated of 20-gauge #304 stainless steel
- Tops and shelves are hemmed and formed for extra perimeter strength
- Finished smooth edges for safety
- 2 Inch ball-bearing swivel casters
- Eight sizes to choose from
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Table Instrument Stainless Steel		
INTB9001	16 Inch × 20 Inch × 34 Inch	1/EA
INTB9002	16 Inch × 30 Inch × 34 Inch	1/EA
INTB9003	18 Inch × 33 Inch × 34 Inch	1/EA
INTB9004	20 Inch × 36 Inch × 34 Inch	1/EA
INTB9005	20 Inch × 48 Inch × 34 Inch	1/EA
INTB9006	24 Inch × 48 Inch × 34 Inch	1/EA
INTB9007	24 Inch × 60 Inch × 34 Inch	1/EA
INTB9008	24 Inch × 72 Inch × 34 Inch	1/EA

Kick Buckets



- Welded stainless steel construction with non-tip mobile base
- Removable 13-quart stainless bucket
- Perimeter rubber bumper to protect walls
- Rubber cushions on ring stand prevent bucket from rattling
- 2 Inch swivel casters for easy mobility
- Dimensions: 13.5 Inch depth × 14 Inch height
- 3-year manufacturer's warranty



Item No.	Description	Quantity
BUKT1101	Bucket Kick With Bumper 2 Inch Casters 13 Quart Stainless Steel 13.5 Inch × 14 Inch	1/EA

Complementary Products

Item No.	Description
H4830YW	Kick Bucket Liner 24 Inch × 30 Inch 1.00 Mil Anti Static Fits Standard 13 Quart Kick Buckets Meets Standard NFPA 56A Capacity 8 Gallon



Waste Receptacles

- 32-quart enamel step-on waste receptacle
- Removable plastic inner liner
- Perfect for exam rooms, restrooms and front office areas
- Step-on foot pedal opens lid
- Meets OSHA standard on blood-borne pathogens
- Steel outer can and lid
- Hazardous waste label included
- Available in white, red, or beige finish
- Weight: 12 lbs.
- Dimensions: 18.25 Inch H x 12.25 Inch W x 11.5 Inch D
- 5-year manufacturer's warranty

Item No.	Description	Quantity
Receptacle Waste Examination Room with Step-On Foot Pedal 32 Quart 18.25 Inch x 12.25 Inch x 11.5 Inch		
RECP2101	Steel White	1/EA
RECP3101	Steel Red	1/EA
RECP4101	Steel Beige	1/EA

LIGHTING-EXAM AND PROCEDURE

Multi-Purpose Mobile Exam Lights

NOT MADE WITH natural rubber latex

- 30 Inch polyvinyl coated steel neck for easy maneuvering
- Patented heat sink and dual shade system won't heat up, even after hours of use
- Shade rotates 1.5 revolutions without wire damage
- Uses standard incandescent (up to 100 watts) or halogen (60 watts) bulbs
- Equipped with 4300 degree Kelvin rating color-corrected, full spectrum halogen bulb
- Easy rolling anti-tip 5-caster base that adjusts from 2-7 foot with a simple twist of the wrist
- Snap-lock coupler allows easy removal of head and neck assembly
- Hospital grade electrical plug
- Available in black, beige or grey
- 5-year manufacturer's warranty



Item No.	Description	Quantity
LGHT7502	Light Examination Multi-Purpose Mobile 30 Inch Flexible Neck	1/EA

Economy Floor Exam Lights

made in the USA **NOT MADE WITH natural rubber latex** **diverse supplier**

- Chrome-plated shaft and pedestal base
- 12 Inch gooseneck for easy positioning
- 360 degree lamp adjustment
- Telescopes 50 Inch to 74 Inch
- Dual shade, 60-watt maximum bulb (not included)
- Hospital grade electrical plug
- 1-year manufacturer's warranty




Item No.	Description	Quantity
LGHT7503	Light Examination 12 Inch Flexible Gooseneck Chrome-Plated with 60 Watt Bulb	1/EA

Procedure Exam Lights

made in the USA **NOT MADE WITH natural rubber latex**

- UL-approved economical 35 watt (4166 foot-candle) halogen light with flexible gooseneck
- High intensity, cool operation
- 33 Inch effective swing radius
- Superior flexibility
- 2 lockable wheels
- No height adjustment
- 3200 degree K color temperature
- 2000-hour bulb life
- Available in wall mount or floor model
- 5-year manufacturer's warranty



Item No.	Description	Quantity
Light Examination Procedure with 35-Watt Halogen Bulb		
LGHT7001	Universal Mount	1/EA
LGHT7002	Floor Stand	1/EA

Please see our full line of Can Liners on pages 100-103 in the Housekeeping Section.

LIGHTING-EXAM AND PROCEDURE



Procedure Exam Lights - K Arm - Halogen



- UL-approved high intensity halogen exam light
- K-Arm technology is highly positionable, with internal springs to ensure the light stays in the desired position without drifting
- Drift-free K-arm with 42 Inches (107 centimeters) arm range
- 3800 degree K Color Temperature
- 99 Color Rendering Index (CRI)
- 45,750 Lux at 24 Inches
- 5-Year manufacturer's warranty

Item No.	Description	Quantity
Light Examination Procedure K-Arm Super Exam 50 with Halogen Bulb		
LGHT7100	Floor Stand	1/EA
LGHT7101	Wall Mount	1/EA
LGHT7102	Rail Mount	1/EA
LGHT7103	Ceiling Mount	1/EA
LGHT7104	Table Mount	1/EA

LIGHTING-EXAM AND PROCEDURE



Procedure Exam Lights - K Arm - LED



- UL-approved high intensity LED exam light
- K-Arm technology is highly positionable, with internal springs to ensure the light stays in the desired position without drifting
- LED source delivers a natural white light and uses 70 percent less energy
- LED sources last 40,000 hours
- 4000 degree K Color Temperature
- 92 Color Rendering Index (CRI)
- 39,000 Lux at 24 Inches
- 5-Year manufacturer's warranty

Item No.	Description	Quantity
Light Examination Procedure K-Arm Super Exam 50 with Light-Emitting Diode Bulb		
LGHT7200	Floor Stand	1/EA
LGHT7201	Wall Mount	1/EA
LGHT7202	Rail Mount	1/EA
LGHT7203	Ceiling Mount	1/EA
LGHT7204	Table Mount	1/EA

K Arm Task Lights



- UL-approved drift free task light
- K-Arm technology is highly maneuverable, with internal springs to ensure the light stays in the desired position without drifting
- 2900 degree K color temperature
- 97 Color Rendering Index (CRI)
- 30,000 Lux at 18 Inches
- Available in floorstand, wall mount, or table mount
- 2-Year manufacturer's warranty



Item No.	Description	Quantity
LGHT7700	Light Task K-Arm Floor Stand 120Volt	1/EA
LGHT7701	Task Light 120 V Halogen K-Arm With 42 Inch Arm Range Wall Mounted UI 153 Certified	1/EA

LIGHTING-EXAM AND PROCEDURE

Procedure Exam Lights-LED

made in the USA NOT MADE WITH natural rubber latex

- The LED source provides a virtually heat-free light so the patient remains comfortable
- Dual handles are easily accessible to position the light head
- 58 percent energy savings over comparable halogen options (at low light level setting) and 18 percent savings when providing double the illumination
- The patented K-arm moves smoothly with no drifting
- Delivers a uniform beam pattern that aids the examination process
- Dual intensity switching allows a choice between two illumination levels
- The low level delivers a 19 percent increase in light level
- Maximum arm reach of 45 Inches (114.30 centimeters)
- Mounting options (floorstand or wall mount) to fit the application
- Delivers consistent light output during power supply fluctuations
- The LED source offers 10 times the operational life of traditional halogen solutions



Item No.	Description	Quantity
Light Procedure LED Nova Exam Light		
LGHT8300	Floor Stand	1/EA
LGHT8301	Wall Mount	1/EA

LIGHTING-MINOR SURGERY AND DIAGNOSTIC

Diagnostic and CoolSpot Lamps

made in the USA NOT MADE WITH natural rubber latex

- UL-approved cool operation lamp with pure white, intense light
- Side friction knob assures a drift-free arm system
- Shadow-free operation by three 50-watt halogen bulbs (5,600 hour average bulb life)
- Floorstand model includes footswitch, locking casters, and cord wrap
- Central, removable, autoclavable SteriHandle™
- 8000 foot candles at 24 Inches
- 3300 degree K color temperature
- Lamp housing and stand are off-white color
- 5-year manufacturer's warranty



Item No.	Description	Quantity
Light Minor Surgery Coolspot II with Three 50-watt Halogen Bulbs		
LGHT8100	Floor Stand	1/EA
LGHT8101	Single Ceiling Mount	1/EA
LGHT8102	Double Ceiling Mount	1/EA
LGHT8103	Wall Mount	1/EA
LGHT8104	Single Head Trolley	1/EA
LGHT8105	Dual Head Trolley	1/EA

LIGHTING-MINOR SURGERY AND DIAGNOSTIC

Minor Surgery Lighting - AIM 50

made in the USA NOT MADE WITH natural rubber latex

- High-intensity of 45,000 lux (4181 fc) at 1 meter
- 3100 K color temperature offers natural and true coloring
- High CRI (Color Rendering Index) of 96
- Uses three quartz halogen bulbs, each rated with a 2000-hour bulb life; bulbs wired in parallel for continuous operation
- Compact 20 Inch (51 Centimeter) diameter lighthead
- Extremely smooth focus mechanism with center sterilizable handle
- 360 degree limitless arm-and-mounting-system rotation around vertical axes (on single ceiling models)
- Cool operation for long hours of use
- 120V standard; 100 and 230/240V models available
- 5-year limited warranty



Item No.	Description	Quantity
Light Surgery Minor Aim 50 with Halogen Bulb		
LGHT8200	Floor Stand	1/EA
LGHT8201	Single Ceiling	1/EA

Minor Surgery Lighting - Halogen

made in the USA NOT MADE WITH natural rubber latex

- High-intensity outpatient surgery or diagnostic light that is available in floorstand, single or double ceiling mounts, and wall-mounted models
- Uses three 35-watt quartz halogen bulbs—each rated with 2000 hour bulb life, wired in parallel for continuous operation
- Shadowless lighting from the Y-shaped light head
- 3500 degree K Color Temperature
- 97 Color Rendering Index (CRI)
- 63,000 Lux at 1 meter
- 63 Inch extension and drift-free positioning
- Floorstand model features base with easy-access footswitch and built-in casters with locking mechanism on back caster
- 5-Year manufacturer's warranty



Item No.	Description	Quantity
Light Surgical Minor Aim 100 with 35-Watt Quartz Halogen Bulb		
LGHT7300	Floor Stand	1/EA
LGHT7301	Single Ceiling	1/EA
LGHT7302	Double Ceiling	1/EA
LGHT7303	Wall Mount	1/EA

LIGHTING-MINOR SURGERY AND DIAGNOSTIC

Minor Surgery Lighting - LED

made in the USA NOT MADE WITH natural rubber latex

- High-intensity outpatient surgery or diagnostic light that is available in floorstand, single or double ceiling mounts, and wall-mounted models
- Multiple reflectors can be tilted using the light handle to adjust size and pattern of the light to the working area
- Shadowless lighting from the Y-shaped light head
- LED source uses 55 percent less energy
- LED sources last 20,000 hours—10 times as long as traditional halogen sources—saving over 1,000 in bulb replacements
- 4300 degree K Color Temperature
- 92 Color Rendering Index (CRI)
- 45,000 Lux at 1 meter; 110,000 Lux at 24 Inch; 10 percent more illumination than halogen
- 5-Year manufacturer's warranty



Item No.	Description	Quantity
Light Surgical Minor Aim with Light-Emitting Diode		
LGHT7400	Floor Stand	1/EA
LGHT7401	Single Ceiling	1/EA
LGHT7402	Double Ceiling	1/EA
LGHT7403	Wall Mount	1/EA

Diagnostic and Outpatient Lamps

made in the USA NOT MADE WITH natural rubber latex

- UL-approved cool operation lamp with pure white, intense light
- Side friction knob assures a drift-free arm system
- Shadow-free operation by three 50-watt halogen bulbs (5,600 hour average bulb life)
- Floorstand model includes footswitch, locking casters, and cord wrap
- Central, removable, autoclavable SteriHandle™
- 8000 foot candles at 24 Inches
- 3300 degree K color temperature
- Lamp housing and stand are off-white color
- 5-year manufacturer's warranty



Item No.	Description	Quantity
Lamp Outpatient Diagnostic 115-Volt		
LGHT7600	Floor Stand	1/EA
LGHT7601	Single Ceiling Mount	1/EA
LGHT7602	Double Ceiling Mount	1/EA
LGHT7603	Wall Mount	1/EA
LGHT8104	Single Head Trolley	1/EA
LGHT7605	Dual Head Trolley	1/EA

BARIATRICS

Bariatric Step Stools

NOT MADE WITH natural rubber latex

- Durable chrome finish with steel top
- Welded reinforced x-frame base
- Optional model with handrail: 36 Inch Height
- Capacity: 600 pounds
- Dimensions:
 - Standard Size: 14.25 Inch Length x 9 Inch Height x 11.25 Inch Width
 - Large Top: 17.25 Inch Length x 9 Inch Height x 12.25 Inch Width



Item No.	Description	Quantity
Stool Step Bariatric 14.25 Inch x 9 Inch x 11.25 Inch Chrome Finish Steel Top 600 Pound Capacity		
BAR3000	Standard Platform	1/EA
BAR300H	With Handrail Standard Platform	1/EA
Stool Step Bariatric 17.25 Inch x 9 Inch x 12.25 Inch Chrome Finish Steel Top 600 Pound Capacity		
BAR300HL	With Handrail Large Platform	1/EA
BAR300L	Large Platform	1/EA

EXAM STOOLS

Economy Exam Stools

made in the USA NOT MADE WITH natural rubber latex

- Full chrome frame
- 1 Inch diameter, S-style all-welded frame
- Soft-feel poly foam padding
- 4 Inch thick seat for extra comfort
- 2 Inch rubber wheel, ball bearing casters
- Round foot ring
- Solid steel machine screw spin height adjustment
- 19.5 Inch diameter base
- Capacity: 250 pounds
- Premium, stain-resistant, woven, knit backed, vinyl upholstery
- 23 standard color choices
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
ST1600F	Stool Examination Four-Leg with Foot Ring 14 Inch Diameter 250 Pound Capacity Chrome	1/EA

EXAM STOOLS

Value Exam Stools



- One-piece, American-made
- Nylon base
- Full size diameter (16 Inch) top and base for comfort and stability
- 5-leg design for greater stability
- Soft-feel poly foam padding
- 4 Inch thick seat for extra comfort
- 2.5 Inch dual wheel hooded casters
- Capacity: 250 pounds
- Easy, no tools required assembly
- Premium, stain-resistant, woven knit-backed, vinyl upholstery
- 23 standard color choices
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
----------	-------------	----------

Stool Examination Value 4 Inch x 16 Inch Round Poly Foam Seat Specify Upholstery 5 Leg 250 Pound Capacity

ST1700	Stool Examination Five-Leg Pneumatic Lift Nylon Base 250 Pound Capacity 16 Inch Diameter 19 Inch to 24.5 Inch Height Range	1/EA
ST1700B	Stool Examination Five-Leg Pneumatic Lift with Contour Backrest Nylon Base 250 Pound Capacity 16 Inch Diameter 19 Inch to 24.5 Inch Height Range	1/EA
ST1700BF	Stool Examination Five-Leg Spin Lift with Foot Ring Fire Rated Nylon Base 250 Pound Capacity 16 Inch Diameter 18.5 Inch to 22.75 Inch Height Range	1/EA
ST1700F	Stool Examination Five-Leg Spin Lift with Foot Ring Nylon Base 250 Pound Capacity 16 Inch Diameter 18.5 Inch to 22.75 Inch Height Range	1/EA
ST1750	Stool Examination Five-Leg Pneumatic Lift Fire Rated Nylon Base 250 Pound Capacity 16 Inch Diameter 19 Inch to 24.5 Inch Height Range	1/EA
ST1750B	Stool Examination Five-Leg Pneumatic Lift with Contour Backrest And Foot Ring Fire Rated Nylon Base 250 Pound Capacity 16 Inch Diameter 19 Inch to 24.5 Inch Height Range	1/EA
ST1750BF	Stool Examination Five-Leg Pneumatic Lift with Contour Backrest And Foot Ring Fire Rated Nylon Base 250 Pound Capacity 16 Inch Diameter 19 Inch to 24.5 Inch Height Range	1/EA
ST1750F	Stool Examination Five-Leg Pneumatic Lift with Contour Backrest And Foot Ring Nylon Base 250 Pound Capacity 16 Inch Diameter 19 Inch to 24.5 Inch Height Range	1/EA

EXAM STOOLS

Deluxe Exam Stools



- 23 Inch premium, die cast, aluminum base with satin finish top
- 5-leg design for greater stability
- Full size diameter top and base for comfort and stability
- Soft-feel poly foam padding
- 4 Inch thick seat for extra comfort
- Chrome hooded dual wheel casters
- Pneumatic height adjustment—large, easy to locate D-shaped height lever
- Capacity: 250 pounds
- Easy, no tools required assembly
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- 23 standard color choices
- 5-year limited warranty



Item No.	Description	Quantity
----------	-------------	----------

Stool Examination Deluxe 4 Inch x 16 Inch Round Poly Foam Seat Specify Upholstery Color 5 Leg Aluminum Base 250 Pound Capacity

ST1800	Aluminum Base Height Adjustable From 17.75 Inch to 23.24 Inch	1/EA
ST1800B	With Backrest Aluminum Base Height Adjustable from 17.75 Inch to 23.24 Inch	1/EA
ST1800BF	With Backrest And Footring Aluminum Base Height Adjustable from 17.75 Inch to 23.24 Inch	1/EA
ST1800F	With Footring Aluminum Base Height Adjustable from 17.75 Inch to 23.24 Inch	1/EA
ST1850	Fire Rated Aluminum Base Height Adjustable from 17.75 Inch to 23.24 Inch	1/EA
ST1850B	With Backrest Fire Rated Aluminum Base Height Adjustable from 17.75 Inch to 23.24 Inch	1/EA
ST1850BF	With Backrest And Footring Fire Rated Aluminum Base Height Adjustable from 17.75 Inch to 23.24 Inch	1/EA
ST1850F	With Footring Fire Rated Aluminum Base Height Adjustable from 17.75 Inch to 23.24 Inch	1/EA

We recognize that esthetics are important, so MediChoice offers matching colors and styles in coordinating pieces for attractively and uniformly furnishing everything from offices to exam rooms.

EXAM STOOLS

Hands Free Exam Stools



- Contemporary styling
- 5-leg cast 24 Inch aluminum base
- Soft Feel poly foam padding
- Soft roll, dual wheel, hooded casters
- Pneumatic cylinder height adjustment
- Foot activated height control
- Capacity: 250 pounds
- Premium, stain resistant, woven, knit-backed vinyl upholstery
- 23 standard color choices
- Hands-free stool (ST2000HF)
 - Seat diameter - 16 Inches
 - Height Range - 19.5 Inches to 24.5 Inches
 - 4 Inch thick round seat
 - Shipping weight: 23 pounds
- Hands-free stool with armrest (ST2000HFA):
 - Seat diameter - 16 Inches
 - Height Range - 19.5 Inches to 24.5 Inches
 - 4 Inch thick round seat
 - Crescent adjustable height armrest with extra firm padding
 - Easy position, depth adjustable, ratchet style, swing away armrest
 - Armrest does not rotate around seat
 - Shipping weight: 36 pounds in 2 boxes
- Hands-free stool with backrest (ST2000HFB):
 - Seat - 18 Inches
 - Height Range - 18 Inches to 23 Inches
 - Comfortable contour seat
 - Adjustable height backrest
 - Shipping weight: 35 pounds in 2 boxes
- 5-year limited warranty

Item No.	Description	Quantity
Stool 4 Inch x 16 Inch Round Poly Foam Seat Specify Upholstery Color 5 Leg 24 Inch Aluminum Base 250 Pound Capacity		
ST2000HF	Hands-Free Foot-Activated Height Control Adjusts from 19.5 Inch to 24.5 Inch	1/EA
ST2000HFA	With Armrest Hands-Free Foot-Activated Height Control Adjusts from 19.5 Inch to 24.5 Inch	1/EA
ST2000HFB	With Backrest Hands-Free Foot-Activated Height Control Adjusts from 18 Inch to 23 Inch	1/EA

EXAM STOOLS

Nylon Base Exam Stools



ST1900B



ST1900

- One-piece American-made nylon base
- 5-leg design for greater stability
- Fast, pneumatic height adjustment
- Large, easy to locate D-Shaped height lever
- Capacity: 250 pounds
- Easy, no tools required assembly
- Premium, stain-resistant, woven, knit-backed, vinyl upholstery
- 23 standard choices
- 4 Inch thick seat for extra comfort
- Depth adjustable backrest (ST1900B-model only)
- Seat Diameter: 16 Inches
- Height Range: 18.75 Inches to 24.25 Inches
- 5-year limited warranty

Item No.	Description	Quantity
Stool Nylon Base 4 Inch x 16 Inch Round Poly Foam Seat Specify Upholstery 5 Leg 250 Pound Capacity		
ST1900	Pneumatic Lift Adjustable from 18.75 Inch to 24.25 Inch	1/EA
ST1900B	With D-Shaped Backrest Pneumatic Lift Adjustable from 18.75 Inch to 24.25 Inch	1/EA

Upholstery Colors

Exam Stools, Treatment Tables, and Blood Drawing Chairs

3AL Alabaster	3AS Allspice	3AB Aubergine	3BK Black	3BG Burgundy
3CG China Green	3CS Clamshell	3CM Cream	3DT Desert Tan	3DG Dove Gray
3EM Emerald	3GP Graphite	3MB Mulberry	3NL Neutral	3PG Purple Gray
3RB Royal Blue	3SP Sapphire	3SB State Blue	3TM Tomato	3VP Viscaya Palm
3WW Wedgewood	3YW Yellow			

Please see our full line of Exam Tables on pages 58-60 in this Section.

EXAM STOOLS

Black Base Exam Stools



- 16 Inch Diameter padded seat covered with seamless, vacuum formed upholstery for easy cleaning and durability
- Pneumatic height adjustment: 17.25 Inches to 22.5 Inches
- Rib-reinforced plastic stool base designed to provide maximum strength
- Standard extra-large twin wheel casters provide easier movement on carpets
- 12 upholstery colors available to complement all MediChoice exam tables
- Capacity: 300 pounds
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Stool Examination Specify Upholstery Color 5 Leg Wheel Caster Black Plastic Base 300 Pound Capacity		
EXAM6000	16 Inch Padded Seat	1/EA
EXAM6050	16 Inch Padded Seat Fire Rated	1/EA

Stainless Contoured Seat Exam Stools



- 14 Inch stainless steel contoured, revolving seat
- Seat adjusts with fast action 1 Inch OD screw
- Height adjustments
 - Non-slip tips: 25 Inches to 32 Inches
 - Casters: 27.5 Inches to 34 Inches
- Capacity: 400 pounds
- Welded footrest
- Choice of slip-resistant floor tips or 2 Inch ball-bearing swivel casters
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Stool Examination Contoured Seat Stainless Steel 4 Leg 400 Pound Capacity		
STOL1851	Non-Slip Floor Tips	1/EA
STOL1852	2 Inch Caster	1/EA

EXAM STOOLS

Stainless Padded Seat Exam Stools



- 14 Inch round cushioned seat available in 12 upholstery colors
- Seat adjusts with fast action 1 Inch OD screw
- Capacity: 400 pounds
- Welded footrest
- Choice of slip-resistant floor tips or 2 Inch ball-bearing swivel casters
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- Height adjustments:
 - Non-slip tips: 27.5 Inches to 34.5 Inches
 - Casters: 30 Inches to 36.5 Inches
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Stool Examination 14 Inch Round Cushioned Seat Specify Upholstery 4 Leg 400 Pound Capacity		
STOL1853	Height Adjustable from 27.5 Inch to 34.5 Inch Non-Slip Floor Tips	1/EA
STOL1854	Height Adjustable from 30 Inch to 36.5 Inch 2 Inch Caster	1/EA
STOL53FR	Fire Rated Height Adjustable from 27.5 Inch to 34.5 Inch Non-Slip Floor Tips	1/EA
STOL54FR	Fire Rated Height Adjustable from 30 Inch to 36.5 Inch 2 Inch Caster	1/EA

Exam Stool Accessories



- For Deluxe, Value, and Economy Stools
- For Deluxe, and Value Stools (Pneumatic Cylinder)
- 5-year limited warranty



Item No.	Description	Quantity
Accessory Exam Stool		
EXTALL	Extra Tall Pneumatic Cylinder	1/EA
GLIDES	Black Nylon Stool Glides	1/EA
SRCASTERS	Soft Roll Stool Casters	1/EA, 4EA/CA

MANUAL EXAM TABLES



Reversible Manual Exam Tables



- Reversible drawers eliminate the need for choosing a right or left handed table and provide flexibility of room layout
- Vacuum-formed, removable, seamless upholstered top allows for easy cleaning, long lasting durability, and helps prevent cross-contamination (Available in 12 colors)
- Easy-glide comfort stirrups have three lateral positions and are infinitely adjustable over the length of the bar
- Pull-out all steel non-slip step (7 Inch high)
- Two side and two front full extension drawers feature quiet roller bearing drawer glides
- Upholstered removable leg rest with plastic 9 Inch x 9 Inch drip pan
- Capacity: 500 pounds
- Dimensions: 27 Inch Width x 58.5 Inch Length x 31.75 Inch Height (77.5 Inch Length with footrest extended)
- Paper roll holder for 18 Inch or 21 Inch rolls
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- 5-year warranty

Item No.	Description	Quantity
Table Examination 500 Pound Capacity 27 Inch x 58.5 Inch x 31.75 Inch		
EXAM5200	Manual Reversible	1/EA
EXAM5250	Manual Reversible Fire Rated	1/EA

Exam Table Accessories

- Optional Drawer Warmers for Manual and Power Examination Tables
- Optional Electrical Outlet

Item No.	Description	Quantity
850	Warmer Drawer Optional Accessory for MediChoice Exam Tables EXAM5000 EXAM5050 EXM701 And EXM701F	1/EA
851	Warmer Drawer Optional Accessory for MediChoice Exam Tables EXAM5200 And EXAM5250	1/EA
852	Outlet Electrical For MediChoice Exam Tables EXM4010 EXM4011 EXM401 EXM401F EXM701 EXM701F EXM5015 EXM5016	1/EA

POWERED EXAM TABLES



Low Access Powered Exam Tables



- Vacuum-formed, removable, seamless upholstered top allows for easy cleaning, long lasting durability, and helps prevent cross-contamination
- Powered height adjustment: 19 Inch to 30 Inch (ADA Compliant)
- Removable stainless steel drip pan
- One touch, easy-glide stirrups with 3 lateral positions and friction lock length adjustment
- Electro-mechanical actuator system
- 8 upholstery and accent colors available
- Capacity: 650 pounds
- Optional drawer warmer
- Concealed paper roll holder (18 Inch or 21 Inch rolls)
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- Articulating back section with gas spring assist
- 10 degree elevated leg rest pad for patient comfort
- Power Hi-Lo
- Pelvic lift
- One foot function control
- Painted steel leg pull out
- Removable upholstered leg cushion
- Front storage drawer
- Dimensions: 27 Inch Width x 58.5 Inch Length x 35 Inch Height (70 Inch Length with footrest extended)
- 3-year manufacturer's warranty

Item No.	Description	Quantity
Table Examination Power Low Access 650 Pound Capacity Steel Frame 27 Inch x 58.5 Inch		
EXM401	High Low Power	1/EA
Table Examination Power Low Access 650 Pound Capacity Steel Frame 27 Inch x 58.5 Inch		
EXM401F	High Low Power Fire Rated	1/EA
Table Examination Power Low Access 650 Pound Capacity Steel Frame 27 Inch x 58.5 Inch		
EXM701	Manual Back	1/EA
Table Examination Power Low Access 650 Pound Capacity Steel Frame 27 Inch x 58.5 Inch		
EXM701F	Manual Back Fire Rated	1/EA

POWERED PODIATRY CHAIRS

Standard Powered Podiatry Chairs



- All steel construction, with adjustable, sliding, removable stainless steel debris tray
- Ultra comfort seamless top with antimicrobial guard (Available in 8 colors)
- Three position memory function foot control
- Fully adjustable leg section, and articulating arms
- Three Functions: Power Hi-Lo; Power Backrest; Power Tilt
- Clear vinyl foot section cover
- Electromechanical actuator system
- Dimensions: 26.5 Inch Width x 74 Inch Length (83 Inch including leg section)—Seat Height 24 Inches to 36 Inches
- Capacity: 375 pounds
- Fire-rated models available to comply with CA, IL, and MA requirements
- 3 year manufacturer's warranty



Item No.	Description	Quantity
Chair Podiatry 375 Pound Capacity 26.6 Inch Width x 83 Inch Length		
EXM5015F	Power with Three Function Programmable Foot Control Fire-Rated	1/EA
EXAM5015	Power with Three Function Programmable Foot Control	1/EA

Low Access Powered Podiatry Chairs



- All steel construction with adjustable, sliding, removable stainless steel debris tray
- Ultra comfort seamless top with antimicrobial guard (8 colors available)
- Three position memory function foot control
- Fully adjustable leg section, and articulating arms
- Clear vinyl foot section cover
- Electromechanical actuator system
- Dimensions 26.5 Inch Width x 74 Inch Length (83 Inch including leg section)—Seat Height 19 Inch to 35 Inch
- Capacity: 600 pounds
- Fire-rated models available to comply with CA, IL, and MA requirements
- 3 position—Power Hi-Lo; Power Backrest; Power Tilt
- 3-Year manufacturer's warranty



Item No.	Description	Quantity
Chair Podiatry 600 Pound Capacity Steel Frame 26.5 Inch x 74 Inch		
EXM5016	Power Low Access with Three Function Programmable Foot Control	1/EA
EXM5016F	Power Low Access with Three Function Programmable Foot Control Fire-Rated	1/EA

POWERED PROCEDURE CHAIRS



Low Access Powered Procedure Chairs



- All steel construction, stainless steel leg section, and removable stainless steel drain pan
- Ultra comfort, seamless top with antimicrobial guard (8 colors available)
- ADA Compliant
- Four function: Power Hi-Lo; Power Back; Power Tilt; Power Leg Section
- Electromechanical actuator system
- Four function hand control (4010 series) or four function programmable foot control (4011 series)
- Dimensions: 27 Inch Width x 85 Inch Length, Seat Height 19 Inch to 35 Inch
- Capacity: 600 pounds
- Trendelenburg position, and power seat tilt or pelvic tilt
- Two position upholstered kneeler/footrest, articulating headrest
- One-touch easy glide stirrups; 3 lateral position and friction lock length adjustment
- Concealed paper roll holder
- Fire-rated models available to comply with CA, IL, and MA requirements
- 3 year manufacturer's warranty

Item No.	Description	Quantity
Chair with Hand Control 600 Pound Capacity Steel Frame 27 Inch x 85 Inch		
EXM4010	Low Access Four Function Power Procedure with Hand Control	1/EA
EXM4010F	Low Access Four Function Power Procedure Fire Rated with Hand Control	1/EA
EXM4011	Low Access Four Function Power Procedure with Foot Control	1/EA
EXM4011F	Low Access Four Function Power Procedure Fire Rated with Foot Control	1/EA

POWERED PROCEDURE CHAIRS



Standard Powered Procedure Chairs



- All steel construction, with stainless steel leg section, and removable stainless drain pan
- Ultra comfort, seamless top with antimicrobial guard
- Easy entry patient access
- Electromechanical actuator system
- Four function hand control (standard) or four function programmable foot control (optional)
- Dimensions: 26.5 Inch Width × 85 Inch Length; Seat Height: 24 Inch to 36 Inch
- Capacity: 400 pound
- Trendelenburg position, and power seat tilt or pelvic tilt
- 2-position upholstered kneeler/footrest
- Articulating head rest
- One touch easy glide stirrups; 3 lateral positions and friction lock length adjustment
- Concealed paper roll holder
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- 3-year manufacturer's warranty

Item No.	Description	Quantity
Chair 400 Pound Capacity Stainless Steel Frame 26.5 Inch × 85 Inch		
EXAM8000	Powered Procedure	1/EA
EXAM8050	Powered Procedure Fire Rated	1/EA

Upholstery Colors*

Exam Tables, Matching Stools and Side Chairs



Adobe Midnight Blue Mocha Onyx



River Rock Sage Sand Grey Steel Blue

BARIATRICS

Bariatric Exam Tables



- Heavy-duty steel and aluminum construction
- Vacuum-formed, seamless upholstered top allows for easy cleaning, long lasting durability, and helps prevent cross-contamination (Available in 12 colors)
- Power high-low mechanism and power back allows for effortless patient positioning
- Adjustable height from 19 Inch to 35 Inch for easy patient or wheelchair accessibility
- Two function hand control
- Capacity: 800 pounds
- Includes paper holder
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- Dimensions: 30 Inch Width × 70 Inch Length (80 Inch with Leg Section Extended)
- 3-year manufacturer's warranty



Item No.	Description	Quantity
Table Examination 800 Pound Capacity 30 Inch × 70 Inch		
EXAM7000	Bariatric Powered	1/EA
EXAM7050	Bariatric Powered Fire Rated	1/EA

Side Chairs



- Premium Naugahyde upholstery available in 8 colors to match MediChoice Exam Tables and Stools
- High density foam padding for added comfort
- Heavy-duty tubular steel construction supports up to 300 pounds
- Wall Saver leg design prevents marring of office walls, and are stackable when not in use
- Available with or without armrests
- Easy to Assemble
- Overall Dimensions: 32 Inch Height × 32 Inch Width × 23.5 Inch Depth
- 3 Year Limited Manufacturer's Warranty



Item No.	Description	Quantity
Chair Side Wall Saver, 300 Pound Capacity 32 Inch × 32 Inch × 23.5 Inch		
SIDE1225	Without Armrest Steel	1/EA
SIDE1226	With Armrest Steel	1/EA

MOBILE WORKSTATIONS

Mobile Workstations

- Ergonomic Kidney Shape with Edge Guard (MWS7400)
 - Work Surface Dimensions: 18 Inch Depth (45 centimeters) × 28 Inch Width (70 centimeters)
 - Height Range: 31 Inch (77.5 centimeters) to 40.5 Inches (101.3 centimeters)
 - Base Dimensions: 24 Inch (60 centimeter) diameter
- Rectangular Shape with Edge Guard (MWS7405)
 - Work Surface Dimensions: 15 Inch Depth (37.5 centimeter) × 26 Inch Width (65 centimeter)
 - Height Range: 31 Inches (77.5 centimeters) to 40.5 Inches (101.3 centimeters)
 - Base Dimensions: 24 Inch (60 centimeter) diameter
- Sturdy overall construction
- 24 Inch diameter base for stability & balance
- Seamless work surface with built-in edge guard
- Pneumatic lift/foot ring operation for height adjustment
- 5 heavy duty casters
- Color Options: Cherry or Maple
- 2 year limited warranty
- Optional Accessories Available



Item No.	Description	Quantity
Workstation Mobile with Edge Guard Available in Maple or Cherry Diameter 2 Inch Wheel Casters 15 Pound Weight Capacity		
MWS7400	Ergonomic Kidney Shape 18 Inch Deep × 28 Inch Width 24 Inch	1/EA
MWS7405	Rectangle Shape 15 Inch Deep × 26 Inch Width 24 Inch	1/EA

Mobile Workstations Accessories

- Optional accessories for Mobile Work Stations
- Only works with MWS7400 and MWS7405



Item No.	Description	Quantity
Accessories Mobile Workstation Fits MediChoice MWS7400 and MWS7405		
MWS180	iPad Tablet Mount Black	1/EA
MWS181	AC Adapter	1/EA
MWS182	Concentric Basket Black Wire	1/EA
MWS183	Center Basket Black Wire	1/EA
MWS184	Locking Cable for Laptop	1/EA
MWS185	Kurly Extension Cord	1/EA

Laminate Colors*
Mobile Workstations



Laminate-7054
Wild Cherry



Laminate-7909
Fusion Maple

MEDICAL CARTS

Pediatric Color-Coded Emergency Carts - Aluminum

NOT MADE WITH natural rubber latex **diverse supplier**

- Featherweight aluminum construction with central lever-locking system for use with breakaway security locks
- Nine drawers with large stylish drawer pulls and metal drawer dividers included for each (Drawer Configuration (9) 3 Inch)
- Includes dual pull out side shelves, and cutaway front plastic top
- Antimicrobial surfaces
- Dimensions: 42 Inch Height × 29 Inch Width × 24.5 Inch Depth
- 5-year manufacturer's warranty



Item No.	Description	Quantity
ACTGK9P	Cart Emergency Pediatric Color-Coded Featherweight Aluminum Central Lever Lock 9 Drawers Antimicrobial Surface 42 Inch × 29 Inch × 24.5 Inch	1/EA

Emergency Carts - Steel

made in the USA NOT MADE WITH natural rubber latex

- Stainless painted steel breakaway locking system with 100 numbered seals included
- Oxygen tank bracket, suction unit shelf, and other cart-mounted accessories that hold critical supplies and equipment are safely mounted within extended bumper
- Unibody construction: one steel sheet for back and side panels with separate undercarriage that extends lifespan 10-15 years
- Painted steel top rail/push handle
- Removable/replaceable plastic top
- Polished chrome drawer pull handles
- Full extension ball-bearing drawer slides
- 5 Inch ball-bearing casters, all full swivel-one directional, one with brake
- Dimensions: 44.5 Inch Height × 47.5 Inch Width × 22 Inch Depth with IV Pole 85 Inch Height
- 12-year manufacturer's warranty



Item No.	Description	Quantity
Cart Emergency Painted Steel Top Rail Push Handle Breakaway Lock 5 Inch Ball-Bearing Premium Caster Red		
MC633	5 Drawers 38.5 Inch × 38 Inch × 22 Inch	1/EA
MC6401	6 Drawers 44.5 Inch × 38 Inch × 22 Inch	1/EA

Accessories

Each cart includes the following:

- Stainless steel 2-prong adjustable I.V. pole
- Plastic cardiac board and brackets
- Drop shelf
- 3 - 6 compartment drawer organizers
- 3 - 7 compartment drawer organizers
- 4 utility hooks
- Oxygen tank holder
- Stainless steel top rail/push handle

MEDICAL CARTS

Deluxe Emergency Carts - Steel



- Stainless painted steel breakaway locking system with 100 numbered seals included
- Oxygen tank bracket, suction unit shelf, and other cart-mounted accessories that hold critical supplies and equipment are safely mounted within extended bumper
- Unibody construction: one steel sheet for back and side panels with separate undercarriage that extends lifespan 10-15 years
- Painted steel top rail/push handle
- Removable/replaceable plastic top
- Polished chrome drawer pull handles
- Optional Hospital grade 6 socket electrical outlet with 15 foot cord available (MC80404A)
- Full extension ball-bearing drawer slides
- 5 Inch ball-bearing casters, all full swivel-one directional, one with brake
- Dimensions: 44.5 Inch Height x 47.5 Inch Width x 22 Inch Depth with IV Pole: 85 Inch Height
- 12 year manufacturer's warranty



Item No.	Description	Quantity
MC6411	Cart Emergency Deluxe Painted Steel Top Rail Push Handle Breakaway Lock 5 Inch Ball-Bearing Premium Caster 6 Drawers 44.5 Inch x 47.5 Inch x 22 Inch Red	1/EA

Accessories

Each cart includes the following:

- Over-Head Shelving
- Stainless Steel Defibrillator Shelf
- Hospital Grade Electrical Outlet with 6 Socket 15 Foot Cord and Cord Wrap
- Extendable IV Pole (extends to 85" high)
- Suction Unit Holder/Shelf
- Cardiac Board (with mounting brackets)
- Oxygen Tank Holder
- Chart Holder (steel)
- 2 Large and 1 Small Utility Baskets
- Locking 5 Quart Sharps Container with Glove Box Holder
- (3) Utility Hooks
- Pack of 50 Break-Away Seals

MEDICAL CARTS

Emergency Carts - Aluminum



- Featherweight aluminum construction with central lever lock for use with breakaway security locks
- Six drawers with large stylish drawer pulls and metal drawer dividers included for each (Drawer Configuration (4) 3 Inch, (1) 6 Inch, (1) 9 Inch)
- Includes dual pull-out side shelves, and cutaway front plastic top
- Also includes cardiac board with mounting bracket, IV pole, and oxygen tank holder
- Antimicrobial surfaces
- Dimensions: 42 Inch Height x 29 Inch Width x 24.5 Inch Depth
- 5-year manufacturer's warranty



Item No.	Description	Quantity
ACTRL6RA	Cart Emergency Featherweight Aluminum Central Lever Lock 6 Drawers Antimicrobial Surface 42 Inch x 29 Inch x 24.5 Inch Red	1/EA

Classic Isolation Carts - Steel



- 3 storage drawers on ball-bearing slides: (1) 6 Inch and (2) 12 Inch
- Each drawer can hold up to 100 pounds
- Four extension ball-bearing drawer slides
- Four 5 Inch ball-bearing premium casters, all full swivel (one tracking, one braking)
- Key lock with 2 keys provided
- Ships in standard yellow, but other colors are available
- Unibody construction: one steel sheet for back and side panels with separate undercarriage that extends lifespan 10-15 years
- Painted steel top rail/push handle
- Removable/replaceable plastic top
- 3 polished chrome drawer pull handles
- Dimensions: 41.5 Inch Height x 32 Inch Width x 22 Inch Depth, coat rack increases height to 62.5 Inch
- 12-year manufacturer's warranty



Item No.	Description	Quantity
MC651	Cart Isolation Painted Steel 3 Drawers Key Lock 41.5 Inch x 32 Inch x 22 Inch Yellow	1/EA

Accessories

Each cart includes the following:

- Coat and hat rack
- Glove box holder
- 4 utility hooks
- Stainless steel top rail/push handle

MEDICAL CARTS

Classic Isolation Carts - Aluminum

NOT MADE WITH natural rubber latex **diverse supplier**

- Featherweight aluminum construction with central key-lock system
- Three drawers with large, stylish drawer pulls and metal drawer dividers included for each (Drawer Configuration (1) 3 Inch, (2) 9 Inch)
- Dual pull out side shelves with cutaway front plastic top
- Antimicrobial surfaces
- Dimensions: 36 Inch Height x 29 Inch Width x 24.5 Inch Depth
- 5-year manufacturer's warranty



Item No.	Description	Quantity
ACSGK3Y	Cart Isolation Featherweight Aluminum 3 Drawers Key Lock 36 Inch x 29 Inch x 24.5 Inch Yellow	1/EA

Mini-Line Isolation Carts - Steel

made in the USA **NOT MADE WITH** natural rubber latex

- 4 storage drawers on ball-bearing slides: (1)-3 Inch and (3)-9 Inch
- Each drawer can hold up to 100 pounds
- Pre-threaded accessory mounted holes
- Key lock with 2 keys provided
- Narrow design with small footprint
- Sturdy steel construction. Available with bumper and caster upgrade.
- Removable/replaceable plastic top
- 4 polished chrome pull handles
- Full extension ball-bearing drawer slides
- 4 single-wheel 3 Inch stem caster set, 2-brake, all fully swivel
- Dimensions: 40 Inch Height x 18 Inch Width x 18 Inch Depth
- 12-year manufacturer's warranty



Item No.	Description	Quantity
MC3154	Cart Isolation Mini Line Painted Steel 4 Drawers 40 Inch x 18 Inch x 18 Inch Yellow	1/EA

MEDICAL CARTS

Anesthesia Carts - Steel

made in the USA **NOT MADE WITH** natural rubber latex

- 6 storage drawers on ball-bearing slides: (4) 3 Inch, (1) 6 Inch, and (1) 12 Inch
- Each drawer can hold up to 100 pounds
- Removable/replaceable plastic top
- 6 chrome drawer pull handles
- Full extension ball-bearing drawer slides
- Four 5 Inch ball-bearing premium casters, all full swivel (one tracking, one braking)
- Available with key locking or keyless entry (push button mechanical combination lock with key lock over-ride that does not require batteries)
- Powder-coat paint for long-lasting, durable finish. Standard color is blue for anesthesia carts, but other colors are available.
- Raised back rail to hold accessories
- Unibody construction: one steel sheet for back and side panels with separate undercarriage that extends lifespan 10-15 years
- 2 aluminum side rails to mount accessories i.e. glove box holder or waste receptacle
- Painted steel top rail/push handle
- Dimensions: 41.5 Inch Height x 32 Inch Width x 22 Inch Depth (55.9 Inch Height with back rail)
- 12-year manufacturer's warranty



Item No.	Description	Quantity
Cart Anesthesia Steel with Aluminum Side Rails		
MC645	6 Drawers Key Lock 32 Inch x 22 Inch Blue	1/EA
MC745E	6 Drawers Upgraded Electronic Lock 32 x 22 Inch Blue	1/EA

Accessories

Each cart includes the following:		3-7 compartment drawer organizers
3-6 compartment drawer organizers	Tape and label dispenser	
5 Compartment tilt bins	Glove box holder	
Drop shelf	2 utility hooks	
Waste receptacle		
Stainless steel top rail/push handle		

MEDICAL CARTS

Anesthesia Carts - Aluminum

NOT MADE WITH
natural rubber latex

- Featherweight aluminum construction with central key lock system
- Six drawers with large stylish drawer pulls and metal drawer dividers included for each (Drawer Configuration (4) 3 Inch, (1) 6 Inch, (1) 9 Inch)
- Includes dual pull-out side shelves, and cutaway front plastic top
- Also comes with waste receptacle, tape dispenser, and base rail system for mounting four compartment tilt bin (also included)
- Antimicrobial surfaces
- Dimensions: 42 Inch Height × 29 Inch Width × 24.5 Inch Depth
- 5 year manufacturer's warranty



Item No.	Description	Quantity
ACTGK6DA	Cart Anesthesia Featherweight Aluminum 6 Drawers Key Lock 42 Inch × 29 Inch × 24.5 Inch Dark Blue	1/EA

Treatment Carts

NOT MADE WITH
natural rubber latex

diverse
supplier

- 2 storage drawers: 3 Inch and 6 Inch to hold compact treatments and 21 Inch storage compartment with locking doors and adjustable shelf to store bulky items
- Each drawer can hold up to 100 pounds
- Four 5 Inch ball-bearing premium casters, all full swivel (one tracking, one braking)
- Key locks on drawers and doors, with 4 keys provided
- Full wrap-around vinyl bumper extends wheel base to increase stability
- Unibody construction: one steel sheet for back and side panels with separate undercarriage that extends lifespan 10-15 years
- Stainless steel top rail/push handle
- Removable/replaceable plastic top
- 4 chrome drawer/door pull handles
- Full extension ball-bearing drawer slides
- Dimensions: 44.5 Inch Height × 32 Inch Width × 22 Inch Depth (does not include w/accessories)



Item No.	Description	Quantity
MC620	Cart Treatment 2 Drawers Key Lock 44.5 Inch × 32 Inch × 22 Inch Stainless Steel Beige	1/EA

Accessories

Each cart includes the following:

- 2 - 3" Adjustable exchange trays
- Drop shelf
- Glove box holder
- Sharps container bracket
- Stainless steel top rail/push handle

MEDICAL CARTS



Monitor Carts - Steel

made in the
USA

NOT MADE WITH
natural rubber latex

- Includes tall over-head shelf system with 2 shelves offering multiple mounting locations
- 6 storage drawers on full ball-bearing sides: (3)-3 Inch, (2)-6 Inch, and (1)-9 Inch
- Unibody construction: one steel sheet for back and side panels with separate undercarriage that extends lifespan 10-15 years
- Removable/replaceable plastic top
- 6 chrome drawer pull handles
- 5 Inch ball-bearing casters, all full swivel—one directional, one with brake
- Durable, hard-baked powder-coat finish
- Optional hospital grade 6 socket electrical outlet with 15 foot cord available (MC80404A)
- Overall cart dimensions: 75 Inch Height × 40 Inch Width × 22 Inch Depth
- 12 year manufacturer's warranty

Item No.	Description	Quantity
MC6651	Cart Monitor Multi-Purpose 6 Drawer with Key Lock 75 Inch × 40 Inch × 22 Inch Stainless Steel Gray	1/EA

MEDICAL CARTS

Cast Carts - Steel

made in the USA NOT MADE WITH natural rubber latex

- 8 storage drawers on ball bearing slides: (1) 3 Inch top drawer with black vinyl mat, (5) 3 Inch drawers with 4 adjustable dividers per drawer, (2) 6 Inch drawers with 4 adjustable dividers per drawer
- Each drawer can hold up to 100 pounds
- Unibody, stainless steel construction: one steel sheet for back and side panels, with separate undercarriage that extends the cart lifespan 10-15 years
- Vinyl padded wrap-around bumper with extended width to accommodate stainless side storage cabinet and stainless steel bucket holder
- Optional top storage compartment accommodates casting supplies
- Four 5 Inch ball-bearing premium casters, all full swivel (one tracking, one braking)
- Carts are ready to put into service with a minimum amount of accessory assembly, delivered with wheels on
- Dimensions: 43 Inch Height x 44 Inch Width x 21 Inch Depth (66 Inch Height with top compartment, 76 Inch Height with pole)
- 12-year manufacturer's warranty



3023MC6026



3023MC6025

Item No.	Description	Quantity
MC6025	8 Drawers 43 Inch x 44 Inch x 21 Inch	1/EA
MC6026	8 Drawers Top Storage 66 Inch x 44 Inch x 21 Inch	1/EA

Accessories

Each cart includes the following:

- Stainless folding side shelf
- Stainless side storage cabinet
- Swivel tourniquet/pump pole
- Stainless steel bucket holder with black vinyl mat
- Aluminum towel bar
- 24" stainless ruler
- Scissors cable

MEDICAL CARTS

Multi-Purpose Phlebotomy Carts - Hybrid Aluminum

NOT MADE WITH natural rubber latex diverse supplier

- Featherweight hybrid aluminum construction with Gatelock key locking system for security
- Four half-width drawers (one each of 3 Inch, 4 Inch, 6 Inch, and 9 Inch) all include removable trays and dividers
- One full width 6 Inch deep drawer
- Also includes: 12 clearview tilt bins, secured sharps container bracket, glove box and paper holders, tilt out waste receptacle, and removable plastic top
- Dimensions: 45 Inch Height x 27 Inch Width x 24 Inch Depth
- 5-year manufacturer's warranty



Item No.	Description	Quantity
ACTWK5W	Cart Medical Phlebotomy Multipurpose Featherweight Aluminum Gatelock Key Lock 5 Drawers 45 Inch x 27 Inch x 24 Inch White	1/EA

Cylinder Carts-Patient Use

NOT MADE WITH natural rubber latex diverse supplier

- Ergonomic, adjustable height, "cane style" handle promotes patient comfort
- Swept-back axle for greater stability
- Chrome finish
- Mounting hook hangs over headboards or fits into IV pole bed mounts
- Accessory hook for mask and extra tubing
- Some assembly required



Item No.	Description	Quantity
MC0XY1	Single 30 Inch x 12 Inch x 8 Inch	1/EA, 2/CA
MC0XY2	Double 30 Inch x 14 Inch x 12 Inch	1/EA, 2/CA

Cylinder Carts-Bulk Transport

NOT MADE WITH natural rubber latex diverse supplier

- Safe transport and storage solutions for multiple E or D medical gas cylinders
- Durable dual-bar reinforced steel construction for years of service
- Clean chrome finish
- Standard 5 Inch heavy-duty casters, 2 locking
- Some assembly required



Item No.	Description	Quantity
MC0XY12	Fits 12 E D Cylinders 39 Inch x 21 Inch x 16 Inch	1/EA
MC0XY6	Fits 6 E D Cylinders 39 Inch x 16 Inch x 11 Inch	1/EA

MEDICAL STORAGE CARTS

Column Carts

- Design incorporates a modular tray system with maximum variability in tray configuration
- 30 percent space savings through flexible dividing of trays and baskets
- Smooth walls for hygienic cleaning
- Simple and Quick divider adjustment
- Panel spacing of 1.25 Inch gives higher density
- Flexible sliding dividers
- Available in: Single, double or triple columns
- Available in: Tambour or Tempered Glass doors
- Locking System Options: Key or Electronic




Item No.	Description	Quantity
Cart Medical Storage Powder-Coated Steel		
Single Column		
MS8128	Single Column Key-Locking Tambour Door 80.5 Inch x 27.8 Inch x 27.9 Inch	1/EA
MS8128DR	Single Column Key-Locking Tempered Glass Door 76 Inch x 27.8 Inch x 29.6 Inch	1/EA
MS8128DREKC	Single Column Electronic-Locking Tempered Glass Door 76 Inch x 27.8 Inch x 29.6 Inch	1/EA
Double Column		
MS8140	Double Column Key-Locking Tambour Door 80.5 Inch x 42 Inch x 27.9 Inch	1/EA
MS8140DR	Double Column Key-Locking Tempered Glass Door 76 Inch x 42 Inch x 27.9 Inch	1/EA
MS8140DREKC	Double Column Electronic-Locking Tempered Glass Door 76 Inch x 42 Inch x 27.9 Inch	1/EA
MS8140SH	Double Column Key-Locking Tambour Door 5 Adjustable Shelves 80.5 Inch x 42 Inch x 27.9 Inch	1/EA
MS8140SPLITB	Double Column Key-Locking Tambour Door Split Column Below and Open Above 80.5 Inch x 42 Inch x 27.9 Inch	1/EA
Triple Column		
MS6140	Three-Quarter Height Double Column Key-Locking Tambour Door 60.5 Inch x 42 Inch x 27.9 Inch	1/EA
MS6140DR	Three-Quarter Height Double Column Key-Locking Tempered Glass Door 57.5 Inch x 42 Inch x 29.6 Inch	1/EA
MS8160	Triple Column Key-Locking Tambour Door 80.5 Inch x 59.6 Inch x 27.9 Inch	1/EA
MS8160R	Triple Column Key-Locking Tambour Door Double-column Shell on Left Single-column on Right 80.5 Inch x 59.6 Inch x 27.9 Inch	1/EA
MS8160LSPLITB	Triple Column Key-Locking Tambour Door Single-column Shell on Left Double-column on Right with Split Double Column Below and Open Above 80.5 Inch x 59.6 Inch x 27.9 Inch	1/EA

MEDICAL STORAGE CARTS

Mobile Medical Storage Cart Accessories

- Accessories used with Mobile Medical Storage Carts Only
- Light weight trays that are highly durable
- Simple and quick divider adjustment
- Accessories help to customize storage cart to users' needs



Item No.	Description	Quantity
MSCATH7	Shelf for Double Wide Medical Storage Cart without Center Panel 7 Catheter Slides 5 Hooks Per Slide	1/EA
MSSH2	Shelf for Full Double Column Medical Storage Cart 5 Adjustable Dividers	1/EA
Shelf for Single Column Medical Storage Cart		
MSCATH3	3 Catheter Slides 5 Hooks Per Slide	1/EA
MSSH1	2 Adjustable Dividers	1/EA
Divider for Medical Storage Cart		
81033	2 Inch Exchange Tray ABS Gray 23.5 Inch x 2 Inch	1/EA
81034	2 Inch Exchange Tray Cross ABS Gray 15.75 Inch x 2 Inch	1/EA
81035	4 Inch Exchange Tray Length ABS Gray 23.625 Inch x 4 Inch	1/EA
81036	4 Inch Exchange Tray Cross ABS Gray 15.75 Inch x 4 Inch	1/EA
81037	8 Inch Exchange Tray ABS Gray 23.5 Inch x 7.75 Inch	1/EA
81038	8 Inch Exchange Tray Cross ABS Gray 15.75 x 8 Inch	1/EA
Exchange Tray for Medical Storage Cart		
81031S	4 Inch Flex Module ABS Gray Pull-out Stoppers 23.625 Inch x 15.75 Inch x 4 Inch	1/EA
81032S	8 Inch Flex Module ABS Gray Pull-out Stoppers 23.625 Inch x 15.75 Inch x 7.875 Inch	1/EA
Divider Accessory Mobile Storage Cart Carbon Steel Wire Powder Painted White		
81073	3 Inch Short Divider Fits Cross-wise	1/EA
81074	3 Inch Long Divider Fits Lengthwise	1/EA
81075	5 Inch Short Divider Fits Cross-wise	1/EA
81076	5 Inch Long Divider Fits Lengthwise	1/EA
81077	8 Inch Short Divider Fits Cross-wise	1/EA
81078	8 Inch Long Divider Fits Lengthwise	1/EA
Basket Accessory Mobile Storage Cart Powder Painted White		
81070	3 Inch Carbon Steel Wire 23.62 Inches x 15.75 Inches x 3 Inches	1/EA
81071	5 Inch Carbon Steel Wire 23.62 Inch x 15.75 Inch x 5 Inch	1/EA
81072	8 Inch Carbon Steel Wire 23.62 Inches x 15.75 Inches x 8 Inches	1/EA
810701	3 Inch Carbon Steel Wire with 1 Short Divider 23.62 Inches x 15.75 Inches x 3 Inches	1/EA
810702	3 Inch Carbon Steel Wire with 1 Long Divider 23.62 Inches x 15.75 Inches x 3 Inches	1/EA
810703	3 Inch Carbon Steel Wire with 1 Long and 1 Short Divider	1/EA
810711	5 Inch Carbon Steel Wire with 1 Short Divider 23.62 Inches x 15.75 Inches x 5 Inches	1/EA
810712	5 Inch Carbon Steel Wire with 1 Long Divider	1/EA
810713	5 Inch Carbon Steel Wire with 1 Long and 1 Short Divider	1/EA
810721	8 Inch Carbon Steel Wire with 1 Short Divider 23.62 Inches x 15.75 Inches x 8 Inches	1/EA
810722	8 Inch Carbon Steel Wire with 1 Long Divider 23.62 Inches x 15.75 Inches x 8 Inches	1/EA
810723	8 Inch Carbon Steel Wire with 1 Long and 1 Short Divider	1/EA

MEDICAL STORAGE CARTS

Catheter Procedure Carts

- Cath Lab in a Cart: Designed especially for Cath Labs and IR department
- Holds all sizes of catheters on full extension holders
- All-inclusive cart that keeps the necessary supplies together in one place
- Cath slides feature a pair of ball bearing slides to withstand the pull force of removing catheters from hooks
- Cath slide shelves can be mounted in a lower position to create storage above for over-sized items
- Key locking



ODRCATH



8128CATH



8140CATH



8DRCATH

Item No.	Description	Quantity
Cart Catheter Procedure Powder-Coated Steel		
Single Column		
MS8128CATH	Single Column with Catheter Slide Shelf 3 Catheter Slides 5 Hooks Per Slide with Tambour Door 80.5 Inch x 27.8 Inch x 27.9 Inch	1/EA
MS8128DRCATH	Single Column with Catheter Slide Shelf 3 Catheter Slides 5 Hooks Per Slide with Tempered Glass Door 76 Inch x 27.8 Inch x 29.6 Inch	1/EA
Double Column		
MS8140CATH	Double Wide with Full Width Catheter Slide Shelf 7 Catheter Slides 5 Hooks Per Slide with Tambour Door 80.5 Inch x 42 Inch x 27.9 Inch	1/EA
MS8140DRCATH	Double Wide with Full Width Catheter Slide Shelf 7 Catheter Slides 5 Hooks Per Slide with Tempered Glass Door 76 Inch x 42 Inch x 29.6 Inch	1/EA

OR AND STRETCHER PADS

Stretcher Replacement Pads



- Various sizes of black pads available to fit all stretcher models
- Covers are radiolucent, waterproof, stain resistant, anti-microbial, anti-bacterial, and washable
- Fire-rated foam fillers conform to CA Technical Bulletin 117
- Standard Pad is constructed of high-density polyurethane foam with a durable, heavy duty Naugahyde cover
- Deluxe Pad is constructed of high-density polyurethane foam with a soft stretchable, urethane-coated, conductive top cover and vinyl sheeting bottom cover
- Additional brand configurations available upon request

Item No.	Description	Quantity
Pad Stretcher for Hausted 31625 2 Piece Black High Density Polyurethane Foam		
MC31623L	Deluxe 2-Way Stretch Cover 3 Inch	1/EA
MC31624L	Deluxe 2-Way Stretch Cover 4 Inch	1/EA
MC31625N	Standard Naugahyde Cover 3 Inch	1/EA
Pad Stretcher for Hausted 31526 Black High Density Polyurethane Foam		
MC31626N	Standard Naugahyde Cover 4 Inch 2 Piece Set	1/EA
MC31663L	Deluxe 2-Way Stretch Cover 29.25 Inch x 76 Inch x 3 Inch	1/EA
Pad Stretcher for Hausted 31661 Black High Density Polyurethane Foam		
MC31663N	Standard Naugahyde Cover 29.5 Inch x 76 Inch x 3 Inch	1/EA
MC31664L	Deluxe 2-Way Stretch Cover 29.25 Inch x 76 Inch x 3 Inch	1/EA
MC31664N	Standard Naugahyde Cover 29.25 Inch x 76 Inch x 3 Inch	1/EA
Pad Stretcher for Hausted 31680 Black High Density Polyurethane Foam		
MC31680N	Standard Naugahyde Cover 24.5 Inch x 72 Inch x 3 Inch	1/EA
MC31689L	Deluxe 2-Way Stretch Cover 24.5 Inch x 72 Inch x 3 Inch	1/EA
Pad Stretcher for Hausted		
MC28764D	Gemini 67177 Deluxe 2-Way Stretch Cover 28 Inch x 76 Inch x 4 Inch	1/EA
MC2469L	24 x 69 x 3 Inch	1/EA
MC75112	75112 Naugahyde Cover 24.5 Inch x 70.5 Inch x 3 Inch	1/EA
MC1006H	Head Piece	1/EA
MC25754	Custom Lectrolite Cover Seamless Foam 4 Inch x 25 Inch x 75 Inch	1/EA
Pad Stretcher for Midmark 900 Black High Density Polyurethane Foam		
MC21921L	Deluxe 2-Way Stretch Cover 24.5 Inch x 71.25 Inch x 2 Inch	1/EA
MC21921N	Naugahyde Cover 24.5 Inch x 71.25 Inch x 2 Inch	1/EA
MC21931L	Deluxe 2-Way Stretch Cover 24.5 Inch x 71.25 Inch x 3 Inch	1/EA
MC21931N	Naugahyde Cover 24.5 Inch x 71.25 Inch x 3 Inch	1/EA
MC21941L	Deluxe 2-Way Stretch Cover 24.5 Inch x 71.25 Inch x 4 Inch	1/EA
MC21941N	Standard Naugahyde Cover 24.5 Inch x 71.25 Inch x 4 Inch	1/EA
Pad Stretcher for Midmark Models 900 1800 1900 Black High Density Polyurethane Foam 29.5 Inch x 77 Inch x 3 Inch		
MC2893L	Deluxe 2-Way Stretch Cover	1/EA
MC2893N	Standard Naugahyde Cover	1/EA
Pad Stretcher for Midmark 540 Black High Density Polyurethane Foam 27 Inch x 77 Inch x 3 Inch		
MC54003	Deluxe 2-Way Stretch Cover	1/EA
Pad Stretcher Black High Density Polyurethane Foam 24 to 27 Inch x 74 to 77 Inch x 3 to 4 Inch		
MC27743L	Deluxe 2-Way Stretch Cover	1/EA
MC27743N	Standard Naugahyde Cover	1/EA

continued on next page

OR AND STETCHER PADS

continued from previous page

Item No.	Description	Quantity
Pad Stretcher Universal Style Black High Density Polyurethane Foam 24 Inch x 76 Inch x 3 Inch		
MC24763L	Deluxe 2-Way Stretch Material	1/EA
MC24763N	Naugahyde Cover	1/EA
Pad Stretcher for Hill-Rom 820 Black High Density Polyurethane Foam 34 Inch x 79.5 Inch x 6 Inch		
MC8206L	Deluxe 2-Way Stretch Cover	1/EA
MC8206N	Standard Naugahyde Cover	1/EA
Pad Stretcher for Hill-Rom 881 Black High Density Polyurethane Foam 27 Inch x 76 Inch x 3 Inch		
MC8813L	Deluxe 2-Way Stretch Cover	1/EA
MC8813N	Standard Naugahyde Cover	1/EA
Pad Stretcher for Hill-Rom 883 Black High Density Polyurethane Foam 27 Inch x 78 Inch x 3 Inch		
MC8833L	Deluxe 2-Way Stretch Cover	1/EA
MC8833N	Standard Naugahyde Cover	1/EA
Pad Stretcher for Hill-Rom		
MC8823	8823 78 x 27 x 4	1/EA
MCP80004	8000 Standard Vinyl Cover Seamless with Cut Sides 4 Inch x 25 Inch x 76 Inch Foam	1/EA
Pad Stretcher for Stryker 946-3 Black High Density Polyurethane Foam 25 Inch x 76 Inch x 3 Inch		
MC9462L	Deluxe 2-Way Stretch Cover	1/EA
MC9463N	Standard Naugahyde Cover	1/EA
Pad Stretcher for Stryker 1060 2 Piece Black High Density Polyurethane Foam 3 Inch		
MC1060L	Deluxe 2-Way Stretch	1/EA
MC1060N	Naugahyde Cover Standard	1/EA
Pad Stretcher for Stryker 1005		
MC1005CC	Comfort Cover 26 Inch x 76 Inch x 3 Inch	1/EA
MC1005N	Black Naugahyde Cover 30 Inch x 76 Inch x 4 Inch Standard High Density Polyurethane Foam	1/EA
Pad Stretcher for Stryker		
MC17103N	Renaissance Stretcher Naugahyde Cover 28 Inch x 75 Inch x 4 Inch	1/EA
MC1007CC	1007 Comfort Cover 30 Inch x 76 Inch x 4 Inch	1/EA
MC1061N	1061 2-Piece Naugahyde Cover	1/EA
MC26744D	Deluxe 2-Way Stretch Cover 24 Inch x 74 Inch x 4 Inch	1/EA
MC24744D	Deluxe 2-Way Stretch Cover 24 Inch x 74 Inch x 4 Inch	1/EA
MC10014N	Model 1001 2-way Stretch Vinyl Cover with Welded Seams 3 Inch x 26 Inch x 76 Inch Foam Pad	1/EA
MC72244L	Lectrolite Standard Vinyl Cover Seamless Foam 4 Inch x 24 Inch x 72 Inch	1/EA
MC1015CC	1015 Comfort Cover 30 Inch x 76 Inch x 3 Inch	1/EA
Pad Stretcher for Skytron		
MC3600DG	3600 3-Piece Antibacterial Washable Gel Pad Cover High-Density Polyurethane Foam	1/EA
MC3501BD	3501B 3-Piece Deluxe 2-Way Cover	1/EA
MC3002D	3002 1-Piece Image Deluxe Cover	1/EA
MC3003D	3003 1-Piece Image Deluxe Cover	1/EA
MC6701D	6701 3-Piece Deluxe 2-Way Stretch Cover	1/EA
Pad Stretcher		
MC30764D	Deluxe 2-Way Stretch Cover 30 Inch x 76 Inch x 4 Inch	1/EA
MC26784D	Rounded Corners Deluxe 2-Way Stretch Cover 26 Inch x 78 Inch x 4 Inch	1/EA
MC78463	For Cyto Liebel HUT-IV-HS Table Standard Vinyl Cover Seamless Foam 3 Inch x 48 Inch x 78 Inch	1/EA
MC1073	For Cardiac Cath Lab Custom Standard Vinyl Cover Seamless Memory Foam 3 Inch x 26 Inch x 107 Inch	1/EA

OR AND STETCHER PADS

Item No.	Description	Quantity
Pad Replacement Stretcher Stryker 26 x 76 x 4 Inch Polyurethane Foam Black		
MC5454N	For Pedigo 5831002 Standard Vinyl Cover Seamless Foam 4 Inch x 31 Inch x 77 Inch	1/EA
MC5453N	Custom For Midmark 545 Naugahyde Cover 3 Inch x 31 Inch x 77 Inch	1/EA
MC27744T	With A Fusion Four Way Stretch Cover Material	1/EA
MC1009N	With Deluxe Stretchable Cover Material	1/EA
MC10612N	Pad Stretcher for Stryker Gynnie 2 Piece Naugahyde Cove 4 Inch Bocy 3 Inch Foot	1/EA
MC25755N	Stretcher Pad for Stryker Naugahyde Cover 25.5 Inch x 75 Inch x 5 Inch	1/EA

OR Table Replacement Pads



- Standard 2 Inch black pad to fit most common OR table models
- Covers are radiolucent, waterproof, stain resistant, anti-microbial, anti-bacterial, and washable
- Fire-rated foam fillers conform to CA Technical Bulletin 117
- Standard Pad is constructed of high-density polyurethane foam
- Deluxe Pad is constructed of high-density polyurethane foam with a soft stretchable, urethane-coated, conductive top cover and vinyl sheeting bottom cover
- Additional brand configurations available upon request

Item No.	Description	Quantity
Pad Table Operating Room for Skytron 3500 3 Piece Black Polyurethane Foam 2 Inch		
MC3500D	Deluxe Soft Stretchable Urethane-Coated Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA
MC3500S	Standard Conductive Vinyl Sheeting Top and Bottom Cover	1/EA
Pad Table Operating Room for Skytron 3600 3 Piece Black Polyurethane Foam 2 Inch		
MC3600D	Deluxe Soft Stretchable Urethane-Coated Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA
MC3600S	Standard Conductive Vinyl Sheeting Top and Bottom Cover	1/EA
Pad Table Operating Room for Shampaine 4900 3 Piece Black Polyurethane Foam 2 Inch		
MC4900D	Deluxe Soft Stretchable Urethane-Coated Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA
MC4900S	Standard Conductive Top Cover Vinyl Sheeting Bottom Cover	1/EA
Pad Table Operating Room for Skytron 3100 2 Piece Black Polyurethane Foam 2 Inch		
MC3100D	Deluxe Soft Stretchable Urethane-Coated Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA
MC3100S	Standard Conductive Top Cover Vinyl Sheeting Bottom Cover	1/EA
Pad Table Operating Room for Amsco 3080/3085 Black Polyurethane Foam 2 Inch		
MC3080D	Deluxe Soft Stretchable Urethane-Coated Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA
MC3080S	Standard Conductive Vinyl Sheeting Top and Bottom Cover	1/EA
Pad Table Operating Room for Skytron 5001/6001/6500/6600 Black Polyurethane Foam 2 Inch		
MC5001D	Deluxe Soft Stretchable Urethane-Coated Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA
MC5001S	Standard Conductive Vinyl Sheeting Top and Bottom Cover	1/EA
Pad Table Operating Room for Amsco 1080/2080 Black Polyurethane Foam 2 Inch		
MC1080D	Deluxe Soft Stretchable Urethane-Coated Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA
MC1080S	Standard Conductive Top Cover and Conductive Vinyl Sheeting Bottom Cover	1/EA

MATTRESS PADS

Mattress Pads

- C-Matt Prevention Hospital Bed Mattress
 - Designed for pressure ulcer prevention and patient comfort
 - Best in class top cover, breathable, bi-directional stretch fabric
 - Fully radio frequency welded seams to prevent fluid ingress
 - Multi-foam core for anatomically dependent pressure redistribution
 - Shear liner for patient comfort
 - Various sizes available
- C-Matt Comfort Ultima Hospital Bed Mattress
 - The Ultimate in patient comfort
 - Conforming viscoelastic memory foam topper
 - Memory foam topper works with the patient's unique body shape to redistribute pressure
 - Breathable, bi-directional stretch top cover
 - Shear liner for patient comfort
 - Various sizes available
- C-Matt Comfort Hospital Bed Mattress
 - Entry level prevention performance
 - Multi-zoned foam core for patient comfort and safety
 - Leak proof top cover
 - Head to foot reversible
 - Various sizes available
- Additional brand configurations available upon request



Item No.	Description	Quantity
Cover Mattress Versa Care Sonically Sealed Blue		
MC11154	Foam	1/EA
MC11155	Air	1/EA
Kit Replacement Mattress Versa Care Includes Blue Cover Viscoelastic Foam Topper Shear Liner		
MC11516	Foam Ca 129	1/EA
MC11518	Air Ca 117	1/EA
Cover Replacement for Total Care Sport Mattress Pads with 14 Hoses Dark Fuscia Includes Urethane Top and Bottom Bravo Silver Technology Sonically Sealed Top and Bottom		
MC11570	P1915eds 36 Inch x 84 Inch x 7 Inch	1/EA
MC11642	Refresher Kit P1915eds Shear Liner Fire Barrier and Foam Topper 36 Inch x 84 Inch x 7 Inch	1/EA
Mattress Bed Comfort Ultima Blue Conforming Viscoelastic Memory Foam Core with Breathable Top Cover Urethane Bi-directional Stretch Radio Frequency Welded Seams with Shear Liner Ca 117		
MC36806C	36 Inch x 80 Inch x 7 Inch	1/EA
MC36846C	36 Inch x 84 Inch x 7 Inch	1/EA
Mattress Bed C-matt Pressure Ulcer Prevention Blue Multi-zone Foam Core with Breathable Top Cover Urethane Bi-directional Stretch Radio Frequency Welded Seams with Shear Liner Ca 117		
MC36806P	36 Inch x 80 Inch x 7 Inch	1/EA
MC84366P	36 Inch x 84 Inch x 7 Inch	1/EA
Mattress Bed Comfort Multi-zoned Foam Core with High Grade Nylon Leak Proof Head to Toe Reversible Blue Top Cover		
MC80366C	36 Inch x 80 Inch x 6 Inch	1/EA
MC84366C	36 Inch x 84 Inch x 7 Inch	1/EA
Mattress Bed C-Matt Pressure Ulcer Prevention Pink Urethane Breathable Top Cover Bi-directional Stretch Radio Frequency Welded Seams		
MC33584P	7 Inch x 35.5 Inch x 84 Inch	1/EA

Cot Replacement Pads

- Recovery 5 cover material
- Chem soft bottom cover
- Fire rated foam fillers conforms to CA 117
- Various sizing options available
- Additional brand configurations available upon request



Item No.	Description	Quantity
MC30705	Cot Mattress Custom 5 Inch x 30 Inch x 70 Inch 5 Soft Vinyl Top Chemsafe Bottom Material 2 Inch 1820 Foam On Top Of 3 Inch 1843 Foam	1/EA

Crib Replacement Pads

- Lectrolite white baby soft cover
- Lectrolite birthday bears cover
- High-density polyurethane foam pad
- Soft, stretchable urethane-coated top cover
- Fire rated foam fillers conforms to CA 117
- Various sizing options available
- Additional brand configurations available upon request



Item No.	Description	Quantity
Mattress Crib Bassinet with Baby Soft White Non-Woven Polyester Cover Foam Pad		
MC12261	12 x 26 x 1.5 Inch	1/EA
Pad Crib Mattress Foam Pad Baby Soft Vinyl Cover Seamless		
MC24564	4 Inch x 24 Inch x 56 Inch	1/EA
MC28583	3 Inch x 28 Inch x 58 Inch	1/EA
MC32574	4 Inch x 32.5 Inch x 57 Inch	1/EA
Pad Crib Mattress Foam Pad Birthday Bears Vinyl Cover Seamless		
MC28584BB	4 Inch x 28 Inch x 58 Inch	1/EA

X-Ray Replacement Pads

- Black conductive cover material
- Fluid proof cover seams
- Radiolucent for artifact free imaging
- Fire rated foam fillers conforms to CA 117
- Various sizing options available
- Additional brand configurations available upon request

Item No.	Description	Quantity
MC28721	X-Ray Table Pad 1 Inch x 30 Inch x 72 Inch Black Vinyl Top with Two Handles at Top And Bottom	1/EA

RECLINERS

Manual Treatment Recliners



- Three position recliner with soft radius lines for an “at home” feeling for the patient room or treatment room
- Gently curved to cradle patient
- Fixed or swing arms (depending on model)
- 3 Inch thermoplastic locking casters
- High-grade polyurethane foam
- 8 vinyl color options of Champion Soft Vinyl and an extensive line of accessory options available
- Easy to use foot release Trendelenburg is standard
- Capacity: 300 pounds
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- 10-year limited manufacturer’s warranty



Item No.	Description	Quantity
----------	-------------	----------

Recliner Treatment Manual 300 Pound Capacity Polyurethane Foam

MCPF74	Fixed Arm	1/EA
MCPF79	Fire Rated Fixed Arm	1/EA
MCPF80	Heat Massage Fixed Arm	1/EA
MCPF81	Heat Massage Fire Rated Fixed Arm	1/EA
MCPF84	Heat Fixed Arm	1/EA
MCPF85	Heat Fire Rated Fixed Arm	1/EA
MCPS70	Swingaway Arm	1/EA
MCPS72	Heat Massage Swingaway Arm	1/EA
MCPS77	Fire Rated Swingaway Arm	1/EA
MCPS78	Heat Massage Fire Rated Swingaway Arm	1/EA
MCPS82	Heat Swingaway Arm	1/EA
MCPS83	Heat Fire Rated Swingaway Arm	1/EA

Complementary Products

Item No.	Description
RA1070	Vinyl Upholstery Upgrade Recliner Factory Option for MediChoice Recliners
RA1080	Table Second Side Factory Option for MediChoice Treatment Chairs
RA1030	Cup Holder Factory Option for MediChoice Treatment Chairs
RA1040	Pillow Adjustable for Head or Lumbar Support Fits All MediChoice Treatment Recliners Factory Add on Option
RA1050	IV Pole Bracket Accessory Fits All MediChoice Treatment Chairs
RA1060	IV Pole With Two Hooks Accessory Fits All MediChoice Treatment Chairs
RA00FB	Table First Side Factory Option for MediChoice Manual Treatment Chairs with Fixed or Swingaway Arms
RA00SB	Table Second Side Factory Option for MediChoice Manual Treatment Chairs with Fixed or Swingaway Arms
RA3020	Push Handle Factory Option for MediChoice Economy Treatment Chairs
RA1020	Bag Hook Fits All MediChoice Treatment Chairs

RECLINERS



Treatment Chairs



- Infinite position recliner made of all steel powder-coated frame construction for greater durability and longer life
- Space saver, zero-wall mechanism provides easy operation and unmatched stability
- Standard models equipped with 5 Inch thermal plastic casters, push handle, fold-away side table, retractable foot tray and swing-away arms option
- Economy model (MC3596 and MCFR3596) equipped with 3 Inch steel casters with front brakes
- 8 standard vinyl colors and an extensive line of accessory options available
- Capacity: 300 pounds
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- Quick release seat for easy cleaning and maintenance
- Z-coil suspension system affords greater comfort
- 10-year limited manufacturer’s warranty

Item No.	Description	Quantity
----------	-------------	----------

Chair Treatment 300 Pound Capacity Steel Frame

MC3547	5 Inch Plastic Caster	1/EA
MC3596	Economy 3 Inch Steel Caster	1/EA
MC3858	Swingaway Arm 5 Inch Plastic Caster	1/EA
MCFR3547	Fire Rated 5 Inch Plastic Caster	1/EA
MCFR3596	Economy Fire Rated 3 Inch Steel Caster	1/EA
MCFR3858	Power Recline Swingaway Arm Fire Rated 5 Inch Plastic Caster	1/EA

RECLINERS



RECLINERS



Standard Manual Recliners



- Quality at an economical price
- User friendly, push-back reclining mechanism, allowing patients to easily adjust to three positions
- Swing open arms are standard for ease of cleaning and lateral patient transfer
- 3 Inch thermoplastic casters with front brakes
- Capacity: 275 pounds
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- Trendelenburg comes standard with a foot-release actuator on both sides
- Fold-away side table and push handle included
- All steel, powder-coated frame construction for durability

Recliners/Transporters



- Infinite positions recliner made of all steel power-coated frame construction for greater durability and longer life
- Space saver, zero-wall mechanism provides ease of operation and unmatched stability
- Swing open arms are standard for ease of cleaning and lateral patient transfer
- 4 Inch dual thermoplastic casters (corrosion resistant and roll with ease)
- 8 standard vinyl color options and an extensive line of accessory options available
- Power control handset and battery back-up are standard (on 1618MC065P model)
- Quick release seat for easy cleaning and maintenance
- Capacity: 300 pounds
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- Z-coil suspension system affords greater comfort

Item No.	Description	Quantity
MC0220	Recliner Clinical Manual 3 Inch Thermoplastic Caster 275 Pound Capacity Steel Frame	1/EA

Item No.	Description	Quantity
----------	-------------	----------

Recliner Clinical Transporter Manual Swingaway Arm 300 Pound Capacity Steel Frame

MC0650	Manual Swingaway Arm	1/EA
MC065P	Power Recline Swingaway Arm	1/EA

Complementary Products

Item No.	Description
RA1020	Bag Hook Fits All MediChoice Treatment Chairs
RA1030	Cup Holder Factory Option for MediChoice Treatment Chairs
RA1040	Pillow Adjustable for Head or Lumbar Support Fits All MediChoice Treatment Recliners Factory Add on Option
RA1050	IV Pole Bracket Accessory Fits All MediChoice Treatment Chairs
RA1060	IV Pole With Two Hooks Accessory Fits All MediChoice Treatment Chairs
RA1070	Vinyl Upholstery Upgrade Recliner Factory Option for MediChoice Recliners
RA1080	Table Second Side Factory Option for MediChoice Treatment Chairs

Complementary Products

Item No.	Description
RA1020	Bag Hook Fits All MediChoice Treatment Chairs
RA1030	Cup Holder Factory Option for MediChoice Treatment Chairs
RA1040	Pillow Adjustable for Head or Lumbar Support Fits All MediChoice Treatment Recliners Factory Add on Option
RA1050	IV Pole Bracket Accessory Fits All MediChoice Treatment Chairs
RA1060	IV Pole with Two Hooks Accessory Fits All MediChoice Treatment Chairs
RA1070	Vinyl Upholstery Upgrade Recliner Factory Option for MediChoice Recliners
RA1080	Table Second Side Factory Option for MediChoice Treatment Chairs
RA1090	Trendelenburg Position Manual Release Factory Option for MediChoice Recliners
RA2010	Electric Heat Recliner Seat And Back Factory Option for MediChoice Manual Recliners
RA2020	Electric Heat Massage Recliner Factory Option for MediChoice Manual Recliners

RECLINERS



Power Recliner Treatment Chairs

made in the USA NOT MADE WITH natural rubber latex

- Infinite position power recliner made of all steel powder-coated frame construction for greater durability and longer life
- Space saver, zero-wall power mechanism provides easy operation and unmatched stability
- Standard models are equipped with 5 Inch thermal plastic casters (2 front locking, 1 rear directional), push handle, fold-away side table, retractable foot tray, and swing-away arms option
- Economy models (MC359P and MCFR359P) equipped with 3 Inch steel casters with front brakes
- Power control handset and battery back-up are standard
- 8 standard vinyl color options and an extensive line of accessory options available
- Capacity: 300 pounds
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- Quick release seat for easy cleaning and maintenance
- Z-coil suspension system affords greater comfort
- 10-year limited manufacturer's warranty

Item No.	Description	Quantity
Chair Treatment 300 Pound Capacity Steel Frame		
MC354P	Power Recline 5 Inch Plastic Caster	1/EA
MC359P	Power Recline Economy 3 Inch Steel Caster	1/EA
MC385P	Power Recline Swingaway Arm 5 Inch Plastic Caster	1/EA
MCFR354P	Power Recline Fire Rated 5 Inch Plastic Caster	1/EA
MCFR359P	Power Recline Economy Fire Rated 3 Inch Steel Caster	1/EA
MCFR385P	Swingaway Arm Fire Rated 5 Inch Plastic Caster	1/EA

BARIATRICS



Bariatric Treatment Chairs

made in the USA NOT MADE WITH natural rubber latex

- Infinite position recliner made of all steel powder-coated frame construction for greater durability and longer life
- Space saver, zero-wall mechanism provides easy operation and unmatched stability
- Foldaway side table, push handle, and retractable foot tray included; optional swing away arms; 5 Inch thermal plastic casters (2 front locking, 1 rear directional)
- 8 standard vinyl colors and an extensive line of accessory options available
- Capacity: 500 pounds
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- 5 Inch thermal plastic casters for superior corrosion resistance and transportability
- Quick release seat for easy cleaning and maintenance
- Z-coil suspension system affords greater comfort
- 10-year limited manufacturer's warranty

Item No.	Description	Quantity
Chair Treatment Bariatric 5 Inch Caster 500 Pound Capacity Steel Frame		
MC3566	Extra-Wide Seat	1/EA
MC3868	Swingaway Arm	1/EA
MCFR3566	Extra-Wide Seat Fire Rated	1/EA
MCFR3868	Power Recline with Swingaway Arm Fire-Rated	1/EA

BARIATRICS



Bariatric Power Recliner Treatment Chairs

made in the **USA** NOT MADE WITH natural rubber latex

- Infinite power position recliner made of all steel powder-coated frame construction for greater durability and longer life
- Space saver, zero-wall power mechanism provides easy operation and unmatched stability
- Foldaway side table, push handle, and retractable foot tray included
- 8 standard vinyl colors and an extensive line of accessory options available
- Capacity: 500 pounds
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- 5 Inch thermal plastic casters for superior corrosion resistance and transportability (2 front locking, 1 rear directional)
- Power control handset and battery back-up are standard
- Quick release seat for cleaning and maintenance
- Z-coil suspension system affords greater comfort
- 10-year limited manufacturers warranty

Item No.	Description	Quantity
Chair Treatment Bariatric 5 Inch Caster 500 Pound Capacity Steel Frame		
MC356P	Power Recline	1/EA
MC386P	Power Recline with Swingaway Arm	1/EA
MCFR356P	Power Recline Fire-Rated	1/EA
MCFR386P	Swingaway Arm Fire Rated	1/EA

BARIATRICS



Bariatric Recliners/Transporters

made in the **USA** NOT MADE WITH natural rubber latex

- Infinite position recliners made of all steel powder-coated frame construction for greater durability and longer life
- Space saver, zero-wall mechanism provides easy operation and unmatched stability
- Quick release seat for easy cleaning and maintenance
- Swing open arms are standard for ease of cleaning and lateral patient transfer
- 4 Inch twin wheel casters (corrosion resistant and roll with ease)
- Eight standard vinyl color options and an extensive line of accessory options available
- Power control handset and battery back-up are standard (on 1618MC066P model)
- Capacity: 500 lbs.
- Covered in medical grade material for easy cleaning and limited exposure to metal parts
- Z-coil suspension system affords greater comfort
- 10-year limited manufacturer's warranty

Item No.	Description	Quantity
MC0660	Manual Bariatric with Swingaway Arm	1/EA
MC066P	Power Recline Bariatric with Swingaway Arm	1/EA

Complementary Products

ITEM NO.	DESCRIPTION
RA1040	Adjustable pillow for head or lumbar support, fits all medichoice treatment recliners factory add on option
RA1050	IV pole bracket accessory fits all medichoice treatment chairs
RA1020	Bag hook fits all medichoice treatment chairs
RA1030	Cup holder factory option for medichoice treatment chair
RA1060	IV pole with two hooks accessory fits all medichoice treatment chairs
RA1070	Vinyl upholstery upgrade recliner factory option for medichoice recliners
RA1080	Table second side factory option for medichoice treatment chair
RA1090	Trendelenburg position manual release factory option for medichoice recliners
RA2010	Electric heat recliner seat and back factory option for medichoice manual recliners
RA2020	Electric heat massage recliner factory option for medichoice manual recliners

SLEEPER CHAIRS

Deluxe Overnighter Sleeper Chairs



- Lever-activated foot rest and actuation release for back rest
- Foot-release Trendelenburg (access from both sides of chair)
- Upright overall depth (with/without handle): 44 Inch/40.5 Inch
- Steel, powder-coated frame
- 3 Inch thermoplastic casters (rear locking)
- Swing arm on patient's left side
- Quick-release seat
- Weight capacity: 400 pounds
- Back height (from floor/seat): 45 Inch/23 Inch
- Overall width (table down/up): 35 Inch/45 Inch
- Seat width: 23.5 Inch



Item No.	Description	Quantity
MC0720	Chair Sleeper Deluxe Overnighter 3 Inch Casters Steel Frame 400 Pound Capacity 45 Inch x 35 Inch x 44 Inch	1/EA

Classic Overnighter Sleeper Chairs



- Made of steel powder-coated frame construction for greater durability and longer life
- Three positions available for maximum comfort: upright, lounge and flat (sleeping)
- Minimum moving parts for longer life and less maintenance
- 3 Inch casters, including rear locking casters for ease of movement
- Capacity: 300 pounds
- No loose parts that will be lost or misplaced
- Rear vinyl guard for wall and vinyl protection
- 10-year limited manufacturer's warranty



Item No.	Description	Quantity
Chair Sleeper 3 Inch Caster 300 Pound Capacity Steel Frame		
MC3526	Overnight Classic	1/EA
MC3526WA	Overnight Medium Oak Arm Cap	1/EA
MCFR3526	Overnight Classic Fire Rated	1/EA

Upholstery Colors†

Treatment Chairs and Sleeper Chairs



Recliner/Treatment Chair Accessories

- These are accessories for the MediChoice treatment chairs
- Add on factory option is added when recliner is produced
- Add on accessory can be added after recliner is produced

Item No.	Description	Quantity
For all MediChoice Recliners Excluding all Clinical Care & Powered Bariatric Recliners		
MC3	Upholstery Fire Rated Finish Ca-133	1/EA
Factory Option for MediChoice Manual Treatment Chairs with Fixed or Swingaway Arms		
RA00FB	Table First Side	1/EA
RA00SB	Table Second Side	1/EA
Fits all MediChoice Treatment Chairs		
RA1010	Arm Cover	1/PR
RA1020	Bag Hook	1/EA
RA1030	Cup Holder Factory Option	1/EA
RA1040	Pillow Adjustable for Head or Lumbar Support Factory Add on Option	1/EA
RA1050	IV Pole Bracket Accessory	1/EA
RA1060	IV Pole with Two Hooks Accessory	1/EA
RA1080	Table Second Side Factory Option	1/EA
Factory Option for MediChoice Recliners		
RA1070	Vinyl Upholstery Upgrade	1/EA
RA1090	Trendelenburg Position Manual Release	1/EA
Factory Option for MediChoice Manual Recliners		
RA2010	Electric Heat Recliner Seat and Back	1/EA
RA2020	Electric Heat Massage Recliner	1/EA
Factory Option for MediChoice Power Recliners		
RA2121	Heat and Massage	1/EA
Factory Option for MediChoice Economy Treatment Chair		
RA3010	Table First Side	1/EA
RA3020	Push Handle	1/EA
RA3030	Tray Foot Retractable	1/EA
RA4020	Controller Dual-Hand	1/EA
RALTBRTC	Tray Lap Curved Black	1/EA

RECLINERS

RECLINERS

Deluxe Clinical Care Recliners



- Recliner mechanism enables caregiver/patient to position the back from seated to TV position to full recline
- Contoured seat and back with integrated footrest enhances patient comfort
- Higher floor-to-back height enables the caregiver to raise the patient up with reduced effort, lessening the chance for a back injury
- Molded side tables on both sides of the recliner conveniently hold prep materials
- Side tables have a smooth surface with rounded edges for easy cleaning and will not trap contaminants
- Available with Heat and Massage, providing three different heat settings for the seat and back, and four massage settings for the upper and lower back
- Choose from 46 standard vinyl upholstery colors with PERMABLOK3 advanced vinyl protection or 10 Silververtex vinyl colors with SILVERGUARD silver ion antimicrobial technology
- Weight Bearing Capacity: 350 pounds evenly distributed
- Optional Drop Arm feature allows for safer lateral transfers for both patient and caregiver
- Trendelenburg Position can be actuated from either side of recliner by a conveniently located foot-pedal
- Built-in arm supports in the back give improved support while in the Trendelenburg Position
- Removable side panels allow for easy access to recliner mechanism for maintenance and cleaning
- Directional wheel lock on right rear caster reduces “crabbing” of recliner and enables the chair to be moved down hallways with less effort
- Four plastic Tente®-swivel casters facilitate transport and maneuverability in confined areas
- Limited lifetime manufacturer’s warranty on frame
- 3-year manufacturer’s warranty on mechanical components
- 2-year limited manufacturer’s warranty on heat and massage
- 1-year limited upholstery manufacturer’s warranty
- Complies with California Technical Bulletin 133 fire rating requirements



Standard Clinical Care Recliners



- Three positions from seated (full upright), TV position, to fully reclined plus Trendelenburg
- Varied-density seat back construction provides excellent support and comfort
- Molded side tables on both sides of the recliner conveniently hold prep materials
- Side tables have a smooth surface with rounded edges for easy cleaning and will not trap contaminants
- Available with Heat and Massage, providing three different heat settings for the seat and back, and four massage settings for the upper and lower back
- Choose from 46 standard vinyl upholstery colors with PERMABLOK3 advanced vinyl protection or 10 Silververtex vinyl with SILVERGUARD silver ion technology
- Removable side panels allow for easy access to recliner mechanism for maintenance and cleaning
- Four plastic Tente®-swivel casters facilitate transport and maneuverability in confined areas
- Directional wheel lock on right rear caster reduces “crabbing” of recliner and enables the chair to be moved down hallways with less effort
- Trigger activated gas cylinder for quick Trendelenburg positioning and independent positioning of back angle
- Weight Bearing Capacity: 290 pounds evenly distributed
- Manufacturer’s limited lifetime warranty on frame
- 3-year warranty on mechanical components
- 2-year warranty on heat and massage unit
- 1-year warranty on upholstered components



Item No.	Description	Quantity
Recliner Clinical Care Fire Rated 350 Pound Capacity		
MC566DH	Deluxe Infinite Position Heat Massage Drop Arm Silverguard Silver Ion Technology	1/EA
MC566DHS	Deluxe Infinite Position Heat Massage Drop Arm Silverguard Silver Ion Technology	1/EA
MC566DS	Deluxe Infinite Position Drop Arm Silverguard Silver Ion Technology	1/EA
MC566GH	Deluxe Infinite Position Heat Massage Silverguard Silver Ion Technology	1/EA
MC566GHS	Deluxe Infinite Position Heat Massage Silverguard Silver Ion Technology	1/EA
MC566GS	Infinite Position Silverguard Silver Ion Technology	1/EA
MCFR566G	Deluxe Infinite Position	1/EA
MCFR566DG	Deluxe Infinite Position Trendelenburg Drop Arm Silverguard Silver Ion Technology	1/EA

Item No.	Description	Quantity
Recliner Clinical Care 290 Pound Capacity		
MC577GH	Standard Heat Massage 3 Position Fire Rated	1/EA
MC577GHS	Standard Trendelenburg Heat Massage 3 Position Fire Rated Silver Upholstery	1/EA
MC577GS	Standard 3 Position Fire Rated Silver Upholstery	1/EA
MCFR577RG	Standard 3 Position Fire Rated	1/EA

RECLINERS



Extra-Wide Clinical Care Recliners

NOT MADE WITH
natural rubber latex

- Dual gas cylinders give smooth power assistance to and from the Trendelenburg position
- Blow-molded side tables conveniently hold prep materials
- Side tables have a smooth surface with rounded edges for easy cleaning and will not trap contaminants
- Extra-wide seat width of 24 Inch designed to comfortably sit and position larger patients
- Available with heat and massage; providing 3 different heat settings for the seat and back, and 4 massage settings for the upper and lower back.
- Choose from 46 standard vinyl upholstery colors with PERMABLOK3 advanced vinyl protection, or 10 Silvertex vinyl colors with SILVERGUARD silver ion technology
- Removable side panels allow for easy access to recliner mechanism for maintenance and cleaning
- Four plastic Tente® swivel casters facilitate transport and maneuverability in confined areas
- Directional wheel lock on right rear caster reduces “crabbing” of recliner and enables the chair to be moved down hallways with less effort
- Ergonomically designed padded armrests for added patient comfort
- Weight Bearing Capacity: 450 pounds evenly distributed
- Limited lifetime manufacturer’s warranty on frame
- 3-year limited manufacturer’s warranty on mechanical components
- 2-year limited manufacturer’s warranty on head and massage units
- 1-year limited manufacturer’s warranty on upholstered components
- Complies with California Technical Bulletin 133 fire rating requirements

Item No.	Description	Quantity
Recliner Clinical Silverguard Silver Ion Technology 450 Pound Capacity		
MC587W	Extra-Wide Infinite Position Fire Rated Trendelenburg	1/EA
MC587WH	Extra-Wide Infinite Position Fire Rated Trendelenburg Heat Massage	1/EA
MC587WHS	Extra-Wide Infinite Position Fire Rated Trendelenburg Heat Massage Silvertex Vinyl	1/EA
MC587WS	Extra-Wide Infinite Position Fire Rated Trendelenburg Silvertex Vinyl	1/EA

BARIATRICS



Powered Bariatric Recliners

made in the USA **NOT MADE WITH**
natural rubber latex

- Heavy-gauge steel, epoxy coated frame recliner with drop arms for easier lateral transfers
- Four-function hand pendant provides motor activation allowing caregivers to safely recline and raise patients with less strain and effort
- Quick-release trigger mechanism on push handle allows chair to swiftly move into Trendelenburg Position
- Fully equipped with rechargeable battery pack for full operation away from a wall outlet
- Available in 4 vinyl upholstery colors with AEGIS Microbe Shield: Blue Ridge, Doe Skin, Jade, and Rosewood
- Four swivel casters facilitate transport and maneuverability in confined areas
- Durable, high impact side tables are easily activated by caregiver and sit flush with chair arm to conveniently hold prep materials
- Legrest and backrest raise and lower independently giving the operator complete control over positioning and comfort of patient
- Complete with attached head pillow for extra support and added comfort
- Extra Wide Seat: 25.25 Inches between arms
- Weight Bearing Capacity: 700 pounds evenly distributed
- Limited manufacturer’s lifetime warranty on frame
- 3-year manufacturer’s warranty on mechanical components
- 2-year limited warranty on heat and massage unit
- 1-year manufacturer’s warranty on upholstered components
- Complies with California Technical Bulletin 133 fire-rating requirements

Item No.	Description	Quantity
MC588W	Recliner Powered Bariatric Fire Rated Drop Arm 700 Pound Capacity Steel Frame	1/EA

Upholstery Colors*

Recliners

					
DOL-6705 Sand	DOL-6714 Burgundy	DOL-6717 Moss	DOL-6724 Jet	DOL-6725 Indigo	DOL-6727 Meteor
					
1305 Convoy	1311 Fence	1315 Sentinel	1318 Steel	1323 Moat	1325 Armor
					
INV-6201 Honey	INV-6203 Sand	INV-6204 Citron	INV-6205 Tide	INV-6206 Harbor	INV-6208 Caramel
					
PAT-427 Bluebridge	PAT-432 Imperial Blue	PAT-454 Royal	PAT-8504 Sand	PAT-8516 Chestnut	PAT-8533 Port
					
PAT-8561 Mulberry	PAT-8567 Dove	PAT-8584 Taupe	PAT-8602 Oatmeal	PAT-8604 Concrete	PAT-8611 Terra Cotta
					
PAT-8615 Berry	PAT-8622 Grey Moss	PAT-8624 Granny Smith	PAT-8626 Eucalyptus	PAT-8628 Sea	REN-9200 Chalk
					
REN-9201 Cobblestone	REN-9202 Periwinkle	REN-9203 Pyrite	REN-9205 Gypsum	REN-9207 Ice Blue	REN-9208 Wineberry
					
REN-9210 Cypress	REN-9212 Blue Jade	REN-9214 Steel Blue	SAV-401 Vintage Gold	SAV-409 Warm Taupe	SAV-9701 Bay Leaf
					
SAV-9711 Tallow	SAV-9715 Eucalyptus	SAV-9722 Desert Sage	SAV-9728 Navy	SAV-9754 Sisal	SAV-9757 Turmeric
					
SAV-9764 Steel Grey	STX-8801 Marine Blue	STX-8805 Jet	STX-8808 Champagne	STX-8809 Taupe	STX-8810 Chestnut
					
STX-8814 Umber	STX-8815 Borscht	STX-8818 Sage	STX-8819 Celery	STX-8822 Storm	TWI-851 Doe Skin
					
TWI-857 Jade	TWI-863 Rosewood				




Bariatric Blood Drawing Chairs with Flip Arms

made in the USA NOT MADE WITH natural rubber latex

- Larger seat with additional depth and an extra center support for added strength
- Adjustable height and depth flip arms, drawer model has 1 flip arm
- Armrests position in place when in use, then flip out of the way for easy exit
- Easy-clean, heavy weight, vinyl upholstery available in 16 color options
- Capacity: 700 pounds
- Overall Dimensions: 48 Inch Width x 32 Inch Depth, drawer model 55 Inch Width x 34 Inch Depth
- Comfortable 2 Inch padded seat, back and armrests
- Adjustable foot levelers
- Chrome armrest with two locking knobs
- Neutral grey powder-coated frame
- Height: 20 Inch
- Arm Height Adjustment: 29 Inch to 34 Inch
- Seat Dimensions: 35 Inch x 21 Inch
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- 5-year limited warranty
- All blood drawing products are non-returnable unless they are unused in unopened, original packaging
- **IMPORTANT NOTE: MediChoice Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.**

Item No.	Description	Quantity
Chair Blood Drawing Bariatric Vinyl Upholstery 700 Pound Capacity		
BLDC4002	With Flip Arm Vinyl Upholstery 48 Inch x 32 Inch	1/EA
BLDC4052	With Drawer Flip Arm Vinyl Upholstery Bariatric 55 Inch x 34 Inch	1/EA
BLDC4502	With Flip Arm Fire Rated Vinyl Upholstery 48 Inch x 32 Inch	1/EA
BLDC4552	With Drawer Flip Arm Fire Rated Vinyl Upholstery 55 Inch x 34 Inch	1/EA



BLOOD DRAWING CHAIRS



Extra-Wide Blood Drawing Chairs with Flip-Arms



- Adjustable height and depth flip arms
- Armrests position in place when in use then flips out of the way for easy exit
- All welded heavy-duty frame with easy-clean durable antimicrobial vinyl upholstery available in 16 colors
- Capacity: 400 pounds
- Accommodates larger patient or child with parent
- Padded and upholstered seat back and armrests
- Flip arms: height adjusts 29 Inch to 39 Inch; depth adjusts 12.5 Inch to 17 Inch
- Dimensions:
 - Overall: 46 Inch Width × 26 Inch Depth
 - Seat: 20 Inch Height × 34 Inch Width × 18 Inch Depth
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements
- 5-year limited warranty
- All blood drawing products are non-returnable unless they are unused in unopened, original packaging
- IMPORTANT NOTE: MediChoice Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended

Item No.	Description	Quantity
Chair Blood Drawing Extra-Wide Vinyl Upholstery 46 Inch × 26 Inch 400 Pound Capacity		
BLDC3002	With Flip Arm	1/EA
BLDCFR3002	With Flip Arm Fire Rated	1/EA

BLOOD DRAWING CHAIRS



Standard Blood Drawing Chairs



- Adjustable height and depth flip arms 27 Inch to 33 Inch
- Armrest locks in place when in use then flips out of the way for easy exit
- One piece, black plastic seat
- All steel, welded frame
- Capacity: 300 pounds
- Dimensions:
 - Overall: 34 Inch Width × 26 Inch Depth
 - Seat: 20 Inch Height × 17 Inch Width × 16 Inch Depth
- Fully adjustable laminate pivot armrests
- 5-year limited warranty
- All blood drawing products are non-returnable unless they are unused in unopened, original packaging
- IMPORTANT NOTE: MediChoice Blood Drawing Chairs are not engineered to restrain a fainting or falling occupant. The chair's occupant should NEVER be left unattended.

Item No.	Description	Quantity
BLDC3001	Chair Blood Drawing 1-Piece Seat Steel Frame Black 34 Inch × 26 Inch 300 Pound Capacity	1/EA

TREATMENT TABLES



Straight Line Treatment Tables

made in the **USA** NOT MADE WITH natural rubber latex

- Upholstered top with 2 Inch foam system padding
- Hardwood legs and adjustable back rest
- Choice of model with or without laminate shelf
- Available in 23 upholstery colors with choice of two wood finishes: dark cherry or natural
- Capacity: 400 pound
- Dimensions: 72 Inch Length × 31 Inch Height; available in 24 Inch Width or 27 Inch Width
- Rounded corners for patient safety
- Tops attached to frame by steel angle
- Stain-resistant upholstery and surfaces
- Paper cutter and dispenser included
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements in 10 upholstery colors
- 5-year limited manufacturer's warranty

Item No.	Description	Quantity
Table Treatment 400 Pound Capacity		
TRTB1901	H-Brace 72 Inch × 31 Inch × 24 Inch	1/EA
TRTB1902	H-Brace 72 Inch × 31 Inch × 27 Inch	1/EA
TRTB1903	With Shelf 72 Inch × 31 Inch × 24 Inch	1/EA
TRTB1904	With Shelf 72 Inch × 31 Inch × 27 Inch	1/EA
TRTBFR1901	H-Brace Fire Rated 72 Inch × 31 Inch × 24 Inch	1/EA
TRTBFR1902	H-Brace Fire Rated 72 Inch × 31 Inch × 27 Inch	1/EA
TRTBFR1903	Fire Rated with Shelf 72 Inch × 31 Inch × 24 Inch	1/EA
TRTBFR1904	Fire Rated with Shelf 72 Inch × 31 Inch × 27 Inch	1/EA

TREATMENT TABLES




Pediatric Treatment Tables with Sliding Doors

made in the **USA** NOT MADE WITH natural rubber latex

- Cabinet style with adjustable back rest and sliding doors on both sides
- All models ship by common carrier only and are fully assembled
- Available in color coordinated door laminate and vinyl: 11 colors of laminate and 23 colors of vinyl to choose from
- Capacity: 400 pounds
- Dimensions: 72 Inch Length × 31 Inch Height × 27 Inch Width
- Easy-clean, scratch-resistant and mar-resistant laminate surfaces
- 2 Inch thick foam padding
- Chip-resistant PVC edge trim
- All stress points reinforced by steel angle
- Optional large colorful graphics
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements in 10 color options
- 5-year limited manufacturer's warranty

Item No.	Description	Quantity
Table Treatment 400 Pound Capacity 72 Inch × 31 Inch × 27 Inch		
PEDT1401	Pediatric with Sliding Door	1/EA
PEDTFR1401	Pediatric with Sliding Door Fire Rated	1/EA



A few of the designs available for Pediatric Treatment Tables

Please see our full line of Exam Table Paper on page 91 in the Exam Room Products Section.

TREATMENT TABLES



Pediatric Treatment Tables with Built-In Digital Scale



- Top rimmed on three sides for safety
- Easy-clean, replaceable center pad and laminate base
- Built-in paper dispenser, paper cutter, and infantometer with folding ends for measuring height
- Digital Doran scale with easy-read display, auto shut-off, and BMI functions
- Available in 14 colors of laminate, 23 colors of vinyl
- Weight capacity: 40 pounds
- Integrated concealed paper dispenser and cutter
- One drawer and two door storage areas
- Comfortable work height of 36 Inches
- Adjustable built-in levelers provide a stable surface
- Optional large colorful graphics
- Fire-rated model available to comply with California, Illinois and Massachusetts requirements in 10 color options
- Dimensions: 43 Inch Length x 36 Inch Height x 24 Inch Width
- 5-year manufacturer's warranty on table; 1-year manufacturer's warranty on scale

Item No.	Description	Quantity
Table Treatment 43 Inch x 36 Inch x 24 Inch		
PEDT1402	Pediatric with Scale	1/EA
PEDTFR1402	Pediatric with Scale Fire Rated	1/EA

Upholstery Colors

Exam Stools, Treatment Tables, and Blood Drawing Chairs



Laminate Colors

Exam Stools, Treatment Tables and Blood Drawing Chairs



OVERBED TABLES

Economy Overbed Tables

NOT MADE WITH
natural rubber latex

- Walnut wood-grain laminated 31 Inch top on .75 Inch particle board core with plastic edge molding
- Chrome plated H-base and column for aesthetics and durability
- Swivel casters provide smooth rolling action
- Table height: 28 Inch to 44 Inch
- Capacity dispersed weight: 30 pounds
- Weight: 22 pounds
- Dimensions:
 - Upper top surface: 31.5 Inch Length × 15 Inch Width
 - Base: 26.2 Inch Length × 15.25 Inch Width
- 1-year manufacturer's limited warranty



Item No.	Description	Quantity
BT89001	Table Overbed Non-Tilt Economy 30 Pound Capacity 31.5 Inch × 15 Inch	1/EA

Heavy-Duty Overbed Tables with Vanity

NOT MADE WITH
natural rubber latex

- Dual top design provides maximum space and includes a sculptured edge for spill containment
- Versatile split tops have a safety lock for secure holding until activation from a single touch lever
- Simple conversion from left to right hand by simply moving the column
- Removable cup holders
- Pneumatic lift height adjustment: 29.5 Inch to 44.5 Inch
- Dimensions:
 - Upper top surface: 18 Inch × 32 Inch
 - Lower top surface: 17.5 Inch × 26.5 Inch
- Powder coated H-base with 2 Inch dual wheel casters for easy movement
- Dual molded urethane ergonomic handles
- Easy set-up by the removal of just one screw
- Hinged mirror with vanity operates from either side
- 1-year manufacturer's warranty



Item No.	Description	Quantity
BT1740	Table Overbed Heavy-Duty with Vanity 18 Inch × 32 Inch	1/EA

UTILITY CARTS

Utility Carts with Wire Baskets

NOT MADE WITH
natural rubber latex **diverse**
supplier

- Strong chrome-plated folding steel frame
- Large 4 Inch rubber ball bearing casters
- Epoxy coated basket holders
- Optional 6 Inch and 12 Inch interchangeable baskets (sold separately)
- Load capacity: 70 pounds
- Shipping weight: 23 pounds cart, 5 pounds wire baskets
- Dimensions: 39 Inch × 19.5 Inch × 30 Inch
- 1-year manufacturer's warranty



Item No.	Description	Quantity
CRTW5401	CRTW5401 Cart Utility Frame Chrome Plated Steel 39 Inch × 9.5 Inch × 30 Inch on Large 4 Inch Rubber Ball-Bearing Casters	1/EA
CRTW5406	CRTW5406 Basket Steel White for Utility Cart CRTW5401 6 Inch	1/EA
CRTW5412	CRTW5412 Basket Steel White for Utility Cart CRTW5401 12 Inch	1/EA

Plastic Utility Carts

NOT MADE WITH
natural rubber latex

- Specifically engineered for healthcare's quick and easy sanitation needs
- Patented leg design for permanent snap-together assembly without tools
- Precision molded, high-density, polyethylene shelves to resist stains and odors
- Shipped unassembled to reduce freight costs
- Capacity: 200 pounds
- 4 Inch all-swivel casters
- Beige
- Dimensions: 17.125 Inch Width × 32.5 Inch Length × 34.875 Inch Height
- Shelf Size: 16 Inch Width × 24 Inch Length
- 1-year manufacturer's warranty



Item No.	Description	Quantity
CRTP5001	Cart Utility High-Density Polyethylene Plastic Beige 4 Inch All Swivel Casters Three-Shelf 17.125 Inch × 32.5 Inch × 34.875 Inch 200 Pound Capacity	1/EA

UTILITY CARTS

Plastic and Aluminum Utility Carts

NOT MADE WITH
natural rubber latex

- Sturdy brushed aluminum legs
- Precision molded, high-density, polyethylene shelves to resist stains and odors
- 4 Inch casters (two with locking brakes)
- Shipped unassembled to reduce freight costs, no tools required for assembly
- Capacity: 300 pounds evenly distributed on 3 shelves
- Actual Cart Weight: 16 pounds
- Quality construction
- Dimensions:
 - Overall: 17.5 Inch Width × 31.75 Inch Length × 38 Inch Height
 - Shelf Size: 16 Inch Width × 25 Inch Length
- 1-year manufacturer's warranty



Item No.	Description	Quantity
C RTP5002	Cart Utility High-Density Polyethylene Aluminum Plastic Black 4 Inch Caster Three-Shelf 17.5 Inch × 31.75 Inch × 38 Inch 300 Pound Capacity Aluminum Plastic Black	1/EA

Stainless Utility Carts

made in the USA **NOT MADE WITH**
natural rubber latex

- 3-shelf all welded stainless steel construction for increased durability and easy sanitation
- Integrated push handle allows maneuverability
- All-swivel casters:
 - Standard duty - 3.5 Inch diameter
 - Medium Duty - 4 Inch diameter
- Featuring new and improved sound deadening material for quieter operation
- Capacity:
 - Standard duty - 300 Pounds
 - Medium Duty - 500 Pounds
- Protective bumpers on legs to safeguard furnishings
- Dimensions:
 - Standard Duty: 16.25 Inch Width × 27.25 Inch Length × 32 Inch Height
 - Medium Duty: 16.75 Inch Width × 27.625 Inch Length × 32.125 Inch Height
 - Shelf Size: 15.5 Inch Width × 24 Inch Length
- 1-year manufacturer's warranty



Item No.	Description	Quantity
C RT S5201	Standard-Duty 3.5 Inch Diameter All-Swivel Caster Three-Shelf 16.25 Inch × 27.25 Inch × 32 Inch 300 Pound Capacity	1/EA
C RT S5202	Medium-Duty 4 Inch Diameter All-Swivel Caster Three-Shelf 16.75 Inch × 27.625 Inch × 32.125 Inch 500 Pound Capacity	1/EA

UTILITY TABLES

Stainless Utility Tables with Shelf and Four-Sided Rail

made in the USA **NOT MADE WITH**
natural rubber latex **diverse supplier**

- Four sided guardrails
- 3 Inch ball-bearing swivel casters
- Top and bottom shelf made of 18-gauge #304 stainless steel with edges bent and hemmed
- Knockdown style
- Dimensions: 20 Inch Wide × 34 Inch Height × 16 Inch Deep
- 3-year manufacturer's warranty



Item No.	Description	Quantity
UTSS2020	Table Utility Stainless Steel with Shelf and Four-Sided Rail 20 Inch × 34 Inch × 16 Inch on 3 Inch Ball-Bearing Swivel Casters	1/EA

Stainless Utility Tables with Shelf, Drawer and Rail

made in the USA **NOT MADE WITH**
natural rubber latex **diverse supplier**

- Constructed of durable #304 stainless steel
- Three-sided guardrail
- Built-in drawer keeps instruments and supplies clean and organized
- Lower shelf for extra storage (edges bent and hemmed for stability and safety)
- Dimensions: 20 Inch Wide × 35.25 Inch Height × 16 Inch Deep
- 2 Inch ball-bearing swivel casters
- Drawer dimensions: 14 Inch Wide × 3.75 Inch Height × 12 Inch Deep
- 3-year manufacturer's warranty



Item No.	Description	Quantity
UTSS2021	Table Utility Stainless Steel with Shelf, Drawer and Rail 20 Inch × 35.25 Inch × 16 Inch on 2 Inch Ball-Bearing Swivel Casters	1/EA

Please see our full line of Cart Covers on page 26 in the Central Supply Section.

COMPUTER ACCESSORIES

Submersible True-Type Keyboards

NOT MADE WITH natural rubber latex **diverse supplier**



- Standard keyboard layout and feel
- Laser etched keys reduce fading
- Antimicrobial: EPA registered imbedded silver ion technology won't wear off
- Available in Wireless or "Plug and Play" USB; no special hardware or software required
- 100% Waterproof and Spill-proof
- Submersible: can be soaked in 1% bleach solution
- Limited lifetime warranty

Item No.	Description	Quantity
Keyboard Computer Submersible Antimicrobial Dishwasher Safe		
SSKB107	True Type USB Quick Connect Black	1/EA, 10EA/CA
SSKB-SSK107	Multimedia True Type USB Quick Connect White	1/EA, 10EA/CA
SSKB-SV107L	Multimedia True Type USB Connect with Long Cable Black	1/EA, 10EA/CA

Silicone Keyboards

NOT MADE WITH natural rubber latex **diverse supplier**



- Silicone keyboards feature 4-point connections to reduce missed key strokes
- Plug-and-play USB connectivity; no special hardware or software required
- Antimicrobial: EPA registered imbedded silver ion technology won't wear off
- 100% Waterproof and Spill-proof; dishwasher safe
- Submersible: can be soaked in 1% bleach solution
- Limited lifetime warranty

Item No.	Description	Quantity
Keyboard Computer Submersible Antimicrobial Dishwasher Safe		
SSKB02	Backlit With Built-In Touchpad Pointing Device All-In-One Silicone Black	1/EA, 20EA/CA
SSKB-S106G2	Backlit All-In-One Silicone Black	1/EA, 20EA/CA
SSKB-SW106G2	Backlit All-In-One Silicone White	1/EA, 20EA/CA
SSKB-SW90PG2	Backlit With Built-In Touchpad Pointing Device All-In-One Silicone White	1/EA, 20EA/CA

COMPUTER ACCESSORIES

Waterproof True-Type Keyboards

NOT MADE WITH natural rubber latex **diverse supplier**



- Standard keyboard layout and feel
- Laser etched keys reduce fading
- Antimicrobial: EPA registered imbedded silver ion technology won't wear off
- Available in Wireless or "Plug and Play" USB; no special hardware or software required
- 100% Waterproof and Spill-proof
- Limited lifetime warranty

Item No.	Description	Quantity
Keyboard Computer Silver Storm Waterproof Antimicrobial Dishwasher Safe True-Type		
SSKB503	Black	1/EA, 10EA/CA
SSKB-STWK503	White	1/EA, 10EA/CA

Computer Mice

NOT MADE WITH natural rubber latex **diverse supplier**



- Wireless or "Plug and Play" USB; no special hardware or software required
- Antimicrobial: EPA registered imbedded silver ion technology won't wear off
- 100% Waterproof and Spill-proof
- Submersible: can be soaked in 1% bleach solution
- Limited lifetime warranty

Item No.	Description	Quantity
Mouse Computer Optical Scroll Wheel Antimicrobial Dishwasher Safe		
SSM-SSM3	Mouse Computer Optical Scroll Wheel Antimicrobial Dishwasher Safe Silicone Coating Black	1/EA, 50EA/CA
SSM-SSWM3	Scroll Wheel Silicone Coating White	1/EA, 100EA/CA
SSM001	Silver Storm Scroll Wheel Black	1/EA, 20EA/CA
SSMW042	Silver Storm Scroll Wheel White	1/EA, 20EA/CA

Ear Curettes

NOT MADE WITH natural rubber latex **diverse supplier**

- Made with a smooth, flexible polypropylene material
- Rounded tips increase patient comfort and reduce risk of injury
- Available in four different tip styles to accommodate a variety of patients, including adult and child sizes
- Color-coded for ease-of-use and product identification
- Single use to help prevent cross-contamination



Item No.	Description	Quantity
Curette Ear Single-Use Disposable Propylene		
CUR1030	Rigid Loop 3 Millimeter Tip Pink	1/EA, 50 EA/BX, 10 BX/CA
CUR1029	Rigid Loop 4 Millimeter Tip Green	1/EA, 50 EA/BX, 10 BX/CA
CUR1028	Flexible Loop 4 Millimeter Tip White	1/EA, 50 EA/BX, 10 BX/CA
CUR1031	Scoop 2 Millimeter Tip Gray For Infant	1/EA, 50 EA/BX, 10 BX/CA

Digital Thermometers, Kits and Sheaths

- Dual-Scale thermometers measure temperature in both Fahrenheit and Celsius: 89.6F-109.4F and 32.0C-43.9C
- Microtechnology ensures precision measurement within 0.2F/0.1C
- Beeper indicates when thermometer is ready for use and when temperature reading is available
- Memory stores most recent measured temperature
- Single-patient use Sheaths provide added infection control and cross contamination protection
- Performs automatic self-test to ensure accuracy
- Operating instructions included
- Automatic shut-off
- Liquid crystal display
- Water resistant
- Mercury-free
- Temperature reading time-range:
 - Oral: 40-90 seconds
 - Rectal: 35-60 seconds
 - Axillary: 60-180 seconds
 - Flexible tip: 10 seconds



Item No.	Description	Quantity
Thermometer Digital		
916812A	Oral Dual-Scale Kit Blue Cap with 20 Sheaths	1/EA, 20 EA/CA
916814A	Oral Dual-Scale Thermometer Blue Cap	1/EA, 25 EA/CA
916820A	Oral Rectal Dual-Scale Neutral Cap	1/EA, 25 EA/CA
916823	Oral Rectal 10-Second Read Blue Cap	1/EA, 25 EA/CA
916815A	Rectal Dual-Scale Red Cap	1/EA, 25 EA/CA
916824	Rectal Dual-Scale 10-Second Read Red Cap	1/EA, 25 EA/CA
916821	Oral Celsius-Scale Blue Top	1/EA, 25 EA/CA
916822	Rectal Celsius-Scale Red Topice	1/EA, 25 EA/CA
206335	Sheath	1/EA, 100 EA/BX, 50 BX/CA

BLOOD PRESSURE

Blood Pressure Cuffs (Reusable)

NOT MADE WITH natural rubber latex



- Two piece style with nylon cuff and removable bladder
- Range guide on outside of cuff will direct clinicians to use a larger or smaller cuff to help prevent the misreading of blood pressure due to incorrect size
- Centered tube bladder helps clinicians achieve optimal arterial compression
- Cuff includes a gauge holder for manual blood pressure use
- Manufactured in compliance with FDA Quality System Regulation 21CFR Part 820; ISO 13485:2003, Medical devices. Meet ANSI/AAMI Standard SP10: 2002

Item No.	Description	Quantity
Cuff Blood Pressure Bladder Nylon Reusable		
M3205SH	1-Tube with Bayonet-Rectus-Hip Connector Child	1/EA, 5 EA/BG
M3204SH	1-Tube with Bayonet-Rectus-Hip Connector Small Adult	1/EA, 5 EA/BG
M3203SH	1-Tube with Bayonet-Rectus-Hip Connector Adult	1/EA, 5 EA/BG
M3203LSH	1-Tube with Bayonet-Rectus-Hip Connector Adult Long	1/EA, 5 EA/BG
M3202SH	1-Tube with Bayonet-Rectus-Hip Connector Large Adult	1/EA, 5 EA/BG
M3202LSH	1-Tube with Bayonet-Rectus-Hip Connector Large Adult Long	1/EA, 5 EA/BG
M3205SC	1-Tube with Screw Connector Child	1/EA, 5 EA/BG
M3204SC	1-Tube with Screw Connector Small Adult	1/EA, 5 EA/BG
M3203SC	1-Tube with Screw Connector Adult	1/EA, 5 EA/BG
M3203LSC	1-Tube with Screw Connector Adult Long	1/EA, 5 EA/BG
M3202SC	1-Tube with Screw Connector Large Adult	1/EA, 5 EA/BG
M3202LSC	1-Tube with Screw Connector Large Adult Long	1/EA, 5 EA/BG
M3205C	2-Tube with Screw Connector Child	1/EA, 5 EA/BG
M3204C	2-Tube with Screw Connector Small Adult	1/EA, 5 EA/BG
M3203C	2-Tube with Screw Connector Adult	1/EA, 5 EA/BG
M3203LC	2-Tube with Screw Connector Adult Long	1/EA, 5 EA/BG
M3202C	2-Tube with Screw Connector Large Adult	1/EA, 5 EA/BG
M3202LC	2-Tube with Screw Connector Large Adult	1/EA, 5 EA/BG
M3205M	2-Tube with Female Subminiature-Locking Connector Child	1/EA, 5 EA/BG
M3204M	2-Tube with Female Subminiature-Locking Connector Small Adult	1/EA, 5 EA/BG
M3203M	2-Tube with Female Subminiature-Locking Connector Adult	1/EA, 5 EA/BG
M3203LM	2-Tube with Female Subminiature-Locking Connector Adult Long	1/EA, 5 EA/BG
M3202M	2-Tube with Female Subminiature-Locking Connector Large Adult	1/EA, 5 EA/BG
M3202LM	2-Tube with Female Subminiature-Locking Connector Large Adult Long	1/EA, 5 EA/BG
M3405	Inflation System 2-Tube with Bulb and Valve Child	1/EA, 5 EA/BG
M3404	Inflation System 2-Tube with Bulb and Valve Small Adult	1/EA, 5 EA/BG
M3403	Inflation System 2-Tube with Bulb and Valve Adult	1/EA, 5 EA/BG
M3403L	Inflation System 2-Tube with Bulb and Valve Adult Long	1/EA, 5 EA/BG
M3402	Inflation System 2-Tube with Bulb and Valve Large Adult	1/EA, 5 EA/BG
M3402L	Inflation System 2-Tube with Bulb and Valve Large Adult Long	1/EA, 5 EA/BG

Size Range Guide: Reusable Cuffs

Size	Circumference Range (cm)
Child	13.0-21.0
Small Adult	20.0-28.0
Adult	27.0-35.0
Adult Long	27.0-35.0
Large Adult	34.0-42.0
Large Adult Long	34.0-42.0

Long sizes have an additional 13.5 centimeters in length

Please see our Nursery Kits on page 128 in the Maternity and Baby Care Section.

BLOOD PRESSURE

REUSABLE BLOOD PRESSURE CUFFS	TUBE		CONNECTOR			BULB & VALVE	SIZE						
	ONE	TWO	BAYONET/RECTUS	SCREW	FEMALE SUBMINIATURE-LOCKING		CHILD	SMALL ADULT	ADULT	ADULT LONG	ADULT LARGE	LARGE ADULT LONG	
1314M3205SH	●		●				●						
1314M3204SH	●		●					●					
1314M3203SH	●		●						●				
1314M3203LSH	●		●							●			
1314M3202SH	●		●								●		
1314M3202LSH	●		●									●	
1314M3205SC	●			●			●						
1314M3204SC	●			●				●					
1314M3203SC	●			●					●				
1314M3203LSC	●			●						●			
1314M3202SC	●			●							●		
1314M3202LSC	●			●								●	
1314M3205C		●		●			●						
1314M3204C		●		●				●					
1314M3203C		●		●					●				
1314M3203LC		●		●						●			
1314M3202C		●		●							●		
1314M3202LC		●		●								●	
1314M3205M		●			●		●						
1314M3204M		●			●			●					
1314M3203M		●			●				●				
1314M3203LM		●			●					●			
1314M3202M		●			●						●		
1314M3202LM		●			●							●	
13140M3405		●				●	●						
13140M3404		●				●		●					
13140M3403		●				●			●				
1314M3403L		●				●				●			
13140M3402		●				●					●		
1314M3402L		●				●						●	

Blood Pressure Cuff Accessories

made in the USA NOT MADE WITH natural rubber latex



- MediChoice offers a wide variety of cuff connection port options as well as cable-to-monitor plugs and adapters for most cuff connection scenarios.
- Manufactured in compliance with FDA Quality System Regulation 21CFR Part 820; ISO 13485:2003, Medical devices. Meet ANSI/AAMI Standard SP10: 2002

Item No.	Description	Quantity
Adapter Blood Pressure Cuff		
M9097	X 2-Tube Male and Female Subminiature to 2-Tube Screw Cuff for Datex-Ohmeda Monitor	1/EA
M9230	Y 1-Tube Male Bayonet to 2-Tube Screw Cuff for Philips HP Agilent Datascope Spacelabs Monitors	1/EA

Item No.	Description	Quantity
M9231	Y 1-Tube Luer Lock to 2-Tube Screw Cuff for Datascope Siemens Monitors	1/EA
M9236	X 2-Tube Female Subminiature to 2-Tube Screw Cuff for Marquette Welch Allyn LXI CVSM Monitor	1/EA
M9269	Y 1-Tube Screw to 2-Tube Screw Cuff for Welch Allyn Spot VSM Monitor	1/EA
M9274	X 2-Tube Bulb and Valve to 2-Tube Screw Cuff for Manual Sphygmomanometers	1/EA
M9314	Y 2-Tube Screw to 1-Tube Bayonet Cuff For Dinamap Monitor	1/EA
M9321	Y 2-Tube Female Subminiature to 1-Tube Bayonet Cuff for Marquette Welch Allyn LXI CVSM Monitors	1/EA
M9333	Y 2-Tube Male and Female Subminiature to 1-Tube Bayonet Cuff For Datex-Ohmeda Monitor	1/EA
M9313	Y 2-Tube Bulb and Valve to 1-Tube Bayonet Cuff for Manual Sphygmomanometers	1/EA
M9345	1-Tube Female Subminiature to 1-Tube Bayonet Cuff for Propaq	1/EA
M9329	1-Tube Screw to 1-Tube Bayonet Cuff for Welch Allyn Spot VSM Monitors	1/EA

BLOOD PRESSURE

Blood Pressure Cuffs (Single-Patient Use)

NOT MADE WITH
natural rubber latex



- Soft fabric cuffs provide patient comfort and are color-coded for easy identification by size
- Single-patient use aids in reducing the risk of cross-contamination and adherence to infection control protocols
- Range guide on outside of cuff directs clinicians to use an appropriately sized cuff to help prevent misreading of blood pressure due to incorrect size
- Centered tube bladder helps clinicians achieve optimal arterial compression
- Cuff includes a gauge holder for manual blood pressure use
- Manufactured in compliance with FDA Quality System Regulation 21CFR Part 820; ISO 13485:2003, Medical devices. Meet ANSI/AAMI Standard SP10: 2002

Item No.	Description	Quantity
Cuff Blood Pressure Bladder Single Patient Use Soft-fabric		
M39100	1-Tube With Bayonet-Rectus-HP Connector Child	1/EA, 20 EA/CA
M39085	1-Tube With Bayonet-Rectus-HP Connector Small Adult	1/EA, 20 EA/CA
M39048	1-Tube With Bayonet-Rectus-HP Connector Adult	1/EA, 20 EA/CA
M39019	1-Tube With Bayonet-Rectus-HP Connector Adult Long	1/EA, 20 EA/CA
M39230	1-Tube With Bayonet-Rectus-HP Connector Large Adult	1/EA, 20 EA/CA
M39060	1-Tube With Bayonet-Rectus-HP Connector Large Adult Long	1/EA, 20 EA/CA
M39091	1-Tube With Screw Connector Child	1/EA, 20 EA/CA
M39092	1-Tube With Screw Connector Small Adult	1/EA, 20 EA/CA
M39073	1-Tube With Screw Connector Adult	1/EA, 20 EA/CA
M39093	1-Tube With Screw Connector Adult Long	1/EA, 20 EA/CA

Item No.	Description	Quantity
M39074	1-Tube With Screw Connector Large Adult	1/EA, 20 EA/CA
M39075	1-Tube With Screw Connector Large Adult Long	1/EA, 20 EA/CA
M39011	2-Tube With Screw Connector Child	1/EA, 20 EA/CA
M39012	2-Tube With Screw Connector Small Adult	1/EA, 20 EA/CA
M39023	2-Tube With Screw Connector Adult	1/EA, 20 EA/CA
M39053	2-Tube With Screw Connector Adult Long	1/EA, 20 EA/CA
M39014	2-Tube With Screw Connector Large Adult	1/EA, 20 EA/CA
M39071	2-Tube With Screw Connector Large Adult Long	1/EA, 20 EA/CA
M39249	2-Tube With Female Subminiature-Locking Connector Child	1/EA, 20 EA/CA
M39213	2-Tube With Female Subminiature-Locking Connector Small Adult	1/EA, 20 EA/CA
M39268	2-Tube With Female Subminiature-Locking Connector Adult	1/EA, 20 EA/CA
M39081	2-Tube With Female Subminiature-Locking Connector Adult Long	1/EA, 20 EA/CA
M39149	2-Tube With Female Subminiature-Locking Connector Large Adult	1/EA, 20 EA/CA
M39018	2-Tube With Female Subminiature-Locking Connector Large Adult Long	1/EA, 20 EA/CA

Size Range Guide: Single-Patient Use Cuffs

Size	Circumference Range (cm)
Child	13.8-21.5
Small Adult	20.5-28.5
Adult	27.5-36.5
Adult Long	27.5-36.5
Large Adult	35.5-46.0
Large Adult Long	35.5-46.0

Long sizes have an additional 13.5 centimeter in length

SINGLE-PATIENT USE BLOOD PRESSURE CUFFS	TUBE		CONNECTOR			SIZE					
	ONE	TWO	BAYONET/RECTUS	SCREW	FEMALE SUBMINIATURE-LOCKING	CHILD	SMALL ADULT	ADULT	ADULT LONG	ADULT LARGE	LARGE ADULT LONG
1314M39100	•		•			•					
1314M39085	•		•				•				
1314M39048	•		•					•			
1314M39019	•		•						•		
1314M39230	•		•							•	
1314M39060	•		•								•
1314M39091	•			•		•					
1314M39092	•			•			•				
1314M39073	•			•				•			
1314M39093	•			•					•		
1314M39074	•			•						•	
1314M39075	•			•							•
1314M39011		•		•		•					
1314M39012		•		•			•				
1314M39023		•		•				•			
1314M39053		•		•					•		
1314M39014		•		•						•	
1314M39071		•		•							•
1314M39249		•			•	•					
1314M39213		•			•		•				
1314M39268		•			•			•			
1314M39081		•			•				•		
1314M39149		•			•					•	
1314M39018		•			•						•

BLOOD PRESSURE

**Sphygmomanometers
Aneroid Deluxe
(Latex-Free)**

**NOT MADE WITH
natural rubber latex** **diverse
supplier**



- Inflation bulb fitted with filter-screen to reduce dust build-up and ensure precise deflation control and reliable service
- Dip-molded inflation bladder withstands heavy use
- Navy-blue nylon cuff with range markings, artery label, gauge holder and hook-and-loop closure
- Black leatherette, zippered carrying case included
- Operating instructions included
- Gauge features:
 - Black enamel, 300mmHg, no-pin-stop, spring-loaded manometer
 - Measures blood pressure from 20mmHg to 300mmHg, certified accurate to +/-3mmHg
- Cuff Sizes:
 - Large Adult: 34.3 to 50.3 centimeters
 - Adult: 25.4 to 40.6 centimeters
 - Child: 18.4 to 26.7 centimeters
- Lifetime calibration warranty on manometer
- Mercury-free

Item No.	Description	Quantity
Sphygmomanometer Aneroid Blood Pressure Deluxe Navy Blue Nylon Cuff With Range Marking Mercury-Free		
994235	Child	1/EA
994205	Adult	1/EA
994215	Large Adult	1/EA

**Sphygmomanometers
Standard (Latex-Free)**

**NOT MADE WITH
natural rubber latex**



- Navy-blue nylon cuff with black bulb and bladder, size guide range markings to help prevent miss-cuffing, artery label, gauge holder, and hook-and-loop closure
- Measures blood pressure from 20mmHg to 300mmHg, certified accurate to +/-3mmHg
- Chrome-plated brass air release valve
- Cuff Sizes:
 - Child: 13 to 19.5 centimeters
 - Adult: 23 to 40 centimeters
 - Large Adult: 34 to 50 centimeters
- Operating instructions included
- Leatherette carrying case with durable nylon zipper
- Meets or exceeds ANSI/AAM/ISO 81060-1
- Mercury-free
- 20-year gauge warranty

Item No.	Description	Quantity
SPG1001	Child	1/EA
SPG1002	Adult	1/EA
SPG1003	Large Adult	1/EA

DOPPLER DEVICES & PROBES

**Doppler Systems
(Obstetrical)**

**made in the
USA** **NOT MADE WITH
natural rubber latex**



- Superior sound quality and sensitivity through new microchip technology
- Available with built-in audio recorder that holds 32 seconds of downloadable doppler sounds
- Units with displays have large screens showing extended heart rate from 50-220 bpm, signal quality, and more
- Interchangeable probes and user-replaceable cables
- Obstetrical probes available in 2 MHz, 2 MHz Waterproof, and 3 MHz
- Built-in battery design with extended life
- 2-year manufacturer's warranty

Item No.	Description	Quantity
Doppler Obstetrical Probe		
250ARSD2	Diagnostic With Display Recharger Audio Recorder 2 Mhz	1/EA
250ARSD3	Diagnostic With Display Recharger Audio Recorder 3 Mhz	1/EA
250ARSDW	Diagnostic With Display Recharger Audio Recorder 2 Mhz Waterproof	1/EA
M150ASD2	Diagnostic Non-Display Audio Recorder 2 Mhz	1/EA
M150ASD3	Diagnostic Non-Display Audio Recorder 3 Mhz	1/EA
M150ASDW	Diagnostic Non-Display Audio Recorder 2 Mhz Waterproof	1/EA
M150RSD2	Diagnostic Non-Display Recharger 2 Mhz	1/EA
M150RSD3	Diagnostic Non-Display Recharger 3 Mhz	1/EA
M150RSDW	Diagnostic Non-Display Recharger 2 Mhz Waterproof	1/EA
M250RSD2	Diagnostic With Display Recharger 2 Mhz	1/EA
M250RSD3	Diagnostic With Display Recharger 3 Mhz	1/EA
M250RSDW	Diagnostic With Display Recharger 2 Mhz Waterproof	1/EA
M350RSD2	Diagnostic Tabletop With Display Recharger 2 Mhz	1/EA
M350RSD3	Diagnostic Tabletop With Display Recharger 3 Mhz	1/EA
M350RSDW	Diagnostic Tabletop With Display Recharger 2 Mhz Waterproof	1/EA
MC150SD2	Diagnostic Non-Display 2 Mhz	1/EA
MC150SD3	Diagnostic Non-Display 3 Mhz	1/EA
MC150SDW	Diagnostic Non-Display 2 Mhz Waterproof	1/EA
MC250SD2	Diagnostic With Display 2 Mhz	1/EA
MC250SD3	Diagnostic With Display 3 Mhz	1/EA
MC250SDW	Diagnostic With Display 2 Mhz Waterproof	1/EA
MCSD2	2 Mhz For Deep Penetration Third Trimester Pregnancy	1/EA
MCSD3	3 Mhz For Early Fetal Heart Detection	1/EA
MCSDE	Transvaginal 5 Mhz Echo Heart	1/EA
MCSDW	Waterproof 2 Mhz Used In Labor And Delivery	1/EA

DOPPLER DEVICES & PROBES

Doppler Systems (Vascular)

made in the USA **NOT MADE WITH natural rubber latex**

- Superior sound quality and sensitivity through new microchip technology
- Built-in audio recorder that holds 32 seconds of doppler sounds, which are downloadable
- Units with displays have large screens showing extended heart rate from 50-220 bpm, signal quality, and more
- Interchangeable probes and user-replaceable cables
- Vascular probes available in 4 MHz, 5 MHz, and 8 MHz
- Built-in battery design with extended life
- 2-year manufacturer's warranty



Item No.	Description	Quantity
Doppler		
250ARSD4	Diagnostic With Display Recharger Audio Recorder 4 Mhz Vascular Probe	1/EA
250ARSD5	Diagnostic With Display Recharger Audio Recorder 5 Mhz Vascular Probe	1/EA
250ARSD8	Diagnostic With Display Recharger Audio Recorder 8 Mhz Vascular Probe	1/EA
M150ASD4	Diagnostic Non-Display Audio Recorder 4 Mhz Vascular Probe	1/EA
M150ASD5	Diagnostic Non-Display Audio Recorder 5 Mhz Vascular Probe	1/EA
M150ASD8	Diagnostic Non-Display Audio Recorder 8 Mhz Vascular Probe	1/EA
M150RSD4	Diagnostic Non-Display Recharger 4 Mhz Vascular Probe	1/EA
M150RSD5	Diagnostic Non-Display Recharger 5 Mhz Vascular Probe	1/EA
M150RSD8	Diagnostic Non-Display Recharger 8 Mhz Vascular Probe	1/EA
M250RSD4	Diagnostic With Display Recharger 4 Mhz Vascular Probe	1/EA
M250RSD5	Diagnostic With Display Recharger 5 Mhz Vascular Probe	1/EA
M250RSD8	Diagnostic With Display Recharger 8 Mhz Vascular Probe	1/EA
M350RSD4	Diagnostic Tabletop With Display Recharger 4 Mhz Vascular Probe	1/EA
M350RSD5	Diagnostic Tabletop With Display Recharger 5 Mhz Vascular Probe	1/EA
M350RSD8	Diagnostic Tabletop With Display Recharger 8 Mhz Vascular Probe	1/EA
MC150SD4	Diagnostic Non-Display 4 Mhz Vascular Probe	1/EA
MC150SD5	Diagnostic 5 Mhz Vascular Probe With Non Display	1/EA
MC150SD8	Diagnostic 8 Mhz Vascular Probe With Non Display	1/EA
MC250SD4	Diagnostic With Display 4 Mhz Vascular Probe	1/EA
MC250SD5	Diagnostic With Display 5 Mhz Vascular Probe	1/EA
MC250SD8	Diagnostic With Display 8 Mhz Vascular Probe	1/EA
MCK225	Stand For M150 and M250 Doppler	1/EA
MCSD4	Probe Vascular 4 Mhz To Ease Location Of Blood Vessel	1/EA
MCSD5	Probe Vascular 5 Mhz To Ease Location Of Blood Vessel	1/EA
MCSD8	Probe Vascular 8 Mhz To Ease Location Of Surface Vessel	1/EA

DOPPLER DEVICES & PROBES

Doppler Systems (ABI)

made in the USA **NOT MADE WITH natural rubber latex**

- The ABI doppler is configured primarily for conducting peripheral vascular exams including ABI exam and for documentation of blood flow
- Sensitive bi-directional dopplers with waveform printout that meets ABI reimbursement criteria
- Automatic calibration and fast waveform printout in 10 seconds
- Printout provided on label easily attaches to report form for reimbursement
- Fast, reimbursable diagnosis of Peripheral Arterial Disease (PAD)
- 1-year manufacturer's warranty



Item No.	Description	Quantity
Doppler Surgical Ankle Brachial Index Display Printer Bi-Directional 8 Mhz Probe		
MCL250AB	With Recharger	1/EA
MCL250AC	With Aneroid, Cuff, & Recharger	1/EA

Doppler Systems (ABI Comprehensive)

made in the USA **NOT MADE WITH natural rubber latex**

- The ABI doppler is configured primarily for conducting peripheral vascular exams including ABI exam and for documentation of blood flow
- Sensitive bi-directional dopplers with waveform printout that meets ABI reimbursement criteria
- Automatic calibration and fast waveform printout in 10 seconds
- Printout provided on label easily attaches to report form for reimbursement
- Fast, reimbursable diagnosis of Peripheral Arterial Disease (PAD)
- 1-year manufacturer's warranty



Item No.	Description	Quantity
MCL500VA	System Doppler Advanced Ankle Brachial Index Display Comprehensive With 9 Blood Pressure Cuffs Included	1/EA


EXAM AID ACCESSORIES

Applicators - General Purpose (Cotton-Tipped)

NOT MADE WITH
natural rubber latex

diverse
supplier

sterile R



- General purpose swab applicator
- Cotton tip with USP-grade cotton fiber
- Available with wood or plastic shaft
- Wood shafts are smooth and splinter-free
- Plastic shaft option made of smooth and sturdy polypropylene
- Fine-grain white birch wooden shaft is tasteless and odorless
- Available in sterile or non-sterile
- Sterile cotton-tipped applicators available individually packaged or 2/pack
- Sterile applicator package printed with ruler on edge for use in measuring size of wounds
- Non-sterile cotton-tipped applicators have autoclave-evident ink (lead-free) band printed on package to aid in identifying sterilized product
- Single-patient use

Item No.	Description	Quantity
Applicator Cotton Tip		
WOD1001	Wood Shaft 3 Inch Non-Sterile	1/EA, 100 EA/PK, 10 PK/BX, 10 BX/CA
WOD1002	Wood Shaft 6 Inch Non-Sterile	1/EA, 100 EA/PK, 10 PK/BX, 10 BX/CA
WOD1003	Wood Shaft 6 Inch Individually Wrapped Sterile	1/PK, 100 PK/BX, 10 BX/CA
WOD1004	Wood Shaft 6 Inch 2-Pack Sterile	2/PK, 100 PK/BX, 10 BX/CA
WOD1005	Plastic Shaft Plastic Polypropylene 6 Inch 2-Pack Sterile	2/PK, 100 PK/BX, 10 BX/CA
WOD2001	Wood Shaft 6 x 0.08 Inch Non-Sterile	1/EA, 1000 EA/BX, 20 BX/CA

Applicators - OB/GYN/Proctoscopic (Rayon-Tipped)

NOT MADE WITH
natural rubber latex

- Extra-absorbent, large rayon tip redesigned with tighter nap to prevent shedding fibers
- Polypropylene applicator sticks provide sturdy, yet flexible, shaft
- Medical-grade rayon fiber
- Hypoallergenic and DEHP-free
- Available in 16-Inch Proctoscopic or 8-Inch OB/GYN style
- OB/GYN applicator available in bulk package
- Autoclavable for use in sterile fields



Item No.	Description	Quantity
Applicator Rayon Tip Plastic Shaft Hypoallergenic DEHP-Free Polypropylene Non-Sterile		
OBP008	OB-GYN 8 Inch	1/EA, 50 EA/BX, 10 BX/CA
OBP008B	OB-GYN Bulk 8 Inch	1/EA, 50 EA/BG, 10 BG/CA
OBP016	Proctoscopic 16 Inch	1/EA, 50 EA/BX, 10 BX/CA

EXAM AID ACCESSORIES

Ultrasound Gel

NOT MADE WITH
natural rubber latex

diverse
supplier



- Bottles designed to feel soft and easy to squeeze; help to minimize repetitive motion that can lead to clinician fatigue
- Both premium and standard gels available in convenient 8 ounce squeeze bottle and 5 Liter container
- Bacteriostatic, water soluble, and perfume-free
- Contains no oils, alcohols, silicones, or fragrances that may harm the transducer
- Provides optimal sound velocity and acoustic impedance during ultrasound procedures
- Proprietary manufacturing method minimizes the presence of air bubbles in the gel for increased image clarity
- Sized to fit standard gel warmers
- Non-sterile, for external use only
- For use only with ultrasound equipment, not for use with automatic or external defibrillators
- Highest quality ultrasound gel made in North America
- Premium Ultrasound Gel
 - Higher viscosity recommended for vascular and echocardiography applications
 - Light glacial blue tint reduces the potential for staining equipment, clothing, and bedding while providing contrast against patient skin
 - Additional thickener inhibits the gel from drying and becoming tacky
 - 20g Individual Packets
 - Single patient packaging is designed to help prevent cross contamination from gel being used on more than one patient
 - Helps reduce the risk of nosocomial infection associated with cleaning, refilling and reusing gel bottles
 - Packets can be warmed in standard gel warmers
- Standard Ultrasound Gel
 - Lower viscosity recommended for abdominal and OB/GYN applications
 - Light aqua blue tint reduces the potential for staining equipment, clothing and bedding while providing contrast against patient skin
 - Non-slip formula makes testing fast and accurate
 - Lower viscosity recommended for general purpose, abdominal and OB/GYN procedures

Item No.	Description	Quantity
Gel Ultrasound		
M50008	Bottles Empty For Standard	1/EA, 12 EA/CA
M500812	Standard 8 Ounce Non-Sterile Fragrance Free Glacial Tint	1/EA, 12 EA/CA
M5005L4	Standard 5 Liter Non-Sterile Fragrance Free Glacial Tint	1/EA, 4 EA/CA
A11010	Bottles Empty For Premium	1/EA, 12 EA/CA
30400400	Premium 20 Gram Non-Sterile Single Patient Packet Fragrance Free Clear	1/EA, 100 EA/BG, 4 BG/CA
10080212	Premium 8 Ounce Non-Sterile Fragrance Free Glacial Tint	1/EA, 12 EA/CA
20040504	Premium 5 Liter Non-Sterile Fragrance Free Glacial Tent	1/EA

EXAM AID ACCESSORIES

Lubricating Jelly

NOT MADE WITH natural rubber latex **diverse supplier** **sterile R**

- Sterile multi-purpose lubricant for most routine medical procedures
- Provides instant lubrication of gloves and instruments for added patient comfort during procedures
- Clear, grease-less, water-soluble formula forms a clear, adhesive film that washes away easily
- Spreads easily over skin as well as most natural and synthetic materials
- Bacteriostatic
- Non-irritating and non-staining
- Alcohol-free
- Four sizes:
 - 2-ounce and 4-ounce tubes
 - 3-gram and 5-gram single-patient use foil packets



Item No.	Description	Quantity
Jelly Lubricating Sterile		
40ZLUB	4 Ounce Tube With Flip Top Lid	1/EA, 12 EA/BX, 6 BX/CA
20ZLUB	2 Ounce Tube With Flip Top Lid	1/EA, 12 EA/BX, 6 BX/CA
3GRAMLUB	3 Gram Packet Single Patient Use	1/EA, 144 EA/BX, 12 BX/CA
5GRAMLUB	5 Gram Packet Single-Patient Use	1/EA, 144 EA/BX, 12 BX/CA

Electrodes (Defibrillation)

NOT MADE WITH natural rubber latex

- Radiolucent to allow for easy, cost-effective use throughout the hospital in less clinically sensitive areas
- Compatible with biphasic defibrillators
 - Standard ISO 9001 (1994)
 - Size: 5.63 x 3.88 Inches
 - Canadian Medical Device Regulations
- Complies with ANSI/AAMI Standard DF39 (1993)
- Wire coating contains PVC
- Non-Sterile
- No adapters required—direct connect to defibrillator cable
- Functionally equivalent to various other brands of defibrillation electrodes
- Product made of vinyl laminated foam, conductive foil, and conductive adhesive hydrogel
- MediChoice adult defibrillation electrodes comply with the following medical device requirements:
 - United States 'Quality System Regulation' (21 CFR 820)
 - European Community's Council Directive 93/42/EEC (MDD) for CE marked products
 - Standard EN 46001 (1996)




Item No.	Description	Quantity
Electrode Defibrillation Electrocardiogram Vinyl Laminated Foam Conductive Foil Adhesive Hydrogel Adult		
MC1710H	Radiolucent Pacing Equal To Phillips Dp2 Amp Dp6 Amp M3713A	2/PK, 10 PK/CA
MC1310P	Non-Radiolucent Pacing Equal To Physio-Control Quik Combo	2/PK, 10 PK/CA
MC1410Z	Radiolucent Pacing Equal To Zoll Stat Padz	2/PK, 10 PK/CA

EXAM AID ACCESSORIES

Electrodes (Resting Tab)

NOT MADE WITH natural rubber latex

- Silver/Silver Chloride (Ag/AgCl) sensing system for stable, reliable electrocardiograph signals
- Conductive adhesive hydrogel for reliable tracings
- High tack adhesive hydrogel for optimal adhesion
- Uniquely formulated hydrogel minimizes residue for faster clean-up
- Non-Sterile



Item No.	Description	Quantity
MC5400	Electrode Electrocardiography Tab Resting Radiolucent	1/EA, 100 EA/BG, 10 BG/BX, 4 BX/CA

Vaginal Specula

NOT MADE WITH natural rubber latex

- Allows examination of vaginal mucosa at 360-degree angle without rotation
- VAG1000 series made from sturdy, clear non-toxic polystyrene
- VAG2000 and VAG3000 series are both made from clear, sturdy, non-toxic polymer blend of polycarbonate and polystyrene
- VAG1000 series has a silent slide and lock mechanism that is self-holding in various open positions and allows for one-handed operation
- Both the VAG2000 and VAG3000 series have built-in locking mechanisms for both elevation and angulation adjustment
- VAG2000 series is compatible with most corded light sources
- VAG3000 series features a fixed, built-in, LED light-source that is powered by 3 alkaline batteries
- VAG3000 series features on/off push-button operation of the light for extended battery life
- VAG3000 series has 30 minute battery life in continuous operation; can be turned on and off as desired during use for extended battery life
- All series provide a wide range of angular and elevational adjustability
- Disposable, single-patient use decreases the risk of contamination
- Non-sterile



Item No.	Description	Quantity
Speculum Vaginal Disposable Clear		
VAG1001	Small	1/EA, 10 EA/BG, 10 BG/CA
VAG1002	Medium	1/EA, 10 EA/BG, 10 BG/CA
VAG1003	Large	1/EA, 10 EA/BG, 10 BG/CA
VAG2001	Small Premium Adjustable Locking	1/EA, 25 EA/BX, 4 BX/CA
VAG2002	Medium Premium Adjustable Locking	1/EA, 25 EA/BX, 4 BX/CA
VAG2003	Large Premium Adjustable Locking	1/EA, 20 EA/BX, 5 BX/CA
VAG3001	With Built-In Led Light Source Small Premium Adjustable Locking	1/EA, 25 EA/BX, 4 BX/CA
VAG3002	With Built-In Led Light Source Medium Premium Adjustable Locking	1/EA, 25 EA/BX, 4 BX/CA
VAG3003	With Built-In Led Light Source Large Premium Adjustable Locking	1/EA, 20 EA/BX, 5 BX/CA

EXAM AID ACCESSORIES

Cotton Balls

NOT MADE WITH natural rubber latex **diverse supplier**

- 100-percent natural cotton
- Soft, absorbent, and uniform in size and weight
- Designed for a variety of uses throughout the hospital
- Packed in recloseable poly bags
- Available in Medium and Large sizes
- Non-sterile



Item No.	Description	Quantity
Ball Cotton Soft Absorbent Non-Sterile		
CTN200L	Large	200/BG, 10 BG/CA
CTN500M	Medium	500/BG, 8 BG/CA

PAPER PRODUCTS

Drape/Stretcher Sheets

made in the USA NOT MADE WITH natural rubber latex

- Durable full sheet construction
- Blankets patient in lightweight comfort
- Perfect for patient draping
- Offered in 2-ply or 3-ply tissue, or tissue/poly/tissue for an added moisture barrier
- Pebble embossed for added softness and comfort



Item No.	Description	Quantity
Sheet Drape		
MC918316	2 Ply Tissue Mauve 40 x 48 Inch	1/EA, 100 EA/CA
MC950453	Stretcher Non-Woven Repels Fluid Blue 40 x 90 Inch	1/EA, 50 EA/CA
M1810824	2 Ply Tissue White 40 x 48 Inch	1/EA, 100 EA/CA
MC18303	Stretcher 2 Ply Tissue White 40 x 60 Inch	1/EA, 100 EA/CA
MC18839	Stretcher Premium 3 Ply Tissue White 40 x 90 Inch	1/EA, 50 EA/CA
M1810837	3 Ply Tissue White 40 x 72 Inch	1/EA, 50 EA/CA
MC18309	Stretcher 3 Ply Tissue White 40 x 72 Inch	1/EA, 50 EA/CA
MC10455	Stretcher Tissue/Poly White 40 x 90 Inch	1/EA, 50 EA/CA
MC18302	2 Ply Tissue White 40 x 48 Inch	1/EA, 100 EA/CA
MC18310	Premium 3 Ply Tissue White 40 x 48 Inch	1/EA, 100 EA/CA
MC80929	Stretcher Tissue/Poly Blue 40 x 90 Inch	1/EA, 50 EA/CA

Please see our full line of Pediatric Treatment Tables on pages 80-81 in the Equipment Section.

PAPER PRODUCTS

Exam Table Paper

made in the USA NOT MADE WITH natural rubber latex

- Gentle patient comfort and durable protection for the exam table
- Available in sturdy, machine-glazed smooth finish or soft, absorbent crepe; poly-backed options available as well
- Available in 14, 18, 21 and 24-Inch widths
- Crepe finish has a softer feel and is more absorbent
- Smooth finish prevents fluid and moisture soak-through
- Non-toxic



Item No.	Description	Quantity
Paper Examination Table Non-Toxic		
MC918701	Pre-Cut Sheet Crepe Finish 18 Inch x 24 Inch	1/EA, 1000 EA/CA
MC186218	Poly-backed 30 Inch Perforated Crepe Finish 18 Inch x 100 Feet	1/RL, 12 RL/CA
MC916143	Premium Crepe Finish Pediatric 14 Inch x 125 Feet	1/RL, 12 RL/CA
MC913142	Premium Smooth 14 Inch x 225 Feet	1/RL, 12 RL/CA
MC11101	Crepe Finish 18 Inch x 125 Feet	1/RL, 12 RL/CA
MC18802	Premium Crepe Finish 18 Inch x 125 Feet	1/RL, 12 RL/CA
MC11102	Smooth Finish 18 Inch x 225 Feet	1/RL, 12 RL/CA
MC11206	Premium Smooth Finish 18 Inch x 225 Feet	1/RL, 12 RL/CA
MC11103	Crepe Finish 21 Inch x 125 Feet	1/RL, 12 RL/CA
MC18804	Premium Crepe Finish 21 Inch x 125 Feet	1/RL, 12 RL/CA
MC916021	Polybacked Crepe Finish Pediatric 21 Inch x 125 Feet	1/RL, 12 RL/CA
MC11104	Smooth Finish 21 Inch x 225 Feet	1/RL, 12 RL/CA
MC11208	Super Premium Smooth 21 Inch x 225 Feet	1/RL, 12 RL/CA
MC914243	Premium Crepe Finish 24 Inch x 125 Feet	1/RL, 12 RL/CA
MC911243	Premium Smooth Finish 24 Inch x 260 Feet	1/RL, 8 RL/CA
MC981418	Premium Crepe Finish Pediatric Printed Bugs & Things 18 Inch x 125 Feet	1/RL, 6 RL/CA
MC981618	Premium Crepe Finish Pediatric Printed Bugs & Things 18 Inch x 125 Feet	1/RL, 6 RL/CA
MC981518	Premium Smooth Finish Pediatric Printed Bugs & Things 18 Inch x 125 Feet	1/RL, 6 RL/CA

Professional Towels

made in the USA NOT MADE WITH natural rubber latex

- Ideal for washing, wiping and draping when absorbency and strength are needed
- Embossed for additional softness, fluid barrier, and strength
- Available in absorbent 2-ply and 3-ply tissue options, with and without polyethylene backing
- Deluxe Waffle Embossed 3-ply tissue with polyethylene backing also available
- Non-toxic



Item No.	Description	Quantity
Towel Professional Multi-purpose		
MC17461	2-Ply Poly Tissue Waffle Embossed White 13 x 18 Inch	1/EA, 500 EA/CA
MC18860	3-Ply Tissue Deluxe Waffle Embossed White 13 x 18 Inch	1/EA, 500 EA/CA
MC11111	3-Ply Poly Tissue Deluxe Waffle Embossed White 13 x 18 Inch	1/EA, 500 EA/CA
MC11109	3-Ply Tissue Diamond Embossed White 13 x 18 Inch	1/EA, 500 EA/CA
MC11110	3-Ply Poly Tissue Waffle Embossed Blue 13 x 18 Inch	1/EA, 500 EA/CA

PAPER PRODUCTS

Procedure Drapes (Sterile)



- Helps ensure clean environment for delicate medical procedures
- Tissue-Poly-Tissue construction provides absorption while preventing fluid pass-through
- Individual packaging protects sterile integrity and helps reduce waste
- Dispenser box helps keep procedure rooms organized
- Available with 3 Inch fenestration to allow access to procedure area



Item No.	Description	Quantity
Drape Procedure Towel Tissue/Poly/Tissue		
MDRP1001	Sterile White 18 x 26 Inch	1/EA, 50 EA/BX, 6 BX/CA
MDRP1002	3 Inch Fenestration Sterile White 18 x 26 Inch	1/EA, 50 EA/BX, 6 BX/CA

EKG Chart Paper



- Consistent trace quality with reliable and smooth printing
- Minimizes electrostatic buildup
- Universal red grid color
- 5-year shelf life of unimaged premium grade thermal paper when stored in original packaging, in darkness, below 77 degrees Fahrenheit (25 degrees Celsius) and 65-percent humidity
- Imaged premium grade thermal paper legible for 5 years
- Product warranty available



Item No.	Description	Quantity
Paper Chart Electrocardiology		
MC9420	3-Channel Z-Fold Thermal Red Grid 8.5 x 11 Inch	300/PD, 8 PD/CA
MC9424	3-Channel Z-Fold Red Grid With Header 216 x 280 Mm	300/PD, 8 PD/CA
MCPCLP12	3-Channel Red Graph 216 x 140 Mm	100/PD, 18 PD/CA
MCMF2428	1-Channel Blank Single Sheet 50 x 120 Mm	100/PK, 5 PK/CA
MCM1707A	3-Channel Z-Fold Red Grid 216 x 280 Mm	200/PD, 8 PD/CA
MCM1910A	Paper Fetal Monitor Electrocardiology Z-Fold Orange Grid 150 x 100 Mm	150/PD, 40 PD/CA
MC2485A	3-Channel Red Grid 216 x 280 Mm	1/EA, 8 EA/CA, 1/PD
MC40477	1-Channel No Grid 50 Mm x 100Feet	1/RL, 10 RL/PK, 10 PK/CA
MSL60196	1-Channel Z-Fold 106 x 120 Mm	230/PD, 100 PD/CA
MC18416	3-Channel Z-Fold Full Stress Red Grid 216 x 280 x 218 Mm	218/PD, 20 PD/CA
MC1100HD	Paper Print Video 20 x 110 Mm	1/RL, 5 RL/BX, 8 BX/CA
MC008038	1-Channel Blank 50 Mm x 100Feet	1/RL, 10 RL/PK, 10 PK/CA
M0007983	Paper Chart Electrocardiology 3-Channel Z-Fold 216 x 280 Mm	200/PD, 10 PD/CA
M0997984	Paper Chart Electrocardiology 3-Channel Z-Fold Red Graph 216 x 280 Mm	200/PD, 10 PD/CA
M0997989	3-Channel Z-Fold Red Graph 216 x 280 Mm	300/PD, 10 PD/CA
M8000300	Recording 3-Channel Fan Fold Red Grid 90 x 90 Mm	200/PD, 50 PD/CA
M9402053	Paper Chart Electrocardiology Recording 3Channel Fan Fold Laser Print Red Grid 260 x 280 Millimeter	500/PD, 5 PD/CA
M997988	3-Channel Z-Fold Blank Red Grid 216 x 140 Mm	300/PD, 10 PD/CA
M4305AAO	Paper Fetal Monitor Electrocardiology Z-Fold Red Grid 152 x 90 Mm	150/PD, 40 PD/CA
M4305CAO	Paper Fetal Monitor Electrocardiology Z-Fold Red Grid With Adhesive Label 152 x 90 Mm	150/PD, 40 PD/CA

STETHOSCOPES

Stethoscopes Disposable (Plastic)



- Lightweight plastic, diaphragm-only blood pressure stethoscope with plastic binaural
- Ideal for single-patient use to prevent cross-contamination and aid with infection control protocol compliance
- PVC tubing in high-visibility yellow or cobalt blue
- Overall length: 30 Inches
- Diaphragm diameter: 1.75 Inches
- Weight: 1.5 ounces
- Individually packaged in a sealed plastic bag
- Disposable



Item No.	Description	Quantity
Stethoscope Blood Pressure Plastic Binaural 1.75 Inch Diaphragm 1.5 Ounce 30 Inch Disposable Single-Patient Use		
994910	Blue	1/EA, 10 EA/BX, 30 BX/CA
994930	Yellow	1/EA, 10 EA/BX, 30 BX/CA

Stethoscopes Disposable (Chrome)



- Lightweight plastic, diaphragm-only blood pressure stethoscope with chrome-plated binaural
- Ideal for single-patient use to prevent cross-contamination and aid with infection control protocol compliance
- PVC tubing in high-visibility yellow or black
- Individually packaged in a sealed plastic bag
- Diaphragm diameter: 1.75 Inches
- Overall length: 30 Inches
- Weight: 2.5 ounces
- Disposable



Item No.	Description	Quantity
Stethoscope Cardiology Disposable Chrome Binaural Lightweight Polyvinyl Chloride Tubing 30 Inch Single-Patient Use		
994700	Black	1/EA, 10 EA/BX, 30 BX/CA
994730	Yellow	1/EA, 10 EA/BX, 30 BX/CA

STETHOSCOPES

Stethoscopes
Dual-Head (Neo-Natal)

NOT MADE WITH
natural rubber latex

diverse
supplier

- Chrome-plated combination chest piece with 20 mm (0.79 Inch) bell and 18 mm (0.71 Inch) diaphragm
- Slip-on, non-chill bell ring
- Adjustable chrome-plated binaural with PVC ear tips
- Flexible black PVC tubing
- Overall length: 30 Inches
- Weight: 4 ounces



Item No.	Description	Quantity
995100	Stethoscope Neonatal Dual-Head Diagnostic 0.98 Inch Diaphragm 4 Ounce Black Pvc Tubing 30 Inch	1/EA, 30 EA/CA

STETHOSCOPES

Stethoscopes
Single-Head (Bowles)

NOT MADE WITH
natural rubber latex

diverse
supplier

- Chrome-plated zinc, diaphragm-only chest piece with traditional, raised-stem design
- Adjustable, chrome-plated brass binaural with PVC ear tips
- Flexible, thick-walled PVC tubing minimizes background interference; black
- Diaphragm diameter: 1.625 Inch
- Overall length: 30 Inches
- Weight: 4 ounces



Item No.	Description	Quantity
994800	Stethoscope Bowles Diagnostic 1.63 Inch Diaphragm 4 Ounce Black Pvc Tubing 30 Inch	1/EA, 30 EA/CA

Stethoscopes
Dual-Head (Infant)

NOT MADE WITH
natural rubber latex

diverse
supplier

- Chrome-plated combination chest-piece with 20mm diaphragm and 18mm bell
- Slip on gray non-chill bell ring for patient comfort
- Adjustable chrome-plated binaural with PVC ear tips
- PVC tubing in pink
- Overall length: 32 Inches
- Weight: 4 ounces



Item No.	Description	Quantity
995250	Stethoscope Cardiology Infant Dual-Head Chrome-Plated Binaural PVC Tubing Pink 32 Inch	1/EA, 30 EA/CA

Stethoscopes
Single-Head
(Light-Weight)

NOT MADE WITH
natural rubber latex

diverse
supplier

- Lightweight aluminum, diaphragm-only stethoscope with nurse-style chest piece
- Adjustable chrome-plated binaural with PVC ear tips
- Flexible PVC tubing in black, red, blue, green, gray, gold, yellow, light blue, teal
- Diaphragm diameter: 1.75 Inches
- Overall length: 30 Inches
- Weight: 4 ounces



Item No.	Description	Quantity
Stethoscope Lightweight Single-Head Cardiology Diaphragm-Only PVC Tubing 30 Inch		
994600	Black	1/EA, 30 EA/CA
994605	Red	1/EA, 30 EA/CA
994610	Blue	1/EA, 30 EA/CA
994615	Green	1/EA, 30 EA/CA
994620	Gray	1/EA, 30 EA/CA
994625	Gold	1/EA, 30 EA/CA
994630	Yellow	1/EA, 30 EA/CA
994635	Teal	1/EA, 30 EA/CA
994640	Light Blue	1/EA, 30 EA/CA

Stethoscopes
Dual-Head (Pediatric)

NOT MADE WITH
natural rubber latex

diverse
supplier

- Lightweight anodized aluminum combination chest piece has a 39 mm (1.5 Inch) bell with non-chill black rim and 26.5 mm (1.0 Inch) diaphragm
- Adjustable chrome-plated brass binaural with PVC ear tips
- Flexible black PVC tubing
- Overall length: 30 Inches
- Weight: 4 ounces



Item No.	Description	Quantity
Stethoscope Dual-Head Pediatric Diagnostic 1.5 Ounce PVC Tubing 30 Inch		
995000	Black	1/EA, 30 EA/CA
995020	Gray	1/EA, 30 EA/CA

STETHOSCOPES



Stethoscopes Dual-Head (Light-Weight)

NOT MADE WITH natural rubber latex **diverse supplier**

- Lightweight, anodized aluminum combination dual-head stethoscope with wide frequency range detection
- Adjustable chrome-plated brass binaural with PVC ear tips
- Flexible PVC tubing in black, red, blue, green, gray, yellow, teal, and pink
- Overall length: 30 Inches; diaphragm diameter: 1.75 Inches; bell diameter 1.85 Inches
- Weight: 4 ounces

Item No.	Description	Quantity
Stethoscope Diagnostic Dual-Head 1.75 Inch Diaphragm 4 Ounce PVC Tubing 30 Inch		
994300	Black	1/EA, 30 EA/CA
994305	Red	1/EA, 30 EA/CA
994310	Blue	1/EA, 30 EA/CA
994315	Green	1/EA, 30 EA/CA
994320	Gray	1/EA, 30 EA/CA
994330	Yellow	1/EA, 30 EA/CA
994335	Teal	1/EA, 30 EA/CA
994350	Pink	1/EA, 30 EA/CA

STETHOSCOPES



Stethoscopes Dual-Head (Adult/Pediatric Professional)

NOT MADE WITH natural rubber latex **diverse supplier**

- Full-range stainless stethoscope with unmatched comfort & superior acoustics
- Comfort sealing eartips and light weight design
- Superior acoustics allow for accurate auscultation of cardiac, pulmonary and vascular sounds
- Stainless steel binaural and thick-walled tubes reduce noise and provide outstanding sound transmission
- Diaphragm is excellent for all sounds and murmurs
- Best for higher frequencies and lung sounds
- Ideal in identifying splitting of sounds, systolic clicks, and ejection sounds
- Bell provides excellent detection and clarity for low and medium frequency heart sounds and murmurs such as a faint gallop or diastolic rumble
- Includes dual non-chill rings, 2 extra pair of soft eartips (large and small), and 1 extra diaphragm
- Flexible PVC tubing in black, burgundy, and navy (adult); Pediatric—black only; Infant—gray only
- Overall length: 28 Inches
- Weight: 130 grams
- Bell and Diaphragm Diameters:
 - Adult: 4.75 cm (bell) and 3.5 cm (diaphragm)
 - Pediatric: 3.9 cm (bell) and 3.0 cm (diaphragm)
 - Infant: 3.1 cm (bell) and 2.25 cm (diaphragm)
- Fully assembled

Item No.	Description	Quantity
Select Dual-Head Stainless Steel Binaural PVC Tubing 28 Inch		
700102	Infant Gray	1/EA
14579	Pediatric Black	1/EA
13579	Adult Black	1/EA
13779	Adult Navy	1/EA
13979	Adult Burgundy	1/EA

STETHOSCOPES

Stethoscopes Dual-Head
(Cardiology Professional)

NOT MADE WITH
natural rubber latex

diverse
supplier

- Full-range cardiology and pulmonology stethoscope
- Can identify extremely low frequencies such as a faint diastolic gallop or rumble, as well as the higher frequencies of the faint early blowing aortic diastolic murmur
- Weight: 155 gm
- Fully assembled
- Diaphragm is excellent for all sounds and murmurs; best for higher frequencies such as faint diastolic murmur of aortic regurgitation
- Ideal in identifying splitting of sounds, systolic clicks and ejection sounds
- Diaphragm is designed to accentuate sounds in the 200 to 500 Hz range where many significant heart sounds occur
- Stainless steel binaural and thick-walled tubes reduce noise and provide outstanding sound transmission
- Bell provides excellent detection and clarity for low and medium frequency heart sounds and murmurs (such as a faint gallop or diastolic rumble)
- Flexible black PVC tubing
- Bell: 4.75 centimeter; Diaphragm: 3.65 centimeter
- Overall length: 25 Inches



Item No.	Description	Quantity
012579	Stethoscope Cardiology Elite Dual-Head Stainless Steel Binaural Black PVC Tubing 25 Inch	1/EA

Stethoscopes Dual-Head
(Sprague-Rappaport)

NOT MADE WITH
natural rubber latex

diverse
supplier

- Threaded chest piece drum of chrome-plated zinc with interchangeable fittings and other accessories for complete diagnostic capability
- Double-tube configuration for better acoustic separation
- Flexible PVC tubing available in black or neon yellow
- Tubing length: 22 Inches
- Overall length: 30 Inches
- Weight: 7.5 ounces
- Fully pre-assembled
- Valve mechanism eliminates acoustic leakage
- Adjustable chrome-plated brass binaural with PVC ear tips
- Accessory kit includes:
 - 1 Adult and 1 Pediatric transparent diaphragms
 - 3 bells (1.25 Inch, 1.0 Inch, and 0.75 Inch)
 - 1 pair mushroom ear tips
 - 1 pair silicon ear tips
 - 1 chest-piece key for easy diaphragm replacement



Item No.	Description	Quantity
Stethoscope Cardiology Sprague-Rappaport With Accessory Kit Pvc Tubing 22 Inch		
994510	Black	1/EA, 30 EA/CA
994535	Neon Yellow	1/EA, 30 EA/CA

SPHYGMOMANOMETERS

Mobile Aneroid
Sphygmomanometers

NOT MADE WITH
natural rubber latex

diverse
supplier

- Height adjustable 5-leg caster stand with range of 40 - 54 Inches from bottom of base to top of unit
- Integral ABS cuff storage compartment
- Manometer angled back 15-degrees for easy viewing
- Large, easy-to-read 6-Inch luminescent dial for high visibility in most lighting conditions
- 300mmHg no-pin stop mechanism withstands the heaviest use
- Inflation system includes black nylon cuff, inflation bladder, bulb and valve, and 8 foot length coiled tubing (when extended)
- Resting basket for cuff when not in use
- Deluxe bulb with filter-screen to reduce dust build-up
- Complete with necessary hardware and assembly instructions
- 2-year inflation system manufacturer's warranty



Item No.	Description	Quantity
Sphygmomanometer Aneroid Mobile		
BPUM3001	Child	1/EA
BPUM4001	Adult	1/EA
BPUM5001	Large Adult	1/EA

Complementary Products:

Item No.	Description
IVWB51B	Basket IV Stand Wire Heavy-Duty Large Epoxy-Coated Steel
IVWB51S	Basket IV Stand Wire Heavy-Duty Large Stainless Steel
IVWB52B	Basket IV Stand Wire Heavy-Duty Small Epoxy-Coated Steel
IVWB52S	Basket IV Stand Wire Heavy-Duty Small Stainless Steel

Wall Mount Aneroid
Sphygmomanometers

NOT MADE WITH
natural rubber latex

diverse
supplier

- Wall bracket swivels 110-degree for easy viewing
- Large, easy-to-read 6-Inch luminescent dial for high visibility in most lighting conditions
- 300 mmHg no-pin stop mechanism withstands the heaviest use
- Inflation system includes black nylon cuff, inflation bladder, bulb and valve, and 8-foot length coiled tubing (when extended)
- Deluxe bulb fitted with filter screen to reduce dust build-up
- Complete with all necessary hardware and assembly instructions
- 2-year inflation system manufacturer's warranty



Item No.	Description	Quantity
Sphygmomanometer Aneroid Wall Mount Blood Pressure Aneroid		
BPUW4002	Adult	1/EA
BPUW5002	Large Adult	1/EA

GLOVES



In today's competitive healthcare climate, medical professionals need to cut costs wherever possible, and high-volume commodity items like gloves are a great place to look for savings. However, careless cost cutting can lead to lower quality and higher incidents of failure. MediChoice exam and high-risk gloves, offer significant savings without sacrificing quality—or comfort.

Glove Selection Guide

Gloves come in many varieties, and each one has a job is best suited for a specific job.

Vinyl exam gloves are a good choice for simpler, low-risk tasks where exposure to blood or other bodily fluids is unlikely. Examples include administering injections, routine suctioning, basic patient examinations, food handling and low-risk housekeeping activities. Vinyl gloves also come in a stretchy version with added polymers that make them feel similar to latex.

Nitrile gloves are ideal for moderate to high-risk situations where exposure to blood, bodily fluids, and other potentially infectious material poses a risk to healthcare workers. Like latex, nitrile gloves form tightly to the hand and can be used for precision tasks without the exposure to latex proteins, which for some clinicians and patients can cause allergic reactions.

Latex, or natural rubber gloves, provides a high level of barrier protection and performs well in moderate to high-risk situations where there is potential for touching blood, bodily fluids, secretions, excretions and items visibly soiled by bodily fluids. They are ideal for wearers that change gloves frequently, wash hands repeatedly, or wear gloves for an extended period of time.

Latex Exam Gloves

- Non-sterile, ambidextrous, powder-free, single-use
- Beaded cuff for added strength and ease in donning
- Textured finish provide strong grip-strength for wet and dry applications
- Contains 50 micrograms or less of total water extractable protein
- Caution: This product contains natural rubber latex which may cause allergic reactions
- Safe use of this glove by or on latex-sensitive individuals has not been established



Item No.	Description	Quantity
Glove Exam Latex 9.5-Inch Cuff Non-Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Contains Natural Rubber Latex		
GLV1000	Extra-Small	1/EA, 100 EA/BX, 10 BX/CA
GLV1001	Small	1/EA, 100 EA/BX, 10 BX/CA
GLV1002	Medium	1/EA, 100 EA/BX, 10 BX/CA
GLV1003	Large	1/EA, 100 EA/BX, 10 BX/CA
GLV1004	Extra-Large	1/EA, 100 EA/BX, 10 BX/CA

NITRILE EXAM GLOVES

ETS Nitrile Exam Gloves

NOT MADE WITH
natural rubber latex

- Fingertip thickness: 6.2mil
- Weight: 4.8g
- Cleared for use with some chemotherapy drugs
- Non-sterile, ambidextrous, powder-free, single-use
- Beaded cuff for added strength and ease in donning
- Textured fingertips provide strong grip-strength for wet and dry applications
- Superior puncture, tear, and chemical resistance with excellent tactile sensitive and dexterity
- Thermal-sensitive compound conforms to hand
- Low-pressure design helps reduce user fatigue
- Made without natural rubber latex for reduced risk of allergic reactions



Item No.	Description	Quantity
Glove Exam Nitrile ETS 9.5-Inch Cuff Non-Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Cleared For Use With Some Chemotherapy Drugs		
GLV2000	Extra-Small	1/EA, 150 EA/BX, 10 BX/CA
GLV2001	Small	1/EA, 150 EA/BX, 10 BX/CA
GLV2002	Medium	1/EA, 150 EA/BX, 10 BX/CA
GLV2003	Large	1/EA, 150 EA/BX, 10 BX/CA
GLV2004	Extra-Large	1/EA, 140 EA/BX, 10 BX/CA

Please see our full line of
Glove Box Holders on page 31
in the Central Supply Section.

NITRILE EXAM GLOVES

Next Generation Nitrile Exam Gloves

NOT MADE WITH natural rubber latex



- Fingertip thickness: 4.7mil
- Weight: 3.5g
- Our most lightweight and flexible nitrile glove
- Cleared for use with some chemotherapy drugs
- Non-sterile, ambidextrous, powder-free, single-use
- Beaded cuff for added strength and ease in donning
- Textured fingertips provide strong grip-strength for wet and dry applications
- Superior puncture, tear, and chemical resistance with excellent tactile sensitivity and dexterity
- Thermal-sensitive compound conforms to hand
- Low-pressure design helps reduce user fatigue
- Made without natural rubber latex for reduced risk of allergic reactions

Item No.	Description	Quantity
----------	-------------	----------

Glove Exam Nitrile XT Next Generation 9.5-Inch Cuff Non-Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Cleared For Use With Some Chemotherapy Drugs

GLV2300	Extra-Small	1/EA, 200 EA/BX, 10 BX/CA
GLV2301	Small	1/EA, 200 EA/BX, 10 BX/CA
GLV2302	Medium	1/EA, 200 EA/BX, 10 BX/CA
GLV2303	Large	1/EA, 200 EA/BX, 10 BX/CA
GLV2304	Extra-Large	1/EA, 180 EA/BX, 10 BX/CA

NITRILE EXAM GLOVES

Stretchy Nitrile Exam Gloves

NOT MADE WITH natural rubber latex



- Fingertip thickness: 5.9mil
- Weight: 4.6g/6.0g
- Available in 9.5-inch and extended 12-inch cuff lengths
- Recommended for administering and/or mixing chemotherapy drugs
- Cleared for use with some chemotherapy drugs
- Non-sterile, ambidextrous, powder-free, single-use
- Beaded cuff for added strength and ease in donning
- Textured fingertips provide strong grip-strength for wet and dry applications
- Superior puncture, tear, and chemical resistance with excellent tactile sensitivity and dexterity
- Thermal-sensitive compound conforms to hand
- Low-pressure design helps reduce user fatigue
- Made without natural rubber latex for reduced risk of allergic reactions

Item No.	Description	Quantity
----------	-------------	----------

Glove Exam Nitrile Non-Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Cleared For Use With Some Chemotherapy Drugs

GLV5000	Extra-Small	1/EA, 100 EA/BX, 10 BX/CA
GLV5001	Small	1/EA, 100 EA/BX, 10 BX/CA
GLV5002	Medium	1/EA, 100 EA/BX, 10 BX/CA
GLV5003	Large	1/EA, 100 EA/BX, 10 BX/CA
GLV5004	Extra-Large	1/EA, 100 EA/BX, 10 BX/CA
GLV7000	Extra-Small 12-Inch Cuff	1/EA, 100 EA/BX, 10 BX/CA
GLV7001	Small 12-Inch Cuff	1/EA, 100 EA/BX, 10 BX/CA
GLV7002	Medium 12-Inch Cuff	1/EA, 100 EA/BX, 10 BX/CA
GLV7003	Large 12-Inch Cuff	1/EA, 100 EA/BX, 10 BX/CA
GLV7004	Extra-Large 12-Inch Cuff	1/EA, 100 EA/BX, 10 BX/CA

VINYL EXAM GLOVES

Vinyl Exam Gloves

NOT MADE WITH natural rubber latex



- Economical alternative for low-risk exposure and tasks
- Non-sterile, ambidextrous, powder-free, single-use
- Beaded cuff for added strength and ease in donning
- Smooth finish provides excellent tactile sensitivity and a comfortable fit
- Made without natural rubber latex for reduced risk of allergic reactions

Item No.	Description	Quantity
----------	-------------	----------

Glove Exam Vinyl 9.5-Inch Cuff Non-Sterile Ambidextrous Powder-Free Beaded Cuff Disposable

GLV3001A	Small	1/EA, 150 EA/BX, 10 BX/CA
GLV3002A	Medium	1/EA, 150 EA/BX, 10 BX/CA
GLV3003A	Large	1/EA, 150 EA/BX, 10 BX/CA
GLV3004A	Extra-Large	1/EA, 130 EA/BX, 10 BX/CA

VINYL EXAM GLOVES

Stretch Vinyl Exam Gloves

NOT MADE WITH natural rubber latex



- Unique polyvinyl is softer and more supple than regular vinyl
- Non-sterile, ambidextrous, powder-free, single-use
- Beaded cuff for added strength and ease in donning
- Smooth finish provides excellent tactile sensitivity and a comfortable fit
- Made without natural rubber latex for reduced risk of allergic reactions

Item No.	Description	Quantity
----------	-------------	----------

Glove Exam Stretch Vinyl 9.5-Inch Cuff Non-Sterile Ambidextrous Powder-Free Beaded Cuff Disposable

GLV4001A	Small	1/EA, 150 EA/BX, 10 BX/CA
GLV4002A	Medium	1/EA, 150 EA/BX, 10 BX/CA
GLV4003A	Large	1/EA, 150 EA/BX, 10 BX/CA
GLV4004A	Extra-Large	1/EA, 130 EA/BX, 10 BX/CA

Sterile Procedure Gloves



- Available in nitrile singles, latex singles and latex pairs
- Sterile ambidextrous, powder-free, single-use
- Packaged in a peel pouch to maintain sterility
- Beaded cuff for added strength and ease in donning
- Textured finish provides strong grip-strength for wet and dry applications
- Caution: The latex gloves contain natural rubber latex which may cause allergic reactions
- Safe use of the latex gloves by or on latex-sensitive individuals has not been established



Item No.	Description	Quantity
Glove Exam Procedure Latex 9.5-Inch Cuff Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Contains Natural Rubber Latex		
MC147214	Small Latex Singles	1/EA, 100 EA/BX, 4 BX/CA
MC147216	Medium Latex Singles	1/EA, 100 EA/BX, 4 BX/CA
MC147218	Large Latex Singles	1/EA, 100 EA/BX, 4 BX/CA
MC147224	Small Latex Pairs	1/PR, 50 PR/BX, 4 BX/CA
MC147226	Medium Latex Pairs	1/PR, 50 PR/BX, 4 BX/CA
MC147228	Large Latex Pairs	1/PR, 50 PR/BX, 4 BX/CA

Item No.	Description	Quantity
Glove Exam Procedure Nitrile 9.5-Inch Cuff Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Not Made With Natural Rubber Latex		
GLV5001S	Small Nitrile Singles	1/EA, 100 EA/BX, 4 BX/CA
GLV5002S	Medium Nitrile Singles	1/EA, 100 EA/BX, 4 BX/CA
GLV5003S	Large Nitrile Singles	1/EA, 100 EA/BX, 4 BX/CA

Complementary Products

Item No.	Description
GP-015	Dispenser Triple Glove Box Clear Polyethylene Terephthalate Plastic 15.75 x 10 x 3.75 Inch Vertical or Horizontal Wall Mounting
GL100-0212	Dispenser Single Glove Box Opaque Quartz Injection Molded ABS Plastic 5.8125 x 10.125 x 4.25 Inch Vertical or Horizontal Wall Mounting
GL100-1214	Dispenser Single Glove Box Semi-Transparent Injection Molded Polycarbonate Plastic 5.8125 x 10.125 x 4.25 Inch Vertical or Horizontal Wall Mounting
GL300-0212	Dispenser Triple Glove Box With Dividers Opaque Quartz Injection Molded ABS Plastic 17.125 x 10.125 x 4.25 Inch Vertical or Horizontal Wall Mounting
GL300-1214	Dispenser Triple Glove Box With Dividers Semi-Transparent Injection Molded Polycarbonate Plastic 17.125 x 10.125 x 4.25 Inch Vertical or Horizontal Wall Mounting

High-Risk Gloves

- Available in 8mil nitrile, 15mil latex and 18mil latex
- 8mil nitrile is appropriate for instrument decontamination and medical waste cleanup
- 15mil latex is an excellent choice for EMS and First Responders
- 18mil latex is ideal for chemotherapy compounding and other related chemotherapy activities
- Extended cuff offers greater protection
- Non-sterile, powder-free, single-use
- Beaded cuff for added strength and ease in donning
- Textured finish provides strong grip-strength for wet and dry applications
- Exceeds ASTM standard for examination gloves for medical and high-risk applications
- Caution: The latex gloves contains natural rubber latex which may cause allergic reactions
- Safe use of the latex gloves by or on latex-sensitive individuals has not been established



Item No.	Description	Quantity
Glove Exam High-Risk Nitrile 8 Mil 12-Inch Cuff Non-Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Cleared For Use With Some Chemotherapy Drugs Not Made With Natural Rubber Latex		
HRA8S	Small Nitrile 8 MM	1/EA, 50 EA/BX, 8 BX/CA
HRA8M	Medium Nitrile 8 MM	1/EA, 50 EA/BX, 8 BX/CA
HRA8L	Large Nitrile 8 MM	1/EA, 50 EA/BX, 8 BX/CA
HRA8XL	Extra-Large Nitrile 8 MM	1/EA, 50 EA/BX, 8 BX/CA

Item No.	Description	Quantity
Glove Exam High-Risk Latex 15 Mil 12-Inch Cuff Non-Sterile Ambidextrous Powder-Free Textured Beaded Cuff Disposable Contains Natural Rubber Latex		
FRA15S	Small Latex 15 MM	1/EA, 50 EA/BX, 10 BX/CA
FRA15M	Medium Latex 15 MM	1/EA, 50 EA/BX, 10 BX/CA
FRA15L	Large Latex 15 MM	1/EA, 50 EA/BX, 10 BX/CA
FRA15XL	Extra-Large Latex 15 MM	1/EA, 50 EA/BX, 10 BX/CA

Item No.	Description	Quantity
Glove Exam High-Risk Latex 18 Mil 12-Inch Cuff Non-Sterile Hand-Specific Powder-Free Textured Beaded Cuff Disposable Cleared For Use With Some Chemotherapy Drugs Contains Natural Rubber Latex		
HSR18S	Small Latex 18 MM	1/PR, 50 PR/BX, 10 BX/CA
HSR18M	Medium Latex 18 MM	1/PR, 50 PR/BX, 10 BX/CA
HSR18L	Large Latex 18 MM	1/PR, 50 PR/BX, 10 BX/CA
HSR18XL	Extra-Large Latex 18 MM	1/PR, 50 PR/BX, 10 BX/CA

Recloseable Bags

NOT MADE WITH
natural rubber latex

- Strong non-recycled polyethylene bag
- Convenient zip-style resealable closure
- Used for collection and transport of patient items
- Available in a wide range of sizes

Item No.	Description	Quantity
Bag Resealable Zip Closure Polyethylene Clear		
ZB03052	3 × 5 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/BX, 10 BX/CA
ZB04062	4 × 6 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/BX, 6 BX/CA
ZB05082	5 × 8 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/BX, 4 BX/CA
ZB06092	6 × 9 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/BX, 4 BX/CA
ZB08102	8 × 10 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/BX, 3 BX/CA
ZB10122	10 × 12 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/CA
ZB10132	10 × 13 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/CA
ZB12152	12 × 15 Inch 2 Mil	1/EA, 100 EA/BG, 10 BG/CA
ZB12154	12 × 15 Inch 4 Mil	1/EA, 100 EA/BG, 5 BG/CA
ZB13184	13 × 18 Inch 4 Mil	1/EA, 100 EA/BG, 5 BG/CA



Heavy-Duty Biohazard Waste Bags

made in the USA **NOT MADE WITH**
natural rubber latex

- Meets federal and state regulations for tear resistance and impact resistance
 - 480 gram Elmendorf tear test (ASTM D1922-94a)
 - 165 gram Dart Impact Resistance Test (ASTM D1709-97)
- Bilingual print (English and Spanish)
- Made from high quality resins



Heavy-Duty Biohazard Waste Bags

Item No.	Type	Capacity	Size	Gauge	Color	Packaging	EA/Case
A4823PRO	Healthcare Printed	8-10 gal.	24 × 23	1.30 mil	Red	FlatPack	500
A4830PRO	Healthcare Printed	8-10 gal.	24 × 30	1.30 mil	Red	FlatPack	500
A4830ZRO	Healthcare Printed	8-10 gal.	24 × 30	3.00 mil	Red	FlatPack	200
A4832PRO	Healthcare Printed	12-16 gal.	24 × 32	1.30 mil	Red	FlatPack	500
B6036HRO	Healthcare Non-Printed	20-30 gal.	30 × 36	0.80 mil	Red	FlatPack	250
A6036PRO	Healthcare Printed	20-30 gal.	30 × 36	1.30 mil	Red	FlatPack	250
A6043PRO	Healthcare Printed	20-30 gal.	30 × 43	1.30 mil	Red	FlatPack	200
A6043PROR01	Healthcare Printed	20-30 gal.	30 × 43	1.30 mil	Red	CorelessRoll	8/25
A6043ARO	Healthcare Printed	20-30 gal.	30 × 43	1.60 mil	Red	FlatPack	200
A6043ZRO	Healthcare Printed	20-30 gal.	30 × 43	3.00 mil	Red	FlatPack	100
Z6143HROP01	Healthcare Printed	20-30 gal.	31 × 43	14.00 mil	Red	FlatPack	250
A6639PRO	Healthcare Printed	33 gal.	33 × 39	1.30 mil	Red	FlatPack	150
A6639ZRO	Healthcare Printed	33 gal.	33 × 39	3.00 mil	Red	FlatPack	100
Z6640HROP01	Healthcare Printed	33 gal.	33 × 40	14.00 mil	Red	FlatPack	250
A8046PRO	Healthcare Printed	40-45 gal.	40 × 46	1.30 mil	Red	FlatPack	200
A8046QRO	Healthcare Printed	40-45 gal.	40 × 46	2.00 mil	Red	FlatPack	100
A8046ZRO	Healthcare Printed	40-45 gal.	40 × 46	3.00 mil	Red	FlatPack	75
Z8048VROP01	Healthcare Printed	40-45 gal.	40 × 48	16.00 mil	Red	FlatPack	250
A7450PRO	Healthcare Printed	44 gal.	37 × 50	1.30 mil	Red	FlatPack	150
A7450ARO	Healthcare Printed	44 gal.	37 × 50	1.60 mil	Red	FlatPack	150
A7450ZRO	Healthcare Printed	44 gal.	37 × 50	3.00 mil	Red	FlatPack	75
C7658PROAA9	Healthcare Printed	60 gal.	38 × 58	1.30 mil	Red	FlatPack	100



Can Liners



- Made from highest quality raw materials
- Material reinforced to enhance physical properties
- Strong “star” seal minimizes leaking and maximizes carrying capacity
- Top-side dual dispensing boxes for easy access
- High-performance standards for dry/wet load capacity, tear and impact specifications
- LLDPE—Linear Low Density Polyethylene. Excellent resistance to puncture and tearing.
- HDPE—High Density Polyethylene. High resistance to puncturing.

Can Liners

Item No.	Type	Capacity	Size	Gauge	Color	Packaging	EA/Case
H08G0	Right Sack	8 gal.	24 × 30	1.00 mil	Gray	FlatPack	500
D4823RBO	LLDPE Liners	8-10 gal.	24 × 23	0.35 mil	Black	FlatPack	1,000
D4823RCO	LLDPE Liners	8-10 gal.	24 × 23	0.35 mil	Clear	FlatPack	1,000
D4823HBO	LLDPE Liners	8-10 gal.	24 × 23	0.55 mil	Black	FlatPack	1,000
D4823HCO	LLDPE Liners	8-10 gal.	24 × 23	0.55 mil	Clear	FlatPack	1,000
H4823MCOR01	LLDPE Liners	8-10 gal.	24 × 23	0.55 mil	Clear	CorelessRoll	10/50
A4824PROR01	Healthcare Printed	8-10 gal.	24 × 24	1.30 mil	Red	CorelessRoll	20/25
Z4824RNOR01	HDPE Liners	8-10 gal.	24 × 24	6.00 mil	Natural	CorelessRoll	20/50
Z4824LNOR01	HDPE Liners	8-10 gal.	24 × 24	8.00 mil	Natural	CorelessRoll	20/50
D4832RBO	LLDPE Liners	12-16 gal.	24 × 32	0.35 mil	Black	FlatPack	1,000
D4832RCO	LLDPE Liners	12-16 gal.	24 × 32	0.35 mil	Clear	FlatPack	1,000
H4832MBO	LLDPE Liners	12-16 gal.	24 × 32	0.50 mil	Black	FlatPack	500
H4832MCO	LLDPE Liners	12-16 gal.	24 × 32	0.50 mil	Clear	FlatPack	500
H4832MCOR01	LLDPE Liners	12-16 gal.	24 × 32	0.50 mil	Clear	CorelessRoll	20/25
H4832MWO	LLDPE Liners	12-16 gal.	24 × 32	0.50 mil	White	FlatPack	500
H4832MWOR01	LLDPE Liners	12-16 gal.	24 × 32	0.50 mil	White	CorelessRoll	20/25
H4832HBO	LLDPE Liners	12-16 gal.	24 × 32	0.70 mil	Black	FlatPack	500
H4832HCO	LLDPE Liners	12-16 gal.	24 × 32	0.70 mil	Clear	FlatPack	500
H4832HCOR01	LLDPE Liners	12-16 gal.	24 × 32	0.70 mil	Clear	CorelessRoll	20/25
C4832PROAA9	Healthcare Printed	12-16 gal.	24 × 32	1.30 mil	Red	FlatPack	500
Z4833RNOR01	HDPE Liners	12-16 gal.	24 × 33	6.00 mil	Natural	CorelessRoll	20/50
Z4833LNOR01	HDPE Liners	12-16 gal.	24 × 33	8.00 mil	Natural	CorelessRoll	20/50
Z4833HNOR01	HDPE Liners	12-16 gal.	24 × 33	13.00 mil	Natural	CorelessRoll	20/25

continued on next page



Please see our full line of Waste Receptacles on page 51 in the Equipment Section.

continued from previous page

Can Liners

Item No.	Type	Capacity	Size	Gauge	Color	Packaging	EA/Case
H5242YC0	LLDPE Liners	20-30 gal.	26 × 42	0.90 mil	Clear	FlatPack	200
D5242TWOR01	LLDPE Liners	20-30 gal.	26 × 42	0.90 mil	White	CorelessRoll	8/25
H5838EXOR01	LLDPE Liners	20-30 gal.	29 × 38	0.80 mil	Blue	CorelessRoll	8/25
H6036LBO	LLDPE Liners	20-30 gal.	30 × 36	0.40 mil	Black	FlatPack	250
H6036MBO	LLDPE Liners	20-30 gal.	30 × 36	0.50 mil	Black	FlatPack	250
H6036MBOR01	LLDPE Liners	20-30 gal.	30 × 36	0.50 mil	Black	CorelessRoll	10/25
H6036MCO	LLDPE Liners	20-30 gal.	30 × 36	0.50 mil	Clear	FlatPack	250
H6036HBO	LLDPE Liners	20-30 gal.	30 × 36	0.65 mil	Black	FlatPack	250
H6036HCO	LLDPE Liners	20-30 gal.	30 × 36	0.65 mil	Clear	FlatPack	250
H6036EWO	LLDPE Liners	20-30 gal.	30 × 36	0.75 mil	White	FlatPack	200
H6036TBO	LLDPE Liners	20-30 gal.	30 × 36	0.90 mil	Black	FlatPack	200
H6036TCO	LLDPE Liners	20-30 gal.	30 × 36	0.90 mil	Clear	FlatPack	200
H6036TCOR01	LLDPE Liners	20-30 gal.	30 × 36	0.90 mil	Clear	CorelessRoll	8/25
H6036TWO	LLDPE Liners	20-30 gal.	30 × 36	0.90 mil	White	FlatPack	200
H6036SGO	LLDPE Liners	20-30 gal.	30 × 36	1.10 mil	Gray	FlatPack	250
Z6037LNO	HDPE Liners	20-30 gal.	30 × 37	8.00 mil	Natural	FlatPack	500
Z6037LNOR01	HDPE Liners	20-30 gal.	30 × 37	8.00 mil	Natural	CorelessRoll	20/25
Z6037MKOR02	HDPE Liners	20-30 gal.	30 × 37	10.00 mil	Black	CorelessRoll	20/25
Z6037MNO	HDPE Liners	20-30 gal.	30 × 37	10.00 mil	Natural	FlatPack	500
Z6037MNOR01	HDPE Liners	20-30 gal.	30 × 37	10.00 mil	Natural	CorelessRoll	20/25
Z6037HNO	HDPE Liners	20-30 gal.	30 × 37	12.00 mil	Natural	FlatPack	500
Z6037HNOR01	HDPE Liners	20-30 gal.	30 × 37	13.00 mil	Natural	CorelessRoll	20/25
Z6037XNO	HDPE Liners	20-30 gal.	30 × 37	16.00 mil	Natural	FlatPack	500
Z6037XNOR01	HDPE Liners	20-30 gal.	30 × 37	16.00 mil	Natural	CorelessRoll	20/25
C6043PROAJ5	Healthcare Printed	20-30 gal.	30 × 43	1.30 mil	Red	FlatPack	200
Z6043YNO	HDPE Liners	20-30 gal.	30 × 43	19.00 mil	Natural	FlatPack	250
H6045HCOR01	LLDPE Liners	20-30 gal.	30 × 45	0.65 mil	Clear	CorelessRoll	10/25
H6045TWO	LLDPE Liners	20-30 gal.	30 × 45	0.90 mil	White	FlatPack	200
H6045TWOR01	LLDPE Liners	20-30 gal.	30 × 45	0.90 mil	White	CorelessRoll	8/25
H5645TCOR01	AccuFit	23 gal.	28 × 45	0.90 mil	Clear	CorelessRoll	8/25
H32CO	Right Sack	32 gal.	33 × 45	1.40 mil	Clear	FlatPack	200
H32WO	Right Sack	32 gal.	33 × 45	1.40 mil	White	FlatPack	200
H6639LCO	LLDPE Liners	33 gal.	33 × 39	0.40 mil	Clear	FlatPack	250
H6639MBO	LLDPE Liners	33 gal.	33 × 39	0.50 mil	Black	FlatPack	250
H6639MCO	LLDPE Liners	33 gal.	33 × 39	0.50 mil	Clear	FlatPack	250
H6639HBO	LLDPE Liners	33 gal.	33 × 39	0.65 mil	Black	FlatPack	250
H6639HCO	LLDPE Liners	33 gal.	33 × 39	0.65 mil	Clear	FlatPack	250
H6639HCOR01	LLDPE Liners	33 gal.	33 × 39	0.65 mil	Clear	CorelessRoll	10/25
H6639EWO	LLDPE Liners	33 gal.	33 × 39	0.75 mil	White	FlatPack	150
H6639TBO	LLDPE Liners	33 gal.	33 × 39	0.90 mil	Black	FlatPack	150
H6639TCO	LLDPE Liners	33 gal.	33 × 39	0.90 mil	Clear	FlatPack	150
H6639TWO	LLDPE Liners	33 gal.	33 × 39	0.90 mil	White	FlatPack	150
H6639SCO	LLDPE Liners	33 gal.	33 × 39	1.10 mil	Clear	FlatPack	250
H6639SGO	LLDPE Liners	33 gal.	33 × 39	1.10 mil	Gray	FlatPack	250
H6639ACOR01	LLDPE Liners	33 gal.	33 × 39	1.50 mil	Clear	CorelessRoll	5/20

continued from previous page

Can Liners

Item No.	Type	Capacity	Size	Gauge	Color	Packaging	EA/Case
Z6640LNO	HDPE Liners	33 gal.	33 × 40	9.00 mil	Natural	FlatPack	500
Z6640MNO	HDPE Liners	33 gal.	33 × 40	12.00 mil	Natural	FlatPack	500
Z6640MNOR01	HDPE Liners	33 gal.	33 × 40	12.00 mil	Natural	CorelessRoll	20/25
Z6640HNO	HDPE Liners	33 gal.	33 × 40	14.00 mil	Natural	FlatPack	250
Z6640HNOR01	HDPE Liners	33 gal.	33 × 40	14.00 mil	Natural	CorelessRoll	10/25
Z6640VNO	HDPE Liners	33 gal.	33 × 40	16.00 mil	Natural	FlatPack	250
Z6640VNOR01	HDPE Liners	33 gal.	33 × 40	16.00 mil	Natural	CorelessRoll	10/25
Z6640XNO	HDPE Liners	33 gal.	33 × 40	17.00 mil	Natural	FlatPack	250
Z6640XNOR01	HDPE Liners	33 gal.	33 × 40	17.00 mil	Natural	CorelessRoll	10/25
Z6640WKO	HDPE Liners	33 gal.	33 × 40	22.00 mil	Black	FlatPack	100
H8046MCO	LLDPE Liners	40-45 gal.	40 × 46	0.55 mil	Clear	FlatPack	250
D8046HBO	LLDPE Liners	40-45 gal.	40 × 46	0.70 mil	Black	FlatPack	250
H8046HBO	LLDPE Liners	40-45 gal.	40 × 46	0.70 mil	Black	FlatPack	125
D8046HCO	LLDPE Liners	40-45 gal.	40 × 46	0.70 mil	Clear	FlatPack	250
H8046HCO	LLDPE Liners	40-45 gal.	40 × 46	0.70 mil	Clear	FlatPack	125
H8046HCOR01	LLDPE Liners	40-45 gal.	40 × 46	0.70 mil	Clear	CorelessRoll	5/25
H8046EWO	LLDPE Liners	40-45 gal.	40 × 46	0.75 mil	White	FlatPack	100
H8046EWOR01	LLDPE Liners	40-45 gal.	40 × 46	0.75 mil	White	CorelessRoll	4/25
H8046TBO	LLDPE Liners	40-45 gal.	40 × 46	0.90 mil	Black	FlatPack	100
H8046TBOR01	LLDPE Liners	40-45 gal.	40 × 46	0.90 mil	Black	CorelessRoll	4/25
H8046TCO	LLDPE Liners	40-45 gal.	40 × 46	0.90 mil	Clear	FlatPack	100
H8046TWO	LLDPE Liners	40-45 gal.	40 × 46	0.90 mil	White	FlatPack	100
H8046SCO	LLDPE Liners	40-45 gal.	40 × 46	1.10 mil	Clear	FlatPack	125
H8046SGO	LLDPE Liners	40-45 gal.	40 × 46	1.10 mil	Gray	FlatPack	125
X8046SKO	Repro Liners	40-45 gal.	40 × 46	1.20 mil	Black	FlatPack	100
X8046PKO	Repro Liners	40-45 gal.	40 × 46	1.30 mil	Black	FlatPack	100
X8046AKO	Repro Liners	40-45 gal.	40 × 46	1.50 mil	Black	FlatPack	100
H8046ACO	LLDPE Liners	40-45 gal.	40 × 46	1.50 mil	Clear	FlatPack	100
H8046QKO	LLDPE Liners	40-45 gal.	40 × 46	2.00 mil	Black	FlatPack	100
X8046QKO	Repro Liners	40-45 gal.	40 × 46	2.00 mil	Black	FlatPack	100
H8046XCO	LLDPE Liners	40-45 gal.	40 × 46	2.00 mil	Clear	FlatPack	100
C8046PROAA9	Healthcare Printed	40-45 gal.	40 × 46	1.30 mil	Red	FlatPack	200
C8046PROAE7	Healthcare Printed	40-45 gal.	40 × 46	1.30 mil	Red	FlatPack	200
C8046PROAG6	Healthcare Printed	40-45 gal.	40 × 46	1.30 mil	Red	FlatPack	100
Z8048MNO	HDPE Liners	40-45 gal.	40 × 48	11.00 mil	Natural	FlatPack	250
Z8048MNOR01	HDPE Liners	40-45 gal.	40 × 48	11.00 mil	Natural	CorelessRoll	10/25
Z8048MNOR03	HDPE Liners	40-45 gal.	40 × 48	12.00 mil	Natural	CorelessRoll	10/25
Z8048HNO	HDPE Liners	40-45 gal.	40 × 48	13.00 mil	Natural	FlatPack	250
Z8048HNOR01	HDPE Liners	40-45 gal.	40 × 48	13.00 mil	Natural	CorelessRoll	10/25
Z8048VNO	HDPE Liners	40-45 gal.	40 × 48	16.00 mil	Natural	FlatPack	250
Z8048VNOR01	HDPE Liners	40-45 gal.	40 × 48	16.00 mil	Natural	CorelessRoll	10/25
Z8048XNO	HDPE Liners	40-45 gal.	40 × 48	17.00 mil	Natural	FlatPack	250
Z8048XNOR01	HDPE Liners	40-45 gal.	40 × 48	17.00 mil	Natural	CorelessRoll	10/25
Z8048WKOR01	HDPE Liners	40-45 gal.	40 × 48	22.00 mil	Black	CorelessRoll	6/25
Z8048UNOR01	HDPE Liners	40-45 gal.	40 × 48	22.00 mil	Natural	CorelessRoll	10/15
Z7446MNOR01	HDPE Liners	44 gal.	37 × 46	12.00 mil	Natural	CorelessRoll	10/25
Z7446HNOR01	HDPE Liners	44 gal.	37 × 46	14.00 mil	Natural	CorelessRoll	10/25
H44CO	Right Sack	44 gal.	37 × 50	1.50 mil	Clear	FlatPack	150
H44GO	Right Sack	44 gal.	37 × 50	1.50 mil	Gray	FlatPack	150

continued from previous page

Can Liners

Item No.	Type	Capacity	Size	Gauge	Color	Packaging	EA/Case
S7646VNO	Tie-Tie	44 gal.	38 × 46	16.00 mil	Natural	CorelessRoll	12/20
H7258MCO	LLDPE Liners	55 gal.	36 × 58	0.55 mil	Clear	FlatPack	200
Z7260XNO	HDPE Liners	55 gal.	36 × 60	17.00 mil	Natural	FlatPack	200
X8646AKO	Repro Liners	56 gal.	43 × 46	1.50 mil	Black	FlatPack	100
D8647HBO	LLDPE Liners	56 gal.	43 × 47	0.70 mil	Black	FlatPack	200
D8647HCO	LLDPE Liners	56 gal.	43 × 47	0.70 mil	Clear	FlatPack	200
H8647HCO	LLDPE Liners	56 gal.	43 × 47	0.70 mil	Clear	FlatPack	100
H8647EWO	LLDPE Liners	56 gal.	43 × 47	0.75 mil	White	FlatPack	100
H8647TCOR01	LLDPE Liners	56 gal.	43 × 47	0.95 mil	Clear	CorelessRoll	4/25
H8647SCO	LLDPE Liners	56 gal.	43 × 47	1.10 mil	Clear	FlatPack	100
H8647SGO	LLDPE Liners	56 gal.	43 × 47	1.10 mil	Gray	FlatPack	100
H8647PCO	LLDPE Liners	56 gal.	43 × 47	1.30 mil	Clear	FlatPack	100
H8648WSO	Silver Sack	56 gal.	43 × 48	1.70 mil	Silver	FlatPack	50
Z8648HNO	HDPE Liners	56 gal.	43 × 48	14.00 mil	Natural	FlatPack	200
Z8648HNOR01	HDPE Liners	56 gal.	43 × 48	14.00 mil	Natural	CorelessRoll	8/25
Z8648VNO	HDPE Liners	56 gal.	43 × 48	16.00 mil	Natural	FlatPack	200
Z8648VNOR01	HDPE Liners	56 gal.	43 × 48	16.00 mil	Natural	CorelessRoll	8/25
Z8648XKO	HDPE Liners	56 gal.	43 × 48	17.00 mil	Black	FlatPack	200
Z8648XNO	HDPE Liners	56 gal.	43 × 48	17.00 mil	Natural	FlatPack	200
Z8648WKO	HDPE Liners	56 gal.	43 × 48	22.00 mil	Black	FlatPack	150
Z8648WKOR01	HDPE Liners	56 gal.	43 × 48	22.00 mil	Black	CorelessRoll	6/25
Z8648WNOR01	HDPE Liners	56 gal.	43 × 48	22.00 mil	Natural	CorelessRoll	6/25
H56CO	Right Sack	56 gal.	44 × 55	1.60 mil	Clear	FlatPack	100
H56GO	Right Sack	56 gal.	44 × 55	1.60 mil	Gray	FlatPack	100
X9250SKO	Repro Liners	56 gal.	46 × 50	1.35 mil	Black	FlatPack	100
H7658MCO	LLDPE Liners	60 gal.	38 × 58	0.55 mil	Clear	FlatPack	200
H7658HBO	LLDPE Liners	60 gal.	38 × 58	0.70 mil	Black	FlatPack	200
H7658HCO	LLDPE Liners	60 gal.	38 × 58	0.70 mil	Clear	FlatPack	200
H7658EWO	LLDPE Liners	60 gal.	38 × 58	0.75 mil	White	FlatPack	100
H7658TBO	LLDPE Liners	60 gal.	38 × 58	0.90 mil	Black	FlatPack	100
H7658TCO	LLDPE Liners	60 gal.	38 × 58	0.90 mil	Clear	FlatPack	100
H7658TWO	LLDPE Liners	60 gal.	38 × 58	0.90 mil	White	FlatPack	100
H7658SCO	LLDPE Liners	60 gal.	38 × 58	1.10 mil	Clear	FlatPack	100
H7658SGO	LLDPE Liners	60 gal.	38 × 58	1.10 mil	Gray	FlatPack	100
H7658SKO	LLDPE Liners	60 gal.	38 × 58	1.25 mil	Black	FlatPack	100
H7658PCO	LLDPE Liners	60 gal.	38 × 58	1.30 mil	Clear	FlatPack	100
X7658AKO	Repro Liners	60 gal.	38 × 58	1.50 mil	Black	FlatPack	100
H7658ACO	LLDPE Liners	60 gal.	38 × 58	1.50 mil	Clear	FlatPack	100
X7658QKO	Repro Liners	60 gal.	38 × 58	2.00 mil	Black	FlatPack	100
Z7660HNO	HDPE Liners	60 gal.	38 × 60	14.00 mil	Natural	FlatPack	200
Z7660VNOR01	HDPE Liners	60 gal.	38 × 60	16.00 mil	Natural	CorelessRoll	8/25
Z7660XEOR01	HDPE Liners	60 gal.	38 × 60	17.00 mil	Green	CorelessRoll	10/20
Z7660XNO	HDPE Liners	60 gal.	38 × 60	17.00 mil	Natural	FlatPack	200
Z7660XNOR01	HDPE Liners	60 gal.	38 × 60	17.00 mil	Natural	CorelessRoll	10/20
Z7660WNOR01	HDPE Liners	60 gal.	38 × 60	22.00 mil	Natural	CorelessRoll	6/25
H0248TKO	LLDPE Liners	65 gal.	51 × 48	1.00 mil	Black	FlatPack	100



Infectious Linen Bags



- Meets federal and state regulations for tear and impact resistance
 - 480 gram Elmendorf tear test (ASTM D1922-94a)
 - 165 gram Dart Impact Resistance Test (ASTM D1709-97)
- Color Coded for Safety and Identification
- Reinforced Seams for Superior Leak Resistance
- Bilingual Printing (English and Spanish)
- Wide range of sizes and thicknesses

Infectious Linen Bags

Item No.	Type	Capacity	Size	Gauge	Color	Packaging	EA/Case
A6043PYO	Healthcare Printed	20-30 gal.	30 × 43	1.30 mil	Yellow	FlatPack	200
B8046SYO	Healthcare N	40-45 gal.	40 × 46	1.20 mil	Yellow	FlatPack	100
A7450PYO	Healthcare Printed	44 gal.	37 × 50	1.30 mil	Yellow	FlatPack	150



Soiled Linen Bags



- Wide Range of Sizes and Thicknesses
- Color Coded Blue or Mint Green for Safety and Identification
- Superior Leak Resistance: Tested for load capacity and seal integrity
- Reinforced Seams
- Bilingual Printing (English & Spanish) where noted

Soiled Linen Bags

Item No.	Type	Capacity	Size	Gauge	Color	Packaging	EA/Case
A6043TXOR01	Healthcare Printed	20-30 gal.	30 × 43	1.00 mil	Blue	CorelessRoll	8/25
A6043YMORH5	Healthcare Printed	20-30 gal.	30 × 43	1.00 mil	miInt	CorelessRoll	8/25
B6043PXO	Healthcare N	20-30 gal.	30 × 43	1.30 mil	Blue	FlatPack	200
A6043PXO	Healthcare Printed	20-30 gal.	30 × 43	1.30 mil	Blue	FlatPack	200
Z6143HXOP01	Healthcare Printed	20-30 gal.	31 × 43	14.00 mil	Blue	FlatPack	250
Z6143HXO	HDPE Liners	20-30 gal.	31 × 43	14.00 mil	Blue	FlatPack	250
Z6640HXOR01	HDPE Liners	33 gal.	33 × 40	14.00 mil	Blue	CorelessRoll	10/25
Z8048YXOH05	Healthcare Printed	40-45 gal.	40 × 48	18.00 mil	Blue	FlatPack	200
A7450PXO	Healthcare Printed	44 gal.	37 × 50	1.30 mil	Blue	FlatPack	150
Z7646VXOPR1	Healthcare Printed	44 gal.	38 × 46	16.00 mil	Blue	CorelessRoll	10/25



Please see our full line of Chrome Hampers on pages 44-45 in the Equipment Section.

Foaming Skin Cleansers



- Multi-use formula—No-Rinse Shampoo, Body Wash & Perineal Wash
- Easy and convenient finger-pump foaming dispenser—no wasteful spillage
- Formulated with aloe vera and vitamins A and E to promote skin health
- pH balanced—enriched with vitamins and antioxidants
- Pleasant fragrance
- OVI21033 contains Cupuacu Butter



Item No.	Description	Quantity
Cleanser Skin Foaming Standard without Rinse		
OVI21033	6 Ounce	1/EA, 12 EA/CA
FC3085	8.5 Ounce	1/EA, 12 EA/CA

Seamless Knit Pants



- Adaptable to hold incontinence liners or pads in place
- Eliminates seam lines under clothing for patient dignity
- Sized from small to bariatric
- Available in two pack, five pack and bulk cases of 100 each
- Constructed of 94% Polyester and 6% Spandex



Item No.	Description	Quantity
Brief Incontinence Seamless Knit Pant 94 Percent Polyester and 6 Percent Spandex		
MBS1001	Yellow And Green Waistband Fits Size Large to Extra Large Bulk Packaged 100 Each Per Case	1/EA, 100 EA/CA
MBS1002	Yellow And Green Waistband Fits Size Large to Extra Large Packaged 2 Each Per Pack with 50 Packs Per Case	2/PK, 50 PK/CA
MBS1005	Yellow And Green Waistband Fits Size Large to Extra Large Packaged 5 Each Per Pack with 20 Packs Per Case	5/PK, 20 PK/CA
MBS1101	Purple And Gray Waistband Fits Size 2XL to 3XL Packaged 100 Each Per Case	1/EA, 100 EA/CA
MBS1102	Purple And Gray Waistband Fits Size 2XL to 3XL Packaged 2 Each Per Pack with 50 Packs Per Case	2/PK, 50 PK/CA
MBS1105	Purple And Gray Waistband Fits Size 2XL to 3XL Packaged 5 Each Per Pack with 20 Packs Per Case	5/PK, 20 PK/CA
MBS1201	White Waistband Fits Size 4XL to Bariatric Bulk Packaged 100 Each Per Case	1/EA, 100 EA/CA
MBS1202	White Waistband Fits Size 4XL to Bariatric Packaged 2 Each Per Pack with 50 Packs Per Case	2/PK, 50 PK/CA
MBS1205	White Waistband Fits Size 4XL to Bariatric Packaged 5 Each Per Pack with 20 Packs Per Case	5/PK, 20 PK/CA
MBS9001	Brown And Blue Waistband Fits Size Small To Medium Bulk Packaged 100 Each Per Case	1/EA, 100 EA/CA
MBS9002	Brown And Blue Waistband Fits Size Small To Medium Packaged 2 Each Per Pack with 50 Packs Per Case	2/PK, 50 PK/CA
MBS9005	Brown And Blue Waistband Fits Size Small To Medium Packaged 5 Each Per Pack with 20 Packs Per Case	5/PK, 20 PK/CA

Underpads



- Super absorbent soft core that quickly and effectively wicks away fluid
- Moisture vapor permeable film allows moisture vapor to transfer away from the patient while protecting bed linens from bodily fluids, leaving the patient dry and comfortable
- Premium Cloth-Like Underpads
 - Strong cloth-like backsheet
 - Resists tearing and ripping while maintaining its integrity
 - 300 pound weight claim
- Premium Strength Underpads
 - Premium strength backsheet assists in the lift and repositioning of patients up to 375 pounds, which helps to reduce the need for drawsheets and reusable underpads



Item No.	Description	Quantity
Underpad Premium Clothlike with Super Absorbent Polymer		
PCL2336	300 Pound Weight Claim 23 Inch × 36 Inch Disposable	1/EA, 12 EA/BG, 6 BG/CA
PCL3036	300 Pound Weight Claim 30 Inch × 36 Inch Disposable	1/EA, 10 EA/BG, 4 BG/CA
Underpad Premium Strength with Super Absorbent Polymer		
PSU2336	375 Pound Weight Claim 23 Inch × 36 Inch Disposable	1/EA, 12 EA/BG, 6 BG/CA
PSU3036	375 Pound Weight Claim 30 Inch × 36 Inch Disposable	1/EA, 10 EA/BG, 4 BG/CA

Odor Eliminators



- Dye Free, Alcohol Free, and Prop 65 Compliant
- Clear formula—can be used on fabrics, skin and hair
- Odor elimination for general patient, necrotic tissue and incontinence odors
- Eliminates odors
- Two formulas: Fresh Air Scented and Unscented
- Available in two convenient sizes: 2 ounce and 8 ounce pump spray bottles



Item No.	Description	Quantity
Eliminator Odor Biological		
OE2240	Scented Pump Spray 2 Ounce	1/EA, 24 EA/CA
OE2241	Unscented Pump Spray 2 Ounce	1/EA, 24 EA/CA
OE8120	Scented Pump Spray 8 Ounce	1/EA, 12 EA/CA
OE8121	Unscented Pump Spray 8 Ounce	1/EA, 12 EA/CA

Please see our full line of Long Torso Seamless Knit Pants on page 130 in the Maternity and Baby care Section.



Washcloths

NOT MADE WITH
natural rubber latex

diverse
supplier

- Soft material is an economical alternative to cloth
- Quick absorption, strong even when wet
- Air-laid washcloths manufactured with cellulose
- Thermabond washcloth manufactured with rayon, polypropylene and viscous material
- First quality material (not factory seconds)
- Available in 10 Inch × 13 Inch and 12 Inch × 13 Inch dimensions



Item No.	Description	Quantity
Washcloth Thermabond Premium Rayon Polypropylene		
101305	10 Inch × 13 Inch	1/EA, 50 EA/BG, 10 BG/CA
Washcloth Airlaid Medium Weight Disposable		
71013	10 Inch × 13 Inch	1/EA, 50 EA/BG, 10 BG/CA
71213	12 Inch × 13 Inch	1/EA, 50 EA/BG, 10 BG/CA

Protective Underwear

made in the
USA

NOT MADE WITH
natural rubber latex

- Underwear-like for better quality of life
- Tear away sides for easy removal
- Cloth-like moisture barrier protects against leaks to keep clothing dry
- Contains super-absorbent polymer for maximum fluid retention while reducing odor and keeping wetness away from skin



Item No.	Description	Quantity
Underwear Protective Pull-On		
PUW1605	Size Small to Medium 34 to 46 Inch	20/BG, 4 BG/CA
PUW1615	Size Large 44 to 54 Inch	18/BG, 4 BG/CA
PUW1625	Size Extra Large 48 to 66 Inch	14/BG, 4 BG/CA

Flushable Cleansing Cloths

made in the
USA

NOT MADE WITH
natural rubber latex

- Flushable, moist wipes for personal cleansing
- 7 Inch × 5 Inch
- Soft and gentle tissues
- Fibers made of 100 percent short cut staple fibers from tree pulp which allows the wipe to disperse like toilet paper
- Safe for sewers and septic systems
- Soft pack with plastic recloseable lid
- Produced with 100% clean wind energy
- Made with artesian glacial aquifer water




Item No.	Description	Quantity
FCC1001	Cloth Flushable Personal Moist Cleansing 7 Inch By 5 Inch Soft Pack Recloseable Lid 48 Wipes in Each Pack Septic Tank Safe 7 Inch × 5 Inch	48/PK, 12 PK/CA

Adult Briefs

made in the
USA

NOT MADE WITH
natural rubber latex

- Heavy Briefs
 - Fastening system features 4 sturdy tape tabs to ensure a secure fit
 - Combines super-absorbent polymer, blue dryness strip and breathable topsheet for added comfort, protection and odor control
 - Diamond and linear embossing provides improved moisture wicking for added dryness
- Ultra Briefs
 - Full fit design provides complete coverage with greater comfort
 - Extra wide crotch for added protection
 - Fastening system features 4 sturdy tape tabs to ensure a secure fit
 - Combines super-absorbent polymer, blue dryness strip, breathable topsheet and extra fluff in the side panel areas for added comfort, protection and odor control
- Cloth-Like Briefs
 - Cloth-like backsheet is ultrasonically bonded for a "quilted" look and feel that is soft and quiet for enhanced comfort and dignity
 - Breathable side panels allow air to pass through, keeping skin dry, cool and comfortable
 - High performance acquisition layer quickly pulls fluid into polymer and fluff core and locks fluid into core to keep skin dry
 - Combines super-absorbent polymer, a high performance blue dryness strip, breathable topsheet and refastenability for added comfort, protection and odor control



Item No.	Description	Quantity
Brief Adult Heavy		
90410	Size Medium Fits Waist or Hip Size 32 to 44 Inches	12/BG, 8 BG/CA
90411	Size Large Fits Waist or Hip Size 45 to 58 Inches	12/BG, 6 BG/CA
90412	Size Extra Large Fits Waist or Hip Size 59 to 64 Inches	10/BG, 6 BG/CA
Brief Adult Ultra with Wetness Indicator		
90413	Size Small Fits Waist or Hip Size 23 to 31 Inches	12/BG, 8 BG/CA
90414	Size Medium Fits Waist or Hip Size 32 to 44 Inches	12/BG, 8 BG/CA
90415	Size Large Fits Waist or Hip Size 45 to 58 Inches	12/BG, 6 BG/CA
90417	Size Extra Large Fits Waist or Hip Size 59 to 64 Inches	10/BG, 6 BG/CA
Brief Adult Cloth-like with Wetness Indicator		
90601	Size 32 to 44 Inches Super Absorbent	10/BG, 10 BG/CA
90602	Size Large Fits Waist or Hip Size 45 to 58 Inches	18/BG, 4 BG/CA
90603	Size Extra Large Fits Waist or Hip Size 59 to 64 Inches	15/BG, 4 BG/CA

Adult Briefs Sizing Chart

Waist or Hip Sizes (in Inches)	
Small	23-31
Medium	32-44
Large	45-58
X-Large	59-64

ADULT WIPES

Flushable Adult Wipes

NOT MADE WITH
natural rubber latex

- Large 9 Inch x 13 Inch pre-moistened disposable, flushable, dispersible wipe
- Single-Use
- Made from Ahlstrom, Hydraspun material
- Hypoallergenic and Alcohol-Free
- Wipes available in Scented, Unscented, or 3.2% Dimethicone skin protectant
- Available in Soft-Pack and Tubs



Item No.	Description	Quantity
Wipe Flushable Unscented Hydraspun Pre-Moistened Disposable		
UFW413	Softpack Large 9 Inch x 13 Inch	24/PK, 24 PK/CA
UFW613	Tub Large 9 Inch x 13 Inch	60/TB, 9 TB/CA
FWD513	Softpack with 3.2 Percent Dimethicone Large 9 x 13 Inch	24/PK, 24 PK/CA

Unscented Adult Wipes

NOT MADE WITH
natural rubber latex

- Large 9 Inch x 13 Inch pre-moistened disposable wipe for incontinence or general cleansing
- Single-Use
- 45-gram Spunlace material
- Hypoallergenic and Alcohol-Free
- Contains Aloe & Lanolin to promote skin health
- Soft-Pack
- Fragrance free



Item No.	Description	Quantity
Wipe Incontinence Adult Unscented Premium Pre-Moistened		
0713MC	50 Each Softpack with Resealable Lid 9 x 13 Inch	50/PK, 12 PK/CA
AW913UBP	24 Each Softpack with Peel Lid 9 x 13 Inch	24/PK, 24 PK/CA

Scented Adult Wipes

NOT MADE WITH
natural rubber latex

- Large 9 Inch x 13 Inch pre-moistened disposable wipe for incontinence or general cleansing
- Single-Use
- 45-gram Spunlace material
- Contains Aloe & Lanolin to promote skin health
- Alcohol-free, hypoallergenic
- Available in single-dispensing soft packs and tubs
- Fresh scent



Item No.	Description	Quantity
Wipe Cleaning Scented Pre-Moistened Premium Adult Spunlace 9 x 13 Inch Disposable		
0913MC	Hard Tub	60/TB, 12 TB/CA
0813MC	Soft Pack Refill	50/PK, 12 PK/CA
AW913BP	Soft Pack with Peel Lid	24/PK, 24 PK/CA

Cleansing Cremes

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Combines superior moisture retention with everyday skin protection
- Formulated with Cupuacu Butter which contains vitamins, antioxidants, and phytosterols
- Trigger-sprayer helps with application process
- Active Ingredient: Dimethicone 1.5%



Item No.	Description	Quantity
OVI21042	Creme Skin Cleansing Boavida 1.5 Percent Dimethicone 8 Ounce	1/EA, 12 EA/CA

Skin Protectant

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Helps prevent minor skin irritation associated with wet skin, urine & stool
- Unique blend of silicones form an effective & persistent shield
- Fast transepidermal penetration—rapidly relieves dry, chapped skin
- Lowers Transepidermal Water Loss (TEWL)—helps skin retain moisture
- Active Ingredient—Dimethicone 2.0%



Item No.	Description	Quantity
OVI21024	Cream Skin Barrier Boavida Canopy 4 Ounce 2 Percent Dimethicone	1/EA, 12 EA/CA

Skin Recovery

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Helps with daily recovery of dry, malnourished skin
- Major moisture absorption properties
- Helps restore the lipid barrier of skin
- Reduces Transepidermal Water Loss (TEWL) from outer skin layer
- CHG & glove compatible
- Active Ingredient: Dimethicone 1.5%



Item No.	Description	Quantity
Cream Skin Recovery Boavida 1.5 Percent Dimethicone		
OVI21002	2 Ounce	1/EA, 24 EA/CA
OVI21004	4 Ounce	1/EA, 12 EA/CA
OVI21009	32 Ounce	1/EA, 6 EA/CA

Skin Treatment

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**

- Treats minor skin irritation associated with wet skin, urine and stool
- Ideal for denuded and/or macerated skin
- Zinc Oxide helps cover and treat affected areas
- Contains menthol to help quickly relieve pain and discomfort
- Formulated with calamine to further soothe irritated skin
- Active Ingredient—Zinc Oxide 20%



Item No.	Description	Quantity
OVI21014	Cream Zinc Oxide Skin Treatment 4 Ounce Boavida Calvida	1/EA, 12 EA/CA

KITS & TRAYS



MediChoice's minor procedure kits and trays have been carefully designed to meet the demands of a wide range of protocols and procedures. Choose your kits based on procedures specific to your healthcare needs, then distribute these standardized kits wherever they are needed. Used throughout an entire facility or system, they can increase efficiency, reduce cost and waste, and most importantly, ensure patients receive a consistent, quality level of care.

- Build your kits to meet your specific requirements
- Choose from a wide range of brand names, including our own MediChoice brand
- A variety of antiseptics, dressings, securement devices and other components are available
- Arrange your components in order of use

Suture Trays

made in the USA NOT MADE WITH natural rubber latex sterile EO

- Stackable tray, wrapped and bagged
- Most commonly used components
- Packaged in typical order of use
- Precision floor grade or economy instruments



Suture Trays	MSUT1002	MSUT1004	MSUT1005
Components			
Packaging			
Hard Tray, 1 Compartment	1		
Hard Tray, 2 Compartment		1	1
Hard Tray w/ Tyvek Lid			
Wrap, CSR, 24" x 24"	1		1
Prep And Dressing			
Gauze, 2" x 2", 8 ply	8		
Gauze, 4" x 4", 8 ply	5		5
Gauze, 4" x 4", 12 ply		15	
Instruments			
Webster, Needle Holder, Satin, 5"	1	1	
Webster, Needle Holder, Smooth, 5"			
Webster, Needle Holder, Fine Tip, 5"			1
Scissors, Iris, Straight, 4½"	1	1	1
Forcep, Adson, Serrated, 4¾"		1	1
Forcep, Adson, Tissue, 4¾", 1 x 2 Teeth	1		
Forcep, Kelly Straight 5½"		1	
Antiseptics			
Lidocaine 2% 10mL		1	
Neosporin		1	
Towels and Drapes			
Drape, Fenestrated, 18" x 26"	1	1	1
Drape, Non-Fenestrated, 18" x 26"			1
Towel, 17" x 19" Poly 2/PP			
Towel, Kaycell, Blue, 19" x 27"	1		
Towel, Kaycell, White, 19" x 27"		2	
Miscellaneous			
Safety Needle 25 GA x ¾"		1	
Safety Needle BD 18GA x 1½"		1	
Safety Needle BD 25GA x 1½"		1	
Syringe 12cc L/L		1	
Medicine Cup, 2 oz.	2		
Ampule Holder 5 mL		1	
Fabric Adhesive Bandage			1
Quantity/Case	20	20	20

Incision and Drainage Trays

NOT MADE WITH natural rubber latex diverse supplier sterile EO

- Packaged in typical order of use
- Most commonly used components
- Stackable tray, wrapped and bagged



Incision and Drainage Trays	INC1001B
Components	
Waste Bag, 12" x 18", White w/Tie	1
Fenestrated Drape, 13" x 18"	1
Towel, 13" x 19", 2-Ply, Tissue/Poly	1
Gauze 3" x 3", 12-Ply	4
PVP Prep Pad	1
Alcohol Prep Pad, Medium	1
Scalpel, #11	1
Forceps, 5", Wireform	1
Scissors, Littauer, 4.5", Wireform	1
Hemostat, Kelly, Serrated, 5.5"	1
Quantity/Case	20

Lumbar Puncture Trays

made in the USA NOT MADE WITH natural rubber latex sterile EO

- Ultra-sharp spinal needles for smooth insertion
- Precisely calibrated manometer with well-defined markings for visual acuity
- Luer slip fitting assures quick and leak-free connections
- Prenumbered vials for accurate specimen identification



Lumbar Puncture Trays	LPT648
Components	
Multi-Compartment Tray with Prep Well	1
Spinal Needle, 20 Gauge x 3.5"	1
Syringe, 3cc with Needle, 25 Gauge x .625"	1
Needle, 20 Gauge x 1.5"	4
Spinal Manometer with Luer Slip, 2-Piece, 550 mm	1
Extension Tube	1
Three-way Stopcock	1
Sponge Sticks	3
Gauze Pads	3
Specimen Vials, with Caps 10 ml	4
Fenestrated Drape with Adhesive Strip	1
Towel	1
Lidocaine Hydrochloride USP, 1% 2 ml with Drug Precaution Sheet	1
CSR Wrap	4
Bandage	1
Patient I.D. Labels	1
Universal Medication Labels	4
Quantity/Case	10

Laceration Trays

made in the USA NOT MADE WITH natural rubber latex sterile EO



- Stackable containers for safe, easy storage
- Convenient easy-to-open packaging
- Packaged in typical order of use

Laceration Trays	LAC1001	LAC1002	LAC1003B	LAC1004	LAC2001	LAC2001A	LAC1010	LAC1011	LAC1012	LAC1013	LAC1015	LAC2004	LAC2006	LAC2007	LAC2018	LAC2020	LAC2021	LAC2022	LAC2023	LAC2024	LAC2025	LAC2026
Components																						
Packaging																						
Hard Tray, 2 Compartment					1	1								1			1				1	
Hard Tray, Rectangular	1	1	1	1				1	1	1	1		1		1		1		1	1	1	1
Hard Tray, 3 Compartment												1				1			1			
Hard Tray w/Tyvek Lid							1															
CSR Wrap	1	1	1	1	1	1		1	1	1	1	1	1	1	1	1		1	1	1	1	1
Prep and Dressing																						
Gauze, 2" x 2", 8-ply			5	5	5	5				5							10					
Gauze, 2" x 2", 12-ply	5							5	5			5	5					5	5	5		5
Gauze, 4" x 4", 8-ply		8	10	10	10	10				10					10			10				10
Gauze, 4" x 4", 12-ply							10	10	10				10	10			10			10		
Gauze, 4" x 4", 16-ply	10																			10		10
Towels/Drapes																						
Drape, Fenestrated, 18" x 26"	1	1		1	1	1			1	1	1	1	1		1		1	1	1	1	1	1
Drape, 17" x 19"		1					1							2								
Drape, 18" x 26", White 3-ply															1							
O.R. Towel, Blue							1	2		2		2	1			1			2			
O.R. Towel, Fenestrated 18" x 26" (3" Fen)			1				1	1								1						
Towel, White 23" x 24"																	1					
Absorbant Towel, Blue						2														2		2
Absorbant Towel, White																	2	3				1
Kaycell Towel, Blue, 19" x 27"	2		2	2	2				2													2
Instruments																						
Scissor, Iris, Straight, 4.5"	1	1	1	1	1	1	1	1	1	1		1	1		1		1	1	1	1	1	1
Scissor, Iris, Curved, 4.5"																1						
Scissor, O.R., 5.5"											1			1								
Forcep, Adson Tissue, 4.75", 1 x 2 Teeth		1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Forcep, Adson Serrated, 4.75"	1						1					1				1	1					1
Forcep, Iris, Straight, 4", 1 x 2 Teeth		1																				
Forcep, Thumb Dressing, Serrated, 5.5"												1							1			
Halstead, Mosquito, Curved, 5"	1	1	1	1	1	1	1	1	1	1		1	1		2	1		1	1	1		1
Baumgartner, Needle Holder, 5.5"														1								
Webster, Needle Holder, Smooth, 5"	1	1			1	1	1	1	1	1	1		1		1	1	1	1		1	1	1
Webster, Needle Holder, Satin, 5"				1																		
Webster, Needle Holder, Large Fingers			1									1							1			
Instrument Tips																						
Needles and Syringes																						
Syringe, Luer Lock, 10cc	1				1	1				1			1		1	1	1	1		1	1	1
Syringe, Luer Lock, 12cc			1	1				1	1													
Syringe, Luer Slip, 10cc			1																			
Syringe, Control, 10cc												1							1			
Syringe, Luer Lock, 20cc																	1				1	
Syringe, Safety-Lok, 1cc w/ 27G x .5" Needle												1							1			
Needle, 18G x 1.5"	1			1	1	1		1	1	1			1		1			1		1	1	1
Needle, 21G x 1.5"															1							
Needle, 22G x 1.5"																	1				1	
Needle, 25G x .625"	1			1	1	1		1	1						1			1		1		1
Needle, 25G x 1.5"											1											
Needle, 27G x .5"					1	1				1							1					1
Needle, 27G x 1.5"	1			1				1	1													
Needle, Safety, 18G x 1.5"			1									1							1			
Needle, Safety, 25G x .625"			1																			
Needle, Safety, 25G x 1"												1	1						1			
Needle, Safety, 27G x .5"			1																			
Needle, 27G x 1.5" Hub A Bev																	1			1		1
Needle, 27G x 1.25" Hub A Bev													1									
Needle, 18G x 1.25" Filter																	1					
Miscellaneous																						
Cotton Tip Applicator, 6", Plastic		2											2		4			4				
Medicine Cups, 2 oz	2	1	2	2	2	2	1	2	2	2			2		2			2		2		2
Ruler						1						1					1		1		1	
Tape Measure			1							1												
Skin Marker																1						
Splashcap																			1			
Supershield												1										
Ampule Breaker																	1					1
CHG Scrub Solution, 4 oz																						1
Label, White .5" x 2"			2																			
Quantity/Case	20	20	20	20	20	20	20	20	20	20	24	20	20	20	20	20	10	20	12	20	10	16

Central Line Dressing Change Trays

NOT MADE WITH natural rubber latex **diverse supplier** sterile EO

- Sterile trays protect antibiotic contents
- Packaged in typical order of use
- Convenient easy-to-open packaging
- Stackable containers for safe, easy storage



Choose MediChoice Transparent Dressings to protect IV sites and surgical incisions. Our dressings conform to hard-to-dress areas, allow constant observation, and can be left in place for several days.



Central Line Dressing Change Trays	CLC1001	CLC1002	CLC1004	CLC1007	CLC1007B	CLC1015	CLC1017	CLC1022	CLC1025	CLC1029	CLC1029A	CLC1035	CLC1058	CLC1090A	CLC1091	CLC1094	CLC1095	CLC4004	CLC3001	CLC1002	
Components																					
Packaging																					
Soft Tray (Blister) w/ Tyvek Lid	1	1	1	1	1	1	1							1	1	1				1	
Foam Tray										1	1										
Hard Tray w/ Tyvek Lid																					1
CSR Wrap 24" x 24"										1	1										1
Prep and Dressing																					
Gauze, 2" x 2", 8-Ply	3	3		1	1		2							2			1				
Gauze, 4" x 4", 8-Ply	2		2						2				2					2	2		
Gauze, 2" x 2", 12-Ply																					1
Gauze, 4" x 4", 12-Ply																					1
Gauze, Split, 2" x 2", Non-Woven		1																			
Gauze, 2" x 2", Non-Woven										4	4				1						
Gauze, 4" x 4", Non-Woven												2									
Tegaderm, 3.5" x 4.5" (#1655)				1	1																
Tegaderm, 4" x 4.75" (#1616)	1	1	1															1			1
Tegaderm, 4" x 4.75" w/Frame (#1626)												1			1						
Tegaderm, 4" x 4.75" PF (#1626W)																	1				
Tegaderm, CHG 3.5" x 4" (#1657)							1	1						1		1					
Tegaderm, CHG IV (#1658)									1												
Tegaderm, Adv 4" x 4.75" (#1688)										1	1		1								
MediChoice Transparent Dressing 4" x 4.5"						1															
Biopatch Protective Disk w/CHG										1	1	1									
Dressing, Non-Adherent, 2" x 2"																					1
Dressing, Non-Adherent, Split, 2" x 2"																					1
Towels/Drapes																					
Drape, 13" x 19", 2-Ply, T/P, Flat	1	1	1				1	1	1			1	1		1		1	1	1	1	1
Drape, 13" x 19.5", 3-Ply, T/P/T				1	1																
Drape, 17" x 19"						1								1		1					
Drape, 18" x 26", White 6" Fen TPT										1	1										
Antiseptics																					
Alcohol Swabstick, 1/Pk									1						1						1
Alcohol Swabsticks, 3/Pk	1	1	1										1						1	1	
Alcohol Prep Pad, Med												3									
Skin Protectant Pad				1	1		1										1				
Skin Protectant Swab, 1/Pk		1												1							
PVP Ointment		1																			
PVP Swabstick, 3/Pk	1	1																			
Chloraprep, 3mL w/ Insert	1		1	1			1	1	1			2		1	1		1	1	1	1	1
Chloraprep, Frepp 1.5mL w/ Insert																					
Chloraprep 5.25mL Triple w/ Insert					1						1					1					
Prevantics Swabstick 1.6mL w/Insert													1								
Prevantics Swabstick 5.1mL w/Insert						1															
Miscellaneous																					
Durapore Tape 2" x 36"										1	1										
Transpore Tape, .75" x 18"						1															
Transpore Tape, 1" x 18"	1	1													1						1
Tape Strip, .5" x 4", 3/Cd			1				1		1	1	1	1	1								
Transparent Tape .75" x 18"																	1				
Cotton Tip Applicator, Wood, 6"	1																				
Dressing Change Label						1						1					1			1	1
Gloves, Walleled, LF/PF, Medium Pair	1	1	1	1	1		1		1							1			1	1	1
Gloves, Nitrile, Walleled, Medium					2					1	1	1		1	1		1			1	1
Gloves, Next Gen Vinyl, Cuffed, Pair								4					1								
Mask, Pleated Earloops, Blue	1	1	1				1			2	2		1	1	1	1	1	1	1	1	1
Mask, Ear Loop, Yellow						2		2	1			2									1
Tape Measure							1														
Hand Sanitizer															1						
Quantity/Case	30	30	30	30	30	30	30	30	30	20	20	30	30	30	30	30	30	30	30	30	20

Soft Tissue Biopsy Trays

NOT MADE WITH natural rubber latex **diverse supplier** sterile EO

- Conveniently packaged tray
- Preenumbered specimen vials with caps
- Single patient use



Bone Marrow Trays

NOT MADE WITH natural rubber latex **diverse supplier** sterile EO

- Contains both Illinois-style Sternal-Iliac and enhanced J-Style Bone Marrow Biopsy/Aspiration needles
- The J-Style needle has a large, easily gripped, ergonomic handle which provides more torque and needle control
- The I-Style aspiration needle provides depth adjustments from 3/8 Inch to 1 1/8 Inch
- Trays are packed in popular "double-deck" configuration
- Single patient use



Soft Tissue Biopsy Trays

STB651

Components

Multi-Compartment Tray with Prep Well	1
Luer Lock Syringe, 20cc	1
Luer Lock Syringe, 5cc	1
Syringe, 3cc with Needle, 25 Gauge x .625"	1
Needle, 20 Gauge x 1.5"	1
Needle, 22 Gauge x 1.5"	1
Sponge Sticks	2
Gauze Pads	4
Frosted Microscope Slides	6
Scalpel	1
Specimen Vials, Pre-numbered, with Caps, 10 mL	2
Fenestrated Drape with Adhesive Strip	1
Towels	2
Lidocaine Hydrochloride USP, 1%, 5ml with Drug Precaution Sheet	1
Ruler, 6"	1
CSR Wrap	1
Bandage	1
Patient I.D. Labels	2
Universal Medication Labels	10
Quantity/Case	10

Bone Marrow Trays

BMA652

Components

Multi-Compartment "Double Deck" Trays	2
J-Style bone Marrow Biopsy Needle, 11 Gauge x 4"	1
I-Style Aspiration Needle, 15 Gauge, Adjustable from .375" to 1.875"	1
Luer Lock Adapter	1
Specimen Tubes Containing Edetate (EDTA) Potassium	2
Frosted Microscope Slides	10
Amber Specimen Vial with Cap	1
Stirring Rod	1
Povidone-Iodine Swabsticks	3
Lidocaine Hydrochloride, USP, 1%, 5 mL with Drug Precaution Sheet	1
Stylet with Guide	1
Scalpel	1
Luer Lock Syringe, 5cc, with Needle, 20 Gauge x 1.5"	1
Luer Lock Syringe, 20cc	1
Needle, 25 Gauge x .625"	1
Needle, 21 Gauge x 1.5"	1
Towel	1
Fenestrated Drape with Adhesive Strip	1
Gauze Pads	5
Bandage	1
CSR Wrap	1
Patient I.D. Label	1
Extra Cap	1
Universal Medication Labels	10
Quantity/Case	10

Please see our full line of Urethral Catheterization and Foley Trays on page 170 in the Urology Section.





Skin Staple Remover Trays

NOT MADE WITH natural rubber latex | diverse supplier | sterile EO

- Sharp instrument every time
- Easy, stable finger-control
- Fine point for easy access to staple
- Convenient easy-to-open packaging
- Packaged in typical order of use

Skin Staple Remover Trays	SSR1001	SSR1002	SSR1003
Components			
Hard Tray w/ Tyvek Lid	1	1	1
Gauze, 3" x 3", 12-Ply	1	1	1
Skin Staple Remover, Tweezer Style	1		1
Skin Staple Remover, Scissor Style		1	
PVP Pred Pad	1	1	
Quantity/Case	50	50	48



Suture Removal Trays

NOT MADE WITH natural rubber latex | diverse supplier | sterile EO

- Sharp instrument every time
- Easy, stable finger-control
- Fine point for easy access to staple
- Convenient easy-to-open packaging
- Packaged in typical order of use

Suture Removal Trays	SUTT001	SUTT002	SUTT003	SUTT004	SUTT005	SUTT006	SUTT007	SUTT008	SUTT009	SUTT010	SUTT011	SUTT012
Components												
Packaging												
Hard Tray w/Tyvek Lid	1	1	1	1	1	1	1	1	1	1	1	1
Prep And Dressing:												
Gauze, 3" x 3", 8-ply									1			
Gauze, 3" x 3", 12-ply	1	1	1	1	1	1	1	1	1	1	1	1
Instruments												
Scissors, Iris, Straight, 4.5"	1	1		1					1	1		
Scissors, Iris, Straight, 4.5", Satin												1
Scissors, Littauer, 4.5"					1							
Scissors, Littauer, 4.5", Wireform			1			1	1					1
Scissors, Littauer, 4.75", Plastic								1				
Forceps, Thumb, 4", Wireform		1					1		1			1
Forcep, Adson, Serrated, 4.75"				1	1	1				1		
Forcep, Adson, Serrated, 4.75", Satin												1
Forceps, Finepoint, 4.5", Plastic	1											
Forceps, Tissue, Plastic						1		1				
Antiseptics												
Alcohol Prep Pad, Large		1	1	1	1							1
PVP Prep Pad		1	1	1	1							
Quantity/Case	50	50	50	50	50	50	50	50	50	50	50	50

Tracheostomy Trays

NOT MADE WITH natural rubber latex **diverse supplier** **sterile EO**

- Highly absorbent pre-split dressing
- Flexible trach brush for easy passage through inner cannula to minimize irritation
- Angle-cut twill ties to aid in easy threading through cannula slots
- Convenient easy-to-open packaging
- Packaged in typical order of use



Tracheostomy Trays	TRAC1001	TRAC1002	TRAC1003	TRAC1004	TRAC1005	TRAC1007	TRAC1008	TRAC1009
Components								
Packaging								
Hard Tray, 3 Compartment w/Tyvek Lid	1	1				1		
Hard Tray, 2 Compartment							1	
Hard Tray, w/Tyvek Lid			1	1	1			1
Prep and Dressing								
Gauze, Split, 4" x 4", 8-Ply						1	2	
Gauze, Split, 4" x 4", 12-Ply	4	4						
Gauze 4" x 4", 12-Ply				3	4			
Gauze 4" x 4", 8-Ply								4
Excilon Gauze						1		
Trach Sponge, 4" x 4"				1	1			
Trach Dressing, Pre-Sewn	1	1						
Dermacea Sponge, 2-pk								1
Miscellaneous								
Cotton Tip Applicator, 4", Plastic				2				
Cotton Tip Applicator, 6", Plastic	3	3	3		2	3	2	5
Forcep, Plastic, Finepoint			1					
Pipe Cleaner	3	3	3	4	4	2		3
Twill Tape, .5" x 36"	1	1	1	1	1	1		
Tracheostomy Brushes	1	1	1	1	1	1	1	1
Towel Blue Folded 9.5" x 13"						1	1	
Saline Solution						1		
Hydrogen Peroxide, Packet, 1 oz		2						2
Gloves, LF/PF, Cuffed Medium, Pair	1	1	1	1	1	1	1	1
Towel, 13.5" x 17", 3 Ply, T/P/T								1
Drape, 13" x 19", 2 Ply, T/P, Flat	1	1						
Drape, 18" x 15", Blue, 2/Pk			2	1	1			
Sterile Saline, 100ml								1
MediChoice Trach Tube Holder, Adult								1
Quantity/Case	20	20	20	20	20	20	20	20

Irrigation Trays with Bulb Syringe

NOT MADE WITH natural rubber latex **diverse supplier** **sterile EO**

- Easy-to-use syringes for irrigation
- Soft pliable bulb with large finger flanges for one-handed operation and optimum control
- Easy-to-read graduations, 60cc capacity



Irrigation Trays with Bulb Syringe	ITRAY130
Components	
60cc Bulb Syringe w/Tip Protector	1
1200cc Basin	1
500cc Grad Container	1
Alcohol Preparation Pad	1
Underpad	1
Patient Identification Label	1
Quantity/Case	20

Irrigation Trays with Piston Syringe

NOT MADE WITH natural rubber latex **diverse supplier** **sterile EO**

- Easy-to-use syringes for irrigation includes a tip protector and a small luer tip adapter
- Ergonomic design with thumb control ring and large finger flanges for one-handed operation and optimum control
- Smooth-motion plunger with three-tiered, leak-proof seal
- Easy-to-read graduations, 60cc capacity



Irrigation Trays with Piston Syringe	PTRAY126	PTRAY136	PTRAY146
Components			
60cc Piston Syringe w/Tip Protector	1	1	1
1000cc Basin	1		1
1200cc Basin		1	
500cc Grad Container	1		1
Alcohol Preparation Pad		1	
Underpad		1	
Patient Identification Label	1	1	1
Luer Adapter	1	1	1
Waterproof CSR Wrap			1
Quantity/Case	20	20	20

Please see our full line of Trach Tube Holders on page 167 in the Respiratory Section.



IV Start Kits

NOT MADE WITH natural rubber latex **diverse supplier** sterile EO



- Packaged in typical order of use
- Convenient easy-to-open packaging
- Stackable containers for safe, easy storage

IV Start Kits	IVS1006	IVS1007	IVS1011	IVS1013	IVS1026	IVS1029	IVS1052C	IVS1090	IVS1020	IVS1043
Components										
Tray										
Soft Tray (Blister) w/Tyvek Lid	1	1	1	1	1	1	1	1	1	1
Dressings										
Gauze, 2" x 2", 8-Ply	2	2	2	2		2			2	
Gauze, 2" x 2", 12-Ply					1					
Gauze, 2" x 2", 4-Ply, Non-Woven							2	2		
Gauze, 3" x 3", 12-Ply										2
Tegaderm Pediatric 2" x 2.25" (#1610)									1	
Tegaderm PF, 2.375" x 2.75" (#1624)		1		1		1				
Tegaderm PF, 2.375" x 2.75" (#1624W)										1
Tegaderm Adv 2.5" x 2.75" (#1683)							1			
Transparent Dressing, 2.375" x 3.75"								1		
OpSite, IV3000 Dressing (#4017D)	1									
Veni-Guard Dressing					1					
Towels/Drapes										
Wrap, 10" x 13", T/P				1						
Miscellaneous										
Transpore Tape, .75" x 18"	1	1	1	1		1	1	1	1	1
Dressing Change Label	1		1	1		1	1		1	1
Tourniquet LF, .75" x 18"	1	1	1	1	1					
Tourniquet LF, 1" x 18"						1		1	1	1
Tourniquet LF, 1" x 18", Orange							1			
Adhesive Bandage			1							
Lidocaine 1% 2mL						1				
1cc Safety Syringe 29G x .5"						1				
Antiseptics										
Alcohol Prep Pads		2	2	2	1	2	1	2	1	1
PVP Pad		1	1	1				1		
PVP-I Swab Ampule .65cc									1	
Iodine Tincture Swab						1				
Idophor Prep Pad					1					
Quantity/Case	50	50	50	50	100	50	100	50	30	100

IV Start Kits with Prevantics

sterile **EO**

- Includes Prevantics™ antiseptic
- Packaged in typical order of use
- Convenient easy-to-open packaging
- Stackable containers for safe, easy storage



IV Start Kits with Prevantics	IVS1001	IVS8003	IVS1030	IVS1040	IVS1056	IVS1062A	IVS1082	IVS2010	IVS2032	IVS2033	IVS2034	IVS2039	IVS2049	IVS2050
Components														
Tray														
Soft Tray (Blister) w/Tyvek Lid	1	1	1	1	1	1	1	1	1	1	1	1	1	1
Dressings														
Gauze, 2" x 2", 8-Ply		2	2			2	2		2					
Gauze, 2" x 2", 4-Ply, Non-Woven	2			2	2			2		2	2	3	2	2
Tegaderm PF, w/Border 2.375" x 2.75" (#1614)				1	1									
Tegaderm PF, 2.375" x 2.75" (#1624)		1				1		1		1	1			
Tegaderm PF, 2.375" x 2.75" (#1624W)							1							
Tegaderm PF w/Strips 2.75" x 3.25 (#1633)			1											
Tegaderm Advanced 2.5" x 2.75" (#1683)									1			1		
MediChoice Transparent Dressing, 2.375" x 2.75"	1													
Dressing Polyskin II													1	
Towels/Drapes														
Wrap, 10" x 13", T/P						1								
Miscellaneous														
Transpore Tape, .75" x 18"	1	1	1		1		1		1	1	1			
Transpore Tape, .75" x 24"				1								1	1	
Transpore Tape, 1" x 18"						1		1						
Transparent Tape, .75" x 18"														1
Durapore Tape, 1" x 18"														
Dressing Change Label	1	1	1	1	1		1	1	1		1		1	1
Tourniquet LF, .75" x 18"				1		1			1		1		1	1
Tourniquet LF, 1" x 18"	1	1	1		1		1					1		
Tourniquet LF, 1" x 18" Orange										1				
Pre-Filled 3cc Syringe Saline 3mL												1		
Ext Set 7 Macrob w/Bond												1		
Antiseptics														
Prevantics Pad, 1.0 mL w/Insert	1	1	1	1	1	1	1	1	1	1	1	2	1	1
Alcohol Prep Pad							1							
Quantity/Case	50	50	100	100	50	50	50	100	50	50	100	50	100	100

CLEAN ROOM ACCESSORIES

Sterile Alcohol



- Contains 70-percent USP-grade Isopropanol
- Compliant with USP
- Eliminates mixing, filtration, processing and QC expense of in-house production of sterile Isopropanol water blends
- Filtered through 0.2 micron in-line filter meets USP grade purified water
- Gamma-irradiated to a sterility assurance level (SAL) of 10⁻⁶
- Prefiltered, packaged, and sterilized for easy introduction into critical environments
- Induction sealed bottle reduces opportunity for leakage during sterilization and shipment
- Double bagged to preserve product integrity
- Lot number printed on bottle assures complete product traceability
- Certificate of processing and sterility included with each case



Item No.	Description	Quantity
IPA01G	Isopropyl Alcohol 70 Percent Gallon Container Sterile	1/EA, 4 EA/CA

LAB ACCESSORIES

Alcohol Dispenser



- One-hand operation with smooth, pump-top dispenser
- Slide-aside lid designed to prevent spills and keep alcohol fumes from escaping
- A must to disinfect repeat use instruments
- Holds up to 250 ml
- Measured pump prevents over-dispensing
- Amber HDPE plastic bottle protects light sensitive liquids



Item No.	Description	Quantity
LAB2805	Dispenser Alcohol Amber Hdpe 250 Milliliter Case Of 12	1/EA, 12 EA/CA

Task Wipe



- 100-percent Rayon
- Durable and soft: won't leave scratches, lint, or residue
- Safe for most delicate surfaces including instruments, lenses, surfaces, and glass slides
- Lint free
- Non-sterile



Item No.	Description	Quantity
MC10509	Wipe Task 10.5 Inch x 9 Inch Non-Woven White Rayon Cleaner Non-Sterile	30/BX, 32 BX/CA

Lab Marker



- Fine tip (1.0mm)
- Ink is quick-drying and permanent
- Use to write on labware, instrument locking tags, sterilization tape, and sterilization wrap
- Can be run through one autoclave cycle without washing off
- Resists cryogenic temperatures down to -150C (-238F)



Item No.	Description	Quantity
MKR8010A	Marker Laboratory 1 Millimeter Fine Tip Black	1/EA, 10 EA/BX, 50 BX/CA

LAB EQUIPMENT

Centrifuges (Bench-Top)



- Six place fixed angle rotor for 2-15 ml tubes
- Quiet operation (58 decibels at 3,400 RPM)
- Fixed speed model at 3,400 RPM for blood only
- Variable speed model 300-4,000 RPM for blood and urine
- 0-30 minute timer with bell tone signal
- Stainless steel inner guard bowl, easy to disinfect
- Brushless, maintenance-free motor
- Lid window for tachometer RPM check
- Suction cup rubber feet hold securely to bench top (contains latex)
- Dimensions: 10.2 x 11.8 x 8 Inches
- Weight:
CEN206F: 15 lbs.,
CEN206V: 16 lbs.
- Power is cut to the motor if the lid is raised during operation
- Adapters available as accessories for 2-3 ml tubes and 3.5-7 ml tubes
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
Centrifuge Benchtop 6-Place 30 Minute Timer 6 x 10 Milliliter or 3 x 15 Milliliter Capacity		
CEN206F	Fixed Speed 3400 RPM	1/EA
CEN206V	Variable Speed 300-4000 RPM	1/EA

Centrifuges (High-Capacity)



- 24 place fixed angle metal rotor for 2-15 ml tubes
- Quiet operation (58 decibels at 3,400 RPM)
- Variable speed 500-4,000 RPM for blood or urine
- Touch-pad control, 2-line digital display RPM & time
- Zero RPM safety lid latch control system
- 60 Minute timer; lid latch opens at end of run
- Stainless steel inner guard bowl for easy disinfection
- Brushless, maintenance-free motor
- Lid window for tachometer RPM check
- Suction cup rubber feet hold securely to bench top (contains natural rubber latex)
- On-board memory remembers last run for quick reset
- Dimensions: 13.75 x 15 x 11 Inches
- Weight: 31 lbs. (unpacked)
- Adapters available as accessories for 2-3 ml tubes and 3.5-7 ml tubes
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
CEN324V	Centrifuge High-Capacity 24-Place Variable Speed 500-4000 RPM 60 Minute Digital Timer	1/EA

LAB EQUIPMENT

Centrifuge (Horizontal-Spin)

diverse supplier

- 6-place horizontal spin rotor for 2-15 ml tubes
- Variable speed 500-4,000 RPM for blood or urine
- Quiet operation (58 decibels at 3,400 RPM)
- Zero RPM safety lid latch control system
- Touch-pad control, 2-line digital display RPM & time
- 60 Minute timer; lid latch opens at end of run
- Stainless steel inner guard bowl for easy disinfection
- Brushless, maintenance-free motor
- Lid window for tachometer RPM check
- Suction cup rubber feet hold securely to bench top (contains natural rubber latex)
- On-board memory, remembers last run for quick reset
- Dimensions: 13.75 x 15 x 11 Inches
- Weight: 32 lbs. (unpacked)
- Adapters available as accessories for 2-3 ml tubes and 3.5-7 ml tubes
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
CEN506V	Centrifuge Horizontal Spin 6-Place Variable Speed 500-4000 RPM 60 Minute Digital Timer	1/EA

Centrifuges (Microhematocrit)

diverse supplier

- 24-place Microhematocrit tube rotor for capillary tubes with replaceable rubber ring and screw-down metal lid
- Fixed speed: 12,000 RPM, RCF of 15,000 g
- 5 minute timer with tone signal at end of run
- Tube reader included
- All-metal housing is easy to clean and disinfect; acid and reagent resistant
- Hematocrit tray made of plastic composite and aluminum for light weight and durability
- Indicator light illuminates when in use
- Safety switch in latch cuts motor power when lid is lifted
- Air-cooled brush AC motor
- Dimensions: 11.6 x 10.3 x 9.3 Inch
- Weight: 17.4 lbs.
- Lid window for tachometer calibration
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
CEN424V	Microhematocrit Centrifuge 24-Place 12000 RPM 5 Minute Timer with Bell Electric Motor Braking and Safety Latch Reader Included 110V	1/EA

LAB EQUIPMENT

Differential Counters

diverse supplier

- Eight (8) keys with totalizer window
- White blood cell maturation count guide
- Signal bell when 100 cells are counted
- All-metal housing and construction is durable, chemical resistant, and easy-to-clean and disinfect.
- Dual reset knobs; resets to zero
- Alternate red and white keys minimize count errors
- Compact size takes up little bench top space
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
CNT808	Counter Differential 8-Key with Totalizer Window Metal Housing Signal Bell Dual Reset Knobs	1/EA

Incubators

diverse supplier

- Incubator with temperature range of ambient to 60 degrees Celsius
- Compact, space-saving size
- Holds up to 20 standard petri dishes (10/shelf) or up to 6 Liter (0.21 cubic feet) total chamber capacity
- Dual analog temperature rheostat controls for coarse and fine adjustment
- LCD digital display, temp shown to 0.1 degrees Celsius
- Superior insulation characteristics of acrylic housing hold constant temperature to within 0.5 degree Celsius
- Clear acrylic door for viewing cultures inside
- Thermometer access port for calibration check
- Bench top weight only 8.5 lbs
- Removable shelf perforated for enhanced air flow
- Chamber-mounted circulating fan for even temperature distribution
- One-year limited manufacturer's warranty
- Available in either 110V or 220V



Item No.	Description	Quantity
INC2010	Incubator 6 Liter (0.21 Cubic Feet) Capacity Ambient to 60 Celsius Clear Door 6.8 x 7.7 x 7 Inch 110-220V	1/EA

LAB EQUIPMENT

Tube Rockers

diverse
supplier

- Select an 8, 15, or 30 Tube Rocker; tubes up to 120mm in length
- 25 Rocking cycles/min for gentle mixing
- Pause button for adding or removing tubes
- Non-slip reversible and autoclavable pad
- Quiet, nearly-silent operation
- Mix blood tubes, injectables, or analyzer controls
- Provides uniformity of specimens
- Mix specimens before analyzer use or for differential count specimen prep
- Compact and affordable
- Replacement pads available
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
Tube Rocker/Mixer with Pause Button 25 Rocking Cycles Per Minute 110V		
MIX308	8 Tubes/Vials	1/EA
MIX315	15 Tubes/Vials	1/EA
MIX330	30 Tubes/Vials	1/EA

Tube Rotators

diverse
supplier

- Variable speed 25-250 RPM platform rotator
- 11 x 11 Inch platform with replaceable rubber pad
- Humidity cover with moisturizing sponge
- 3-Position switch for Constant-On, Timed Mixing, and Off
- Tie-down flask spring included
- Digital time and speed display
- Digital timer: 0-99 hours, 59 min.
- Indicator bell tone signals end of run
- Holds printed well slides, 96-well immunoassay microtiter plates, RPR and VDRL slides
- Compact, durable, easy to clean and maintain
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
MIX600V	Rotator Digital 11 x 11 Inch Mixing Plate with Humidity Cover Variable Speed 25-250 RPM Digital Display 120V	1/EA

Vortex Mixers

diverse
supplier

- Variable speed up to 3000 RPM
- Dual-mode switch to easily go from Auto-On to Touch mode
- Vortex-mix batched frozen specimens
- Suction cup rubber feet anchor to bench top
- Power-On LED indicator light
- Soft rubber mixing cup
- Very quiet operation
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
MIX720V	Mixer Tube Vortex Variable Speed Up To 3000 RPM	1/EA

MICROSCOPES & ACCESSORIES

Microscopes (Economy)

diverse
supplier

- Economical Binocular Lab Microscope with 10x eyepieces
- Achromatic 4X, 10X, 40X, and 100X (Oil Immersion) objectives
- LED Illumination built-in; no blue filter required
- Graduated mechanical stage with graduated coaxial coarse and fine focus controls
- Iris diaphragm and intensity rheostat for illumination control
- 4-Position reverse nosepiece for stage access
- Focusable NA 1.25 Abbe Condenser
- Stage mounted slide holder
- Dual voltage power supply 110V & 220V
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
MIC1000	Microscope Laboratory Binocular Economy Achromat or Plan 4X 10X 40X 100X Objectives and 10X Eyepieces Integrated Mechanical Stage 1.25 Abbe Condenser Led Illuminator	1/EA

Microscopes (Binocular)

diverse
supplier

- Binocular lab microscope with Achromat or Plan objectives
- Achromat or our best Plan 4X, 10X, 40X, and 100X (oil) objectives
- 10x High-Eyepoint eyepieces with dual diopter control
- Users do not need to remove glasses
- LED Illumination with intensity rheostat light control
- Dual slide holder on graduated mechanical stage
- Focusable NA 1.25 Abbe condenser with Iris diaphragm
- Eyepieces inclined at 30 degree angle
- Graduated coaxial coarse and fine focus controls with ball bearing movement for smooth, precise adjustments
- Dual voltage power supply, 110V & 220V
- Plastic dust cover included
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
Microscope Laboratory Binocular 4X 10X 40X 100X Objectives and 10X Eyepieces Integrated Mechanical Stage 1.25 Abbe Condenser LED Illuminator 110V-220V		
MIC380	Achromat	1/EA
MIC380P	Plan	1/EA

Complementary Products

Item No.	Description
MIC591K	Kit Optical Cleaning

MICROSCOPES & ACCESSORIES

Microscopes (Advanced)



- Advanced binocular or camera-ready trinocular (T)
- Our best plan 4X, 10X, 40X, and 100X (oil) objectives
- Five hole nosepiece, place for an extra objective
- 10x High-Eyepoint eyepieces with dual diopter control
- Trinocular version accepts a camera for EMR applications
- LED Illumination with intensity rheostat light control
- Koehler field iris for expert illumination control
- Dual slide holder on graduated mechanical stage
- Focusable NA 1.25 Abbe condenser with Iris diaphragm
- Eyepieces inclined at 30-degree angle
- Optional phase contrast objectives and condenser
- Graduated coaxial coarse and fine focus controls
- Variety of camera options depending on application
- Universal power supply 110V to 240V
- 5-year limited manufacturer's warranty



Item No.	Description	Quantity
Microscope Laboratory with Seidentopf Head Paired with 10X Eyepiece Set of Plan Objectives 4X 10X 40X 100X Integrated Mechanical Stage Led Koehler Illumination with Intensity Adjustment 1.25 Abbe Condenser Power Input 110V-240V		
MIC804	Binocular	1/EA
MIC804T	Trinocular	1/EA

Microscope Optical Cleaning Kits



- All-inclusive for optics maintenance
- Blow brush to remove dust from lens surfaces
- Lens cleaning solution for coated optics
- Premium lens cleaning tissue
- Soft, lint-free lens cloth
- Place mats for eye pieces
- Hard-sided storage case
- Detailed cleaning instructions included




Item No.	Description	Quantity
MIC591K	Kit Optical Cleaning	1/EA

MICROSCOPES & ACCESSORIES

Slides and Cover Glasses



- Exact size tolerances and absence of ragged edges
- Available in a wide variety of sizes and types to meet all laboratory needs
- All cover glass is cut, sorted, weighed and individually selected for uniform quality, thickness, flatness, and size
- Microscope Slides:
 - Precision-sized for accurate uniformity
 - Each pack is wrapped to protect against humidity and moisture
- Cover Glass Thickness:
 - No. 1 Thickness: 0.13 - 0.17 millimeter
 - No. 1.5 Thickness: 0.16 - 0.19 millimeter
 - No. 2 Thickness: 0.17 - 0.25 millimeter
- Frosted slides feature sandblasted end for marking specimens
- Superfrost slides feature raised white ink-frosted end for marking



Item No.	Description	Quantity
Glass Slide Plain		
MC300006	Economy Chinese Green 25 × 75 × 1 Millimeter	72/BX, 20 BX/CA
OMC630	Standard Swiss White 25 × 75 × 1 Millimeter	72/BX, 20 BX/CA
MC2950	Premium Swiss White 25 × 75 × 1 Millimeter	72/BX, 20 BX/CA
MC2950F	Premium Swiss White 3 × 1 Inch 1 Millimeter	72/BX, 20 BX/CA
Glass Slide Frosted		
MC311006	Economy Chinese Green 25 × 75 × 1 Millimeter	72/BX, 20 BX/CA
MC2951	Premium Swiss White 25 × 75 × 1 Millimeter	72/BX, 20 BX/CA
MC2951F	Premium Swiss White 3 × 1 Inch 1 Millimeter	72/BX, 20 BX/CA
MC2955	Premium Swiss White 25 × 75 × 1.2 Millimeter	72/BX, 20 BX/CA
Glass Slide Superfrost Premium		
MC4951	Swiss White 25 × 75 × 1 Millimeter	72/BX, 20 BX/CA
MC4951F	Swiss White 3 × 1 Inch 1 Millimeter	72/BX, 20 BX/CA
MC4955F	Swiss White 3 × 1 Inch 1.2 Millimeter	72/BX, 20 BX/CA
Glass Cover 1 Ounce Box		
MC22221	22 × 22 Millimeter Size 1	1/OZ, 10 OZ/CA
MC222215	22 × 22 Millimeter Size 1.5	1/OZ, 10 OZ/CA
MC22222	22 × 22 Millimeter Size 2	1/OZ, 10 OZ/CA
MC22401	22 × 40 Millimeter Size 1	1/OZ, 10 OZ/CA
MC22501	22 × 50 Millimeter Size 1	1/OZ, 10 OZ/CA
MC24301	24 × 30 Millimeter Size 1	1/OZ, 10 OZ/CA
MC24501	24 × 50 Millimeter Size 1	1/OZ, 10 OZ/CA
MC245015	24 × 50 Millimeter Size 1.5	1/OZ, 10 OZ/CA
MC24601	24 × 60 Millimeter Size 1	1/OZ, 10 OZ/CA

PHLEBOTOMY ACCESSORIES

Tourniquets (Latex)

- Easy to apply and remove with minimal curling
- 100-percent natural latex with low powder content
- Variety of packaging options available
- High tensile and tear strength
- Single-patient use



Item No.	Description	Quantity
Tourniquet 1 x 18 Inch Contains Natural Rubber Latex Ivory		
TRN6001	Rolled and Banded 250/bag	1/EA, 250 EA/BG, 10 BG/CA
TRN6002	Flat Pack 100/box	1/EA, 100 EA/BX, 10 BX/CA
TRN6003	Flat Pack 250/box	1/EA, 250 EA/BX, 4 BX/CA

Tourniquets (Non-Latex)

NOT MADE WITH natural rubber latex

- Easy to apply and remove with minimal curl to stay on patient's arm and increase patient comfort
- Orange tourniquets provide a visual cue to clinicians that aid in reducing the risk of leaving tourniquets on non-communicative and sensory-impaired patients
- Super premium and premium products have added thickness, elasticity, and tensile strength for added quality and patient comfort
- Rolled and banded tourniquets aid in infection control initiatives; bands are not made with natural rubber latex
- Blue and orange tourniquets are powder-free while white tourniquets have low-powder
- Range of styles and packaging options to meet all needs, including pediatric and bariatric
- Made from strong polyisoprene material
- High tensile and tear strength
- Single patient use



Item No.	Description	Quantity
Tourniquet Polyisoprene		
PC6006	White Textured Super Premium Flat Pack 1 x 18 Inch	1/EA, 250 EA/BX, 4 BX/CA
PC6006OR	Orange Textured Super Premium Flat Pack 1 x 18 Inch	1/EA, 250 EA/BX, 4 BX/CA
TRN1181	Orange Textured Premium Rolled and Banded 1 x 18 Inch	1/EA, 250 EA/BX, 4 BX/CA
TRN1182	Orange Textured Premium Bulk Bag 1 x 18 Inch	1/EA, 250 EA/BX, 4 BX/CA
TRN1183	Blue Textured Premium Rolled and Banded 1 x 18 Inch	1/EA, 100 EA/BX, 10 BX/CA
TRN1184	Blue Textured Premium Flat Pack 1 x 18 Inch	1/EA, 250 EA/BX, 4 BX/CA
TRN1185	Blue Smooth Perforated Roll 1 x 18 Inch	1/EA, 25 EA/BX, 20 BX/CA
TRN1186	Orange Smooth Flat Pack 1 x 18 Inch	1/EA, 25 EA/BX, 40 BX/CA
TRN1187	Orange Smooth Pediatric Bulk Bag 0.375 x 12 Inch	1/EA, 300 EA/BX, 12 BX/CA
TRN1188	Orange Smooth Bariatric Rolled and Banded 1.5 x 32 Inch	1/EA, 40 EA/BX, 12 BX/CA
TRN1190	Orange Textured Perforated Roll 1 Inch x 18 Inches	1/EA, 25 EA/BX, 40 BX/CA
TRN1191	Blue Textured Perforated Roll	1/EA, 25 EA/BX, 40 BX/CA
TRN1192	Blue Textured Flat Pack 1 Inch x 18 Inches	1/EA, 250 EA/BX, 4 BX/CA

PHLEBOTOMY ACCESSORIES

Safety Lancets

NOT MADE WITH natural rubber latex sterile R

- Meets or exceeds OSHA requirement for healthcare safety and compliance
- Patented Comfort Zone Technology combined with side-activation trigger designed to deliver a near pain-free sampling experience for the patient
- Wide-range of needle gauges (18G, 21G, 23G, 28G, and 30G) provides maximum comfort and appropriate blood volumes
- Easy-to-identify colors allow the user to quickly identify and select the correct gauge lancet
- Manufactured in England for exceptional quality and consistency



Item No.	Description	Quantity
Lancet Blood Sampling Safety Unistik 3		
MC5040	Dual-Use 18 Gauge x 1.8 Millimeter Depth Pink	1/EA, 200 EA/BX, 20 BX/CA
MC5060	High Flow 21 Gauge x 2.0 Millimeter Depth Dark Gray	1/EA, 200 EA/BX, 20 BX/CA
MC5030	Normal 23 Gauge x 1.8 Millimeter Depth Medium Flow Light Gray	1/EA, 200 EA/BX, 20 BX/CA
MC5020	Comfort 28 Gauge x 1.8 Millimeter Depth Low Flow White	1/EA, 200 EA/BX, 20 BX/CA
MC5070	Lowest Flow 30 Gauge x 1.5 Millimeter Depth Purple	1/EA, 200 EA/BX, 20 BX/CA

Item No.	Clinical Application	Color
MC5020	For low flow blood glucose testing, HbA1c	Cream with white cap
MC5030	For medium flow blood glucose testing, cholesterol, HbA1c	Cream with light gray cap
MC5060	For high flow blood glucose, cholesterol, HGB, PT/INR and other like tests	Cream with dark gray cap
MC5040	For high flow fingersticks and low flow heelsticks	Cream with pink cap
MC5070	For lowest flow fingerstick blood glucose testing	Cream with purple cap

Phlebotomy Trays

NOT MADE WITH natural rubber latex diverse supplier

- Phlebotomy tray with removable Tube-CUBE™ holders
- Includes 4 Tube-CUBE™ holders: two 9-place for 8-18mm tubes, one 4-place for 26-30mm (50ml) tubes, and one single-place holder for up to 66mm tube
- Includes two hand-Squeeze™ drawing aides
- Pull-out storage tray on bottom for small accessories (e.g.: needles, bandages, etc.)
- Made of durable, easy to clean and disinfect LDPE
- Light and sturdy construction make it easy to carry when full
- Ability to batch samples by route, floor, or time
- Single-place holder perfect for alcohol dispenser, blood-gas ice cup, etc.
- Holds numerous accessories such as butterflies, IV starter, etc. (accessories not included)
- Can accommodate sharps container (not included)



Item No.	Description	Quantity
LAB4870	Tray Phlebotomy with Tube Cube Holder 9.5 x 7 x 15 Inch	1/EA

Item No.	Description
LAB2805	Dispenser Alcohol Amber HDPE 250 Milliliter Case Of 12 Not Made with Natural Rubber Latex MediChoice

PHLEBOTOMY ACCESSORIES

Tube Organizers

NOT MADE WITH natural rubber latex **diverse supplier**

- Three tier tube organizer and dispenser
- Organize by tube type, size, or cap
- Transparent acrylic allows easy visibility to inventory
- Can be wall mounted or placed directly on bench top
- Nine compartments hold up to 210 tubes
- Includes 3 microshelves for short tubes
- Constructed from durable 0.1875-inch acrylic
- Bonus top shelf for accessories, glove box, or other boxed supplies



Item No.	Description	Quantity
LAB5230	Tube Rack Dispenser 3-Tier 9-Compartment Clear Acrylic	1/EA

Tube Racks

NOT MADE WITH natural rubber latex **diverse supplier**

- Each rack holds up to 50 tubes
- Alpha numeric cross reference for tube tracking and batching
- Tightly join racks together by T-joint for batch control
- Select 10-13 millimeter or 13-16 millimeter tube hole sizes
- Made of durable, autoclavable polypropylene
- Will not become brittle in low temperature refrigerators
- Easy-to-clean and disinfect
- Easy-to-label by date, time, shift, and batch
- Packaged in sets of two racks each



Item No.	Description	Quantity
Rack Storage T-Racks Expandable Polypropylene 8.75 x 3.75 x 5.25 Inch		
LAB5450	50 Place for 13-16 Millimeter Tubes	2/PK
LAB5451	50 Place 10-13 Millimeter Tubes	2/PK

POINT OF CARE TESTING

Fecal Occult Test (gFOB)

made in the USA NOT MADE WITH natural rubber latex

- Qualitative method for detecting occult blood, which may be indicative of gastrointestinal disease
- Test results within 60 seconds
- Includes an enhancer making tests more sensitive and readable than standard Guaiac tests
- Detects hemoglobin as low as 0.32 milligrams per 100 milliliters
- 36-month shelf life
- CLIA waived



Item No.	Description	Quantity
----------	-------------	----------

Kit Fecal Occult Blood Test

MCHS100	100 Two Test Slide and 100-Applicator Sticks 2 x 10 Milliliter Developer	100/BX, 20 BX/CA
MCHS1000	1000 Two Test Single Slide and 1000 Applicator Sticks 20 x 10 Milliliter Developer	1000/BX, 3 BX/CA
MCHS34	34 Triple Slide and Applicator Sticks Developer	34/BX, 10 BX/CA
MCHSDV8	8 x 10 Milliliter Developer	1/EA, 8 EA/BX
MCHSPP50	1 Triple Collection Slide and 3 Applicator Sticks with Patient Mailer 50-Tests	50/BX, 16 BX/CA
MHSME-100	Foil Mailing Envelopes 100 Postal Approved for Fecal Specimen	1/EA, 100 EA/BX
HSDVER8	8 x 10 Milliliter Developer Hema Screen Enhanced Readability	1/EA, 8 EA/BX
HSER100	100 Single Slide and 100 Applicator Sticks 2 x 10 Milliliter Developing Solution with 2 Apertures to Sample 2 Different Sites Hema Screen Enhanced Readability	100/BX, 10 BX/CA
HSP50ER	1-Triple Slide 3-Applicator Sticks 3 x 10 Milliliter Developing Solution with Foil Mailing Envelope Hema-Screen Enhanced Readability 50-Tests Temp Controlled 59-86F	50/BX, 16 BX/CA
HSTAT50	50 Single 2-Test Slide and 2 Applicator Sticks Developer Built Into The Slide Hema Screen Stat	50/BX, 10 BX/CA

Fecal Occult Test (iFOB)

made in the USA NOT MADE WITH natural rubber latex

- Qualitative immunoassay test for the determination of human hemoglobin in patient fecal specimens.
- The hema-screen™ SPECIFIC features an innovative sampling method that utilizes one card for collection of two (2) fecal specimens and uses only one (1) immunochemical specimen preparation tube/test cassette system for analysis
- Meets the American Cancer Society and the U.S. Task Force for Colorectal Cancer (CRC) Foundation recommendation of multiple-day testing
- The sampling tab technology replaces unpleasant liquid buffer tubes and eliminates human error due to patient oversampling
- Specific for human hemoglobin, eliminating interactions from red meat, vitamin C, and other factors which may lead to false positives test results
- Excellent 30-day stability
- 24-month shelf life from date of manufacture
- CLIA waived



Item No.	Description	Quantity
----------	-------------	----------

Fecal Occult Blood Test

HSSP-25	Kit Immunochemical Hema-Screen Specific Rapid with Collection Cards and Mailers 25-Tests	25/BX, 10 BX/CA
HSSPCAS	Kit Immunochemical Hema-Screen Specific Rapid Tubes and Cassettes 25-Tests	25/BX, 10 BX/CA
HSSPCON	Control Immunochemical Positive Negative 3 Milliliter	1/EA, 4 EA/CA
HSSPENV	Kit Immunochemical Hema-Screen Specific Patient Mailer with 20 Envelopes	20/BX, 10 BX/CA

POINT OF CARE TESTING

Pregnancy Test (hCG Cards/Cassettes)

made in the USA **NOT MADE WITH natural rubber latex**

- Qualitative immunoassay test for the rapid detection of hCG in human urine or serum for early stage detection of pregnancy
- 99-Percent sensitivity, specificity, and accuracy
- 20 mIU/ml sensitivity in urine specimens
- 25 mIU/ml sensitivity in serum specimens
- Results in 3 minutes (urine); 5 minutes (serum)
- Each one-step pregnancy test is individually sealed in a foil pouch and may be stored at room temperature
- Room temperature storage
- 24-month shelf life
- CLIA waived for urine specimen only



Item No.	Description	Quantity
Kit Human Chorionic Gonadotropin Pregnancy Test		
MC001U	Urine 25-Test	1/EA, 25 EA/BX, 10 BX/CA
MC02SU	Urine/Serum 25-Test	1/EA, 25 EA/BX, 10 BX/CA
HCGCON2	Control Human Chorionic Gonadotropin Refrigerated Positive Negative 2.5 Milliliter	1/EA

Strep A Test

NOT MADE WITH natural rubber latex

- A qualitative immunoassay dipstick test device for the detection of Strep A antigen from direct throat swab
- Detection level is 1.0 x 10⁶ organisms/swab
- 97-Percent Sensitivity, 95-Percent Specificity
- High clinical sensitivity at low colonization levels
- Quick results in 5 minutes
- External positive and negative controls included
- CLIA waived



Item No.	Description	Quantity
IDSTREP-25	Test Immunoassay Dipstick Strep A 25-Tests	25/BX, 10 BX/CA

Contents:	
25 Dipsticks	25 Sterile swabs
25 Disposable extraction test tubes	1 Package insert
1 Reagent A (10ml)	1 Reagent B (10ml)
1 Positive Control	1 Negative Control
1 Procedure card	1 Workstation

POINT OF CARE TESTING

Mono Test

made in the USA **NOT MADE WITH natural rubber latex**

- Qualitative immunoassay for the detection of Infectious Mononucleosis Heterophile Antibodies in whole blood, serum, or plasma
- Single reagent procedure provides easy to read results in 3-8 minutes
- No age limitations
- Overall accuracy 99-percent
- Superior 99.9-percent sensitivity for earlier IM detection
- Specificity of 98.8-percent
- External positive and negative controls included
- CLIA waived for whole blood



Item No.	Description	Quantity
IDM-25	Kit Immunoassay Test Mono for Blood Serum Plasma 25-Tests	25/BX, 10 BX/CA

Contents:	
25 Cassettes	1 Vial Positive Control
25 Disposable Pipettes	1 Vial Negative Control
1 Vial Developing Buffer	1 Product Insert

H Pylori Test

NOT MADE WITH natural rubber latex

- Rapid immunoassay for the qualitative detection of anti-Helicobacter pylori antibody (IgA/ IgG/ IgM) in human whole blood, serum, or plasma
- Sensitivity of 95.9-percent with 91.6-percent overall agreement to biopsy histology
- Easy to use, single-reagent procedure
- Results available in 10 minutes
- CLIA waived for whole blood
- Convenient finger stick sample collection
- External positive and negative controls included



Item No.	Description	Quantity
IHP-10	Kit Helicobacter Pylori Test for Blood Serum Plasma 10-Tests	10/BX, 10 BX/CA

SPECIMEN COLLECTION

Specimen Containers (24 Hour)

NOT MADE WITH natural rubber latex

- Amber container with white screw-on, leak-resistant lid
- Container made of durable polypropylene and cap made of polyethylene
- 3000 ml (0.8 gallon) capacity with imprinted graduation on container
- Integrated handle aids in container transport
- Non-sterile



Item No.	Description	Quantity
SC2400	Container Specimen Amber Urine Collection with Screw Lid 3000 Milliliter Polypropylene Non-Sterile	1/EA, 40 EA/CA

SPECIMEN COLLECTION

Specimen Containers (Red Lab)

NOT MADE WITH natural rubber latex

- Red container with white snap-on lid
- Leak-resistant lid snaps on for easy closure and odor containment
- Lid imprinted with areas to write the cup's contents, patient name, room number, and date
- Non-Sterile
- Single-patient use



Item No.	Description	Quantity
SCR800	Container Specimen Red with Imprinted White Snap Lid Single-Patient Use 8 Ounce 240 Milliliter	1/EA, 25 EA/BG, 10 BG/CA

Specimen Containers (Tamper Proof)

NOT MADE WITH natural rubber latex

- Specimen containers available in two options: sterile fluid pathway or in fully-sterile O.R. peel pouch
- SPC400 has tamper-evident label while SPC405 and SPC405OR have tamper-evident ring
- Patient label included on all products
- Leak-resistant closure
- Graduated in ounces and milliliters on two sides to make it easy to read and measure specimens
- Premium quality polypropylene and polyethylene make the containers durable, crack-resistant, and optically clear
- Textured flat cap surface provides additional writing surface in stackable design format
- SPC405 and SPC405OR:
 - Tri-sensory cap design with visual, audible, and tactile cues ensures proper cap application
 - Audible ratchet sound when lid is tightened
 - Stress-flex cap design promotes increased sealing as pressure builds up within the container
 - 95 kPa compliant for transport in pneumatic tube systems
 - Lot number on each container allows full traceability of container and contents



Item No.	Description	Quantity
Container Specimen Patient Label Tamper-Evident Seal Single-Patient Use		
SPC400	Wide-Mouth Leak-Proof Gray Lid 4 Ounce 120 Milliliter Sterile Fluid Pathway	1/EA, 100 EA/CA
SPC405	Wide-Mouth Leak-Proof Orange Lid 90 Milliliter Sterile Fluid Pathway Pneumatic Tube Approved	1/EA, 100 EA/BG, 4 BG/CA
SPC405OR	Wide-Mouth Leak-Proof Orange Lid 90 Milliliter Sterile In Operating Room Peel Pouch Pneumatic Tube Approved	1/EA, 100 EA/CA

SPECIMEN COLLECTION

Specimen Containers (Urinalysis)

NOT MADE WITH natural rubber latex

- Non-sterile urinalysis specimen containers for single-patient use
- NSC400:
 - 4 oz. (120 ml) polypropylene container with gray plastic screw-on lid
 - Pre-printed label for recording patient information on each cup
 - Embossed graduations in both ounce and milliliter
- NSC650:
 - 6.5 oz. (190 ml) polypropylene container with pre-printed paper lid for recording patient information
 - Cup includes built-in pouring spout to aid sample transfer in lab
 - Frosted panel with embossed graduations



Item No.	Description	Quantity
Container Specimen Urinalysis Non-Sterile Single-Patient Use Polypropylene Clear		
NSC400	Gray Lid 4 Ounce 120 Milliliter	1/EA, 25 EA/SL, 20 SL/CA
NSC650	Frosted Panel with Pre-Printed Paper Lid 6.5 Ounce 190 Milliliter	1/EA, 25 EA/SL, 20 SL/CA

Specimen Bags

NOT MADE WITH natural rubber latex

- Available in a variety of standard sizes in clear or colored film with easy to use zip closure, slide closure, or tamper-evident adhesive closure systems
- Bags with zip closures provide leak-resistant dependability
- All bags include extra outside pocket for paperwork
- Biohazard symbol provides immediate identification that meets federal regulations and OSHA requirements
- Slide and adhesive closures suitable for transporting bodily fluids and tissues (not leak proof)
- Tamper-evident (adhesive closure) bags available for sensitive specimens
- 2 Mil polyethylene material provides sturdy and reliable construction
- Zip closure bags approved for use in pneumatic systems



Item No.	Description	Quantity
Bag Specimen Biohazard Printed Polyethylene with Outside Pocket 2 Mil		
SB0609	Zip-Closure Clear 6 x 9 Inch Leak-Resistant	1/EA, 100 EA/PK, 10 PK/CA
SB0610T	Tamper-Evident Adhesive Clear 6 x 9 Inch	1/EA, 100 EA/PK, 10 PK/CA
SB0810	Zip-Closure Clear 8 x 10 Inch Leak-Resistant	1/EA, 100 EA/PK, 10 PK/CA
MC4510	Zip-Closure Tamper-Evident Adhesive Clear 7 x 9 Inch	1/EA, 100 EA/PK, 10 PK/CA
MZ0609BE	Zip-Closure Blue 6 x 9 Inch Leak-Resistant	1/EA, 100 EA/PK, 10 PK/CA
MC4510	Slide-Closure Clear 7 x 9 Inch Leak-Resistant	1/EA, 100 EA/PK, 10 PK/CA
MC4999	Zip-Closure Clear 9 x 12 Inch Leak-Resistant	1/EA, 100 EA/PK, 10 PK/CA
MC5993	Tamper-Evident Adhesive Clear 10 x 10 Inch	1/EA, 100 EA/PK, 10 PK/CA
MC600STB	Stat Tamper-Evident Adhesive Clear 6 x 10 Inch	1/EA, 100 EA/PK, 10 PK/CA

Item No.	Description	Quantity
Bag Specimen Zip-Closure Biohazard Non-Printed Polyethylene with Outside Pocket 2 Mil		
MC4993	Clear 8 x 8 Inch	1/EA, 100 EA/PK, 10 PK/CA
MC4996	Clear 6 x 9 Inch	1/EA, 100 EA/PK, 10 PK/CA

Baby Shampoo



- No-tear formula
- Rinses quickly and easily
- Gentle cleansing
- pH balanced
- Convenient flip-top cap for easy dispensing
- Compares to national brand tear-free formula
- Quaternium-15 Free



Item No.	Description	Quantity
Shampoo and Baby Bath Tear-Free		
PC8122	Quaternium-15 Free 2 Ounce Baby Fresh Scent	1/EA, 144 EA/CA

Baby Lotion



- Quaternium-15 Free
- Gentle pH balanced formula with conditioners and emollients that hydrate sensitive skin
- Baby fresh scent
- Does not contain animal fat
- Convenient flip-top cap for easy dispensing
- Absorbs quickly and easily and does not leave skin greasy
- Available in 2 ounce and 4 ounce



Item No.	Description	Quantity
Lotion Baby with Flip-Top Cap Q15-Free		
BL4020	2 Ounce	1/EA, 144 EA/CA
BL4040	4 Ounce	1/EA, 96 EA/CA

BABY POWDER

Talc Free Baby Powder



- Soothing, talc-free powder
- Contains food-grade corn starch
- Fresh, clean scent
- Twist-top cap



Item No.	Description	Quantity
Powder Baby Talc-Free Cornstarch		
PC1015	2 Ounce	1/EA, 96 EA/CA
PC1035	4 Ounce	1/EA, 48 EA/CA

Baby Wipes



- Tear-resistant, soft cloth-like material
- Individually treated wipes provide even distribution of moisture
- Available in tub or resealable soft pack
- Baby-powder scent or fragrance-free
- Hypo-allergenic, alcohol free



Item No.	Description	Quantity
Wipes Baby Scented 6.5 Inch x 8.5 Inch Soft Cloth		
88201	Softpack 40 Count	40/PK, 24 PK/CA
88203	Tub 40 Count	40/TB, 24 TB/CA
88205	Softpack with Flip Top Cap 80 Count	80/PK, 12 PK/CA
Wipes Baby Unscented 6.5 Inch x 8.5 Inch Soft Cloth		
88202	Softpack 40 Count	40/PK, 24 PK/CA
88204	Tub 40 Count	40/TB, 24 TB/CA
88206	Softpack Aloe Lanolin 80 Count	80/PK, 12 PK/CA

Ear and Ulcer Syringes



- Used for aspiration and irrigation
- Available sizes include PVC (1 ounce, 2 ounce and 3 ounce) and TPE (2 ounce and 3 ounce) sizes
- Typically use 1 ounce for infants, 2 ounce for pediatrics and 3 ounce for adults
- Ribbed bulb for secure grip and better control
- Smooth tapered tip reduces tissue trauma
- Red ear syringes are more sturdy for increased suction power
- Flat bottom allows syringe to stand upright
- Individually packaged



Item No.	Description	Quantity
Syringe Bulb Ear-Ulcer Polyvinyl Chloride		
EAR301	1 Ounce Infant Green	1/EA, 50 EA/CA
EAR303	3 Ounce Adult Green	1/EA, 50 EA/CA
EAR402	2 Ounce Pediatric Blue	1/EA, 50 EA/CA
Syringe Bulb Ear-Ulcer Thermoplastic Elastomer		
EAR502	2 Ounce Pediatric Red	1/EA, 50 EA/CA
EAR503	3 Ounce Adult Red	1/EA, 50 EA/CA

Infant Heel Warmers

made in the **USA** **NOT MADE WITH** natural rubber latex

- Easy-to-activate, disposable packets provide instant warmth
- Hypoallergenic lift and stick tape holds warmer in place and allows for easy adjustment
- Food-grade, non-toxic heating agents
- Won't harm infants' skin or eyes
- Patented cloth activator disc is easy to activate with less bag ruptures
- Heat to approximately 104 degrees F from room temperature; warmth lasts up to 10 minutes



Item No.	Description	Quantity
69922	Warmer Heel Infant with Heel Strap Clear Liquid with Adhesive Band 3.5 x 5.5 Inch	1/EA, 25 EA/BX, 4 BX/CA

Infant Heel Lancets

NOT MADE WITH natural rubber latex

- Designed to quickly and comfortably sample neonatal blood with minimal pain and discomfort
- Ergonomic design offers multiple grip options
- Meets or exceeds OSHA and NCCLS guidelines for healthcare safety and compliance
- Easy-to-identify colors and size indicators on side of device (Premie/Full Term) make choosing the correct heelstick easy for the user
- Quick release sterility tab helps prevent accidental incisions
- Alignment guides allow for precise placement



Item No.	Description	Quantity
Heelstick Infant Single Use Sterile		
MC5080	Premie Incision Depth 0.85 Millimeter and Length 1.75 Millimeter	1/EA, 50 EA/BX, 6 BX/CA
MC5090	Full-Term Incision Depth 1.00 Millimeter and Length 2.50 Millimeter	1/EA, 50 EA/BX, 6 BX/CA

Fetal Monitoring Belts

made in the **USA** **NOT MADE WITH** natural rubber latex **diverse supplier**

- Soft woven polyester belts stretch easily
- Available in buttonhole with buttons, button hole without buttons, Velcro style or mesh
- Buttonhole style is 2.275 inch wide x 48 inch long, buttonholes placed every .5 inch
- Velcro style is 1.5 inch x 42 inch featuring a 1 inch x 3 inch rounded ends for added safety
- Packaged in pairs containing one pink belt and one blue belt; Mesh belts are both white
- Disposable
- Single Patient Use



Item No.	Description	Quantity
Belt Fetal Monitoring Polyester Mesh Material		
222362	White 2 Inch Wide x 48 Inch Long	2/PK, 50 PK/CA
223240	Blue and Pink Velcro Style Closure 1.5 Inch Wide x 42 Inch Long	2/PK, 50 PK/CA
224220	Blue and Pink Button Closure 2.375 Inch Wide x 48 Inch Long	2/PK, 50 PK/CA
224221	Pink and White Button Hole with Buttons Closure 2.375 Inch Wide x 48 Inch Long	2/PK, 50 PK/CA

ABDOMINAL FETAL MONITORING BANDS

Abdominal Fetal Monitoring Bands

made in the **USA** **NOT MADE WITH** natural rubber latex **diverse supplier**

- Seamless tubular stockinette style bands
- Tubular style holds monitoring devices in place
- Soft woven material for added comfort
- Flexible and easy to use
- Great alternative to fetal monitoring belts and straps



Item No.	Description	Quantity
Band Monitoring Abdominal Fetal Stockinette Style Woven		
FMB10PRE	Pre-Cut 10 Inch Wide x 16 Inch Long	1/EA, 50 EA/CA
FMB10RL	11 Yard Roll 10 Inches Wide	1/RL
FMB12PRE	Pre-Cut 12 Inch Wide x 16 Inch Long	1/EA, 50 EA/CA
FMB12RL	11 Yard Roll 12 Inches Wide	1/RL
FMB14PRE	Pre-Cut 14 Inch Wide x 16 Inch Long	1/EA, 50 EA/CA
FMB14RL	11 Yard Roll 14 Inches Wide	1/RL
FMB16PRE	Pre-Cut 16 Inch Wide x 16 Inch Long	1/EA, 50 EA/CA
FMB16RL	11 Yard Roll 16 Inches Wide	1/RL

Please see our full line of Safety Lancets on page 122 in the Laboratory Section.

Baby Comb and Brush

NOT MADE WITH
natural rubber latex



- Baby Brush and Comb Set
 - Extra-soft nylon bristle brush
 - Fine-tooth plastic comb
 - Gender neutral yellow color
 - Individually packaged in poly bag
- Baby Comb
 - Two-sided comb, fine teeth on one side and extra-fine teeth on the other

Item No.	Description	Quantity
BACBRUSH	Brush and Comb Set Baby Extra Soft Nylon Bristle Brush	1/EA, 24 EA/BX, 6 BX/CA
BACOMB	Baby Comb Fine Tooth 2-Sided	1/EA, 12 EA/PK, 12 PK/CA

Baby Bonnets

made in the USA

NOT MADE WITH
natural rubber latex

diverse supplier



- Knitted polyester yarn construction for softness, comfort and stretch
- Assists in reducing the infant's loss of body temperature after delivery
- Pink, blue and white stripe design helps reduce inventory SKUs
- Available in single-ply and double-ply

Item No.	Description	Quantity
Bonnet Baby Tri-Color with Pink Blue and White Polyester Cotton Blend 4 Inch Wide x 6 Inch Long		
99955	Double Ply	1/EA, 25 EA/BG, 8 BG/CA
99956	Single Ply	1/EA, 25 EA/BG, 8 BG/CA

Baby Booties

made in the USA

NOT MADE WITH
natural rubber latex



- Size Range: 0-3 Months
- Made with fibers of Polyester and Spandex for warmth and comfort
- Terrycloth inside for additional softness against baby's skin
- Adequate stretch for ease of fit, but the booty is not likely to fall off
- Each pair packaged in individual polybags, 48 pairs in each case

Item No.	Description	Quantity
Bootie Baby Individually Bagged Pair 95% Polyester and 5% Spandex		
BBT100	Blue	1/PR, 48 PR/CA
BBT200	Pink	1/PR, 48 PR/CA
BBT300	White	1/PR, 48 PR/CA

NURSERY KITS



MediChoice Nursery Kits are an inexpensive way to prepare a new mom with a few basic baby grooming supplies that can be used while in the hospital and at home. The kits contain a bulb syringe, a soft pack of unscented baby wipes, a fine-tooth comb for gentle grooming, a digital thermometer, and a 2-ounce bottle of tearless baby shampoo.

Customization Options:

As an added convenience, MediChoice offers customizable Nursery Kits to help improve the patient care experience and strengthen your healthcare facility's brand name. The Kits can be tailored to meet your specific needs.

- Items in the kit come pre-packed by MediChoice
- Contents can be customized to patient's needs
- Kits can be designed with your name, logo and colors

Nursery Kits

NOT MADE WITH
natural rubber latex

diverse supplier

sterile EO

- 60 cubic centimeters, EO sterile bulb syringe
- 80 count, unscented soft pack baby wipes
- Fine tooth comb
- Digital Thermometer
- 2 ounce, tearless baby shampoo



Item No.	Description	Quantity
MBKT101	Kit Basic Nursery NICU Components Bulb Syringe Baby Wipes Baby Comb Oral Digital Thermometer and 2 Ounce Baby Shampoo	1/EA, 12 EA/CA



Please see our full line of Quiet Kits on page 13 in the Admissions Section.

Umbilical Cord Clamp Cutters

NOT MADE WITH natural rubber latex **diverse supplier**



- Comfortable and stable handles
- Cuts easily through the plastic cord clamp
- Individually packaged
- Safe
- Single use only
- Non-sterile

Item No.	Description	Quantity
88100	Clipper Umbilical Cord Clamp Disposable Non-Sterile Single Use	1/EA, 6 EA/BX, 10 BX/CA

Umbilical Cord Clamps

NOT MADE WITH natural rubber latex **diverse supplier** **sterile EO**



- Double-grip surface ensures firm, continuous pressure and prevents seeping
- Secure snap closure
- Hinge guard prevents umbilical cord from slipping into hinge area
- Individually packaged in easy-to-open soft pouches
- Includes illustrated instructions
- Non-toxic, molded plastic

Item No.	Description	Quantity
CLMP101	Clamp Umbilical Cord Double Grip Sterile	1/EA, 6 EA/BX, 5 BX/CA

Baby Scale Liners

made in the USA **NOT MADE WITH** natural rubber latex



- Two ply embossed tissue for comfort, softness and absorbency
- Poly backing provides added strength and prevents fluid leak-through
- Special sizing options means scale is fully covered which protects both baby and scale

Item No.	Description	Quantity
Liner Baby Scale Two Ply Embossed Tissue with Poly Backing White		
MBSL1322	13 Inch x 22 Inch	1/EA, 250 EA/CA
MBSL1325	13 Inch x 25 Inch	1/EA, 250 EA/CA

Footprinters

NOT MADE WITH natural rubber latex **diverse supplier**



- Molded plastic printers with rigid backing for quick and easy impression
- High-quality, washable, non-toxic ink for clear, accurate print
- Neatly stacked in convenient dispenser box

Item No.	Description	Quantity
88102	Footprinter Impression Rigid Disposable	1/EA, 25 EA/BX, 4 BX/CA

PERINEAL WASHES

Perineal Wash Spray

made in the USA **NOT MADE WITH** natural rubber latex **diverse supplier**



- No rinse formula cleans, conditions and moisturizes delicate skin
- pH balanced to better match sensitive skin
- Refreshing, floral fragrance
- Pump can be used at any angle, including upside-down
- Can also be used for ostomy care

Item No.	Description	Quantity
Spray Perineal Wash without Rinse		
PW2004	4 Ounce	1/EA, 48 EA/CA
PW2085	8.5 Ounce	1/EA, 48 EA/CA
PPB1080	Irrigation Bottle 8 Ounce	1/EA, 50 EA/CA

Perineal Cold Packs

NOT MADE WITH natural rubber latex



- Extra-absorbent OB pad with adhesive strip to keep pack in place
- Supplies uniform therapeutic temperatures for postpartum relief
- Helps reduce discomfort and swelling
- Designed for quick, easy activation for instant relief
- Degree Range: Cools to approximately 33° F/.5° C
- Maintains a safe therapeutic temperature for up to one hour
- Single use

Item No.	Description	Quantity
69902	Cold Pack Therapy Perineal and Extra-Absorbent OB Pad with Adhesive Strip 4 Inch x 11 Inch Disposable	1/EA, 24 EA/CA

Please see our full line of Baby Scales on page 37-38 in the Equipment Section.

Obstetric Pads



- Secure adhesive strip on all designs
- Unscented
- Moisture barrier for protection and comfort
- Peach pad is contoured for optimal absorption, fit and performance
- Overnight pad provides superior absorbency and extra protection when lying down
- Super Pads are individually wrapped for greater sanitation and easy dispensing
- Super Pads have eight-channel design for increased leakage protection



Item No.	Description	Quantity
Pad Obstetric with Adhesive Strip		
OB0308MX	Maxi 2.75 Inch x 8 Inch	12/BG, 24 BG/CA
OB03110L	Overnight Long 3 Inch x 11 Inch	12/BG, 24 BG/CA
Pad Obstetric Individually Wrapped with Adhesive Strip		
OB0613PC	Contoured Peach 7 Inch x 14 Inch Individually Wrapped	1/EA, 20 EA/BG, 10 BG/CA
OB0310SL	2.75 Inch x 10 Inch Individually Wrapped	1/EA, 24 EA/BG, 12 BG/CA

Mesh Briefs



- Soft, breathable mesh
- Holds pads in place
- Available in sizes Large to 4XL
- Color coded stitching and waistband indicates size
- Available in bulk, 2 each in a bag and 5 each in a bag



Item No.	Description	Quantity
Briefs Mesh		
MB4002	Large Brown Waistband	2/PK, 50 PK/CA
MB4005	Large Brown Waistband	5/PK, 20 PK/CA
MB4100	Large Brown Waistband	1/EA, 100 EA/CA
MB5002	Extra Large Green Waistband	2/PK, 50 PK/CA
MB5005	Extra Large Green Waistband	5/PK, 20 PK/CA
MB5100	Extra Large Green Waistband	1/EA, 100 EA/CA
MB6002	2X Large Purple Waistband	2/PK, 50 PK/CA
MB6005	2X Large Purple Waistband	5/PK, 20 PK/CA
MB6100	2X Large Purple Waistband	1/EA, 100 EA/CA
MB7002	3X Large Light Blue Waistband	2/PK, 50 PK/CA
MB7005	3X Large Light Blue Waistband	5/PK, 20 PK/CA
MB7100	3X Large Light Blue Waistband	1/EA, 100 EA/CA
MB8002	4X Large White Waistband	2/PK, 50 PK/CA
MB8005	4X Large White Waistband	5/PK, 20 PK/CA
MB8100	4X Large White Waistband	1/EA, 100 EA/CA


Long Torso Seamless Knit Pants



- High waist long torso stretch brief
- Finished leg seams
- Adaptable to hold incontinence liners or pads in place
- Eliminates seam lines under clothing for patient dignity
- Sized from small to bariatric pairs, two pack, five pack and bulk cases of 100 each
- Constructed of 94% Polyester and 6% Spandex
- Great for use in Postpartum Care, Gynecological Needs, Incontinence Care and Postoperative Care



Item No.	Description	Quantity
Brief Seamless Long Torso Knit Pant 94% Polyester and 6% Spandex		
MBLT1001	Brown and Blue Waistband Fits Size Small to Medium	1/EA, 100 EA/CA
MBLT1002	Brown and Blue Waistband Fits Size Small to Medium	2/PK, 50 PK/CA
MBLT1005	Brown and Blue Waistband Fits Size Small to Medium	5/PK, 20 PK/CA
MBLT1101	Yellow and Green Waistband Fits Size Large to Extra Large	1/EA, 100 EA/CA
MBLT1102	Yellow and Green Waistband Fits Size Large to Extra Large	2/PK, 50 PK/CA
MBLT1105	Yellow and Green Waistband Fits Size Large to Extra Large	5/PK, 20 PK/CA
MBLT1201	Purple and Gray Waistband Fits Size 2XL to 3XL	1/EA, 100 EA/CA
MBLT1202	Purple and Gray Waistband Fits Size 2XL to 3XL	2/PK, 50 PK/CA
MBLT1205	Purple and Gray Waistband Fits Size 2XL to 3XL	5/PK, 20 PK/CA
MBLT1301	White Waistband Fits Size 4XL to Bariatric	1/EA, 100 EA/CA
MBLT1302	White Waistband Fits Size 4XL to Bariatric	2/PK, 50 PK/CA
MBLT1305	White Waistband Fits Size 4XL to Bariatric	5/PK, 20 PK/CA



BATHROOM SAFETY AIDS

Raised Toilet Seats

NOT MADE WITH natural rubber latex

- White, heavy-duty plastic construction is lightweight, strong, and durable
- Front clamping mechanism locks securely to the toilet and prevents shifting of the seat
- Wide, contoured seating surface is comfortable
- Available with push-button, removeable armrests (RTS200) for secure grip, comfort, and support
- 300-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
Seat Raised Toilet Front Clamping 300-pound Evenly Distributed Weight Capacity Heavy-duty Plastic Seat		

RTS100	16.5 Inch Wide × 16.5 Inch Deep × 4.5 Inch High	1/EA, 3 EA/CA
RTS200	20.75 Inch Overall Width 18.5 Inch Wide Between Arms 16.5 Inch Deep × 4 Inch High	1/EA, 2 EA/CA

Adjustable Height Raised Toilet Seat

NOT MADE WITH natural rubber latex

- Features four clip-on brackets that fit most toilets
- Adjustable angle and height maximize comfort and reduce bending and patient effort
- Solid molded plastic seat ensures durability
- VersaGuard brackets will not scrape or mar toilets
- 300-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
RTS300	Seat Raised Toilet Adjustable Height 14.75 Inch Wide × 15.5 Inch Deep Adjustable Height from 3.75 Inch to 5.75 Inch Versaguard Bracket	1/EA

Bath Seats

- Lightweight, durable, and rust-resistant anodized aluminum frame
- Platinum grey plastic seat available with optional backrest (BTH231)
- Seat height adjustable in 1-inch increments
- Delivered assembled with tool-free back attachment
- Skid-resistant leg tips contain natural rubber latex
- 350-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
Overall 20 Inch Wide × 18 Inch Deep Seat Area 20 Inch Wide × 12 Inch Deep Height Adjustable from 14 to 20 Inch		
BTH230	Seat Bath	1/EA, 4 EA/CA
BTH231	Seat Bath with Backrest	1/EA, 4 EA/CA

BATHROOM SAFETY AIDS

Bariatric Bath Seats

- Durable, blow-molded plastic seat available with optional backrest (BTH251)
- Lightweight, durable, and rust-resistant anodized aluminum frame
- Seat height adjustable in 1-inch increments
- Platinum grey color
- Delivered assembled with tool-free back attachment
- Skid-resistant leg tips contain natural rubber latex
- 600-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
----------	-------------	----------

Seat Bath Bariatric Contains Natural Rubber Latex		
BTH250	Without Backrest Overall 18 Inch Wide × 16.75 Inch Deep Seat Area 17.75 Inch Wide × 12 Inch Deep Height Adjustable from 14 to 17 Inch	1/EA, 2 EA/CA
BTH251	With Backrest Overall 18 Inch Wide × 20.5 Inch Deep Seat Area 17 Inch Wide × 16.25 Inch Deep Height Adjustable from 16 to 20 Inch	1/EA, 2 EA/CA

Transfer Bench with Backrest

- Provides assistance to those who have difficulty stepping over bathtub walls
- Lightweight, durable, rust-resistant anodized aluminum frame
- Drainage holes in seat minimize water build-up
- Side support arm for lateral transfer
- Adjustable seat height in half-inch increments
- Reversible back attachment to accommodate any bathroom
- Platinum grey color
- Comes assembled with tool-free arm and back attachments
- Skid-resistant leg tips contain natural rubber latex
- 350-pound capacity with even distribution of weight
- Lifetime limited manufacturer's warranty



Item No.	Description	Quantity
BTH441	Bench Transfer with Backrest Overall 29.25 Inch Wide × 22 Inch Deep Seat Area 26.5 Inch Wide × 16 Inch Deep Height Adjustable from 18.25 to 22.25 Inch 350-pound Evenly Distributed Weight Capacity Platinum Grey Plastic Seat Anodized Aluminum Frame Contains Natural Rubber Latex	1/EA, 2 EA/CA

BATHROOM SAFETY AIDS

Bariatric Transfer Bench with Backrest

- Durable blow-molded plastic bench and backrest with extra deep seat for added stability
- Drainage holes in seat minimize water buildup
- Lightweight, durable, and rust-resistant anodized aluminum frame
- Seat height adjustable from 16 inches to 20 inches in 1-inch increments
- Platinum grey color
- Delivered assembled with tool-free back attachment
- Back attachment reverses to accommodate any bathroom
- Skid-resistant leg tips contain natural rubber latex
- 600-pound capacity with even distribution of weight
- Lifetime limited warranty



Item No.	Description	Quantity
BTH451	Bench Transfer Bariatric with Backrest Overall 29.5 Inch Wide x 20.5 Inch Deep Seat Area 26.75 Inch Wide x 16.25 Inch Deep Height Adjustable from 16 to 20 Inch 600-pound Evenly Distributed Weight Capacity Platinum Grey Plastic Seat Anodized Aluminum Frame Contains Natural Rubber Latex	1/EA, 2 EA/CA

Shower Transport Chairs

NOT MADE WITH
natural rubber latex

- Designed for safe transport to and from the shower
- Rust-resistant, anodized, non-corrosive aluminum frame
- Comfortable nylon backrest removes for easy cleaning
- Cushioned armrests enhance user comfort
- 5-inch swivel casters with rear wheel safety locks facilitate maneuverability in tight spaces
- Rolls easily over most toilets
- 250-pound capacity with even distribution of weight
- Limited lifetime manufacturer's warranty



Item No.	Description	Quantity
BTH310	Shower Chair Transport Aluminum Frame Removable Nylon Backrest Overall 21.5 Wide x 24 Deep x 39 High Seat Area 14 Wide x 16.5 Deep x 22 High 5 Inch Locking Casters	1/EA

BATHROOM SAFETY AIDS

PVC Shower Commode Chairs

NOT MADE WITH
natural rubber latex

- Easy to assemble and break down without tools for convenient transport
- Constructed of the highest quality, healthcare-grade PVC
- Deluxe elongated, extra heavy-duty, contoured open front seat
- Dual reinforcement in the caster fittings for extra support and threaded stem casters for additional safety
- Quick drying, removable royal blue mesh back support
- Ergonomic curved arms provide support without sharp edges and include anti-slip handgrips for added safety
- Clamp Joint Bracket Locking System is made of high-grade plastic and stainless steel
- 300-pound capacity with even distribution of weight
- One-year limited warranty on PVC frame; six-month limited warranty on all other materials



Item No.	Description	Quantity
BTH810	Shower Chair PVC Removable Blue Mesh Back Overall 18 Inch Wide x 22 Inch Deep x 40 Inch High Seat Area 18 Inch Wide x 18 Inch Deep x 21 Inch High	1/EA

COMMODOE CHAIRS

Folding Steel Commodes

NOT MADE WITH
natural rubber latex

- Versatile, 3-in-1 design allows commode to be used as a bedside commode, a toilet safety frame, or a raised toilet seat for a standard toilet
- Folds flat to less than 4 inches for convenient storage and travel
- Welded steel frame provides strength and stability
- Sturdy plastic seat snaps on and off easily
- Height adjusts from 17 inches to 23 inches in 1-inch increments
- Comes complete with commode pail, cover, and splash shield
- Tool-free assembly
- 300-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
COMM7001	Commode Folding Welded Steel Frame Plastic Snap-On Seat Height Adjusts from 17 Inch to 23 Inch in 1-inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield	1/EA, 4 EA/CA

Complementary Products

Item No.	Description
COMM7007	Commode Pail Set 7-Quart Capacity with Cover Gray Fits Most Standard Commodes
COMM7008	Commode Pail Set Autoclavable 134 Degrees Celsius for 3 Minutes Polypropylene Plastic 7-Quart Capacity with Cover and Handle White Fits Most Standard Commodes
COMM7009	Commode Pail Set 12-Quart Capacity with Cover and Handle Gray Fits Most Standard Commodes



Please see additional Fall Management products, including our full line of Fall Management Slippers on page 17 in the Apparel-Patient Section.

COMMODOE CHAIRS

Painted Steel Commodes

NOT MADE WITH natural rubber latex

- Removable back bar facilitates use as a bedside commode, raised toilet seat, or toilet safety frame
- Welded, powder-coated steel frame provides strength and stability and is rust resistant and easy to clean
- Back bar detaches to accommodate use as toilet safety frame
- Sturdy plastic seat snaps on and off easily
- Seat height adjusts from 15.5 inches to 23.5 inches in 1-inch increments
- Comes complete with commode pail, cover, and splash shield
- Contoured arm rests enhance user comfort
- 300-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
COMM6001	Commode Painted Welded Steel Frame Plastic Snap-On Seat Height Adjusts from 15.5 Inch to 23.5 Inch in 1-Inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield	1/EA, 4 EA/CA

Complementary Products

Item No.	Description
COMM7007	Commode Pail Set 7-Quart Capacity with Cover Gray Fits Most Standard Commodes
COMM7008	Commode Pail Set Autoclavable 134 Degrees Celsius for 3 Minutes Polypropylene Plastic 7-Quart Capacity with Cover and Handle White Fits Most Standard Commodes
COMM7009	Commode Pail Set 12-Quart Capacity with Cover and Handle Gray Fits Most Standard Commodes

COMMODOE CHAIRS

Drop Arm Commodes

NOT MADE WITH natural rubber latex

- Arms drop independently below the seat level for safe lateral transfers
- Easy-to-release mechanism can be operated by patients lacking manual dexterity
- Arms lock automatically when in the raised position for added safety
- Welded, powder-coated, steel frame is durable, rust resistant, and easy to clean
- Removable back bar facilitates use as a bedside commode, raised toilet seat, or toilet safety frame
- Cushioned armrests enhance user comfort
- Height adjusts from 18.75 inches to 22.75 inches in 1-inch increments
- Comes complete with commode pail, cover, and splash shield
- 300-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
COMM6002	Commode Drop-Arm Welded Powder-Coated Steel Frame Plastic Snap-On Seat Height Adjusts from 18.75 Inch to 22.75 Inch in 1-Inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield	1/EA

Complementary Products

Item No.	Description
COMM7007	Commode Pail Set 7-Quart Capacity with Cover Gray Fits Most Standard Commodes
COMM7008	Commode Pail Set Autoclavable 134 Degrees Celsius for 3 Minutes Polypropylene Plastic 7-Quart Capacity with Cover and Handle White Fits Most Standard Commodes
COMM7009	Commode Pail Set 12-Quart Capacity with Cover and Handle Gray Fits Most Standard Commodes

Heavy-Duty Aluminum Commodes

NOT MADE WITH natural rubber latex

- Versatile, 3-in-1 design allows commode to be used as a bedside commode, a toilet safety frame, or a raised toilet seat for a standard toilet
- Aluminum frame is lightweight and easy to transport
- Sturdy plastic seat snaps on and off easily
- Height adjusts from 18 inches to 23 inches in 1-inch increments
- Comes complete with commode pail, cover, and splash shield
- Contoured arm rests enhance user comfort
- 400-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
COMM7002	Commode Heavy-Duty Aluminum Plastic Snap-on Seat Height Adjusts from 18 Inch to 23 Inch in 1-inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield	1/EA, 2 EA/CA

Complementary Products

Item No.	Description
COMM7007	Commode Pail Set 7-Quart Capacity with Cover Gray Fits Most Standard Commodes
COMM7008	Commode Pail Set Autoclavable 134 Degrees Celsius for 3 Minutes Polypropylene Plastic 7-Quart Capacity with Cover and Handle White Fits Most Standard Commodes
COMM7009	Commode Pail Set 12-Quart Capacity with Cover and Handle Gray Fits Most Standard Commodes

Bariatric Drop Arm Commodes

NOT MADE WITH natural rubber latex

- Arms drop independently below the seat level for safe lateral transfers
- Easy-to-release mechanism can be operated by patients lacking manual dexterity
- Arms lock automatically when in the raised position for added safety
- Welded-steel construction and extra-wide base provide added stability
- Removable back bar facilitates use as a bedside commode, raised toilet seat, or toilet safety frame
- Extra-wide platform seat with contoured and indented front is comfortable for long periods of sitting
- Cushioned armrests enhance user comfort
- Height adjusts from 19.5 inches to 23 inches in 1-inch increments
- Comes complete with commode pail, cover, and splash shield
- 600-pound capacity with even distribution of weight
- Limited lifetime warranty



Item No.	Description	Quantity
BAR706	Commode Bariatric Drop-Arm Welded Steel Height Adjusts from 19.5 inch to 23 inch in 1-inch Increments 23 inch Seat Width Includes Pail with Cover and Splash Shield	1/EA

Complementary Products

Item No.	Description
COMM7007	Commode Pail Set 7-Quart Capacity with Cover Gray Fits Most Standard Commodes
COMM7008	Commode Pail Set Autoclavable 134 Degrees Celsius for 3 Minutes Polypropylene Plastic 7-Quart Capacity with Cover and Handle White Fits Most Standard Commodes
COMM7009	Commode Pail Set 12-Quart Capacity with Cover and Handle Gray Fits Most Standard Commodes

COMMODE CHAIRS

Commode Pail Sets

NOT MADE WITH
natural rubber latex



- Designed to fit all standard commodes as well as other leading manufacturers' commodes
- Cover included
- Available in 7- or 12-quart capacity
- 12-quart and autoclavable pails include easy-grip handles
- Autoclavable pail (COMM7008) designed to withstand autoclave temperatures up to 134 degrees Celsius for 3 minutes

Item No.	Description	Quantity
Commode Pail Sets Fit Most Standard Commodes		
COMM7007	7-Quart Capacity with Cover Gray	1/EA, 6 EA/CA
COMM7008	7-Quart Capacity with Cover and Handle White Autoclavable	1/EA, 6 EA/CA
COMM7009	12-Quart Capacity with Cover and Handle Gray	1/EA, 6 EA/CA

Complementary Products

Item No.	Description
COMM6001	Commode Painted Welded Steel Frame Plastic Snap-on Seat Height Adjusts from 15.5 Inch to 23.5 Inch in 1-inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield 300-Pound Evenly Distributed Weight Capacity
COMM6002	Commode Drop-Arm Welded Powder-coated Steel Frame Plastic Snap-on Seat Height Adjusts from 18.75 Inch to 22.75 Inch in 1-inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield 300-Pound Evenly Distributed Weight Capacity
COMM7001	Commode Folding Welded Steel Frame Plastic Snap-on Seat Height Adjusts from 17 Inch to 23 Inch in 1-inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield 300-Pound Evenly Distributed Weight Capacity
COMM7002	Commode Heavy-Duty Aluminum Plastic Snap-on Seat Height Adjusts from 18 Inch to 23 Inch in 1-inch Increments 14-Inch Seat Width Includes Pail with Cover and Splash Shield 400-Pound Evenly Distributed Weight Capacity
BAR706	Commode Bariatric Drop-Arm Welded Steel Height Adjusts from 19.5 Inch to 23 Inch in 1-inch Increments 23 Inch Seat Width Includes Pail with Cover and Splash Shield 600-Pound Evenly Distributed Weight Capacity

CRUTCHES

Aluminum Crutches

NOT MADE WITH
natural rubber latex



- Anodized aluminum construction for durability and strength
- Solid composite cuff secures I-beam foot piece to crutch frame for added stability and safety
- Easy-to-read height scale and simple push-button mechanism adjusts crutch height in 1-inch increments
- Handgrip position easily adjusted with wing nuts
- Comfortable and durable underarm cushions and handgrip pads
- Skid-resistant vinyl tips provide excellent traction
- 300-pound capacity at a 45-degree angle

Item No.	Description	Quantity
Crutch Aluminum Push-button Height Adjustment		
77708A	Adult Accommodates 5 feet 2 inch to 5 feet 10 inch in 1-inch Increments	1/PR, 8 PR/CA
77709A	Child Accommodates 4 feet to 4 feet 6 inch in 1-inch Increments	1/PR, 8 PR/CA
77710A	Tall Accommodates 5 feet 10 inch to 6 feet 6 inch in 1-inch Increments	1/PR, 8 PR/CA
77711A	Youth Accommodates 4 feet 6 inch to 5 feet 2 inch in 1-inch Increments	1/PR, 8 PR/CA

Accessories

Item No.	Description
MAA1001	Cushion Underarm Crutch Pair Fits MediChoice Aluminum Crutches 77708A, 77709A, 77710A, and 77711A
MAA1002	Handgrip Crutch Pair Fits MediChoice Aluminum Crutches 77708A, 77709A, 77710A, and 77711A
MAA1003	Tip Crutch Pair Fits MediChoice Aluminum Crutches 77708A, 77709A, 77710A, and 77711A

Bariatric Steel Crutches

NOT MADE WITH
natural rubber latex

- Heavy-duty steel construction for durability and strength
- Easy-to-read height scale and simple push-button mechanism adjusts crutch height in 1-inch increments
- Handgrip position easily adjusted with wing nuts
- Comfortable and durable underarm cushions and handgrip pads
- Skid-resistant vinyl tips provide excellent traction
- 600-pound capacity with even distribution of weight
- Limited lifetime manufacturer's warranty



Item No.	Description	Quantity
Crutch Bariatric Steel Push-button Height Adjustment		
BAR206R	Adult Accommodates 5 feet 2 inch to 5 feet 10 inch in 1-inch Increments	1/PR
BAR206T	Tall Accommodates 5 feet 10 inch to 6 feet 6 inch in 1-inch Increments	1/PR

WALKERS

Single-Release Folding Walkers

NOT MADE WITH
natural rubber latex

- Scratch-resistant, 1-inch aluminum tubing provides maximum strength while remaining lightweight
- One-button mechanism allows quick and easy, single-handed folding and unfolding
- Folds to depth of 4 inches for convenient storage
- Rugged steel cross-bars add strength without interfering with patient stride
- Height adjustable from 32 inches to 36 inches
- 250-pound weight capacity



Item No.	Description	Quantity
77741	Walker Folding Single-Release 1-inch Aluminum 22 Inch Wide × 18.5 Inch Deep Height Adjustable from 32 to 36 Inch 250-Pound Weight Capacity	1/EA, 4 EA/CA

Accessories

Item No.	Description
MAA2001	Tips Walker Pair Fits Walkers 77741, 77742, 77743, and WLK4500
MAA2002	Wheels Walker 5-inch Pair Fits Walkers 77741, 77742, 77743, and WLK4500
MAA2003	Glide Caps Walker Pair Fits Walkers 77741, 77742, 77743, and WLK4500
MAA2004	Ski Glides Walker Pair Fits Walkers 77741, 77742, 77743, and WLK4500

Dual-Release Folding Walkers

NOT MADE WITH
natural rubber latex

- Scratch-resistant, 1-inch aluminum tubing provides maximum strength while remaining lightweight
- Rugged steel cross-bars add strength without interfering with patient stride
- Dual-release folding mechanism allows frame to be folded easily to a depth of 4 inches and is ideal for users with limited dexterity
- Contoured, handgrips provide additional comfort and security
- Available with or without wheels
- 5-inch wheels featured on WLK4500 and WLK3205 are ideal for navigating outdoor terrain
- 3-inch wheels featured on 77743 provide a low profile for use over flat surfaces, especially indoors
- 300-pound weight capacity



Item No.	Description	Quantity
Walker Folding Dual-Release 300-Pound Weight Capacity		
77742	1-inch Aluminum 22.64 Inch Wide × 18.7 Inch Deep Height Adjustable from 32 to 39 Inch Contoured PVC Handgrips	1/EA, 4 EA/CA
77743	With 3-inch Front Wheels 1-inch Aluminum 23 Inch Wide × 18.7 Inch Deep Height Adjustable from 32 to 39 Inch Contoured PVC Handgrips	1/EA, 4 EA/CA
WLK4500	Adult with 5-inch Front Wheels 25.39 Inch Wide × 18.7 Inch Deep Height Adjustable from 32 to 39 Inch	1/EA, 4 EA/CA
WLK3200	Youth 1-inch Aluminum 22 Inch Wide × 18.5 Inch Deep Height Adjustable from 25.5 to 32.5 Inch	1/EA, 4 EA/CA
WLK3205	Youth with 5-inch Front Wheels 1-inch Aluminum 22 Inch Wide × 18.5 Inch Deep Height Adjustable from 25.5 to 32.5 Inch	1/EA, 4 EA/CA

Accessories

Item No.	Description
MAA2001	Tips Walker Pair Fits Walkers 77741, 77742, 77743, and WLK4500
MAA2002	Wheels Walker 5-inch Pair Fits Walkers 77741, 77742, 77743, and WLK4500
MAA2003	Glide Caps Walker Pair Fits Walkers 77741, 77742, 77743, and WLK4500
MAA2004	Ski Glides Walker Pair Fits Walkers 77741, 77742, 77743, and WLK4500

WALKERS

Bariatric Dual-Release Folding Walkers

NOT MADE WITH
natural rubber latex

- Durable, 1-inch anodized aluminum tubing provides maximum strength while remaining lightweight
- Specialized design provides more width and depth for bariatric patients
- Raised H-frame facilitates over-toilet positioning
- Dual-release folding mechanism allows frame to be folded easily and is ideal for users with limited dexterity
- Contoured, textured handgrips provide additional comfort and security
- Height adjustable from 31 inches to 38 inches
- Suitable for patients ranging from 5 feet 4 inches to 6 feet 2 inches
- Available with or without 5-inch wheels
- 600-pound weight capacity
- Limited lifetime warranty



Item No.	Description	Quantity
Walker Bariatric Folding Dual-Release 1-inch Anodized Aluminum		
BAR806	26.125 Inch Wide × 18 Inch Deep Height Adjustable from 31 to 38 Inch	1/EA, 2 EA/CA
BAR806W	With 5-inch Front. 26.125 Inch Wide × 18 Inch Deep Height Adjustable from 31 to 38 Inch	1/EA, 2 EA/CA

Lightweight Rigid Walkers

NOT MADE WITH
natural rubber latex

- Non-folding design provides extra rigidity for use where space is not an issue
- Scratch-resistant, 1-inch aluminum tubing provides maximum strength while remaining lightweight
- Angled legs provide maximum patient support in all directions and help resist tipping
- Angled side brace helps support patients as they stand up and sit down
- Fits easily around standard toilet fixtures
- Height adjustable from 31 inches to 38 inches
- 400-pound weight capacity
- Limited lifetime warranty



Item No.	Description	Quantity
WLK1000	Walker Rigid Lightweight 1-inch Aluminum 22.5 Inch Wide × 20.5 Inch Deep Height Adjustable from 31 to 38 Inch	1/EA, 4 EA/CA

Mobility Aid Accessories

NOT MADE WITH
natural rubber latex

- Complementary accessories for MediChoice Aluminum Crutches and selected MediChoice Walkers
- Crutch accessories fit MediChoice Aluminum Crutches 77708A, 77709A, 77710A, and 77711A
- Walker accessories fit MediChoice Walkers 77741, 77742, 77743, and WLK4500
- Walker glide caps and ski glides make movement with a wheeled walker easier and smoother for patients
- All accessories provided in pairs



Item No.	Description	Quantity
MAA1001	Cushion Underarm Crutch	1/PR, 18 PR/CA
MAA1002	Handgrip Crutch	1/PR, 18 PR/CA
MAA1003	Tip Crutch	1/PR, 18 PR/CA
MAA2001	Tips Walker	1/PR, 20 PR/CA
MAA2002	Wheels Walker 5 Inch	1/PR, 6 PR/CA
MAA2003	Glide Caps Walker	1/PR, 20 PR/CA
MAA2004	Ski Glides Walker	1/PR, 20 PR/CA

CANES

Aluminum Canes

NOT MADE WITH
natural rubber latex

- Made of lightweight, high-strength aluminum with a scratch-resistant finish
- Easy-to-use, one-button height adjustment in 1-inch increments— from 30 to 39 Inches
- Locking ring safely secures height-adjusted cane and prevents rattling
- Contoured PVC handle provides a comfortable grip
- Slip-resistant tip enhances safety and stability
- 250-pound capacity in an upright position



Item No.	Description	Quantity
77740	Cane Aluminum 0.875-inch Diameter Contoured PVC Handle Height Adjustable in 1-inch Increments from 30 to 39 Inch Slip-resistant Tip 250-Pound Capacity Pound Capacity at Upright Position	1/EA, 6 EA/CA

Adjustable Offset Canes

NOT MADE WITH
natural rubber latex

- Offset handle places the user's weight directly over the base for maximum balance and control
- Attractive finish resists marring and scratching
- Easy push-button height adjustment with locking collet for added security
- Soft nitrile grip contoured for maximum comfort
- Available in black or chrome finish
- Height adjusts in 1-inch increments from 31 to 39 Inches
- Convenient wrist strap
- Sculpted grip designed to reduce the chance of slippage even when wet
- 250-pound capacity in an upright position
- Limited lifetime manufacturer's warranty



Item No.	Description	Quantity
Cane Aluminum. Offset Contoured Nitrile Handle Height Adjustable in 1-inch Increments from 31 to 39 Inch Slip-resistant Tip		
CANE942	Black Finish	1/EA, 6 EA/CA
CANE943	Chrome Finish	1/EA, 6 EA/CA

CANES

Aluminum Quad Canes

NOT MADE WITH
natural rubber latex

- Offset handle places the user's weight directly over the base for maximum balance and control
- Low profile, four-point base further enhances stability
- Easy-to-use, one-button height adjustment in 1-inch increments—from 30 to 39 Inches
- Locking ring safely secures height-adjusted cane and prevents rattling
- Tool-free adjustment for left- or right-handed use
- Attractive, scratch-resistant chrome finish
- 300-pound capacity in an upright position
- Limited lifetime manufacturer's warranty



Item No.	Description	Quantity
CANE618	Cane Quad Aluminum with Large 12 x 8 Inch Base Offset Handle Adjusts for Left or Right Hand 300-Pound Capacity at Upright Position	1/EA, 4 EA/CA

Solid Maple Wooden Canes

NOT MADE WITH
natural rubber latex

- Made of strong, maple wood with an attractive walnut finish
- 36-inch, non-adjustable length
- 1-inch diameter CANE501 has a 260-pound capacity in an upright position
- 0.875-inch diameter CANE578 has a 240-pound capacity in an upright position
- Contoured handle provides a comfortable grip
- Slip-resistant tip enhances safety and stability
- 1-year limited manufacturer's warranty



Item No.	Description	Quantity
Cane Solid Maple Wooden 36 Inch Long Contoured Handle Slip-resistant Tip		
CANE501	Walnut Finish 1-inch Diameter	1/EA, 6 EA/CA
CANE578	Walnut Finish 0.875-inch Diameter	1/EA, 6 EA/CA

WHEELCHAIRS

Standard Wheelchairs

NOT MADE WITH
natural rubber latex

- Chrome-plated carbon steel frame built for durability
- Dual axle allows seat to be lowered for hemiplegic patients
- Simply lift up on the seat to fold the wheelchair for convenient storage
- Composite mag wheels are molded with contours to enhance strength and durability
- Available with swingaway footrests or elevating legrests with swingaway capability
- Available with removable desk-length or full-length armrests or with fixed, full-length armrests
- Flame-resistant vinyl upholstery available in black, blue, or maroon
- Full-size chart pocket on backrest
- 250-pound weight capacity (300-pound capacity for 20-inch seat width)
- 3-year manufacturer's warranty on frame on crossbrace



Item No.	Description	Quantity
Wheelchair Standard Flame-resistant Vinyl Upholstery Chrome-plated Carbon Steel Frame Adjustable Height for Foot Propulsion		
WC2001	18 inch Wide x 16 inch Deep Seat Black Fixed Footrest Fixed Full Padded Arm	1/EA, 1/CA
WC2001BL	18 inch Wide x 16 inch Deep Seat Blue Fixed Footrest Fixed Full Padded Arm	1/EA
WC2001M	18 inch Wide x 16 inch Deep Seat Maroon Fixed Footrest Fixed Full Padded Arm	1/EA
WC2002	18 inch Wide x 16 inch Deep Seat Black Swingaway Footrest Fixed Full Padded Arm	1/EA, 1/CA
WC2002BL	18 inch Wide x 16 inch Deep Seat Blue Swingaway Footrest Fixed Full Padded Arm	1/EA
WC2002M	18 inch Wide x 16 inch Deep Seat Maroon Swingaway Footrest Fixed Full Padded Arm	1/EA
WC2003	18 inch Wide x 16 inch Deep Seat Black Swingaway Footrest Removable Desk Padded Arm	1/EA, 1/CA
WC2003BL	18 inch Wide x 16 inch Deep Seat Blue Swingaway Footrest Removable Desk Padded Arm	1/EA
WC2003M	18 inch Wide x 16 inch Deep Seat Maroon Swingaway Footrest Removable Desk Padded Arm	1/EA
WC2004	18 inch Wide x 16 inch Deep Seat Black Elevating Legrest Removable Desk Padded Arm	1/EA, 1/CA
WC2004BL	18 inch Wide x 16 inch Deep Seat Blue Elevating Legrest Removable Desk Padded Arm	1/EA
WC2004M	18 inch Wide x 16 inch Deep Seat Maroon Elevating Legrest Removable Desk Padded Arm	1/EA
220RADPS	20 inch Wide x 16 inch Deep Seat Black Swingaway Footrest Removable Desk Padded Arm	1/EA

Accessories

Item No.	Description
WCA02	Holder Oxygen Tank Wheelchair Fits All MediChoice Standard Heavy-Duty and Lightweight Wheelchairs
WCAIVO	Holder Combination Oxygen Tank and IV Fits All MediChoice Standard Wheelchairs
WCAIVP	Pole IV Telescoping Fits All MediChoice Standard Heavy-Duty and Lightweight Wheelchairs
WCAATH	Anti-Theft Device Anti-Fold Fits All MediChoice Standard Heavy-Duty Lightweight and Reclining Wheelchairs
WCAOATH	Anti-Theft Device Overhead Fits All MediChoice Standard Heavy-Duty and Lightweight Wheelchairs
GEYATL	Anti-Tipper Rear Wheelchair Adjustable Position Fits All MediChoice Standard Heavy-Duty and Lightweight Wheelchairs
WCAAT	Anti-Tipper Rear Wheelchair Fixed Position Fits All MediChoice Standard Heavy-Duty and Lightweight Wheelchairs

WHEELCHAIRS

Heavy-Duty Wheelchairs

NOT MADE WITH
natural rubber latex

- Heavy-duty reinforced chromed steel frame supports weights up to 450 pounds (22-inch and 24-inch seat) and 300 pounds (20-inch seat)
- Simply lift up on the seat to fold the wheelchair for convenient storage
- Composite mag wheels are molded with contours to enhance strength and durability
- Durable, double-embossed, flame-resistant vinyl upholstery available in black, blue, or maroon
- Oversized chart pocket on backrest
- Available with aluminum swingaway footrests or elevating legrests with padded calf and swingaway capability
- Available with removable desk-length or full-length armrests
- 3-year manufacturer's warranty on frame and crossbrace



Item No.	Description	Quantity
Wheelchair Heavy Duty Flame-resistant Double-embossed Vinyl Upholstery Reinforced Chrome Steel Frame		
WC2005	20 Inch Wide x 16 Inch Deep Seat Black Swingaway Footrest Removable Full Padded Arm	1/EA, 1/CA
WC2005BL	20 Inch Wide x 16 Inch Deep Seat Blue Swingaway Footrest Removable Full Padded Arm	1/EA
WC2005M	20 Inch Wide x 16 Inch Deep Seat Maroon Swingaway Footrest Removable Full Padded Arm	1/EA
WC2006	20 Inch Wide x 16 Inch Deep Seat Black Elevating Legrest Removable Full Padded Arm	1/EA, 1/CA
WC2006BL	20 Inch Wide x 16 Inch Deep Seat Blue Elevating Legrest Removable Full Padded Arm	1/EA
WC2006M	20 Inch Wide x 16 Inch Deep Seat Maroon Elevating Legrest Removable Full Padded Arm	1/EA
WC2007	22 Inch Wide x 18 Inch Deep Seat Black Swingaway Footrest Removable Desk Padded Arm	1/EA, 1/CA
WC2007BL	22 Inch Wide x 18 Inch Deep Seat Blue Swingaway Footrest Removable Desk Padded Arm	1/EA
WC2007M	22 Inch Wide x 18 Inch Deep Seat Maroon Swingaway Footrest Removable Desk Padded Arm	1/EA
WC2008	22 Inch Wide x 18 Inch Deep Seat Black Elevating Legrest with Upholstered Calf Pads Removable Desk Padded Arm	1/EA, 1/CA
WC2008BL	22 Inch Wide x 18 Inch Deep Seat Blue Elevating Legrest with Upholstered Calf Pads Removable Desk Padded Arm	1/EA
WC2008M	22 Inch Wide x 18 Inch Deep Seat Maroon Elevating Legrest with Upholstered Calf Pads Removable Desk Padded Arm	1/EA
WC2009	24 Inch Wide x 18 Inch Deep Seat Black Swingaway Footrest Removable Desk Padded Arm	1/EA, 1/CA
WC2009BL	24 Inch Wide x 18 Inch Deep Seat Blue Swingaway Footrest Removable Desk Padded Arm	1/EA
WC2009M	24 Inch Wide x 18 Inch Deep Seat Maroon Swingaway Footrest Removable Desk Padded Arm	1/EA
WC2010	24 Inch Wide x 18 Inch Deep Seat Black Elevating Legrest with Upholstered Calf Pads Removable Desk Padded Arm	1/EA, 1/CA
WC2010BL	24 Inch Wide x 18 Inch Deep Seat Blue Elevating Legrest with Upholstered Calf Pads Removable Desk Padded Arm	1/EA
WC2010M	24 Inch Wide x 18 Inch Deep Seat Maroon Elevating Legrest with Upholstered Calf Pads Removable Desk Padded Arm	1/EA
WC2110	24 Inch Wide x 18 Inch Deep Seat Black Elevating Legrest with Upholstered Calf Pads Removable Full Padded Arm	1/EA, 1/CA
WC2110BL	24 Inch Wide x 18 Inch Deep Seat Blue Elevating Legrest with Upholstered Calf Pads Removable Full Padded Arm	1/EA
WC2110M	24 Inch Wide x 18 Inch Deep Seat Maroon Elevating Legrest with Upholstered Calf Pads Removable Full Padded Arm	1/EA

WHEELCHAIRS

Bariatric Folding Wheelchairs

NOT MADE WITH
natural rubber latex

- Reinforced welds and gussets on frame
- Black, flame-resistant, reinforced nylon upholstery
- Simply lift up on the seat to fold the wheelchair for convenient storage
- Detachable desk armrests can be reversed for easier propelling
- Available with swingaway footrests or elevating legrests
- Available with 26 inch or 30 inch wide seats
- 650-pound weight capacity
- 8 inch × 2 inch solid front casters
- 24-inch solid rear wheels
- Limited lifetime manufacturer's warranty on frames and crossbraces



Item No.	Description	Quantity
Wheelchair Bariatric Folding Removable Desk Arm Black Reinforced Nylon Upholstery		
WCBAR26E	Elevating Legrest 26 Inch Wide × 19.5 Deep Seat 35.5 Inch Overall Width	1/EA
WCBAR26S	Swing Away Footrest 26 Inch Wide × 19.5 Deep Seat 35.5 Inch Overall Width	1/EA
WCBAR30E	Elevating Legrest 30 Inch Wide × 19.5 Deep Seat 40 Inch Overall Width	1/EA
WCBAR30S	Swing Away Footrest 30 Inch Wide × 19.5 Deep Seat 40 Inch Overall Width	1/EA

Transport Wheelchairs

- Compact and lightweight, this economical wheelchair makes caregiver-assisted transport safe, convenient, and comfortable
- Folding frame and back allow the wheelchair to become extremely compact for storage in small places
- Chrome-plated carbon steel frame built for durability
- Strong, 8-inch wheels in the front and back provide maneuverability and a smooth ride
- Comfortable, durable, and easy-care black padded nylon upholstery
- Wheels contain natural rubber latex
- 250-pound weight capacity
- 3-year manufacturer's warranty on frame and crossbrace



Item No.	Description	Quantity
WC1001	Wheelchair Transport Folding Frame and Back 19 Inch Wide × 16 Inch Deep Seat 8-inch Wheels Front and Back Black Padded Nylon Upholstery Automotive Style Seatbelt Swingaway Footrest Fixed Armrest with Full Arm Pads Chrome-plated Carbon Steel Frame	1/EA, 1/CA

WHEELCHAIRS

Reclining Wheelchairs

NOT MADE WITH
natural rubber latex

- Reclining backrest offers patient comfort and positioning
- Full back and headrest support
- Chrome-plated carbon steel frame built for durability
- Reclines from 90 to 160 degrees
- Rear anti-tippers standard
- Dual axle allows seat to be lowered for hemiplegic patients
- Removable 10-inch headrest
- Simply lift up on the seat to fold the wheelchair for convenient storage
- Padded nylon upholstery is comfortable and durable
- Elevating legrest models have height-adjustable, full padded armrests
- Swingaway footrest models have removable desk armrests
- 250-pound weight capacity
- 3-year manufacturer's warranty on frame on crossbrace



Item No.	Description	Quantity
Wheelchair Reclining Full Back and Headrest Support Removable 10-inch Headrest Reclines 90 to 160 Degrees Black Padded Nylon Upholstery Chrome-plated Carbon Steel Frame		
4R18RADPS	18 Inch × 17 Inch Seat Swingaway Footrest Removable Desk Padded Arm	1/EA
4R18HAFPE	18 Inch × 17 Inch Seat Elevating Legrest Removable Full Padded Arm	1/EA
4R20HAFPE	20 Inch × 17 Inch Seat Elevating Legrest Height Adjustable Full Padded Arm	1/EA

Lightweight Wheelchairs

NOT MADE WITH
natural rubber latex

- Weighs only 35 pounds
- Chrome-plated carbon steel frame built for durability
- Scratch-resistant composite hand rim
- Dual axle allows seat to be lowered for hemiplegic patients
- Simply lift up on the seat to fold the wheelchair for convenient storage
- Comfortable padded nylon upholstery
- 250-pound weight capacity
- 3-year manufacturer's warranty on frame and crossbrace



Item No.	Description	Quantity
318RADPE	Wheelchair Lightweight 18 Inch Wide × 16 Inch Deep Seat Black Padded Nylon Upholstery Elevating Legrest Removable Desk Padded Arm Chrome-plated Carbon Steel Frame	1/EA

Wheelchair Accessories

NOT MADE WITH
natural rubber latex

- Complementary products to be used with MediChoice Wheelchairs
- Some replacement parts also included
- Additional accessories/parts can be added as needed if available



Item No.	Description	Quantity
Oxygen Tank Holders and IV Poles <i>Fits All MediChoice Standard, Heavy-Duty, and Lightweight Wheelchairs</i>		
WCAO2	Oxygen Tank Holder	1/EA
WCAIVP	Telescoping IV Pole	1/EA
<i>Fits MediChoice Standard Wheelchairs Only</i>		
WCAIVO	Combination Oxygen Tank and IV Pole	1/EA
Rear Anti-Tippers <i>Fits All MediChoice Standard, Heavy-Duty and Lightweight Wheelchairs</i>		
GEYATL	Adjustable Position	1/EA
WCAAT	Fixed Position	1/PR
Anti-Theft/Anti-Fold Devices <i>Fits All MediChoice Standard, Heavy-Duty, Lightweight, and Reclining Wheelchairs</i>		
WCAATH	Anti-Theft/Anti-Fold	1/EA
<i>Fits All MediChoice Standard, Heavy-Duty and Lightweight Wheelchairs</i>		
WCAOATH	Overhead Anti-Theft Device	1/EA
Replacement Upholstery <i>Black Vinyl Seat and Back—Includes All Necessary Hardware</i>		
GEY2ORRUP	For MediChoice 20-inch Wheelchairs—220RADPS WC2005 WC200620	1/EA
GEY218UP	For MediChoice 18-inch Wheelchairs—WC2001 WC2002 WC2003 WC200418	1/EA
GEY22HDUP	For MediChoice 22-inch Wheelchairs—WC2007 WC200822	1/EA
GEY24HDUP	For MediChoice 24-inch Wheelchairs—WC2009 WC2010 WC211024	1/EA
Swingaway Footrest Assemblies <i>Fits All MediChoice Heavy-Duty Wheelchairs</i>		
GEYFRA	Footrest Assembly with Aluminum Footplate	1/PR
<i>Fits All MediChoice Standard, Lightweight, and Reclining Wheelchairs</i>		
GEYFRC	Footrest Assembly with Composite Footplate	1/KT
WCAFRCHL	Footrest Assembly with Composite Footplate and Heel Loop	1/KT
Elevating Legrest Assemblies with Black Vinyl Calf Pads <i>Fits All MediChoice Heavy-Duty Wheelchairs</i>		
GEYELRA	Elevating Legrest Assembly with Aluminum Footplate	1/PR
<i>Fits All MediChoice Standard, Lightweight, and Reclining Wheelchairs</i>		
WCAELRC	Elevating Legrest Assembly with Composite Footplate	1/PR
<i>Fits All MediChoice Standard and Reclining Wheelchairs with 2-inch Seat Extender</i>		
WCAELRLP	Elevating Legrest Assembly with Composite Footplate	1/PR
Articulating Legrest Assemblies with Black Vinyl Calf Pads <i>Fits All MediChoice Standard, Lightweight, and Reclining Wheelchairs</i>		
GEYALR	Articulating Legrest Assembly with Composite Footplate	1/EA

Item No.	Description	Quantity
Armrest Assemblies with Black Vinyl Pads <i>Fits All MediChoice Standard, Heavy-Duty, Lightweight, and Reclining Wheelchairs</i>		
WCAHADP	Desk Arm Adjustable Height	1/PR
WCARADP	Desk Arm Fixed Height	1/PR
WCAHAFF	Full Arm Adjustable Height	1/PR
WCARAFP	Full Arm Fixed Height	1/PR
Replacement Armpads (Armpads Only) <i>Fits All MediChoice Transport Wheelchairs</i>		
ECTC-ARMPAD	Full Black Nylon	1/PR
<i>Fits All MediChoice Standard, Heavy-Duty, Lightweight, and Reclining Wheelchairs</i>		
GEY2AP	Full Black Vinyl	1/EA
GEY2APD	Desk Black Vinyl	1/EA
Wheelchair Cushions, Foam Base		
300	Deep Contoured with Fluid Tripad Stretch Cover 16 × 16 × 4.5 Inch, 250-Pound Capacity	1/EA
305MJ	Mild Contoured with Moisture-Resistant Cover and No-Slip Bottom 16 × 16 × 2 Inch, 250-Pound Capacity	1/EA
2466	Mild Contoured with Moisture-Resistant Cover and No-Slip Bottom 18 × 16 × 2 Inch, 250-Pound Capacity	1/EA
JBASICP1618	High-density Contoured with Water-Resistant, Dartex-Coated Top Cover 16 × 18 × 3 Inch, 300-Pound Capacity	1/EA
JBASICP1816	High-density Contoured with Water-Resistant, Dartex-Coated Top Cover 18 × 16 × 3 Inch, 300-Pound Capacity	1/EA
JBASICP2016	High-density Contoured with Water-Resistant, Dartex-Coated Top Cover 20 × 16 × 3 Inch, 300-Pound Capacity	1/EA
JBASICP2018	High-density Contoured with Water-Resistant, Dartex-Coated Top Cover 20 × 18 × 3 Inch, 300-Pound Capacity	1/EA
JBASICP2020	High-density Contoured with Water-Resistant, Dartex-Coated Top Cover 20 × 20 × 3 Inch, 300-Pound Capacity	1/EA
JBASICP2218	High-density Contoured with Water-Resistant, Dartex-Coated Top Cover 22 × 18 × 3 Inch, 500-Pound Capacity	1/EA
Front Caster Assemblies <i>Including 8 Inch Casters, Fork Assembly, and all Necessary Hardware</i>		
ECTCCASTER	MediChoice Transport Wheelchairs	1/PR
GEY2CA	Black-Poly, Scratch-Resistant Casters for MediChoice 18-20 Inch Standard, Heavy-Duty, Lightweight, and Reclining Wheelchairs	1/PR
GEY2HDCA	Pneumatic Casters for MediChoice 22-24 Inch Heavy-Duty Wheelchairs	1/PR
Rear Wheel Assemblies <i>Including 24 × 1 Inch Rear Rim and Tire</i>		
GEY1WA	Composite Rim with Low-Poly Tires for MediChoice Standard and Lightweight Wheelchairs	1/EA
GEY2HDWA	Chrome Rim with Full-Poly Tires for MediChoice Heavy-Duty Wheelchairs	1/EA
Wheel Lock Assemblies <i>For MediChoice Standard, Heavy-Duty, and Lightweight Wheelchairs</i>		
GEY1WL	Fixed Armrest	1/EA
GEY2WL	Push-To-Lock Style for Wheelchairs without Elevating Legrests	1/PR
GEY2WL1	Push-To-Lock Style for Wheelchairs with Elevating Legrests	1/PR

Please see our full line of Wheelchair Scales on page 39 in the Equipment Section.

Please see our full line of Wheelchair Cushions on page 168 in the Skin Care Section.

Yankauers

NOT MADE WITH natural rubber latex **sterile EO**

- PVC suction instrument
- Transparent plastic allows visual check for fluids
- Available with regular open tip or bulb tip with or without a control vent
- Curve in the Yankauer allows for ease of site use and minimizes hand and arm activity
- Non-conductive



Item No.	Description	Quantity
Yankauer Suction Bulb Tip		
YANK1001	Non-Vented Polyvinyl Chloride Blue Sterile	1/EA, 50 EA/CA
YANK1002	Vented Polyvinyl Chloride Blue Sterile	1/EA, 50 EA/CA
YANK3001	Non-Vented Rigid Polyvinyl Chloride Purple Sterile	1/EA, 50 EA/CA
Yankauers Suction Open Tip		
YANK2001	Non-Vented Polyvinyl Chloride Blue Sterile	1/EA, 50 EA/CA
YANK2002	Vented Polyvinyl Chloride Blue Sterile	1/EA, 50 EA/CA

Decanting Devices

NOT MADE WITH natural rubber latex **sterile EO**

- Used for aseptic decanting of solutions from containers
- Vial decanter has air inlets to ensure smooth pouring
- Made of durable ABS plastic




Item No.	Description	Quantity
BD10102	Decanter Bag for Aseptic Transfer of Fluids from Flexible Containers Acrylonitrile-Butadiene-Styrene Plastic Sterile	1/EA, 50 EA/CA
VD10106	Decanter Vial for Aseptic Transfer of Fluids from Small Vials Acrylonitrile-Butadiene-Styrene Plastic Sterile	1/EA, 50 EA/CA
TD10108	Device Transfer for Aseptic Transfer of Fluids from One Container to Another Sterile	1/EA, 50 EA/CA

Electrosurgical Tip Cleaners

made in the USA NOT MADE WITH natural rubber latex **sterile R**

- Foam pad measures 1¾ Inch × 1¾ Inch with abrasive material on one side, strong adhesive on the other
- Cleans both mono-polar and bi-polar probes
- Adheres to both disposable and non-disposable drapes
- Radiopaque strip in adhesive side for x-ray detection
- Individually packaged



Item No.	Description	Quantity
TIP1001	Cleaner Tip Electrosurgical Pencil 1.75 Inch × 1.75 Inch Sterile	1/EA, 100 EA/CA

Suction Tubing

NOT MADE WITH natural rubber latex **sterile EO**

- Durable and flexible tubing withstands heavy-duty suctioning
- Available in 6', 10', 12' and 20' lengths
- Soft female color-coded connector on each end
- Male adapter included
- Non-conductive, sterile PVC
- Suction power is 22 Inch of mercury



Item No.	Description	Quantity
Tubing Suction 2 Connector Non-Conductive Heavy-Duty Durable Flexible with Male Adapter Polyvinyl Chloride Sterile		
ST1001	6 Feet × ⅜ Inch	1/EA, 50 EA/CA
ST1002	10 Feet × ⅜ Inch	1/EA, 50 EA/CA
ST1003	12 Feet × ⅜ Inch	1/EA, 20 EA/CA
ST1004	20 Feet × ⅜ Inch	1/EA, 20 EA/CA
ST2001	6 Feet × ¼ Inch	1/EA, 50 EA/CA
ST2002	10 Feet × ¼ Inch	1/EA, 50 EA/CA
ST2003	12 Feet × ¼ Inch	1/EA, 20 EA/CA
ST2004	20 Feet × ¼ Inch	1/EA, 20 EA/CA
ST3001	6 Feet × ⅝ Inch	1/EA, 50 EA/CA
ST3002	10 Feet × ⅝ Inch	1/EA, 50 EA/CA
ST3003	12 Feet × ⅝ Inch	1/EA, 20 EA/CA
ST3004	20 Feet × ⅝ Inch	1/EA, 20 EA/CA

Solidifiers

made in the USA NOT MADE WITH natural rubber latex

- Completely and quickly solidifies liquid medical waste, fluids and saline solutions
- Solidified liquid medical waste cannot spill or splash which reduces the possibility of employee exposure to potentially infectious waste
- Solidified liquid medical waste is in compliance with DOT regulations
- Contains no potentially harmful chemicals such as Chlorine, Phenol, Gluteraldehyde or Iodine
- Available in convenient unit dose bottles and in a bulk bottle
- Spout caps provide easy access to smaller ports on suction canister lids



Item No.	Description	Quantity
Solidifier Medical Waste		
SLD500	500 Cubic Centimeter Screw Cap	1/EA, 180 EA/CA
SLD1200	1200 Cubic Centimeter Screw Cap	1/EA, 180 EA/CA
SLD1500	1500 Cubic Centimeter Screw Cap	1/EA, 180 EA/CA
SLD2000	2000 Cubic Centimeter Screw Cap	1/EA, 130 EA/CA
SLD2500	2500 Cubic Centimeter Screw Cap	1/EA, 130 EA/CA
SLD3000	3000 Cubic Centimeter Screw Cap	1/EA, 90 EA/CA
SLD14000	14000 Cubic Centimeter Spout Cap	1/EA, 24 EA/CA
SLD21200	1200 Cubic Centimeter Spout Cap	1/EA, 154 EA/CA
SLD21500	1500 Cubic Centimeter Spout Cap	1/EA, 154 EA/CA
SLD22000	2000 Cubic Centimeter Spout Cap	1/EA, 108 EA/CA
SLD23000	3000 Cubic Centimeter Spout Cap	1/EA, 108 EA/CA

Insufflation Tubing

NOT MADE WITH natural rubber latex **sterile EO**

- High-flow 10' PVC tubing
- Available Luer-to-Luer or Luer-to-Quick-Connect
- 1 micron hydrophobic ULPA filter
- DEHP-free
- Pliable tubing/no memory
- Most are packaged sterile and in an individual peel pouch
- Item 1314INSF02NS is non-sterile and packaged in bulk
- Storz and Stryker compatible heated insufflation tubing is also available



Item No.	Description	Quantity
Tubing Insufflation 0.1 Micron Filter DEHP-Free Polyvinyl Chloride		
INSF02NS	Quick-Connect 10 Feet Non-Sterile	1/EA, 50 EA/CA
INSF1001	Luer 10 Feet Sterile	1/EA, 60 EA/CA
INSF1002	Quick-Connect 10 Feet Sterile	1/EA, 60 EA/CA
INSF1003	Luer Funnel Connect Stryker 10 Feet Sterile	1/EA, 60 EA/CA
INSF1004	Luer Round Connect Storz 10 Feet Sterile	1/EA, 60 EA/CA
INSF2004	Luer Round Connect Heated Storz 10 Feet Sterile	1/EA, 60 EA/CA
INSF2003	Luer Heated Stryker 10 Feet Sterile	1/EA, 5 EA/BX, 4 EA/CA

Biopsy Punches

NOT MADE WITH natural rubber latex **diverse supplier** **sterile EO**

- Seamless stainless steel blade with protective cap on punch blade
- Machined blade minimizes tissue tears
- High-impact polystyrene hollow handle for easy sample collection
- Punch size is printed on handle
- Package is color coded by size for easy identification
- Single use



Item No.	Description	Quantity
Punch Biopsy Dermal Stainless Steel Polystyrene Sterile		
DP0200	2 Millimeter	1/EA, 25 EA/BX
DP0300	3 Millimeter	1/EA, 25 EA/BX
DP0350	3.5 Millimeter	1/EA, 25 EA/BX
DP0400	4 Millimeter	1/EA, 25 EA/BX
DP0500	5 Millimeter	1/EA, 25 EA/BX
DP0600	6 Millimeter	1/EA, 25 EA/BX
DP0800	8 Millimeter	1/EA, 25 EA/BX

Disposable Safety Scalpels

NOT MADE WITH natural rubber latex **diverse supplier**

- Disposable Safety Scalpels with Moveable Sheath
 - Highest quality stainless steel Sheffield blade, made in England
 - Stable, fixed stainless steel blade with moveable sheath to prevent accidental blade exposure
 - OSHA-compliant
 - Protective, transparent sheath can be opened, closed and locked with a single hand; audible double-click when sheath is in locked position
 - Sturdy, "easy-grip" acrylo-nitrile handle fits comfortably in either hand
 - Individually packaged
- Disposable Safety Scalpels with Retractable Blade
 - Retractable, Swiss stainless steel blade has a razor sharp edge for superior cutting and clean incisions
 - Ribbed, ambidextrous handle is top activated providing a well-balanced grip
 - Non-slip ergonomic handle with metric scale has a one-piece design for improved integrity
 - Enclosed front shield enhances blade stability and security
 - Audible click mechanism indicates blade and final lock position for safe handling and disposal
 - Fully extended blade and transparent blade holder enables clear line of sight and full incision depth
 - Color coded slider button on handle for easy blade size identification



Item No.	Description	Quantity
Scalpel Surgical Safety Stainless Steel Disposable Sterile		
55610	Fixed Blade with Movable Sheath Size 10	1/EA, 10 EA/BX
55611	Fixed Blade with Movable Sheath Size 11	1/EA, 10 EA/BX
55615	Fixed Blade with Movable Sheath Size 15	1/EA, 10 EA/BX
SCPL10	Retractable Blade Size 10	1/EA, 10 EA/BX
SCPL11	Retractable Blade Size 11	1/EA, 10 EA/BX
SCPL15	Retractable Blade Size 15	1/EA, 10 EA/BX

Eye Protectors

NOT MADE WITH natural rubber latex **sterile EO**



- Protects the eyes from blood, fluid, and foreign materials or unintentional contact
- Offers additional protection for the patient when under general anesthesia
- Rigid clear plastic protective eye covers
- Hypo-allergenic self adhesive foam cushion
- Easy peel backing with courtesy tabs for fast accurate positioning over the eyes
- Single Patient Use

Item No.	Description	Quantity
Protector Anesthesia Eye Foam Single-Patient Use Plastic Clear Sterile		
EP1500	Child	1/EA, 25 EA/CA
EP2500	Adult	1/EA, 25 EA/CA
EP3500	Double Foam Adult	1/EA, 25 EA/CA



Esmark Bandages

diverse supplier **sterile R**

- Smooth, comfortable material that ensures patient protection
- Consistent elongation for uniform compression and easy application
- Safely provides maximum pressure to extremities for exsanguination before surgery
- Conveniently packaged in sterile pouch

Item No.	Description	Quantity
Bandage Compression Esmark Sterile		
55501	12 Feet x 3 Inch Contains Natural Rubber Latex	1/EA, 25 EA/CA
55502	9 Feet x 3 Inch Contains Natural Rubber Latex	1/EA, 20 EA/CA
55503	12 Feet x 4 Inch Contains Natural Rubber Latex	1/EA, 20 EA/CA
55504	9 Feet x 4 Inch Contains Natural Rubber Latex	1/EA, 20 EA/CA
55505	12 Feet x 6 Inch Contains Natural Rubber Latex	1/EA, 20 EA/CA
55506	9 Feet x 6 Inch Contains Natural Rubber Latex	1/EA, 20 EA/CA
55511	Synthetic 12 Feet x 3 Inch	1/EA, 20 EA/CA
55512	Synthetic 9 Feet x 3 Inch	1/EA, 20 EA/CA
55513	Synthetic 12 Feet x 4 Inch	1/EA, 20 EA/CA
55514	Synthetic 9 Feet x 4 Inch	1/EA, 20 EA/CA
55515	Synthetic 12 Feet x 6 Inch	1/EA, 20 EA/CA
55516	Synthetic 9 Feet x 6 Inch	1/EA, 20 EA/CA

Kick Bucket Liners

NOT MADE WITH natural rubber latex



- Anti-static liner
- Fits standard 13 quart kick buckets
- Meets Standard NFPA 56A
- Capacity: 8 gallons
- Dimensions: 24 Inch x 30 Inch
- Gauge: 1.00 Mil

Item No.	Description	Quantity
H4830YW	Kick Bucket Liner 24 Inch x 30 Inch	1/EA, 125 EA/BX, 4 BX/CA

Complementary Products

Item No.	Description
BUKT1101	Bucket Kick with Bumper 2 Inch Casters 13 Quart Stainless Steel 13.5 Inch x 14 Inch



Disposable OR Towels

NOT MADE WITH natural rubber latex **sterile R**

- Provide a cost-effective alternative for reducing the risks associated with reusables
- Highly-absorbent, 100% woven cotton, pre-washed and de-linted.
- All towels measure 16 Inch by 24 Inch
- Available in three colors and a variety of packaging configurations
- Most are packaged sterile, pre-folded, poly-wrapped and in convenient peel pouch
- Item #ORT400B, ORT100G and ORT100GX are non-sterile and packaged in bulk

Item No.	Description	Quantity
Towel Operating Room Pre-washed De-linted Cotton		
ORT02B	Blue 2 Pack Sterile	2/PK, 40 PK/CA
ORT02BX	Blue Xray Detectable 2 Pack Sterile	2/PK, 40 PK/CA
ORT02W	White 2 Pack Sterile	2/PK, 40 PK/CA
ORT02WX	White Xray Detectable 2 Pack Sterile	2/PK, 40 PK/CA
ORT04B	Blue 4 Pack Sterile	4/PK, 20 PK/CA
ORT04BX	Blue Xray Detectable 4 Pack Sterile	4/PK, 20 PK/CA
ORT04G	Green 4 Pack Sterile	4/PK, 20 PK/CA
ORT04W	White 4 Pack Sterile	4/PK, 20 PK/CA
ORT04WX	White Xray Detectable 4 Pack Sterile	4/PK, 20 PK/CA
ORT06B	Blue 6 Pack Sterile	6/PK, 12 PK/CA
ORT06BX	Blue Xray Detectable 6 Pack Sterile	6/PK, 12 PK/CA
ORT06GX	Green Xray Detectable 6 Pack Sterile	6/PK, 12 PK/CA
ORT08B	Blue 8 Pack Sterile	8/PK, 10 PK/CA
ORT100G	Green Bulk Non-Sterile	1/EA, 100 EA/CA
ORT100GX	Green Xray Detectable Bulk Non-Sterile	1/EA, 100 EA/CA
ORT10B	Blue 10 Pack Sterile	10/PK, 8 PK/CA
ORT12B	Blue 12 Pack Sterile	12/PK, 6 PK/CA
ORT400B	Blue Bulk Non-Sterile	1/EA, 400 EA/CA

Scope Warmers

NOT MADE WITH natural rubber latex **sterile EO**

- Simple one step activation for instant heat
- Heats scope in approximately one minute. Scope stays warm 45 minutes to an hour
- Eliminates fogging during laparoscopic procedures
- Accommodates all scope sizes and lengths
- Latex and DEHP Free



Item No.	Description	Quantity
DLW510	Disposable Laparoscope Warmer	1/EA, 10 EA/BX, 3 BX/CA

Cordless Cauteries

NOT MADE WITH natural rubber latex **sterile EO**

- Battery-operated, single-use instrument available in fine tip or loop tip
- Fast heat-up and cool-down; Heats to 2200° F
- ABS plastic casing; plastic switch integrated with brass conductor strip to keep it in place
- Guaranteed five-year shelf life (manufacture date on product; expiration date on packaging)
- Operating instructions on product package
- Conveniently packaged in sterile pouches; 10 per shelf pack



Item No.	Description	Quantity
Cautery Cordless Hi Temp Sterile		
055401	Fine Tip	1/EA, 10 EA/BX
055403	Loop Tip	1/EA, 10 EA/BX

Hypo Hyperthermia Blankets

made in the USA NOT MADE WITH natural rubber latex

- Made from soft, thin, light weight polyurethane material for better performance and patient comfort
- Unique pinwheel channel design which promotes excellent water flow and improved temperature therapy
- White sides of pads and blankets can be wetted to provide moist therapy
- Superior draping characteristics help prevent sliding or shifting during use
- Low internal pressures (usually 3-5 psi) on all pads and blankets, resulting in reducing contact point pressure, helping to prevent skin reddening or breakdown
- Smaller pads are also available for heat therapy (or for cooling dependent upon machine capability)



Item No.	Description	Quantity
Blanket Therapeutic Hypo-Hyperthermia Polyurethane		
HTB030	Pediatric 24 x 30 Inch	1/EA, 10 EA/CA
HTB060	Adult 24 x 60 Inch	1/EA, 10 EA/CA
HTBW30	Pediatric 24 x 30 Inch Reusable	1/EA, 10 EA/CA
HTBW60	Adult 24 x 60 Inch Reusable	1/EA, 10 EA/CA

Tubing Connectors

NOT MADE WITH natural rubber latex **sterile EO**

- Sturdy PVC connectors provide tight seal
- Easy to connect to a variety of tubing sizes
- Available in sterile or non-sterile



Item No.	Description	Quantity
Connector Tubing Tight Seal Sturdy Polyvinyl Chloride		
CON1001	Sims Sterile	1/EA, 100 EA/BX, 10 BX/CA
CON1002	Sims Non-Sterile	1/EA, 100 EA/BX, 10 BX/CA
CON2001	5-In-1 Sterile	1/EA, 100 EA/BX, 10 BX/CA
CON2002	5-In-1 Non-Sterile	1/EA, 100 EA/BX, 10 BX/CA

Anti-Fog Solution

made in the USA NOT MADE WITH natural rubber latex **sterile EO**

- Non-abrasive
- Dries in seconds
- Inhibits condensation and clouding for uninterrupted procedures
- Individual sterile packages
- Nontoxic, laboratory-tested, safe to use in body cavities
- Radiopaque, adhesive-backed sponge



Item No.	Description	Quantity
FOG1001	Anti-Fog Solution, Endoscopic, Non-Toxic, Non-Abrasive, 6 Cubic Centimeter with Radiopaque Adhesive-Backed Sponge	1/EA, 48 EA/CA

Vessel Loops

made in the USA NOT MADE WITH natural rubber latex **sterile EO**

- Made from high-grade medical elastomer, stronger than traditional silicone loops
- Excellent atraumatic surface quality, smoother and less abrasive than silicone loops
- 18 Inch length gives you two to three inches more than conventional loops for more convenient fastening
- Radiopaque (visible to X-ray) and color-coded; Red - Arteries, Blue - Veins, White - Nerves, Yellow - Ureters
- Available in maxi and mini sets



Item No.	Description	Quantity
Loop Vessel 18 Inch Elastomer Sterile		
LOOP1001	Maxi Set Blue 2 Per Card	2/PK, 50 PK/BX
LOOP1002	Maxi Set Blue 2 Per Pack	2/PK, 50 PK/BX
LOOP1003	Maxi Set Red 2 Per Card	2/PK, 50 PK/BX
LOOP1004	Maxi Set Red 2 Per Pack	2/PK, 50 PK/BX
LOOP1005	Maxi Set White 2 Per Card	2/PK, 50 PK/BX
LOOP1006	Maxi Set White 2 Per Pack	2/PK, 50 PK/BX
LOOP1007	Mini Set Yellow 2 Per Card	2/PK, 50 PK/BX
LOOP1008	Maxi Set Yellow 2 Per Pack	2/PK, 50 PK/BX
LOOP2001	Mini Set Blue 2 Per Card	2/PK, 50 PK/BX
LOOP2002	Mini Set Blue 2 Per Pack	2/PK, 50 PK/BX
LOOP2003	Mini Set Red 2 Per Card	2/PK, 50 PK/BX
LOOP2004	Mini Set Red 2 Per Pack	2/PK, 50 PK/BX
LOOP2005	Mini Set White 2 Per Card	2/PK, 50 PK/BX
LOOP2006	Mini Set White 2 Per Pack	2/PK, 50 PK/BX
LOOP2007	Mini Set Yellow 2 Per Card	2/PK, 50 PK/BX
LOOP2008	Mini Set Yellow 2 Per Pack	2/PK, 50 PK/BX

Utility Markers



- Non-smearing
- Fast drying
- Will not fade or run
- Ink Color: Black



Item No.	Description	Quantity
Marker Utility Black		
MKR5025	Non-Sterile	1/EA, 25 EA/BX
MKR5109	With 9 Blank Labels Sterile	1/EA, 50 EA/BX

Laparotomy Sponges



- Sterile contents with consistent 5 count packs
- X-ray detectable loops or rings
- Pre-washed, delinted
- 100% cotton, highly-absorbent
- Available in soft packs or rigid trays for delivery to sterile field
- Blue loops & rings provide ease of counting and retraction from the body cavity




Item No.	Description	Quantity
Sponge Laparotomy Pre-Washed X-Ray Detectable Sterile		
LAP1212	12 Inch x 12 Inch 5 Per Pack	50/PK, 40 PK/CA
LAP1818	18 Inch x 18 Inch 5 Per Pack	50/PK, 40 PK/CA
LAP1818R	18 Inch x 18 Inch with Rings 5 Per Pack	50/PK, 40 PK/CA
LAP1818T	18 Inch x 18 Inch Rigid Tray 5 Per Pack	50/PK, 40 PK/CA
LAP418	4 Inch x 18 Inch 5 Per Pack	50/PK, 40 PK/CA
LAP836	8 Inch x 36 Inch 5 Per Pack	50/PK, 40 PK/CA

Electrosurgical Pencils/Tips



- Pencils are available in push button or rocker switch configurations
- Pencils include 10 foot cord with universal 3-prong connector that will fit all competitors' generators
- Universal design allows for interchangeable tip electrodes
- Ergonomic lightweight design reduces the stress and tension on your hands during lengthy procedures
- Packaged in a convenient insulated safety holster



Item No.	Description	Quantity
----------	-------------	----------

Electrode Electrosurgical Blade Polytetrafluoroethylene Coated Sterile

ES25020	2.5 Inch	1/EA, 50 EA/CA
ES25021	Extended Insulation 2.5 Inch	1/EA, 50 EA/CA
ES27520	2.75 Inch	1/EA, 50 EA/CA
ES27521	Extended Insulation 2.75 Inch	1/EA, 50 EA/CA
ES65020	6.5 Inch	1/EA, 50 EA/CA
ES65021	Extended Insulation 6.5 Inch	1/EA, 50 EA/CA
ES65024	Stainless Steel 6.5 Inch	1/EA, 50 EA/CA

Electrode Electrosurgical Needle Polytetrafluoroethylene Coated Sterile

ES27530	2.75 Inch	1/EA, 50 EA/CA
ES27531	Extended Insulation 2.75 Inch	1/EA, 50 EA/CA
ES65031	Extended Insulation 6.5 Inch	1/EA, 50 EA/CA

Pencil Electrosurgical with Holster 10 Foot Cord with Universal 3-Prong Connector Sterile

ES5004	Button Switch Stainless Steel Blade Electrode 2.5 Inch	1/EA, 50 EA/CA
ES5104	Rocker Switch Stainless Steel Blade Electrode 2.5 Inch	1/EA, 50 EA/CA
ES8000	Button Switch Polytetrafluoroethylene Blade Electrode 2.5 Inch	1/EA, 50 EA/CA
ES8100	Rocker Switch Polytetrafluoroethylene Blade Electrode 2.5 Inch	1/EA, 50 EA/CA

Light Handle Covers



- Eliminate labor-intensive processing costs
- Reduce contamination risk
- Increase productivity
- Reduce setup time



Item No.	Description	Quantity
----------	-------------	----------

Cover Light Handle Sterile

LHC1001	Flexible 1 Per Pack	1/PK, 160 PK/CA
LHC1002	Flexible 2 Per Pack	2/PK, 80 PK/CA
LHC1003	Flexible 3 Per Pack	3/PK, 64 PK/CA
LHC3001	Rigid 1 Per Pack	1/PK, 64 PK/CA
LHC3002	Rigid 2 Per Pack	2/PK, 40 PK/CA

Sponge Counters

NOT MADE WITH natural rubber latex

- Clear visibility for monitoring sponge usage
- Accurate for assessing absorbed fluid
- Meet AORN recommendations for correct count
- Available with clear or blue back panel
- Attach easily to IV stand or wall
- 10 4 Inch x 4 Inch pockets that may be expanded to form five 4 Inch x 8 Inch pockets



Item No.	Description	Quantity
Bag Sponge Counter 10 or 5 Pockets		
SPCT1001	Clear Back	1/EA, 50 EA/BX, 5 BX/CA
SPCT1002	Blue Back	1/EA, 50 EA/BX, 5 BX/CA

Skin Adhesive

NOT MADE WITH natural rubber latex sterile EO

- LiquiBand Exceed
 - Pure 2-octyl cyanoacrylate formulation for maximum durability and flexibility
 - Porous felt applicator tip for even, consistent application without squirting and dripping
 - Unique “no clog” applicator tip
 - Elliptical shaped applicator allowing narrow or wide application strips
 - 0.8 mL volume per applicator
 - Easy to use winged applicator allowing for safe and easy activation and helps provide control throughout application
 - High patient satisfaction with atraumatic wound closure, low exothermic reaction and a water resistant barrier
- LiquiBand Flow Control
 - Fast-drying n-butyl formulation
 - Flow control tip provides greater control, safer application



Item No.	Description	Quantity
Adhesive Skin LiquiBand Sterile		
LBX001MC	Exceed Octyl 0.8 Gram	1/EA, 10 EA/BX, 6 BX/CA
LFC002MC	Flow Control N-Butyl 0.5 Gram	1/EA, 12 EA/BX, 6 BX/CA

Instrument Magnetic Drapes

NOT MADE WITH natural rubber latex sterile EO

- Flexible drapes easily conform to the shape of the surface
- Raised magnets reduce risk of instruments falling off the sterile field, allowing for easy removal and transfer of instruments
- Non-slip foam backing helps prevent shifting during procedures
- Non-magnetic transfer tray may be used for hands-free passing of sharp instruments, minimizing risk of sharps-related injuries
- Constructed of durable medical grade silicone



Item No.	Description	Quantity
Drape Magnetic Pad Foam Backing Medical Grade Silicone Sterile		
MAG1620	16 x 20	1/EA, 20 EA/CA
MAG1016	10 x 16	1/EA, 40 EA/CA
MAG1016S	10 x 16 with Stay 1 Neutral Zone	1/EA, 40 EA/CA

Surgical Blades

NOT MADE WITH natural rubber latex diverse supplier sterile R

- Computerized grinding and proprietary polishing processes ensure consistent, clean, sharp cutting edges
- Two VCI strips ensure product and packaging integrity (reduces risk of blade piercing foil packs)
- Produced from high quality, high durability Swiss Carbon Steel and Stainless Steel
- Individually packaged in easy-to-peel foil
- Color-coded packaging for quick and easy blade identification
- Manufactured in FDA inspected facilities and conform to ISO 7740 and ISO 13485



Item No.	Description	Quantity
3001T-10	Stainless Steel #10	1/EA, 100 EA/BX
3001T-11	Stainless Steel #11	1/EA, 100 EA/BX
3001T-15	Stainless Steel #15	1/EA, 100 EA/BX
2001T-10	Carbon Steel #10	1/EA, 100 EA/BX
2001T-11	Carbon Steel #11	1/EA, 100 EA/BX
2001T-15	Carbon Steel #15	1/EA, 100 EA/BX

SKIN MARKERS

Dual Tip Sterile Skin Markers

made in the USA NOT MADE WITH natural rubber latex sterile R

- Assists with Joint Commission's Universal Protocol for Preventing Wrong Site, Wrong Procedure, Wrong Person Surgery
- Regular tip on one end, fine tip on the other end
- Delivers an accurate and reliable mark to the incision site
- Gentian violet ink does not rub off during prep
- Barrel imprinted with ruler



Item No.	Description	Quantity
MKR7132	Dual-Tip Skin Marker, Regular & Fine Tip	1/EA, 100 EA/CA

Please see our full line of IV Stands on page 45 in the Equipment Section.

SKIN MARKERS

Standard Sterile Skin Markers



- Assists with Joint Commission's Universal Protocol for Preventing Wrong Site, Wrong Procedure, Wrong Person Surgery
- Tapered pen that minimizes surgeon hand fatigue
- Delivers an accurate and reliable mark to the incision site
- Gentian violet ink does not rub off during prep
- Conveniently packaged with ruler and/or labels
- Single patient use



Item No.	Description	Quantity
Skin Marker, Sterile		
77733	With Ruler	1/EA, 25 EA/BX, 4 BX/CA
77734	With Ruler and Labels	1/EA, 25EA/BX, 4 BX/CA

Skin and Utility Markers with Pre-printed Labels



- 38 pre-printed medication labels designed for quick identification in the operating room
- Epinephrine and Adrenaline labels printed in red
- 6 generic medication labels & 4 blank labels
- Black utility marker for waterproof and smearproof medication labeling
- Duo-tip skin marker designed for variable line thickness
- 6 Inch flexible ruler
- Assists in compliance with Joint Commission's National Patient Safety goal for medication labeling



Item No.	Description	Quantity
SURL048	Skin and Utility Marker with 6 Inch Flexible Ruler and 48 Pre-Printed Labels	1/EA, 25 EA/BX, 12 BX/CA

Correct Site Non-Sterile Skin Markers



- Assists with Joint Commission's Universal Protocol for Preventing Wrong Site, Wrong Procedure, Wrong Person Surgery
- Mini size skin markers are designed to reduce cost and biowaste generation and encourage disposal after single patient use
- Full size barrel diameter for hand comfort and better grip
- Mini skin markers contain Gentian violet ink and proof-of-use packaging
- Skin Markers available in individual poly sleeve or dispenser box



Item No.	Description	Quantity
Marker Skin Correct Site Non-Sterile Single-Patient Use		
MKR2100	Mini Gentian Violet Dispenser Box of 100 Pens	1/EA, 100 EA/BX, 12 BX/CA
MKR2250	With Sleeve Mini Gentian Violet	1/EA, 250 EA/BX, 5 BX/CA

SURGICAL BRAS

Surgical Vests



- Breast wrap supports ribs, sternum, and breast tissue
- Immobilizes side and breast tissue, eliminating the need for splints and taping
- Holds dressings in place without tape
- Elastic back panel provides comfort
- Velcro® brand hook and loop fasteners provides adjustable compression
- Adjustable padded straps for support and comfort



Item No.	Description	Quantity
Vest Surgical with Hook-Loop Fastener		
SVST102	Medium 36-38 Inch	1/EA
SVST103	Large 38-40 Inch	1/EA
SVST104	Extra Large 40-42 Inch	1/EA
SVST105	2X Large 42-44 Inch	1/EA
SVST106	3X Large 44-46 Inch	1/EA
SVST107	4X Large 46-48 Inch	1/EA
SVST108	5X Large 48-50 Inch	1/EA

Premium Surgical Bras



- Constructed of cool comfortable cotton/Spandex
- Front closure allows for easy access and application
- 2 Inch wide closure in front allows for adjustable fit and compression within each size
- Padded and adjustable shoulder straps to provide maximum comfort to tender tissue
- Shoulder strap is completely removable in front for easy application
- Removable drain bulb holder under each cup to keep drainage tube and bulb out of the way of surgeon and nursing team
- Hidden bust seam



Item No.	Description	Quantity
Bra Surgical Premium Cotton Spandex Hook and Loop Front Closure		
BRA101	Small 34-36 Inch	1/EA
BRA102	Medium 36-38 Inch	1/EA
BRA103	Large 38-40 Inch	1/EA
BRA104	Extra Large 40-42 Inch	1/EA
BRA105	2X Large 42-44 Inch	1/EA
BRA106	3X Large 44-46 Inch	1/EA
BRA107	4X Large 46-48 Inch	1/EA

SURGICAL BRAS

Standard Surgical Bras

NOT MADE WITH
natural rubber latex

- Full coverage therapeutic breast support provides light compression
- Hook and eye front closure for ease in dressing
- Each size designed to fit a range of B/C/D cups
- Constructed of comfortable, cotton/Spandex material
- Available in white and black
- Velcro® brand hook and loop fasteners on straps permit easy access to surgical site



Item No.	Description	Quantity
Bra Surgical Standard Cotton Spandex Hook and Eye Front Closure		
BRA232B	Black Size 32	1/EA
BRA232W	White Size 32	1/EA
BRA234B	Black Size 34	1/EA
BRA234W	White Size 34	1/EA
BRA236B	Black Size 36	1/EA
BRA236W	White Size 36	1/EA
BRA238B	Black Size 38	1/EA
BRA238W	White Size 38	1/EA
BRA240B	Black Size 40	1/EA
BRA240W	White Size 40	1/EA
BRA242B	Black Size 42	1/EA
BRA242W	White Size 42	1/EA

FOAM POSITIONING PRODUCTS

Foam Head Positioners

made in the USA **NOT MADE WITH**
natural rubber latex

- Molded from high-density, disposable foam
- Provide support, stability, elevation and safe alignment of head and neck
- Slotted head cradle is designed to accommodate an endotracheal tube
- Donut cut-out is left to allow for multiple uses
- Foam Density: 1.35 pounds



Item No.	Description	Quantity
Positioner Head 1.35 Pound Foam		
P06290	Donut Cutout 9 Inch	1/EA, 36 EA/CA
P06310	Cradle	1/EA, 24 EA/CA
P06320	Cradle Slotted	1/EA, 24 EA/CA

Abduction Pillows

- Molded from high-density, disposable foam
- Provide support and stability for proper healing after hip surgery
- Concave sides allow for snug fit
- Soft foam straps secure the patient
- Vacuum packed for easy storage
- Foam Density: 1.35 pounds



Item No.	Description	Quantity
Pillow Abduction Hip Contour 1.35 Pound Foam Disposable		
P06431	Small 18 x 12 x 6 Inch	1/EA, 6 EA/CA
P06432	Medium 22 x 15 x 6 Inch	1/EA, 6 EA/CA
P06433	Large 25 x 18 x 6 Inch	1/EA, 6 EA/CA

FOAM POSITIONING PRODUCTS

Foam Heel and Elbow Protectors

made in the USA **NOT MADE WITH**
natural rubber latex

- Molded from high-density, disposable foam
- Provide support during surgical procedures
- Reduces the risk of pressure sores and nerve damage
- Stability and protection for heels and elbows
- Foam Density: 1.35 pounds



Item No.	Description	Quantity
Protector Convoluted Foam		
P06110	Ulnar	1/PR, 36 PR/CA
P06120	Heel Pressure-Relieving	1/PR, 36 PR/CA

Foam Positioning Wedges

made in the USA **NOT MADE WITH**
natural rubber latex

- Molded from high density disposable foam
- Provides support and stability for safe positioning of patients
- Reduces risk of pressure sores and nerve damage
- Vacuum packed for easy storage




Item No.	Description	Quantity
P06545	Foam Positioning Wedge, Multi-Purpose, 45 Degree, 23 x 11 x 11 Inch	1/EA, 8 EA/CA

Please see our full line of Turn and Positioning Slider Sheets on page 169 in the Skin Care Section.

FOAM POSITIONING PRODUCTS

Arm Board Overlays

NOT MADE WITH natural rubber latex



- Ideal cushioning; better than pillows or towels
- Promotes air circulation
- Standard 16 Inch x 8 Inch size
- 1.5 Inch thick foam
- Meets California Bureau of Home Furnishings Bulletin 117

Item No.	Description	Quantity
PAD16815	Foam Armboard Overlay Pad, 16 x 18 Inch, 1.5 Inch thick	1/EA, 12 EA/CA

Cradle Positioners

made in the USA NOT MADE WITH natural rubber latex



- Molded from high density disposable foam
- Provides support and stability for extremities
- Reduces risk of pressure sores and nerve damage
- Vacuum packed for easy storage

Item No.	Description	Quantity
P06610	Foam Cradle Positioner, Multipurpose	1/EA, 12 EA/CA

STOCKINETTES

Tubular Stockinettes

NOT MADE WITH natural rubber latex



- Soft fabric provides comfortable wrap under casts
- Provide consistent pressure for wound dressings
- Packaged in 25-yard rolls
- Off-white color reduces glare from O.R. lighting

Item No.	Description	Quantity
----------	-------------	----------

Stockinette Orthopedic Tubular 1-Ply Cotton Non-sterile		
STK0225	2 Inch x 25 Yard	1/EA, 12 EA/CA
STK1325	3 Inch x 25 Yard	1/EA, 12 EA/CA
STK1425	4 Inch x 25 Yard	1/EA, 10 EA/CA
STK1525	5 Inch x 25 Yard	1/EA, 6 EA/CA
STK1625	6 Inch x 25 Yard	1/EA, 6 EA/CA
STK1825	8 Inch x 25 Yard	1/EA, 4 EA/CA
STK1025	10 Inch x 25 Yard	1/EA, 4 EA/CA
STK1225	12 Inch x 25 Yard	1/EA, 2 EA/CA

STOCKINETTES

Impervious Stockinettes

NOT MADE WITH natural rubber latex sterile EO



- Helps prevent contamination during procedures
- Moisture-repellent barrier helps guard against harmful bacteria
- Packaged in a peel pouch with low-memory CSR wrap
- Also available with pull-tabs for quick applications
- Textured synthetic liner with seamless, impervious elastomeric sleeve
- Packaged in peel pouch with CSR wrap

Item No.	Description	Quantity
----------	-------------	----------

Stockinette Orthopedic Impervious Moisture-Repellent Barrier Synthetic Sterile		
STK6630	6 x 30 Inch	1/EA, 18 EA/CA
STK6948	9 x 48 Inch	1/EA, 12 EA/CA
STK6248	12 x 48 Inch	1/EA, 10 EA/CA
STK5948	With Pull Tab 9 x 48 Inch	1/EA, 12 EA/CA
STK1248	With Pull Tab 12 x 48 Inch	1/EA, 10 EA/CA

Sterile 1-Ply Stockinettes

NOT MADE WITH natural rubber latex sterile EO



- Pre-cut in convenient sizes
- Sewn ends prevent fraying of fabric
- Off-white color reduces glare from O.R. lighting
- Rolled and packaged in a peel pouch with low-memory CSR wrap
- Fabric blend is 50/50 polyester/cotton offering benefits of both cotton and polyester yarns; the breathability of cotton and the natural wicking of polyester

Item No.	Description	Quantity
----------	-------------	----------

Stockinette Tubular Pre-Cut 1-Ply Sewn End Cotton Polyester off-White		
STK2336	3 x 36 Inch	1/EA, 40 EA/CA
STK2436	4 x 36 Inch	1/EA, 40 EA/CA
STK2448	4 x 48 Inch	1/EA, 36 EA/CA
STK2460	4 x 60 Inch	1/EA, 30 EA/CA
STK2648	6 x 48 Inch	1/EA, 30 EA/CA
STK2660	6 x 60 Inch	1/EA, 30 EA/CA
STK2672	6 x 72 Inch	1/EA, 24 EA/CA

Please see our full line of Mattress Overlays on page 168 in the Skin Care Section.



STOCKINETTES



Sterile 2-Ply Stockinettes

NOT MADE WITH natural rubber latex **sterile EO**

- Pre-cut in convenient sizes
- Sewn ends prevent fraying of fabric
- Off-white color reduces glare from O.R. lighting
- Rolled and packaged in a peel pouch with low-memory CSR wrap
- Fabric blend is 50/50 polyester/cotton offering benefits of both cotton and polyester yarns; the breathability of cotton and the natural wicking of polyester

Item No.	Description	Quantity
Stockinette Tubular Pre-Cut 2-Ply Sewn End Cotton Polyester off-White Sterile		
STK3336	3 × 36 Inch	1/EA, 36 EA/CA
STK3436	4 × 36 Inch	1/EA, 24 EA/CA
STK3448	4 × 48 Inch	1/EA, 20 EA/CA
STK3460	4 × 60 Inch	1/EA, 18 EA/CA
STK3648	6 × 48 Inch	1/EA, 18 EA/CA
STK3660	6 × 60 Inch	1/EA, 12 EA/CA
STK3672	6 × 72 Inch	1/EA, 12 EA/CA
STK3860	8 × 60 Inch	1/EA, 12 EA/CA
STK3872	8 × 72 Inch	1/EA, 12 EA/CA

Bias-Cut Stockinettes

NOT MADE WITH natural rubber latex **sterile EO**

- Soft, comfortable cotton blend offering benefits of both cotton and polyester yarns; the breathability of cotton and the natural wicking of polyester
- Off-white color reduces glare from O.R. lighting
- Single layers conform to shape of limbs and other body parts
- Sterile available in four widths, 3 Inch to 8 Inch
- Sterile packaged in a peel pouch with CSR wrap

Item No.	Description	Quantity
Stockinette Orthopedic Bias-Cut 1-Ply Blend Cotton Polyester Off-White Sterile		
STK4034	3 Inch × 4 Yard	1/EA, 20 EA/CA
STK4044	4 Inch × 4 Yard	1/EA, 20 EA/CA
STK4045	4 Inch × 5 Yard	1/EA, 20 EA/CA
STK4064	6 Inch × 4 Yard	1/EA, 18 EA/CA
STK4084	8 Inch × 4 Yard	1/EA, 15 EA/CA

Gait Belts

NOT MADE WITH
natural rubber latex

- Provides a secure hand hold for therapists and caregivers while working with patients
- Heavy-duty cotton webbing treated with bacteriostatic agent to provide long-term protection even after repeated washings
- Strong, durable nickel-plated metal buckle offers maximum hold while reducing wear on the cotton belt
- Available in natural color or pastel stripes
- 2 inches wide and available in 54-, 60-, or 72-inch lengths



Item No.	Description	Quantity
Gait Belt Heavy-Duty Cotton Webbing Treated with Bacteriostatic Agent Nickel-Plated Metal Buckle		
GAIT541	Pastel Stripe 2 x 54 Inch	1/EA, 12 EA/CA
GAIT542	Natural 2 x 54 Inch	1/EA, 12 EA/CA
GAIT601	Pastel Stripe 2 x 60 Inch	1/EA, 12 EA/CA
GAIT602	Natural 2 x 60 Inch	1/EA, 12 EA/CA
GAIT721	Pastel Stripe 2 x 72 Inch	1/EA, 12 EA/CA
GAIT722	Natural 2 x 72 Inch	1/EA, 12 EA/CA

Patient Restraints

NOT MADE WITH
natural rubber latex

diverse
supplier

- Provides a soft, safe, and secure fit around patient wrists
- Ambidextrous design with one size that fits most adults
- Made of machine washable cotton, polyester, and polyester fiberfill
- Optional straps included with all MediChoice Padded Mitts except MTT2811



Item No.	Description	Quantity
Mitts One Size Fits All Pair Cotton Polyester and Polyester Fiberfill		
MTT2810	Soft Hand No Finger Separators Breathable Cover Polystyrene Bead Fill Hook-and-Loop Wrist Closure	1/PR, 12 PR/CA
MTT2811	Easy-View Flap to View Fingers No Finger Separators Double-Padded Hook-and-Loop Wrist Closure	1/PR, 12 PR/CA
MTT2813	Assured Control with Mesh Back 5 Finger Separators Quick-Release Buckle Wrist Closure Optional 48 Inch Strap	1/PR, 12 PR/CA
MTT2814	Essential Secure with Mesh Back No Finger Separators Hook-and-Loop Wrist Closure Optional 48.5 Inch Strap	1/PR, 12 PR/CA
MTT2816	Finger Control with Mesh Back 5 Finger Separators Hook-and-Loop Wrist Closure Optional 48.5 Inch Strap	1/PR, 12 PR/CA
MTT2819	Double Essential Secure with Mesh Back No Finger Separators Double-Padded Hook-and-Loop Wrist Closure Optional 48.75 Inch Strap	1/PR, 12 PR/CA

Heel and Elbow Protectors

NOT MADE WITH
natural rubber latex

diverse
supplier

- Padding shields sensitive skin areas and minimizes friction and shear on elbows and heels
- Knitted flexible stretch weave conforms around the heel or elbow to help hold product in place
- No irritating straps, ties, or seams
- Sewn-in contour foam pad eliminates rolling
- Color-coded sizing
- Machine washable



Item No.	Description	Quantity
Protector Heel and Elbow Padded Knit 72 Percent Acrylic 20 Percent Spandex 8 Percent Nylon Fabric		
EHP1001	Small Yellow	1/PR, 12 PR/CA
EHP1002	Medium Blue	1/PR, 12 PR/CA
EHP1003	Large White	1/PR, 12 PR/CA
EHP1004	Extra-Large Green	1/PR, 12 PR/CA
EHP1005	2X Large Beige	1/PR, 12 PR/CA

Abdominal Binders

NOT MADE WITH
natural rubber latex

- Knit elastic panels provide post-operative support for abdominal muscles
- Hook-and-loop closure provides a more secure, more precise, and more comfortable fit
- Stitching between the horizontal panels allows the binder to conform to the body contours and stay in position
- Plush inner surface and butted seams help ensure patient comfort
- Durable stitching allows cut-away for drainage tubing
- Binders can be combined to accommodate patient girth



Item No.	Description	Quantity
Binder Abdominal Universal Knit Elastic Hook-and-Loop Closure		
ABB1050	10 Inch Wide Stretches to Fit 26 to 50 Inch	1/EA, 50 EA/CA
ABB3P45	Tri-panel 9 Inch Wide Stretches to Fit Up to 45 Inch	1/EA, 50 EA/CA
ABB3P62	Tri-panel 9 Inch Wide Stretches to Fit 46 to 62 Inch	1/EA, 50 EA/CA
ABB3P72	Tri-panel 9 Inch Wide Stretches to Fit 55 to 72 Inch	1/EA, 50 EA/CA
ABB4P45	Four-panel 12 Inch Wide Stretches to Fit Up to 45 Inch	1/EA, 50 EA/CA
ABB4P62	Four-panel 12 Inch Wide Stretches to Fit 46 to 62 Inch	1/EA, 50 EA/CA
ABB4P72	Four-panel 12 Inch Wide Stretches to Fit 55 to 72 Inch	1/EA, 50 EA/CA
ABB4P96	Four-panel 12 Inch Wide Stretches to Fit 72 to 96 Inch	1/EA, 50 EA/CA

Please see additional Patient Safety products, including Bathroom Safety Aids, on pages 131-132 in the Mobility Aids Section.

CERVICAL COLLARS

Two Piece Cervical Collars

NOT MADE WITH natural rubber latex

- Reduces rotation and slightly hyperextends the cervical spine
- Provides rigid immobilization
- Excellent for emergency tracheotomies and quick access
- Molded and shaped to conform to any patient
- Constructed of lightweight, closed-cell foam material
- Hook-and-loop closure facilitates easy application and removal



Item No.	Description	Quantity
----------	-------------	----------

Collars Cervical Two Piece 1.25 Inch High (Measured from the Chin to the Top of the Sternum)

OSG2209	Infant 6-8 Inch Neck Circumference	1/EA, 12 EA/CA
---------	------------------------------------	----------------

Collars Cervical Two Piece 1.75 Inch High (Measured from the Chin to the Top of the Sternum)

OSG2208	Pediatric 8-11 Inch Neck Circumference	1/EA, 12 EA/CA
---------	--	----------------

Collars Cervical Two Piece 2.25 Inch High (Measured from the Chin to the Top of the Sternum)

OSG2201	Small 10-13 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2202	Medium 13-16 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2203	Large 16-19 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2204	Extra-Large 19 Inch or Greater Neck Circumference	1/EA, 12 EA/CA

Collars Cervical Two Piece 3.25 Inch High (Measured from the Chin to the Top of the Sternum)

OSG2301	Small 10-13 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2302	Medium 13-16 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2303	Large 16-19 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2304	Extra-Large 19 Inch or Greater Neck Circumference	1/EA, 12 EA/CA

Collars Cervical Two Piece 4.25 Inch High (Measured from the Chin to the Top of the Sternum)

OSG2401	Small 10-13 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2402	Medium 13-16 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2403	Large 16-19 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2404	Extra-Large 19 Inch or Greater Neck Circumference	1/EA, 12 EA/CA

Collars Cervical Two Piece 5.25 Inch High (Measured from the Chin to the Top of the Sternum)

OSG2501	Small 10-13 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2502	Medium 13-16 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2503	Large 16-19 Inch Neck Circumference	1/EA, 12 EA/CA
OSG2504	Extra-Large 19 Inch or Greater Neck Circumference	1/EA, 12 EA/CA

CERVICAL COLLARS

Foam Cervical Collars

NOT MADE WITH natural rubber latex

- Styled to give uniform support
- Hook-and-loop closure
- Breathable and washable stockinette cover
- Choice of 1 inch medium- or firm-density foam



Item No.	Description	Quantity
----------	-------------	----------

Collar Cervical Serpentine Covered with Stockinette Hook-Loop Closure 1-inch Medium-Density Foam

OSG2101	Small 3.5 Inch High × 18.25 Inch Long	1/EA
OSG2102	Medium 4 Inch High × 20.25 Inch Long	1/EA
OSG2103	Large 4.5 Inch High × 20.5 Inch Long	1/EA

Collar Cervical Serpentine Covered with Stockinette Hook-Loop Closure 1-inch Firm-Density Foam

OSG2111	Small 3.5 Inch High × 18.25 Inch Long	1/EA
OSG2112	Medium 4 Inch High × 20.25 Inch Long	1/EA
OSG2113	Large 4.5 Inch High × 22.5 Inch Long	1/EA

Arm Slings

NOT MADE WITH natural rubber latex

- Easy to secure arm sling for immobilization and support of arm
- Breathable polyester blend
- Adjustable, wide shoulder strap with hook and loop closure for perfect fit
- Neck pad offers optimal patient comfort
- Universal sling can be folded to proper length to accommodate most patients



Item No.	Description	Quantity
----------	-------------	----------

Sling Arm Envelope Style Breathable Navy Polyester with Shoulder Pad Adjustable Hook-and-Loop Closure

ASL2102	Extra-Small Child 7 × 11 Inch	1/EA, 24 EA/CA
ASL2103	Small 7.5 × 12.5 Inch	1/EA, 24 EA/CA
ASL2104	Medium 8 × 15 Inch	1/EA, 24 EA/CA
ASL2105	Large 8.5 × 17.5 Inch	1/EA, 24 EA/CA
ASL2106	Extra-Large 9 × 20 Inch	1/EA, 24 EA/CA
ASL2150	Universal 7 × 18 Inch	1/EA, 24 EA/CA

SHOULDER IMMOBILIZERS

Universal Shoulder Immobilizers

NOT MADE WITH natural rubber latex made in the USA

- Our most popular shoulder immobilizer
- Forearm cuff, humeral cuff and shoulder strap may be placed in desired position
- Web shoulder strap with buckle adjustment
- Available in straight or contoured
- Foam laminate construction with hook-and-loop closure
- Fits rib circumferences between 24 and 50 inches
- Fits either left or right arm



Item No.	Description	Quantity
Immobilizer Shoulder Web Shoulder Strap with Buckle Adjustment Universal Foam Laminate 24-50 Inch		
OSG304C	Contoured	1/EA
OSG304S	Straight	1/EA

Sling and Swathe

NOT MADE WITH natural rubber latex made in the USA

- Secures arm against body
- May be used as an arm sling if swathe is not desired
- Fits either left or right arm
- Soft foam construction
- Hook-and-loop closure



Item No.	Description	Quantity
OSG302U	Immobilizer Shoulder Sling Swathe Soft Foam Hook-and-Loop Closure Universal Adult 54 Inch	1/EA, 6 EA/CA

Velpeau Shoulder Immobilizers

NOT MADE WITH natural rubber latex

- Breathable, cool cotton material
- Removable foam swathe strap
- O-Ring on strap for easy adjustment
- Hook-and-loop closure
- Fits either left or right arm



Item No.	Description	Quantity
Immobilizer Shoulder Velpeau Hook-and-Loop Closure with Waist Strap Cotton		
OSG3031	Small 38 Inch Shoulder Strap 36 Inch Waist Strap	1/EA
OSG3032	Medium 45 Inch Shoulder Strap 39 Inch Waist Strap	1/EA
OSG3033	Large 48 Inch Shoulder Strap 43 Inch Waist Strap	1/EA
OSG3034	Extra-Large 52 Inch Shoulder Strap 46 Inch Waist Strap	1/EA

KNEE IMMOBILIZERS

Universal Knee Immobilizers

NOT MADE WITH natural rubber latex made in the USA

- Provides support and stabilization for the leg and knee
- Posterior stays feature a center I-beam construction and are contoured and reinforced to provide rigidity and a custom fit
- Two movable stays for medial and lateral placement
- Soft outer foam laminate material adds comfort and shows less staining
- Trimmable, universal design fits most patients
- Hook-and-loop straps for easy closure



Item No.	Description	Quantity
Immobilizer Knee Hook-and-Loop Closure Provides Support and Stabilization for the Leg and Knee Universal Foam Laminate		
OSG5012	12 Inch	1/EA, 6 EA/CA
OSG5014	14 Inch	1/EA, 12 EA/CA
OSG5016	16 Inch	1/EA, 12 EA/CA
OSG5018	18 Inch	1/EA, 12 EA/CA
OSG5020	20 Inch	1/EA, 12 EA/CA
OSG5022	22 Inch	1/EA, 12 EA/CA
OSG5024	24 Inch	1/EA, 12 EA/CA
OSG5026	26 Inch	1/EA, 12 EA/CA
OSG5028	28 Inch	1/EA, 12 EA/CA

Sizing Chart

Reorder No.	Size	Length
OSG5012	Universal	12 Inches
OSG5014	Universal	14 Inches
OSG5016	Universal	16 Inches
OSG5018	Universal	18 Inches
OSG5020	Universal	20 Inches
OSG5022	Universal	22 Inches
OSG5024	Universal	24 Inches
OSG5026	Universal	26 Inches
OSG5028	Universal	28 Inches

Measure 2-3 inches above the lateral malleolus (the outer bump of the ankle) and 6-8 inches below the femoral head (hip joint) to determine the required length

KNEE IMMOBILIZERS

Universal Knee Immobilizers with Compression Straps

NOT MADE WITH natural rubber latex made in the USA

- Elastic hook and loop straps for quick and easy application
- Two movable stays for medial and lateral placement
- Reinforced posterior stays for extra rigidity
- Contoured posterior stays
- Universal sizing will fit most patients



Item No.	Description	Quantity
Immobilizer Knee with Elastic Compression Strap Contoured Universal		
OSG5116	16 Inch	1/EA, 12 EA/CA
OSG5118	18 Inch	1/EA, 12 EA/CA
OSG5120	20 Inch	1/EA, 12 EA/CA
OSG5122	22 Inch	1/EA, 12 EA/CA
OSG5124	24 Inch	1/EA, 12 EA/CA

Sizing Chart

Reorder No.	Posterior Stays	Size	Length
OSG5116	Contoured	Universal	16 Inches
OSG5118	Contoured	Universal	18 Inches
OSG5120	Contoured	Universal	20 Inches
OSG5122	Contoured	Universal	22 Inches
OSG5124	Contoured	Universal	24 Inches

Measure 2-3 inches above the lateral malleolus (the outer bump of the ankle) and 6-8 inches below the femoral head (hip joint) to determine the required length

KNEE IMMOBILIZERS

Tri-Panel Universal Knee Immobilizers

NOT MADE WITH natural rubber latex

- Tri-panel design for more anatomical fit
- Movable side panels for medial and lateral stay positioning
- Two medial and lateral stays
- Three contoured posterior stays for better fit
- Foam laminate material with tricot lining
- Hook-and-loop closure



Item No.	Description	Quantity
Immobilizer Knee Tripanel with Hook-and-Loop Closure Universal		
OSG5212	12 Inch	1/EA, 6 EA/CA
OSG5214	14 Inch	1/EA, 12 EA/CA
OSG5216	16 Inch	1/EA, 12 EA/CA
OSG5218	18 Inch	1/EA, 12 EA/CA
OSG5220	20 Inch	1/EA, 12 EA/CA
OSG5222	22 Inch	1/EA, 12 EA/CA
OSG5224	24 Inch	1/EA, 12 EA/CA

Sizing Chart

Reorder No.	Size	Length
OSG5212	Universal	12 Inches
OSG5214	Universal	14 Inches
OSG5216	Universal	16 Inches
OSG5218	Universal	18 Inches
OSG5220	Universal	20 Inches
OSG5222	Universal	22 Inches
OSG5224	Universal	24 Inches

Measure 2-3 inches above the lateral malleolus (the outer bump of the ankle) and 6-8 inches below the femoral head (hip joint) to determine the required length

ANKLE BRACES

Ankle Stirrup Braces with Foam and Terry Liners

NOT MADE WITH natural rubber latex

- Designed to help support and protect unstable or injured ankles
- Terry-covered foam liner helps provide uniform compression and offers additional patient comfort
- Standard, 10.5-inch and small, 8-inch braces designed for acute injuries while Trainer, 9-inch brace is used primarily for chronic indications
- Easy application and proper fit with adjustable heel pads, swivel straps, and vertical side straps
- Re-engineered shell is light, strong, and conforms to the ankle with minimal bulk in the shoe
- Fits either left or right ankle



Item No.	Description	Quantity
Brace Ankle Stirrup with Terry-Covered Foam Liner Fits Right or Left Ankle		
BRC4001	Standard 10.5 Inch	1/EA
BRC4002	Trainer 9 Inch	1/EA
BRC4003	Small 8 Inch	1/EA

ANKLE BRACES

Ankle Stirrup Braces with Gel and Air

NOT MADE WITH natural rubber latex

- Designed to help support and protect unstable or injured ankles
- Dual chamber design combines lightweight gel, air, and closed cell-foam materials for maximum support, compression, and cold therapy
- Standard, 10.5-inch brace designed for acute injuries while Trainer, 9-inch brace is used primarily for chronic indications
- Easy application and proper fit with adjustable heel pads, swivel straps, and vertical side straps
- Semi-rigid shell is light, strong, and conforms to the ankle with minimal bulk in the shoe
- Fits either left or right ankle



Item No.	Description	Quantity
Brace Ankle Stirrup Lightweight with Gel and Air Bladder Fits Right or Left Ankle		
BRC3001	Standard 10.5 Inch	1/EA, 12 EA/CA
BRC3002	Trainer 9 Inch	1/EA, 12 EA/CA

Ankle Stirrup Braces with Air

NOT MADE WITH natural rubber latex

- Designed to help support and protect unstable or injured ankles
- Pre-inflated air bladder combines air and foam for maximum support and comfort
- Valve and tubing provided to adjust bladder to desired compression
- Standard, 10.5-inch brace designed for acute injuries while Trainer, 9-inch brace is used primarily for chronic indications
- Easy application and proper fit with adjustable heel pads, swivel straps, and vertical side straps
- Semi-rigid shell is light, strong, and conforms to the ankle with minimal bulk in the shoe
- Fits either left or right ankle



Item No.	Description	Quantity
Brace Ankle Stirrup Lightweight with Air Bladder Fits Right or Left Ankle		
BRC2001	Standard 10.5 Inch	1/EA, 12 EA/CA
BRC2002	Trainer 9 Inch	1/EA, 12 EA/CA

Ankle Stirrup Rehabilitation Kit

NOT MADE WITH natural rubber latex

- Complete system for the functional management of the injured ankle from the initial injury through rehabilitation
- Standard MediChoice Ankle Brace with Air provides support and protection for the injured ankle
- Elastic bandage provides circumferential compression to control swelling
- Flexible cold pack provides cold therapy to help control swelling and pain
- Exercise band provides for ankle rehabilitation exercises
- Application and exercise instructions provide necessary patient information for brace application and ankle rehabilitation



Item No.	Description	Quantity
BRC2050	Kit Rehabilitation Ankle Stirrup Includes Standard Ankle Brace with Air Elastic Bandage Flexible Cold Pack Exercise Band and Instructions	1/EA

ANKLE BRACES

Figure Eight Ankle Braces

NOT MADE WITH
natural rubber latex

- Designed to support and protect ankles
- Figure eight configuration eliminates taping
- Thin and strong ballistic-type nylon provides maximum support and durability
- Knit tongue is breathable and comfortable
- No seam on the foot for greater comfort
- Elastic panel at the Achilles area for added comfort
- Fits either left or right ankle



Item No.	Description	Quantity
Brace Ankle Figure Eight Strap with Knit Tongue Black Nylon		
BRC6001	Small 11-12 Inch	1/EA, 4 EA/CA
BRC6002	Medium 12-13 Inch	1/EA, 4 EA/CA
BRC6003	Large 13-14 Inch	1/EA, 4 EA/CA
BRC6004	Extra-Large 14-15 Inch	1/EA, 4 EA/CA

Sizing Chart

Reorder No.	Length	Size	Male Shoe Size	Female Shoe Size
BRC6001	11-12 Inches	Small	6-7	7-9
BRC6002	12-13 Inches	Medium	8-10	10-11
BRC6003	13-14 Inches	Large	11-12	12-13
BRC6004	14-15 Inches	Extra Large	13-14	14+

ANKLE WALKERS

Low Profile Fixed Ankle Walkers

NOT MADE WITH
natural rubber latex

- High-top walker designed for the treatment of stable fractures and ankle sprains
- Deluxe foam liner designed for greater patient comfort
- Wider footbed provides a more comfortable fit
- Lightweight design includes contoured, metal reinforced uprights for greater strength and a more anatomical fit at the calf
- Low profile, continuous arching rocker design promotes a smooth, energy-conserving gait
- Cushioned heel for greater energy absorption
- Padded strap included for optional padding at ankle



Item No.	Description	Quantity
Walker Ankle Fixed High-Top Deluxe Foam Liner Metal-Reinforced Uprights		
OSG6010	Extra-Small	1/EA, 6 EA/CA
OSG6011	Small	1/EA, 6 EA/CA
OSG6012	Medium	1/EA, 6 EA/CA
OSG6013	Large	1/EA, 6 EA/CA
OSG6014	Extra-Large	1/EA, 6 EA/CA

Sizing Chart

Reorder No.	Size	Male Shoe Size	Female Shoe Size
OSG6010	Extra-Small	3½	4½
OSG6011	Small	5-8½	4-7
OSG6012	Medium	9-11	7½-10
OSG6013	Large	11½-13	10½-12
OSG6014	Extra Large	13½+	12½+

Air Shell Ankle Walkers

NOT MADE WITH
natural rubber latex

- Three-piece design includes posterior and anterior shells for maximum immobilization
- Dual air chambers positioned medial/lateral at the ankle can be inflated for ideal compression and a custom fit
- Lightweight low profile design and rocker sole facilitate a natural gait
- Deluxe liner includes a detachable toe cover
- Wider foot base design for greater patient comfort and support
- Pair of walker socks included



Item No.	Description	Quantity
Walker Ankle Air Shell Lightweight Low Profile High-Top with Socks		
OSG6021	Small	1/EA, 6 EA/CA
OSG6022	Medium	1/EA, 6 EA/CA
OSG6023	Extra-Large	1/EA, 6 EA/CA

Sizing Chart

Reorder No.	Size	Male Shoe Size	Female Shoe Size
OSG6021	Small	4-7	5-8
OSG6022	Medium	7½-11	8½-12
SG6023	Large/XL	11½+	12½+

← Please see our full line of Mobility Aids, including Crutches, Canes, and Walkers, on pages 134-136 in the Mobility Aids Section.

NIGHT SPLINTS

Plantar Fasciitis Night Splints

NOT MADE WITH natural rubber latex

- Designed to provide relief from pain and discomfort of plantar fasciitis and Achilles tendonitis
- Allows the patient to apply a consistent 90 degree angle stretch to the plantar fascia by holding the foot in gentle dorsiflexion
- Soft, brushed nylon and terry lining provides greater patient comfort and protection against irritation and skin breakdown
- Removable, foam wedge insert provides additional 5 inches of stretch on the plantar fascia if needed
- Hook-and-loop closure
- Fits either left or right foot



Item No.	Description	Quantity
----------	-------------	----------

Splint Night Plantar Fasciitis Hook-and-Closure Nylon and Terry Lining Removable Foam Wedge Insert

OSG7001	Small	1/EA
OSG7002	Medium	1/EA
OSG7003	Large	1/EA

Sizing Chart

Reorder No.	Size	Male Shoe Size	Female Shoe Size
OSG7001	Small	3-5	5-7
OSG7002	Medium	5½-6½	7½-8½
OSG7003	Large	7-8	9-10

Plantar Fasciitis Stretch Splints

NOT MADE WITH natural rubber latex

- Designed to provide relief from pain and discomfort of plantar fasciitis and Achilles tendonitis
- Allows the patient to apply consistent stretch to the plantar fascia by holding the foot in gentle dorsiflexion
- Greater than 90 degree stretch possible
- Lightweight design for patient comfort while sleeping
- Slip-resistant sole for greater patient safety
- Fits either left or right foot



Item No.	Description	Quantity
----------	-------------	----------

Splint Night Stretch Plantar Fasciitis Small Hook-and-Loop Closure

OSG7011	Small	1/EA
OSG7012	Medium	1/EA
OSG7013	Large	1/EA

Sizing Chart

Reorder No.	Size	Male Shoe Size	Female Shoe Size
OSG7011	Small	Up to 7	Up to 7½
OSG7012	Medium	7½-10	8-10
OSG7013	Large	10½+	11+

POST OP SHOES

Classic Post Op Shoes

NOT MADE WITH natural rubber latex

- Canvas, foam-padded upper
- Hook-and-loop closure
- Soft inner sole bonded to wood
- Slip-resistant, ribbed sole



Item No.	Description	Quantity
----------	-------------	----------

Shoe Post Op Classic Black Canvas Foam-Padded Hook-and-Loop Closure

OSG6328	Pediatric	1/EA, 12 EA/CA
OSG6321	Women Small	1/EA, 12 EA/CA
OSG6322	Women Medium	1/EA, 12 EA/CA
OSG6323	Women Large	1/EA, 12 EA/CA
OSG6331	Men Small	1/EA, 12 EA/CA
OSG6332	Men Medium	1/EA, 12 EA/CA
OSG6333	Men Large	1/EA, 12 EA/CA
OSG6334	Men Extra-Large	1/EA, 12 EA/CA

Sizing Chart

Reorder No.	Size	Shoe Size
OSG6328	Pediatric	12-1
OSG6321	Women Small	4-6
OSG6322	Women Medium	6½-8
OSG6323	Women Large	8½-10
OSG6331	Men Small	6-8
OSG6332	Men Medium	8½-10
OSG6333	Men Large	10½-12
OSG6334	Men X-Large	12½-14

Alcohol Prep Pads

NOT MADE WITH
natural rubber latex



- For injection site preparation and immediate disinfection of the skin or wound area
- Saturated with 70-percent isopropyl alcohol
- Manufactured according to current FDA regulations, GMP and USP requirements
- Meets highest bacteriological standards
- All pads 2-ply, packaged in individual foil pouches, and sterile
- Available in large and medium sizes
- Single patient use
- Pad Size:
 - Medium 1.19 × 2.63 Inch
 - Large: 1.75 × 3.5 Inch

Item No.	Description	Quantity
Pad Prep Alcohol Saturated With 70 Percent Isopropyl Alcohol 2-Ply Sterile Single-Patient Use		
APP102A	Medium 1.19 Inch × 2.63 Inch	1/EA, 200 EA/BX, 20 BX/CA
APP104A	Large 1.75 Inch × 3.5 Inch	1/EA, 100 EA/BX, 10 BX/CA

TOWELETTES

BZK Towelettes

NOT MADE WITH
natural rubber latex



- For use whenever antiseptic cleansing/refreshing is necessary
- Ideal for general cleansing of hands, face, and skin areas
- Alcohol-free formula ensures stingless cleansing for minor cuts, scrapes and burns
- Saturated with 0.13-percent Benzalkonium Chloride (BZK).
- Cleans and cools simultaneously
- Non-toxic and non-irritating formula mild enough for everyday use

Item No.	Description	Quantity
BZK100	Towelette Benzalkonium Chloride Antiseptic Non-Toxic Non-Irritating 5 × 8 Inch Sterile	1/EA, 100 EA/BX, 10 BX/CA

Cleansing Towelettes

NOT MADE WITH
natural rubber latex



- Cleanses and refreshes on-the-go without the need for soap and water
- Intended for gentle cleansing of the face, hands, and body as needed
- Non-toxic and non-irritating antiseptic ideal for personal hygiene, gentle enough for every-day use
- Saturated with 5-percent Isopropyl Alcohol and 0.5-percent Benzalkonium Chloride
- Individually packaged
- Towelette size is 5 × 7.75 Inches

Item No.	Description	Quantity
CLT100	Towelette Cleansing 5 × 7.75 Inch Sterile	1/EA, 100 EA/BX, 10 BX/CA

Obstetrical Towelettes

NOT MADE WITH
natural rubber latex



- For personal hygiene at every napkin change and whenever cleansing/refreshing is needed
- Cleans and cools simultaneously
- Aids in bactericidal cleansing for patient use
- Saturated with Benzalkonium Chloride (0.13-percent), SD Alcohol-40 (1-percent), Propylene Glycol (1-percent), Camphor (0.05-percent), Sodium Hydroxide (0.05-percent), Purified Water (97.77-percent)
- 5 × 8 Inch towelettes come packed in individual foil pouches
- Single-patient use

Item No.	Description	Quantity
OBT100	Towelette Obstetrical Saturated Single-Patient Use 5 × 8 Inch Sterile	1/EA, 100 EA/BX, 10 BX/CA

Nail Polish Remover Pads

NOT MADE WITH
natural rubber latex



- For the removal of nail polish when respiratory evaluation is necessary
- One pad sufficient for up to ten nails
- Individual foil packets keep pads moist
- Single patient use
- Contains Methyl Acetate, Butoxydiglycol, Sodium Citrate, Denatonium Benzoate, Tocopheryl Acetate, and distilled water

Item No.	Description	Quantity
NPR100	Pad Remover Nail Polish Acetone-Free Single-Patient Use 1.25 × 2.63 Inch	1/EA, 100 EA/BX, 10 BX/CA

Adhesive Tape Remover Wipes

NOT MADE WITH
natural rubber latex



- Aids in the removal of tape and other adhesive residues from skin
- Non-irritating and non-sensitizing formula
- Easy to use in a convenient, individually-wrapped foil package
- Single-patient use

Item No.	Description	Quantity
ARW050	Wipe Adhesive Tape Remover 30 × 60 Millimeter	1/EA, 50 EA/BX, 20 BX/CA

Lemon Glycerin Swabs

NOT MADE WITH
natural rubber latex

- Gentle and soothing; ideal for NPO and ventilator patients
- Foil packets keep swabs moist and aid in the prevention of cross contamination
- Glycerin stimulates salivary glands and flavoring stimulates taste buds
- Each packet contains three 4-Inch swabs
- Contains Glycerin (U.S.P.), Citric Acid (U.S.P.), and Lemon Flavoring



Item No.	Description	Quantity
LGS100	Swabsticks Lemon Glycerin 3-Pack 4-Inch	3/PK, 25 PK/BX, 10 BX/CA

Oral Swabettes

made in the USA **NOT MADE WITH**
natural rubber latex

- Rigid foam head helps gently loosen debris and mucus while stimulating oral tissue
- Gentle enough for fragile and tender oral tissues
- Extra-long 6-Inch handle to help avoid cross contamination
- Available in unflavored or fresh mint flavor
- Plastic stem won't absorb oral secretions
- Variety of packaging options: bulk, 10/pack, and 20/pack



Item No.	Description	Quantity
Swabette Oral Rigid Foam Head Plastic 6 Inch Stem		
OSB101	Unflavored Unwrapped Pink	1/EA, 250 EA/BX, 4 BX/CA
OSB102	Unflavored Unwrapped Pink	20/PK, 25 PK/BX, 2 BX/CA
OSB103	Unflavored Individually Wrapped Bulk Pink	1/EA, 10 EA/PK, 25 PK/BX, 4 BX/CA
OSB201	Mint Flavored Individually Wrapped Bulk Green	1/EA, 250 EA/BX, 4 BX/CA
OSB202	Mint Flavored Unwrapped Green	20/PK, 25 PK/BX, 2 BX/CA
OSB203	Mint Flavored Individually Wrapped Green	1/EA, 10 EA/PK, 25 PK/BX, 4 BX/CA

Alcohol Swabsticks

- Use as a topical antiseptic
- Ideal for use prior to injections
- Effective for cleaning minor cuts, scrapes, and burns to reduce risk of infection
- Non-woven swabs saturated with 70-percent isopropyl alcohol, USP



Item No.	Description	Quantity
Swabsticks Alcohol Sterile		
APS1000	1 Per Pack	1/PK, 50 PK/BX, 20 BX/CA
APS3000	3 Per Pack	3/PK, 25 PK/BX, 20 BX/CA

Instant Hand Sanitizers

made in the USA **NOT MADE WITH**
natural rubber latex **diverse supplier**

- Kills 99.9-percent of germs without the need for water or rinsing
- Enriched with Aloe Vera, Vitamin E, and other emollients to condition and moisturize skin during repeated use
- Available in pleasant citrus fragrance or fragrance-free
- Gel contains 62-percent ethyl alcohol
- Meets CDC standards for proper hand washing protocols
- Hypoallergenic



Item No.	Description	Quantity
Sanitizer Hand Gel 62 Percent Ethyl Alcohol		
88300	Citrus-Scented Flip-Top Cap 4 Ounce	1/EA, 12 EA/BX, 8 BX/CA
088300FF	Fragrance-Free Flip-Top Cap 4 Ounce	1/EA, 12 EA/BX, 8 BX/CA
88303	Citrus-Scented Pump Dispenser 8 Ounce	1/EA, 12 EA/BX, 4 BX/CA
088303FF	Fragrance-Free Pump Dispenser 8 Ounce	1/EA, 12 EA/BX, 4 BX/CA

Tongue Depressors

NOT MADE WITH
natural rubber latex **sterile EO**

- Fine-grain white birch is tasteless and odorless
- Smooth, rounded edges with splinter-free surface and high tensile strength
- Available in non-sterile or individually-packaged sterile
- Available in Standard size (6 × 0.69 Inch) and Junior size (5.5 × 0.63 Inch)
- Can also be used for spreading ointments and stirring/mixing compounds and solutions



Item No.	Description	Quantity
Blade Depressor Tongue Non-Sterile		
W0D3001	Junior 5.5 Inch	1/EA, 500 EA/BX, 10 BX/CA
W0D3002	Senior 6 Inch 500 EA/BX	1/EA, 500 EA/BX, 10 BX/CA
W0D3005	Senior 6 Inch 100 EA/BX	1/EA, 100 EA/BX, 12 BX/CA
Blade Depressor Tongue Sterile		
W0D3004	Junior 5.5 Inch	1/EA, 100 EA/BX, 10 BX/CA
W0D3003	Senior 6 Inch	1/EA, 100 EA/BX, 10 BX/CA



Penlights

NOT MADE WITH natural rubber latex

- Battery-operated, lightweight penlight with pocket clip switch for activation
- Pinpoint lens-end lamp focuses light beam for greater illumination
- Imprinted pupil gauge on side of barrel for easy measurement of pupillary reflex
- ABS casing with aluminum conductor strip
- Battery shelf life of 2 years from date of manufacture
- Measures 5 Inches in length by 0.5 Inches in diameter
- Resealable plastic package includes six pens
- 1-year guarantee



Item No.	Description	Quantity
PEN3051	Penlight Surgical Battery Operated with Pupil Gauge Pocket Clip Switch 5 x 0.5 Inch Disposable	1/EA, 6 EA/PK, 6 PK/BX, 10 BX/CA

REFILLABLE ICE BAGS

Refillable Ice Bags with Clip Closure

NOT MADE WITH natural rubber latex

- Wide mouth opening for easy filling
- Clip is surrounded by outer layers to prevent irritation of the skin and virtually eliminate leakage
- Clip sleeve, not tape, secures clip to bag
- Conforms to joints without restricting movement
- Sturdy ties offer secure donning and maximum maneuverability of ice bag
- Three-layer design of soft spunlace outer layer, absorbent middle layer, and impervious plastic bladder
- Plastic bladder covered by dual outer layers drastically reduces condensation
- Bladder attached to outer covering prevents it from slipping down into the bag



Item No.	Description	Quantity
Bag Ice Refillable Wide Mouth Opening Clip Closure With Ties 3-Layer Design		

IB0613	Small 6 x 13 Inch	1/EA, 25 EA/BX, 4 BX/CA
IB0713U	Universal 7 x 13 Inch	1/EA, 25 EA/BX, 4 BX/CA
IB0813	Large 8 x 13 Inch	1/EA, 25 EA/BX, 4 BX/CA

Instant Cold Packs

NOT MADE WITH natural rubber latex

- For temporary relief of minor pain and swelling caused by sprains, minor burns, contusions, oral pain, and whenever cold therapy is indicated
- A simple squeeze and shake activates the pack
- Safety protocol recommends that the pack be wrapped in towel or provided cloth prior to use
- Available in 3 sizes: Large 6 x 9 Inch, Medium 6 x 6.5 Inch, and Junior 5 x 7 Inch
- Medium size now available in soft flocked packaging material for added patient comfort
- Non-woven covers included with CWP903 and CWP904
- Single-patient use
- Disposable



Item No.	Description	Quantity
Pack Cold Instant Disposable Single-Patient Use		
CWP903	Junior 5 x 7 Inch	1/EA, 48 EA/CA
CWP908	Medium Flocked 6 x 6.5 Inch	1/EA, 24 EA/CA
CWP904	Large with Cover 6 x 9 Inch	1/EA, 24 EA/CA

Reusable Hot or Cold Gel Packs

NOT MADE WITH natural rubber latex

- Heat in the microwave for warm therapy; freeze for cold therapy
- Soft-freeze gel remains flexible to contour to patient's body
- Packs include separate spunlace covers to protect patient skin
- Fast relief for bruises, swelling, muscle spasms/pain, headache, stiff/sore joint, and muscle tension and cramps
- Non-toxic, non-irritating and non-caustic polymer gel
- Single-patient use



Item No.	Description	Quantity
CWP905	Pack Gel Hot-Cold Reusable Single-Patient Use 6 x 10 Inch	1/EA, 24 EA/CA

Instant Warm Packs

made in the USA **NOT MADE WITH natural rubber latex** **diverse supplier**

- Soothing heat therapy for minor muscle aches, joint pain, stiffness and muscle spasms or whenever warm therapy is needed
- A simple squeeze and quick shake activates the pack
- Warms to 115 degrees and provides lasting heat for approximately 20 minutes
- Available in flock covered pack or laminated pack with protective sleeve (safety protocol recommends that all packs be wrapped in towel or provided cloth sleeve prior to use)
- Convenient single use disposable compress
- DEHP-free



Item No.	Description	Quantity
Pack Warm Instant 115° Fahrenheit Maximum Temperature Disposable		
CWP901	With Sleeve 6 x 9 Inch	1/EA, 24 EA/CA
CWP951	Flocked 6 x 9 Inch	1/EA, 24 EA/CA

Heat Therapy Pads

made in the USA NOT MADE WITH natural rubber latex diverse supplier

- Made from soft, thin, light weight polyurethane material for better performance and patient comfort
- Low internal pressures (usually 3-5 psi) on pads result in reducing contact point pressure to help prevent skin reddening and breakdown
- Unique pinwheel channel design promotes excellent water flow and improved temperature therapy
- Superior draping characteristic reduces sliding and shifting during use
- Thinner polyurethane will stay malleable even when used for cooling therapy
- Pads can be used with both hot and cold therapy pumps
- Disposable and reusable pads available



Item No.	Description	Quantity
Pad Heat Therapy Polyurethane		
HTP003	Neck 4 x 19 Inch Single Patient Use	1/EA, 20 EA/CA
HTP017	Small 12 x 17 Inch Single Patient Use	1/EA, 20 EA/CA
HTP020	Medium 14 x 20 Inch Single Patient Use	1/EA, 20 EA/CA
HTP024	Large 18 x 24 Inch Single Patient Use	1/EA, 20 EA/CA
HTPW17	Small 12 x 18.5 Inch Reusable	1/EA, 20 EA/CA
HTPW20	Medium 14 x 20 Inch Reusable	1/EA, 20 EA/CA
HTPW24	Large 18 x 24 Inch Reusable	1/EA, 20 EA/CA

Heat Therapy Pumps

made in the USA NOT MADE WITH natural rubber latex diverse supplier

- Bright and easy to read display, especially in dim patient rooms, providing continuous real-time temperature settings
- Digital controller with proprietary software manages prescribed set-point temperature within 1-degree Fahrenheit
- Three temperature safety limits guard against overheating with high limit switches that can be tested without disassembling; Biomed friendly
- Uses ordinary tap water; distilled water not required
- Limited lifetime warranty with use of MediChoice pads
- Mercury free



Item No.	Description	Quantity
HTP500	Pump Heat Therapy Mercury-Free	1/EA

Enema Bag Kits

NOT MADE WITH natural rubber latex

- Integrated funnel built into enema bag for easy filling
- 60-Inch PVC tubing has a pre-lubricated tapered tip with rounded eyes for added patient comfort and ease of insertion
- Castile soap packet included is moisture proof (0.3 fl. oz.)
- Underpad provided for added convenience and spill containment
- Plastic shut-off clamp for simple operation
- Volume capacity: 1500 cc



Item No.	Description	Quantity
ENBG1500	Kit Enema Bag with 60 Inch Pvc Flush Tube 1 Underpad 1 Castile Soap 1500 Milliliter with Slide Clamp	1/EA, 50 EA/CA

Multipurpose Tube Holders

made in the USA NOT MADE WITH natural rubber latex

- Provides standardized system for securing catheters and medical tubing
- Helps reduce medical tube trauma by limiting tube movement
- Hypoallergenic tape reduces skin irritation caused by other taping methods
- Aids in the reduction of catheter restriction by applying uniform pressure
- Provides comfort to patients with indwelling catheters
- Increases patient comfort and reduces cost of frequent re-taping
- Permits repeated access to tubing without restricting flow
- Can remain in place comfortably for up to one week; ideal for patients with fragile skin
- Easy to adjust fit for individual patient needs
- May be used with or without inner adhesive tab



Item No.	Description	Quantity
MTH1002	Holder Tube Multipurpose Non-Slip Loop Non-Sterile	1/EA, 50/BX, 400/CA

Alcohol and Peroxide

NOT MADE WITH natural rubber latex diverse supplier

- First-aid antiseptic
- Reduces the risk of infection in cuts, scrapes, and burns
- Alcohol may also be used as a rubbing compound
- Peroxide can be used as an oral debriding agent
- Peroxide also aids in the removal of phlegm and mucus associated with oral wounds
- USP Grade



Item No.	Description	Quantity
----------	-------------	----------

Solution First Aid Antiseptic 70-Percent Alcohol Isopropyl		
IA7004	4 Ounce	1/EA, 60 EA/CA
IA7016	16 Ounce	1/EA, 12 EA/CA

Solution First Aid Antiseptic 3-Percent Hydrogen Peroxide		
HP0304	4 Ounce	1/EA, 24 EA/CA
HP0308	8 Ounce	1/EA, 12 EA/CA
HP0316	16 Ounce	1/EA, 12 EA/CA

Complimentary Products

Item No.	Description
LAB2805	Alcohol Dispenser Amber HDPE 250 Milliliter Case Of 12

PILL CRUSHERS & SPLITTERS

Multi-Purpose Pill Crushers

NOT MADE WITH natural rubber latex

- Convenient all-in-one design splits, crushes, and stores pills
- Integrated large compartment can be used as a drinking cup
- Crushes tablets to a fine powder in seconds with a simple twist
- Compact size for easy transport
- Stainless steel blade integrated under protective plastic cap for safety
- Multiple compartments to store pills
- Made of heavy duty plastic; dishwasher safe up to 140 degrees Fahrenheit




Item No.	Description	Quantity
PLC1220	Crusher Pill Multi-Purpose Splitter Plastic	1/EA, 20 EA/CA

Heavy-Duty Pill Crushing System

NOT MADE WITH natural rubber latex

- Heavy-duty steel construction for durability
- Crushes multiple pills at once with minimal effort
- Easy to use, one-hand operation
- Quiet operation will not disturb patient
- Integrated pouch dispenser tray with 50 pouches
- Aids in accurate dosage for patient
- Unique design of crush plates ensures consistent performance
- Dimensions: 11 Inch length x 3 Inch width x 4 Inch height
- Weight: 4 pounds



Item No.	Description	Quantity
PLC5000	System Pill Crushing Heavy-Duty with 50 Pill Pouches Steel	1/EA, 6 EA/CA
PLC5800	Pouch Pill Straight Seam	50/BG, 20 BG/BX, 8 BX/CA

Post Mortem Bags and Kits

NOT MADE WITH natural rubber latex

- Made of strong high grade plastic with bound edges
- Rustproof zipper opens full length of bag for maximum access
- Each bag includes three ID tags and three contamination tags
- Post Mortem kits includes: 3 ID tags, 1 chin strap, 1 cellulose pad, 2-60 Inch ties, and 3-36 Inch ties; personal effects bag included with PM700K and PM800K
- Shroud kits include 54 x 108 Inch white opaque linen plastic sheet, 44 Inch chin strap, 1 underpad, 2-60 Inch ties, 3-36 Inch ties, and 3 white ID tags with string ties
- Thickness: Standard bags 7 Mil, Heavy Duty bags 9 Mil, and Disaster Bag 21.6 Mil (Mil equivalent to 0.0254 mm)
- Conforms to OSHA's universal precautions regarding containment of body fluids and protection from blood borne pathogens



Item No.	Description	Quantity
----------	-------------	----------

Bag Postmortem

PM100B	Disaster Includes 36 x 94 Inch Bag with 6 Nylon Handles 4 Yellow Warning Tags 1 Red Biohazard Tag 3 ID Tags 3 Contamination Tags 1 Zip Lock Biohazard Bag Black Max Weight Capacity 600 Pounds 21.6 Mil	1/EA, 10 EA/CA
PM200B	Infant Includes 18 x 30 Inch Bag 3 ID Tags 3 Contamination Tags Max Weight Capacity 300 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM300B	Pediatric Includes 30 x 48 Inch Bag 3 ID Tags 3 Contamination Tags Max Weight Capacity 300 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM400B	Adult Includes 36 x 90 Inch Bag with Straight Zip 3 ID Tags 3 Contamination Tags Max Weight Capacity 350 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM500B	Adult Includes 36 x 90 Inch Bag with Curved Zip 3 ID Tags 3 Contamination Tags Max Weight Capacity 350 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM800B	XL Includes 54 x 90 Inch XL Bag with Straight Zip 3 ID Tags 3 Contamination Tags 1 Chin Strap 3 60 Inch Ties 3 36 Inch Ties Max Weight Capacity 350 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM900BXL	Adult XL Black with Locking Double Zipper Includes: 47.5 x 99 Inch Bag Locking Double Zipper 3 ID Tags 3 Disease Warning Tags 4 Yellow Contamination Tags 1 Red Bio-Hazard Tag 1 Personal Effects Zip Closure Bag Max Weight Capacity 1000 Lbs Mil-21.6	1/EA, 5 EA/CA
PM900BDZ	Adult with Double Locking Zipper Includes: 36.5 x 92.5 Inch Adult Bag 1 Locking Double Zipper 3 ID Tags 3 Disease Warning Tags 4 Yellow Contamination Tags 1 Red Bio-Hazard Tag 1 Personal Effects Zip Closure Bag Max Weight Capacity 350 Lbs Mil-5.5	1/EA, 10 EA/CA

Kit Postmortem

PM400K	Adult Includes 36 x 90 Inch Bag with Straight Zip 3 ID Tags 3 Contamination Tags 1 Chin Strap 1 Underpad 2 60 Inch Ties 3 36 Inch Ties Max Weight Capacity 350 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM401K	Adult Includes 36 x 90 Inch Bag with Straight Zip 3 ID Tags 3 Contamination Tags 172 x 108 Inch Plastic Sheet 1 Chin Strap 1 Underpad 2 60 Inch Ties 3 36 Inch Ties 1 Polybag With Tie Max Weight Capacity 350 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM500K	Adult Includes 36 x 90 Bag with Curved Zip 3 ID Tags 3 Contamination Tags 1 Chin Strap 1 Underpad 2 60 Inch Ties 3 36 Inch Ties Max Weight Capacity 350 Pounds 7.0 Mil	1/EA, 10 EA/CA
PM600S	Shroud Adult Includes 54 x 108 Inch Shroud Sheet 3 ID Tags 3 Contamination Tags 1 Chin Strap 1 Underpad 2 60 Inch Ties 3 36 Inch Ties	1/EA, 10 EA/CA
PM700K	Adult Heavy Wrap Includes 36 x 90 Inch Heavy Duty Bag 3 ID Tags 3 Contamination Tags 1 Chin Strap 2 Underpads 2 60 Inch Ties 3 36 Inch Ties 1 Set Precaution Labels 1 Clear Polybag with Tie Max Weight Capacity 500 Pounds 9.0 Mil	1/EA, 10 EA/CA
PM800K	Bariatric Includes 72 x 90 Inch XL Bag with Curved Zip 3 ID Tags 3 Contamination Tags 1 Chin Strap 2 Underpads 2 60 Inch Ties 3 36 Inch Ties 1 Polybag with Tie 4 Precaution Labels Max Weight Capacity 500 Pounds 9.0 Mil	1/EA, 10 EA/CA
PM800S	Shroud XL Adult Includes 72 x 108 Inch XL Shroud Sheet 3 ID Tags 3 Contamination Tags 1 Chin Strap 1 Underpad 2 60 Inch Ties 3 36 Inch Ties	1/EA, 10 EA/CA

Please see our full line of Soufflé and Medicine Cups on pages 6-7 in the Admissions Section.

CABLES & LEADS

Cables and Leads

NOT MADE WITH
natural rubber latex

- No special adaptors required, Plug & Play
- Superior patient comfort
- Wide variety of connectors available
- Multi-manufacturer equivalents
- Available in adult, pediatric, infant, and neonate sizes
- Leads available in snap and grabber connections
- One stop cable shopping for most OEM patient monitoring equipment
- CE-Certified product
- Will not void any equipment warranties
- Durable with a 12-month warranty



Supplier Name	Supplier Item #	Description	MediChoice Item #
Bard	153624HM	Philips Temp Cable, 2 Pin Round Plug to Dual Connector, 3.0m	CL15101
BCI	3044	BCI, Adult SpO2 Sensor, 1.1m	CL17101
BCI	3178	Pediatric SoftTip SpO2 sensor, DB9 Connector, 1.1m	CL17102
Drager	3375834	Drager Compatible SpO2 Adaptor Cable, Nellcor, 2.2m	CL9101
Drager	MP00953-02	Drager Compatible NIBP Hose, 2.2m	CL9501
Drager	5956441	Drager compatible 3 Lead Set, Grabber, 0.9m	CL9201
Drager	5956458	Drager compatible 5 Lead Set, Grabber, 0.9m	CL9202
Drager	3368391	Drager Compatible Multimed Cable, 2.5m	CL9401
Edwards Lifesciences	896083021	IBP Cable, Philips Style TruWave PX1800, 4.0m	CL14101
GE	2016560-001	GE CAM 14 Coiled Patient Cable, 1.3m	CL8407
GE	2021141-001	GE Compatible One Piece 3 Lead Set, Grabber	CL8601
GE	2021197-001	IBP Adaptor Cable, 11 Pin to Edwards, 4.0m	CL8801
GE	2056813-010	Clip Set, 10 Pack, Banana Style, 4mm	CL8701
GE	2006644-001	GE compatible Nellcor SpO2 Cable, 3.0m	CL8101
GE	E9004GE	GE Compatible Nellcor SpO2 Cable, 3.0m	CL8101
GE	2017002-001	SpO2 Adaptor Cable, Masimo, 2.2m, *Ref*	CL8104
GE	2017002-003	SpO2 Adaptor Cable, Masimo, 2.2m, *Ref*	CL8104
GE	2027263-002	SpO2 Adaptor Cable, Masimo Ref , 2.2m *REF*	CL8107
GE	2021406-001	GE compatible Nellcor OXI SpO2 Cable, 3.0m	CL8102
GE	2025350-001	GE compatible Nellcor OXI SpO2 Cable, 3.0m	CL8102
GE	2008773-001	Nellcor OXI SpO2 Adapter Cable, 3.0m	CL4104
GE	2027263-001	Masimo compatible SpO2 Adaptor Cable, 2.2m *REF*	CL6101
GE	8841	NIBP Hose, Dual Tube, 2.2m, Screw Connector, Dinamap 8100	CL8506
GE	E9001TJ	GE Compatible NIBP Hose, 2.5m	CL8501
GE	9461-217	GE Compatible NIBP Hose, 2.5m	CL8502
GE	2017008-003	GE Compatible NIBP Hose, 2.5m	CL8502
GE	2017008-001	GE Compatible NIBP Hose, 2.5m	CL8503

Supplier Name	Supplier Item #	Description	MediChoice Item #
GE	414873-001	GE Compatible NIBP Hose, 2.5m	CL8503
GE	107365	GE Compatible NIBP Hose, 2.5m	CL8504
GE	107363	GE Compatible NIBP Hose, 2.5m	CL8505
GE	107368	GE Compatible Neonatal/Infant NIBP Hose, 2.5m	CL8507
GE	88845	GE Compatible NIBP Hose	CL8508
GE	9461-203	GE Compatible NIBP Hose	CL8509
GE	414874-001	GE Compatible NIBP Hose	CL8510
GE	412682-001	GE compatible 3 Lead Set, Grabber, 0.9m	CL8201
GE	411203-001	GE Compatible 3 Lead Set, Grabber, 0.9m	CL8201
GE	412682-002	GE Compatible 3 Lead Set, Grabber, 0.9m	CL8201
GE	412682-001	GE Compatible 3 Lead Set, Grabber, 0.9m	CL8201
GE	411203-001	GE Compatible 3 Lead Set, Snap, 0.9m	CL8214
GE	412681-001	GE compatible 5 Lead Set, Grabber, 0.9m	CL8202
GE	414556-001	GE compatible 5 Lead Set, Grabber, 0.9m	CL8202
GE	E9003CJ	GE compatible 5 Lead Set, Grabber, 0.9m	CL8202
GE	411200-001	GE compatible 5 Lead Set, Snap, 0.9m	CL8203
GE	411202-001	GE compatible 5 Lead Set, Snap, 0.9m	CL8203
GE	E9008KB	GE compatible 5 Lead Set, Snap, 0.9m	CL8203
GE	414556-002	GE compatible 5 Lead Set, Grabber, 1.3m	CL8204
GE	416467-001	GE compatible 5 Lead Set, Grabber, 1.3m	CL8204
GE	394111-005	GE compatible 5 Lead Set, Apex Pro FH, Grabber, 0.9m	CL8205
GE	394111-011	GE compatible 5 Lead Set, Apex Pro FH, Grabber, 0.9m	CL8205
GE	394111-006	GE compatible 5 Lead Set, Apex Pro FH, Snap, 0.9m	CL8206
GE	394111-012	GE compatible 5 Lead Set, Apex Pro FH, Snap, 0.9m	CL8206
GE	414556-005	GE compatible 5 Lead Set, Grabber, Mixed Length	CL8207
GE	412681-005	GE compatible 5 Lead Set, Grabber, Mixed Length	CL8207
GE	411200-005	GE compatible 5 Lead Set, Snap, Mixed Length	CL8208

continued on next page

continued from previous page

Supplier Name	Supplier Item #	Description	MediChoice Item #
GE	411202-005	GE compatible 5 Lead Set, Snap, Mixed Length	CL8208
GE	421930-001	GE compatible 6 Lead Set, Grabber, 0.9m	CL8219
GE	394111-009	6 Lead Set, Grabber, 0.9m, ApexPro FH	CL8220
GE	412682-002	GE Compatible 3 Lead Set, Grabber, 51"	CL8221
GE	38401817	GE compatible 10 Lead Set, Banana, Mixed Length	CL8210
GE	420101-002	GE compatible 10 Lead Set, Banana, Mixed Length	CL8210
GE	416467-003	GE compatible 5 V-Lead Set, Grabber, 0.9m	CL8211
GE	900716-001	GE compatible 3 lead, DIN, Mini-clip, 0.9m	CL8212
GE	TS-F4-GE	GE compatible Clip Style SpO2 Sensor, 3.0m	CL8301
GE	TS-E4-GE	GE compatible Ear Clip Style SpO2 Sensor, 3.0m	CL8302
GE	2027258-001	Masimo Compatible SpO2 Sensor, Adult, 1.1m *REF*	CL5302
GE	420907-001	Nonin Compatible, Adult SpO2 Sensor, 1.1m	CL10101
GE	2017004-001	3 Lead ECG Trunk, Neonatal, AHA, 3.6m, New GE #412994-001	CL8406
GE	2017003-001	GE compatible 5 Lead Trunk Cable, 2.5M	CL8402
GE	412931-001	GE compatible 5 Lead Trunk Cable, 2.5M	CL8402
GE	2017003-001	GE compatible 5 Lead Trunk Cable, 3.6M	CL8402
GE	22341809	GE compatible 10 Lead EKG Trunk Cable, 2.5M	CL8403
GE	2017006-001	GE compatible 12 Lead Trunk Cable, 2.5M	CL8404
GE	2016998-001	Dual Temp Cable, Used with YSI400 Reusable Probes, 0.3m	CL8405
GE	2001925-003	GE CAM 14 Single Universal Lead, 40"	CL8408
GE	416468-001	GE compatible 5 V-Lead Set, Grabber, 1.3m	CL8409
GE	900179-201	CAM 14/AM4/AM5 Electrode Connector Set, Banana, 14/Set	CL8410
GE	None	GE Fixed One Piece ECG cable, 5 leads with grabbers	CL8411
GE	900178-001	CAM 14/AM4/AM5 Electrode Connector Set, Banana, 10/Set	CL8412
GE	900177-001	GE AM4 Lead Set, 10 Lead, Mixed Length	CL8413
GE	420101-001	GE CAM 14 Lead Set, 14 Lead, Mixed Length	CL8414
GE	2021701-001	Temp adaptor cable, YSI400, 3.0m	CL8802
GE	2021196-001	IBP Adaptor Cable, 11 Pin to Medex/Abbott, 4.0m	CL8803
GE/Datex	545317	GE/Datex Compatible 3 Lead Set, Grabber, 0.9m	CL7202
GE/Datex	545303	GE/Datex Compatible 5 Lead ECG Trunk Cable, 3.0m	CL7402
GE/Datex	545318	GE/Datex Compatible 5 Lead Set, Grabber, 0.9m	CL7203
GE/Datex	OXY-C3	GE/Datex compatible SpO2 Adaptor Cable, 2.2m	CL7101
GE/Datex	OXY-SL3	GE/Datex compatible SpO2 Adaptor Cable, 2.2m	CL7102
GE/Datex	OXY-OL3	GE/Datex compatible SpO2 Adaptor Cable, 2.2m	CL7103
GE/Datex	545327	GE/Datex compatible 3 Lead Set, Snap, 0.9m	CL7201
GE/Datex	OXY-F-DB	GE/Datex compatible SpO2 Sensor, 1.1m	CL7301
GE/Datex	OXY-F1-H	GE/Datex compatible SpO2 Sensor, 1.1m	CL7302
GE/Datex	OXY-F-UN	GE/Datex compatible SpO2 Sensor, 1.1m	CL7303
GE/Datex	OXY-E-UN	GE/Datex compatible SpO2 Ear Clip Sensor, 1.1m	CL7304
GE/Datex	OXY-F4-N	GE/Datex compatible SpO2 Sensor, 3.0m	CL7305
GE/Datex	OXY-F4-H	GE/Datex compatible SpO2 Sensor, 3.0m	CL7306
GE/Datex	OXY-E-DB	GE/Datex compatible SpO2 Ear Clip Sensor, 1.1m	CL7307
GE/Datex	545302	GE/Datex compatible 3 Lead ECG Trunk Cable, 3.0m	CL7401
GE/Datex	165641	Temp adaptor cable, YSI400, 3.0m	CL7501
GE/Datex	TS-N3	TruSignal Interconnect Cable with Datex Connector, 2.2m	CL7502
Masimo	LNC-10	Masimo compatible SpO2 Adaptor Cable, 2.2m *REF*	CL6101
Masimo	2281	Philips Compatible SpO2 Adaptor Cable, Masimo, 2.2m *REF*	CL6102
Masimo	LNC MP10	Masimo Compatible SpO2 Adaptor Cable, 2.2m *REF*	CL6102

Supplier Name	Supplier Item #	Description	MediChoice Item #
Masimo	1269	Masimo Compatible SpO2 Sensor, Adult, 1.1m *REF*	CL5301
Masimo	LNOP DCI	Masimo Compatible SpO2 Sensor, Adult, 1.1m *REF*	CL5301
Masimo	LNOP DCIP	Masimo Compatible SpO2 Sensor, Pediatric, 1.1m *REF*	CL6302
Masimo	1863	Masimo Compatible SpO2 Sensor, Adult, 1.1m *REF*	CL5302
Masimo	LNCS DCI	Masimo Compatible SpO2 Sensor, Adult, 1.1m *REF*	CL5302
Masimo	LNCS DCIP	Masimo Compatible SpO2 Sensor, Pediatric, 1.1m *REF*	CL6304
Masimo	LNCS TC-I	Masimo Compatible LNCS TC-I SpO2 Ear Sensor, 1.1m *REF*	CL6305
Masimo	1969	Masimo Compatible Finger Sensor, 3.0m *REF*	CL5304
Medtronic	3006218-02	One Piece 3 Lead Trunk and Leads, Snap	CL12101
Mindray/Datascope	0012-00-0516-02	Datascope Comp SpO2 Adaptor Cable, 2.2m	CL5205
Mindray/Datascope	0600-00-0026-01	Nellcor Compatible Finger Sensor, 1.1m	CL4105
Mindray/Datascope	0600-00-0026-02	Datascope one piece SpO2 Trunk/Sensor, 3.0m	CL4106
Mindray/Datascope	0012-00-0991-01	Datascope compatible 3 Lead Trunk Cable, DIN, 3.0m	CL5402
Mindray/Datascope	0012-00-0622-01	Datascope compatible 3 Lead Set, Snap, 0.9m	CL5206
Mindray/Datascope	0012-00-1464	Nellcor OXI SpO2 Adapter Cable, 3.0m	CL4104
Mindray/Datascope	0012-00-1261-08	Datascope compatible 3 lead Set, Snap, 0.9m	CL5201
Mindray/Datascope	0012-00-1261-02	Datascope compatible 5 lead Set, Snap, 0.9m	CL5202
Mindray/Datascope	0012-00-1262-01	Datascope compatible 5 lead Set, Grabber, 0.9m	CL5209
Mindray/Datascope	0012-00-1262-02	Datascope compatible 5 lead Set, Grabber, 0.9m	CL5209
Mindray/Datascope	0012-00-1262-03	Datascope compatible 5 lead Set, Grabber, 0.9m	CL5209
Mindray/Datascope	0010-30-42727	Datascope compatible 5 lead Set, Grabber, 0.9m	CL5209
Mindray/Datascope	0012-00-1514-02	Datascope compatible 5 lead Telemetry Set, Grabber, 0.9m	CL5203
Mindray/Datascope	0012-00-1503-02	Datascope compatible 5 lead Telemetry Set, Snap, 0.9m	CL5204
Datascope	0012-00-1262-09	Datascope 3-Lead Set, Grabber, 0.9m	CL5207
Mindray/Datascope	0600-00-0047	Masimo Compatible Finger Sensor, 1.1m *REF*	CL5301
Mindray/Datascope	0600-00-0126	Masimo Compatible SpO2 Sensor, Adult, 1.1m *REF*	CL5302
Mindray/Datascope	0012-00-1099-01	Masimo Compatible LNOP SpO2 Cable, 2.2m *REF*	CL5303
Mindray/Datascope	0012-00-1255-01	Datascope compatible 3/4 ECG Trunk, 6 Pin, 3.0m	CL5401
Mindray/Datascope	0683-04-004	Hose Nibp 2.5m Datascope Mindray Compatible	CL5403
Mortara	9325-001-50	Alligator Clip, 10 Pack, Banana Style, 4mm	CL13101
Mortara	9293-032-50	One Piece 10 Lead Trunk and Leads, Banana	CL13102
Nellcor	Nellcor DEC-4	Nellcor compatible SpO2 Extension Cable, 1.1m	CL4101
Nellcor	Nellcor DEC-8	Nellcor compatible SpO2 Extension Cable, 2.2m	CL4102
Nellcor	SCP10	Nellcor compatible sensor cable, non-OXI, 3.0m	CL4103
Nellcor	DOC-10	Nellcor OXI compatible SpO2 Adaptor Cable, 3.0m	CL4104
Nellcor	D-YS	Nellcor Compatible Multisite Sensor, OXI	CL4301
Nellcor	DS-100A	Nellcor OXI Compatible Finger Sensor, 1.1m	CL4302
Nihon Kohden	BR-903PA	Nihon Kohden compatible 3 Lead Set, Grabber, 0.9m	CL3204

continued on next page

CABLES & LEADS

continued from previous page

Supplier Name	Supplier Item #	Description	MediChoice Item #
Nihon Kohden	JL-900P	Nihon Kohden compatible SpO2 Adaptor Cable, 3.0m	CL3101
Nihon Kohden	JL-650P	Nihon Kohden Compatible SpO2 Adapter Cable, 3.0m	CL3102
Nihon Kohden	YN-901P	Nihon Kohden compatible Adult/Pediatric NIBP Hose, 2.5m	CL3501
Nihon Kohden	YN-921P	Nihon Kohden compatible Neonatal NIBP Hose	CL3502
Nihon Kohden	BR-913PA	Nihon Kohden compatible 3 Lead Set, Snap, 0.9m	CL3201
Nihon Kohden	BR-906PA	Nihon Kohden compatible 6 Lead Set, Grabber, 0.9m	CL3202
Nihon Kohden	BR-916PA	Nihon Kohden compatible 6 Lead Set, Snap, 0.9m	CL3203
Nihon Kohden	JC-906PA	Nihon Kohden compatible 6 Lead Trunk Cable, 3.0m	CL3401
Nonin	8000AA	Nonin Compatible, Adult SpO2 Sensor, 1.1m	CL10101
Philips	989803129231	Alligator Clip, 10 Pack, Banana Style, 3mm	CL2701
Philips	989803129161	Philips Compatible 10 Lead Set, Banana, Mixed Length	CL2215
Philips	989803143181	Philips Compatible One Piece 3 Lead Set, Grabber	CL2601
Philips	989803103941	Philips Compatible 5 Lead ECG Trunk Cable, 3.0m	CL2406
Philips	989803104521	Philips Compatible 5 Lead Set, Snap, Mixed Length	CL2212
Philips	M1733A	Philips Compatible One Piece 3 Lead Set, Snap	CL2602
Philips	M1943A	Philips Compatible SpO2 Adapter Cable, 1.1m	CL2108
Philips	M1943AL	Philips Compatible SpO2 Adapter Cable, 3.0m	CL2102
Philips	M1941A	Philips Compatible SpO2 Adapter Cable, 2.2m	CL2103
Philips	M1943NL	Philips Compatible SpO2 Adapter Cable, OXI, 3.0m	CL2104
Philips	M1940A	Philips Compatible SpO2 Adapter Cable, 2.2m	CL2105
Philips	M1900B	Philips Compatible SpO2 Adapter Cable, 2.2m	CL2106
Philips	1821	Philips Compatible SpO2 Adaptor Cable, Masimo MP12, 2.2m *REF*	CL2107
Philips	451261000761	Philips Compatible SpO2 Adaptor Cable, Masimo MP12, 2.2m *REF*	CL2107
Philips	M1020-61100	Philips Compatible SpO2 Adaptor Cable, Masimo MP12, 2.2m *REF*	CL2107
Philips	989803148221	Philips Compatible SpO2 Adaptor Cable, Masimo LNC-MP10, 2.2m *REF*	CL6102
Philips	M1598B	Philips Compatible NIBP Hose, 2.5m	CL2501
Philips	M1599B	Philips Compatible NIBP Hose, 2.5m	CL2501
Philips	M1596B	Philips Compatible Neonatal NIBP Hose, 2.5m	CL2502
Philips	M1597B	Philips Compatible Neonatal NIBP Hose, 2.5m	CL2502
Philips	M1671A	Philips Compatible 3 Lead Set, Grabber, 0.9m	CL2201
Philips	M1673A	Philips Compatible 3 Lead Set, Snap, 0.9m	CL2202
Philips	989803152051	Philips Compatible 5 Lead Telemetry Set, Grabber, 0.9m	CL2203
Philips	989803152071	Philips compatible 5 Lead Telemetry Set, Snap	CL2204
Philips	989803128651	Philips Compatible SpO2 Adapter Cable, 3.0m	CL2102
Philips	M1968A	Philips compatible 5 lead Set, Grabber, Mixed Length	CL2205
Philips	989803125841	5 Lead Set AAMI, grabber, mixed length	CL2205
Philips	M1644A	Philips Compatible 5 Lead Set, Snap, Mixed Length	CL2206
Philips	M1976A	Philips Compatible 5 Lead Limb Set, Grabber, 0.9m	CL2207
Philips	989803125881	Philips compatible 5-Lead Set, Grabber, 0.9m	CL2207
Philips	M1602A	Philips Compatible 5 Lead Limb Set, Snap, 0.9m	CL2208
Philips	M1603A	Philips Compatible 3 Lead Set, Grabber, 0.9m	CL2209
Philips	M1605A	Philips Compatible 3 Lead Set, Snap, 0.9m	CL2210
Philips	M1623A	Philips Compatible 5 Lead Set, Grabber, 0.9m	CL2211

Supplier Name	Supplier Item #	Description	MediChoice Item #
Philips	M1625A	Philips Compatible 5 Lead Set, Snap, Mixed Length	CL2212
Philips	M2593A	Philips Compatible Telemetry 5-Lead Set, Grabber, 0.9m	CL2213
Philips	M2592A	Philips Compatible Telemetry 5-Lead Set, Snap, 0.9m	CL2214
Philips	989803106331	Philips Telemetry 5-Lead Set, Snap, 0.9M with Shield	CL2214
Philips	M2590A	Philips Telemetry 3-Lead Set, Grabber, 0.9m	CL2216
Philips	989803143201	Philips one piece 5 lead cable, Grabber, 3.0m	CL2217
Philips	M1192A	Philips Compatible Pediatric SoftTip SpO2 sensor, 8-Pin, 3.0m	CL2301
Philips	M1196A	Philips Compatible SpO2 Sensor, 3.0m	CL2302
Philips	M1193A	Philips Compatible Neonate Wrap SpO2 Sensor, 3.0m	CL2303
Philips	M1194A	Philips Compatible Ear Clip SpO2 Sensor, 1.1m	CL2304
Philips	M1191A	Philips Compatible Adult SoftTip SpO2 Sensor, 3.0m	CL2305
Philips	M1191AL	Philips Compatible Adult SoftTip SpO2 Sensor, 3.0m	CL2305
Philips	M1191B	Philips Compatible Adult SoftTip SpO2 Sensor, 3.0m	CL2305
Philips	M1190A	Philips Compatible SpO2 Sensor, 3.0m	CL2306
Philips	M1669A	Philips Compatible 3 Lead ECG Trunk Cable, 3.0m	CL2401
Philips	M1668A	Philips Compatible 5 Lead ECG Trunk Cable, 3.0m	CL2402
Philips	M1663A	Philips Compatible 5+5 Lead ECG Trunk Cable, 3.0m	CL2403
Philips	M1949A	Philips Compatible 5+5 Lead ECG Trunk Cable, 3.0m	CL2403
Philips	989803125831	Philips compatible 5+5 ECG Trunk Cable, 12 Pin, 2.7m	CL2403
Philips	M1580A	Philips Compatible 3 Lead ECG Trunk Cable, 3.0m	CL2404
Philips	M1500A	Philips Compatible 3 Lead ECG Trunk Cable, 3.0m	CL2405
Philips	989803103811	Philips Compatible 3 Lead ECG Trunk Cable, 3.0m	CL2405
Philips	M1520A	Philips Compatible 5 Lead ECG Trunk Cable, 3.0m	CL2406
Nova Med	11-CA-401-12	Philips Temp Cable, 2 Pin Round Plug to Dual Connector, 3.0m	CL2407
Spacelabs	700-0030-00	Spacelabs compatible SpO2 Adapter Cable	CL1101
Spacelabs	700-0906-01	Spacelabs compatible SpO2 Adapter Cable, M Tech, 2.2m *REF*	CL1102
Spacelabs	700-0792-00	Spacelabs compatible Nellcor OXI SpO2 cable, 3.0m	CL1103
Spacelabs	700-0007-00	Spacelabs compatible 3 Lead Set, Snap, 0.9m	CL1201
Spacelabs	700-0007-08	Spacelabs compatible 5 Lead Set, Snap, 0.9m	CL1202
Spacelabs	012-0605-00	Spacelabs compatible 5 Lead Set, DIN, Snap, 0.9m	CL1203
Spacelabs	700-0006-35	Spacelabs compatible 3 Lead Set, Grabber, 0.9m	CL1204
Spacelabs	012-0286-00	Spacelabs compatible 5 Lead Set, DIN, Grabber, 0.9m	CL1205
Spacelabs	012-0284-00	Spacelabs compatible 3 Lead Set, Grabber, DIN, Unshielded, 0.9m	CL1206
Spacelabs	700-0008-08	Spacelabs compatible 3 Lead Trunk Cable, 3.0m	CL1401
Spacelabs	700-0008-06	Spacelabs compatible 5 Lead Trunk Cable, 3.0m	CL1402
Welch Allyn	008-0315-00	One Piece 3 Lead Trunk and Leads, Snap	CL16101
Zoll	8000-0025-02	Zoll Compatible One Piece 3 Lead Trunk and Leads, Snap	CL11101

Specimen Traps

NOT MADE WITH natural rubber latex **sterile EO**

- Transparent 40cc or 70cc capacity traps include 1cc graduations, screw-on cap, removable male connector and either a pressure sensitive label or area on the trap to record information
- Transparent 20cc traps include graduations, TOUCH-TROL vacuum breaker, whistle tip suction catheter with Slide-Tex finish and hydrophobic filter
- Convenient second cap for transport to lab
- Manufactured in durable polyurethane with PVC tubing



Item No.	Description	Quantity
Trap Specimen Mucous Polyurethane Sterile		
MST2008	Vacuum Breaker with 8 French Catheter 20 Cubic Centimeter	1/EA, 50 EA/CA
MST2010	Vacuum Breaker with 10 French Catheter 20 Cubic Centimeter	1/EA, 50 EA/CA
MST4000	40 Cubic Centimeter	1/EA, 50 EA/CA
MST7000	70 Cubic Centimeter	1/EA, 50 EA/CA

Suction Catheters

NOT MADE WITH natural rubber latex **sterile EO**

- Flexible PVC suction catheters
- Clear control vent for controlled suctioning
- Available in coiled or straight style
- Beveled whistle tip and offset eyes promote maximum evacuation and helps reduce mucosal trauma, Pediatric sizes (6Fr, 8Fr, and 10Fr) have rounded or "DeLee" tip
- Adult sizes (12Fr, 14Fr, 16Fr, and 18Fr) are 22 inch long and graduated for measuring proper depth insertion
- Pediatric sizes (6Fr, 8Fr, and 10Fr) are 16 inch long, and graduated for measuring proper depth insertion



Item No.	Description	Quantity
Catheter Suction with Control Vent Polyvinyl Chloride Sterile		
Straight		
RSC06D	DeLee Tip 6Fr	1/EA, 50 EA/CA
RSC08D	DeLee Tip 8Fr	1/EA, 50 EA/CA
RSC10D	DeLee Tip 10Fr	1/EA, 50 EA/CA
RSC12	Whistle Tip 12Fr	1/EA, 50 EA/CA
RSC14	Whistle Tip 14Fr	1/EA, 50 EA/CA
RSC16	Whistle Tip 16Fr	1/EA, 50 EA/CA
RSC18	Whistle Tip 18Fr	1/EA, 50 EA/CA
Coiled		
RSC06CD	DeLee Tip 6Fr	1/EA, 50 EA/CA
RSC08CD	DeLee Tip 8Fr	1/EA, 50 EA/CA
RSC10CD	DeLee Tip 10Fr	1/EA, 50 EA/CA
RSC12C	Whistle Tip 12Fr	1/EA, 50 EA/CA
RSC14C	Whistle Tip 14Fr	1/EA, 50 EA/CA

Suction Catheter Kits

NOT MADE WITH natural rubber latex **sterile EO**

- Flexible PVC suction catheters
- Clear control vent for controlled suctioning
- Beveled whistle tip and offset eyes promote maximum evacuation and help reduce mucosal trauma. Pediatric sizes (6Fr, 8Fr and 10Fr) have a rounded or "DeLee" tip
- Adult sizes (12Fr and 14Fr) have whistle tip and are 22 inch long and graduated for measuring proper depth insertion
- Pediatric sizes (6Fr, 8Fr and 10Fr) are 16 inch long, and graduated for measuring proper depth insertion
- Kits include pre-cuffed powder-free vinyl gloves and pop-up solution container in a "walleled" jacket



Item No.	Description	Quantity
Kit Suction Catheter French 1 Pair Pre-Cuffed Powder-Free Vinyl Gloves 1 Pop-Up Solution Container in a Walleled Jacket Sterile		
RSC06KD	DeLee Tip 6Fr	1/EA, 50 EA/CA
RSC08KD	DeLee Tip 8Fr	1/EA, 50 EA/CA
RSC10KD	DeLee Tip 10Fr	1/EA, 50 EA/CA
RSC12K	Whistle Tip 12Fr	1/EA, 50 EA/CA
RSC14K	Whistle Tip 14Fr	1/EA, 50 EA/CA

Suction Catheter Trays

NOT MADE WITH natural rubber latex **sterile EO**

- Flexible PVC suction catheters
- Clear control vent for controlled suctioning
- Beveled whistle tip and offset eyes promote maximum evacuation and help reduce mucosal trauma. Pediatric sizes (6Fr, 8Fr and 10Fr) have rounded or "DeLee" tip
- Adult sizes (12Fr and 14Fr) are 22 inch long and graduated for measuring proper depth insertion
- Pediatric sizes (6Fr, 8Fr and 10Fr) are 16 inch long, and graduated for measuring proper depth insertion
- Includes pre-cuffed powder-free vinyl gloves, 100ML bottle of saline solution in rigid basin



Item No.	Description	Quantity
Tray Suction Catheter Polyvinyl Chloride 1 Pair Pre-cuffed Powder-free Vinyl Gloves 1 Bottle 100 Milliliter Saline Rigid Basin Sterile		
RSC06TSD	DeLee Tip 6Fr	1/EA, 36 EA/CA
RSC08TSD	DeLee Tip 8Fr	1/EA, 36 EA/CA
RSC10TSD	DeLee Tip 10Fr	1/EA, 36 EA/CA
RSC12TS	Beveled Whistle Tip 12Fr	1/EA, 36 EA/CA
RSC14TS	Beveled Whistle Tip 14Fr	1/EA, 36 EA/CA

Nebulizers

NOT MADE WITH
natural rubber latex

- Easy to assemble and setup to use with aerosol masks
- Clear cup for observing medicine
- Flat sturdy cup allows nebulizers to stand up for filling; helps prevent spillage; 5 mL/cc volume
- Graduated volume markings on cup
- Delivers fine, dense mist at any angle up to 90 degrees
- Leak-proof, secure seal



Item No.	Description	Quantity
----------	-------------	----------

Nebulizer Respiratory 7 Foot Kink-Resistant Tubing Leak-Proof Secure Seal

RSP0710	Tee Mouthpiece Flexhose	1/EA, 50 EA/CA
RSP0720	With Adult Mask	1/EA, 50 EA/CA

Oxygen Masks

NOT MADE WITH
natural rubber latex

- Swivel snout accommodates multiple patient positions
- Transparent mask with comfortable anatomical form enables clear observation of patient
- Elastic strap with end caps are adjustable and easy to thread
- DEHP-free



Item No.	Description	Quantity
----------	-------------	----------

Mask Oxygen 7 Foot Kink-Resistant Tubing

RSP1010P	Medium Concentration Pediatric	1/EA, 50 EA/CA
RSP1010A	Medium Concentration Adult Universal Connector	1/EA, 50 EA/CA
RSP1020A	Non-Rebreather with Safety Vent Adult Universal Connector	1/EA, 50 EA/CA
RSP1021A	Non-Rebreather without Safety Vent Adult Universal Connector	1/EA, 50 EA/CA

Trach Mask

NOT MADE WITH
natural rubber latex

- Transparent mask with comfortable anatomical form
- Easy to setup and apply
- Adjustable elastic braid helps secure seal
- Enables clear observation of patient



Item No.	Description	Quantity
----------	-------------	----------

RSP2010	Mask Tracheostomy without Tubing Adult	1/EA, 50 EA/CA
---------	--	----------------

Face Tents

NOT MADE WITH
natural rubber latex

- Transparent mask with comfortable anatomical form enables clear observation of patient
- Elastic strap with end caps are adjustable and easy to thread
- DEHP-free



Item No.	Description	Quantity
----------	-------------	----------

Tent Face Adult

RSP3010	Without Tubing	1/EA, 50 EA/CA
RSP3020	With 7 Foot Aerosol Tubing	1/EA, 50 EA/CA

Aerosol Masks

NOT MADE WITH
natural rubber latex

- Swivel snout accommodates multiple patient positions
- Transparent mask with comfortable anatomical form enables clear observation of patient
- Elastic strap with end caps are adjustable and easy to thread
- DEHP Free



Item No.	Description	Quantity
----------	-------------	----------

Mask Aerosol Vented

RSP5010	Adult	1/EA, 50 EA/CA
RSP5010P	Pediatric	1/EA, 50 EA/CA

CO2 Cannulas

NOT MADE WITH
natural rubber latex

- Provides oxygen while monitoring end-tidal CO2
- Precision sampling of end-tidal CO2
- Soft tubing relieves pressure behind ears thus reducing issues with excoriated skin
- Reduction in labor costs associated with having to apply ointment, creams or lotion behind ears
- Soft anatomic prong design helps ensure patient comfort
- Allows patient to speak, eat and drink



Item No.	Description	Quantity
----------	-------------	----------

Cannula Carbon Dioxide Oxygen Sampling Comfort Soft Plus Adult

RSP6120A	7 Foot Kink Resistant Tubing with Female Luer Universal Connector	1/EA, 25 EA/CA
RSP6110A	7 Foot Kink Resistant Tubing with Male Luer Universal Connector	1/EA, 25 EA/CA
RSP6121A	7 Foot Kink Resistant Tubing with Female Luer 2 Inch Pigtail Universal Connector	1/EA, 25 EA/CA
RSP6122A	10 Foot Kink Resistant Tubing with Female Luer Universal Connector	1/EA, 25 EA/CA
RSP6112A	10 Foot Kink Resistant Tubing with Male Luer Universal Connector	1/EA, 25 EA/CA

Oxygen Cannulas

NOT MADE WITH
natural rubber latex

- Soft, anatomic, over-the-ear design helps ensure patient comfort and fit
- Allows patient to speak, eat and drink
- Soft tubing relieves pressure behind ears thus reducing issues with excoriated skin
- Reduction in labor costs associated with having to apply ointment, creams or lotion behind ears
- Available in a variety of configurations
- DEHP-free



Item No.	Description	Quantity
Cannula Nasal with Kink-Resistant Tubing		
Soft Tip with Kink-Resistant Tubing		
RSP7010	7Ft Tubing Adult	1/EA, 50 EA/CA
RSP7020A	14Ft Tubing Universal Connector Adult	1/EA, 50 EA/CA
RSP7030A	25Ft Tubing Universal Connector Adult	1/EA, 25 EA/CA
Straight Tip with Kink-Resistant Tubing		
RSP8010A	7Ft Tubing Universal Connector Adult	1/EA, 50 EA/CA
RSP8020A	14Ft Tubing Universal Connector Adult	1/EA, 50 EA/CA
Comfort Soft Plus with Kink-Resistant Tubing		
RSP7105A	7Ft Tubing Universal Connector Pediatric	1/EA, 50 EA/CA
RSP7110	7Ft Tubing Adult	1/EA, 50 EA/CA
RSP7110A	7Ft Tubing Universal Connector Adult	1/EA, 50 EA/CA
RSP7120A	14Ft Tubing Universal Connector Adult	1/EA, 50 EA/CA
RSP7130	25Ft Tubing Adult	1/EA, 25 EA/CA

Oxygen Tubing

NOT MADE WITH
natural rubber latex

- Clear tubing for a variety of respiratory applications
- Kink-resistant, six channel, inner lumen design
- DEHP-free



Item No.	Description	Quantity
Tubing Oxygen Kink-Resistant		
RSP9000	7 Foot	1/EA, 50 EA/CA
RSP9010	With Barbed Connector 7 Foot	1/EA, 50 EA/CA
RSP9030	With Barbed Connector 25 Foot	1/EA, 25 EA/CA

Trach Tube Holders

made in the USA **NOT MADE WITH**
natural rubber latex

- Two-piece adjustable neck band provides a secure fit for most patients
- Neck band includes elastic segment for greater comfort and less constriction
- Soft neck band lining aids in patient comfort and keeping neck area drier
- Hook-and-loop closure attaches quickly and easily to most tracheostomy tubes
- Tagless for patient comfort



Item No.	Description	Quantity
Holder Tube Tracheostomy Hook and Loop Closure		
TTH1001	Adult 18 Inch	1/EA, 10 EA/BX
TTH2001	Infant	1/EA, 10 EA/BX
TTH3001	Bariatric 28 Inch	1/EA, 5 EA/BX

Sodium Chloride Solution

made in the USA **NOT MADE WITH**
natural rubber latex **sterile A**

- Convenient twist-off cap allows easy dispensing
- Accurate inverted milliliter graduations
- Safe for IPPB, medication nebulization, and tracheal lavage
- Compatible with closed suction catheters
- No preservatives



Item No.	Description	Quantity
Solution 9 Percent Sodium Chloride Sterile		
UD9003	3 Milliliter	1/EA, 100 EA/BX, 10 BX/CA
UD9005	5 Milliliter	1/EA, 100 EA/BX, 10 BX/CA
UD9015	15 Milliliter	1/EA, 48 EA/BX, 3 BX/CA

Respiratory Set Up Bags

NOT MADE WITH
natural rubber latex

- Durable and economical bag is 12 inch x 15 inch x 4 inch
- Keeps respiratory setups dust-free and organized
- Clear bag makes it easy to see what is inside
- Clinicians can write patient information directly on bag with pen
- Hang from respirator or at patient bedside
- Low-density polyethylene is recyclable



Item No.	Description	Quantity
RSP0010	Bag Respiratory Patient Setup	1/EA, 250 EA/CA

Please see our full line of Tracheostomy Trays on page 114 in the Kits and Trays Section.

MATTRESS OVERLAYS

Synthetic Mattress Overlays

made in the USA NOT MADE WITH natural rubber latex

- Absorbs and transfers moisture to reduce skin breakdown
- Helps circulate air under body and distribute body weight
- Machine washable
- Shed resistant
- 99% Eco-friendly
- Available in 3 sizes




Item No.	Description	Quantity
Overlay Mattress Bed Pad Decubitus Synthetic		
PAD2430	24 × 30 × 1 Inch	1/EA, 10 EA/CA
PAD3040	30 × 40 × 1 Inch	1/EA, 10 EA/CA
PAD3060	30 × 60 × 1 Inch	1/EA, 10 EA/CA

Pressure Relieving Foam Mattress Overlays

made in the USA NOT MADE WITH natural rubber latex

- Distributes body weight evenly for pressure redistribution
- Quality convoluted foam promotes air circulation for enhanced comfort
- Available in multiple sizes and thicknesses
- Superior craftsmanship for lasting value




Item No.	Description	Quantity
Overlay Mattress Standard Eggcrate Foam 33 × 72 Inch		
PAD12000	2 Inch Peak × 0.375 Inch Base	1/EA, 18 EA/CA
PAD12500	2.5 Inch Peak × 0.375 Inch Base	1/EA, 8 EA/CA
PAD12750	2.75 Inch Peak × 0.50 Inch Base	1/EA, 8 EA/CA
PAD13000	3 Inch Peak × 1 Inch Base	1/EA, 8 EA/CA
PAD13130	3 Inch Peak × 0.25 Inch Base	1/EA, 8 EA/CA
PAD14000	4 Inch Peak × 1 Inch Base	1/EA, 6 EA/CA

WHEELCHAIR CUSHIONS

Pressure Relieving Foam Seat Cushions

made in the USA NOT MADE WITH natural rubber latex

- Provides maximum weight distribution for enhanced pressure management
- Foam sections adapt to individual needs based on pressure points in order to optimize pressure relief
- Unique cubic formation provides pressure redistribution while simultaneously minimizing shear forces
- Construction provides cooling channels throughout the cushion
- Available with an incontinent proof cover
- 6-month manufacturer's warranty



Item No.	Description	Quantity
Cushion Seat Pressure Relieving Foam		
WCC1816M	Medium-Density 18 × 16 × 4 Inch 300 Pound Capacity	1/EA, 6 EA/CA
WC1816MS	Medium-Density with Cover 18 × 16 × 4 Inch 300 Pound Capacity	1/EA, 6 EA/CA
WC1816HS	High-Density with Cover 18 × 16 × 3 Inch 250 Pound Capacity	1/EA, 6 EA/CA

WHEELCHAIR CUSHIONS

Air Lift Seat Cushions

made in the USA NOT MADE WITH natural rubber latex

- Provides pressure relief and redistribution for decubitus prevention
- Light, comfortable and economical cushion can be used inside or outside
- Can be used under any temperature conditions
- Easy to inflate or deflate, shipped pre-inflated
- Wipe-clean, antibacterial, heat-sealed vinyl
- Comes with safety straps
- Vinyl meets flammability requirements of CA 117
- 1 year manufacturer's warranty on standard size cushions and hand pump; 6-month warranty on bariatric cushion



Item No.	Description	Quantity
Cushion Wheelchair Seat Air Lift Pre-Inflated		
WCC1717	Vinyl Tan 17 × 17 Inch 350 Pound Capacity	1/EA, 12 EA/CA
WCC1919	Vinyl Tan 19 × 19 Inch 350 Pound Capacity	1/EA, 12 EA/CA
WCC2228	Vinyl Green Bariatric 22 × 28 Inch 700 Pound Capacity	1/EA, 6 EA/CA
Optional Pump		
751580M	Pump Air Handheld Optional For Wheelchair Cushion Universal	1/EA

Air Lock Seat Cushions

made in the USA NOT MADE WITH natural rubber latex

- Provides pressure relief and redistribution for decubitus prevention
- Soft and comfortable air flotation surface
- Wipe-clean, antibacterial, heat-sealed vinyl
- Can be used under any temperature conditions
- Shipped pre-inflated
- Comes with safety straps
- Vinyl meets flammability requirements of CA 117
- 1 year manufacturer's warranty on standard size cushion; 6-month warranty on bariatric size cushion



Item No.	Description	Quantity
Cushion Wheelchair Seat Air Lock Pre-inflated No Valve		
ALSC1919	Vinyl Tan 19 × 19 Inch 350 Pound Capacity	1/EA, 12 EA/CA
ALSC2228	Bariatric Vinyl Green 22 × 28 inch 700 Pound Capacity	1/EA, 6 EA/CA

Please see our full line of Wheelchairs on pages 137-138 in the Mobility Aids Section.

POSITIONING AIDS

Turn and Positioning Bed Systems

made in the USA NOT MADE WITH natural rubber latex

- Safe transferring and positioning system with interchangeable transfer sheet for caregiver
- 30° wedge for optimum coccyx off-loading position provides pressure relief
- Slider sheet for easy patient transfer and reduction of friction and shearing
- Slider sheets connect with bolsters (Velcro) to stabilize system
- Multiple attached handles aid with safer transferring and minimize caregiver back injury
- Wedges are covered with a easy wipe-clean low shear cover
- Choice of two lengths: 16" and 24"
- Can be set up with multiple wedges
- Waterproof barrier for fluid control
- "High grab" upper surface for placement of bed pads (pads not provided)



Item No.	Description	Quantity
Bed System 30 Degree		
TPS6200	Standard 24 Inch Wedge with Blue Nylon Slider Sheet with Four Handles 42 Inch Width x 50 Inch Length 375 Pound Capacity	1/EA, 5 EA/CA
TPS6250	Standard 24 Inch Wedge with White and Tan Mesh Slider Sheet with Continuous Handles 42 Inch Width x 50 Inch Length 250 Pound Capacity	1/EA, 5 EA/CA
TPS6260	Premium Two 16 Inch Wedges with Blue Nylon Slider Sheet with Continuous Handles 38 Inch Width x 50 Inch Length 375 Pound Capacity	1/EA, 5 EA/CA
TPS6261	Patient Transfer Blue Nylon Slider Sheet with Continuous Handles 38 Inch Width x 50 Inch Length 375 Pound Capacity	1/EA, 5 EA/CA

PRESSURE RELIEVING HEEL PROTECTORS

Ultimate Off-Loading Heel Boots

made in the USA NOT MADE WITH natural rubber latex

- Universally sized heel protector designed for the treatment and prevention of pressure ulcers
- Off-loading, zero pressure boot with supersoft poly-fill
- Open heel design effectively eliminates heel pressure by suspending the heel
- Ultra-soft, open weave inner fabric with pillow style cushioning to minimize skin shearing and friction
- Extended height prevents sheets/blankets from contact with toes
- Reduces risk of plantar flexion (foot drop)
- Easy to apply with adjustable straps for custom fit
- Wipe clean outer fabric
- DVT compatible
- 6-month manufacturer's warranty



Item No.	Description	Quantity
Protector Heel Pressure-Relieving Ultimate		
HPB6300	Ultimate Off-Loading	1/EA, 8 EA/CA
HPB6355	Ultimate Off-Loading with Bilateral Anti-Rotation Wedge	1/EA, 8 EA/CA

PRESSURE RELIEVING HEEL PROTECTORS

Premium Heel Protectors

made in the USA NOT MADE WITH natural rubber latex

- Universally sized, open heel design effectively eliminates heel pressure by suspending the heel in the air
- Comfortable heel off-loading at an economical price
- Super soft covering that reduces friction and shear
- Ultra-soft, microfiber lining minimizes friction
- Anti-foot drop construction controls dorsal flexion
- Easy to apply with two versatile straps that can be positioned in multiple configurations
- Removable gel pack to reduce pressure under the ankle, control heat buildup, and redistribute force
- DVT compatible and meets CA 117 flammability standards
- 6-month manufacturer's warranty



Cutaway of Off-Loading Heel (Cutaway not an actual product.)




Item No.	Description	Quantity
Protector Heel Pressure-Relieving Premium		
HPB305A	Flocked Vinyl Fabric with Four Air Holes	
HPB305G	Microfiber with Gel Pack	1/EA, 6/BG, 24/CA
HPB6382	Microfiber without Gel Pack	1/EA, 16/CA
HPB6382W	With Wedge Stabilizer Bi-Lateral Anti-Rotation	1/EA, 8/CA

Standard Heel Protectors

made in the USA NOT MADE WITH natural rubber latex

- Universally sized, standard foam open heel protector off-loads the heel for the prevention and treatment of decubitus heel ulcers
- Promotes faster healing of pressure ulcers and soft foam helps prevent skin tears
- Heel insert covered with low friction fabric for reduction of sweating and friction burns
- Includes auxiliary foam pad for extending the length or height
- Two soft, comfortable cloth tie straps
- Multiple ventilation holes for cooling affect
- Available in either smooth or convoluted foam
- Reduces risk of foot drop
- DVT compatible and meets CA 117 flammability standards
- 6-month manufacturer's warranty



Item No.	Description	Quantity
Protector Heel Pressure-Relieving Standard		
HPB970C	With Convoluted Foam	1/EA, 6 EA/BG, 4 BG/CA
HPB987S	With Smooth Foam	1/EA, 6 EA/BG, 4 BG/CA

Heel Pillows

NOT MADE WITH natural rubber latex

- Provides pressure relief for the feet, heels and ankles
- Breathable cotton/polyester cover provides maximum cool comfort
- Microfiber fill provides maximum heel cushioning for protection against skin breakdown
- One size fits all with convenient hook and loop closure
- Long-lasting and washable



Item No.	Description	Quantity
HP2005	Pillow Heel Universal Pair One Size Fits All Cotton Polyester Cover with Microfiber Fill Hook and Loop Fasteners	1/PR, 6 PR/CA

Calculi Strainer

NOT MADE WITH
natural rubber latex

- Plastic funnel with a fine-mesh filter for trapping kidney stones
- Durable and rinseable polypropylene plastic allows for patient use at home
- Polyester mesh screen filters the smallest stones while allowing unobstructed flow through the strainer
- Fits securely into 24-hour specimen collectors
- Built-in tab for easier sanitary handling



Item No.	Description	Quantity
URST2001	Strainer Calculi Specimen Collection Polypropylene Reusable	1/EA, 25 EA/BG, 4 BG/CA

Irrigation Syringes

NOT MADE WITH
natural rubber latex **sterile EO**

- Designed for one-handed operation and optimal control
- Quick and easy direction of flow, fluid visualization, and ease of measurement
- Luer adapter tip allows for use with smaller-lumen tubing
- Piston syringe double-flanged plunger ensures adequate vacuum
- Bulb syringe is designed to stand on end, making it excellent for use in clinical settings



Item No.	Description	Quantity
BULB011P	Bulb Syringe	1/EA, 50 EA/CA
PIST015P	Piston Syringe	1/EA, 50 EA/CA

Catheter Leg Straps

NOT MADE WITH
natural rubber latex

- Hook-and-loop closure attaches anywhere on strap allowing secure and comfortable fit for most leg sizes
- Plush elastic 2-inch leg band is soft and comfortable against the skin with evenly-distributed compression
- Catheter securement loop is easy to use and holds the catheter firmly in place to prevent bag detachment or patient injury
- Now available with silicone, no-slip backing (CLS1400S) for added patient security
- Universal size fits most patients and catheters




Item No.	Description	Quantity
CLS1400	Universal Catheter Securement	1/EA, 10 EA/BX
CLS1400S	With Silicone No-Slip Backing Universal Catheter Securement	1/EA, 10 EA/BX

Leg Bags

NOT MADE WITH
natural rubber latex **sterile EO**

- Soft vinyl construction for durability and patient comfort
- Universal connections attach to indwelling catheters or external catheters
- Flip-style drainage valve opens and closes with a simple flip, providing visual confirmation the drain is either open or closed
- Pre-attached, woven, adjustable, button-hole leg straps hold bag securely and can be easily cleaned
- Available with or without 18-inch extension tubing for added flexibility
- Offered in medium, 600 mL (19 oz.), and large, 900 mL (32 oz.) capacity bags



Item No.	Description	Quantity
SLB0019	23 Inch Leg Strap Medium 600 Milliliter Vinyl	1/EA, 48 EA/CA
SLB0019T	18 Inch Extension Tubing 23 Inch Leg Strap Medium 600 Milliliter Vinyl	1/EA, 48 EA/CA
SLB0032	30 Inch Leg Strap Large 900 Milliliter Vinyl	1/EA, 48 EA/CA
SLB0032T	18 Inch Extension Tubing 30 Inch Leg Strap Large 900 Milliliter Vinyl	1/EA, 48 EA/CA

Foley Insertion Trays

NOT MADE WITH
natural rubber latex **sterile EO**

- Sterile compact kit provides everything required to insert a Foley catheter with ease
- Allows choice of drain bag and catheter size to meet the needs of specific patients
- Catheter not included



Item No.	Description	Quantity
INTRAY80	10 Cubic Centimeter Pre-Filled Syringe of Sterile Water	1/EA, 20 EA/CA
INTRAY90	30 Cubic Centimeter Pre-Filled Syringe of Sterile Water	1/EA, 20 EA/CA

Urethral Catheterization Trays

sterile EO

- Sterile compact kit provides all necessary components needed for a complete procedure
- Available with either a 14FR plastic catheter or a 14FR red rubber catheter



Item No.	Description	Quantity
UCT870	With 14Fr Plastic Catheter Not Made with Natural Rubber Latex	1/EA, 20 EA/CA
UCT871	With 14Fr Red Rubber Catheter Contains Natural Rubber Latex	1/EA, 20 EA/CA

Pediatric Urine Collector

NOT MADE WITH
natural rubber latex **sterile EO**

- Versatile, unisex design provides a comfortable fit for both infants and toddlers
- Adhesive surface holds securely, removes easily, and folds in half to send a sealed specimen to the lab
- PVC construction with graduations to 200ml



Item No.	Description	Quantity
PUC300	Collector Urine Pediatric 200 Milliliter Polyvinyl Chloride	1/EA, 500 EA/CA



Please see our full line of Kits and Trays on pages 108-117 in the Kits and Trays Section.

ADHESIVE BANDAGES

Sheer Adhesive Bandages

NOT MADE WITH natural rubber latex sterile EO

- Lofty, absorbent, non-adherent pad will not stick to the wound
- Ultrasonic, micro perforations increase breathability to speed up the healing process.
- Durable, long-lasting adhesion and protection
- Tri-fold wrapper for easy opening



Item No.	Description	Quantity
Bandage Adhesive Sheer Non-Adherent Pad Sterile		
ABN2075A	0.75 x 3 Inch	1/EA, 100 EA/BX, 12 BX/CA
ABN2100A	1 x 3 Inch	1/EA, 100 EA/BX, 12 BX/CA
ABN2240	2 x 4 Inch	1/EA, 50 EA/BX, 24 BX/CA
ABN2087	0.875 Inch Round Spot	1/EA, 100 EA/BX, 24 BX/CA

Flexible Fabric Adhesive Bandages

NOT MADE WITH natural rubber latex sterile EO

- Lofty, absorbent, non-adherent pad will not stick to the wound
- Unique cotton/rayon blend is pliable and comfortable
- Breathable to speed up the healing process
- Durable, long-lasting adhesion and protection
- Tri-fold wrapper for easy opening



Item No.	Description	Quantity
Bandage Adhesive Flexible Fabric Non-Adherent Pad Cotton Rayon Sterile		
ABN3075A	0.75 x 3 Inch	1/EA, 100 EA/BX, 12 BX/CA
ABN3100A	1 x 3 Inch	1/EA, 100 EA/BX, 12 BX/CA
ABN3235	2 x 3.5 Inch	1/EA, 50 EA/BX, 12 BX/CA

Soft-Flex Adhesive Bandages

NOT MADE WITH natural rubber latex sterile EO

- Lofty, absorbent, non-adherent pad will not stick to the wound
- Durable cotton/rayon blend
- Island pad seals on all sides
- Breathable to speed up the healing process
- Pliable and comfortable
- Durable, long-lasting adhesion and protection
- Tri-fold wrapper for easy opening



Item No.	Description	Quantity
Bandage Adhesive Soft-Flex Non-Adherent Pad Cotton Rayon Sterile		
ABN5075	0.75 x 3 Inch	1/EA, 100 EA/BX, 12 BX/CA
ABN5100	1 x 3 Inch	1/EA, 100 EA/BX, 12 BX/CA
ABN5087	0.875 Inch Round Spot	1/EA, 100 EA/BX, 12 BX/CA
ABN5125	1.25 Inch Oval	1/EA, 100 EA/BX, 12 BX/CA

ADHESIVE BANDAGES

Children's Designer Adhesive Bandages

NOT MADE WITH natural rubber latex sterile EO

- Lofty, absorbent, non-adherent pad will not stick to the wound
- Ultrasonic, micro perforations increase breathability to speed up the healing process
- Patented Stat Strip® one step opening design
- Durable, long-lasting adhesion and protection
- Two fun designs: Assorted Neon and Swimming Fish Character



Item No.	Description	Quantity
Bandage Adhesive Pediatric Non-Adherent Pad Sterile		
ABN4075D	0.75 x 3 Inch Character	1/EA, 100 EA/BX, 12 BX/CA
ABN4075N	0.75 x 3 Inch Assorted Neon	1/EA, 100 EA/BX, 12 BX/CA

First Aid Ointments

made in the USA

- First aid antibiotic topical ointment
- Helps prevent infection in minor cuts, scrapes and burns
- USP Formulation
- Active ingredient in Bacitracin ointment: Bacitracin 500 I.U./Gram
- Active ingredient in Triple Antibiotic Ointment: Bacitracin 400 I.U./gram, Neomycin 3.5 milligrams, Polymyxin B Sulfate 5,000 I.U./gram
- Bacitracin is comparable to the active ingredient in Baciquent
- Triple Antibiotic Ointment is comparable to the active ingredients in Neosporin
- Inactive ingredient: petrolatum
- Available in individual single use packets: 144 Packets/Box, 12 Boxes/Case



Item No.	Description	Quantity
Ointment Antibiotic First Aid Non-Sterile		
OINT109	Bacitracin Ointment 0.9 Gram Packet	1/EA, 144 EA/BX, 12 BX/CA
OINT110	Bacitracin Ointment 1 Ounce Tube	1/EA, 12 EA/CA
OINT209	Triple Antibiotic Ointment 0.9 Gram Packet	1/EA, 144 EA/BX, 12 BX/CA

Skin Closure Strips

NOT MADE WITH
natural rubber latex **sterile EO**



- Reinforced sterile skin closure strips to secure, close and support small cuts and surgical incisions
- Supports skin following suture or staple removal
- Hypoallergenic adhesive is strong, yet comfortable and breathable
- Ideal for body areas that do not bend or flex and where minimal swelling is expected after procedures
- Coated release liner for easy application
- Holds firmly in place
- Translucent to aid with site visualization
- Removes easily without trauma to the skin
- Provides improved cosmetic results

Item No.	Description	Quantity
Strip Skin Closure Reinforced Adhesive		
SCS1040	0.125 × 3 Inch	5/PK, 50 PK/BX, 4 BX/CA
SCS1041	0.25 × 3 Inch	3/PK, 50 PK/BX, 4 BX/CA
SCS1042	0.25 × 1.5 Inch	6/PK, 50 PK/BX, 4 BX/CA
SCS1046	0.25 × 4 Inch	10/PK, 50 PK/BX, 4 BX/CA
SCS1047	0.5 × 4 Inch	6/PK, 50 PK/BX, 4 BX/CA
SCS1048	1 × 5 Inch	4/PK, 50 PK/BX, 2 BX/CA

Cohesive Bandages

NOT MADE WITH
natural rubber latex **diverse**
supplies **sterile EO**



- Self-adherent wrap for even, consistent compression in retaining dressings, protecting wounds, immobilizing injuries or securing devices to patients
- Adheres only to itself and provides consistent compression to help minimize injured areas
- Slightly tacky feel eliminates slipping or loosening
- Fabric provides breathability and patient comfort
- Easy edge to quickly locate starting point of bandage
- Tears easily from the roll
- Available in sterile and non-sterile packaging
- Color options available; random assortment of colors
- Textured surface enhances cohesive action for more uniform application and security
- Cohesive technology eliminates hardening or casting up

Item No.	Description	Quantity
Bandage Cohesive Self-Adherent Non-Sterile		
CBN1001	1 Inch × 5 Yard Tan	1/EA, 12 EA/BX, 6 BX/CA
CBN1002	2 Inch × 5 Yard Tan	1/EA, 12 EA/BX, 6 BX/CA
CBN1003	3 Inch × 5 Yard Tan	1/EA, 12 EA/BX, 6 BX/CA
CBN1004	4 Inch × 5 Yard Tan	1/EA, 12 EA/BX, 6 BX/CA
CBN1006	6 Inch × 5 Yard Tan	1/EA, 12 EA/BX, 4 BX/CA
CBN102BL	2 Inch × 5 Yard Blue	1/EA, 36 EA/BX, 2 BX/CA
CBN102CP	2 Inch × 5 Yard Color Pack	1/EA, 36 EA/BX, 2 BX/CA
CBN115CP	1.5 Inch × 5 Yard Color Pack	1/EA, 36 EA/BX, 2 BX/CA
Bandage Cohesive Self-Adherent Sterile		
CBN1103	3 Inch × 5 Yard Tan	1/EA, 36 EA/CA
CBN1104	4 Inch × 5 Yard Tan	1/EA, 36 EA/CA
CBN1106	6 Inch × 5 Yard Tan	1/EA, 36 EA/CA

Tubular Net Bandages

NOT MADE WITH
natural rubber latex



- Retainer net assists in wound management, providing easy examination and dressing changes without the use of tape
- Open weave design allows air flow and permits freedom of movement
- Elastic and polyester net bandage cuts easily and stretches to fit various body parts and shapes
- Multi-use
- Ravel-free

Item No.	Description	Quantity
Bandage Tubular Net Stretch Elastic Polyester		
NB5501	Size 1, 25 Yards Fits Fingers, Toes, Wrists	1/RL, 1 RL/BX, 10 BX/CA
NETBND01	Size 1, 25 Yards Fits Fingers, Toes, Wrists	1/RL, 1 RL/BX, 10 BX/CA
NB5502	Size 2, 25 Yards Fits Small Hand, Arm, Leg, Foot	1/RL, 1 RL/BX, 10 BX/CA
NETBND02	Size 2, 25 Yards Fits Small Hand, Arm, Leg, Foot	1/RL, 1 RL/BX, 10 BX/CA
NB5503	Size 3, 25 Yards Fits Medium Hand, Arm, Leg, Foot	1/RL, 1 RL/BX, 10 BX/CA
NETBND03	Size 3, 25 Yards Fits Medium Hand, Arm, Leg, Foot	1/RL, 1 RL/BX, 10 BX/CA
NB5504	Size 4, 25 Yards Fits Large Hand, Arm, Leg, Foot	1/RL, 1 RL/BX, 10 BX/CA
NETBND04	Size 4, 25 Yards Fits Large Hand, Arm, Leg, Foot	1/RL, 1 RL/BX, 10 BX/CA
NB5505	Size 5, 25 Yards Fits Small Head, Shoulder, Thigh	1/RL, 1 RL/BX, 10 BX/CA
NETBND05	Size 5, 25 Yards Fits Small Head, Shoulder, Thigh	1/RL, 1 RL/BX, 10 BX/CA
NB5506	Size 6, 25 Yards Fits Large Head, Shoulder, Thigh	1/RL, 1 RL/BX, 10 BX/CA
NETBND06	Size 6, 25 Yards Fits Large Head, Shoulder, Thigh	1/RL, 1 RL/BX, 10 BX/CA
NB5507	Size 7, 25 Yards Fits Small Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NETBND07	Size 7, 25 Yards Fits Small Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NB5508	Size 8, 25 Yards Fits Medium Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NETBND08	Size 8, 25 Yards Fits Medium Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NB5509	Size 9, 25 Yards Fits Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NETBND09	Size 9, 25 Yards Fits Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NB5510	Size 10, 25 Yards Fits Extra-Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NETBND10	Size 10, 25 Yards Fits Extra-Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NB5511	Size 11, 25 Yards Fits Double Extra-Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NETBND11	Size 11, 25 Yards Fits Double Extra-Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NB5512	Size 12, 25 Yards Fits Triple Extra-Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NETBND12	Size 12, 25 Yards Fits Triple Extra-Large Chest, Back, Perineum, Axilla	1/RL, 1 RL/BX, 10 BX/CA
NETBND22	Size 22, 25 Yards Custom Size, Stretches to 80 Inches in Circumference	1/RL, 1 RL/BX, 10 BX/CA
NET512P	Size 5, Pre-cut, 12 Inch	1/EA, 50 EA/CA

ELASTIC BANDAGES

Woven Elastic Bandages

NOT MADE WITH natural rubber latex diverse supplier sterile EO



- Cool and comfortable cotton bandage absorbs more moisture than synthetic bandages
- Standard economical, premium and premium plus bandages made of medium to heavy gauge yarns
- Controlled stretch for consistent support and compression
- Premium and orthopedic bandages resistant to deterioration from ointments and medications
- Washable orthopedic bandages made of heavy gauge yarns and cotton covered spandex
- Excellent strength and consistent elasticity
- Lays flat and easily conforms to all body contours
- Metal clip closure; double length bandage is self-closure

Item No.	Description	Quantity
Bandage Compression Clip Closure Standard Elastic Woven Yarn Non-Sterile		
REB1002	2 Inch x 4.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB1003	3 Inch x 4.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB1004	4 Inch x 4.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB1006	6 Inch x 4.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
Bandage Compression Clip Closure Premium Elastic Woven Yarn Non-Sterile		
REB2002	2 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2003	3 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2004	4 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2006	6 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
Bandage Compression Clip Closure Premium Elastic Woven Yarn Sterile		
REB2102	2 Inch x 5 Yard	1/EA, 20 EA/CA
REB2103	3 Inch x 5 Yard	1/EA, 20 EA/CA
REB2104	4 Inch x 5 Yard	1/EA, 20 EA/CA
REB2106	6 Inch x 5 Yard	1/EA, 20 EA/CA
Bandage Compression Clip Closure Premium Plus Elastic Woven Yarn Non-Sterile		
REB3002	2 Inch x 5.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB3003	3 Inch x 5.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB3004	4 Inch x 5.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB3006	6 Inch x 5.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
Bandage Compression Clip Closure Premium Plus Elastic Woven Yarn Sterile		
REB3102	2 Inch x 5.5 Yard	1/EA, 10 EA/CA
Bandage Compression Clip Closure Orthopedic Elastic Woven Yarn Non-Sterile		
REB4002	2 Inch x 5 Yard	1/EA, 10 EA/BX, 5 BX/CA
REB4004	4 Inch x 5 Yard	1/EA, 10 EA/BX, 5 BX/CA
REB4006	6 Inch x 5 Yard	1/EA, 10 EA/BX, 5 BX/CA
Bandage Compression Self-Closure Premium Elastic Woven Yarn Non-Sterile		
REB6011	Double Length 6 Inch x 11 Yard	1/EA, 12 EA/BX, 2 BX/CA

ELASTIC BANDAGES

Knit Elastic Bandages

NOT MADE WITH natural rubber latex diverse supplier sterile EO



- Cool and comfortable cotton bandage
- Suitable for general support and compression
- Length is stretched length
- Color-coded for easy size identification
- Ideal secondary bandage
- Knitted construction provides superior absorbency and breathability
- Heavy gauge yarn and cotton covered spandex for superior tensile strength and compression
- Resistant to deterioration from ointments and medications
- Scalloped edges prevent rolling
- Double self-closure allows swift and easy application during surgical procedures
- Stretch enables bandage to better conform to applied surface
- Available with dual self-closure flaps or metal clip closures

Item No.	Description	Quantity
Bandage Compression Dual Self-Closure Standard Elastic Knit Cotton Non-Sterile		
REB2012	2 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2013	3 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2014	4 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2016	6 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
Bandage Compression Clip Closure Standard Elastic Knit Cotton Non-Sterile		
REB2022	2 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2023	3 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2024	4 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB2026	6 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
Bandage Compression Dual Self-Closure Premium Elastic Knit Cotton Non-Sterile		
REB3012	2 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB3013	3 Inch x 5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB3014	4 Inch x 5.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
REB3016	6 Inch x 5.5 Yard	1/EA, 10 EA/BX, 10 BX/CA
Bandage Compression Dual Self-Closure Premium Elastic Knit Cotton Sterile		
REB3112	2 Inch x 5 Yard	1/EA, 20 EA/CA
REB3113	3 Inch x 5 Yard	1/EA, 20 EA/CA
REB3114	4 Inch x 5.5 Yard	1/EA, 20 EA/CA
REB3116	6 Inch x 5.5 Yard	1/EA, 20 EA/CA

COMPRESSION BANDAGE SYSTEMS

3-Layer Reduced Compression Bandage System

NOT MADE WITH natural rubber latex **diverse supplier**




- 3-Layer system for optimum compression and absorbency in the treatment of venous leg ulcers of varying stages
 - Layer 1: Non-woven absorbent padding for comfort and protection of bony joints
 - Layer 2: Breathable cotton/poly crepe conforming bandage anchors and smoothes padded layer while absorbing exudate
 - Layer 3: Self-adhering cohesive bandage secures system together
- Pressure: 17-25 mm/Hg
- Three individually packaged layers in a resealable pouch
- Adhesive tape strips and detailed instruction guide included

Item No.	Description	Quantity
BND3000	Bandage Compression Reduced 3-Layer System	1/EA, 8 EA/CA

3-Layer Compression Bandage System with Trak™ Bandage

diverse supplier



- 3-Layer system for optimum compression and absorbency in the treatment of venous leg ulcers of varying stages and other related conditions
 - Layer 1: Non-woven absorbent padding for comfort and protection of bony joints
 - Layer 2: Patented Trak™ compression bandage provides a calculated guide for clinicians to deliver consistent compression with every application
 - Layer 3: Self-adhering cohesive bandage secures system together
- Pressure 35-40 mm Hg
- Three individually packaged layers in a resealable pouch
- Adhesive tape strips and detailed instruction guide included

Item No.	Description	Quantity
BND3040	Bandage Compression With Trak 3-Layer System Contains Natural Rubber Latex	1/EA, 8 EA/CA

COMPRESSION BANDAGE SYSTEMS

4-Layer Compression Bandage System

NOT MADE WITH natural rubber latex **diverse supplier**



- Four Layer Bandage Kit is designed for optimum compression and absorbency in the treatment of venous leg ulcers and related conditions
 - Layer 1: Non-woven absorbent padding for comfort and protection of bony joints
 - Layer 2: Breathable cotton/poly crepe conforming bandage anchors and smoothes padded layer while absorbing exudate
 - Layer 3: Light compression bandage conforms to any contour while providing half of the compression of the system
 - Layer 4: Self-adhering cohesive bandage secures system together
- Pressure: 35-40 mm/Hg
- Four individually packaged layers in a resealable pouch
- Adhesive tape strips and detailed instruction guide included

Item No.	Description	Quantity
BND9001	Bandage Compression Reduced 4-Layer System	1/EA, 8 EA/CA

Paste Bandages

NOT MADE WITH natural rubber latex **diverse supplier**



- 100% cotton gauze bandage used in the compression treatment of ulcers, edema and strains of the extremities
- Non-raveling gauze evenly impregnated with a non-hardening zinc oxide paste
- Conforms to patient's leg without gaps or wrinkles
- No gelatin or preservatives
- Available with calamine lotion to ease itching

Item No.	Description	Quantity
Bandage Compression Paste Cotton Unna Boot		
99903	3 Inch × 10 Yard	1/EA, 12 EA/CA
99904	4 Inch × 10 Yard	1/EA, 12 EA/CA
Bandage Compression Paste Cotton With Calamine Unna Boot		
99903C	3 Inch × 10 Yard	1/EA, 12 EA/CA
99904C	4 Inch × 10 Yard	1/EA, 12 EA/CA

Tubular Elastic Support Bandages



- Tubular elastic, cotton/spandex bandage for support and compression
- Ideal for effectively managing sprains and strains, soft-tissue injuries, general edema, and post-cast removal
- Easy to apply and stays in place without pins, clips or tape
- Delivers consistent support with a single application
- Convenient alternative to flat wrap bandages



Item No.	Description	Quantity	Item No.	Description	Quantity
Bandage Tubular Elastic Support Natural			Bandage Tubular Elastic Support Cotton Spandex Beige		
TSB5150	Size A 1½ Inch × 11 Yard Infant Feet Arm	1/EA	TSB7250	Size B 2½ Inch × 11 Yard Small Hand Arm	1/EA
TSB5250	Size B 2½ Inch × 11 Yard Small Hand Arm	1/EA	TSB9250	Size B 2½ Inch × 36 Inch Small Hand Arm	1/EA, 12 EA/CA
TSB5275	Size C 2¾ Inch × 11 Yard Medium Arm Small Ankle	1/EA	TSB7275	Size C 2¾ Inch × 11 Yard Medium Arm Small Ankle	1/EA
TSB5300	Size D 3 Inch × 11 Yard Large Arm Medium Ankle Small Knee	1/EA	TSB9275	Size C 2¾ Inch × 36 Inch Medium Arm Small Ankle	1/EA, 12 EA/CA
TSB5350	Size E 3½ Inch × 11 Yard Large Ankle Medium Knee Small Thigh	1/EA	TSB7300	Size D 3 Inch × 11 Yard Large Arm Medium Ankle Small Knee	1/EA
TSB5400	Size F 4 Inch × 11 Yard Large Knee Medium Thigh	1/EA	TSB9300	Size D 3 Inch × 36 Inch Large Arm Medium Ankle Small Knee	1/EA, 12 EA/CA
TSB5450	Size G 4½ Inch × 11 Yard Large Thigh	1/EA	TSB7350	Size E 3½ Inch × 11 Yard Large Ankle Medium Knee Small Thigh	1/EA
TSB5550	Size I 5½ Inch × 11 Yard Large Knee XL Thigh XS Trunk	1/EA	TSB9350	Size E 3½ Inch × 36 Inch Large Ankle Medium Knee Small Thigh	1/EA, 12 EA/CA
TSB5675	Size J 6¾ Inch × 11 Yard Small Trunk	1/EA	TSB7400	Size F 4 Inch × 11 Yard Large Knee Medium Thigh	1/EA
TSB5825	Size K 8¼ Inch × 11 Yard Medium Trunk	1/EA	TSB9400	Size F 4 Inch × 36 Inch Large Knee Medium Thigh	1/EA, 12 EA/CA
TSB5013	Size L 13 Inch × 11 Yard Large Trunk	1/EA	TSB7450	Size G 4½ Inch × 11 Yard Large Thigh	1/EA

MediChoice® Tubular Elastic Support Bandage Compression Sizing Chart				Low Pressure 5-10 mmHG		Medium Pressure 10-20 mmHG		High Pressure 20-30 mmHG	
Product #	Color	Size Code	Size	Low Pressure Circumference		Medium Pressure Circumference		High Pressure Circumference	
				Centimeters	Inches	Centimeters	Inches	Centimeters	Inches
TSB5150	Natural	Size A	1.50"	19-22 cm	7.5-8.7 in	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in
TSB5250	Natural	Size B	2.50"	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in
TSB5275	Natural	Size C	2.75"	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in
TSB5300	Natural	Size D	3.00"	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in	35-40 cm	13.8-17.3 in
TSB5350	Natural	Size E	3.50"	30-35 cm	11.8-13.8 in	35-44 cm	13.8-17.3 in	40-48 cm	17.3-18.9 in
TSB5400	Natural	Size F	4.00"	35-44 cm	13.8-17.3 in	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in
TSB5450	Natural	Size G	4.50"	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in	64-72 cm	25.2-28.3 in
TSB5550	Natural	Size I	5.50"	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in	64-72 cm	25.2-28.3 in
TSB5675	Natural	Size J	6.75"	64-72 cm	25.2-28.3 in	72-86 cm	28.3-33.8 in	86-96 cm	33.8-37.8 in
TSB5825	Natural	Size K	8.25"	72-86 cm	28.3-33.8 in	86-96 cm	33.8-37.8 in	96-108 cm	37.8-42.5 in
TSB5013	Natural	Size L	8.25"	72-86 cm	28.3-33.8 in	86-96 cm	33.8-37.8 in	96-108 cm	37.8-42.5 in
TSB7250	Beige	Size B	2.50"	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in
TSB9250	Beige	Size B	2.50"	22-24 cm	8.7-9.4 in	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in
TSB7275	Beige	Size C	2.75"	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in
TSB9275	Beige	Size C	2.75"	24-27 cm	9.4-10.6 in	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in
TSB7300	Beige	Size D	3.00"	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in	35-40 cm	13.8-17.3 in
TSB9300	Beige	Size D	3.00"	27-30 cm	10.6-11.8 in	30-35 cm	11.8-13.8 in	35-40 cm	13.8-17.3 in
TSB7350	Beige	Size E	3.50"	30-35 cm	11.8-13.8 in	35-44 cm	13.8-17.3 in	40-48 cm	17.3-18.9 in
TSB9350	Beige	Size E	3.50"	30-35 cm	11.8-13.8 in	35-44 cm	13.8-17.3 in	40-48 cm	17.3-18.9 in
TSB7400	Beige	Size F	4.00"	35-44 cm	13.8-17.3 in	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in
TSB9400	Beige	Size F	4.00"	35-44 cm	13.8-17.3 in	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in
TSB7450	Beige	Size G	4.50"	44-58 cm	17.3-22.8 in	58-64 cm	22.8-25.2 in	64-72 cm	25.2-28.3 in

ARM SLEEVES

Arm Sleeves



- Seamless tubular protective sleeves with finished hems
- Gently protects fragile skin against tears, scrapes and abrasions
- Color-coded for easy size identification: Small to X-Large
- Promotes warmth
- 93% durable soft cotton, 7% Spandex



Item No.	Description	Quantity
Sleeve Arm Protective		
ARM0001	Small, Red	1/PR, 12 PR/CA
ARM0002	Medium, Green	1/PR, 12 PR/CA
ARM0003	Large, Yellow	1/PR, 12 PR/CA
ARM0004	Extra-Large, Blue	1/PR, 12 PR/CA

WOUND PREP SOLUTIONS

Saline Wound Flush



- First aid wound cleansing and irrigation wash for gentle and effective cleansing of minor cuts and burns
- Gentle, pressurized stream for precise wound flushing and removal of dirt, debris, and foreign material; 4 to 13 psi
- Drug and preservative free with a two year shelf life
- Unique can and spray design allows for even flow at any angle
- Environmentally friendly, non-aerosol can uses a bag-on-valve technology
- Improved cap design that is snap-secure
- No Sting
- Sterile



Item No.	Description	Quantity
SWF071	Wound Flush Saline Spray 7.1 Ounce Sterile	1/EA, 12 EA/CA

Povidone Iodine Prep Solutions



- Effective for skin and mucous membrane preparation prior to surgical procedures
- Contains 10% weight by volume Povidone iodine, USP solution equivalent to 1% available iodine
- Broad-spectrum antimicrobial with highly effective surfactants
- No sting and non-irritating to skin, wounds and mucosa
- Easy to use flip-top cap with secure butterfly seal to indicate first use and prevent leaks prior to use



Item No.	Description	Quantity
Solution Skin Prep Povidone Iodine Antimicrobial Non-Irritating		
PVP202	2 Ounce	1/EA, 48 EA/CA
PVP204	4 Ounce	1/EA, 36 EA/CA
PVP208	8 Ounce	1/EA, 24 EA/CA
PVP216	16 Ounce	1/EA, 12 EA/CA
PVP232	32 Ounce	1/EA, 12 EA/CA

WOUND PREP SOLUTIONS

Povidone Iodine Scrub Solutions



- Surgical scrub solution for hand scrubbing and topical degerming of the skin
- Contains 7.5% weight by volume Povidone iodine, USP solution equivalent to 0.75% available iodine
- Broad-spectrum antimicrobial with highly effective surfactants
- No sting, non-irritating to skin, wounds and mucosa
- Easy to use flip-top cap with secure butterfly seal to indicate first use and prevent leaks prior to use



Item No.	Description	Quantity
Solution Surgical Scrub Povidone Iodine		
PVP504	4 Ounce	1/EA, 36 EA/CA
PVP516	16 Ounce	1/EA, 12 EA/CA
PVP532	32 Ounce	1/EA, 12 EA/CA

Povidone Iodine Swabsticks Sterile



- Saturated with a 10% USP Povidone-Iodine solution (equivalent to 1% available iodine).
- Used for general antiseptic and germicidal cleaning
- No sting and non-irritating to skin
- Individually packaged, sterile single or triple packs
- Aluminum paper pouch packaging with easy-tear opening



Item No.	Description	Quantity
Skin Prep Antimicrobial Swabsticks Povidone Iodine Sterile		
PVPS100	Single Packs	1/PK, 50 PK/BX, 20 BX/CA
PVPS300	Triple Packs	3/PK, 25 PK/BX, 20 BX/CA

Transparent Film Dressings



- Easy to control, frame delivery system with unique pull tab for easy frame removal
- Clear film protects skin, IV sites and surgical incisions while allowing visibility to assess the wound site
- Conforms easily to areas that are difficult to apply dressings
- Can be left securely in place for several days
- Excellent moisture vapor permeability, yet impermeable to bacteria and liquids
- Removes easily and is non-irritating with repeated use



Item No.	Description	Quantity
Dressing Wound Transparent Film Sterile		
TD1002	2.375 Inch x 2.75 Inch	1/EA, 100 EA/BX, 4 BX/CA
TD1004	4 Inch x 4.5 Inch	1/EA, 50 EA/BX, 4 BX/CA

HYDROCOLLOIDS

Thin Hydrocolloids

NOT MADE WITH natural rubber latex sterile EO



- Promotes moist wound healing by facilitating tissue granulation, epithelization and autolysis while minimizing skin maceration
- Dressings are lower profile and thinner than other hydrocolloids, yet have the capacity to manage fluid in excess of traditional thin hydrocolloids
- Super thin dressing is only 0.1mm in thickness
- Thin Absorbent Skin Adhesive (TASA) technology creates superior performance and fluid handling capacity and breathability
- Completely transparent, allowing for easy wound site visualization; minimizes the need for frequent dressing changes, therefore extended wear time
- Low-friction polyurethane cover helps prevent skin shear
- Skin-friendly dressing formulated to be gentle to periwound skin

Item No.	Description	Quantity
Dressing Wound Hydrocolloid Thin Absorbent		
AWC9501	2 Inch x 2 Inch	1/EA, 10 EA/BX, 44 BX/CA
AWC9502	3.9 Inch x 3.9 Inch	1/EA, 10 EA/BX, 24 BX/CA
AWC9503	5.9 Inch x 5.9 Inch	1/EA, 5 EA/BX, 17 BX/CA

Standard Hydrocolloids

NOT MADE WITH natural rubber latex sterile EO



- Promotes moist wound healing by facilitating tissue granulation, epithelization and autolysis while minimizing skin maceration
- High fluid handling capacity in a thin and conformable dressing provides long-term wear and fewer dressing changes
- Flexible design allows the dressing to easily conform to the body, ultimately enhancing patient comfort
- Rounded corners help to prevent edge lift
- Formulated to be gentle to periwound skin

Item No.	Description	Quantity
Dressing Wound Hydrocolloid Standard		
AWC9511	2 Inch x 2 Inch	1/EA, 10 EA/BX, 44 BX/CA
AWC9512	3.9 Inch x 3.9 Inch	1/EA, 10 EA/BX, 24 BX/CA
AWC9513	5.9 Inch x 5.9 Inch	1/EA, 5 EA/BX, 17 BX/CA
AWC9514	7.9 Inch x 7.9 Inch	1/EA, 5 EA/BX, 15 BX/CA

HYDROCOLLOIDS

Odor Control Hydrocolloids

NOT MADE WITH natural rubber latex sterile EO



- Promotes moist wound healing by facilitating tissue granulation, epithelization and autolysis while minimizing skin maceration
- Thin border with integrated high adhesion formulation ensures a low profile and enhanced dressing securement
- Central island with odor-absorbing technology provides exceptional fluid handling capacity with the added benefit of odor management
- Two layers for maximum absorption (0.6 MM center)
- Can remain in place several days; fewer dressing changes and less disturbance to the wound bed
- Rounded corners help to prevent edge lift
- Flexible design allows the dressings to easily conform to the body, ultimately enhancing patient comfort
- Sacral dressing is designed specifically for use in the highly contoured sacral area

Item No.	Description	Quantity
Dressing Wound Hydrocolloid Bordered With Odor Control Sterile		
AWC9101	2 Inch x 2 Inch	1/EA, 10 EA/BX, 36 EA/CA
AWC9102	3.9 Inch x 3.9 Inch	1/EA, 10 EA/BX, 20 EA/CA
AWC9103	5.9 Inch x 5.9 Inch	1/EA, 5 EA/BX, 17 EA/CA
Dressing Wound Hydrocolloid Bordered Sacral With Odor Control Sterile		
AWC9105	3.3 Inch x 3.9 Inch	1/EA, 10 EA/BX, 20 EA/CA

Wound Cleansers

made in the USA NOT MADE WITH natural rubber latex



- No rinse, non-irritating formulation gently cleanses wounds to aid in dressing changes with minimal disruption to patient comfort
- Non-ionic surfactant gently loosens and removes wound debris, protein, and microorganisms
 - Maintains a moist wound environment
 - Adjustable spray nozzle settings: spray, stream, or off
 - Indications: Stage I-IV Pressure ulcers, stasis ulcers and diabetic ulcers, partial and full thickness wounds, infected and non-infected wounds, post-surgical wounds, first and second-degree burns, cuts, abrasions, and minor skin irritations

Item No.	Description	Quantity
Wound Cleanser Trigger Spray Bottle Sterile		
SWF208	8 Ounce	1/EA, 6 EA/CA
SWF212	12 Ounce	1/EA, 6 EA/CA

Hydrogels

made in the USA NOT MADE WITH natural rubber latex sterile EO



- Protects against abrasions, friction, desiccation, and contamination
- Covers and helps create a moist environment for wound healing
- Clear gel permits visual examination of wound
- Non-irritating and non-cytotoxic
- Indications: Stage I-IV pressure ulcers, partial and full thickness wounds, leg ulcers, surgical wounds, first and second-degree burns, lacerations, abrasions, and skin tears

Item No.	Description	Quantity
Wound Gel Hydrogel Sterile		
WG1001	Hydrogel Bellows Bottle 1 Ounce	1/EA, 10 EA/BX, 3 BX/CA

NUMERIC

300.....	139	77895.....	23	994610.....	93	089RP30.....	10
850-852.....	58	81033-81038.....	66	994615.....	93	089RP40.....	10
2466.....	139	81070-81078.....	66	994620.....	93	089RP40W.....	10
12579.....	95	88100.....	129	994625.....	93	089RP50.....	10
13579.....	94	88102.....	129	994630.....	93	089RP60.....	10
13779.....	94	88201-88206.....	126	994635.....	93	089RP65.....	10
13979.....	94	88300.....	158	994640.....	93	0913MC.....	107
14579.....	94	88303.....	158	994700.....	92	OMC630.....	121
30500.....	15	90410-90415.....	106	994730.....	92	13140M3402-13140M3405.....	85
31000.....	15	90417.....	106	994800.....	93	1314M3202C.....	85
33000.....	15	90601-90603.....	106	994910.....	92	1314M3202LC.....	85
33302-33305.....	7	99903-99904.....	174	994930.....	92	1314M3202LM.....	85
55501-55506.....	142	99955-99956.....	128	995000.....	93	1314M3202LSC.....	85
55511-55516.....	142	101305.....	106	995020.....	93	1314M3202LSH.....	85
55562.....	24	123526.....	15	995100.....	93	1314M3202M.....	85
55582-55583.....	24	123528.....	15	995250.....	93	1314M3202SC.....	85
55610-55611.....	141	123530.....	15	1060300.....	6	1314M3202SH.....	85
55615.....	141	123532.....	15	10080212.....	89	1314M3203C.....	85
55882-55883.....	24	123534.....	15	20040504.....	89	1314M3203LC.....	85
66600-66605.....	17	123536.....	15	30400400.....	89	1314M3203LM.....	85
66607-66611.....	17	123626.....	15	03304L.....	21	1314M3203LSC.....	85
66699.....	17	123628.....	15	03304M.....	21	1314M3203LSH.....	85
69902.....	129	123630.....	15	03304S.....	21	1314M3203M.....	85
69922.....	127	123632.....	15	03304X.....	21	1314M3203SC.....	85
71013.....	106	123634.....	15	03504BL.....	22	1314M3203SH.....	85
71213.....	106	123636.....	15	03504BM.....	22	1314M3204C.....	85
77685.....	23	181101-181104.....	16	03504BS.....	22	1314M3204M.....	85
77695.....	23	206335.....	84	03504BX.....	22	1314M3204SC.....	85
77701-77704.....	20	222362.....	127	03504L.....	22	1314M3204SC.....	85
77706-77707.....	20	223240.....	127	03504M.....	22	1314M3204SH.....	85
77721-77722.....	23	224220-224221.....	127	03504S.....	22	1314M3205C.....	85
77724-77729.....	20, 23	230500.....	15	03504X.....	22	1314M3205M.....	85
77733-77734.....	146	231000.....	15	03504XS.....	22	1314M3205SC.....	85
77740.....	136	700102.....	94	05503L.....	21	1314M3205SH.....	85
77741-77743.....	135	810701-810703.....	66	05503X.....	21	1314M3402L.....	85
77751-77753.....	24	810711-810713.....	66	055401.....	143	1314M3403L.....	85
77755.....	24	810721-810723.....	66	055403.....	143	1314M39011-1314M39012.....	86
77757.....	25	916812A.....	84	055803L.....	21	1314M39014.....	86
77760.....	23	916814A.....	84	055803X.....	21	1314M39018-1314M39019.....	86
77761.....	24	916815A.....	84	0713MC.....	107	1314M39023.....	86
77765.....	24	916820A.....	84	0813MC.....	107	1314M39048.....	86
77772.....	25	916821-916824.....	84	088300FF.....	158	1314M39053.....	86
77776.....	25	994205.....	87	088303FF.....	158	1314M39060.....	86
77779.....	25	994215.....	87	089DP10.....	9	1314M39071.....	86
77780.....	24	994235.....	87	089DP20.....	9	1314M39073-1314M39075.....	86
77785.....	23	994300.....	94	089DP30.....	9	1314M39081.....	86
77795.....	23	994305.....	94	089DP40.....	9	1314M39085.....	86
77857.....	25	994310.....	94	089DP50.....	9	1314M39085.....	86
77860.....	23	994315.....	94	089DP55.....	9	1314M39091-1314M39093.....	86
77872.....	25	994320.....	94	089DP60.....	9	1314M39100.....	86
77873.....	24	994330.....	94	089DP65.....	9	1314M39149.....	86
77876.....	25	994335.....	94	089DP65.....	9	1314M39213.....	86
77879.....	25	994350.....	94	089EC100.....	9	1314M39230.....	86
77880.....	24	994350.....	94	089EC10A.....	9	1314M39249.....	86
77885.....	23	994510.....	95	089LP20.....	10	1314M39268.....	86
		994535.....	95	089LP30.....	10	14-182.....	30
		994600.....	93	089RP10.....	10	2001T-10-2001T-11.....	145
		994605.....	93	089RP20.....	10		

2001T-15.....	145	77779XL.....	25	ABN2087.....	171	BAR206R.....	134
220RADPS.....	137	77780XL.....	24	ABN2100A.....	171	BAR206T.....	134
250ARSD2-250ARSD3.....	87	77785XL.....	23	ABN2240.....	171	BAR3000.....	54
250ARSD4-250ARSD5.....	88	77795XL.....	23	ABN3075A.....	171	BAR300H.....	54
250ARSD8.....	88	77857XL.....	25	ABN3100A.....	171	BAR300HL.....	54
250ARSDW.....	87	77860XL.....	23	ABN3235.....	171	BAR300L.....	54
20ZLUB.....	90	77872XL.....	25	ABN4075D.....	171	BAR706.....	133
3001T-10-3001T-11.....	145	77873XL.....	24	ABN4075N.....	171	BAR806.....	135
3001T-15.....	145	77876XL.....	25	ABN5075.....	171	BAR806W.....	135
305MJ.....	139	77879XL.....	25	ABN5087.....	171	BASN2050.....	48
318RADPE.....	138	77880XL.....	24	ABN5100.....	171	BAST1001-BAST1002.....	47
3304XX.....	21	77895XL.....	23	ABN5125.....	171	BB-001.....	29
33302BL.....	7	81031S.....	66	ACSGK3Y.....	63	BB-034.....	29
3504BXX.....	22	81032S.....	66	ACTGK6DA.....	64	BBT100.....	128
3504BXXX.....	22	916812A.....	84	ACTGK9P.....	61	BBT200.....	128
3504XX.....	22	916814A.....	84	ACTRL6RA.....	62	BBT300.....	128
3504XXX.....	22	916815A.....	84	ACTWK5W.....	65	BD10102.....	140
3GRAMLUB.....	90	916820A.....	84	ADKT101-ADKT102.....	13	BG001-0111.....	29
40ZLUB.....	90	99903C.....	174	ADPTRMC.....	40	BG001-0212.....	29
4R18HAFPE.....	138	99904C.....	174	ALSC1919.....	168	BG002-0111.....	29
4R18RADPS.....	138			ALSC2228.....	168	BG002-0212.....	29
4R20HAFPE.....	138			AM1000-AM1005.....	14	BG008-0111.....	30
5503BL.....	21	A		APP102A.....	157	BK001-0300.....	30
5503BX.....	21	A11010.....	89	APP104A.....	157	BK002-0300.....	30
5503BXX.....	21	A4823PRO.....	99	APS1000.....	158	BL4020.....	126
5503BXXX.....	21	A4824PROR01.....	100	APS3000.....	158	BL4040.....	126
5503BXXXX.....	21	A4830PRO.....	99	ARM0001-ARM0004.....	176	BLDC3001-BLDC3002.....	78
5503XX.....	21	A4830ZRO.....	99	ARW050.....	157	BLDC4002.....	77
5503XXX.....	21	A4832PRO.....	99	ASL2102-ASL2106.....	151	BLDC4052.....	77
5503XXX.....	21	A6036PRO.....	99	ASL2150.....	151	BLDC4502.....	77
5503XXX.....	21	A6043ARO.....	99	AW913BP.....	107	BLDC4552.....	77
5503XXX.....	21	A6043PRO.....	99	AW913UBP.....	107	BLDCFR3002.....	78
55562XL.....	24	A6043PROR01.....	99	AWC9101-AWC9103.....	177	BMA652.....	112
55582XL.....	24	A6043PXO.....	104	AWC9105.....	177	BND3000.....	174
55583XL.....	24	A6043PYO.....	104	AWC9501-AWC9503.....	177	BND3040.....	174
55803XX.....	21	A6043TXOR01.....	104	AWC9511-AWC9514.....	177	BND9001.....	174
55803XXX.....	21	A6043YMORH5.....	104			BP-007.....	29
55882XL.....	24	A6043ZRO.....	99	B		BP-011-BP-012.....	29
55883XL.....	24	A6639PRO.....	99	B6036HRO.....	99	BP-021-BP-023.....	29
5GRAMLUB.....	90	A6639ZRO.....	99	B6043PXO.....	104	BP-090.....	29
751580M.....	168	A7450ARO.....	99	B8046SYO.....	104	BPUM3001.....	95
77695XL.....	23	A7450PRO.....	99	BA12100.....	14	BPUM4001.....	95
77708A.....	134	A7450PXO.....	104	BA1240.....	14	BPUM5001.....	95
77709A.....	134	A7450PYO.....	104	BA1250.....	14	BPUW4002.....	95
77710A.....	134	A7450ZRO.....	99	BA1260.....	14	BPUW5002.....	95
77711A.....	134	A7779.....	20	BA1270.....	14	BRA101-BRA107.....	146
77751XL.....	24	A8046PRO.....	99	BA1280.....	14	BRA232B.....	147
77752XL.....	24	A8046QRO.....	99	BA1290.....	14	BRA232W.....	147
77753XL.....	24	A8046ZRO.....	99	BA50100.....	14	BRA234B.....	147
77755XL.....	24	ABB1050.....	150	BA5040.....	14	BRA234W.....	147
77757XL.....	25	ABB3P45.....	150	BA5050.....	14	BRA236B.....	147
77760L.....	23	ABB3P62.....	150	BA5060.....	14	BRA236W.....	147
77760XL.....	23	ABB3P72.....	150	BA5070.....	14	BRA238B.....	147
77761XL.....	24	ABB4P45.....	150	BA5080.....	14	BRA238W.....	147
77765L.....	24	ABB4P62.....	150	BA5090.....	14	BRA240B.....	147
77765XL.....	24	ABB4P72.....	150	BACBRUSH.....	128	BRA240W.....	147
77765XL.....	24	ABB4P96.....	150	BACOMB.....	128		
77772XL.....	25	ABN2075A.....	171				
77776XL.....	25						

Catalog Number Index

BRA242B.....	147	CEN324V.....	118	CL8601.....	162	D4823RBO.....	100
BRA242W.....	147	CEN424V.....	119	CL8701.....	162	D4823RCO.....	100
BRC2001–BRC2002.....	154	CEN506V.....	119	CL8801.....	162	D4832RBO.....	100
BRC2050.....	154	CL002-0111.....	30	CL8802–CL8803.....	163	D4832RCO.....	100
BRC3001–BRC3002.....	154	CL003-0111.....	30	CL9101.....	162	D5242TWOR01.....	101
BRC4001–BRC4003.....	154	CL10101.....	163, 164	CL9201–CL9202.....	162	D8046HBO.....	102
BRC6001–BRC6004.....	155	CL1101–CL1103.....	164	CL9401.....	162	D8046HCO.....	102
BT1740.....	81	CL11101.....	164	CL9501.....	162	D8647HBO.....	103
BT89001.....	81	CL1201–CL1206.....	164	CLC1001–CLC1002.....	111	D8647HCO.....	103
BTH230–BTH231.....	131	CL12101.....	163	CLC1004.....	111	DC5000–DC5001.....	12
BTH250–BTH251.....	131	CL13101–CL13102.....	163	CLC1007.....	111	DLW510.....	143
BTH310.....	132	CL1401–CL1402.....	164	CLC1007B.....	111	DP0200.....	141
BTH441.....	131	CL14101.....	162	CLC1015.....	111	DP0300.....	141
BTH451.....	132	CL15101.....	162	CLC1017.....	111	DP0350.....	141
BTH810.....	132	CL16101.....	164	CLC1022.....	111	DP0400.....	141
BUKT1101.....	50	CL17101–CL17102.....	162	CLC1025.....	111	DP0500.....	141
BULB011P.....	170	CL2102–CL2108.....	164	CLC1029.....	111	DP0600.....	141
BW100-0212.....	30	CL2201–CL2217.....	164	CLC1029A.....	111	DP0800.....	141
BW3024G.....	49	CL2301–CL2306.....	164	CLC1035.....	111	DS0500.....	6
BW3024S.....	49	CL2401–CL2407.....	164	CLC1058.....	111	DSP01L.....	21
BW3036G.....	49	CL2501–CL2502.....	164	CLC1090A.....	111	DSP01M.....	21
BW3036S.....	49	CL2601–CL2602.....	164	CLC1091.....	111	DSP1XL.....	21
BW3072DG.....	50	CL2701.....	164	CLC1094–CLC1095.....	111	DSP2XL.....	21
BW3072DS.....	50	CL3101–CL3102.....	164	CLC3001.....	111	DST01L.....	21
BW3072G.....	50	CL3201–CL3203.....	164	CLC4004.....	111	DST01M.....	21
BW3072S.....	50	CL3204.....	163	CLMP101.....	129	DST1XL.....	21
BW3504.....	50	CL3401.....	164	CLS1400.....	170	DST2XL.....	21
BW3530ST.....	49	CL3501–CL3502.....	164	CLS1400S.....	170	E	
BW3534ST.....	49	CL4101–CL4103.....	163	CLT100.....	157	EAR301.....	126
BWX01-0412.....	30	CL4104.....	162, 163	CNT808.....	119	EAR303.....	126
BZK100.....	157	CL4105–CL4106.....	163	COMM6001–COMM6002.....	133	EAR402.....	126
		CL4301–CL4302.....	163	COMM7001.....	132	EAR502–EAR503.....	126
		CL5201–CL5207.....	163	COMM7002.....	133	EBD500.....	13
		CL5209.....	163	COMM7007–COMM7009.....	134	ECTC-ARMPAD.....	139
C4832PROAA9.....	100	CL5301–CL5304.....	163	CON1001–CON1002.....	143	ECTCCASTER.....	139
C6043PROAJ5.....	101	CL5401–CL5403.....	163	CON2001–CON2002.....	143	ED-097.....	29
C7658PROAA9.....	99	CL6101.....	162, 163	CP-071–CP-073.....	29	ED-750.....	30
C8046PROAA9.....	102	CL6102.....	163, 164	CP-075.....	29	ED-760.....	30
C8046PROAE7.....	102	CL6302.....	163	CP-076.....	30	EDSP01L.....	21
C8046PROAG6.....	102	CL6304–CL6305.....	163	CRT202.....	41	EDSP01M.....	21
CANE501.....	136	CL7101–CL7103.....	163	CRTP5001.....	81	EDSP1XL.....	21
CANE578.....	136	CL7201–CL7203.....	163	CRTP5002.....	82	EDSP2XL.....	21
CANE618.....	136	CL7301–CL7307.....	163	CRTS5201–CRTS5202.....	82	EFS200–EFS201.....	34
CANE942–CANE943.....	136	CL7401–CL7402.....	163	CRTW5401.....	81	EGC1000.....	12
CART1750.....	48	CL7501–CL7502.....	163	CRTW5406.....	81	EGC1010.....	12
CBN1001–CBN1004.....	172	CL8101–CL8102.....	162	CRTW5412.....	81	EGS100.....	34
CBN1006.....	172	CL8104.....	162	CTN200L.....	91	EHP1001–EHP1005.....	150
CBN102BL.....	172	CL8107.....	162	CTN500M.....	91	ENBG1500.....	160
CBN102CP.....	172	CL8201–CL8207.....	162	CUR1028–CUR1031.....	84	ENZ001.....	34
CBN1103–CBN1104.....	172	CL8208.....	162, 163	CWP901.....	159	ENZ005.....	34
CBN1106.....	172	CL8210–CL8212.....	163	CWP903–CWP905.....	159	ENZ015.....	34
CBN115CP.....	172	CL8214.....	162	CWP908.....	159	ENZ501.....	34
CC3848F.....	26	CL8219–CL8221.....	163	CWP951.....	159	EP-105.....	29
CC3848R.....	26	CL8301–CL8302.....	163			EP1500.....	142
CC4548F.....	26	CL8402–CL8414.....	163			EP2500.....	142
CC4548R.....	26	CL8501–CL8510.....	162			EP3500.....	142
CEN206F.....	118			D4823HBO.....	100		
CEN206V.....	118			D4823HCO.....	100		

C

D

ES25020–ES25021	144	FG6002BB	27	GB-074	31	GLV4004A	97
ES27520–ES27521	144	FG6004–FG6005	27	GB-144	31	GLV5000–GLV5001	97
ES27530–ES27531	144	FM001-0213	30	GD-011	29	GLV5001S	98
ES5004	144	FM002-0213	29	GD-012	31	GLV5002	97
ES5104	144	FM100-0212	33	GD-042	31	GLV5002S	98
ES65020–ES65021	144	FM200-0212	33	GD-111	31	GLV5003	97
ES65024	144	FM300-0212	33	GEY1WA	139	GLV5003S	98
ES65031	144	FMB10PRE	127	GEY1WL	139	GLV5004	97
ES8000	144	FMB10RL	127	GEY20RRUP	139	GLV7000–GLV7004	97
ES8100	144	FMB12PRE	127	GEY218UP	139	GN001-0212	30
EXAM5015	59	FMB12RL	127	GEY22HDUP	139	GN100-0212	30
EXAM5200	58	FMB14PRE	127	GEY24HDUP	139	GP-013–GP-015	31
EXAM5250	58	FMB14RL	127	GEY2AP	139	GP-020	31
EXAM6000	57	FMB16PRE	127	GEY2APD	139	GP-025	31
EXAM6050	57	FMB16RL	127	GEY2CA	139	GP-042	31
EXAM7000	60	FOG1001	143	GEY2HDCA	139	GP-061	31
EXAM7050	60	FP-009–FP-010	33	GEY2HDWA	139	GP-072–GP-073	31
EXAM8000	60	FP-015	33	GEY2WL	139	GP-106–GP-107	31
EXAM8050	60	FP-017	33	GEY2WL1	139	GP-120	31
EXM401	58	FP-021–FP-022	33	GEYALR	139	GP-310	31
EXM4010	59	FP-028	33	GEYATL	139	GP-320	31
EXM4010F	59	FP-038	33	GEYELRA	139	GP-330	31
EXM4011	59	FP-084	33	GEYFRA	139	GP-340	31
EXM4011F	59	FP-119–FP-122	33	GEYFRC	139	GS-004–GS-006	31
EXM401F	58	FPB040	18	GL003-0413	31	GS-060	31
EXM5015F	59	FPB060	18	GL010-0223	31	GS-073	31
EXM5016	59	FRA15L	98	GL010-0233	31	GS-108–GS-109	31
EXM5016F	59	FRA15M	98	GL011-0613	31	GS-122–GS-124	31
EXM701	58	FRA15S	98	GL020-0212	31	GW-055	31
EXM701F	58	FRA15XL	98	GL020-0223	31		
EXTALL	57	FS1100F	20	GL020-0233	31		
	F	FT0040	13	GL022-0613	31	H	
FB-001–FB-005	33	FT0060	13	GL030-0223	31	H0248TKO	103
FB-007	30	FT3500	48	GL030-0233	31	H08GO	100
FB-010	33	FT3500H	48	GL033-0613	31	H32CO	101
FB-021	33	FWD513	107	GL040-0212	31	H32WO	101
FB-037	33		G	GL040-0223	31	H44CO	102
FB-039–FB-041	33	GA10100	15	GL040-0233	31	H44GO	102
FB-063–FB-064	33	GA1050	15	GL042-0212	31	H4823MCOR01	100
FB-090–FB-091	33	GA1060	15	GL044-0613	31	H4830YW	142
FC3085	105	GA1070	15	GL100-0212	31	H4832HBO	100
FCC1001	106	GA1080	15	GL100-0212	31	H4832HCO	100
FD-005	33	GA1090	15	GL100-1214	31	H4832HCOR01	100
FD-010	33	GA50100	15	GL100214	31	H4832MBO	100
FD-062	33	GA5050	15	GL300-0212	31	H4832MCO	100
FD-068	33	GA5060	15	GL300-1214	31	H4832MCOR01	100
FD-111–FD-112	33	GA5070	15	GLIDES	57	H4832MWO	100
FG1001–FG1002	26	GA5080	15	GLV1000–GLV1004	96	H5242YCO	101
FG14921–FG14922	27	GA5090	15	GLV2000–GLV2004	96	H5645TCOR01	101
FG2001–FG2010	27	GAIT541–GAIT542	150	GLV2300–GLV2304	97	H56CO	103
FG3001–FG3002	27	GAIT601–GAIT602	150	GLV3001A	97	H56GO	103
FG4001–FG4002	27	GAIT721–GAIT722	150	GLV3002A	97	H5838EXOR01	101
FG4101	27	GB-001–GB-003	31	GLV3003A	97	H6036EWO	101
FG5001–FG5002	27	GB-027	31	GLV3004A	97	H6036HBO	101
FG5104	27	GB-059	31	GLV4001A	97	H6036HCO	101
FG6001–FG6002	27	GB-067–GB-068	31	GLV4002A	97	H6036LBO	101
				GLV4003A	97	H6036MBO	101

Catalog Number Index

H6036MBOR01.....	101	H8647SCO.....	103	HTP024.....	160	IVS1030.....	117
H6036MCO.....	101	H8647SGO.....	103	HTP500.....	160	IVS1040.....	117
H6036SGO.....	101	H8647TCOR01.....	103	HTPW17.....	160	IVS1041.....	116
H6036TBO.....	101	H8648WSO.....	103	HTPW20.....	160	IVS1043.....	115, 116
H6036TCO.....	101	HAMP1000.....	45	HTPW24.....	160	IVS1052C.....	115
H6036TCOR01.....	101	HAMP8001–HAMP8003.....	44			IVS1056.....	117
H6036TWO.....	101	HAMP8004.....	45	IA7004.....	160	IVS1062A.....	117
H6045HCOR01.....	101	HAMP9000.....	44	IA7016.....	160	IVS1063.....	116
H6045TWO.....	101	HCGCON2.....	124	IB0613.....	159	IVS1069.....	116
H6045TWOR01.....	101	HDS100.....	13	IB0713U.....	159	IVS1071B.....	116
H6639ACOR01.....	101	HP-010.....	29	IB0813.....	159	IVS1072.....	116
H6639EWO.....	101	HPO304.....	160	IDM-25.....	124	IVS1080.....	116
H6639HBO.....	101	HPO308.....	160	IDSTREP-25.....	124	IVS1082.....	117
H6639HCO.....	101	HPO316.....	160	IHP-10.....	124	IVS1084.....	116
H6639HCOR01.....	101	HP2005.....	169	INC1001B.....	109	IVS1086–IVS1087.....	116
H6639LCO.....	101	HPB305A.....	169	INC2010.....	119	IVS1090.....	115
H6639MBO.....	101	HPB305G.....	169	INL001.....	34	IVS1093.....	116
H6639MCO.....	101	HPB6300.....	169	INL005.....	34	IVS1095–IVS1096.....	116
H6639SCO.....	101	HPB6355.....	169	INL015.....	34	IVS2000.....	116
H6639SGO.....	101	HPB6382.....	169	INL101.....	34	IVS2003.....	116
H6639TBO.....	101	HPB6382W.....	169	INL110.....	34	IVS2005.....	116
H6639TCO.....	101	HPB970C.....	169	INSF02NS.....	141	IVS2009.....	116
H6639TWO.....	101	HPB987S.....	169	INSF1001–INSF1004.....	141	IVS2010.....	117
H7258MCO.....	103	HR1000.....	26	INSF2003–INSF2004.....	141	IVS2011.....	116
H7658ACO.....	103	HR3100.....	26	INTB9001–INTB9008.....	50	IVS2016.....	116
H7658EWO.....	103	HR4100.....	26	INTRAY80.....	170	IVS2030–IVS2031.....	116
H7658HBO.....	103	HR6100.....	26	INTRAY90.....	170	IVS2032–IVS2034.....	117
H7658HCO.....	103	HRA8L.....	98	IPA01G.....	118	IVS2036.....	116
H7658MCO.....	103	HRA8M.....	98	ISP001.....	34	IVS2039.....	117
H7658PCO.....	103	HRA8S.....	98	ISP080.....	34	IVS2040–IVS2042.....	116
H7658SCO.....	103	HRA8XL.....	98	IST1001.....	26	IVS2045–IVS2046.....	116
H7658SGO.....	103	HS0610.....	26	ITRAY130.....	114	IVS2048.....	116
H7658SKO.....	103	HS0615.....	26	IVOS50.....	46	IVS2049–IVS2050.....	117
H7658TBO.....	103	HS0813.....	26	IVBH56.....	46	IVS3003.....	116
H7658TCO.....	103	HS1215.....	26	IVBW48.....	46	IVS3005.....	116
H7658TWO.....	103	HSDVER8.....	123	IVBW485.....	46	IVS8003.....	117
H8046ACO.....	102	HSER100.....	123	IVCH44.....	46	IVSH43.....	46
H8046EWO.....	102	HSHR18L.....	98	IVCR1010–IVCR1013.....	45	IVSH57.....	46
H8046EWOR01.....	102	HSHR18M.....	98	IVCR24H.....	45	IVSS1050–IVSS1051.....	45
H8046HBO.....	102	HSHR18S.....	98	IVCR44H.....	45	IVST46.....	46
H8046HCO.....	102	HSHR18XL.....	98	IVOH55.....	46	IVUC41.....	46
H8046HCOR01.....	102	HSP50ER.....	123	IVPF42.....	46	IVWB51B.....	47
H8046MCO.....	102	HSSP-25.....	123	IVS1001.....	117	IVWB51S.....	47
H8046QKO.....	102	HSSPCAS.....	123	IVS1003–IVS1005.....	116	IVWB52B.....	47
H8046SCO.....	102	HSSPCON.....	123	IVS1006 IVS1007.....	115	IVWB52S.....	47
H8046SGO.....	102	HSSPENV.....	123	IVS1009.....	116		
H8046TBO.....	102	HSTAT50.....	123	IVS1011.....	115	J	
H8046TBOR01.....	102	HTB030.....	143	IVS1013.....	115	JBASICP1618.....	139
H8046TCO.....	102	HTB060.....	143	IVS1017.....	116	JBASICP1816.....	139
H8046TWO.....	102	HTBW30.....	143	IVS1020.....	115	JBASICP2016.....	139
H8046XCO.....	102	HTBW60.....	143	IVS1020A.....	116	JBASICP2018.....	139
H8647EWO.....	103	HTP003.....	160	IVS1021.....	116	JBASICP2020.....	139
H8647HCO.....	103	HTP017.....	160	IVS1023.....	116	JBASICP2218.....	139
H8647PCO.....	103	HTP020.....	160	IVS1024A.....	116		
				IVS1026.....	115	K	
				IVS1029.....	115	KS010-0412.....	33
						KS010-0420.....	33
						KS010-0434.....	33

Catalog Number Index

MC1080S.....	68	MC2451BL.....	6	MC3082WH.....	6	MC5060.....	122
MC1100HD.....	92	MC2451CR.....	6	MC3083BL.....	6	MC5070.....	122
MC11101–MC11104.....	91	MC2451WH.....	6	MC3100D.....	68	MC5080.....	127
MC11105–MC11108.....	16	MC2451YL.....	6	MC3100S.....	68	MC5090.....	127
MC11109–MC11111.....	91	MC2451YL12.....	6	MC311006.....	121	MC5400.....	90
MC11115.....	16	MC2453WH.....	6	MC3154.....	63	MC54003.....	67
MC11154–MC11155.....	69	MC24564.....	69	MC31623L.....	67	MC5453N.....	68
MC11206.....	91	MC24601.....	121	MC31624L.....	67	MC5454N.....	68
MC11208.....	91	MC2469L.....	67	MC31625N.....	67	MC566DH.....	75
MC11516.....	69	MC24744D.....	68	MC31626N.....	67	MC566DHS.....	75
MC11518.....	69	MC24763L.....	68	MC31663N.....	67	MC566DS.....	75
MC11570.....	69	MC24763N.....	68	MC31664N.....	67	MC566GH.....	75
MC11642.....	69	MC2481BL.....	6	MC31663L.....	67	MC566GHS.....	75
MC12261.....	69	MC2481CR.....	6	MC31664L.....	67	MC566GS.....	75
MC1310P.....	90	MC2481WH.....	6	MC31680N.....	67	MC577GH.....	75
MC1410Z.....	90	MC2483BL.....	6	MC31689L.....	67	MC577GHS.....	75
MC14182.....	30	MC2485A.....	92	MC32574.....	69	MC577GS.....	75
MC147214.....	98	MC250SD2–MC250SD3.....	87	MC33584P.....	69	MC587W.....	76
MC147216.....	98	MC250SD4–MC250SD5.....	88	MC3500D.....	68	MC587WH.....	76
MC147218.....	98	MC250SD8.....	88	MC3500S.....	68	MC587WHS.....	76
MC147224.....	98	MC250SDW.....	87	MC3501BD.....	68	MC587WS.....	76
MC147226.....	98	MC25754.....	67	MC3526.....	74	MC588W.....	76
MC147228.....	98	MC25755N.....	68	MC3526WA.....	74	MC5993.....	125
MC150SD2–MC150SD3.....	87	MC26744D.....	68	MC3547.....	70	MC600STB.....	125
MC150SD4–MC150SD5.....	88	MC26784D.....	68	MC354P.....	72	MC6025–MC6026.....	65
MC150SD8.....	88	MC27743L.....	67	MC3566.....	72	MC620.....	64
MC150SDW.....	87	MC27743N.....	67	MC356P.....	73	MC633.....	61
MC17103N.....	68	MC27744T.....	68	MC3596.....	70	MC6411.....	62
MC1710H.....	90	MC28583.....	69	MC359P.....	72	MC645.....	63
MC17461.....	91	MC28584BB.....	69	MC3600D.....	68	MC651.....	62
MC18156.....	16	MC28721.....	69	MC3600DG.....	68	MC6651.....	64
MC18302–MC18303.....	91	MC28764D.....	67	MC3600S.....	68	MC6701D.....	68
MC18309–MC18310.....	91	MC2893L.....	67	MC36806C.....	69	MC72244L.....	68
MC18416.....	92	MC2893N.....	67	MC36806P.....	69	MC745E.....	63
MC186218.....	91	MC2950–MC2951.....	121	MC36846C.....	69	MC75112.....	67
MC18802.....	91	MC2950F.....	121	MC3858.....	70	MC78463.....	68
MC18804.....	91	MC2951F.....	121	MC385P.....	72	MC80366C.....	69
MC18838.....	16	MC2955.....	121	MC3868.....	72	MC80831.....	16
MC18839.....	91	MC3.....	74	MC386P.....	73	MC80929.....	91
MC18842.....	16	MC300006.....	121	MC40477.....	92	MC8206L.....	68
MC18844.....	16	MC3002D.....	68	MC4305.....	26	MC8206N.....	68
MC18846.....	16	MC3003D.....	68	MC4510.....	125	MC84366C.....	69
MC18860.....	91	MC3051BL.....	6	MC4807.....	26	MC84366P.....	69
MC18964.....	16	MC3051CR.....	6	MC4900D.....	68	MC8813L.....	68
MC21921L.....	67	MC3051DB.....	6	MC4900S.....	68	MC8813N.....	68
MC21921N.....	67	MC3051GR.....	6	MC4951.....	121	MC8823.....	68
MC21931L.....	67	MC3051OR.....	6	MC4951F.....	121	MC8833L.....	68
MC21931N.....	67	MC3051RD.....	6	MC4955F.....	121	MC8833N.....	68
MC21941L.....	67	MC3051YL.....	6	MC4993.....	125	MC901–MC903.....	18
MC21941N.....	67	MC3053BL.....	6	MC4996.....	125	MC910516.....	16
MC22221–MC22222.....	121	MC30705.....	69	MC4999.....	125	MC910536.....	16
MC222215.....	121	MC30764D.....	68	MC500.....	18	MC911243.....	91
MC22401.....	121	MC3080D.....	68	MC5001D.....	68	MC913142.....	91
MC22501.....	121	MC3080S.....	68	MC5001S.....	68	MC914243.....	91
MC24301.....	121	MC3081BL.....	6	MC5020.....	122	MC916021.....	91
MC24501.....	121	MC3081CR.....	6	MC5030.....	122	MC916143.....	91
MC245015.....	121	MC3081GR.....	6	MC5040.....	122		

Catalog Number Index

OVI21002.....	107	PC9080.....	8	PM401K.....	161	PVP232.....	176
OVI21004.....	107	PCL2336.....	105	PM500B.....	161	PVP504.....	176
OVI21009.....	107	PCL3036.....	105	PM500K.....	161	PVP516.....	176
OVI21014.....	107	PD100-0212.....	30	PM600S.....	161	PVP532.....	176
OVI21024.....	107	PD200-0212.....	30	PM700K.....	161	PVPS100.....	176
OVI21033.....	105	PD300-0212.....	30	PM800B.....	161	PVPS300.....	176
OVI21042.....	107	PEDT1401.....	79	PM800K.....	161	PW2004.....	129
	P	PEDT1402.....	80	PM800S.....	161	PW2085.....	129
P06110.....	147	PEDTFR1401.....	79	PM900BDZ.....	161		R
P06120.....	147	PEDTFR1402.....	80	PM900BXL.....	161	RA00FB.....	74
P06290.....	147	PEN3051.....	159	PP4200-PP4201.....	11	RA00SB.....	74
P06310.....	147	PEW050.....	20	PP4200D.....	11	RA1010.....	74
P06320.....	147	PEW100.....	20	PPB032C.....	10	RA1020.....	74
P06431-P06433.....	147	PEW250.....	20	PPB1010.....	10	RA1030.....	74
P06545.....	147	PEZ001.....	34	PPB1080.....	129	RA1040.....	74
P06610.....	148	PEZ005.....	34	PPB1110.....	10	RA1050.....	74
PAD12000.....	168	PEZ015.....	34	PPB1210.....	10	RA1060.....	74
PAD12500.....	168	PEZ101.....	34	PPB2010.....	11	RA1070.....	74
PAD12750.....	168	PEZ110.....	34	PPB2020.....	11	RA1080.....	74
PAD13000.....	168	PEZ301.....	34	PPB2110.....	11	RA1090.....	74
PAD13130.....	168	PG4-0220.....	22	PPB2120.....	11	RA2010.....	74
PAD14000.....	168	PG4-1001.....	22	PPB2122.....	11	RA2020.....	74
PAD16815.....	148	PG4-1073.....	22	PPB3010.....	11	RA2121.....	74
PAD2430.....	168	PG4-1092.....	22	PPB3110.....	11	RA3010.....	74
PAD3040.....	168	PG4-1211.....	22	PPB4010.....	11	RA3020.....	74
PAD3060.....	168	PG4-1331.....	22	PPB4110.....	11	RA3030.....	74
PBBX100.....	12	PG4-1563.....	22	PPB6100.....	11	RA4020.....	74
PC1015.....	126	PG4-1573.....	22	PS004-0111.....	31	RALTBRTC.....	74
PC102FF.....	8	PG4-1680.....	22	PS004-0212.....	31	RC001-0111.....	30
PC1035.....	126	PG4-2001.....	22	PS007-0212.....	29	RC010111.....	30
PC104FF.....	8	PG4-2073.....	22	PS012-0512.....	29	REB1002-REB1004.....	173
PC2000.....	7	PG4-2121.....	22	PS014-0212.....	29	REB1006.....	173
PC2515.....	8	PG4-2161.....	22	PS015-0212.....	29	REB2002-REB2004.....	173
PC2600.....	9	PG4-2321.....	22	PS016-0111.....	29	REB2006.....	173
PC2610.....	9	PG4-2331.....	22	PS016-0212.....	29	REB2012-REB2014.....	173
PC28202-PC28203.....	8	PG4-2333.....	22	PS016-0223.....	29	REB2016.....	173
PC29022-PC29023.....	8	PG4-2561.....	22	PS016-0233.....	29	REB2022-REB2024.....	173
PC3000.....	7	PG4-3363.....	22	PS016-0512.....	29	REB2026.....	173
PC4020.....	8	PG4-4021.....	23	PS017-0212.....	29	REB2102-REB2104.....	173
PC4035.....	8	PG4-4051.....	22	PS020212.....	29	REB2106.....	173
PC5015.....	8	PG4-4073.....	23	PS030212.....	29	REB3002-REB3004.....	173
PC6002.....	8	PG4-4241.....	23	PS060512.....	29	REB3006.....	173
PC6004.....	8	PG4-4331.....	23	PSU2336.....	105	REB3012-REB3014.....	173
PC6006.....	122	PG4-5033.....	22	PSU3036.....	105	REB3016.....	173
PC6006OR.....	122	PG4-5063.....	22	PTB100.....	7	REB3102.....	173
PC7000.....	7	PG4-5073.....	22	PTRAY126.....	114	REB3112-REB3114.....	173
PC7181.....	12	PG4-6403.....	22	PTRAY136.....	114	REB3116.....	173
PC7190-PC7191.....	12	PIST015P.....	170	PTRAY146.....	114	REB4002.....	173
PC7200-PC7201.....	12	PLC1220.....	161	PUC300.....	170	REB4004.....	173
PC7202C1.....	12	PLC5000.....	161	PUW1605.....	106	REB4006.....	173
PC7306.....	9	PLC5800.....	161	PUW1615.....	106	REB6011.....	173
PC7385.....	9	PM100B.....	161	PUW1625.....	106	RECP2101.....	51
PC7500.....	8	PM200B.....	161	PVP202.....	176	RECP3101.....	51
PC8020.....	7	PM300B.....	161	PVP204.....	176	RECP4101.....	51
PC8122.....	126	PM400B.....	161	PVP208.....	176	RG-006.....	29
PC9020.....	8	PM400K.....	161	PVP216.....	176	RG-008.....	33

RP-011.....	29	RSP8010A.....	167	SLD2000.....	140	SSKB-SSK107.....	83
RP-030.....	29	RSP8020A.....	167	SLD21200.....	140	SSKB-STWK503.....	83
RS001-0212.....	33	RSP9000.....	167	SLD21500.....	140	SSKB-SV107L.....	83
RS001-0223.....	33	RSP9010.....	167	SLD22000.....	140	SSKB-SW106G2.....	83
RS001-0233.....	33	RSP9030.....	167	SLD23000.....	140	SSKB-SW90PG2.....	83
RS001-0512.....	33	RTS100.....	131	SLD2500.....	140	SSKB02.....	83
RS002-0212.....	33	RTS200.....	131	SLD3000.....	140	SSKB107.....	83
RS004-0512.....	33	RTS300.....	131	SLD500.....	140	SSKB503.....	83
RS005-0512.....	33	RZR2000.....	9	SLP01UB.....	17	SSM-SSM3.....	83
RSC06CD.....	165			SLP01UY.....	17	SSM-SSWM3.....	83
RSC06D.....	165	S		SLP1001-SLP1006.....	17	SSM001.....	83
RSC06KD.....	165	S7646VNO.....	103	SLP12DG.....	17	SSMW042.....	83
RSC06TSD.....	165	SB0609.....	125	SLP12DO.....	17	SSR1001-SSR1003.....	113
RSC08CD.....	165	SB0610T.....	125	SLP12DR.....	17	ST1001-ST1004.....	140
RSC08D.....	165	SB0810.....	125	SLP12DY.....	17	ST1600F.....	54
RSC08KD.....	165	SC2400.....	124	SLP14DG.....	17	ST1700.....	55
RSC08TSD.....	165	SCAL101.....	44	SLP14DO.....	17	ST1700B.....	55
RSC10CD.....	165	SCAL1284.....	37	SLP14DP.....	17	ST1700BF.....	55
RSC10D.....	165	SCAL1374.....	38	SLP14DR.....	17	ST1700F.....	55
RSC10KD.....	165	SCAL1601-SCAL1602.....	38	SLP14DY.....	17	ST1750.....	55
RSC10TSD.....	165	SCAL1603.....	37	SLP16DG.....	17	ST1750B.....	55
RSC12.....	165	SCAL1604.....	39	SLP16DO.....	17	ST1750BF.....	55
RSC12C.....	165	SCAL1607.....	39	SLP16DR.....	17	ST1750F.....	55
RSC12K.....	165	SCAL1610.....	39	SLP16DY.....	17	ST1800.....	55
RSC12TS.....	165	SCAL1644.....	41	SLP18DG.....	17	ST1800B.....	55
RSC14.....	165	SCAL1676.....	39	SLP18DO.....	17	ST1800BF.....	55
RSC14C.....	165	SCAL1763.....	39	SLP18DR.....	17	ST1800F.....	55
RSC14K.....	165	SCAL17MC.....	41	SLP18DY.....	17	ST1850.....	55
RSC14TS.....	165	SCAL1813.....	41	SLP2001-SLP2007.....	17	ST1850B.....	55
RSC16.....	165	SCAL1856.....	40	SLP4004-SLP4005.....	17	ST1850BF.....	55
RSC18.....	165	SCAL1882.....	40	SLP6DY.....	17	ST1850F.....	55
RSP0010.....	167	SCAL20MC.....	41	SLP7DG.....	17	ST1900.....	56
RSP0710.....	166	SCAL31MC.....	40	SLP7DO.....	17	ST1900B.....	56
RSP0720.....	166	SCAL38CH.....	42	SLP7DR.....	17	ST2000HF.....	56
RSP1010A.....	166	SCAL38MC.....	42	SLP7DY.....	17	ST2000HFA.....	56
RSP1010P.....	166	SCAL41MC.....	43	SLP8DG.....	17	ST2000HFB.....	56
RSP1020A.....	166	SCAL42MC.....	43	SLP8DO.....	17	ST2001-ST2004.....	140
RSP1021A.....	166	SCAL45MC.....	43	SLP8DR.....	17	ST3001-ST3004.....	140
RSP2010.....	166	SCAL46MC.....	43	SLP8DY.....	17	STB651.....	112
RSP3010.....	166	SCAL49MC.....	43	SLPH3001-SLPH3005.....	18	STEP1750-STEP1752.....	48
RSP3020.....	166	SCAL81DH.....	40	SN003-0212.....	33	STK0225.....	148
RSP5010.....	166	SCAL88DH.....	42	SP-003.....	30	STK1025.....	148
RSP5010P.....	166	SCPL10-SCPL11.....	141	SP-573.....	29	STK1225.....	148
RSP6110A.....	166	SCPL15.....	141	SPC400.....	125	STK1248.....	148
RSP6112A.....	166	SCR800.....	125	SPC405.....	125	STK1325.....	148
RSP6120A.....	166	SCS1040-SCS1042.....	172	SPC4050R.....	125	STK1425.....	148
RSP6121A.....	166	SCS1046-SCS1048.....	172	SPCT1001-SPCT1002.....	145	STK1525.....	148
RSP6122A.....	166	SD-010.....	30	SPG1001-SPG1003.....	87	STK1625.....	148
RSP7010.....	167	SIDE1225-SIDE1226.....	60	SRCASTERS.....	57	STK1825.....	148
RSP7020A.....	167	SJ0502-SJ0506.....	21	SS1215.....	26	STK2336.....	148
RSP7030A.....	167	SJ0602-SJ0606.....	21	SS1218.....	26	STK2436.....	148
RSP7105A.....	167	SLB0019.....	170	SS3509.....	26	STK2448.....	148
RSP7110.....	167	SLB0019T.....	170	SS5210.....	26	STK2460.....	148
RSP7110A.....	167	SLB0032.....	170	SS5214.....	26	STK2648.....	148
RSP7120A.....	167	SLB0032T.....	170	SS7513.....	26	STK2660.....	148
RSP7130.....	167	SLD1200.....	140	SSKB-S106G2.....	83	STK2672.....	148
		SLD14000.....	140				
		SLD1500.....	140				

Carts, Cylinder Patient Use.....	65	Clamps, Universal.....	46
Carts, Deluxe Emergency-Steel.....	62	Classic Isolation Carts - Aluminum.....	63
Carts, Emergency Pediatric Color-Coded, Aluminum.....	61	Classic Isolation Carts - Steel.....	62
Carts, Emergency-Aluminum.....	62	Classic Overnighter Sleeper Chairs.....	74
Carts, Emergency-Steel.....	61	Classic Post Op Shoes.....	156
Carts, Mini-Line Isolation-Steel.....	63	Cleaner, Enzyme Pre-Soak.....	34
Carts, Monitor-Steel.....	64	Cleaners, Cart.....	35
Carts, Multi-Purpose Phlebotomy-Hybrid Aluminum.....	65	Cleaners, Electrosurgical Tip.....	140
Carts, Plastic and Aluminum Utility.....	82	Cleaners, Neutral PH Instrument.....	35
Carts, Plastic Utility.....	81	Cleaners, Specialty Products.....	34
Carts, Stainless Utility.....	82	Cleaning Kits, Microscope Optical.....	121
Carts, Utility with Wire Baskets.....	81	Cleaning Wipes, Baby.....	126
Cases, Eye Class.....	12	Cleanser, Denture Tablets.....	9
Cases, Pillow.....	10	Cleansers, Foaming Skin.....	105
Cast Carts - Steel.....	65	Cleansers, Skin.....	177
Catheter Leg Bags.....	170	Cleansing Cloths, Flushable.....	106
Catheter Leg Straps.....	170	Cleansing Cremes.....	107
Catheter Procedure Carts.....	67	Cleansing Towelettes.....	157
Catheter, Suction Kits.....	165	Clear Polyethylene Telemetry Pouches.....	18
Catheter, Suction Trays.....	165	Clinical Care Recliners, Deluxe.....	75
Catheterization Trays, Urethral.....	170	Clinical Care Recliners, Extra-Wide.....	76
Catheters, Suction.....	165	Clinical Care Recliners, Standard.....	75
Cauteries, Cordless.....	143	Clippers.....	9
Central Line Dressing Trays.....	111	Closure Strips, Skin.....	172
Centrifuges (Bench-Top).....	118	CO2 Cannulas.....	166
Centrifuges (High-Capacity).....	118	Coats, Premium Disposable Lab.....	22
Centrifuges (Horizontal-Spin).....	119	Coats, Standard Disposable Lab.....	21
Centrifuges (Microhematocrit).....	119	Cohesive Bandages.....	172
Cervical Collars, Foam.....	151	Cold Packs, Instant.....	159
Cervical Collars, Two Piece.....	151	Cold Packs, Perineal.....	129
Chair Scales, Convertible.....	42	Cold Packs, Refillable Ice Bags with Clip Closure.....	159
Chair Scales, Digital.....	39	Cold Reusable Gel Packs.....	159
Chairs, Bariatric Blood Drawwing with Flip Arms.....	77	Collars, Foam Cervical.....	151
Chairs, Bariatric Recline Treatment.....	73	Collars, Two Piece Cervical.....	151
Chairs, Bariatric Treatment.....	72	Collection Containers, 24 Hour.....	124
Chairs, Bath.....	131	Collection Containers, Tamper Proof.....	125
Chairs, Classic Overnighter Sleeper.....	74	Collection Containers, Urinalysis.....	125
Chairs, Drop Arm Commodes.....	133	Collectors, Pediatric Urine.....	170
Chairs, Extra-Wide Blood Drawing with Flip Arms.....	78	Color Coded Guedel Airways.....	15
Chairs, Low Access Powered Podiatry.....	59	Column Carts.....	66
Chairs, Low Access Powered Procedure.....	59	Column Scales, High Capacity Digital with BMI.....	39
Chairs, Power Recline Treatment.....	72	Column Scales, High Capacity Digital with BMI and Digital Height Rod.....	39
Chairs, PVC Shower Commode.....	132	Comb, Baby.....	128
Chairs, Shower Transport.....	132	Comb/Hairbrush/Toothbrush.....	7
Chairs, Side.....	60	Combination Bariatric Scales.....	38
Chairs, Standard Blood Drawing.....	78	Comfort Pack.....	13
Chairs, Standard Powered Podiatry.....	59	Commode Chairs, PVC Shower.....	132
Chairs, Standard Powered Procedure.....	60	Commode Pail Sets.....	134
Chairs, Treatment.....	70	Commodes, Bariatric Drop Arm.....	133
Children's Designer Adhesive Bandages.....	171	Commodes, Drop Arm.....	133
Chrome Disposable Stethoscopes.....	92	Commodes, Folding Steel.....	132
Chrome Double Post Mayo Stands.....	47	Commodes, Heavy-Duty Aluminum.....	133
Chrome Hampers, 18 Round.....	45	Commodes, Painted Steel.....	133
Chrome Hampers, 18 Round Tilt-Lid.....	44	Commodes, Toilet Seat Raised.....	131
Chrome Hampers, Adjustable Folding.....	44	Compression Arm Sleeves.....	176
Chrome Hampers, Linen with Optional Labels.....	44	Compression Bandage System with Trak™ Bandage, 3-Layer.....	174
Chrome Hampers, Square Tilt-Lid.....	44	Compression Bandage System, 3-Layer Reduced.....	174
Chrome IV Stands.....	45	Compression Bandage System, 4-Layer.....	174
Chrome Linen Hampers with Optional Labels.....	44	Computer Mice.....	83
Chrome Step Stools.....	48	Computer Workstation Accessories, Mobile.....	61
Clamps, Towel.....	27	Computer Workstations, Mobile.....	61
Clamps, Umbilical Cord.....	129	Connectors, Tubing.....	143

Convertible Chair Scales.....	42	Diagnostic and Outpatient Lamps.....	54
CoolSpot and Diagnostic Lamps.....	53	Diaper and Organ Scales.....	40
Cord Clamp Cutters, Umbilical.....	129	Diapers, Adult.....	106
Cord Clamps, Umbilical.....	129	Differential Counters.....	119
Cordless Cauteries.....	143	Digital BMI Scales with Handrail.....	40
Correct Site Non-Sterile Skin Markers.....	146	Digital Chair Scales.....	39
Cot Replacement Pads.....	69	Digital Column Scales, High Capacity with BMI and Digital Height Rod.....	39
Cotton Balls.....	91	Digital Column Scales, High Capacity with BMI.....	39
Cotton Blend Bias Cut Stockinettes.....	149	Digital Floor Scales with Cable Remote Control.....	40
Cotton Stockinette, Tubular.....	148	Digital Portable Scales.....	41
Cotton Tipped Applicators, General Purpose.....	89	Digital Scales with BMI, Digital Height Rod and Wireless.....	37
Counter Top Blanket Warmers.....	49	Digital Scales with Handrail.....	41
Counters, Differential.....	119	Digital Take Home Scales.....	44
Counters, Sponge.....	145	Digital Thermometers Kits and Sheaths.....	84
Cover Glasses and Slides.....	121	Dispensers, Alcohol.....	118
Cover Gowns.....	23	Dispensers, Apparel.....	30
Cover Gowns, Light-Weight.....	25	Dispensers, Emesis Bags.....	11
Cover Gowns, Mid-Weight.....	25	Dispensers, Face Mask.....	33
Cover Slips and Slides.....	121	Dispensers, Speciality and Accessories.....	30
Coveralls, Premium Disposable.....	21	Disposable Aprons.....	20
Coveralls, Standard Disposable.....	21	Disposable Blood Pressure Cuffs.....	86
Covers, Cart.....	26	Disposable Boot Covers.....	23
Covers, Light Handle.....	144	Disposable Bouffant Caps.....	20
Cradle Positioners.....	148	Disposable Cover Gowns, Light-Weight.....	25
Cream Denture Adhesive.....	9	Disposable Cover Gowns, Mid-Weight.....	25
Cremes, Cleansing.....	107	Disposable Exam Capes.....	16
Crib Replacement Pads.....	69	Disposable Exam Gowns.....	16
Crusher, Multi-Purpose Pill.....	161	Disposable Exam Shorts.....	16
Crushing System, Pill Heavy-Duty.....	161	Disposable Gowns, AAMI Rated Material.....	25
Crutches, Aluminum.....	134	Disposable Gowns, Fluid-Resistant.....	23
Crutches, Bariatric Steel.....	134	Disposable Gowns, Polyethylene.....	24
Cuff Accessories, Blood Pressure.....	85	Disposable Isolation Gowns, AAMI Rated.....	24
Cuffs, Blood Pressure (Reusable).....	84	Disposable Isolation Gowns, Non-AAMI.....	24
Cuffs, Blood Pressure (Single-Patient Use).....	86	Disposable Lab Coats, Premium.....	22
Cups, Denture.....	12	Disposable Lab Coats, Standard.....	21
Cups, Medicine.....	7	Disposable OR Towels.....	142
Curettes, Ear.....	84	Disposable Pillows.....	9
Cushions, Air Lift Seat.....	168	Disposable Premium Coveralls.....	21
Cushions, Air Lock Seat.....	168	Disposable Procedure Masks.....	22
Cushions, Pressure Relieving Foam Seat.....	168	Disposable Razors.....	9
Custom Minor Procedure Kits and Trays.....	108	Disposable Safety Scalpels.....	141
Cutters, Umbilical Cord Clamp.....	129	Disposable Scrubs.....	21
Cylinder Carts - Bulk Transport.....	65	Disposable Shoe Covers.....	23
Cylinder Carts - Patient Use.....	65	Disposable Shoe Covers, SMS.....	23
		Disposable Staff Jackets.....	21
		Disposable Standard Coveralls.....	21
		Disposable Stethoscopes, Chrome.....	92
		Disposable Stethoscopes, Plastic.....	92
		Disposable Surgeon Caps.....	20
		Disposable Towels.....	91
		Disposable Underwear, Protective.....	106
		Disposable Utensils.....	7
		Disposable Washcloths.....	106
		Doppler Systems (ABI Comprehensive).....	88
		Doppler Systems (ABI).....	88
		Doppler Systems (Obstetrical).....	87
		Doppler Systems (Vascular).....	88
		Double Door Blanket Warmers.....	50
		Drainage and Incision Trays.....	109
		Drape/Stretcher Sheets.....	91
		Drapes, Magnetic Instrument.....	145
		Drapes, Procedure, Sterile.....	92

D

Daily Moisturizing Lotion.....	8
Decanting Devices.....	140
Defibrillation Electrodes.....	90
Defogger.....	143
Deluxe Clinical Care Recliners.....	75
Deluxe Emergency Carts - Steel.....	62
Deluxe Exam Stools.....	55
Deluxe Overnighter Sleeper Chairs.....	74
Denture Cups.....	12
Denture Products.....	9
Deodorant, Roll-On.....	8
Depressors, Tongue.....	158
Dermal Biopsy Punches.....	141
Devices, Decanting.....	140
Diagnostic and CoolSpot Lamps.....	53

Dressing Trays, Central Line.....	111	Exam Gloves, Latex.....	96
Dressing, Hydrogels.....	177	Exam Gloves, Next Generation Nitrile.....	97
Dressing, Odor Control Hydrocolloids.....	177	Exam Gloves, Stretch Vinyl.....	97
Dressing, Standard Hydrocolloids.....	177	Exam Gloves, Stretchy Nitrile.....	97
Dressing, Thin Hydrocolloids.....	177	Exam Gloves, Vinyl.....	97
Dressings, Transparent Film.....	176	Exam Gowns.....	16
Drinking Straws, Flexible.....	6	Exam Lights, Diagnostic and Outpatient.....	54
Drop Arm Commodes.....	133	Exam Lights, Economy Floor.....	51
Drop Arm Commodes, Bariatric.....	133	Exam Lights, Multi-Purpose Mobile.....	51
Dual Tip Sterile Skin Markers.....	145	Exam Lights, Procedure.....	51
Dual-Head Professional Stethoscopes, Adult/Pediatric.....	94	Exam Lights, Procedure - K Arm - Halogen.....	52
Dual-Head Stethoscopes, Cardiology Professional.....	95	Exam Lights, Procedure - K Arm - LED.....	52
Dual-Head Stethoscopes, Infant.....	93	Exam lights, Procedure LED.....	53
Dual-Head Stethoscopes, Light-Weight.....	94	Exam Procedure Drapes, Sterile.....	92
Dual-Head Stethoscopes, Neo-Natal.....	93	Exam Shorts, Disposable.....	16
Dual-Head Stethoscopes, Pediatric.....	93	Exam Stool Accessories.....	57
Dual-Head Stethoscopes, Sprague-Rappaport.....	95	Exam Stools, Black Base.....	57
Dual-Release Folding Walkers.....	135	Exam Stools, Deluxe.....	55
E			
Ear and Ulcer Syringes.....	126	Exam Stools, Economy.....	54
Ear Curettes.....	84	Exam Stools, Hands Free.....	56
Ear Plugs.....	13	Exam Stools, Nylon Base.....	56
Earbuds and Headsets.....	13	Exam Stools, Stainless Contoured Seat.....	57
ECG Chart Paper.....	92	Exam Stools, Stainless Padded Seat.....	57
ECG Defibrillation Electrodes.....	90	Exam Stools, Value.....	55
ECG Leads and Cables.....	162	Exam Table Accessories.....	58
ECG Resting Tab Electrodes.....	90	Exam Table Paper.....	91
Economy Exam Stools.....	54	Exam Tables, Bariatric.....	60
Economy Floor Exam Lights.....	51	Exam Tables, Low Access Powered.....	58
Economy Microscopes.....	120	Exam Tables, Reversible Manual.....	58
Economy Overbed Tables.....	81	Extension Cables.....	162
Eggcrate.....	148	Extra-Wide Blood Drawing Chairs with Flip-Arms.....	78
EKG Chart Paper.....	92	Extra-Wide Clinical Care Recliners.....	76
EKG Defibrillation Electrodes.....	90	Eye Glass Cases.....	12
EKG Resting Tab Electrodes.....	90	Eye Level Scales.....	38
Elastic Bandages, Woven.....	173	Eye Protectors.....	142
Elastic Bandages, Knit.....	173	Eyewear, Protective.....	20
Elastic Support Bandages, Tubular.....	175	F	
Elbow and Heel Protectors.....	150	Face Mask Dispensers.....	33
Elbow and Heel, Foam Protectors.....	147	Face Mask, Oxygen.....	166
Electrocardiogram Chart Paper.....	92	Face Tents.....	166
Electrodes (Defibrillation).....	90	Faceshields.....	20
Electrodes (Resting Tab).....	90	Facial Tissues.....	13
Electrosurgical Pencils/Tips.....	144	Fall Management Blanket.....	18
Electrosurgical Tip Cleaners.....	140	Fall Management Slippers.....	17
Eliminators, Odor.....	105	Fascitis Night Splints, Plantar.....	156
Emergency Carts - Aluminum.....	62	Fascitis Stretch Splints, Plantar.....	156
Emergency Carts - Steel.....	61	Fecal Occult Test (gFOB).....	123
Emergency Carts, Deluxe-Steel.....	62	Fecal Occult Test (iFOB).....	123
Emergency Carts, Pediatric Color-Coded, Aluminum.....	61	Fetal Monitoring Bands, Abdominal.....	127
Emesis Bags and Dispensers.....	11	Fetal Monitoring Belts.....	127
Emesis Basins.....	11	Figure Eight Ankle Braces.....	155
Enema Bag Kits.....	160	Film Dressings, Transparent.....	176
Enzymatic Cleaner, Pre-Soak.....	34	Fingernail Clippers.....	9
Enzyme Cleaner Pre-Soak.....	34	First Aid Ointments.....	171
Enzyme Spray Pre-Cleaner.....	34	Fixed Angle Rotor Bench-Top Centrifuges.....	118
Esmark Bandages.....	142	Fixed Angle Rotor High-Capacity Centrifuges.....	118
ETS Nitrile Exam Gloves.....	96	Fixed Ankle Walkers, Low Profile.....	155
Exam Capes.....	16	Flexible Drinking Straws.....	6
Exam Drape/Strecher Sheets.....	91	Flexible Fabric Adhesive Bandages.....	171
Exam Gloves, ETS Nitrile.....	96	Floor Exam Lights, Economy.....	51
Exam Gloves, High-Risk.....	98	Floor Scales, Digital with Cable Remote Control.....	40
		Fluid-Resistant Gowns.....	23

Flushable Adult Wipes.....	107	Hand and Body Lotion.....	8
Flushable Cleansing Cloths.....	106	Hand Sanitizers, Instant.....	158
Foam Cervical Collars.....	151	Hand-Held Aneroid Sphygmomanometers, Deluxe.....	87
Foam Head Positioners.....	147	Hand-Held Aneroid Sphygmomanometers, Latex-Free.....	87
Foam Heel and Elbow Protectors.....	147	Handle Covers, Light.....	144
Foam Liner, Carafe.....	11	Handles, Steering.....	46
Foam Mattress Overlays, Pressure Relieving.....	168	Handrail Scales, Portable.....	42
Foam Padding.....	148	Hands Free Exam Stools.....	56
Foam Positioning Wedges.....	147	Hard Sole Slippers.....	18
Foam Seat Cushions, Pressure Relieving.....	168	Head Cover.....	20
Foam Telemetry Pouches.....	18	Head Positioners, Foam.....	147
Foam, Abduction Pillows.....	147	Headsets.....	13
Foam, Cradle Positioners.....	148	Heat Therapy Pads.....	160
Foam, Heel and Elbow Protectors.....	150	Heat Therapy Pumps.....	160
Foaming Skin Cleansers.....	105	Heavy Base Chrome IV Stands.....	45
Folding Bariatric Wheelchairs.....	138	Heavy-Duty Aluminum Commodes.....	133
Folding Steel Commodes.....	132	Heavy-Duty Biohazard Waste Bags.....	99
Folding Walkers, Bariatric Dual Release.....	135	Heavy-Duty Overbed Tables with Vanity.....	81
Folding Walkers, Dual-Release.....	135	Heavy-Duty Pill Crushing System.....	161
Folding Walkers, Single Release.....	135	Heavy-Duty Wheelchairs.....	137
Folding Wheelchair Scales with Handrail.....	43	Heavy-Duty Wire Baskets.....	47
Foley Insertion Trays.....	170	Heel and Elbow Protectors.....	150
Footprinters.....	129	Heel and Elbow Protectors, Foam.....	147
Forceps.....	26	Heel Boots, Ultimate Off-Loading.....	169
Fork.....	7	Heel Pillows.....	169
Four-Layer Compression Bandage System.....	174	Heel Protectors, Premium.....	169
Fracture Bedpan.....	10	Heel Protectors, Standard.....	169
Friction Mats, Instrument.....	145	Heel Stick Devices, Infant.....	127
G		Heel Warmers, Infant.....	127
Gait Belts.....	150	Height Rods.....	37
Gel packs, Reusable Hot or Cold.....	159	Helicobacter Pylori Test.....	124
Gel, Ultrasound.....	89	Hemostats.....	27
General Purpose Cotton-Tipped Applicators.....	89	High Capacity Digital Column Scales with BMI.....	39
gFOB Test.....	123	High Capacity Digital Column Scales with BMI and Digital Height Rod.....	39
Glove Box Holders.....	31	High-Capacity Centrifuges.....	118
Gloves, ETS Nitrile Exam.....	96	High-Risk Gloves.....	98
Gloves, Latex Exam.....	96	Holder, Apparel.....	30
Gloves, Sterile Procedure.....	98	Holder, Glove Box.....	31
Gloves, Stretch Vinyl Exam.....	97	Holder, Multipurpose Tube.....	160
Gloves, Stretchy Nitrile Exam.....	97	Holder, Needle.....	27
Gloves, Vinyl Exam.....	97	Holder, Oxygen Cylinder.....	46
Glue, Skin Adhesive.....	145	Holder, Trach Tube.....	167
Glycerin Swabs, Lemon.....	158	Hooks, Organizer.....	46
Gowns, AAMI Rated Material.....	25	Horizontal-Spin Centrifuges.....	119
Gowns, Cover.....	23	Hot Packs, Instant.....	159
Gowns, Exam.....	16	Hot Reusable Gel Packs.....	159
Gowns, Fluid-Resistant.....	23	Human Chorionic Gonadotropin Test, Cards/Cassettes.....	124
Gowns, Polyethylene.....	24	Hydrocolloids, Odor Control.....	177
Graduates, Triangular.....	7	Hydrocolloids, Standard.....	177
Guedel Airways, Color-Coded.....	15	Hydrocolloids, Thin.....	177
H		Hydrogels.....	177
H Pylori Test.....	124	Hygiene Stations, Respiratory.....	33
Hair Shampoo, Baby.....	126	Hyperthermia, Hypo Blankets.....	143
Hairbrush.....	7	Hypo Hyperthermia Blankets.....	143
Halogen, Minor Surgery Lighting.....	53	I	
Halogen, Procedure Exam Lights - K Arm.....	52	ID Cards.....	6
Hamper Accessories.....	45	iFOB Test.....	123
Hampers, 18 Round Chrome.....	45	Immobilizers with Compression Straps, Universal Knee.....	153
Hampers, 18 Round Tilt-Lid Chrome.....	44	Immobilizers, Tri-Panel Universal Knee.....	153
Hampers, Adjustable Folding Chrome.....	44	Immobilizers, Universal Knee.....	152
Hampers, Chrome Linen with Optional Labels.....	44	Immobilizers, Universal Shoulder.....	152
Hampers, Chrome Square Tilt-Lid.....	44	Immobilizers, Velpeau Shoulder.....	152

Impervious Stockinettes.....	148	Kits, Nursery.....	128
In-Bed Patient Positioning.....	169	Kits, Post Mortem.....	161
Incision and Drainage Trays.....	109	Kits, Quiet.....	13
Incubators.....	119	Kits, Suction Catheter.....	165
Infant and Baby Scales.....	38	Knee Immobilizers, Tri-Panel Universal.....	153
Infant Heel Lancets.....	127	Knee Immobilizers, Universal.....	152
Infant Heel Warmers.....	127	Knee Immobilizers, Universal Knee with Compression Strap.....	153
Infant Measuring Mats.....	37	Knife.....	7
Infant Scale Liner.....	129	Knit Elastic Bandages.....	173
Infant Scales, Neo-Natal and Intensive Care.....	37	Knit Pants, Seamless.....	105
Infant Stethoscopes, Dual-Head.....	93		
Infectious Linen Bags.....	104	L	
Infusion Pump Supports.....	46	Lab Coats, Premium Disposable.....	22
Insertion Trays, Foley.....	170	Lab Coats, Standard Disposable.....	21
Instant Cold Packs.....	159	Lab Markers.....	118
Instant Hand Sanitizers.....	158	Lab Microscopes, Binocular.....	120
Instant Warm Packs.....	159	Lab Specimen Containers, Red.....	125
Instrument Cleaners, Neutral pH.....	35	Laceration Trays.....	110
Instrument Lubricants.....	34	Lamps, Diagnostic and CoolSpot.....	53
Instrument Magnetic Drapes.....	145	Lamps, Diagnostic and Outpatinet.....	54
Instrument Stands, Mobile.....	47	Lancets, Safety.....	122
Instrument Tables.....	50	Lap Sponges.....	144
Instruments, Forceps.....	26	Laparotomy Sponges.....	144
Instruments, Hemostats.....	27	Latex Exam Gloves.....	96
Instruments, Needle Holders.....	27	Latex High Risk Exam Gloves.....	98
Instruments, Scissors.....	27	LED, Minor Surgery Lighting.....	54
Insufflation Tubing.....	141	LED, Procedure Exam Lights - K Arm.....	52
Intensive Care Infant Scales, Neo-Natal.....	37	Leg Bags.....	170
Irrigation Bulb, Ear and Ulcer Syringes.....	126	Leg Straps, Catheter.....	170
Irrigation Syringes.....	170	Lemon Glycerin Swabs.....	158
Irrigation Trays with Bulb Syringe.....	114	Light Handle Covers.....	144
Irrigation Trays with Piston Syringe.....	114	Light-Weight Cover Gowns.....	25
Isolation Carts, Classic-Aluminum.....	63	Light-Weight Stethoscopes, Dual Head.....	94
Isolation Carts, Classic-Steel.....	62	Light-Weight Stethoscopes, Single-Head.....	93
Isolation Carts, Mini-Line-Steel.....	63	Lighting, Minor Surgery AIM 50.....	53
Isolation Cover Gowns, Yellow.....	23	Lights, Diagnostic and Outpatient.....	54
Isolation Gowns, AAMI Rated.....	24	Lights, Economy Floor Exam.....	51
Isolation Gowns, Non-AAMI.....	24	Lights, K Arm Task.....	52
Isolation Masks.....	22	Lights, Minor Surgery - Halogen.....	53
IV Pole Outlet Strips.....	46	Lights, Minor Surgery - LED.....	54
IV Stands, Chrome.....	45	Lights, Multi-Purpose Mobile Exam.....	51
IV Stands, Heavy Base Chrome.....	45	Lights, Procedure Exam.....	51
IV Stands, Stainless.....	45	Lights, Procedure Exam - K Arm - Halogen.....	52
IV Start Kits.....	115	Lights, Procedure Exam - K Arm - LED.....	52
IV Start Kits with ChloroPrep.....	116	Lights, Procedure Exam LED.....	53
IV Start Kits With Preventics.....	117	Lightweight Rigid Walkers.....	135
		Lightweight Wheelchairs.....	138
J		Limited Reusable Pillows.....	10
Jackets, Disposable Staff.....	21	Linen Hampers, Chrome with Optional Labels.....	44
Jelly, Lubricating.....	90	Linen, Infectious Bags.....	104
		Linen, Soiled Bags.....	104
K		Liners, Baby Scale.....	129
K Arm Task Lights.....	52	Liners, Can.....	100
Keyboards, Silicone.....	83	Liners, Kick Bucket.....	142
Keyboards, Submersible True-Type.....	83	Lip Balm.....	8
Keyboards, Waterproof True-Type.....	83	LiquiBand Skin Adhesive.....	145
Kick Bucket Liners.....	142	Long Torso Seamless Knit Pants.....	130
Kick Buckets.....	50	Loops, Vessel.....	143
Kidney Dish, Emesis Basins.....	11	Lotion, Baby.....	126
Kidney Stone Collector, Calculi Strainer.....	170	Lotion, Daily Moisturizing.....	8
Kits, Ankle Stirrup Rehabilitation.....	154	Lotion, Hand and Body.....	8
Kits, Enema Bag.....	160	Low Access Powered Exam Tables.....	58
Kits, IV Start.....	115	Low Access Powered Podiatry Chairs.....	59
Kits, IV Start with ChloroPrep.....	116		

Low Access Powered Procedure Chairs.....	59	Mid-Weight Cover Gowns.....	25
Low Profile Fixed Ankle Walkers.....	155	Mini-Line Isolation Carts - Steel.....	63
Lubricants, Instrument.....	34	Minor Procedure Custom Kits and Trays.....	108
Lubricating Jelly.....	90	Minor Surgery Lighting - AIM 50.....	53
Lumbar Puncture Trays.....	109	Minor Surgery Lighting - Halogen.....	53
		Minor Surgery Lighting - LED.....	54
		Mitts, Padded.....	150
M		Mixers, Vortex.....	120
Magnetic Drapes, Instrument.....	145	Mobile ABI Comprehensive Doppler Systems.....	88
Male Urinals.....	11	Mobile Aneroid Sphygmomanometers.....	95
Manual Recliners, Standard.....	71	Mobile Instrument Stands.....	47
Manual Treatment Recliners.....	70	Mobile Workstations.....	61
Markers, Correct Site Non-Sterile Skin.....	146	Mobile Workstations Accessories.....	61
Markers, Dual Tip Sterile Skin.....	145	Mobility Aid Accessories.....	136
Markers, Lab.....	118	Moisturizer, Baby Lotion.....	126
Markers, Skin and Utility with Pre-printed Labels.....	146	Monitor Carts - Steel.....	64
Markers, Standard Sterile Skin.....	146	Monitoring Bands, Abdominal Fetal.....	127
Markers, Utility.....	144	Monitoring Belts, Fetal.....	127
Masks, Aerosol.....	166	Mono Test.....	124
Masks, Anesthesia.....	14	Mononucleosis Test.....	124
Masks, Face Tents.....	166	Mouse, Computer.....	83
Masks, Oxygen.....	166	Mouthwash.....	8
Masks, Procedure.....	22	Mucous Specimen Traps.....	165
Masks, Procedure with Visor.....	22	Mugs, Patient.....	10
Masks, Surgical Tie.....	22	Multi-Purpose Mobile Exam Lights.....	51
Masks, Surgical with Anti-Fog.....	23	Multi-Purpose Phlebotomy Carts - Hybrid Aluminum.....	65
Masks, Surgical with Visor.....	22	Multi-Purpose Pill Crushers.....	161
Masks, Trach.....	166	Multipurpose Tube Holders.....	160
Maternity Pad.....	130		
Mats, Baby Measuring.....	37	N	
Mats, Magnetic Instrument.....	145	Nail Clippers.....	9
Mattress Overlays, Pressure Relieving Foam.....	168	Nail Polish Remover Pads.....	157
Mattress Pads.....	69	Nasopharyngeal Airways.....	15
Mattress, Synthetic Overlays.....	168	Nebulizers.....	166
Mayo Stands, Chrome Double Post.....	47	Needle Holders.....	27
Mayo Stands, Stainless.....	47	Neo-Natal and Intensive Care Infant Scales.....	37
Measuring Mats, Baby.....	37	Neo-Natal Scales with Handles 1 Gram Accuracy.....	41
Measuring Rods.....	37	Neo-Natal Stethoscopes, Dual-Head.....	93
Mechanical Physician Scales.....	38	Net Bandages, Tubular.....	172
Medical Cart Accessories.....	66	Neutral PH Instrument Cleaners.....	35
Medical Carts, Anesthesia - Aluminum.....	64	Next Generation Nitrile Exam Gloves.....	97
Medical Carts, Anesthesia - Steel.....	63	NIBP Hoses.....	162
Medical Carts, Cast-Steel.....	65	Night Splints, Plantar Fasciitis.....	156
Medical Carts, Classic Isolation - Aluminum.....	63	Nitrile Exam Gloves, ETS.....	96
Medical Carts, Classic Isolation - Steel.....	62	Nitrile Exam Gloves, Next Generation.....	97
Medical Carts, Cylinder Bulk Transport.....	65	Nitrile Exam Gloves, Stretchy.....	97
Medical Carts, Cylinder Patient Use.....	65	Nitrile High Risk Exam Gloves.....	98
Medical Carts, Deluxe Emergency - Steel.....	62	No-Rinse Shampoo and Body Bath.....	8
Medical Carts, Emergency - Pediatric Color-Coded, Aluminum.....	61	Non-AAMI Isolation Gowns.....	24
Medical Carts, Emergency - Aluminum.....	62	Non-Woven Telemetry Pouches with Neck and Waist Straps.....	18
Medical Carts, Emergency - Steel.....	61	Non-Woven Telemetry Pouches with Neck Straps.....	18
Medical carts, Mini-Line Isolation - Steel.....	63	Nursery Kits.....	128
Medical Carts, Monitor - Steel.....	64	Nylon Base Exam Stools.....	56
Medical Storage Cart Accessories.....	66		
Medical Storage Carts, Column.....	66	O	
MediChoice Storage Carts, Catheter Procedure.....	67	OB Doppler Systems.....	87
Medicine Cups.....	7	OB/GYN/Proctoscopic Rayon Tipped Applicators.....	89
Mesh Briefs.....	130	Obstetric Pads.....	130
Mice, Computer.....	83	Obstetrical Doppler Systems.....	87
Microhematocrit Centrifuges.....	119	Obstetrical Towelettes.....	157
Microscope Optical Cleaning Kits.....	121	Odor Control Hydrocolloids.....	177
Microscopes (Advanced).....	121	Odor Eliminators.....	105
Microscopes (Binocular).....	120	Offset Adjustable Canes.....	136
Microscopes (Economy).....	120	Optical Cleaning Kits, Microscope.....	121

Saline Wound Flush	176	Single-Head Stethoscopes, Bowles	93
Sanitizers, Instant Hand	158	Single-Head Stethoscopes, Light-Weight	93
Scale Liners, Baby	129	Skin Adhesive	145
Scales, High Capacity Digital Column with BMI and Digital Height Rod	39	Skin and Utility Markers with Pre-printed Labels	146
Scales with BMI, Digital Height Rod and Wireless	37	Skin Cleansers	177
Scales, Combination Baraitric	38	Skin Cleansers, Foaming	105
Scales, Convertible Chair	42	Skin Cleansing Cremes	107
Scales, Digital Chair	39	Skin Cleansing Towelettes	157
Scales, Digital Floor with Cable Remote Control	40	Skin Closure Strips	172
Scales, Digital Portable	41	Skin Markers, Correct Site Non-Sterile	146
Scales, Digital with Handrail	41	Skin Markers, Dual Tip Sterile	145
Scales, Digital BMI with Handrail	40	Skin Markers, Standard Sterile	146
Scales, Folding Wheelchair with Handrail	43	Skin Protectant	107
Scales, High Capacity Digital Column with BMI	39	Skin Recovery	107
Scales, Infant and Baby	38	Skin Staple Remover Trays	113
Scales, Mechanical Physician	38	Skin Treatment	107
Scales, Neo-Natal and Intensive Care Infant Scales	37	Sleeper Chairs, Classic Overnighter	74
Scales, Neo-Natal with Handles 1 Gram Accuracy	41	Sleeper Chairs, Deluxe Overnighter	74
Scales, Portable Bariatric with BMI	41	Sleeves, Arm	176
Scales, Portable Handrail	42	Slides and Cover Glasses	121
Scales, Take Home	44	Sling and Swathes	152
Scales, Tilt and Roll Wheelchair with Handrail	43	Slings, Arm	151
Scales, Wheelchair	39	Slippers, Fall Management	17
Scalpel Blades	145	Slippers, Hard Sole	18
Scalpels, Disposable Safety	141	Slippers, Smooth Side Out	17
Scented Adult Wipes	107	Slippers, Textured Side Out	17
Scissors	27	Slippers, Universal	17
Scissors, Utility	26	Smooth Side Out Slippers	17
Scope Warmers	143	SMS Shoe Covers	23
Scrub Solutions, Povidone Iodine	176	Soap Dishes and Wash Basins	11
Scrubs, Disposable	21	Socks, Baby	128
Seamless Kint pants, Long Torso	130	Sodium Chloride Solution	167
Seamless Knit Pants	105	Soft Tissue Biopsy Trays	112
Seat Cushions, Air Lock	168	Soft-Flex Adhesive Bandages	171
Seat Cushions, Pressure Relieving Foam	168	Soiled Linen Bags	104
Seat, Air Lift Cushions	168	Solid Maple Wooden Canes	136
Seats, Adjustable Height Raised Toilet	131	Solidifiers	140
Seats, Bariatric Bath	131	Solution, Anti-Fog	143
Seats, Bath	131	Solution, Sodium Chloride	167
Seats, Raised Toilet	131	Souffle Cups	6
Self-Adhering Bandages, Cohesive	172	Specialty Dispensers and Accessories	30
Sensors	162	Specialty Products	34
Set Up Bags, Respiratory	167	Specimen Bags	125
Sets, Commode Pail	134	Specimen Containers (24 Hour)	124
Shampoo and Body Bath	7	Specimen Containers (Red Lab)	125
Shampoo and Body Bath, No-Rinse	8	Specimen Containers (Tamper Proof)	125
Shampoo and Body Wash	8	Specimen Containers (Urinalysis)	125
Shampoo, Baby	126	Specimen Pans	7
Shaving Cream	8	Specimen Traps	165
Sheer Adhesive Bandages	171	Specula, Vaginal	90
Shirt, Disposable Scrubs	21	Sphygmomanometers Aneroid (Latex-Free)	87
Shoe Covers	23	Sphygmomanometers Aneroid Deluxe (Latex-Free)	87
Shoes, Classic Post Op	156	Sphygmomanometers, Mobile Aneroid	95
Shorts, Exam Disposable	16	Sphygmomanometers, Wall Mount Aneroid	95
Shoulder Immobilizers, Universal	152	Splints, Plantar Fascitis Night	156
Shoulder mmobilizers, Velpeau	152	Splints, Plantar Fascitis Stretch	156
Shower Chairs, PVC Commode	132	Sponge Counters	145
Shower Transport Chairs	132	Spoon	7
Side Chairs	60	Sprague-Rappaport Stethoscopes, Dual-Head	95
Silicone Keyboards	83	Spray, Perineal Wash	129
Single Door Blanket Warmers	50	Square Tilt-Lid Chrome Hampers	44
Single Release Folding Walkers	135	Stability Weights	46

Stackable Commode Bedpan	10	Stethoscopes Disposable (Plastic)	92
Stackable, Stainless Step Stools	48	Stethoscopes Dual-Head (Adult/Pediatric Professional)	94
Staff Cover Gowns, Light-Weight	25	Stethoscopes Dual-Head (Cardiology Professional)	95
Staff Cover Gowns, Mid-Weight	25	Stethoscopes Dual-Head (Infant)	93
Staff Gowns, AAMI Rated Material	25	Stethoscopes Dual-Head (Light-Weight)	94
Staff Gowns, Cover	23	Stethoscopes Dual-Head (Neo-Natal)	93
Staff Gowns, Fluid-Resistant	23	Stethoscopes Dual-Head (Pediatric)	93
Staff Gowns, Polyethylene	24	Stethoscopes Dual-Head (Sprague-Rappaport)	95
Staff Isolation Gowns, AAMI Rated	24	Stethoscopes Single-Head (Bowles)	93
Staff Isolation Gowns, Non-AAMI	24	Stethoscopes Single-Head (Light-Weight)	93
Staff, Disposable Jackets	21	Stirrup, Ankle Braces with Air	154
Stainless Contoured Seat Exam Stools	57	Stirrup, Ankle Braces with Foam and Terry Liners	154
Stainless IV Stands	45	Stirrup, Ankle Braces with Gel and Air	154
Stainless Mayo Stands	47	Stirrup, Ankle Rehabilitation Kits	154
Stainless Padded Seat Exam Stools	57	Stockinettes, Bias-Cut	149
Stainless Stackable Step Stools	48	Stockinettes, Impervious	148
Stainless Step Stools	48	Stockinettes, Sterile 1-Ply	148
Stainless Utility Carts	82	Stockinettes, Sterile 2-Ply	149
Stainless Utility Tables with Shelf and Four-Sided Rail	82	Stockinettes, Tubular	148
Stainless Utility Tables with Shelf, Drawer and Rail	82	Stool Accessories, Exam	57
Standard Blood Drawing Chairs	78	Stools, Bariatric Step	54
Standard Clinical Care Recliners	75	Stools, Black Base Exam	57
Standard Disposable Coveralls	21	Stools, Chrome Step	48
Standard Disposable Lab Coats	21	Stools, Deluxe Exam	55
Standard Heel Protectors	169	Stools, Economy Exam	54
Standard Hydrocolloids	177	Stools, Hands Free Exam	56
Standard Manual Recliners	71	Stools, Nylon Base Exam	56
Standard Powered Podiatry Chairs	59	Stools, Stainless Contoured Seat Exam	57
Standard Powered Procedure Chairs	60	Stools, Stainless Padded Seat Exam	57
Standard Shoe Covers	23	Stools, Stainless Stackable Step	48
Standard Sterile Skin Markers	146	Stools, Stainless Step	48
Standard Surgical Bras	147	Stools, Value Exam	55
Standard Wheelchairs	137	Storage Cart Accessories, Medical	66
Stands, Basin	47	Straight Line Treatment Tables	79
Stands, Chrome Double Post Mayo	47	Strainers, Calculi	170
Stands, Chrome IV	45	Straps, Catheter Leg	170
Stands, Heavy Base Chrome IV	45	Straws, Flexible Drinking	6
Stands, IV Stainless	45	Strep A Test	124
Stands, Mobile Instrument	47	Stretch Splints, Plantar Fasciitis	156
Stands, Stainless Mayo	47	Stretch Vinyl Exam Gloves	97
Stations, Respiratory Hygiene	33	Stretcher Replacement Pads	67
Steel, Anesthesia Carts	63	Stretcher/Drape Sheet	91
Steel, Bariatric Crutches	134	Stretchy Nitrile Exam Gloves	97
Steel, Cast Carts	65	Strips, IV Pole Outlet	46
Steel, Emergency Carts	61	Submersible True-Type Keyboards	83
Steel, Folding Commodes	132	Suction Catheter Kits	165
Steel, Painted Commodes	133	Suction Catheter Trays	165
Steel, Deluxe Emergency Carts	62	Suction Catheters	165
Steel, Mini-Line Isolation Carts	63	Suction Tubing	140
Steel, Monitor Carts	64	Suction Yankauers	140
Steering Handles	46	Support Bandages, Tubular Elastic	175
Step Stools, Bariatric	54	Supports, Abdominal Binders	150
Step Stools, Chrome	48	Supports, Infusion Pump	46
Step Stools, Stainless	48	Surgeon Caps	20
Step Stools, Stainless Stackable	48	Surgeon Hoods	20
Sterile 1-Ply Stockinettes	148	Surgery Lighting, Minor AIM 50	53
Sterile 2-Ply Stockinettes	149	Surgical Blades	145
Sterile Alcohol	118	Surgical Bras, Premium	146
Sterile Impervious Stockinette	148	Surgical Bras, Standard	147
Sterile Procedure Gloves	98	Surgical Masks with Anti-Fog	23
Sterilization Pouches	26	Surgical Masks with Visor	22
Stethoscopes Disposable (Chrome)	92	Surgical Tie Masks	22

Surgical Vests	146	Topical Cleansing Towelettes	157
Suture Removal Trays	113	Topical Ointments, First Aid	171
Suture Trays	109	Topical Povidone Iodine Solutions, Prep	176
Swabettes, Oral	158	Topical Povidone Iodine Solutions, Scrub	176
Swabs, Alcohol	158	Topical Povidone Iodine Swabsticks, Sterile	176
Swabs, Cotton	91	Tops, Disposable Scrubs	21
Swabs, Lemon Glycerin	158	Tourniquets (Latex)	122
Swabs, Oral	158	Tourniquets (Non-Latex)	122
Swabsticks, Alcohol	158	Towel Clamps	27
Swathe and Sling	152	Towelettes, BZK	157
Synthetic Impervious Stockinette	148	Towelettes, Cleansing	157
Synthetic Mattress Overlays	168	Towelettes, Obstetrical	157
Syringes, Ear and Ulcer	126	Towels, Disposable OR	142
Syringes, Irrigation	170	Towels, Professional	91
T			
Table-Top Obstetrical Doppler Systems	87	Trach Masks	166
Table-Top Vascular Doppler Systems	88	Trach Tube Holders	167
Table, Exam Accessories	58	Tracheostomy Trays	114
Table, Exam Paper	91	Transfer Bench with Backrest	131
Table, OR Replacement Pads	68	Transfer Bench, Bariatric with Backrest	132
Tables, Bariatric Exam	60	Transparent Film Dressings	176
Tables, Economy Overbed	81	Transport Chairs, Shower	132
Tables, Heavy-Duty Overbed with Vanity	81	Transport Wheelchairs	138
Tables, Instrument	50	Transporters/Recliners	71
Tables, Low Access Powered Exam	58	Transporters/Recliners, Bariatric	73
Tables, Pediatric Treatment with Built-in Digital Scale	80	Traps, Specimen	165
Tables, Pediatric Treatment with Sliding Doors	79	Trash Can Bags	100
Tables, Reversible Manual Exam	58	Trash Cans	51
Tables, Stainless Utility with Shelf and Four-Sided Rail	82	Trays, Bone Marrow	112
Tables, Stainless Utility with Shelf, Drawer and Rail	82	Trays, Central Line Dressing	111
Tables, Straight Line Treatment	79	Trays, Foley Insertion	170
Tablets, Denture Cleanser	9	Trays, Incision and Drainage	109
Take Home Scales	44	Trays, Irrigation with Bulb Syringe	114
Talc Free Baby Powder	126	Trays, Irrigation with Piston Syringe	114
Tamper Proof Specimen Containers	125	Trays, Laceration	110
Tape Remover Wipes, Adhesive	157	Trays, Lumbar Puncture	109
Task Lights, K Arm	52	Trays, Phlebotomy	122
Task Wipes	118	Trays, Skin Staple Remover	113
Telemetry Pouches, Clear Polyethylene	18	Trays, Soft Tissue Biopsy	112
Telemetry Pouches, Foam	18	Trays, Suction Catheter	165
Telemetry Pouches, Non-Woven with Neck and Waist Straps	18	Trays, Suture	109
Telemetry Pouches, Non-Woven with Neck Straps	18	Trays, Suture Removal	113
Tents, Face	166	Trays, Tracheostomy	114
Textured Side Out Slippers	17	Trays, Urethral Catheterization	170
Therapy Pads, Heat	160	Treatment Chairs	70
Therapy Pumps, Heat	160	Treatment Chairs, Bariatric	72
Thermometers, Digital with Sheaths	84	Treatment Chairs, Bariatric Power Recline	73
Thin Hydrocolloids	177	Treatment Chairs, Power Recline	72, 76
Three-Layer Compression Bandage System with Trak™ Bandage	174	Treatment Recliners, Manual	70
Three-Layer Reduced Compression Bandage System	174	Treatment Tables, Pediatric with Built-in Digital Scale	80
Tilt and Roll Wheelchair Scales with Handrail	43	Treatment Tables, Pediatric with Sliding Doors	79
Tip Cleaners, Electrosurgical	140	Treatment Tables, Straight Line	79
Tips/Pencils, Electrosurgical	144	Treatment, Skin	107
Tissues, Facial	13	Tri-Panel Universal Knee Immobilizers	153
Toenail Clippers	9	Triangular Graduates	7
Toilet Seats, Adjustable Height Raised	131	Trinocular Microscopes, Advanced	121
Toilet Seats, Raised	131	Triple Antibiotic Ointment, First Aid	171
Tongue Depressors	158	Tube Holders, Multipurpose	160
Toothbrush	7	Tube Holders, Trach	167
Toothpaste	9	Tube Organizers	123
Topical Cleansers	177	Tube Racks	123
Topical Cleansing Creams	107	Tube Rockers	120
		Tube Rotators	120

Tube Rotor Microhematocrit Centrifuges	119
Tubing Connectors	143
Tubing, Insufflation	141
Tubing, Oxygen	167
Tubing, Suction	140
Tubular Bandage, Stockinette	148
Tubular Bandage, Stockinette 1-Ply Sterile	149
Tubular Bandage, Stockinette 2-Ply Sterile	148
Tubular Elastic Support Bandages	175
Tubular Net Bandages	172
Tubular Stockinettes	148
Tumblers and Carafes	11
Turn and Positioning	169
Twenty-Four Hour Specimen Containers	124
Two Piece Cervical Collars	151

U

Ulcer and Ear Syringes	126
Ultimate Off-Loading Heel Boots	169
Ultrasound Gel	89
Umbilical Cord Clamp Cutters	129
Umbilical Cord Clamps	129
Under Counter Blanket Warmers	49
Underpads	105
Underwear, Protective	106
Underwear, Seamless Knit	105
Underwear, Seamless Knit, Long Torso	130
Universal Clamps	46
Universal Knee Immobilizers	152
Universal Knee Immobilizers with Compression Straps	153
Universal Shoulder Immobilizers	152
Universal Slippers	17
Universal Surgeon Hood	20
Unna Boot	174
Unscented Adult Wipes	107
Urethral Catheterization Trays	170
Urinals, Male	11
Urinalysis Specimen Containers	125
Urinary Leg Bags	170
Urine Collectors, Pediatric	170
Urine Containers, Specimen	125
Urine Specimen Containers, 24 Hour	124
Utility and Skin Markers with Pre-printed Labels	146
Utility Carts with Wire Baskets	81
Utility Carts, Plastic	81
Utility Carts, Plastic and Aluminum	82
Utility Carts, Stainless	82
Utility Markers	144
Utility Scissors	26
Utility Tables, Stainless with Shelf and Four-Sided Rail	82
Utility Tables, Stainless with Shelf, Drawer and Rail	82

V

Vaginal Specula	90
Value Exam Stools	55
Vascular Doppler Systems	88
Velpeau Shoulder Immobilizers	152
Vessel Loops	143
Vests, Surgical	146
Vinyl Exam Gloves	97
Vinyl, Exam Gloves Stretch	97
Vomit Bags and Dispensers	11
Vomit Basins	11
Vortex Mixers	120

W

Walker Accessories	136
Walkers, Air Shell Ankle	155
Walkers, Bariatric Dual Release Folding	135
Walkers, Dual-Release Folding	135
Walkers, Lightweight Rigid	135
Walkers, Low Profile Fixed Ankle	155
Walkers, Single Release Folding	135
Wall Mount Aneroid Sphygmomanometers	95
Warm Packs, Instant	159
Warmers, Counter Top Blanket	49
Warmers, Double Door Blanket	50
Warmers, Infant Heel	127
Warmers, Scope	143
Warmers, Single Door Blanket	50
Warmers, Under Counter Blanket	49
Wash Basins and Soap Dishes	11
Wash, Perineal Spray	129
Washcloths	106
Waste Bags, Heavy-Duty Biohazard	99
Waste Disposal, Solidifiers	140
Waste Receptacles	51
Waterproof True-Type Keyboards	83
Wax Remover	84
Wedges, Foam Positioning	147
Weights, Stability	46
Wheelchair Accessories	139
Wheelchair Scales	39
Wheelchair Scales, Folding with Handrail	43
Wheelchair Scales, Tilt and Roll with Handrail	43
Wheelchair Seat Cushions, Air Lock	168
Wheelchairs, Bariatric Folding	138
Wheelchairs, Heavy-Duty	137
Wheelchairs, Reclining	138
Wheelchairs, Standard	137
Wheelchairs, Transport	138
Wheelchairs, Lightweight	138
Wipes, Baby	126
Wipes, Flusable Cleansing Cloths	106
Wipes, Flushable Adult	107
Wipes, Scented Adult	107
Wipes, Task	118
Wipes, Unscented Adult	107
Wire Baskets, Heavy-Duty	47
Wooden Canes	136
Workstation Accessories, Mobile	61
Workstations, Mobile	61
Wound Cleansers	177
Wound Closure Strips	172
Wound Closure, Skin Adhesive	145
Wound Dressing, Hydrogels	177
Wound Dressing, Odor Control Hydrocolloids	177
Wound Dressing, Standard Hydrocolloids	177
Wound Dressing, Thin Hydrocolloids	177
Wound Flush, Saline	176
Wound Protectant, Hydrogels	177
Woven Elastic Bandages	173
Wraps, Cohesive	172

X

X-Ray Replacement Pads	69
------------------------------	----

Y

Yankauers	140
Yellow Cover Gowns	23

A

Abdominal Fetal Monitoring Bands.....	127
Adjustable Folding Chrome Hampers.....	44
Adult/Pediatric Professional Stethoscopes, Dual-Head.....	94
Advanced Microscopes.....	121
Alcohol and Peroxide.....	160
Alcohol Dispensers.....	118
Aluminum Emergency Carts.....	62
Aneroid Sphygmomanometers, Deluxe.....	87
Aneroid Sphygmomanometers, Mobile.....	95
Aneroid Sphygmomanometers, Wall Mount.....	95
Antiperspirant Deodorant, Roll On.....	8
Applicators - General Purpose (Cotton-Tipped).....	89
Arm Sleeves.....	176

B

Baby Bassinets.....	48
Baby Bonnets.....	128
Baby Shampoo.....	126
Bandage System with Trak™ Bandage, 3-Layer Compression.....	174
Bandage System, 3-Layer Reduced Compression.....	174
Bandage System, 4-Layer Compression.....	174
Bandages, Cohesive.....	172
Bandages, Esmark.....	142
Bandages, Knit Elastic.....	173
Bandages, Paste.....	174
Bandages, Tubular Elastic Support.....	175
Bandages, Tubular Net.....	172
Bandages, Woven Elastic.....	173
Bands, Abdominal Fetal Monitoring.....	127
Bariatric Exam Tables.....	60
Basin Stands.....	47
Basins, Emesis.....	11
Basins, Wash.....	11
Bassinets.....	48
Bathroom Take Home Scales.....	44
Beanies, Baby.....	128
Bedpans.....	10
Belts, Fetal Monitoring.....	127
Bench-Top Centrifuges.....	118
Bench-Top Horizontal-Spin Centrifuges.....	119
Binocular Microscopes.....	120
Binocular Microscopes, Advanced.....	121
Binocular Microscopes, Economy.....	120
Biopsy Punches.....	141
Biopsy Trays, Soft Tissue.....	112
Black Base Exam Stools.....	57
Blanket Warmers, Counter Top.....	49
Blanket Warmers, Double Door.....	50
Blanket Warmers, Single Door.....	50
Blanket Warmers, Under Counter.....	49
Blanket, Fall Management.....	18

Body Bath and Shampoo.....	7
Body Bath and Shampoo, No Rinse.....	8
Body Lotion, and Hand.....	8
Bone Marrow Trays.....	112
Bonnets, Baby.....	128
Bowles Stethoscopes, Single-Head.....	93
Briefs, Mesh.....	130
Buckets, Kick.....	50
Bulk Transport, Cylinder Carts.....	65

C

Carafes and Tumblers.....	11
Cart Cleaners.....	35
Carts, Classic Isolation - Aluminum.....	63
Carts, Cylinder Bulk Transport.....	65
Carts, Cylinder Patient Use.....	65
Carts, Emergency Pediatric Color-Coded, Aluminum.....	61
Carts, Emergency - Aluminum.....	62
Carts, Emergency - Steel.....	61
Carts, Multi-Purpose Phlebotomy - Hybrid Aluminum.....	65
Carts, Utility with Wire Baskets.....	81
Central Line Dressing Trays.....	111
Centrifuges (Bench-Top).....	118
Centrifuges (High-Capacity).....	118
Centrifuges (Horizontal-Spin).....	119
Centrifuges (Microhematocrit).....	119
Chairs, Low Access Powered Podiatry.....	59
Chairs, Low Access Powered Procedure.....	59
Chairs, Side.....	60
Chairs, Standard Powered Podiatry.....	59
Chairs, Standard Powered Procedure.....	60
Chrome Disposable Stethoscopes.....	92
Chrome Hampers, 18 Round.....	45
Chrome Hampers, 18 Round Tilt-Lid.....	44
Chrome Hampers, Adjustable Folding.....	44
Clamps, Umbilical Cord.....	129
Classic Isolation Carts - Aluminum.....	63
Cleaner, Enzyme Pre-Soak.....	34
Cleaners, Cart.....	35
Cleaners, Neutral PH Instrument.....	35
Cleaners, Specialty Products.....	34
Cleaning Kits, Microscope Optical.....	121
Cleansers, Foaming Skin.....	105
Cleansing Creams.....	107
Clear Polyethylene Telemetry Pouches.....	18
Clippers.....	9
Cohesive Bandages.....	172
Comfort Pack.....	13
Compression Arm Sleeves.....	176
Compression Bandage System with Trak™ Bandage, 3-Layer.....	174
Compression Bandage System, 3-Layer Reduced.....	174
Compression Bandage System, 4-Layer.....	174

Computer Mice	83	Enzyme Cleaner Pre-Soak	34
Cord Clamp Cutters, Umbilical	129	Enzyme Spray Pre-Cleaner	34
Cord Clamps, Umbilical	129	Esmark Bandages	142
Correct Site Non-Sterile Skin Markers	146	Exam Lights, Economy Floor	51
Cotton Balls	91	Exam Stools, Black Base	57
Cotton Tipped Applicators, General Purpose	89	Exam Stools, Stainless Contoured Seat	57
Counter Top Blanket Warmers	49	Exam Stools, Stainless Padded Seat	57
Counters, Differential	119	Exam Tables, Bariatric	60
Cover Gowns	23	Exam Tables, Low Access Powered	58
Cremes, Cleansing	107	Exam Tables, Reversible Manual	58
Cups, Medicine	7		
Curettes, Ear	84	F	
Cutters, Umbilical Cord Clamp	129	Fall Management Blanket	18
Cylinder Carts - Bulk Transport	65	Fall Management Slippers	17
Cylinder Carts - Patient Use	65	Fetal Monitoring Bands, Abdominal	127
		Fetal Monitoring Belts	127
D		Film Dressings, Transparent	176
Daily Moisturizing Lotion	8	Fingernail Clippers	9
Deodorant, Roll-On	8	Fixed Angle Rotor Bench-Top Centrifuges	118
Dermal Biopsy Punches	141	Fixed Angle Rotor High-Capacity Centrifuges	118
Differential Counters	119	Floor Exam Lights, Economy	51
Digital Take Home Scales	44	Foam Liner, Carafe	11
Dispensers, Alcohol	118	Foam Telemetry Pouches	18
Disposable Safety Scalpels	141	Foam, Heel and Elbow Protectors	150
Disposable Stethoscopes, Chrome	92	Foaming Skin Cleansers	105
Disposable Stethoscopes, Plastic	92	Footprinters	129
Disposable Washcloths	106	Four-Layer Compression Bandage System	174
Double Door Blanket Warmers	50	Fracture Bedpan	10
Dressing Trays, Central Line	111		
Dressings, Transparent Film	176	G	
Dual-Head Professional Stethoscopes, Adult/Pediatric	94	Gel, Ultrasound	89
Dual-Head Stethoscopes, Infant	93	General Purpose Cotton Tipped Applicators	89
Dual-Head Stethoscopes, Light-Weight	94	Gowns, Cover	23
Dual-Head Stethoscopes, Neo-Natal	93		
Dual-Head Stethoscopes, Pediatric	93	H	
		Hair Shampoo, Baby	126
E		Hampers, 18 Round Chrome	45
Ear Curettes	84	Hampers, 18 Round Tilt-Lid Chrome	44
ECG Chart Paper	92	Hampers, Adjustable Folding Chrome	44
Economy Floor Exam Lights	51	Hand and Body Lotion	8
Economy Microscopes	120	Hand Sanitizers, Instant	158
EKG Chart Paper	92	Hand-Held Aneroid Sphygmomanometers, Deluxe	87
Elastic Bandages, Women	173	Heat Therapy Pads	160
Elastic Bandages, Knit	173	Heat Therapy Pumps	160
Elastic Support Bandages, Tubular	175	Heel and Elbow Protectors	150
Elbow and Heel Protectors	150	High-Capacity Centrifuges	118
Electrocardiogram Chart Paper	92	Horizontal-Spin Centrifuges	119
Eliminators, Odor	105	Hot Packs, Instant	159
Emergency Carts - Aluminum	62		
Emergency Carts - Steel	61	I	
Emergency Carts, Pediatric Color-Coded, Aluminum	61	ID Cards	6
Emesis Basins	11	Incubators	119
Enzymatic Cleaner, Pre-Soak	34	Infant Stethoscopes, Dual-Head	93
		Instant Hand Sanitizers	158
		Instant Warm Packs	159

Instrument Cleaners, Neutral PH	35	Mesh Briefs	130
Instrument Lubricants.....	34	Mice, Computer	83
Instrument Tables	50	Microhematocrit Centrifuges.....	119
Irrigation Trays with Bulb Syringe.....	114	Microscope Optical Cleaning Kits	121
Irrigation Trays with Piston Syringe	114	Microscopes (Advanced).....	121
Isolation Carts, Classic - Aluminum	63	Microscopes (Binocular).....	120
Isolation Cover Gowns, Yellow	23	Microscopes (Economy).....	120
IV Start Kits.....	115	Mitts, Padded.....	150
IV Start Kits with ChloraPrep	116	Mixers, Vortex	120
	J	Mobile Aneroid Sphygmomanometers.....	95
Jelly, Lubricating.....	90	Monitoring Bands, Abdominal Fetal	127
	K	Monitoring Belts, Fetal	127
Keyboards, Silicone.....	83	Mouse, Computer	83
Keyboards, Submersible True-Type	83	Mouthwash.....	8
Keyboards, Waterproof True-Type.....	83	Mugs, Patient.....	10
Kick Buckets.....	50	Multi-Purpose Phlebotomy Carts - Hybrid Aluminum	65
Kidney Dish, Emesis Basins	11		N
Kits, IV Start	115	Nail Clippers	9
Kits, IV Start with ChloraPrep	116	Neo-Natal Stethoscopes, Dual-Head	93
Kits, Nursery.....	128	Net Bandages, Tubular	172
Kits, Quiet.....	13	Neutral PH Instrument Cleaners	35
Knit Elastic Bandages	173	No-Rinse Shampoo and Body Bath	8
Knit Pants, Seamless	105	Non-Woven Telemetry Pouches with Neck and Waist Straps.....	18
	L	Non-Woven Telemetry Pouches with Neck Straps.....	18
Lab Markers	118	Nursery Kits.....	128
Lab Microscopes, Binocular	120		O
Light-Weight Stethoscopes, Dual-Head	94	Obstetric Pads	130
Light-Weight Stethoscopes, Single-Head.....	93	Odor Eliminators	105
Lights, Economy Floor Exam	51	Optical Cleaning Kits, Microscope	121
Lip Balm	8	OR Table Replacement Pads.....	68
Long Torso Seamless Knit Pants.....	130	Organizers, Tube	123
Lotion, Daily Moisturizing.....	8		P
Lotion, Hand and Body	8	Packs, Instant Warm	159
Low Access Powered Exam Tables.....	58	Padded Mitts	150
Low Access Powered Podiatry Chairs.....	59	Pads, Heat Therapy.....	160
Low Access Powered Procedure Chairs.....	59	Pads, OB.....	130
Lubricants, Instrument.....	34	Pads, OR Table Replacement.....	68
Lubricating Jelly	90	Pads, Stretcher Replacement.....	67
	M	Pail, Kick Buckets	50
Male Urinals	11	Pans, Specimen	7
Markers, Correct Site Non-Sterile Skin.....	146	Panties, Mesh	130
Markers, Lab	118	Panties, Seamless Knit.....	105
Maternity Pad	130	Panties, Seamless Knit, Long Torso	130
Mayo Stands, Stainless.....	47	Paper Cups	6
Medical Carts, Classic Isolation - Aluminum.....	63	Paste Bandages	174
Medical Carts, Cylinder Bulk Transport	65	Patient Blanket, Fall Management	18
Medical Carts, Cylinder Patient Use.....	65	Patient Identification Cards	6
Medical Carts, Emergency Pediatric Color-Coded, Aluminum.....	61	Patient Mugs	10
Medical Carts, Emergency - Aluminum.....	62	Patient Restraints, Padded Mitts.....	150
Medical Carts, Emergency - Steel	61	Patient Slippers, Fall Management	17
Medicine Cups.....	7	Patient Slippers, Smooth Side Out.....	17

Patient Slippers, Textured Side Out.....	17	Rotators, Tube.....	120
Patient Slippers, Universal.....	17	Round Chrome Hampers.....	45
Patient Use Cylinder Carts.....	65	Round Tilt-Lid Chrome Hampers	44
Pediatric Color-Coded Emergency Carts - Aluminum	61		
Pediatric Stethoscopes, Dual-Head	93	S	
Pediatric/Adult Professional Stethoscopes, Dual-Head.....	94	Safety Mitts, Padded	150
Perineal Wash Spray	129	Safety Scalpels, Disposable	141
Peroxide and Alcohol.....	160	Saline Wound Flush	176
Phlebotomy Trays.....	122	Sanitizers, Instant Hand	158
Plastic Cups.....	7, 11	Scales, Take Home.....	44
Plastic Emesis Basins	11	Scalpel Blades	145
Plastic Pitcher.....	11	Scalpels, Disposable Safety	141
Plastic Stethoscopes, Disposable	92	Seamless Kint pants, Long Torso.....	130
Podiatry Chairs, Low Access Powered	59	Seamless Knit Pants	105
Podiatry Chairs, Standard Powered	59	Self-Adhering Bandages, Cohesive.....	172
Pontoon Bedpan	10	Shampoo and Body Bath.....	7
Portion Cups	6	Shampoo and Body Bath, No-Rinse.....	8
Pouches, Non-Woven Telemetry with Neck and Waist Straps.....	18	Shampoo, Baby.....	126
Pouches, Non-Woven Telemetry with Neck Straps.....	18	Side Chairs	60
Pouches, Telemetry Clear Polyethylene	18	Silicone Keyboards	83
Pouches, Telemetry Foam	18	Single Door Blanket Warmers.....	50
Powered Exam Tables, Low Access	58	Single-Head Stethoscopes, Bowles	93
Powered Podiatry Chairs, Low Access.....	59	Single-Head Stethoscopes, Light-Weight.....	93
Powered Podiatry Chairs, Standard	59	Skin Cleansers, Foaming.....	105
Powered Procedure Chairs, Low Access.....	59	Skin Cleansing Creams	107
Powered Procedure Chairs, Standard	60	Skin Markers, Correct Site Non-Sterile	146
Pre-Cleaner, Enzyme Spray	34	Skin Protectant	107
Pre-Soak, Enzyme Cleaner.....	34	Skin Recovery.....	107
Procedure Chairs, Low Access Powered.....	59	Skin Staple Remover Trays.....	113
Procedure Chairs, Standard Powered	60	Skin Treatment	107
Professional Dual-Head Stethoscopes, Adult/Pediatric.....	94	Sleeves, Arm.....	176
Protectant, Skin	107	Slippers, Fall Management.....	17
Protectors, Heel and Elbow	150	Slippers, Smooth Side Out	17
Pumps, Heat Therapy.....	160	Slippers, Textured Side Out	17
Punches, Biopsy	141	Slippers, Universal	17
		Smooth Side Out Slippers	17
Q		Soap Dishes and Wash Basins	11
Quiet Kits.....	13	Soft Tissue Biopsy Trays.....	112
		Souffle Cups	6
R		Specialty Products	34
Racks, Tube	123	Specimen Pans.....	7
Recovery, Skin	107	Sphygmomanometers Aneroid Deluxe (Latex-Free)	87
Removal Trays, Suture	113	Sphygmomanometers, Mobile Aneroid	95
Remover Trays, Skin Staple.....	113	Sphygmomanometers, Wall Mount Aneroid	95
Replacement Pads, OR Table.....	68	Spray, Perineal Wash.....	129
Replacement Pads, Stretcher	67	Stackable Commode Bedpan	10
Retention Net Bandages	172	Stackable, Stainless Step Stools	48
Reusable Surgical Blades.....	145	Staff Gowns, Cover	23
Reversible Manual Exam Tables	58	Stainless Contoured Seat Exam Stools	57
Rinse Aids.....	35	Stainless Mayo Stands.....	47
Rinse, Mouthwash	8	Stainless Padded Seat Exam Stools	57
Rockers, Tube	120	Stainless Stackable Step Stools	48
Roll-on Deodorant.....	8		

Stainless Step Stools	48	Trays, Central Line Dressing.....	111
Stainless Utility Tables with Shelf and Four-Sided Rail	82	Trays, Irrigation with Bulb Syringe.....	114
Stainless Utility Tables with Shelf, Drawer and Rail.....	82	Trays, Irrigation with Piston Syringe	114
Standard Powered Podiatry Chairs	59	Trays, Phlebotomy.....	122
Standard Powered Procedure Chairs.....	60	Trays, Skin Staple Remover.....	113
Stands, Basin.....	47	Trays, Soft Tissue Biopsy.....	112
Stands, Stainless Mayo.....	47	Trays, Suture Removal	113
Steel Emergency Carts.....	61	Trays, Tracheostomy	114
Step Stools, Stainless	48	Treatment, Skin	107
Step Stools, Stainless Stackable	48	Trinocular Microscopes, Advanced	121
Stethoscopes Disposable (Chrome).....	92	Tube Organizers.....	123
Stethoscopes Disposable (Plastic).....	92	Tube Racks.....	123
Stethoscopes Dual-Head (Adult/Pediatric Professional)	94	Tube Rockers.....	120
Stethoscopes Dual-Head (Infant).....	93	Tube Rotators	120
Stethoscopes Dual-Head (Light-Weight).....	94	Tube Rotor Microhematocrit Centrifuges	119
Stethoscopes Dual-Head (Neo-Natal).....	93	Tubular Elastic Support Bandages.....	175
Stethoscopes Dual-Head (Pediatric)	93	Tubular Net Bandages.....	172
Stethoscopes Single-Head (Bowles).....	93	Tumblers and Carafes.....	11
Stethoscopes Single-Head (Light-Weight)	93		
Stools, Black Base Exam.....	57	U	
Stools, Stainless Contoured Seat Exam	57	Ultrasound Gel.....	89
Stools, Stainless Padded Seat Exam.....	57	Umbilical Cord Clamp Cutters	129
Stools, Stainless Stackable Step	48	Umbilical Cord Clamps.....	129
Stools, Stainless Step	48	Under Counter Blanket Warmers	49
Stretcher Replacement Pads.....	67	Underwear, Seamless Knit	105
Submersible True-Type Keyboards	83	Underwear, Seamless Knit, Long Torso.....	130
Support Bandages, Tubular Elastic	175	Universal Slippers.....	17
Surgical Blades.....	145	Urinals, Male.....	11
Suture Removal Trays	113	Utility Carts with Wire Baskets	81
Swabs, Cotton	89	Utility Tables, Stainless with Shelf and Four-Sided Rail	82
		Utility Tables, Stainless with Shelf, Drawer and Rail.....	82
T			
Table, OR Replacement Pads.....	68	V	
Tables, Bariatric Exam	60	Vomit Basins.....	11
Tables, Instrument	50	Vortex Mixers.....	120
Tables, Low Access Powered Exam.....	58		
Tables, Reversible Manual Exam	58	W	
Tables, Stainless Utility with Shelf and Four-Sided Rail	82	Wall Mount Aneroid Sphygmomanometers	95
Tables, Stainless Utility with Shelf, Drawer and Rail.....	82	Warm Packs, Instant	159
Take Home Scales	44	Warmers, Counter Top Blanket	49
Telemetry Pouches, Clear Polyethylene	18	Warmers, Double Door Blanket	50
Telemetry Pouches, Foam	18	Warmers, Single Door Blanket.....	50
Telemetry Pouches, Non-Woven with Neck and Waist Straps.....	18	Warmers, Under Counter Blanket	49
Telemetry Pouches, Non-woven with Neck Straps.....	18	Wash Basins and Soap Dishes	11
Textured Side Out Slippers.....	17	Wash, Perineal Spray	129
Therapy Pads, Heat.....	160	Washcloths	106
Therapy Pumps, Heat.....	160	Waterproof True-Type Keyboards	83
Three-Layer Compression Bandage System with Trak™ Bandage	174	Wax Remover	34
Three-Layer Reduced Compression Bandage System	174	Wound Flush, Saline	176
Toenail Clippers	9	Woven Elastic Bandages.....	173
Topical Cleansing Creams.....	107	Wraps, Cohesive	172
Tracheostomy Trays	114		
Transparent Film Dressings.....	176	Y	
Trays, Bone Marrow	112	Yellow Cover Gowns.....	23

Sustainability + Owens & Minor

As a leading healthcare logistics company, we are committed to environmentally responsible business practices. It's our core business strategy to aggregate multiple, disparate supply chains into a single, efficient delivery channel. This means we ship with more complete truckloads and ultimately reduce the number of direct, single-vendor delivery trucks required to serve our customers' daily needs.

Our Portland, Oregon distribution center models our commitment. It reuses 90% of packing materials—redistributing used pallets, recycling 15,000 in 2015—and recycles more than 175 tons of cardboard and nearly 10 tons of plastics and paper each year.

Because we manage approximately 7 million square feet of warehouse and office space and a private fleet of 360 delivery vehicles, we have an opportunity to positively impact the environment through sustainability programs that focus on efficiency in fleet management, recycling, and energy consumption reduction.

We are proud that the Owens & Minor fleet was awarded the highest grade for energy efficiency by the EPA's SmartWay Transport Partnership.

Supplier Diversity + Owens & Minor = A MediChoice Advantage

Owens & Minor has an established commitment to advancing supplier diversity. Within the healthcare industry, we have earned a highly regarded position due to the tangible impact and contributions derived from our relationships with minority-, women-, and veteran-owned business enterprises. Investing in supply chain diversity is a business imperative that benefits our provider customers, the suppliers and the community. We believe a thriving community of diverse suppliers generates healthy competition and innovation, and can optimize the supply chain and drive out cost, which positively impacts shareholder value.

The MediChoice brand has grown as a result of our dedication to providing high-quality products to the point-of-need in the most streamlined, cost-effective way possible. And we do all this with superior service levels. Likewise, our commitment to supplier diversity in our sourcing relationships has grown with it. Currently, products sourced from diverse suppliers account for nearly one-fourth of MediChoice sales. By choosing MediChoice, hospitals realize significant savings while advancing their organization's supplier diversity goals.

This Diversity Supplier indicator  can be seen throughout the catalog.

01 HOME OFFICE

9120 Lockwood Blvd.
Mechanicsville, VA 23116
(804) 723-7000 • FAX (804) 723-7100

69 ALBUQUERQUE LDC

8000 Bluewater Rd. NW, Ste. H
Albuquerque, NM 87121
(505) 352-1300 • FAX (505) 352-3800

70 ALLENTOWN MDC

7437 Industrial Blvd.
Allentown, PA 18106-9378
(610) 706-3149 • FAX (610) 391-9480

30 ATLANTA MDC

3375 Cobb International Blvd.
Kennesaw, GA 30152-4399
(770) 426-3700 • FAX (770) 426-3730

57 AUSTIN LDC

2120 Grand Ave. Pkwy. South, Ste. 100
Austin, TX 78728
(512) 252-8450 • FAX (512) 252-8542

89 BALTIMORE MDC

1425 Magellan Rd.
Hanover, MD 21076
(410) 865-7700 • FAX (410) 691-2364

71 BOSTON MDC

20 Freedom Way
Franklin, MA 02038-2584
(508) 520-2000 • FAX (508) 520-2080

20 CHARLOTTE MDC

202 Commerce Blvd.
Kings Mountain, NC 28086
(704) 730-8117 • FAX (704) 730-1514

98 CHICAGO RDC

437 Tower Blvd.
Carol Stream, IL 60188
(630) 588-2800 • FAX (630) 588-2318

75 CINCINNATI MDC

1605 Worldwide Blvd.
Hebron, KY 41048
(859) 282-7910 • FAX (859) 282-7912

58 CLEVELAND SLC

1160 Marquette St.
Cleveland, OH 44114
(216) 416-8344 • FAX (216) 416-8403

41 DALLAS MDC

550 Lakeside Pkwy., Ste. 200
Flower Mound, TX 75028
(972) 538-5500 • FAX (972) 393-7627

53 DENVER LDC

11440 East 56th Ave., Ste. 150
Denver, CO 80239-2228
(303) 576-4200 • FAX (303) 375-1859

48 DES MOINES SLC

4488 112th St.
Urbandale, IA 50322
(515) 276-0206 • FAX (515) 251-8367

44 DETROIT MDC

36663 Van Born Rd., Ste. 390
Romulus, MI 48174-4160
(734) 728-2020 • FAX (734) 728-2222

45 FT. LAUDERDALE LDC

14599 NW 8th St.
Sunrise, FL 33325-6248
(954) 845-2200 • FAX (954) 845-2250

60 GAINESVILLE SLC

4807 NE 63rd Ave.
Gainesville, FL 32609
(352) 264-9939 • FAX (352) 264-2289

49 HARLINGEN LDC

Valley Industrial Center #1
1805 N. Loop 499, Ste. 130
Harlingen, TX 78550
(956) 428-6401 • FAX (956) 428-8857

83 HONOLULU LDC

91-240 Kalaeloa Blvd.
Kapolei, HI 96707
(808) 682-2444 • FAX (808) 682-0555

51 HOUSTON MDC

2700 Brittmoore
Houston, TX 77043-1402
(713) 939-1464 • FAX (713) 939-3325

81 INDIANAPOLIS LDC

6585 East 30th St., Ste. A
Indianapolis, IN 46219-1101
(317) 568-4500 • FAX (317) 549-1801

56 JACKSON LDC

326 U.S. Hwy. 49 South
Richland, MS 39218-9439
(601) 933-6280 • FAX (601) 933-6292

59 JACKSONVILLE LDC

8489 Westside Industrial Dr.
Jacksonville, FL 32219
(904) 596-4600 • FAX (904) 783-8376

85 KANSAS CITY LDC

1800 Southern Rd.
Kansas City, MO 64120-1144
(816) 483-6613 • FAX (816) 231-5902

64 KNOXVILLE LDC

3551 Workman Rd.
Knoxville, TN 37921
(865) 588-0676 • FAX (865) 588-0705

03 LOUISVILLE MDC

6201 Global Distribution Way
Louisville, KY 40228
(502) 671-7550 • FAX (502) 491-3955

67 MEMPHIS LDC

5755 FedEx Lane, Ste 100
Olive Branch, MS 38654
(662) 893-7231 • FAX (662) 893-7819

66 MINNEAPOLIS MDC

2151 County Rd. H2W
Mounds View, MN 55112-4729
(763) 785-9100 • FAX (763) 785-2615

62 NASHVILLE LDC

3308 Briley Park Blvd. South
Nashville, TN 37207
(615) 650-9200 • FAX (615) 650-9311

68 NEW ORLEANS LDC

310 James Dr. East, Ste. 100
Saint Rose, LA 70087
(504) 472-5783 • FAX (504) 472-5750

77 OMAHA LDC

10288 South 134th St.
Omaha, NE 68138
(402) 455-4444 • FAX (402) 896-0884

37 PHILADELPHIA MDC

1220 Forest Pkwy.
West Deptford, NJ 08066
(856) 423-9900 • FAX (856) 423-7387

80 PHOENIX LDC

8313 W. Pierce St., Ste. 100
Tolleson, AZ 85353
(602) 269-7121 • FAX (602) 278-7883

78 PITTSBURGH LDC

134 Industrial Park Rd., Ste. 1000
Greensburg, PA 15601-5539
(724) 837-5700 • FAX (724) 830-7025

82 PORTLAND LDC

27695 SW 95th Ave.
Wilsonville, OR 97070-8215
(503) 685-9800 • FAX (503) 682-0228

84 RALEIGH LDC

4200 Global St.
Raleigh, NC 27610
(919) 231-6367 • FAX (919) 231-6352

87 RICHMOND MDC

11900 North Lakeridge Pkwy.
Ashland, VA 23005
(804) 752-5200 • FAX (804) 798-6728

90 SALT LAKE CITY LDC

2475 South 3200 West, Ste. A
West Valley City, UT 84119-1236
(801) 972-8188 • FAX (801) 973-0734

52 SAN DIEGO LDC

1424 Corporate Center Drive,
Suites 130 & 140
San Diego, CA 92154
(619) 661-4001 • FAX (619) 661-4040

50 SAN FRANCISCO MDC

18520 Stanford Rd.
Tracy, CA 95377
(209) 833-4600 • FAX (209) 833-4639

92 SEATTLE LDC

4301 West Valley Hwy. East, Ste. 300
Sumner, WA 98390
(253) 863-1855 • FAX (253) 863-1865

65 SOUTHERN CALIFORNIA RDC

5125 Ontario Mills Pkwy.
Ontario, CA 91764
(909) 944-2100 • FAX (909) 944-2198

14 ST. LOUIS LDC

15 Gateway Commerce Ctr. Dr. West
Edwardsville, IL 62025
(618) 931-4396 • FAX (618) 931-4661

16 TULSA LDC

8831 East Pine, Ste. 100
Tulsa, OK 74115-5801
(918) 838-8012 • FAX (918) 838-2688

PACKAGING FACILITIES

Asheville (MAI)

25 Heywood Rd.
Arden, NC 28704
(828) 338-7540 • FAX (828) 681-8828

Williamsburg (MAI)

9000 Westmont Dr.
Toano, VA 23168
(757) 566-3510 • FAX (757) 566-0212

OTHER FACILITIES

301 California Sales Office

(West Coast Sales Office)
452-4 Sespe Ave.
Fillmore, CA 93015
(805) 524-0243 • FAX (805) 524-0220

54 O&M Stanford

820 Quarry Rd., RM H0540
Palo Alto, CA 94304
(650) 736-1966 • FAX (650) 736-9536

KEY

LDC = Local Distribution Center
MDC = Market Distribution Center
SLC = Strategic Logistics Center
MAI = Medical Action, Inc
RDC = Regional Distribution Center

This listing is printed quarterly and is subject to change. A more up-to-date version may be available at: www.owens-minor.com

This listing is for general information purposes. If shipping and billing addresses are required, please contact the appropriate distribution center.

NOTES



About MediChoice®

MediChoice is Owens & Minor's brand of medical-surgical products and durable medical equipment that is expanding to meet the changing needs of healthcare. The MediChoice product portfolio reflects Owens & Minor's longstanding commitment to building strong and diverse supply chain relationships with leading manufacturer partners. Products we source from women-, minority-, and veteran-owned businesses make up over one-fourth of the MediChoice product line.

By offering a single brand for thousands of regularly purchased items, MediChoice can increase efficiency and support product standardization initiatives. With MediChoice, the savings can be immediate. Our competitively priced products save hospitals an average of 10 percent on the most regularly purchased supplies. With near-perfect fill rates, MediChoice customers get premium service on common items while controlling their supply spend.

Let MediChoice help you meet your patient-care needs, and solve your budget concerns.

About Owens & Minor

OWENS & MINOR, INC (NYSE: OMI) is a leading healthcare logistics company dedicated to Connecting the World of Medical Products to the Point of Care™ by providing vital supply chain services to healthcare providers and manufacturers of healthcare products. Owens & Minor provides logistics services across the spectrum of medical products—from disposable medical supplies to devices and implants.

With logistics platforms strategically located in the United States and Europe, Owens & Minor serves markets where three-quarters of global healthcare spending occurs. Owens & Minor's customers span the healthcare market from independent hospitals to large, integrated healthcare networks, as well as group purchasing organizations, healthcare products manufacturers, and the federal government. A FORTUNE 500 company, Owens & Minor is headquartered in Richmond, Virginia, and has annualized revenues exceeding \$9 billion. For more information about Owens & Minor, visit the company website at www.owens-minor.com.

Contact Information

For more information on these or other MediChoice products, contact your local Owens & Minor representative or visit our website at medichoice.com.

